# MONOSPOT<sup>®</sup>

**REF** 776060

 **20 tests**

**IVD** For in vitro diagnostic use only

ENGLISH p 2

A qualitative and semiquantitative differential slide test for the Heterophile antibody of infectious mononucleosis

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### SUMMARY AND EXPLANATION

Certain, apparently unique heterophile antibodies appear in the serum of patients during the infectious mononucleosis process. Measurements of these antibodies are performed by indirect procedures. One type of heterophile antibody was originally described by Forssman and bears his name. The Forssman antibody is one of a group of heterophile antibodies which has the ability to react with antigens apparently unrelated to those that stimulated their production. [1] The heterophile antibody associated with infectious mononucleosis is not the same as the Forssman antibody, since the IM antigen is a distinct entity. [2] A differential absorption of the patient's serum or plasma is therefore necessary to distinguish the specific heterophile antibody of infectious mononucleosis from other heterophile antibodies. [3]

The basic principle of the absorption steps in MONOSPOT are comparable to those originally suggested by Davidsohn in his sheep agglutinin test. [1,2] For detection of the specific heterophile antibody associated with IM, MONOSPOT utilises fresh stabilised horse erythrocytes (red blood cells) which have been shown to be more sensitive than sheep erythrocytes or formalinised horse erythrocytes. [4] Infectious diseases such as influenza, rubella, hepatitis and others may cause clinical symptoms that mimic infectious mononucleosis and present problems in diagnosis. [5] The final diagnosis of infectious mononucleosis depends upon clinical, haematological and serologic findings. [6]

Although usually seen in adolescents and young adults, infectious mononucleosis has been reported in children as young as one and one-half years, [7] and also in patients over the age of 60. A few documented examples of recurrence of infectious mononucleosis have been reported, however, they have been extremely rare. [8,9] It has been shown that a non-specific respiratory infection may cause a significant heterophile antibody response in patients who have had infectious mononucleosis, months or years previously. [9]

The specific heterophile antibody of IM is usually demonstrable by the sixth to tenth day after the onset of illness. It is at its highest level during the second to third week, and may persist for at least six weeks. [10] The level of antibody activity is not correlated with the severity of the disease or the degree of lymphocytosis. [11]

### PRINCIPLE OF THE PROCEDURE

The procedure is based on agglutination of horse red blood cells by the heterophile antibody of infectious mononucleosis. Since horse red blood cells contain both Forssman and IM antigens, a differential absorption of the patient's serum is necessary to distinguish the specific heterophile antibody of IM from those of the Forssman type. This is accomplished by absorbing the serum or plasma with both guinea pig kidney and beef erythrocyte stroma. Guinea pig kidney contains only the Forssman antigen while beef erythrocytes contain only the antigen associated with infectious mononucleosis. Therefore, guinea pig kidney will absorb only heterophile antibodies of the Forssman type while beef erythrocytes will absorb only the heterophile antibody of infectious mononucleosis. Agglutination of horse red blood cells by the absorbed patient specimen is indicative of a positive reaction for the IM heterophile antibody (see Results).

### REAGENTS

1. **Reagent I** (Guinea Pig Antigen) is suspension of guinea pig kidney antigen preserved with sodium azide (1%).
2. **Reagent II** (Beef Erythrocyte Stroma) is a suspension of beef erythrocyte stroma antigen preserved with sodium azide (1%).
3. **Indicator Cells** (Horse Erythrocytes) is a suspension of stabilised horse red blood cells preserved with chloramphenicol (1:3,000), neomycin sulfate (1:10,000) and gentamycin (1:20,000).
4. **Positive Control Serum** (Human) is human serum containing the heterophile antibody of infectious mononucleosis, preserved with sodium azide (0.1%).
5. **Negative Control Serum** (Human) is human serum not containing the heterophile antibody of infectious mononucleosis, preserved with sodium azide (0.1%).

For in vitro diagnostic use.

CAUTION: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

### REAGENT PREPARATION

The reagents are ready to use. Prior to use, however, Reagent I, Reagent II and Indicator Cells should be shaken well to provide a homogeneous suspension.

### STORAGE INSTRUCTIONS

When properly stored, the kit is stable through the expiration date printed on the package label. The expiration date of the kit corresponds to the shortest dated reagent. Contents of the entire kit should be discarded on the expiration date. Reagents must be stored at 2 to 8°C when not in use. DO NOT FREEZE REAGENTS. Deterioration of Indicator Cells will be signified by hemolysis. Proper reactivity may be verified by a positive control test (see Procedure).

## **SPECIMEN COLLECTION AND PREPARATION**

No special preparation of the patient is required prior to specimen collection. Serum, or plasma mixed with anticoagulants including EDTA, sodium oxalate, potassium oxalate, sodium citrate, ACD solution or heparin, may be used. Inactivation of the serum is not necessary, however, inactivated serum may be used. Serum or plasma samples should be clear and particle-free. The heterophile antibody of infectious mononucleosis is relatively stable and uncontaminated specimens or serum or plasma may be stored at 2 to 8°C for several days following collection before testing. If prolonged storage is desired, the patient specimen should be frozen.

A sample for testing with MONOSPOT may be obtained from finger stick blood using the procedure outlined below:

1. Using standard heparinized or non-heparinized capillary pipettes (75 mm length, 1.1 mm to 1.2 mm inside diameter, 85µL volume) obtain a minimum of four capillary pipettes full of blood from a finger puncture.
2. Prepare capillary pipettes and centrifuge according to standard procedure for hematocrit testing.
3. Break capillary pipettes at interface between serum or plasma and cells and use serum or plasma from two capillary pipettes as patient sample for each side of the slide according to the directions for use. Sufficient patient specimen for the test (0.05 mL for each side of the slide) is recovered by this procedure if a hematocrit of less than 50% is obtained. If the hematocrit should be greater than 50%, additional capillary pipettes may be necessary to yield the required specimen volume.

## **PROCEDURE**

### **Materials provided**

Reagent I (Guinea Pig Kidney) 1 mL

Reagent II (Beef Erythrocyte Stroma) 1 mL

Indicator Cells (Horse Erythrocytes) 20 tests

Positive Control Serum (Human) 0.5 mL

Negative Control Serum (Human) 0.5 mL

Glass Slide - Microcapillary Pipettes, 20 mcl capacity, for delivery of Indicator Cells - Rubber Bulbs for use with Microcapillary Pipettes - Plastic Pipettes for delivery of serum or plasma samples - Wooden Applicators for mixing.

All materials necessary to perform the tests are included in the package - no additional material is needed.

### **Directions for use – Qualitative method**

1. Place the slide on a flat surface under a direct light source.
2. Invert the vial of Indicator Cells several times to uniformly suspend the cells. Using one of the disposable microcapillary pipettes provided, deliver 10 mcl of Indicator Cells to one corner of both squares on the slide.  
To use the disposable microcapillary pipettes:
  - a. Insert the end of the pipette marked with a heavy black line one-quarter inch into the neck of the rubber bulb.
  - b. Hold the rubber bulb between the thumb and third finger.
  - c. Tilt the vial of Indicator Cells, insert the pipette, and allow the pipette to fill by capillary action to the top mark (20 mcl). Do not draw cells into the bulb.
  - d. To deliver 10 mcl of Indicator Cells to the slide, place the index finger over the hole in the top of the bulb and squeeze gently until the level of cells in the pipette reaches the first mark. Touch pipette tip to a corner of square "I" to release the cells.
  - e. Deliver the remaining 10 mcl of Indicator cells from the pipette to a corner of square "II".
3. Place one drop of thoroughly mixed Reagent I in the center of square "I".
4. Place one drop of thoroughly mixed Reagent II in the center of square "II".
5. Using a disposable plastic pipette, add one drop of the serum or plasma being tested to the center of each square on the slide.
6. Avoiding the Indicator Cells, mix Reagent I and Reagent II with the serum or plasma in each square at least ten times using a clean wooden applicator for each square. Then, using no more than ten stirring motions, blend in the Indicator Cells over the entire surface of each square. Start a timer upon completion of the final mixing.
7. DO NOT PICK UP OR MOVE THE SLIDE DURING THE REACTION PERIOD.
8. Observe for agglutination for no longer than one minute after final mixing.

Both Positive and Negative Control Serum are provided in each package. These control sera should be used to check the reagents upon arrival in the laboratory and periodically during the dating period.

### **RESULTS - Qualitative method**

Proper performance of the test will indicate the presence or absence of the heterophile antibody of infectious mononucleosis by the presence or absence of agglutination, as indicated below.

#### **POSITIVE TEST**

If the agglutination pattern is stronger on the left side of the slide ( Square I), the test is positive.

#### **NEGATIVE TEST**

If the agglutination pattern is stronger on the right side of the slide ( Square II), the test is negative.

If no agglutination appears on the either square of the slide ( I and II), or if agglutination is equal on both squares of the slide ( I and II), the test is negative.

Agglutination which appears after the one-minute reading time, or which appears when the slide is moved or rocked, is not to be interpreted as a positive result for the heterophile antibody of infectious mononucleosis.

A positive test result is indicative of the heterophile antibody specific for infectious mononucleosis. A titration procedure, as described below, may be performed to provide a semi-quantitative indication of the level of infectious mononucleosis heterophile antibody.

**TITRATION PROCEDURE - Directions for use – Semi-quantitative method**

1. Prepare the following dilutions of serum known to be positive for the heterophile antibody of infectious mononucleosis.

TUBE	DILUTION	COMPOSITION
1	1:2	0.5 mL of patients serum + 0.5 mL physiological saline; mix
2	1:4	0.5 mL of tube 1 + 0.5 mL physiological saline; mix
3	1:8	0.5 mL of tube 2 + 0.5 mL physiological saline; mix
4	1:16	0.5 mL of tube 3 + 0.5 mL physiological saline; mix
5	1:32	0.5 mL of tube 4 + 0.5 mL physiological saline; mix

2. Place a titration slide on a flat surface under a direct light source.
3. Invert the vial of Indicator Cells several times to uniformly suspend the cells. Using one of the disposable microcapillary pipettes provided, deliver 10 µl of Indicator Cells to one corner of each square.
4. Place one drop of thoroughly mixed Reagent I in the center of each of the five squares.
5. Using a disposable plastic pipette, place one drop of each of the dilutions under test in the center of its corresponding square ( Tube 1 to square 1 , etc.).
6. Avoiding the Indicator Cells, mix the reagent and diluted specimen in each square at least ten times with a clean wooden applicator. Then, using no more than ten stirring motions, blend in the Indicator Cells over the entire surface of each square. Start a timer upon completion of the final mixing.
7. DO NOT PICK UP OR MOVE SLIDE DURING THE REACTION PERIOD.
8. Observe for agglutination for no longer than one minute after final mixing.

**RESULTS Semi-quantitative method**

1. The highest dilution in which visible agglutination occurs is the end point. (If agglutination is present in all dilution, extend doubling dilutions from tube 5 and repeat the procedure as outlined above.)
2. Record the titration value obtained. The titer with MONOSPOT cannot be compared with titration values obtained with other slide or tube test procedures since there are significant variations in sensitivity. However, the relative magnitude of titers performed on specimens with different techniques is proportional to the concentration of IM heterophile antibody present.
3. Although the titration value is not indicative of the severity of the disease, sequential testing of the patient with MONOSPOT may provide information of value to the clinician.

**LIMITATIONS OF THE PROCEDURE**

To obtain accurate results, only clear, particle-free serum or plasma samples should be used. False positive results with MONOSPOT have been reported.[12-14] Since the clinical symptoms of IM are similar to many other diseases, it is often difficult to disprove the theoretical possibility of a concomitant infection.[15] The simultaneous occurrence of infectious mononucleosis and hepatitis has been reported.[16] A result which may be interpreted as false positive may possibly be due to residual infectious mononucleosis antibody present after clinical symptoms have subsided.[5,10] Seronegative infectious mononucleosis has also been reported.[5] Because of a delayed heterophile antibody response, it is possible that clinical and haematological symptoms of IM will appear before serologic confirmation is possible. [17]

**PERFORMANCE CHARACTERISTICS**

When the test with MONOSPOT is performed in accordance with the directions for use, accuracy greater than 99% is obtainable.[18-22] The result must be combined with clinical and haematological features for diagnosis of infectious mononucleosis. MONOSPOT has been shown to be highly specific for the heterophile antibody of infectious mononucleosis through the use of differential absorption.[2,23] MONOSPOT utilises fresh, stabilised horse red blood cells, which have been shown to be more sensitive than sheep red blood cells or formalinised horse red blood cells.[6,23-25] These performance characteristics were gathered from a large number of papers which have appeared in the scientific literature, in addition to testing performed by the manufacturer.

Monospot	Infectious Mononucleosis [26]		
	Pos	Neg	Total
Pos	94	2	96
Neg	6	118	104
Total	100	120	220

Sensitivity = 94/100 x 100% = 94% - Specificity = 118/120 x 100% = 98%  
 VPP = 97.9% - VPN = 95.1% - Accuracy = 96.3%

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## REFERENCES

1. Davidsohn I., Walker P.H. The nature of the heterophilic antibodies in infectious mononucleosis. I. Am. J. Clin. Pathol. 1935; 5: 455.
2. Davidsohn I. Serologic diagnosis of infectious mononucleosis. JAMA 1937; 108: 289.
3. Lee C.L., Davidsohn I., Panczyszyn O. Horse agglutinins in infectious mononucleosis II - the spot test. Am. J. Clin. Pathol. 1968; 49: 12.
4. Lee C.L., Davidsohn I., Slaby R. Horse agglutinins in infectious mononucleosis. Am. J. Clin. Pathol. 1968; 49: 3.
5. Carter R.L., Penman H.G. Infectious mononucleosis. Blackwell Scientific Publications 1968.
6. Lee C.L., Davidsohn I. Spot test for infectious mononucleosis. Exhibit, Joint Annual Meeting ASCP and CAP 1967.
7. Starling K.A., Fernbach D.J. Infectious mononucleosis in the preschool child. Letters to the Journal, JAMA 1968; 203: 810.
8. Bender C.E. Recurrent mononucleosis. JAMA 1962; 182: 954.
9. Hoagland R.J. Resurgent heterophil-antibody reaction after mononucleosis. N. Engl. J. Med. 1963; 269: 1307.
10. Lee C.L., Davidsohn I. Serologic tests for infectious mononucleosis. Comm on Continuing Educ. Council on Immunohematology ASCP 1972.
11. Wintrobe M.M. Infectious mononucleosis. In Harrison T.R.: Principles of Internal Medicine. The Blackiston Company 1950.
12. Sadoff L., Goldsmith O. False positive infectious mononucleosis spot test in pancreatic carcinoma. JAMA 1971; 218: 1297.
13. Philips G.M. False-positive MONOSPOT test result in rubella. Letters to the Journal, JAMA 1972; 222: 585.
14. Horwitz C.A., Polesky H., Stillman T., Ward P.C.J., Henle G., Henle W. Persistent haemagglutination for infectious mononucleosis in rheumatoid arthritis. Br. Med. J. 1973; 1: 591.
15. Tri T.B., Herbst K. False-positive MONOSPOT or rubella and mononucleosis. Letters to the Journal, JAMA 1973; 255: 63.
16. Madhavan T. hepatitis and infectious mononucleosis. Letters to the Journal, JAMA 1973; 225: 314.
17. Bæhner R.L., Shuler S.E. Infectious mononucleosis in childhood. Clin. Pediatr. (Phila) 1967; 6: 393.
18. Davidsohn I., Henry J.B. Todd-Sanford clinical diagnosis by laboratory methods. 14th edition. W.B. Saunders Company 1969.
19. Wahren B. Diagnosis of infectious mononucleosis by the MONOSPOT test. Am. J. Clin. Pathol. 1969; 52: 303.
20. Galloway E. Comparison of three slide tests for infectious mononucleosis with Davidsohn's presumptive and differential heterophile test. Can. J. Med. Tech. Dec. 1969; 31: 197.
21. Seitanidis B. A comparison of the MONOSPOT with the Paul-Bunnell test in infectious mononucleosis and other diseases. J. Clin. Pathol. 1969; 22: 321.
22. Evans D.M.D., Sanerkin N.G., Lewis J. Facilitating the laboratory diagnosis of infectious mononucleosis. Am. J. Clin. Pathol. 1969; 52: 702.
23. Lee C.L. Spot test for infectious mononucleosis. Bull. of Pathol. 1968; 9: 78.
24. Barrett A.M. The serological diagnosis of glandular fever (infectious mononucleosis); a new technique. J. Hyg. (Camb) 1941; 41: 330.
25. Wilkinson P.C., Carmichael D.S. Immunochemical characterization and serologic behavior of antibodies against red cells in infectious mononucleosis. J. Lab. Clin. Med. 1964; 64: 529.
26. Tilton RC et al. Comparative evaluation of three commercial tests for detection of heterophile antibody in patients with infectious mononucleosis. J Clin Microbiol 1988;26:275-278.

 **Meridian Bioscience Europe Srl**  
Via dell' Industria, 7  
I-20020 Villa Cortese, Milano - Italy  
Tel: +39 (0331) 43 36 36  
Fax: +39 (0331) 43 36 16  
E-mail: info@mdeur.com

**Meridian Bioscience Europe sa/nv**  
Rue de l'Industrie, 7  
B-1400 Nivelles - Belgium  
Tel: +32 (67) 89 59 59  
Fax: +32 (67) 89 59 58  
E-mail: info@mdeur.be

**Meridian Bioscience Europe France**  
Le Quadra, Promenade des Anglais, 455  
F-06299 Nice Cedex 3 - France  
Tel: +33 (4) 93 18 72 10  
Fax: +33 (4) 93 18 72 11  
E-mail: info@meridianbioscience.fr

**Meridian Bioscience Europe bv**  
Halderheiweg, 6  
NL-5282 SN Boxtel - The Netherlands  
Tel: +31 (411) 62 11 66  
Fax: +31 (411) 62 48 41  
E-mail: meridian@wxs.nl

**Meridian Bioscience, Inc.**  
Corporate Office  
3471 River Hills Drive  
Cincinnati, Ohio 45244, USA  
Tel: +1 (513) 271 37 00  
Fax: +1 (513) 271 37 62  
E-mail: mbi@meridianbioscience.com

Website [www.mdeur.com](http://www.mdeur.com)

Expiry date / Date de péremption/ Scadenza/ Fecha de caducidad/ Haltbarkeitsdatum

Lot number / Numéro de lot / Numero di lotto/ Número de lote / Chargenbezeichnung

For in vitro diagnostic use / Usage in vitro/ Per uso diagnostico in vitro Uso diagnóstico in vitro / In vitro Diagnosticum

This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro/ Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro/ Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro/ Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG

Catalogue number / Référence article/Codice dell'articolo/ Número de catálogo / Artikelnummer

Please read pack insert / Lire attentivement le mode d'emploi / Leggere il foglietto informativo/ Leer instrucciones de uso / Bitte Packungsbeilage beachten

Manufactured by / Fabriqué par/ Prodotto da / Fabricado por / Hergestellt von

Contains sufficient for <n> tests / Contenu suffisant pour "n"tests/ Contenido suficiente per "n" prove/ Contenido suficiente para <n> ensayos / Enthält ausreichend für <n> Tests

Conservare a/ Store at / Conserver entre/ Conservar a/ Lagerung bei

Rev. 2/04