# Fungal Immunodiffusion Reagents

**REF** 100201, 100401, 100405,  
100601, 100801, 101009,  
101012, 103001, 103101

**IVD**

Rx Only

**INTENDED USE**  
The Meridian Bioscience Immunodiffusion Reagents are standardized, purified preparations for the in vitro determination of precipitating antibodies to two systemic fungal pathogens: *Histoplasma capsulatum* (*H. capsulatum*) and *Coccidioides immitis* (*C. immitis*). Reagents are available for detection of either *Coccidioides immitis* "TP" (early) or "F" (late) antibodies. These reagents are optimized for use in the Ouchterlony double diffusion method.

**SUMMARY AND EXPLANATION OF THE TEST**  
*H. capsulatum* and *C. immitis* are causative agents of deep-seated mycoses. These fungi present a diagnostic challenge to the microbiology laboratory and the physician.

Radiographically, lesions produced by the systemic fungi may be difficult to distinguish from each other, tuberculous lesions, neoplasms, (benign and malignant) cancerous tissue and each other. Symptoms are often unremarkable and can mimic various pneumonias, sarcoidosis, cancer and other maladies. Culturally and histologically the organisms can be difficult to demonstrate, even after repeated attempts.<sup>3-7</sup>

Frequently serology offers the only evidence available to guide treatment, suggest prognosis or lead to the selection of more definitive diagnostic techniques, such as intensive culture or biopsy. In addition, quantitative serology, such as complement fixation testing, can provide important information on the effects of chemotherapy.<sup>18</sup>


Fungal immunodiffusion tests are simple and reliable tools for the evaluation of suspected mycoses. Anticomplementary sera may be tested with these procedures. The tests also provide specificity data on reactions obtained by the complement fixation method. No expensive equipment is required and the procedures are simple enough to be performed by any laboratory, thus providing excellent screening methods.

**BIOLOGICAL PRINCIPLES**  
The Meridian Bioscience Fungal Immunodiffusion Reagents are based on the principle of double diffusion as described by Oudin and Ouchterlony. An antibody and its homologous soluble antigen are placed in separate wells cut in an agarose diffusion medium and allowed to diffuse outward. Between the two wells, a concentration gradient of each of the reaction components is established ranging from antigen excess closest to the antigen well, to antibody excess closest to the antibody well. A visible line of precipitate forms at the point of equivalence.<sup>10-12</sup>

The well patterns of the Meridian immunodiffusion plate are arranged to provide each patient test well with a reference known system so that identity reactions are readily apparent. (See figure 1)


Antigens or antibodies may be tested for "identity" by placing a test well of the substance in question adjacent to the wells of a known system. If the antigen antibody complexes are identical, the precipitin lines form an unbroken line of identity with the known system. Partial and nonidentity reactions are also possible. (See figure 2)

Figure 1. Arrangement of sera and antigens in immunodiffusion well pattern.

Control positive serum, wells 1 & 4. Patient sera, wells 2, 3, 5 & 6. Antigen, well 7.

Figure 2. Identification of Immunodiffusion Bands

A partial identity reaction occurs when certain components of the antigens (or antibodies) are identical and others are not. The "spur" represents the components which are unrelated. A nonidentity reaction will occur when the antigen-antibody complexes are different. The resulting "X" or crossed reaction indicates that two unrelated complexes are present.

**REAGENTS/MATERIALS PROVIDED**

- Antigens:**
  - Histoplasmin: Catalog #100201, A purified mycelia cultural filtrate from the growth of *Histoplasma capsulatum* containing the "H" and "M" antigens in 1.0 mL size. The preparation contains 100 units/mL\* of each antigen. Lipids, media components, and most other antigenic components have been removed. (Note: The "H" band appears closest to the antibody well while the "M" band is closest to the antigen well.)
  - Coccidioidin "F": Catalog #100401 (1 mL) / Catalog #100405 (5 mL), Coccidioides "F" antigens are available in 1.0 mL and 5.0 mL sizes. The available antigens consist of cultural filtrates containing the "F" and "TP" antigens of *Coccidioides immitis* at a concentration of 100 units/mL.\*
  - Coccidioidin "TP": Catalog #103001, Coccidioides "TP" antigen consisting of purified cultural filtrates containing the "TP" antigen of *Coccidioides immitis* in a 1.0 mL size. The preparation contains 100 units/mL of antigen\* These antigens contain thimerosal at a final concentration of 1:10,000 as a preservative. Care should be exercised, however, to avoid microbial contamination.

\*These units are internal Meridian standards designed only to ensure lot to lot reproducibility of antigen concentration.

- Control Sera:** The following lyophilized control sera are provided in 1.0 mL size. Control sera contain either goat, rabbit or human serum with antibodies directed against their respective antigens. These control sera are preserved with 0.1% sodium azide. Care should be exercised to avoid microbial contamination.
  - Histoplasmosis Positive Control Catalog #100601: Contains specific antibodies directed against the "H" and "M" antigens of *H. capsulatum*. (Note: the "H" band appears near the antibody well, while the "M" band is closest to the antigen well.)
  - Coccidioidomycosis "F" Positive Control Catalog #100801: Contains specific antibodies directed against the "F" antigen of *C. immitis*.
  - Coccidioidomycosis "TP" Positive Control Catalog #103101: Contains specific antibodies directed against the "TP" antigen of *C. immitis*.

- Immunodiffusion Plates:** The immunodiffusion plates consist of 0.9% agarose in glycine/phosphate buffer at pH 7.4 + 0.2.
  - 1 Series Plates (12 plates/pack) Catalog #101012: Agarose gel plate with wells to test four patients per plate.
  - 4 Series Plates (9 plates/pack) Catalog #101009: Agarose gel plate with wells to test 16 patients per plate.

**Note:** 1 Series plates must be used for *C. immitis* "TP" testing. 1 Series or 4 Series Plates may be used for *C. immitis* "F" and *H. capsulatum* testing.

**Note:** Recent advances in plate formulation have eliminated the need for Immunodiffusion Band Intensifying Fluid.

**Note:** Due to the standardization necessary to the production of high quality Meridian Bioscience fungal serology reagents, performance of tests with antigens, control sera and plates with materials other than those produced by Meridian Bioscience, Inc. cannot be guaranteed. The user assumes full responsibility for any modification to the procedures published herein.

**MATERIALS AND EQUIPMENT REQUIRED BUT NOT SUPPLIED**

- Reagent quality water.
- Moist chamber; any convenient container such as a petri dish, plastic box or glass jar with a tight fitting cover containing moist filter paper or paper toweling is satisfactory provided the ID plates are stationary, level and remain hydrated during incubation.
- Reading light. Dark field plate readers are commercially available; however, a satisfactory system can be devised using a high intensity lamp in which the front of the bulb shield reflector is covered with black construction paper with a small hole (1-2 cm in diameter) cut for illumination. Alternatively, the plate may be held at a 45 degree angle to almost any bright indirect light source for adequate visualization.
- Pipette, capillary tube, or equivalent

- PRECAUTIONS**
- All reagents are for in vitro diagnostic use only.
  - When handling blood specimens, adequate measures should be taken to prevent the dissemination of etiological agents potentially present in the specimen.
  - The control sera are preserved with sodium azide. Accumulation of this chemical in metal plumbing fixtures may represent an explosion hazard. It is therefore recommended that excess control serum simply be discarded in an appropriate waste receptacle or flushed down a drain with copious amounts of water.

**HAZARD and PRECAUTIONARY STATEMENTS**

There are no known hazards associated with this product: 100201, 100401, 100405, 103001.

REF 100601 & 100801 & 103101	<p><b>Signal Word:</b> Danger</p> <p><b>Hazardous Statements:</b> H300 – Fatal if swallowed H400 – Very toxic to aquatic life H410 – Very toxic to aquatic life with long lasting effects EUH032 – Contact with acids liberates very toxic gas <b>Precautionary Statement – EU (§28, 1272/2008)</b> P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor/ physician</p>
 Control Serum	

REF 101009 & 101012	<p><b>Signal Word:</b> Warning</p> <p><b>Hazardous Statements:</b> H302 - Harmful if swallowed</p>
 Immunodiffusion Plates	

**SHELF LIFE AND STORAGE**

Fungal antigens are best stored at refrigerated temperature (2-8 C). Repeated freezing and thawing is detrimental to these antigens. Extreme heat (more than 60 C) is also to be avoided. When stored at 2-8 C the antigens are stable until the expiration date indicated on each vial.

Positive control sera are stable in the lyophilized state until the stated expiration as indicated on each vial when stored at 2-8 C. Once reconstituted, the sera can be stored at 2-8 C for one month. If the sera will not be used within one month, it should be aliquoted and frozen at -20 C. At -20 C a minimum nine month life may be expected. Repeated freezing and thawing should be avoided. When the positive control sera are in use, the period at room temperature should be kept as short as possible.

Immunodiffusion plates should be stored unopened at 2-8 C until the stated expiration date on the package. If plate packages are opened and the plates must be stored before use, they must be kept in a tightly sealed, humid container to avoid drying. A resealable plastic container with a small moist sponge or towel is ideal.

**SPECIMEN COLLECTION AND PREPARATION**

- For optimal results, serum is aseptically obtained from the patient.
- If a delay is encountered in specimen processing, refrigeration for up to 72 hours is permissible. Specimens may be stored up to six months at -20 C with no loss of activity, provided they are not repeatedly thawed and refrozen. Specimens in transit between laboratories should be maintained at 2-8 C for optimal results. Specimens may be preserved with 1:10,000 thimerosal or 0.1% sodium azide if necessary.
- Greater sensitivity can be achieved when performing the Cocci "TP" antibody test by concentrating the patient's serum approximately 8:1.<sup>13-16</sup> This can be accomplished by pervaporation or use of a 125000 molecular weight cutoff membrane filter (e.g. Amicon™ B-125 or equivalent). While various strategies are available for pervaporation, the serum may be simply placed in a short length of dialysis tubing and air dried or gently warmed with a hand dryer until the correct volume reduction is obtained.
- In a limited clinical study, a 15% increase in sensitivity was obtained when specimens were concentrated by pervaporation and allowed to prediffuse for two hours prior to addition of antigen. (Pappagianis, et al., personal communication)

**TEST PROCEDURE**



**A. Reconstitution of positive control serum.**

- Using a 1.0 mL pipet or tuberculin syringe add 0.95 mL of reagent quality water to the control serum vials.
- Gently "flick" the bottom of the vials to mix.
- Allow the vials to stand at room temperature for about 1 hour, occasionally inverting to mix. (Do not shake)
- If, on its first use, the positive control serum produces a negative result, it may be that sufficient reconstitution time was not allowed. In this case repeat the test ensuring that sufficient reconstitution time is allowed.
- Aliquot and freeze any control serum which will not be consumed within one month. (See **SHELF LIFE AND STORAGE**)

**B. Test procedure**

- Please refer to figures 1, 3 and 4 for proper arrangement of the test well pattern.
- Fill the control serum wells (1 & 4) of an appropriate ID plate for each fungal analyte to be tested.

**Figure 3.**

<p><b>FUNGAL IMMUNODIFFUSION ANALYSIS</b> (SAMPLE FORM)</p> <p>Plate No. <u>5</u>      Date <u>8-23-17</u></p> <p>Coccidioides "TP"</p> 	<p><b>FUNGAL IMMUNODIFFUSION ANALYSIS</b> (SAMPLE FORM)</p> <p>Plate No. <u>5</u>      Date <u>8-23-17</u></p> <p>Histoplasma</p> 
---	---

- Record the name, date, and/or lab number of the first patient on line 2 of the left hand column (figure 4) of the recording form and fill well 2 with the patient's serum. Using a fresh capillary tube/pipette tip for each patient, repeat this procedure for additional sera in wells 3, 5 and 6.
- Fill the center well with appropriate fungal antigen.
- Number and date the ID plate and the analysis form.
- Place the ID plate in a moist chamber and incubate at room temperature for 24 hours. Note: bands may be evident in as little as three to six hours for some sera. It will still be necessary to confirm the results of the other sera later.
- After 24 hours read and record the ID bands on the analysis form (figure 4). (See Reading the Test) An interim report should be issued at this point if no identity reactions are observed. Positive results should be reported immediately.
- An additional 48 hour incubation period is recommended to confirm a negative result. A final report is made at the conclusion of this period.

**Figure 4.**

Well Number	Histoplasmin Reading	Coccidioidin Reading
1. Control Serum	Histo Control	Cocci Control
2. CARTER, R. 3287	H and M BANDS	NEGATIVE
3. LESTER, C. 2519	NEGATIVE	NEGATIVE
4. Control Serum	Histo Control	Cocci Control
5. BARNES, T. 1362	NEGATIVE	NEGATIVE
6. HASKAMP, D. 1154	NEGATIVE	1 - BAND
7. Antigen	Histoplasmin	Coccidioidin
(Name, Date, Lab No., etc.)	After 24 hour incubation, draw bands observed on the well diagram above and record the number of bands on the blanks provided. NOTE: The control bands have been drawn on the well diagram above. If these bands are not observed on the plates after 24 hour incubation, the test is not valid and must be repeated.	

#### READING THE TEST

Any bright indirect light source may be used to aid in seeing the precipitin bands. A dark background is preferred. For instance, the plate may be held next to the shade of a high intensity lamp whose light beam is directed straight downward, over a black countertop. A bright direct light source may also be used if the plate is held in the light beam at a 45 degree angle against a dark background. Light boxes, providing indirect light against a dark background and prepared specifically for reading immunodiffusion plates, are also commercially available.

Particular attention should be paid to the orientation of bands produced by the patient serum in relation to the control bands. A smooth junction of the bands indicates an identity reaction. When patient antibody levels are low, only a slight bend in the control band toward a position in front of the patient well may be apparent. Should no reaction be observed between the positive control and the antigen, the test must be repeated.

#### INTERPRETATION OF RESULTS

A band of identity with the positive control indicates patient antibody against the fungal antigen in question. Partial identity reactions are also regarded as positive for antibody. Non-identity reactions are regarded as a negative test.<sup>11</sup>

#### QUALITY CONTROL

**This test should be performed per applicable local, state, or federal regulations or accrediting agencies.**

User quality control is performed on each test run by the inclusion of a positive control serum for each antigen. Positive control bands should be present at 24 hours. If these bands are not observed on the plates after the 24 hour incubation, the test is not valid and must be repeated. If a negative control is desired, a known negative patient may be used as negative control. The results obtained with the controls should be recorded in an appropriate log book to maintain high quality testing and comply with regulatory requirements.

**If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Service Department at 1-800-343-3858 (US) or your local distributor.**

#### EXPECTED VALUES

In general, an identity reaction against a given antigen is indicative of active or recent past infection. Partial identity reactions are also important indicators of probable disease, especially when no identity reaction is seen among other agents. Non-identity reactions may be apparent when the disease state is caused by a mycotic agent other than the one tested. A non-identity reaction, however, is a negative test for the organism in question.

If activity against any of the fungal agents is observed, a vigorous attempt should be made to demonstrate the organism culturally for confirmation.

#### Histoplasmosis

The clinically significant antigens for *H. capsulatum* are designated "H" and "M" antigens. (Note: The "H" band appears closest to the antibody well while the "M" band is closest to the antigen well.) The "M" band has been found in about 63% of patients with active histoplasmosis, while the "H" band is apparent in only 27%. The "H" band is rarely observed in the absence of antibodies against the "M" antigen. The "M" band may appear in patients who have recently recovered from histoplasmosis as well as in the serum of previous histoplasmosis patients who have been recently skin tested.

Since the test may be negative in as many as 10% of culturally demonstrated cases, the absence of either an "H" or an "M" band does not rule out histoplasmosis.<sup>1, 21</sup>

#### Coccidiomycosis

Antibodies directed against the "F" antigen indicate active or recent past (up to 1 year) infection. The ID test is usually positive within four weeks after infection and remains positive throughout the clinically active disease.

Antibodies directed against the "TP" antigen indicate recent, active infection. These are primarily IgM antibodies which appear within several days of infection. This ID test is usually positive within 1 week after infection and reverts to negative within a few months.

Latex and complement fixation testing may provide important additional information regarding patient status.<sup>4-6</sup>

Crossreactions may be seen in patients harboring other systemic fungi (especially *H. capsulatum*) so care must be exercised when reading for identity reactions<sup>9, 18, 21</sup> (See **LIMITATIONS OF THE PROCEDURE**). Partial identity reactions are also important indicators of probable disease. Nonidentity reactions may be apparent when the disease state is caused by a mycotic agent other than the one tested. A nonidentity reaction, however, is a negative test for the organism in question.

If activity against *Coccidioides immitis* is observed, a vigorous attempt should be made to demonstrate the organism culturally for confirmation.

#### LIMITATIONS OF THE PROCEDURE

The high rate of negative serologic tests observed among certain culturally demonstrable cases limits the predictive value of a negative test.<sup>20</sup>

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**REFERENCES**

- Bauman, DS. and CD. Smith. Comparison of Immunodiffusion and complement fixation tests in the diagnosis of histoplasmosis. J. Clin. Microbiol. 1975;2:77-80.
- Crowle, AJ. *Immunodiffusion*, Second Edition. Academic Press. New York. 1973; Chapters 1, 2 and 4.
- Heiner, DC. Diagnosis of histoplasmosis using precipitin reactions in agar gel. Pediatrics 1958;22:616-627.
- Huppert, M. and JW. Bailey. The use of immunodiffusion tests in coccidioidomycosis 1. The accuracy and reproducibility of the immunodiffusion test which correlates with complement fixation. Amer. J. Clin. Path. 1965;44:364-368.
- Huppert, M. and JW. Bailey. The use of immunodiffusion tests in coccidioidomycosis 11. An immunodiffusion test as a substitute for the tube precipitin test. Amer. J. Clin. Path. 1965;44:369-373.
- Kaufman, L. Value of immunodiffusion tests in the diagnosis of systemic mycotic diseases. Ann. Clin. Lab. Sci. 1973;3:141-146.
- Kaufman, L. Practical applications of immunologic procedures for the diagnosis of opportunistic fungus infections. In: *Proceedings of the Second International Conference on Opportunistic Fungal Infections*. CC. Thomas Publishing, Springfield. 1975;pp 47-57.
- Kaufman, L. and MJ. Clark. Value of the concomitant use of complement fixation and immunodiffusion tests in the diagnosis of coccidioidomycosis. Appl. Microbiol. 1974;28:641-43.
- Ouchterlony, O. Antigen-antibody reactions in gels and the practical application of this phenomenon in the laboratory diagnosis of diphtheria. Med. Diss. Stockholm, 1949.
- Ouchterlony, O. *Handbook of Immunodiffusion and Immunolectrophoresis*. Ann Arbor Publishers Inc., Ann Arbor, 1968.
- Oudin, J. L'analyse, immunochimique qualitative. Methode par diffusion des antigenes au sein de l'immunoserum precipitant gelose. Premiere parse. Ann Inst. Pasteur 1948;75:30-52.
- Palmer, DF. et al. *Serodiagnosis of Mycotic Diseases*. CC. Thomas Publishing, Springfield. 1977;pp. 7-18.
- Pappagianis, D. Serology and Serodiagnosis of Coccidioidomycosis. In: *Coccidioidomycosis*. David A. Stearns, M.D., Editor, Plenum Publishing.
- Pappagianis, D. Immunoserologic Procedures in the Diagnosis and Management of Pulmonary Mycoses. In: *Handbook of Fungus Diseases*. H. Einstein, Editor. Amer. College of Chest Physicians. 1981.
- Pappagianis, D. Antibody in Cerebrospinal Fluid in Non-meningitic Coccidioidomycosis. Sabouraudia 1972;10:173-179.
- Pappagianis, D., NJ. Lindsey, CE. Smith, and M. Saito. Antibodies in Human Coccidioidomycosis: Immunolectrophoretic Properties. Proceedings of the Society for Experimental Biology and Medicine 1985;118:118-122.
- Pappagianis, D. False-positive Reactions of Cerebrospinal Fluid and Diluted Sera with the Coccidioidin Latex-agglutination Test. Am. J. Clin. Path. 1976; 66:916-921.
- Smith, CE. Coccidioidomycosis. Med. Clin. North Amer. 1943; 27:790-807.
- Smith, CE., MT. Salito, RR. Beard, KM. Kepp, RW. Clark, and BV. Eddie. Serological tests in the diagnosis of coccidioidomycosis. Am. J. Hyg. 1950;52:1-2.
- Tosh, FE. Problems in the serological diagnosis of opportunistic fungus diseases. In: *Proceedings of the Second International Conference on Opportunistic Fungal Infections*. C. C. Thomas Publishing, Springfield, 1975;pp. 36-43.
- Kaufman L. Value of immunodiffusion tests in the diagnosis of systemic mycotic diseases. Ann Clin Lab Sci 1973; 3:14146.

SN10000


REV. 09/23

 <p>Manufactured By</p>	<p><b>Meridian Bioscience, Inc.</b>  3471 River Hills Drive  Cincinnati, OHIO - 45244 USA  <a href="http://www.meridianbioscience.com">www.meridianbioscience.com</a></p> <p><u>Contacts:</u>  Main Telephone (+1) 513.271.3700  Customer Service/Orders 800.543.1980  Technical Support Center 800.343.3858  Information Fax: 513.272.5432  Ordering Fax: 513.271.0124  E-mail: <a href="mailto:info@meridianbioscience.com">info@meridianbioscience.com</a></p>
	<p><b>Meridian Bioscience Europe, SRL</b>  Via Dell'Industria 7, 20035 Villa Cortese  (Milano) ITALY  <a href="http://www.meridianbioscience.com">www.meridianbioscience.com</a></p> <p><u>Contacts:</u>  Main Telephone (+39) 0331.433636  E-mail: <a href="mailto:info@meridianbioscience.eu">info@meridianbioscience.eu</a>  Technical Support: <a href="mailto:MBE-TechService@meridianbioscience.eu">MBE-TechService@meridianbioscience.eu</a>  Customer Service/Orders:</p> <ul style="list-style-type: none"> <li>• For Italian Customers:  <a href="mailto:ordini@meridianbioscience.com">ordini@meridianbioscience.com</a></li> <li>• For Distributors / International Customers:  <a href="mailto:Export.CustomerService@meridianbioscience.eu">Export.CustomerService@meridianbioscience.eu</a></li> </ul>

**INTERNATIONAL SYMBOL USAGE**

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.