Fungal Immunodiffusion System

REF 100201, 100601

IVD In vitro diagnostic medical device

Rx Only

INTENDED USE

The Meridian Bioscience Fungal Immunodiffusion Reagents are standardized, purified preparations for the *in vitro* determination of precipitating antibodies to two systemic fungal pathogens: *Histoplasma capsulatum* (*H. cap.*) and *Coccidioides immitis* (*C. imm.*). These reagents are optimized for use in the Ouchterlony double diffusion method.

SUMMARY AND EXPLANATION OF THE TEST

H. cap. and *C. imm.* are causative agents of deep-seated mycoses. These fungi present a diagnostic challenge to the microbiology laboratory and the physician.

Radiographically, lesions produced by the systemic fungi may be difficult to distinguish from tuberculous lesions, neoplasms, cancerous tissue and each other. Symptoms are often unremarkable and may mimic various pneumonias, sarcoidosis, cancer and other maladies. Culturally and histologically the organisms may be difficult to demonstrate even after repeated attempts.³⁻⁵

Frequently serology offers the only evidence available to guide treatment, suggest prognosis or lead to the selection of more definitive diagnostic techniques such as intensive culture or biopsy. In addition, quantitative serology, such as complement fixation testing, can provide important information on the effects of chemotherapy.¹¹


The fungal immunodiffusion test is a simple and reliable tool for the evaluation of suspected mycoses. Anticomplementary sera may be tested using this technique. The test also provides specificity data on reactions obtained by the complement fixation method. No expensive equipment is required and the technique is simple enough to be performed by any laboratory, thus providing an excellent screening method.

BIOLOGICAL PRINCIPLES

The Meridian Bioscience Fungal Immunodiffusion System is based on the principle of double diffusion as described by Oudin and Ouchterlony. An antibody and its homologous soluble antigen are placed in separate wells cut in an agarose diffusion medium and allowed to diffuse outward. Between the two wells, a concentration gradient of each of the reaction components is established ranging from antigen excess closest to the antibody well, to antibody excess closest to the antigen well. A visible line of precipitate forms at the point of equivalence.⁹⁻¹⁰

Antigens or antibodies may be tested for "identity" by placing a test well of the substance in question adjacent to the wells of a known system. If the antigen-antibody complexes are identical, the precipitin lines form an unbroken line of identity with the known system. Partial and nonidentity reactions are also possible. (See Figure 1)


Figure 1

A partial identity reaction occurs when certain components of the antigens (or antibodies) are identical and others are not. The "spur" represents the components which are unrelated. A nonidentity reaction will occur when the antigen-antibody complexes are different. The resulting "X" or crossed reaction indicates that two unrelated complexes are present.

The well patterns of the Meridian immunodiffusion plate are arranged to provide each patient test well with a reference known system so that identity reactions are readily apparent. (See Figure 2)

Figure 2: Arrangements of Sera and Antigens in Immunodiffusion Well Pattern.

Positive control serum, wells 1 & 4; Patient sera, wells 2, 3, 5, & 6; Antigen, well 7.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

A. **Antigens.** The following antigens are available in a 0.2 mL size as a part of the Fungal Immunodiffusion System (Catalog #100201, 100601) or separately in a 1.0 mL size. Coccidioidin antigen is also available in a 5.0 mL size.

- Histoplasmin. A purified mycelia cultural filtrate from the growth of *Histoplasma capsulatum* containing the "H" and "M" antigens. The preparation contains 100 units/mL* of each antigen. Lipids, media components, and most other antigenic components have been removed. (Note: The "H" band appears closest to the antibody well while the "M" band is closest to the antigen well.)
- Coccidioidin. A purified cultural filtrate containing the "F" antigen of *Coccidioides immitis* at a concentration of 100 units/mL*.

These antigens contain thimerosal at a final concentration of 0.01% as a preservative. Care should be exercised, however, to avoid microbial contamination.

*These units are internal Meridian standards designed only to insure lot to lot reproducibility of antigen concentration.

B. **Control Sera.** The following lyophilized control sera are provided in 0.4 mL size as a part of the Fungal Immunodiffusion System or are available separately in the 1.0 mL size. All control sera are raised in hyper-immune (not infected) goats against purified antigen preparations.

- Histoplasmosis positive control. Contains specific antibodies directed against the "H" and "M" antigens of *H. capsulatum*. (Note: The "H" band appears near the antibody well, while the "M" band is closest to the antigen well.)
- Coccidioidomycosis positive control. Contains specific antibodies directed against the "F" antigen of *C. immitis*.

These positive control sera are preserved with 0.1% sodium azide. Nevertheless, care should be exercised to avoid microbial contamination.

NOTE: Due to the standardization necessary to the production of high quality fungal serology reagents, performance of antigens, control sera and plates with materials other than those produced by Meridian Bioscience, Inc. cannot be guaranteed. The user assumes full responsibility for any modifications to the procedures published herein.

MATERIALS NOT PROVIDED

- Reagent quality water.
- Moist chamber. Any convenient container such as a petri dish, plastic box or glass jar with a tight fitting cover containing moist filter paper or paper toweling is satisfactory provided the ID plates are stationary, level and remain hydrated during incubation.

- Reading light. Darkfield plate readers are commercially available; however, a satisfactory system can be devised using a high-intensity lamp in which the front of the bulb shield-reflector is covered with black construction paper with a small hole (1-2 cm in diameter) cut for illumination. Alternatively, the plate may be held at a 45 degree to almost any bright indirect light source for adequate visualization.
- Immunodiffusion plates (Meridian Catalog Number 101009).

PRECAUTIONS

- All reagents are for *in vitro* diagnostic use only.
- When handling blood specimens, adequate measures should be taken to prevent the dissemination of etiologic agents potentially present in the specimen.
- The control sera are preserved with sodium azide. Accumulation of this chemical in metal plumbing fixtures may represent an explosion hazard. It is therefore recommended that excess control serum simply be discarded in an appropriate waste receptacle or flushed down a drain with copious amounts of water.

HAZARD and PRECAUTIONARY STATEMENTS

Refer to the SDS, available at www.meridianbioscience.com for Hazards and Precautionary Statements.

SHELF LIFE AND STORAGE

Fungal antigens are best stored at refrigerator temperature 2-8 C. Repeated freezing and thawing is detrimental to these antigens. When stored at 2-8 C the antigens are stable until the expiration as indicated on each bottle.

Positive control sera are stable in the lyophilized state until the expiration as indicated on each bottle when stored at 2-8 C. Once reconstituted it is suggested that they be aliquotted and frozen if they will not be consumed within one month. At -20 C a minimum nine month life may be expected. Repeated freezing and thawing should be avoided. When the positive control sera are in use, the period at room temperature should be kept as short as possible.

SPECIMEN COLLECTION AND PREPARATION

For optimal results, sterile serum is obtained from the blood of a patient. If a failure is encountered in specimen processing, refrigeration for up to 72 hours is permissible. Specimens may be stored up to six months at -20 C with no loss of activity provided they are not repeatedly thawed and refrozen. Specimens in transit between laboratories should be maintained at 2-8 C for optimal results. Specimens may be preserved with 0.01% thimerosal or 0.1% sodium azide if necessary.


TEST PROCEDURE

- Reconstitution of positive control sera**
 - Using a 1.0 mL pipet or tuberculin syringe add 0.38 mL (for 0.4 mL vials) or 0.95 mL (for 1.0 mL vials) of reagent quality water to the control serum vials.
 - Gently "flick" the bottom of the vials to mix.
 - Allow the vials to stand at room temperature for about one hour, occasionally inverting to mix. (Do not shake)
 - If, upon its first use, the serum appears inactive, it may be that sufficient standing time was not allowed. In this case repeat the test.
 - Aliquot and freeze any control serum which will not be consumed within one month. (See **SHELF LIFE AND STORAGE**)
- Test Procedure**
 - Please refer to figures 2, 3 and 4 for proper designation of the test well patterns.
 - Using the provided capillary pipets (Fungal Immunodiffusion System only), fill the control serum wells of an ID plate with the appropriate control sera as follows:
 - Histoplasma* control serum - Series I, wells 1 & 4.
 - Coccidioides* control serum - Series III, wells 1 & 4.

Figure 3

FUNGAL IMMUNODIFFUSION ANALYSIS (SAMPLE FORM)

Plate No. 5 Date 12-5-01

- Record the name, date, and/or lab number of the first patient on line two of the left hand column (Figure 4) of the recording form and, using a capillary tube, fill the well two of each of series I, II, III and IV with the patient's serum. Using a fresh capillary tube for each patient, repeat this procedure for additional sera in wells three, five and six of each series.
- Using a fresh capillary tube, fill the center well or series I with Histoplasmin antigen and Coccidioidin antigen (series III).
- Number and date the ID plate and the analysis form.
- Place the ID plate in a moist chamber and incubate at room temperature for 24 hours. NOTE: Bands may be evident in as little as six hours for some antigens. It will still be necessary to confirm the results of the other antigens later.
- After 24 hours read and record the ID bands on the analysis form (Figure 4). (See **Reading the Test**) An interim report should be issued at this point if no identity reactions are observed. Positive results should be reported immediately.

Figure 4

Well Number	SERIES I (Histoplasmin) Reading	SERIES III (Coccidioidin) Reading
1.	<u>Control Serum</u>	<u>Histo Control</u>
2.	<u>CARTER, R. 3287</u>	<u>H and M BANDS</u>
3.	<u>LESTER, C. 2519</u>	<u>NEGATIVE</u>
4.	<u>Control Serum</u>	<u>Histo Control</u>
5.	<u>BARNES, T. 1362</u>	<u>NEGATIVE</u>
6.	<u>HASKAMP, D. 1154</u>	<u>NEGATIVE</u>
7.	<u>Antigen</u>	<u>Coccidioidin</u>
(Name, Date, Lab No., etc.)	After 24 hour incubation, draw bands observed on the well diagram above and record the number of bands on the blanks provided. NOTE: The control bands have been drawn on the well diagram above. If these bands are not observed on the plates after 24 hour incubation, the test is not valid and must be repeated.	

- An additional 48 hour incubation period is recommended to confirm a negative result with any reagent. A final report is made at the conclusion of this period.

INTERPRETATION OF RESULTS

READING THE TEST

Any bright indirect light source may be used to aid in seeing the precipitin bands. A dark background is preferred. For instance, the plate may be held next to the shade of a high intensity lamp whose light beam is directed straight downward, over a black countertop.

A bright direct light source may also be used if the plate is held in the light beam at a 45 degree angle against a dark background.

Light boxes, providing indirect light against a dark background, prepared specifically for reading immunodiffusion plates are also commercially available.

Particular attention should be paid to the orientation of bands produced by the patient serum in relation to the control bands. A smooth junction of the bands is indicative of an identity reaction. If the control band is seen to bend toward a position in front of the patient well, it is indicative of patient antibody at a low titer.

Should no reaction be observed between the positive control and the antigen the test must be repeated.

Such reactions may be tested by flooding with 1 M sodium citrate for five minutes. Reactions which disappear after this treatment are due to C reactive protein and may be disregarded.

INTERPRETATION

A band of identity with a known positive control is indicative of patient antibody against the antigen in question. Partial identity reactions are regarded as positive for antibody against the antigen only if no other identity reaction is present on the plate. Non-identity reactions are regarded as a negative test.¹¹

REPORTING OF TEST RESULTS

If identity or partial identity reactions are observed the test may be reported as "Positive for antibodies to (name of antigen)" or "Identity reaction with (name of antigen)". Non-identity reactions and sera which are unreactive may be reported as "Negative for antibodies to (name of antigen) by immunodiffusion".

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies. User quality control is performed on each test run by the inclusion of a positive control serum for each antigen. A negative control may usually be found among the various patients tested.

The results obtained with the controls should be recorded in an appropriate log book to maintain high quality testing and comply with regulatory requirements.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

EXPECTED VALUES

In general, an identity reaction against a given antigen is indicative of active or recent past infection. Partial identity reactions are also important indicators of probable disease, especially when no identity reaction is observed among other agents. Non-identity reactions may be apparent when the disease state is caused by a mycotic agent other than the one tested. A non-identity reaction, however, is a negative test for the organism in question.

If activity against any of the fungal agents is observed, a vigorous attempt should be made to demonstrate the organism culturally for confirmation.

Histoplasmosis

The clinically significant antigens for *H. cap* are designated "H" and "M" antigens. (Note: The "H" band appears closest to the antibody well while the "M" band is closest to the antigen well.) The "M" band has been found in about 63% of patients with active histoplasmosis, while the "H" band is apparent in only 27%. The "H" band is rarely observed in the absence of antibodies against the "M" antigen. The "M" band may appear in patients who have recently recovered from histoplasmosis as well as in the serum of previous histoplasmosis patients who have been recently skin tested.

Since the test may be negative in as many as 10% of culturally demonstrated cases, the absence of either an "H" or an "M" band does not rule out histoplasmosis.^{1, 5}

Coccidioidomycosis

Antibodies directed against the "F" antigen are indicative of active or recent past (up to one year) infection. The ID test is usually positive within four weeks after infection and remains positive throughout the clinically active disease. Latex and complement fixation testing may provide important additional information regarding patient status.⁴

Cross-reactions may be seen in patients harboring other systemic fungi (especially *H. cap*.) so care must be exercised when reading for identity reactions.^{7, 11, 13}

LIMITATIONS OF THE PROCEDURE

The high rate of negative serologic tests observed among culturally demonstrable cases limits the predictive value of a negative test.

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REFERENCES

- Bauman DS and Smith CD. Comparison of immunodiffusion and complement fixation tests in the diagnosis of histoplasmosis. *J Clin Microbiol* 1975; 2:77-80.
- Crowle AJ. Immunodiffusion, Second Edition. Academic Press, New York, Chapters 1, 2 and 4. 1973.
- Heiner DC. Diagnosis of histoplasmosis using precipitin reactions in agar gel. *Pediatrics*, 1958; 22:616-627.
- Huppert M and Bailey JW. The use of immunodiffusion tests in coccidioidomycosis. *Amer J Clin Path* 1965; 44:364-368.
- Kaufman L. Value of immunodiffusion tests in the diagnosis of systemic mycotic diseases. *Ann Clin Lab Sci* 1973; 3:14146.
- Kaufman L. Practical applications of immunologic procedures for the diagnosis of opportunistic fungus infections. In: Proceedings of the Second International Conference on Opportunistic Fungal Infections. C. C. Thomas Publishing, Springfield, 1975; pp. 47-57.
- Kaufman L and Clark MJ. Value of the concomitant use of complement fixation and immunodiffusion tests in the diagnosis of coccidioidomycosis. *Appl. Microbiol.*, 1974; 28:641-43.
- Ouchterlony O. Antigen-antibody reactions in gels and the practical application of this phenomenon in the laboratory diagnosis of diphtheria. *Med. Diss.*, Stockholm. 1949.
- Ouchterlony O. Handbook of Immunodiffusion and Immunoelectrophoresis. Ann Arbor Publishers, Inc., Ann Arbor. 1968.
- Oudin J. L'analyse, immunochimique qualitative. Methode par diffusion des antigenes au sein de l'immunserum precipitant gelose. *Premier parse.* Ann Inst. Pasteur, 1948; 75:30-52.
- Palmer DF, et al., Serodiagnosis of Mycotic Diseases. C. C. Thomas Publishing, Springfield, 1977; pp 78.
- Smith CE. Coccidioidomycosis. *Med. Clin. North Amer.*, 1943; 27:790-807.
- Smith CE, Saito MT, Beard RR, Kepp KM, Clark RW and Eddie BV. Serological tests in the diagnosis of coccidioidomycosis. *Am J Hyg* 1950; 52:1-2.
- Tosh FE. Problems in the serological diagnosis of opportunistic fungus diseases. In: Proceedings of the Second International Conference on Opportunistic Fungal Infections. C. C. Thomas Publishing, Springfield, 1975; pp. 36-43.

SN10400

REV. 04/16

 Manufactured By	Meridian Bioscience, Inc. USA/Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 Telephone: 513.271.3700 Orders/Customer Service: 800.543.1980 Technical Support Center: 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124
 Authorized Representative	Meridian Bioscience Europe S. r. L Via dell'Industria, 7 20020 Villa Cortese, Milano ITALY Tel: +39 0331 43 36 36 Fax: +39 0331 43 36 16 Email: info@meridianbioscience.eu WEB: www.meridianbioscience.eu

Meridian Bioscience Europe s.a./n.v.
 2 Avenue du Japon - 1420 Braine l'Alleud
 BELGIUM
 Tel: +32 (0) 67 89 59 59
 Fax: +32 (0) 67 89 59 58
 Email: info.bn@meridianbioscience.eu

Meridian Bioscience Europe France
 34 rue de Ponthieu - 75008 Paris
 FRANCE
 Tel: +33 (0) 1 42 56 04 40
 Fax: +33 (0) 9 70 06 62 10
 Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
 Postbus 301 - 5460 AH Veghel
 NETHERLANDS
 Tel: +31 (0) 411 62 11 66
 Fax: +31 (0) 411 62 48 41
 Email: info.bn@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione. Sistema pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenzubereitung, in dem sich Probenverdünnungspuffer befindet
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Numero de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contents sufficient for n tests / Contenuto sufficiente per n saggi / Contenu suffisant pour n tests / Contenido suficiente para n ensayos / Inhalt ausreichend für n Proben	SOLN STOP	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Coniugato enzimatico / Conjugado enzimático / Enzymkonjugat
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontroltest
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	REAG	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Date di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweis
CONJ	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
SUBS	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
R. Only	Caution: US Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner. / La legge Federale negli Stati Uniti prevede che questo test sia venduto solo a un operatore autorizzato ad operare nella salute. / Attention: La vente ou la prescription de ce dispositif par un médecin agréé est soumise à des restrictions par la loi fédérale des US. / Cuidado: Los leyes federales de los EE.UU. limitan la venta de dispositivos por o por el orden de licenciados en la práctica del cuidado de la salud. / Vorzicht: Der Verkauf oder die Verschreibung dieses Geräts durch einen zugelassenen Arzt unterliegt den Beschränkungen des US Bundesgesetzes.	DET REAG	Detection Reagent / Reagente Diretto / Réactif de Détection / Reactivo de Detección / Nachweis Reagenz

For technical assistance, call Technical Support Services at (800) 343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at (800) 543-1980.