COCCIDIOIDES "F" ANTIGEN FOR FUNGAL IMMUNODIFFUSION

REF

100401, 100405,
100801, 103001, 103101

IVD

In vitro diagnostic medical device

Rx Only

INTENDED USE

Meridian Bioscience Coccidioides Immunodiffusion Reagents are standardized, purified preparations for the in vitro determination of precipitating antibodies to *Coccidioides immitis* (*C. imm.*). Reagents are available for detection of either "TP" (early) or "F" (late) antibodies. These reagents are optimized for use in the Ouchterlony double diffusion method.

SUMMARY AND EXPLANATION OF THE TEST

Coccidioides immitis is one of the causative agents of deep-seated mycoses. Though somewhat localized in the southwestern United States, modern transportation has increased the likelihood of infection for individuals visiting the region.

The disease presents a diagnostic challenge to the physician and laboratorian. Radiographically, lesions produced by the systemic fungi may be difficult to distinguish from each other, tuberculous lesions, neoplasms, (benign and malignant). Symptoms are often unremarkable and can mimic various pneumonias, sarcoidosis, cancer and other maladies. Culturally and histologically the organisms can be difficult to demonstrate, even after repeated attempts.³

Frequently serology offers the only evidence available to guide treatment, suggest prognosis or lead to the selection of more definitive diagnostic techniques such as intensive culture or biopsy. In addition, quantitative serology, such as complement fixation testing, can provide important information on the effects of chemotherapy.¹⁸

The Coccidioides Immunodiffusion Systems are simple and reliable tools for the evaluation of suspected mycoses. Anticomplementary sera may be tested with these procedures. The tests also provide specificity data on reactions obtained by the complement fixation method. No expensive equipment is required and the procedures are simple enough to be performed by any laboratory, thus providing excellent screening methods.


Use of the immunodiffusion, coccidioides "F" (IDCF) test alone for the serodiagnosis of coccidioidomycosis is not recommended since, in general, only IgG antibody is detected. When used together with the Coccidioides Latex Agglutination System for *C. imm* antibodies, 93% of active coccidioidomycosis patients were positive in one or both tests (CL or ID). In addition, IgM antibody against *Coccidioides immitis* may be detected using the Meridian Coccidioides IDTP antigen and control serum set.^{4,6,12}

BIOLOGICAL PRINCIPLES

The Meridian Bioscience Coccidioides Immunodiffusion System is based on the principle of double diffusion as described by Oudin and Ouchterlony. An antibody and its homologous soluble antigen are placed in separate wells cut in an agarose diffusion medium and allowed to diffuse outward. Between the two wells, a concentration gradient of each of the reaction components is established ranging from antigen excess closest to the antibody well, to antibody excess closest to the antigen well. A visible line of precipitate forms at the point of equivalence.¹⁰⁻¹²

Antigens or antibodies may be tested for "identity" by placing a test well of the substance in question adjacent to the wells of a known system. If the antigen antibody complexes are identical, the precipitin lines form an unbroken line of identity with the known system. Partial and nonidentity reactions are also possible. (See figure 1)


Figure 1. Identification of Immunodiffusion Bands

A partial identity reaction occurs when certain components of the antigens (or antibodies) are identical and others are not. The "spur" represents the components which are unrelated. A nonidentity reaction will occur when the antigen-antibody complexes are different. The resulting "X" or crossed reaction indicates that two unrelated complexes are present.

The well patterns of the Meridian immunodiffusion plate are arranged to provide each patient test well with a reference known system so that identity reactions are readily apparent. (see figure 2)

Figure 2. Arrangement of sera and antigens in immunodiffusion well pattern.

Early antibody, (IgM) is detected by the "TP" antigen. Late antibody, (IgG) is detected by the "F" antigen. The results obtained with the "TP" antigen parallel that of tube precipitin and latex tests, while the results obtained with the "F" antigen are similar to those observed with complement fixation.^{14, 16, 18}

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

- Antigen:** Coccidioides antigens are available in 1.0 mL and 5.0 mL sizes. The available antigens consist of purified cultural filtrates containing either the "F" or "TP" antigen of *Coccidioides immitis* at a concentration of 100 units/mL. The antigens contain thimerosal at a final concentration of 1:10,000 as a preservative. Care should be exercised, however, to avoid microbial contamination. *These units are internal Meridian standards designed only to ensure lot to lot reproducibility of antigen concentration.
- Control Sera:** Anti-Coccidioides lyophilized control serum is provided in a 1.0 mL size. Control serum contains goat, rabbit or human serum with antibodies directed against either the "F" or "TP" antigen of *C. immitis* (see PRECAUTIONS). These positive control sera are preserved with 0.1% sodium azide. Care should be exercised to avoid microbial contamination.
- Immunodiffusion Plates:** Twelve ID plates are provided. When four sera per plate are tested, sufficient materials are provided to test 48 patient sera. ID plates may be obtained by ordering Meridian Catalog #101012. The immunodiffusion plates consist of 0.9% agarose in glycine/phosphate buffer at pH 7.4 + 0.2.

Note: Recent advances in plate formulation have eliminated the need for Immunodiffusion Band Intensifying Fluid.

Note: Due to the standardization necessary to the production of high quality Meridian Bioscience fungal serology reagents, performance of tests with antigens, control sera and plates with materials other than those produced by Meridian Bioscience, Inc. cannot be guaranteed. The user assumes full responsibility for any modification to the procedures published herein.

MATERIALS AND EQUIPMENT REQUIRED BUT NOT SUPPLIED

- Reagent quality water.

- Moist chamber; any convenient container such as a petri dish, plastic box or glass jar with a tight fitting cover containing moist filter paper or paper toweling is satisfactory provided the ID plates are stationary, level and remain hydrated during incubation.
- Reading light. Dark field plate readers are commercially available; however, a satisfactory system can be devised using a high intensity lamp in which the front of the bulb shield reflector is covered with black construction paper with a small hole (1-2 cm in diameter) cut for illumination. Alternatively, the plate may be held at a 45 degree angle to almost any bright indirect light source for adequate visualization.

PRECAUTIONS

- All reagents are for in vitro diagnostic use only.
- Patient specimens and control serum may contain infectious agents. The control serum may contain human serum which was screened for HBsAg and antibody to HIV and found to be negative. However, no test can offer complete assurance that human blood will not transmit HIV, hepatitis or other infectious agents. Therefore, patient specimens, control serum, used ID plates and other potentially hazardous material should be handled at Biosafety Level 2 as recommended in the CDC/NIH manual "Biosafety in Microbiology and Biomedical Laboratories".
- The control sera are preserved with sodium azide. Accumulation of this chemical in metal plumbing fixtures may represent an explosion hazard. It is therefore recommended that excess control serum simply be discarded in an appropriate waste receptacle or flushed down a drain with copious amounts of water.

HAZARD and PRECAUTIONARY STATEMENTS

Refer to the SDS, available at www.meridianbioscience.com for Hazard and Precautionary Statements.

SHELF LIFE AND STORAGE

Coccidioides antigens are best stored at refrigerator temperature (2-8 C). Repeated freezing and thawing is detrimental to these antigens. Extreme heat (more than 60 C) is also to be avoided. When stored at 2-8 C the antigens are stable at least 24 months from the date of manufacture.

Positive control sera are stable in the lyophilized state for a minimum of 24 months from the date of manufacture when stored at 2-8 C. Once reconstituted it is suggested that they be aliquotted and frozen if they will not be consumed within one month. At -20 C a minimum nine month life may be expected. Repeated freezing and thawing should be avoided. When the positive control sera are in use, the period at room temperature should be kept as short as possible.

Immunodiffusion plates are individually packaged and if unopened and stored at 2-8 C, an 18 month minimum shelf life may be expected. If plate packages are opened and the plates must be stored before use, they must be placed in a tightly sealed, humid container to avoid drying. A resealable plastic container with a small moist sponge or towel is ideal.

SPECIMEN COLLECTION AND PREPARATION

- For optimal results, sterile serum is obtained from the patient's blood.
- If a delay is encountered in specimen processing, refrigeration for up to 72 hours is permissible. Specimens may be stored up to six months at -20 C with no loss of activity, provided they are not repeatedly thawed and refrozen. Specimens in transit between laboratories should be maintained at 2-8 C for optimal results. Specimens may be preserved with 1:10,000 thimerosal or 0.1% sodium azide if necessary.
- Greater sensitivity can be achieved when performing the "TP" antibody test by concentrating the patient's serum approximately 8:1.¹³⁻¹⁶ This can be accomplished by pervaporation or use of a 125000 molecular weight cutoff membrane filter (e.g. Amicon™ B-125 or equivalent). While various strategies are available for pervaporation, the serum may be simply placed in a short length of dialysis tubing and air dried or gently warmed with a hand dryer until the correct volume reduction is obtained.
- In a limited clinical study, a 15% increase in sensitivity was obtained when specimens were concentrated by pervaporation and allowed to prediffuse for two hours prior to addition of antigen. (Pappagianis, et al., personal communication)

TEST PROCEDURE

A. Reconstitution of positive control serum.

- Using a 1.0 mL pipet or tuberculin syringe add 0.95 mL of reagent quality water to the control serum vials.
- Gently "flick" the bottom of the vials to mix.
- Allow the vials to stand at room temperature for about 1 hour, occasionally inverting to mix. (Do not shake)
- If, on its first use, the serum appears inactive, it may be that sufficient standing time was now allowed. In this case repeat the test.
- Aliquot and freeze any control serum which will not be consumed within one month. (See SHELF LIFE AND STORAGE)

B. Test procedure

- Please refer to figures 2, 3 and 4 for proper designation of the test well pattern.
- Using the provided capillary pipets (Coccidioides Immunodiffusion System only), fill the control serum wells (1 & 4) of an ID plate.

Figure 3.

FUNGAL IMMUNODIFFUSION ANALYSIS
(SAMPLE FORM)

Plate No. 5 Date 8-23-02

TEST SERIES I

- Record the name, date, and/or lab number of the first patient on line 2 of the left hand column (figure 4) of the recording form (provided in the Coccidioides Immunodiffusion System only) and, using a capillary tube, fill well 2 with the patient's serum. Using a fresh capillary tube for each patient, repeat this procedure for additional sera in wells 3, 5 and 6.
- Using a fresh capillary tube, fill the center well with coccidioides antigen.
- Number and date the ID plate and the analysis form.
- Place the ID plate in a moist chamber and incubate at room temperature for 24 hours. Note: bands may be evident in as little as three to six hours for some sera. It will still be necessary to confirm the results of the other sera later.
- After 24 hours read and record the ID bands on the analysis form (figure 4). (See Reading the Test) An interim report should be issued at this point if no identity reactions are observed. Positive results should be reported immediately.

Figure 4.

Well Number	TEST
1. Control Serum	Cocci Control
2. CARTER, R. 3287	NEGATIVE
3. LESTER, C. 2519	NEGATIVE
4. Control Serum	Cocci Control
5. BARNES, T. 1362	NEGATIVE
6. HASKAMP, D. 1154	I-BAND
7. Antigen	Coccidioidin

(Name, Date, Lab No., etc.)

After 24 hour incubation, draw bands observed on the well diagram above and record the number of bands on the blanks provided.
NOTE: The control bands have been drawn on the well diagram above. If these bands are not observed on the plates after 24 hour incubation, the test is not valid and must be repeated.

8. An additional 48 hour incubation period is recommended to confirm a negative result. A final report is made at the conclusion of this period.

READING THE TEST

Any bright indirect light source may be used to aid in seeing the precipitin bands. A dark background is preferred. For instance, the plate may be held next to the shade of a high intensity lamp whose light beam is directed straight downward, over a black countertop. A bright direct light source may also be used if the plate is held in the light beam at a 45 degree angle against a dark background. Light boxes, providing indirect light against a dark background and prepared specifically for reading immunodiffusion plates, are also commercially available.

Particular attention should be paid to the orientation of bands produced by the patient serum in relation to the control bands. A smooth junction of the bands indicates an identity reaction. When patient antibody levels are low, only a slight bend in the control band toward a position in front of the patient well may be apparent. Should no reaction be observed between the positive control and the antigen, the test must be repeated.

INTERPRETATION OF RESULTS

A band of identity with the positive control indicates patient antibody against *Coccidioides immitis*. Partial identity reactions are also regarded as positive for antibody. Non-identity reactions are regarded as a negative test.¹¹

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies. User quality control is performed on each test run by the inclusion of a positive control serum for each antigen. A negative control may usually be found among the various patients tested. The results obtained with the controls should be recorded in an appropriate log book to maintain high quality testing and comply with regulatory requirements. **If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Service Department at 1-800-343-3858 (US) or your local distributor.**

EXPECTED VALUES

Antibodies directed against the "F" antigen indicate active or recent past (up to 1 year) infection. The ID test is usually positive within four weeks after infection and remains positive throughout the clinically active disease.

Antibodies against the "TP" antigen indicate recent, active infection. These are primarily IgM antibodies which appear within several days of infection. This ID test is usually positive within 1 week after infection and reverts to negative within a few months.

Latex and complement fixation testing may provide important additional information regarding patient status.⁴⁻⁶

Cross-reactions may be seen in patients harboring other systemic fungi (especially *H. cap.*) so care must be exercised when reading for identity reactions^{9, 18, 21} (See **LIMITATIONS OF THE PROCEDURE**). Partial identity reactions are also important indicators of probable disease. Nonidentity reactions may be apparent when the disease state is caused by a mycotic agent other than the one tested. A nonidentity reaction, however, is a negative test for the organism in question.

If activity against *Coccidioides immitis* is observed, a vigorous attempt should be made to demonstrate the organism culturally for confirmation.

LIMITATIONS OF THE PROCEDURE

The high rate of negative serologic tests observed among certain culturally demonstrable cases limits the predictive value of a negative test.²⁰ Cross-reactions may be observed in patients harboring infections caused by other systemic fungi. While these may be easily distinguished from the identity lines produced against coccidioidin, their true identity will not be apparent in a system limited to one antigen.

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REFERENCES

- Bauman, DS. and CD. Smith. Comparison of Immunodiffusion and complement fixation tests in the diagnosis of histoplasmosis. *J. Clin. Microbiol.* 1975;2:77-80.
- Crowle, AJ. *Immunodiffusion*. Second Edition. Academic Press. New York. 1973; Chapters 1, 2 and 4.
- Heiner, LC. Diagnosis of histoplasmosis using precipitin reactions in agar gel. *Pediatrics* 1958;22:616-627.
- Huppert, M. and JW. Bailey. The use of immunodiffusion tests in coccidioidomycosis 1. The accuracy and reproducibility of the immunodiffusion test which correlates with complement fixation. *Amer. J. Clin. Path.* 1965;44:364-368.
- Huppert, M. and JW. Bailey. The use of immunodiffusion tests in coccidioidomycosis 11. An immunodiffusion test as a substitute for the tube precipitin test. *Amer. J. Clin. Path.* 1965;44:369-373.
- Kaufman, L. Value of immunodiffusion tests in the diagnosis of systemic mycotic diseases. *Ann. Clin. Lab. Sci.* 1973;3:141-146.
- Kaufman, L. Practical applications of immunologic procedures for the diagnosis of opportunistic fungus infections. In: *Proceedings of the Second International Conference on Opportunistic Fungal Infections*. CC. Thomas Publishing, Springfield. 1975; pp 47-57.
- Kaufman, L. and MJ. Clark. Value of the concomitant use of complement fixation and immunodiffusion tests in the diagnosis of coccidioidomycosis. *Appl. Microbiol.* 1974;28:641-43.
- Ouchterlony, O. Antigen-antibody reactions in gels and the practical application of this phenomenon in the laboratory diagnosis of diphtheria. *Med. Diss. Stockholm*, 1949.
- Ouchterlony, O. *Handbook of Immunodiffusion and Immunelectrophoresis*. Ann Arbor Publishers Inc., Ann Arbor, 1968.
- Oudin, J. L'analyse immunochimique qualitative. Methode par diffusion des antigenes au sein de l'immunserum precipitant gelose. *Premiere parse. Ann Inst. Pasteur* 1948;75:30-52.
- Palmer, DF. et al. *Serodiagnosis of Mycotic Diseases*. CC. Thomas Publishing, Springfield. 1977; pp. 7-18.
- Pappagianis, D. Serology and Serodiagnosis of Coccidioidomycosis. In: *Coccidioidomycosis*. David A. Stearns, M.D., Editor, Plenum Publishing.
- Pappagianis, D. Immunoserologic Procedures in the Diagnosis and Management of Pulmonary Mycoses. In: *Handbook of Fungus Diseases*. H. Einstein, Editor. Amer. College of Chest Physicians. 1981.
- Pappagianis, D. Antibody in Cerebrospinal Fluid in Non-meningitic Coccidioidomycosis. *Sabouraudia* 1972;10:173-179.
- Pappagianis, D., NJ. Lindsey, CE. Smith, and M. Saito. Antibodies in Human Coccidioidomycosis: Immunoelectrophoretic Properties. *Proceedings of the Society for Experimental Biology and Medicine* 1985;118:118-122.
- Pappagianis, D. False-positive Reactions of Cerebrospinal Fluid and Diluted Sera with the Coccidioid Latex-agglutination Test. *Am. J. Clin. Path.* 1976; 68:916-921.
- Smith, CE. Coccidioidomycosis. *Med. Clin. North Amer.* 1943; 27:790-807.
- Smith, CE., MT. Salito, RR. Beard, KM. Kepp, RW. Clark, and BV. Eddie. Serological tests in the diagnosis of coccidioidomycosis. *Am. J. Hyg.* 1950;52:1-2.
- Tosh, FE. Problems in the serological diagnosis of opportunistic fungus diseases. In: *Proceedings of the Second International Conference on Opportunistic Fungal Infections*. C. C. Thomas Publishing, Springfield, 1975; pp. 36-43.

SN10000

REV. 04/16

<p>Manufactured By</p>	<p>Meridian Bioscience, Inc. USA/Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 Telephone: 513.271.3700 Orders/Customer Service: 800.543.1980 Technical Support Center: 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124</p>
	<p>Meridian Bioscience Europe S. r. l. Via dell'Industria, 7 20020 Villa Cortese, Milano ITALY Tel: +39 0331 43 36 36 Fax: +39 0331 43 36 16 Email: info@meridianbioscience.eu WEB: www.meridianbioscience.eu</p>

Meridian Bioscience Europe s.a./n.v.
2 Avenue du Japon - 1420 Braine l'Alleud
BELGIUM
Tel: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
Email: info.bn@meridianbioscience.eu

Meridian Bioscience Europe France
34 rue de Ponthieu - 75008 Paris
FRANCE
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 9 70 06 62 10
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
Postbus 301 - 5460 AH Veghel
NETHERLANDS
Tel: +31 (0) 411 62 11 66
Fax: +31 (0) 411 62 48 41
Email: info.bn@meridianbioscience.eu

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EC.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünngungspuffer befindet
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Numero de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <=> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" essais / Contenido suficiente para "n" ensayos / Inhalt ausreichend für <=> Prüfungen	SOLN STOP	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stoplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Coniugato enzimatico / Conjugé enzymatique / Conjugado enzimático / Enzymkonjugat
SN	Serial number / Numero di serie / Numéro de série / Numero de serie / Seriennummer	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	REAG	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
CONJ	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillon / Diluyente de muestra / Probenverdünnungspuffer
SUBS	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
R. Only	Caution: US Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner. / La Legge Federale negli Stati Uniti prevede che questo test sia venduto solo a un operatore abilitato ad operare nella salute. / Attention: La vente ou la prescription de ce dispositif par un médecin agréé est soumise à des restrictions par la loi fédérale des US. / Cuidado: Los dispositivos de los EU limitan la venta de dispositivos por o por la orden de licenciados en la práctica del cuidado de la salud. / Vorsicht: Der Verkauf oder die Verschreibung dieses Geräts durch einen zugelassenen Arzt unterliegt den Beschränkungen des US Bundesgesetzes.	DET REAG	Detection Reagent / Reagente Diretta / Réactif de Détection / Reactivo de Detección / Nachweis Reagenz

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de símbolos, Zeichenerklärung)