merifluor®

Pneumocystis

Direct Immunofluorescent Detection Procedure for *Pneumocystis carinii* in Respiratory Tract Specimens

REF 222030

IVD

Rx Only

INTENDED USE

MERIFLUOR Pneumocystis is an in vitro procedure for the detection of *Pneumocystis carinii* cysts and trophozoites in direct smears of respiratory tract specimens.

SUMMARY AND EXPLANATION OF THE TEST

Pneumocystis carinii is a unicellular fungus which causes life threatening pneumonia in severely immuno-compromised individuals.^{1,2} The available evidence indicates that most children experience mild or subclinical infection with *P. carinii* and that the organism may persist in a latent state until reactivation occurs following immuno-suppression.^{1,3} Until recently, the populations at risk have been premature infants, patients with hematopoietic malignancies or immunoglobulin defects, and patients receiving immunosuppressive treatment such as cortico-steroids for cancer or organ transplantation.^{1,2,4,5,6} Since 1981, the inherent susceptibility of individuals with acquired immune deficiency syndrome (AIDS) has significantly increased the prevalence of pneumocystis pneumonia. Pneumocystis pneumonia develops in 60-85% of AIDS patients and is by far the most common opportunistic infection in this group.^{1,7,8}

Three developmental forms of *P. carinii* are discernible with light microscopy: the trophozoite, the sporozoite, and the cyst.^{1,9} Trophozoites, 2-5 µm in diameter, are pleomorphic and tend to cluster together in masses. Trophozoites undergo maturation to become cysts. The cyst, 5-8 µm in diameter, is the most easily recognized form of *P. carinii*. Cysts contain up to eight crescent or pleomorphic-shaped sporozoites, 1-2 µm in diameter, and tend to cluster in masses within an extracellular matrix material.

Currently, a clinical diagnosis of *P. carinii* pneumonia is based on radiologic evidence and direct microscopic observation of respiratory specimens using nonspecific techniques.^{1,9} Nonspecific staining methods, such as methenamine silver nitrate and toluidine blue, stain the cyst wall but have no affinity for trophozoites or sporozoites.^{1,9} Alternatively, the Wright Giemsa and Gram staining methods stain trophozoites and sporozoites, but fail to stain cyst walls.^{1,9} Due to the nonspecific nature of these stains proper identification of the organism has been difficult.^{10,13} Detection of *P. carinii* in bronchoalveolar lavage fluids, induced-sputum and impression smears of lung tissue has been accomplished by immunofluorescent assay.^{11,12,13}

MERIFLUOR Pneumocystis provides a simple, highly specific, direct immunofluorescent procedure for the identification of *P. carinii* cysts and trophozoites in respiratory tract specimens.

BIOLOGICAL PRINCIPLES

MERIFLUOR Pneumocystis utilizes the principle of direct immunofluorescence. The Detection Reagent contains FITC labeled monoclonal antibodies directed against cell wall and matrix antigens of *P. carinii* cysts, sporozoites and trophozoites. Prepared slides of respiratory specimens are treated with the Detection Reagent. The monoclonal antibodies attach to pneumocystis antigens present in the specimen. The slides are rinsed to remove unbound antibodies, mounted with MERIFLUOR Mounting Medium, and examined for bright apple-green fluorescence and characteristic morphology of *P. carinii* cysts using a fluorescence microscope. Sporozoites, trophozoites and the extracellular matrix should also fluoresce. Background material present in the specimen is counterstained orange to red.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. Detection Reagent - FITC labeled anti-*P. carinii* monoclonal antibodies in a buffered solution containing protein, Evans Blue and 0.07% sodium azide.
2. Mounting Medium - Buffered glycerol containing an antifluorescent and 0.05% sodium azide.
3. Microscope Slides
4. Transfer pipettes

MATERIALS NOT PROVIDED

1. *P. carinii* positive control slides of human origin
2. *P. carinii* negative control slide
3. Centrifuge
4. Humidity chamber
5. Wash bottle
6. Microscope slide coverslips
7. Acetone, reagent grade (store in tightly closed container to prevent absorption of water).
8. 35-37 °C incubator
9. Deionized water or tap water
10. Pipettes
11. Dithiothreitol
12. Fluorescence microscope, equipped with a filter system for fluorescein isothiocyanate (FITC) with the following parameters: Excitation wavelength - 490-500, Barrier filter - 510-530.

PRECAUTIONS

1. All reagents and components are for in vitro diagnostic use only.
2. Do not mix reagents from different kit lot numbers.
3. Do not use kit components beyond the indicated expiration date on the kit label.
4. Protect the Detection Reagent from exposure to light.
5. All reagents should be inverted gently before use.
6. Clean microscope slides prior to use with 95% ethanol.
7. Patient specimens may contain HIV or other infectious agents and should be handled by properly trained personnel and disposed of as potential biohazards. Wear disposable gloves while handling specimens and performing the test procedure.
8. Do not allow the specimen to dry after the Detection Reagent is added.
9. If specimen material is not seen upon scanning the slide wells, loss of sample may have occurred. This is usually due to an overly vigorous wash procedure, insufficient fixation of the specimen, or insufficient cleaning of the slide.
10. The Detection Reagent contains Evans Blue which is a potential carcinogen. Avoid contact with skin. Wash area thoroughly if contact occurs.

WARNING: The reagents in this kit contain sodium azide, which is a skin irritant. Avoid skin contact with the kit components. Disposal of reagents containing sodium azide into lead or copper plumbing can result in the formation of explosive metal azides. This can be avoided by flushing with a large volume of water.

HAZARD and PRECAUTIONARY STATEMENTS

 ! / Counting Medium	<p>Signal Word Danger</p> <p>Hazard Statements H317 – May cause an allergic skin reaction H341 – Suspected of causing genetic defects H350 – May cause cancer Contains Formaldehyde</p> <p>Precautionary Statements - EU (§28, 1272/2008) P280 – Wear eye protection/face protection P201 – Obtain special instructions before use P281 – Use personal protective equipment as required P308 + P313 – If exposed or concerned: Get medical advice/attention P321 – See SDS Section 4 or Section 11 for additional medical treatment information</p>
--	--

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-8 °C and return the kit promptly to the refrigerator after each use. **CAUTION: DO NOT FREEZE.**

SPECIMEN COLLECTION AND PREPARATION

Bronchoalveolar lavage (**BAL**), bronchial wash (**BW**) and induced sputum (**IS**) specimens should be collected using standard procedures. Pretreatment of specimens with fixatives such as formaldehyde, paraformaldehyde, or ethanol is not recommended, as this can cause a decrease or loss of fluorescence.

The following chart summarizes the time frame in which specimens are suitable for processing with MERIFLUOR Pneumocystis:

BAL/BW	Unpreserved	About 1 month at 4 C Stain within 8 hours or freeze
	Acetone fixed on slide	
Induced Sputa	Digested or non-digested	1 day at 4 C 1-2 months
	Digested & frozen	
	Acetone fixed on slide	Stain within 8 hours or freeze
Tissue	Unfixed	Freeze immediately or do touch prep and

A. Preparation of BW and BAL Specimens

1. Centrifugation Method
 - a. Specimens should be concentrated to increase the recovery of *P. carinii*. Mix the specimen thoroughly, transfer a minimum of 5 mL to a centrifuge tube and centrifuge at 1,800 xg for 10 minutes. Remove and discard all but 0.5 mL of the supernatant. Thoroughly resuspend the pellet in the remaining 0.5 mL of fluid.
 - b. Remove 25-50 µL of the resuspended specimen with a transfer pipette (marking closest to tip is 50 µL) and place it on the appropriate well of a microscope slide.
 - c. Spread the specimen over the well using the side of the transfer pipette.
 - d. Continue immediately with staining (see PROCEDURE).
2. Shandon Cytospin™ 2 Method
 - a. The concentration procedure may be omitted by centrifuging 0.5 mL-0.6 mL of unconcentrated **BAL** or **BW** at 900 rpm for 5 minutes using 1 white and 1 tan filter. The marking closest to the bulb of the transfer pipette (supplied with kit) is 0.3 mL.
 - b. Continue immediately with staining (see PROCEDURE).

B. Preparation of IS Specimens

Induced sputum specimens can be processed by preparing smears directly from sputa. However, increased recovery can be obtained by concentration after liquefaction of the sputa with dithiothreitol.^{11, 13}

1. Direct smears from **IS**.
 - a. With an applicator stick, gather all the opaque, white or cream-colored flecks of material present in the specimen to the periphery of the specimen.
 - b. Remove the flecks with a transfer pipette and place them on the appropriate well of a microscope slide.
 - c. Spread these flecks and their associated material over the well using the side of the transfer pipette.
 - d. Continue immediately with staining (see PROCEDURE).

C. Preparation of Biopsy Specimens

Transbronchial, open lung, or other biopsy specimens must not be formalin-fixed. Transport specimen in a container on gauze slightly moistened with saline. The tissue should not be immersed in saline nor allowed to dry as either condition can result in poor touch preparations.

1. Prepare a freshly cut surface on a fragment of tissue. Touch the cut surface of the specimen to an appropriately labeled microscope slide. Make several non-overlapping imprints within the well, avoiding smearing. The use of several cut surfaces will increase the opportunity to detect *P. carinii* organisms.
2. While imprints are still moist on the slide, fix by adding 1-2 drops of acetone and allow to air dry completely. Proceed with step 2 of PROCEDURE.

PREPARATION OF POSITIVE CONTROL

A positive control should be included each time the test is performed. Bronchoalveolar lavage or bronchial wash specimens previously confirmed as positive for *P. carinii* and which have not been treated with fixatives such as formaldehyde, paraformaldehyde or ethanol, may be used. Prepare control slides as described in SPECIMEN PREPARATION. Slides which have been acetone fixed may be frozen at -20 C for up to one month. It is not necessary to repeat acetone fixation on thawed slides. Control slides may also be obtained commercially. Positive control specimens must exhibit the staining characteristics of a positive specimen as described in INTERPRETATION OF RESULTS.

TEST PROCEDURE

1. Allow specimen to air dry, then fix by adding 1-2 drops of acetone to specimen on slide and allow to air dry completely. **NOTE: Specimens must be tested within 8 hours.**
2. Use a separate transfer pipette to place 50 µL (marking closest to tip) of Detection Reagent in each specimen and control well. With the transfer pipette spread the reagent over the entire well taking care not to disturb the specimen.
3. Incubate the slide in a humidified chamber for 30 minutes at 35-37 C.
4. Tap the slide gently on a paper towel to remove excess reagent. Use a wash bottle to rinse the slide with a gentle stream of water until the excess Detection Reagent has been removed. **NOTE: Do not submerge the slide during rinsing.** Exercise care during the rinse process to avoid disturbing the specimen or causing cross contamination of the specimens.
5. Remove excess water by tapping the slide on a clean paper towel. **NOTE: Do not allow the slide to dry.**
6. Add one or two drops of mounting medium to each slide and apply a coverslip.
7. Examine each well thoroughly using a fluorescent microscope at 200-300X magnification. Protect slide from light. **NOTE: The slide should be read immediately. Otherwise store the slide at 2-8 C in the dark and read within 24 hours.**

INTERPRETATION OF RESULTS

1. Positive test result: Any specimen which contains two typical cysts exhibiting apple-green fluorescence of characteristic morphology should be considered positive for the presence of *P. carinii*. Generally cysts, each 5-8 µm in diameter, are found together with trophozoites in clusters. Clusters can be variable in size and may appear with or without "honeycomb" like structures. Some cysts fluoresce evenly throughout their structure whereas other cysts may fluoresce mainly on their periphery and produce a "honeycomb" appearance within the clusters.^{12, 13} If no cysts are observed and symptoms persist, another specimen should be obtained and examined.
2. Negative test result: A specimen which fails to exhibit the characteristic fluorescence described above.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

1. A positive control specimen of human origin should be obtained to verify the proper performance of the kit.

2. At the time of each use the kit components should be visually examined for obvious signs of contamination, freezing or clumping.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

EXPECTED VALUES

The MERIFLUOR Pneumocystis direct immuno-fluorescent test will identify the presence of *P. carinii* in respiratory specimens. The incidence of *P. carinii* infection varies between hospitals and patient populations. *Pneumocystis pneumonia* develops in 60-85% of untreated AIDS patients but is likely to be lower in populations receiving prophylactic antimicrobial therapy.

LIMITATIONS OF THE PROCEDURE

MERIFLUOR Pneumocystis is an in vitro diagnostic test to aid in the diagnosis of *P. carinii* infection. The presence of *P. carinii* is associated with pneumonia in immunosuppressed individuals. In addition, the presence of cysts in a specimen does not preclude the existence of other microorganisms or another underlying condition as the causative agent of a patient's symptoms. The final diagnosis should be made by a clinician familiar with the history of the patient, the results of a physical examination, and the laboratory data. NOTE: Meridian Bioscience has not evaluated the suitability of commercially marketed liquefaction reagents for treating sputa that will be tested by MERIFLUOR Pneumocystis. It is the users responsibility for validating the compatibility of specimens treated with such reagents with the reagents supplied as part of the MERIFLUOR assay.

SPECIFIC PERFORMANCE CHARACTERISTICS

Induced sputa and bronchoalveolar lavages were collected at three hospitals from patients suspected of having *Pneumocystis carinii* pneumonia. All specimens were stained with the MERIFLUOR Pneumocystis test kit and with a histologic stain: Wright-Giemsa (DiffQuik™) or Gomori methenamine silver. Specimens from two sites were also stained with a competitor's direct immunofluorescent test kit. The results of the 65 specimens included in the analysis are shown in Table 1.

Table 1

MERIFLUOR Pneumocystis	Histologic Stain		Competitor's Assay	
	+	-	+	-
	44	1*	31	0
	2	18	0	12
Sensitivity	96%			
Specificity	95%			
Agreement with histologic stain	95%			
Agreement with competitor assay	100%			

*A bronchoalveolar lavage obtained from this patient 24 hours later was positive by both histologic and immunofluorescent stains, indicating an increased sensitivity for the MERIFLUOR Pneumocystis test.

The following organisms from culture did not fluoresce when tested with MERIFLUOR Pneumocystis:

Aspergillus flavus, *Aspergillus fumigatus*, *Aspergillus niger*, *Blastomyces dermatitidis*, *Candida albicans*, *Candida guilliermondii*, *Candida krusei*, *Candida pseudotropicalis*, *Candida tropicalis*, *Coccidioides immitis*, *Cryptococcus neoformans*, *Cryptosporidium* spp., *Histoplasma capsulatum*, *Legionella pneumophila*, *Nocardia asteroides*, *Saccharomyces cerevisiae*, *Toxoplasma gondii*

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REFERENCES

1. MacFarlane, JT and RG Finch, Pneumocystis carinii pneumonia 1985, Thorax 40: 561-570.
2. Walzer, PD, DP Perl, DJ Krogstad, PG Rawson, and MG Schultz. Pneumocystis carinii pneumonia in the United States. Epidemiologic, diagnostic and clinical features. Ann.Intern.Med 1974; 80:83-93.
3. Stagno, S., LL Pifer, WT Hughes, DM Brasfield and RE Tiller. Pneumocystis carinii pneumonitis in young immunocompetent infants. Pediatrics 1980; 66:56-62.
4. Hughes, WT, RA Price, H Kim, TP Coburn, D Girsby and S Feldman. Pneumocystis carinii pneumonitis in children with malignancies. J.Pediatr.1973; 82:404-415.
5. Pifer, LL, HB Niell, BJ Morrison, JD Counce, Jr, JM Freeman, DR Woods and CL Neely. Pneumocystis carinii antigenemia in adults with malignancy, infection, or pulmonary disease. J.Clin.Microbiol 1984; 20:887-890.
6. Hardy, AM, CP Wajszczuk, AF Suffredini, TR Hakala and M Ho. Pneumocystis carinii pneumonia in renal-transplant recipients treated with cyclosporine and steroids. J.Infect.Dis. 1984; 149:143-147.
7. Brenner, M, FP Ognibene, EE Lack, JT Simmons, AF Suffredini, HC Lane, AS Fauci, JE Parrillo, JH Shelhamer and H Masur. Prognostic factors and life expectancy of patients and acquired immunodeficiency syndrome and Pneumocystis carinii pneumonia. Am.Rev.Respir.Dis. 1987; 136:1199-1206.
8. Posner, D and FA Kahn. Respiratory infection in AIDS: A comprehensive care approach. J.Respir.Dis. 1987; 8:83-94.
9. Hughes, WT In Pneumocystis carinii Pneumonitis. Vol 182 CRC Press Inc., Boca Ratón, FL 1987; p 9-32.
10. Milder, JE PD Walzer, JD Coonrod and ME Rudledge. Comparison of histological and immunological techniques for detection of Pneumocystis carinii rat bronchial lavage fluid. J.Clin.Microbiol, 1980; 11:409-417.
11. Lim, SK, WC Eveland and RJ Porter. Direct fluorescent antibody method for the diagnosis of Pneumocystis carinii pneumonitis from sputa or tracheal aspirates from humans. Appl.Microbiol. 1974; 27:144-149.
12. Gill, VJ, G Evans, F Stock, JE Parrillo, H Masur and JA Kovacs. Detection of Pneumocystis carinii by fluorescent antibody stain using a combination of three monoclonal antibodies. J.Clin.Microbiol. 1987;25:1837-1840.
13. Kovacs, JA, V Ng, H Masur, G Leoung, WK Hadley, G Evans, HC Lane, FP Ognibene, J Shelhamer, JE Parrillo and VJ Gill. Diagnosis of Pneumocystis carinii: Improved detection in sputum using monoclonal antibodies. N.Eng.J.Med. 1988; 318:589-593.

SN11060

REV. 07/20

Manufactured By

Meridian Bioscience, Inc.
Corporate Office
3471 River Hills Drive
Cincinnati, Ohio 45244 USA
Telephone: 513.271.3700
Orders/Customer Service:
800.543.1980
Technical Support Center:
800.343.3858
Information Fax: 513.272.5432
Ordering Fax: 513.271.0124

EC REP

Authorized Representative

Meridian Bioscience Europe S. r. L
Via dell'Industria, 7
20035 Villa Cortese, Milano
ITALY
Tel: +39 0331 43 36 36
Fax: +39 0331 43 36 16
Email: info@meridianbioscience.eu
WEB: www.meridianbioscience.com/eu

Meridian Bioscience Europe s.a./n.v.
2 Avenue du Japon - 1420 Braine l'Alleud
BELGIUM
Tel: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
Email: info.bn@meridianbioscience.eu

Meridian Bioscience Europe France
34 rue de Ponthieu - 75008 Paris
FRANCE
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 9 70 06 62 10
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
Postbus 301 - 5460 AH Veghel
NETHERLANDS
Tel: +31 (0) 411 62 11 66
Fax: +31 (0) 411 62 48 41
Email: info.bn@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simblos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
	In vitro diagnostic medical device / Dispositivo medico-diagnetico in vitro / Dispositivo medical de diagnostic in vitro / Dispositivo medico para diagnosticos in vitro / In-vitro-Diagnostikum	EC REP	Authorized representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandatario dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EG on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnetici in vitro / Ce produit est conforme aux exigences de la Directive 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / Sistema zur Probenverarbeitung, in dem sich Probenverdünnungspuffer befindet
	Catalogue number / Numero di catalogo / Référence du catalogue / Numero de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Eingefrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabbricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Solvente per valutazione delle prestazioni / Réactif IVD réservé à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Content sufficient for n tests / Contenuto sufficiente per n esami / Contenu suffisant pour n tests / Contenido suficiente para n ensayos / Inhalt ausreichend für n Prüflingen	SOLN STOP	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopflüssigkeit
	Temperature limitation / Limite di temperatura / Limites de température / Límite de temperatura / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Coniugato enzimatico / Conjugué enzymatique / Conjugado enzimático / Enzymkonjugat
	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / testgerät	REAG	Reagent / Reagents / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentration 20X / Soluzione di lavaggio concentrazione 20X / Solution de lavage concentrée 20X / Solución tampon de lavado 20X / 20fach konzentriertes Waschkonzentrat
	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig	DET REAG	Detection Reagent / Reagente Directo / Réactif de Détection / Reactivo de Detección / Nachweis Reagenz
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / ne pas utiliser si le paquet est endommagé / No use si el paquete esta dañado / Nicht verwenden, wenn die Verpackung beschädigt ist	TUBE	Empty Tube / Provaletta vuota / Tube vide / Tubo vacío / Leeres Gefäß

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.