# MERIFLUOR<sup>®</sup>

## CRYPTOSPORIDIUM/GIARDIA

Direct immunofluorescent detection procedure for the simultaneous detection of *Cryptosporidium* oocysts and *Giardia* cysts in fecal material

REF 250050

IVD In vitro diagnostic medical device

Rx Only

### INTENDED USE

**MERIFLUOR *Cryptosporidium/Giardia* (MERIFLUOR C/G)** is an in vitro direct immunofluorescent detection procedure for the simultaneous detection of *Cryptosporidium* oocysts and *Giardia* cysts in fecal material.

### SUMMARY AND EXPLANATION OF THE TEST

#### **Cryptosporidium:**

Protozoans of the genus *Cryptosporidium* are encountered worldwide and produce a self-limiting gastroenteritis in immunocompetent individuals. The primary symptoms include explosive, watery diarrhea, accompanied by vomiting, abdominal cramping, and low grade fever lasting two to fourteen days.<sup>1,2</sup> In immunocompromised patients, symptoms are more severe and persistent and may result in mortality particularly in patients with Acquired Immune Deficiency Syndrome.<sup>1,2</sup> Attempts to find effective chemotherapeutic agents have been unsuccessful thus far.

*Cryptosporidium* spp. have been found in a variety of domestic animals as well as poultry, calves, lambs, and goats. The first isolation of *Cryptosporidium* from a human was reported in 1976. Although several species have been described, little host specificity has been indicated.<sup>1,2</sup> The organism is a small, 2-6 micron, coccidian. The parasite undergoes its life cycle within the microvillous border of the intestine, although its occurrence in other sites has been reported.<sup>1,2</sup> The infective oocysts are excreted in feces from the host. Fecal-oral contamination is the most likely route of transmission.<sup>1,2</sup>

Early diagnosis of cryptosporidiosis using fecal material may eliminate the need for invasive procedures and facilitate supportive measures. Although various methods for the detection of oocysts in feces have been described, modified acid-fast staining is the most common procedure performed.<sup>1</sup> The acid-fast nature of the oocyst makes it difficult to differentiate from other noncryptosporidial acid-fast organisms of comparable size.

#### **Giardia:**

Giardiasis is one of the most common intestinal parasitic infections in the world.<sup>3</sup> Contaminated food, untreated surface-water polluted with cyst-containing animal feces in conjunction with inadequate filtration in water treatment facilities are the primary sources.<sup>3,4</sup> In addition, public health authorities include giardiasis as a sexually transmitted disease.<sup>3</sup>

Organisms of the genus *Giardia* were first observed in 1675 by Antony van Leeuwenhoek.<sup>3</sup> *Giardia* sp. are protozoan flagellates which exist in two forms: 1) a mobile and actively feeding trophozoite (9-21 µm long and 5-15 µm wide) found on the surface of the small intestine; 2) an infective oval-shaped cyst form which is smaller (8-12 µm long and 7-10 µm wide) and contains two to four nuclei.<sup>3,4</sup>

Clinical manifestations of *Giardia* infection range from asymptomatic cyst-passing with only mild gastrointestinal complaints to severely acute cases involving explosive, watery and foul-smelling diarrhea with abdominal cramps as well as bloating, flatulence, anorexia and nausea.<sup>3,4</sup> Though some cases present a self-limiting condition, symptoms of an acute or subacute infection may persist for several days to several months.<sup>3,4</sup> Therapeutic drug treatment using quinacrine hydrochloride, metronidazole, and furazolidone has been recommended for asymptomatic cyst-passers as well as symptomatic patients.<sup>3,4</sup>

Clinical diagnosis of giardiasis is routinely accomplished by direct microscopic examination of stool or stool concentrates for the presence of cysts and/or trophozoites. Conventional staining methods yield a sensitivity rate of 50-70%; lower rates have been observed when cysts were present at low levels.<sup>4,5</sup> Examination of duodenal and jejunal fluid aspirates and small intestine biopsy specimens have been used to improve conventional sensitivity rates but they are invasive and costly. Other methods include serodiagnosis by immunofluorescence and enzyme immunoassay procedures, however, neither technique indicates an active *Giardia* infection.<sup>4,5</sup>

An association between *Cryptosporidium* and *Giardia* infection has not been demonstrated. However numerous incidents of coinfection with both organisms have been documented.<sup>6</sup>

The **MERIFLUOR C/G** direct immunofluorescent detection procedure rapidly and accurately detects *Cryptosporidium* oocysts and *Giardia* cysts in clinical specimens. The specificity of monoclonal antibodies to *Cryptosporidium* and *Giardia* is combined with the sensitivity of the direct immunofluorescent assay format. The immunofluorescent technique is simple to perform and has been demonstrated to have a higher sensitivity than traditional staining procedures.<sup>7,8</sup>

### BIOLOGICAL PRINCIPLES

**MERIFLUOR C/G** utilizes the principle of direct immuno-fluorescence. The Detection Reagent contains a mixture of FITC labeled monoclonal antibodies directed against cell wall antigens of *Cryptosporidium* oocysts and *Giardia* cysts. A prepared fecal specimen is treated with the Detection Reagent and a Counterstain. The monoclonal antibodies attach to *Cryptosporidium* or *Giardia* antigens present in the specimen. The slides are rinsed to remove unbound antibodies, a coverslip is affixed with Mounting Medium, and the slides are examined for apple green color and characteristic morphology of *Cryptosporidium* oocysts and *Giardia* cysts using a fluorescent microscope. Background material in the specimen is counterstained dull orange to red.

### REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. **Detection Reagent:** FITC labeled anti-*Cryptosporidium* and anti-*Giardia* monoclonal antibodies in a buffered solution containing a protein stabilizer and 0.1% sodium azide.
2. **Counterstain:** Eriochrome Black solution.
3. **20X Wash Buffer:** Concentrated wash buffer with a preservative.
4. **Positive Control:** Formalinized stool preparation of *Cryptosporidium* oocysts and *Giardia* cysts containing 0.09% thimerosal.
5. **Negative Control:** Formalinized stool containing 0.09% thimerosal.
6. **Mounting Medium:** Buffered glycerol containing formalin, an anti-quincher and 0.05% sodium azide.
7. **Transfer Loops**
8. **Treated Slides**

\*NOTE: These reagents contain formaldehyde, which is a potential cancer hazard.

### MATERIALS NOT PROVIDED

1. Distilled or deionized water
2. Wash bottle
3. Humidity chamber
4. Microscope slide coverslips (22 x 50 mm), #1 thickness
5. Fluorescent microscope equipped with a filter system for fluorescein isothiocyanate (FITC) with the following parameters: Excitation wavelength: 490-500 nm, Barrier filter: 510-530 nm.
6. Applicator sticks

### PRECAUTIONS

1. All reagents and components are for in vitro diagnostic use only.
2. Do not mix reagents from different kit lot numbers.
3. Do not use kit components beyond the indicated expiration date on the kit label.
4. Protect the Detection Reagent and Counterstain from exposure to light.
5. All reagents should be mixed thoroughly before use.
6. Replace the reagent vial caps on their respective vials.
7. When distributing specimens or reagents onto the slide wells use care not to scratch the treated surface with the loop or applicator stick.

8. The Mounting Medium, Positive Control, and Negative Control contain formaldehyde, which is an irritant. Avoid contact with skin or clothing. If contact occurs, thoroughly flush area with water.
9. The Positive Control contains formalin killed *Cryptosporidium* oocysts and *Giardia* cysts, however, it should be handled as a potentially infectious specimen.
10. The positive and negative controls contain formaldehyde which is a potential cancer hazard.
11. Patient specimens may contain HIV or other infectious agents and should be handled by properly trained personnel and disposed of as potential biohazards. Wear disposable gloves while handling patient specimens and performing the test procedure.
12. Specimens preserved in PVA or MF/MIF are not suitable for use with **MERIFLUOR C/G**. Formalin, SAF and EcoFix® preserved specimens may be used.
13. If stool material is not seen upon scanning the slide wells, loss of sample may have occurred. This is usually due to an overly vigorous wash procedure or insufficient specimen drying time.

**WARNING:** Some reagents in this kit contain sodium azide which is a skin irritant and harmful if swallowed. Avoid skin contact with reagents. Disposal of reagents containing sodium azide into lead or copper plumbing can result in the formation of explosive metal azides. This can be avoided by flushing with a large volume of water during such disposal.

### HAZARD AND PRECAUTIONARY STATEMENTS

Refer to the SDS, available at [www.meridianbioscience.com](http://www.meridianbioscience.com) for Hazards and Precautionary Statements.

### SHELF LIFE AND STORAGE

The expiration date of the **MERIFLUOR C/G** kit is indicated on the kit label. To maintain optimal performance, the kit should be stored at 2-8 C and returned promptly to the refrigerator after each use. **CAUTION: DO NOT FREEZE.**

### REAGENT PREPARATION

1. Bring the entire kit to room temperature before use.
2. Prepare 1X wash buffer as needed. For example: 2.5 mL of 20X wash buffer + 47.5 mL of distilled or deionized water. Unused 1X wash buffer can be stored at room temperature up to three months.

### SPECIMEN COLLECTION AND PREPARATION

1. Stool specimens preserved in 10% formalin, sodium acetate-acetic acid-formalin (SAF) or EcoFix should be used. The preservative-stool mixture must be allowed to stand for a minimum of 30 minutes at RT to ensure adequate fixation. **MERIFLUOR C/G** is designed for use with the Para-Pak® line of parasitology products. Collect the specimen as directed by the appropriate Para-Pak package insert supplied with the collection system. Complimentary Para-Pak Formalin or SAF package inserts can be supplied by Meridian Bioscience, Inc., Technical Support Center. Please call 1-800-343-3858.
2. To increase the probability of detection in stools with low numbers of cysts or oocysts, the specimen may be concentrated prior to processing (see **SPECIFIC PERFORMANCE CHARACTERISTICS**). The Para-Pak® Macro-CON® and Para-Pak® CON-Trate® Systems are recommended for this purpose. The removal of fecal lipids by ethyl acetate may be particularly beneficial if the fecal specimen is mucoid. Dilipidation will facilitate the adherence of fecal material to the treated microscope slide.

### TEST PROCEDURE

1. Use a transfer loop to transfer a drop of fecal material to a treated slide well. Spread the specimen over the entire well. Do not scratch the treated surface of the slide.
2. Use a new transfer loop to transfer a drop of Positive Control to a separate treated slide well. Spread the Positive Control over the entire well. Do not scratch the treated surface of the slide.
3. Use a new transfer loop to transfer a drop of Negative Control to a separate treated slide well.
4. Allow the slides to air dry completely at room temperature (usually requires 30 minutes).
5. Place one drop of Detection Reagent in each well.
6. Place one drop of Counterstain in each well.
7. Mix the reagents with an applicator stick and spread over the entire well. Do not scratch the treated surface of the slide.
8. Incubate the slides in a humidified chamber for 30 minutes at room temperature. **Note: Protect from light.**
9. Use a wash bottle to rinse the slides with a gentle stream of 1X Wash Buffer until excess Detection Reagent and Counterstain is removed. **Note: Do not submerge the slides during rinsing.** Avoid disturbing the specimen or causing cross contamination of the specimens.
10. Remove excess buffer by tapping the long edge of the slide on a clean paper towel. **Note: Do not allow slide to dry.**
11. Add one drop of Mounting Medium to each well and apply a coverslip.
12. Scan each well thoroughly using 100-200X magnification. The presence of *Cryptosporidium* oocysts should be confirmed at higher magnification.

### INTERPRETATION OF RESULTS

#### **Control Reactions:**

1. *Cryptosporidium* oocysts are round to slightly oval in shape, 2-6 µm in diameter. The oocyst wall will stain bright apple green. A suture line may also be visible.
2. *Giardia* cysts are oval-shaped organisms, 8-12 µm long. The cyst wall will stain bright apple green.
3. Background material should counterstain dull orange to red.
4. The negative control well should not exhibit apple green fluorescence of characteristic cyst or oocyst morphology.
5. Some background coloration may sometimes result as an artifact of slide preparation. This background fluorescence is easily discernable by its lack of vibrant, apple green color or characteristic morphology of the cyst.

#### **Patient Specimens:**

1. ***Cryptosporidium* positive test result:** Any stool specimen exhibiting one or more oocysts with an apple green color and characteristic morphology should be considered positive for the presence of *Cryptosporidium* sp.
2. ***Giardia* positive test result:** Any stool specimen exhibiting one or more cysts with an apple green color and characteristic morphology should be considered positive for the presence of *Giardia*.
3. **Negative test result:** A stool specimen with no apple green fluorescence should be considered negative for the presence of *Cryptosporidium* oocysts and *Giardia* cysts. Background fluorescent debris which does not exhibit vibrant apple green color or characteristic morphology should also be considered negative.

### QUALITY CONTROL

**This test should be performed per applicable local, state, or federal regulations or accrediting agencies.**

To provide quality control of the reagents contained in the kit, the Positive and Negative Controls should be evaluated each time a patient specimen is tested. **If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact the Meridian Technical Services Department at 1-800-343-3858 (US) or your local distributor.**

At the time of each use, the kit components should be visually examined for obvious signs of contamination, freezing, or leakage.

It is suggested that the results of each quality control check be recorded in an appropriate log book to maintain high quality testing procedures and comply with regulatory agencies.

### EXPECTED VALUES

The **MERIFLUOR C/G** direct immunofluorescent test will identify the presence of *Cryptosporidium* oocysts and *Giardia* cysts in stool specimens. Typically, patients in either active disease state shed oocysts or cysts, respectively, in moderate to large numbers. Some patients may continue to excrete oocysts or cysts in small numbers after complete clinical recovery.<sup>2,4</sup> Variations in the rate of positivity may be expected due to age, location, and general health conditions of the patient population under study.<sup>1,3</sup>

### LIMITATIONS OF THE PROCEDURE

**MERIFLUOR C/G** is an in vitro diagnostic procedure to aid in the diagnosis of *Cryptosporidium* and *Giardia* infection. While the presence of *Cryptosporidium* or *Giardia* is often associated with diarrhea and vomiting, shedding of oocysts or cysts by asymptomatic individuals has been observed.<sup>1,4</sup> In addition, the presence of oocysts or cysts in a stool specimen does not preclude the existence of other microorganisms or another underlying condition as the causative agent of a patient's symptoms. For these reasons appropriate concurrent testing for other etiologic agents should be considered.

## SPECIFIC PERFORMANCE CHARACTERISTICS

Stool specimens preserved in 10% Formalin were received from several reference labs and tested in a blind study. Specimens were concentrated by the Para-Pak® CON-Trate® System and stained with **MERIFLUOR C/G (DFA), MERIFLUOR Cryptosporidium** (IFA), **MERIFLUOR Giardia** (IFA), hot modified acid-fast, and iodine.

### 1. Cryptosporidium

The following tables compare the results of the various staining procedures for *Cryptosporidium* oocysts in 115 stools. For specimens containing numerous oocysts, positive results were obtained by the three relevant staining procedures. For the 15 specimens with discrepant results, low numbers of oocysts were observed by the positive staining method. To resolve the discrepant results, these specimens were stained twice more by each method.

DFA	Acid Fast		Resolved	
	+	-	+	-
+	24	13	32	5
-	2	76	1	77

Sensitivity	92%	97%
Specificity	85%	94%
Correlation	87%	95%

### 2. Giardia

The following tables compare the results of the various staining procedures for *Giardia* cysts in 115 concentrated stool specimens. For the six specimens with discrepant results, low numbers of cysts were observed by the positive staining method. To resolve the discrepant results, these specimens were stained by the trichrome staining method using a paired PVA preserved specimen.

DFA	Iodine		Resolved	
	+	-	+	-
+	38	4	42	0
-	2	71	0	73

Sensitivity	95%	100%
Specificity	95%	100%
Correlation	95%	100%

### 3. Nonconcentrated Stools

Specimens were stained with **MERIFLUOR C/G** direct from the formalin vial or after concentration. Two of the 32 specimens positive for *Cryptosporidium* in concentrated stool were negative in the nonconcentrated sample. In both instances few oocysts were observed in the concentrated sample. There was total agreement with 40 *Giardia* positive specimens.

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## REFERENCES

- Fayer R and Ungar BLP. *Cryptosporidium* spp. and Cryptosporidiosis. Microbiol. Rev. 1986; 50:458-483.
- Janoff EN and Reller LB. *Cryptosporidium* species, a protein protozoan. J Clin Microbiol. 1987;25:967-975.
- Stevens David P. Selective primary health care: strategies for control of disease in the developing world. XIX. Giardiasis. Rev. Inf. Dis. 1985;7:530-535.
- Dupont L and Sullivan PS. Giardiasis: the clinical spectrum, diagnosis and therapy. Ped Inf Dis. 1986;6:S131-138.
- Ungar BLP, Yolken RH, Nash TE and Quinn TC. Enzyme-linked immunosorbent assay for the detection of *giardia lamblia* in fecal specimens. J Inf Dis. 1984; 149:90-97.
- Jokipii L, Pohjola S and Jokipii AMM. *Cryptosporidiosis* associated with traveling and giardiasis. Gastroenterol. 1985;89:838-842.
- Garcia LS, Brewer TC and Bruckner DA. Fluorescence detection of *cryptosporidium* oocysts in human fecal specimens by using monoclonal antibodies. J Clin Microbiol. 1987;25:119-121.
- Data on file, Meridian Bioscience, Inc. Cincinnati, OH 45244.

SN11220

REV. 04/16

Manufactured By

Meridian Bioscience, Inc.  
USA/Corporate Office  
3471 River Hills Drive  
Cincinnati, Ohio 45244  
Telephone: 513.271.3700  
Orders/Customer Service:  
800.543.1980  
Technical Support Center:  
800.343.3858  
Information Fax: 513.272.5432  
Ordering Fax: 513.271.0124

---

Meridian Bioscience Europe S. r. L.  
Via dell'Industria, 7  
20020 Villa Cortese, Milano  
ITALY  
Tel: +39 0331 43 36 36  
Fax: +39 0331 43 36 16  
Email: info@meridianbioscience.eu  
WEB: www.meridianbioscience.eu

Authorized Representative

Meridian Bioscience Europe s.a./n.v.  
2 Avenue du Japon - 1420 Braine l'Alleud  
BELGIUM  
Tel: +32 (0) 67 89 59 59  
Fax: +32 (0) 67 89 59 58  
Email: info.bn@meridianbioscience.eu

Meridian Bioscience Europe France  
34 rue de Ponthieu - 75008 Paris  
FRANCE  
Tel: +33 (0) 1 42 56 04 40  
Fax: +33 (0) 9 70 06 62 10  
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.  
Postbus 301 - 5460 AH Veghel  
NETHERLANDS  
Tel: +31 (0) 411 62 11 66  
Fax: +31 (0) 411 62 48 41  
Email: info.bn@meridianbioscience.eu

## INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	<b>CONTROL +</b>	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
<b>LOT</b>	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	<b>CONTROL -</b>	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
<b>IVD</b>	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	<b>EC REP</b>	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
<b>CE</b>	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostico in vitro / Ce produit répond aux exigences de la Directive 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostika 98/79/EG.	<b>SMP PREP DIL SPE</b>	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet
<b>REF</b>	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	<b>BUF RXN</b>	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	<b>SOLN STOP</b>	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	<b>CONJ ENZ</b>	Enzyme Conjugate / Coniugato enzimatico / Conjugé enzymatique / Conjugado enzimático / Enzymkonjugat
<b>SN</b>	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	<b>CONTROL</b>	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
<b>TEST</b>	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	<b>REAG</b>	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsverstadium	<b>BUF WASH</b>	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Washpuffer
<b>BUF</b>	Buffer / Soluzione tamponne / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
<b>CONJ</b>	Conjugate / Coniugato / Conjugat / Conjugado / Konjugat	<b>DIL SPE</b>	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillon / Diluyente de muestra / Probenverdünnungspuffer
<b>SUBS</b>	Substrate / Substrato / Substrat / Substrato / Substrat	<b>BUF WASH 20X</b>	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
<b>R. Only</b>	Caution: US Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner. / La legge Federale negli Stati Uniti prevede che questo test sia venduto solo a un operatore abilitato ad operare nella salute. / Attention: la vente ou la prescription de ce dispositif par un médecin agréé est soumise à des restrictions par la loi fédérale des US. / Cuidado: las leyes federales de los EU limitan la venta de dispositivos por o por la orden de licenciados en la práctica del curso de la salud. / Vorsicht: Der Verkauf oder die Verschreibung dieses Geräts durch einen zugelassenen Arzt unterliegt den Beschränkungen des US-Bundesgesetzes.	<b>DET REAG</b>	Detection Reagent / Reagente Diretto / Réactif de Détection / Reactivo de Detección / Nachweis Reagenz

For technical assistance, call Technical Support Services at (800) 343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at (800) 543-1980.