Cryptosporidium / Giardia

Direct immunofluorescent detection procedure for the simultaneous detection
of *Cryptosporidium* oocysts and *Giardia* cysts in fecal material

REF 250050

IVD

Rx Only

INTENDED USE

MERIFLUOR *Cryptosporidium/Giardia* (MERIFLUOR C/G) is an in vitro direct immunofluorescent detection procedure for the simultaneous detection of *Cryptosporidium* oocysts and *Giardia* cysts in fecal material.

SUMMARY AND EXPLANATION OF THE TEST

Cryptosporidium:

Protozoans of the genus *Cryptosporidium* are encountered worldwide and produce a self-limiting gastroenteritis in immunocompetent individuals. The primary symptoms include explosive, watery diarrhea, accompanied by vomiting, abdominal cramping, and low grade fever lasting two to fourteen days.^{1,2} In immunocompromised patients, symptoms are more severe and persistent and may result in mortality particularly in patients with Acquired Immune Deficiency Syndrome.^{1,2} Attempts to find effective chemotherapeutic agents have been unsuccessful thus far.

Cryptosporidium spp. have been found in a variety of domestic animals as well as poultry, calves, lambs, and goats. The first isolation of *Cryptosporidium* from a human was reported in 1976. Although several species have been described, little host specificity has been indicated.^{1,2} The organism is a small, 2-6 micron, coccidian. The parasite undergoes its life cycle within the microvillous border of the intestine, although its occurrence in other sites has been reported.^{1,2} The infective oocysts are excreted in feces from the host. Fecal-oral contamination is the most likely route of transmission.^{1,2}

Early diagnosis of cryptosporidiosis using fecal material may eliminate the need for invasive procedures and facilitate supportive measures. Although various methods for the detection of oocysts in feces have been described, modified acid-fast staining is the most common procedure performed.¹ The acid-fast nature of the oocyst makes it difficult to differentiate from other noncryptosporidial acid-fast organisms of comparable size.

Giardia:

Giardiasis is one of the most common intestinal parasitic infections in the world.³ Contaminated food, untreated surface-water polluted with cyst-containing animal feces in conjunction with inadequate filtration in water treatment facilities are the primary sources.^{3,4} In addition, public health authorities include giardiasis as a sexually transmitted disease.³

Organisms of the genus *Giardia* were first observed in 1675 by Antony van Leeuwenhoek.³ *Giardia* sp. are protozoan flagellates which exist in two forms: 1) a mobile and actively feeding trophozoite (9-21 µm long and 5-15 µm wide) found on the surface of the small intestine; 2) an infective oval-shaped cyst form which is smaller (8-12 µm long and 7-10 µm wide) and contains two to four nuclei.^{3,4}

Clinical manifestations of *Giardia* infection range from asymptomatic cyst-passing with only mild gastrointestinal complaints to severely acute cases involving explosive, watery and foul-smelling diarrhea with abdominal cramps as well as bloating, flatulence, anorexia and nausea.^{3,4} Though some cases present a self-limiting condition, symptoms of an acute or subacute infection may persist for several days to several months.^{3,4} Therapeutic drug treatment using quinacrine hydrochloride, metronidazole, and furazolidone has been recommended for asymptomatic cyst-passers as well as symptomatic patients.^{3,4}

Clinical diagnosis of giardiasis is routinely accomplished by direct microscopic examination of stool or stool concentrates for the presence of cysts and/or trophozoites. Conventional staining methods yield a sensitivity rate of 50-70%: lower rates have been observed when cysts were present at low levels.^{4,5} Examination of duodenal and jejunal fluid aspirates and small intestine biopsy specimens have been used to improve conventional sensitivity rates but they are invasive and costly. Other methods include serodiagnosis by immunofluorescence and enzyme immunoassay procedures, however, neither technique indicates an active *Giardia* infection.^{4,5}

An association between *Cryptosporidium* and *Giardia* infection has not been demonstrated. However numerous incidents of coinfection with both organisms have been documented.⁶

The MERIFLUOR C/G direct immunofluorescent detection procedure rapidly and accurately detects *Cryptosporidium* oocysts and *Giardia* cysts in clinical specimens. The specificity of monoclonal antibodies to *Cryptosporidium* and *Giardia* is combined with the sensitivity of the direct immunofluorescent assay format. The immunofluorescent technique is simple to perform and has been demonstrated to have a higher sensitivity than traditional staining procedures.^{7,8}

BIOLOGICAL PRINCIPLES

MERIFLUOR C/G utilizes the principle of direct immuno-fluorescence. The Detection Reagent contains a mixture of FITC labeled monoclonal antibodies directed against cell wall antigens of *Cryptosporidium* oocysts and *Giardia* cysts. A prepared fecal specimen is treated with the Detection Reagent and a Counterstain. The monoclonal antibodies attach to *Cryptosporidium* or *Giardia* antigens present in the specimen. The slides are rinsed to remove unbound antibodies, a coverslip is affixed with Mounting Medium, and the slides are examined for apple green color and characteristic morphology of *Cryptosporidium* oocysts and *Giardia* cysts using a fluorescent microscope. Background material in the specimen is counterstained dull orange to red.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. Detection Reagent: FITC labeled anti-*Cryptosporidium* and anti-*Giardia* monoclonal antibodies in a buffered solution containing a protein stabilizer and 0.1% sodium azide.
2. Counterstain: Eriochrome Black solution.
3. 20X Wash Buffer: Concentrated wash buffer with a preservative.
4. Positive Control*: Formalinized stool preparation of *Cryptosporidium* oocysts and *Giardia* cysts containing 0.09% thimerosal.
5. Negative Control*: Formalinized stool containing 0.09% thimerosal.
6. Mounting Medium: Buffered glycerol containing formalin, an antifluorescent and 0.05% sodium azide.
7. Transfer Loops
8. Treated Slides

*NOTE: These reagents contain formaldehyde, which is a potential cancer hazard.

MATERIALS NOT PROVIDED

1. Distilled or deionized water
2. Wash bottle
3. Humidity chamber
4. Microscope slide coverslips (22 x 50 mm), #1 thickness
5. Fluorescent microscope equipped with a filter system for fluorescein isothiocyanate (FITC) with the following parameters: Excitation wavelength: 490-500 nm, Barrier filter: 510-530 nm.
6. Applicator sticks

PRECAUTIONS

1. All reagents and components are for in vitro diagnostic use only.
2. Do not mix reagents from different kit lot numbers.
3. Do not use kit components beyond the indicated expiration date on the kit label.
4. Protect the Detection Reagent and Counterstain from exposure to light.
5. All reagents should be mixed thoroughly before use.
6. Replace the reagent vial caps on their respective vials.
7. When distributing specimens or reagents onto the slide wells use care not to scratch the treated surface with the loop or applicator stick.
8. The Mounting Medium, Positive Control, and Negative Control contain formaldehyde, which is an irritant. Avoid contact with skin or clothing. If contact occurs, thoroughly flush area with water.
9. The Positive Control contains formalin killed *Cryptosporidium* oocysts and *Giardia* cysts, however, it should be handled as a potentially infectious specimen.
10. The positive and negative controls contain formaldehyde which is a potential cancer hazard.
11. Patient specimens may contain HIV or other infectious agents and should be handled by properly trained personnel and disposed of as potential biohazards. Wear disposable gloves while handling patient specimens and performing the test procedure.
12. Specimens preserved in PVA or MF/MIF are not suitable for use with MERIFLUOR C/G. Formalin, SAF and EcoFix® preserved specimens may be used.
13. If stool material is not seen upon scanning the slide wells, loss of sample may have occurred. This is usually due to an overly vigorous wash procedure or insufficient specimen drying time.

WARNING: Some reagents in this kit contain sodium azide which is a skin irritant and harmful if swallowed. Avoid skin contact with reagents. Disposal of reagents containing sodium azide into lead or copper plumbing can result in the formation of explosive metal azides. This can be avoided by flushing with a large volume of water during such disposal.

HAZARD and PRECAUTIONARY STATEMENTS

 Detection Antibody	Signal Word Danger Hazard Statements H301 – Toxic if swallowed Precautionary Statements - EU (§28, 1272/2008) P301 + P310 – IF SWALLOWED: Immediately call a POISON CENTER or doctor/ physician. P321 - Specific treatment (see supplemental first aid instructions with this material)
  Counterstain	Signal Word Danger Hazard Statements H302 - Harmful if swallowed H319 - Causes serious eye irritation H402 - Harmful to aquatic life H411 - Toxic to aquatic life with long lasting effects H360D - May damage the unborn child H226 - Flammable liquid and vapor Precautionary Statements - EU (§28, 1272/2008) P201 - Obtain special instructions before use P281 - Use personal protective equipment as required P308 + P313 - If exposed or concerned: Get medical advice/ attention P202 - Do not handle until all safety precautions have been read and understood P370 + P378 - In case of fire: Use dry sand, dry chemical or alcohol-resistant foam for extinction P210 - Keep away from heat/sparks/open flames/hot surfaces. - No smoking
  Negative Control Positive Control	Signal Word Danger Hazard Statements H301 - Harmful if swallowed H311 – Toxic in contact with skin H317 - May cause an allergic skin reaction H330 - Fatal if inhaled H341 – Suspected of causing genetic defects H350 – May cause cancer Contains Formaldehyde, Methyl alcohol Precautionary Statements - EU (§28, 1272/2008) P301 + P310 – IF SWALLOWED: Immediately call a POISON CENTER or doctor/ physician P280 – Wear protective gloves/ protective clothing/ eye protection/ face protection P310 – Immediately call a POISON CENTER or doctor/ physician P304 + P340 – IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing P260 – Do not breathe dust/fume/gas/mist/vapors/spray P403 + P233 – Store in a well-ventilated place. Keep container tightly closed P280 – Wear eye protection/ face protection P321 – See SDS Section 4 or Section 11 for specific medical treatment information P201 – Obtain special instructions before use P281 – Use personal protective equipment as required P308 + P313 - If exposed or concerned: Get medical advice/ attention P320 – See SDS Section 4 for any urgent treatment information
 Wash Buffer	Signal Word Danger Hazard Statements H302 – Harmful if swallowed H311 – Toxic in contact with skin Precautionary Statements - EU (§28, 1272/2008) P280 – Wear protective gloves/ protective clothing/ eye protection/ face protection
 Mounting Medium	Signal Word Danger Hazard Statements H317 – May cause an allergic skin reaction H341 – Suspected of causing genetic defects H350 – May cause cancer Contains Formaldehyde Precautionary Statements - EU (§28, 1272/2008) P280 – Wear protective gloves/ protective clothing/ eye protection/ face protection P201 – Obtain special instructions before use P281 – Use personal protective equipment as required P308 + P313 - If exposed or concerned: Get medical advice/ attention P321 – See SDS Section 4 or Section 11 for additional medical treatment information

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-8 C and return the kit promptly to the refrigerator after each use. **CAUTION: DO NOT FREEZE.**

REAGENT PREPARATION

1. Bring the entire kit to room temperature before use.
2. Prepare 1X wash buffer as needed. For example: 2.5 mL of 20X wash buffer + 47.5 mL of distilled or deionized water. Unused 1X wash buffer can be stored at room temperature up to three months.

SPECIMEN COLLECTION AND PREPARATION

1. Stool specimens preserved in 10% formalin, sodium acetate-acetic acid-formalin (SAF) or EcoFix should be used. The preservative-stool mixture must be allowed to stand for a minimum of 30 minutes at RT to ensure adequate fixation. **MERIFLUOR C/G** is designed for use with the Para-Pak® line of parasitology products. Collect the specimen as directed by the appropriate Para-Pak package insert supplied with the collection system. Complimentary Para-Pak Formalin or SAF package inserts can be supplied by Meridian Bioscience, Inc., Technical Support Center. Please call 1-800-343-3858.
2. To increase the probability of detection in stools with low numbers of cysts or oocysts, the specimen may be concentrated prior to processing (see **SPECIFIC PERFORMANCE CHARACTERISTICS**). The Para-Pak® Macro-CON-G and Para-Pak® CON-Trate® Systems are recommended for this purpose. The removal of fecal lipids by ethyl acetate may be particularly beneficial if the fecal specimen is mucoid. Dilipidation will facilitate the adherence of fecal material to the treated microscope slide.

TEST PROCEDURE

1. Use a transfer loop to transfer a drop of fecal material to a treated slide well. Spread the specimen over the entire well. Do not scratch the treated surface of the slide.
2. Use a new transfer loop to transfer a drop of Positive Control to a separate treated slide well. Spread the Positive Control over the entire well. Do not scratch the treated surface of the slide.
3. Use a new transfer loop to transfer a drop of Negative Control to a separate treated slide well.
4. Allow the slides to air dry completely at room temperature (usually requires 30 minutes).
5. Place one drop of Detection Reagent in each well.
6. Place one drop of Counterstain in each well.
7. Mix the reagents with an applicator stick and spread over the entire well. Do not scratch the treated surface of the slide.
8. Incubate the slides in a humidified chamber for 30 minutes at room temperature. **Note: Protect from light.**
9. Use a wash bottle to rinse the slides with a gentle stream of 1X Wash Buffer until excess Detection Reagent and Counterstain is removed. **Note: Do not submerge the slides during rinsing.** Avoid disturbing the specimen or causing cross contamination of the specimens.
10. Remove excess buffer by tapping the long edge of the slide on a clean paper towel. **Note: Do not allow slide to dry.**
11. Add one drop of Mounting Medium to each well and apply a coverslip.
12. Scan each well thoroughly using 100-200X magnification. The presence of *Cryptosporidium* oocysts should be confirmed at higher magnification.

INTERPRETATION OF RESULTS

Control Reactions:

1. *Cryptosporidium* oocysts are round to slightly oval in shape, 2-6 µm in diameter. The oocyst wall will stain bright apple green. A suture line may also be visible.
2. *Giardia* cysts are oval-shaped organisms, 8-12 µm long. The cyst wall will stain bright apple green.
3. Background material should counterstain dull orange to red.
4. The negative control well should not exhibit apple green fluorescence of characteristic cyst or oocyst morphology.
5. Some background coloration may sometimes result as an artifact of slide preparation. This background fluorescence is easily discernable by its lack of vibrant, apple green color or characteristic morphology of the cyst.

Patient Specimens:

1. ***Cryptosporidium* positive test result:** Any stool specimen exhibiting one or more oocysts with an apple green color and characteristic morphology should be considered positive for the presence of *Cryptosporidium* sp.
2. ***Giardia* positive test result:** Any stool specimen exhibiting one or more cysts with an apple green color and characteristic morphology should be considered positive for the presence of *Giardia*.
3. **Negative test result:** A stool specimen with no apple green fluorescence should be considered negative for the presence of *Cryptosporidium* oocysts and *Giardia* cysts. Background fluorescent debris which does not exhibit vibrant apple green color or characteristic morphology should also be considered negative.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

To provide quality control of the reagents contained in the kit, the Positive and Negative Controls should be evaluated each time a patient specimen is tested. If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact the Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

At the time of each use, the kit components should be visually examined for obvious signs of contamination, freezing, or leakage.

It is suggested that the results of each quality control check be recorded in an appropriate log book to maintain high quality testing procedures and comply with regulatory agencies.

EXPECTED VALUES

The MERIFLUOR C/G direct immunofluorescent test will identify the presence of *Cryptosporidium* oocysts and *Giardia* cysts in stool specimens. Typically, patients in either active disease state shed oocysts or cysts, respectively, in moderate to large numbers. Some patients may continue to excrete oocysts or cysts in small numbers after complete clinical recovery.^{2,4} Variations in the rate of positivity may be expected due to age, location, and general health conditions of the patient population under study.^{1,3}

LIMITATIONS OF THE PROCEDURE

MERIFLUOR C/G is an in vitro diagnostic procedure to aid in the diagnosis of *Cryptosporidium* and *Giardia* infection. While the presence of *Cryptosporidium* or *Giardia* is often associated with diarrhea and vomiting, shedding of oocysts or cysts by asymptomatic individuals has been observed.^{1,4} In addition, the presence of oocysts or cysts in a stool specimen does not preclude the existence of other microorganisms or another underlying condition as the causative agent of a patient's symptoms. For these reasons appropriate concurrent testing for other etiologic agents should be considered.

SPECIFIC PERFORMANCE CHARACTERISTICS

Stool specimens preserved in 10% Formalin were received from several reference labs and tested in a blind study. Specimens were concentrated by the Para-Pak® CON-Trate® System and stained with MERIFLUOR C/G (DFA), MERIFLUOR Cryptosporidium (IFA), MERIFLUOR Giardia (IFA), hot modified acid-fast, and iodine.

1. *Cryptosporidium*

The following tables compare the results of the various staining procedures for *Cryptosporidium* oocysts in 115 stools. For specimens containing numerous oocysts, positive results were obtained by the three relevant staining procedures. For the 15 specimens with discrepant results, low numbers of oocysts were observed by the positive staining method. To resolve the discrepant results, these specimens were stained twice more by each method.

DFA	Acid Fast				Resolved			
	+	-	+	-	+	-	+	-
+	24	13	32	5				
-	2	76	1	77				

Sensitivity 92%

97%

Specificity 85%

94%

Correlation 87%

95%

2. *Giardia*

The following tables compare the results of the various staining procedures for *Giardia* cysts in 115 concentrated stool specimens. For the six specimens with discrepant results, low numbers of cysts were observed by the positive staining method. To resolve the discrepant results, these specimens were stained by the trichrome staining method using a paired PVA preserved specimen.

DFA	Iodine				Resolved			
	+	-	+	-	+	-	+	-
+	38	4	42	0				
-	2	71	0	73				

Sensitivity 95%

100%

Specificity 95%

100%

Correlation 95%

100%

3. Nonconcentrated Stools

Specimens were stained with MERIFLUOR C/G direct from the formalin vial or after concentration. Two of the 32 specimens positive for *Cryptosporidium* in concentrated stool were negative in the nonconcentrated sample. There was total agreement with 40 *Giardia* positive specimens.

CROSSREACTIVITY

The following organisms did not fluoresce when tested with MERIFLUOR C/G:³

Blastocystis hominis, *Candida albicans*, *Chilomastix mesnili*, *Dientamoeba fragilis*, *Endolimax nana*, *Entamoeba coli*, *Entamoeba hartmanni*, *Entamoeba histolytica*, *Iodamoeba butschlii*, *Strongyloides stercoralis*, *Trichomonas hominis*, *Trichomonas vaginalis*

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REFERENCES

1. Fayer R and Ungar BLP. *Cryptosporidium* spp. and *Cryptosporidiosis*. *Microbiol. Rev.* 1986; 50:458-483.
2. Janoff EN and Reller LB. *Cryptosporidium* species, a protein protozoan. *J Clin Microbiol.* 1987;25:967-975.
3. Stevens David P. Selective primary health care: strategies for control of disease in the developing world. XIX. Giardiasis. *Rev. Inf. Dis.* 1985;7:530-535.
4. Dupont L and Sullivan PS. Giardiasis: the clinical spectrum, diagnosis and therapy. *Ped Inf Dis.* 1986;86:S131-138.
5. Ungar BLP, Yolken RH, Nash TE and Quinn TC. Enzyme-linked immunosorbent assay for the detection of *giardia lamblia* in fecal specimens. *J Inf Dis.* 1984; 149:90-97.
6. Jokipii L, Pohjola S and Jokipii AMM. *Cryptosporidiosis* associated with traveling and giardiasis. *Gastroenterol.* 1985;89:838-842.
7. Garcia LS, Brewer TC and Bruckner DA. Fluorescence detection of *cryptosporidium* oocysts in human fecal specimens by using monoclonal antibodies. *J Clin Microbiol.* 1987;25:119-121.
8. Data on file, Meridian Bioscience, Inc. Cincinnati, OH 45244.

SN11220

REV. 10/20

 Manufactured By	Meridian Bioscience, Inc. Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 USA Telephone: 513.271.3700 Orders/Customer Service: 800.543.1980 Technical Support Center: 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124	Meridian Bioscience Europe s.a./n.v. 2 Avenue du Japon - 1420 Braine l'Alleud BELGIUM Tel: +32 (0) 67 89 59 59 Fax: +32 (0) 67 89 59 58 Email: info.bn@meridianbioscience.eu
 Authorized Representative	Meridian Bioscience Europe S. r. L. Via dell'Industria, 7 20035 Villa Cortese, Milano ITALY Tel: +39 0331 43 36 36 Fax: +39 0331 43 36 16 Email: info@meridianbioscience.eu WEB: www.meridianbioscience.com/eu	Meridian Bioscience Europe France 34 rue de Ponthieu - 75008 Paris FRANCE Tel: +33 (0) 1 42 56 04 40 Fax: +33 (0) 9 70 06 62 10 Email: info.fr@meridianbioscience.eu
Australian Sponsor Emergo Australia Level 20, Tower II Darling Park 201 Sussex Street Sydney, NSW 2000 Australia		Meridian Bioscience Europe B.V. Zekeringstraat 17 A 1014BM Amsterdam NETHERLANDS Tel: +31 (0) 411 62 11 66 Fax: +31 (0) 411 62 48 41 Email: Info.bn@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
	In vitro diagnostic medical device / Dispositivo medico-diagnóstico in vitro / Dispositivo médico de diagnóstico in vitro / Dispositivo medico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EG on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Proben vorbereitung, in dem sich Proben verdünnungs puffer befindet
	Catalogue number / Número de catálogo / Référence du catalogue / Numero de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Eingrenzen
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampon de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Solamente para evaluación de las prestaciones / Réactifs IVD réservés à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" test / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	SOLN STOP	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limite di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Coniugato enzimatico / Conjugue enzymatique / Conjunto enzimático / enzymkeimjag
	Serial number / Número de serie / Numéro de série / Número de serie / Seriennummer	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / testgerät	REAG	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
	Conjugate / Coniugato / Conjugué / Conjunto / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillon / Diluyente de muestra / Probenverdünnungspuffer
	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentration 20X / Soluzione diluente 20X / Solution de lavage concentrée 20X / Solución tamponnée de lavado 20X / 20fach konzentriertes Waschkonzentrat
	Prescription Use Only / Per l'uso su prescrizione medica / Uniquelement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig	DET REAG	Detection Reagent / Reagente Directo / Réactif de Détection / Reactivo de Detección / Nachweis Reagenz
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / ne pas utiliser si le paquet est endommagé / No use si el paquete está dañado / Nicht verwenden, wenn die Verpackung beschädigt ist	TUBE	Empty Tube / Provaletta vuota / Tube vide / Tube vacío / Leeres Gefäß

For technical assistance, call Technical Support Services at (800) 343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at (800) 543-1980.