Para-Pak® Systems
for use with LV-PVA and 10% Formalin
(Low Viscosity PVA Fixative)
(Patent #5624554)

REF 300112, 330284, 900112

IVD

Rx Only

INTENDED USE

Para-Pak based systems provide standardized procedures for the routine collection, transportation, preservation, and examination of stool specimens for intestinal parasites. Kit systems are designed for easy use by individuals not trained in microbiological procedures and afford an excellent means of minimizing the adverse effects of delay in specimen transportation. Bulk products, for preservation of fresh specimens, are also available.

SUMMARY AND EXPLANATION OF THE TEST

Diagnosis of intestinal parasitic disease is confirmed by recovery and identification of helminth eggs and larvae, or protozoan trophozoites and cysts in the clinical parasitology laboratory. Timely collection and transportation of "fresh" stool specimens to the laboratory cannot always be insured. Workload conditions and priorities in clinical laboratories frequently do not permit immediate examination of "fresh" specimens. Procedures such as incubation, refrigeration or freezing of stool specimens will not guarantee the recovery of all diagnostic stages of all parasites.^{5, 6, 9, 10, 14, 15}

Proper use of the Para-Pak systems thus assures the parasitologist that diagnostic stages of intestinal parasites will be preserved, if present, in the fecal material.

BIOLOGICAL PRINCIPLES

Polyvinyl alcohol is a water soluble plastic which when combined with Schaudin's fixative provides an excellent preservative-fixative for protozoan trophozoites. The resulting permanent slide lends itself to commonly used staining procedures such as iron-hematoxylin or Wheatley's trichrome.^{2, 4, 5, 16}

Para-Pak low viscosity polyvinyl alcohol, fixative has certain advantages over previous formulations. The low viscosity mixture allows easier miscibility of the specimen, better fixation due to faster penetration, faster drying and superior film clarity.

Para-Pak 10% Formalin preserved specimens may be examined directly and concentrated for recovery of eggs, larvae, and protozoan cysts.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

Each kit consists of one vial containing 15 mL LV-PVA-fixative and one vial containing 15 mL 10% formalin preservative; each with built-in collection spoons. Simple directions for patients and nursing personnel are also provided.

Catalog #300112

10% Buffered Neutral Formalin (item #9004)

LV-PVA (item #9001)

Catalog #330284

10% Buffered Neutral Formalin (item #9004)

LV-PVA (item #9001)

C&S Vials (item #9006)

Catalog #900112

LV-PVA (item #9001)

MATERIALS NOT PROVIDED

1. Ethyl acetate (suggested) or diethyl ether (optional)
2. Zinc sulfate solution (specific gravity = 1.18)
3. Physiological saline
4. Cotton tipped applicator sticks
5. Microscope slides and coverslips
6. Centrifuge
7. Microscope
8. Transfer pipettes

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Avoid contact of fixative solution with the skin and eyes. Should contact occur, flush with running water. If irritation should develop, see a physician.
3. Fixative solutions are poisonous. If ingested, dilute by drinking milk or water. Then call local poison center or physician immediately.
4. LV-PVA fixative is corrosive for metals. Keep out of contact with metals. LV-PVA is flammable.
5. Due to the infectious nature of unpreserved stools, care and handwashing should be employed when the specimen is collected and manipulated.
6. Formalin is a potential cancer hazard.
7. Any serious incident that has occurred in relation to the device should be reported to Meridian Bioscience, Inc., 3471 River Hills Drive, Cincinnati, Ohio 45244 USA or Technical Support Center 800.343-3858.
8. **IMPORTANT:** See SDS for additional safety and hazard information.

Para-Pak 10% Formalin

Signal Word

Danger

Hazard Statements

H301 - Toxic if swallowed
 H311 - Toxic in contact with skin
 H317 - May cause an allergic skin reaction
 H341 - Suspected of causing genetic defects
 H350 - May cause cancer
 H402 - Harmful to aquatic life
 H330 - Fatal if inhaled
 H315 - Causes skin irritation
 H318 - Causes serious eye damage
 H370 - Causes damage to organs

Precautionary Statements - EU (§28, 1272/2008)

P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor/ physician
 P321 - Specific treatment (see supplemental first aid instructions with this material)
 P280 - Wear protective gloves/ protective clothing
 P403 + P233 - Store in a well-ventilated place. Keep container tightly closed
 P280 - Wear eye protection/ face protection
 P321 - See SDS Section 4 or Section 11 for specific medical treatment information
 P201 - Obtain special instructions before use
 P281 - Use personal protective equipment as required
 P308 + P313 - IF exposed or concerned: Get medical advice/ attention
 P202 - Do not handle until all safety precautions have been read and understood
 P264 - Wash face, hands and any exposed skin thoroughly after handling
 P270 - Do not eat, drink or smoke when using this product
 P260 - Do not breathe dust/fume/gas/mist/vapors/spray
 P271 - Use only outdoors or in a well-ventilated area
 P284 - Wear respiratory protection
 P272 - Contaminated work clothing should not be allowed out of the workplace
 P307 + P311 - IF exposed: Call a POISON CENTER or doctor/ physician
 P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
 P310 - Immediately call a POISON CENTER or doctor/ physician
 P302 + P352 - IF ON SKIN: Wash with plenty of soap and water
 P312 - Call a POISON CENTER or doctor/ physician if you feel unwell
 P361 - Remove/Take off immediately all contaminated clothing
 P332 + P313 - If skin irritation occurs: Get medical advice/ attention
 P363 - Wash contaminated clothing before reuse
 P304 + P340 - IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing
 P310 - Immediately call a POISON CENTER or doctor/ physician
 P405 - Store locked up
 P403 + P233 - Store in a well-ventilated place. Keep container tightly closed
 P501 - Dispose of contents/ container to an approved waste disposal plant.

Para-Pak LV-PVA

Signal Word

Danger

Hazard statements

H301 - Toxic if swallowed
 H311 - Toxic in contact with skin
 H314 - Causes severe skin burns and eye damage
 H331 - Toxic if inhaled
 H336 - May cause drowsiness or dizziness
 H341 - Suspected of causing genetic defects
 H361 - Suspected of damaging fertility or the unborn child
 H373 - May cause damage to organs through prolonged or repeated exposure
 H411 - Toxic to aquatic life with long lasting effects
 H225 - Highly flammable liquid and vapor
 Contains Isopropyl alcohol, Acetic acid, Mercuric chloride, Hydrochloric acid

Precautionary Statements - EU (§28, 1272/2008)

P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor/ physician
 P321 - Specific treatment (see supplemental first aid instructions with this material)
 P280 - Wear protective gloves/ protective clothing
 P322 - Specific measures (see .? on this label)
 P321 - See SDS Section 4 or Section 11 for specific medical treatment information
 P403 + P233 - Store in a well-ventilated place. Keep container tightly closed
 P260 - Do not breathe dust/fume/gas/mist/vapors/spray
 P280 - Wear protective gloves/ protective clothing/ eye protection/ face protection
 P303 + P361 + P353 - IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/shower
 P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
 P310 - Immediately call a POISON CENTER or doctor/ physician
 P280 - Wear eye protection/ face protection
 P281 - Use personal protective equipment as required
 P201 - Obtain special instructions before use
 P202 - Do not handle until all safety precautions have been read and understood
 P308 + P313 - IF exposed or concerned: Get medical advice/ attention
 P314 - Get medical advice/ attention if you feel unwell
 P501 - Dispose of contents/ container to an approved waste disposal plant.
 P210 - Keep away from heat/sparks/open flames/hot surfaces. - No smoking
 P370 + P378 - In case of fire: Use dry sand, dry chemical or alcohol-resistant foam for extinction

SHELF LIFE AND STORAGE

Expiration dating for Para-Pak System is indicated on the outer packaging. Store at room temperature (15-30 C). Excessive heat and cold should be avoided.

SPECIMEN COLLECTION AND PREPARATION

1. The patient should be cautioned against the use of antacids, barium, bismuth, antidiarrheal medication, or oily laxatives prior to collection of the specimen.
2. To assure recover of parasitic elements which are passed intermittently and in fluctuating numbers, three specimens spaced a few days apart must be examined. In the case of hospitalized patients, it is suggested that all fecal passages be collected for a designated length of time to avoid prolonging the hospital stay.^{5, 10}
3. The specimen is ideally passed into a bedpan but must not be contaminated with urine. Alternatively, a large plastic bag or "saran wrap" may be placed over the toilet seat opening and the specimen passed into the bag. A thoroughly cleaned and dried milk carton, cut so as to remove the upper two thirds of the carton may also be used. It will be easier to collect the specimen if the water supply to the toilet is shut off and water drained from the bowl by "flushing twice".
4. An appropriate (i.e., bloody, slimy, watery) area of stool should be selected and sampled with the collection spoons provided in the caps of the containers. Sufficient stool is added to each container to bring the liquid level up to the "Fill to Here" line. This will result in approximately 5 mL of sample. To insure adequate sampling of the formed stool, material should be removed from the sides, ends, and middle of the bolus.
5. Agitate each specimen with the spoon along the sides of the container, tighten the cap and shake firmly to insure that the specimen is adequately mixed. When mixing is completed, the specimen should appear homogenous.
6. Return the vials to their container, seal the container, and label appropriately.

TEST PROCEDURE

The Para-Pak Systems lend themselves to a wide variety of procedures in common use. The following discussion is not exhaustive and alternatives may be found in the literature cited. While variations may exist from lab to lab, a thorough examination should include at least four steps:

1. Gross examination:
Record the presence of blood, worms, mucus, or proglottids.
2. Direct microscopic examination from the 10% formalin preserved specimen:
 - a. Place a clean glass slide on a sheet of newsprint.
 - b. Add a drop of saline (iodine may be substituted if desired) to the slide.
 - c. Add a representative sample of formalin preserved specimen to the drop of saline and mix thoroughly with the collecting spoon. The newsprint must be just legible through the slide.
 - d. Place a double width coverslip on the suspension and examine immediately.
3. Permanent slides for staining with iron-hematoxylin, Wheatley's Trichome (Meridian Catalog #400101), etc.:
 - a. Pour some of the LV-PVA fixed material into a paper towel and allow to stand for 3 minutes. This will absorb out excess LV-PVA and is considered critical to obtain the best possible staining.⁹
 - b. Using an applicator stick or brush, spread (avoid smearing) some of the specimen from the paper towel onto one or more clean glass slides. For best adherence, spread the material to the edge of the slide.
 - c. The slides are dried overnight at room temperature or for several hours in a 37 C incubator or slide warmer. Accelerated drying may cause some morphological distortion. The slides must be dried thoroughly to avoid washing the film off during staining.
4. Concentration procedures:

One or more concentration procedures should be employed. No one concentration procedure works equally well for all parasites,^{1, 5, 14} however, two in common usage that lend themselves well to use with the Para-Pak Systems are:

A. Formalin-ether (ethyl acetate)^{7, 18} sedimentation:

1. Mix the 10% Formalin or LV-PVA specimen thoroughly. The specimen is now ready for processing with the Para-Pak Macro-Con or Para-Pak CON-Trate Stool Concentration Systems. See the appropriate package insert for further directions.
If Para-Pak Macro-Con or Para-Pak CON-Trate® are not available, a sufficient quantity of specimen must be strained into a 15 mL conical centrifuge tube through one layer of narrow mesh or two layers of wide mesh gauze to provide the 1.0 mL of sediment. This amount will vary with the size and density of the specimen.
2. Add approximately 6-8 mL of 10% Formalin (or saline), mix thoroughly and allow to stand five minutes.
3. Add 3 mL of ethyl acetate or ether then stopper and shake the tube vigorously for at least 30 seconds. Carefully remove the stopper.
4. Centrifuge at 500 xg for 10 minutes (1800-2200 rpm for most tabletop centrifuges).
5. Four layers will be apparent:
 - a. Top layer: ethyl acetate or ether
 - b. Second layer: plug of debris
 - c. Third layer: formalin
 - d. Bottom layer: sediment
6. After ringing the plug of debris from the sides of the tube with an applicator stick, carefully decant the top 3 layers. While keeping the tube inverted, a cotton swab may be used to remove debris sticking to the sides of the tube. This is particularly important for obtaining suitable results with ethyl acetate and avoids solvent bubbles in the wet mount.
7. Add a few drops of physiological saline or 10% Formalin to resuspend the remaining sediment. If the resulting slides are too dense (newsprint should be legible through them) more saline or formalin may be added.
8. Iodine and saline mounts are suggested for microscopic examination.
NOTE: Rarely, a specimen may be extremely thick or mucoid, necessitating a preliminary wash. The wash, in physiological saline or tap water, may be implemented between Steps 1e and 2. The washed material should be centrifuged as in Step 4.

B. Zinc Sulfate flotation:

1. Thoroughly mix a representative portion of the 10% formalin stool suspension or **fresh** unpreserved specimen in a 15 mL centrifuge tube and q.s. with tap water to approximately 10-12 mL. The amount of specimen to use will vary with its size and density.
2. Centrifuge one minute at 1000-1200 xg.
3. If sediment is about one milliliter in volume, decant the supernatant fluid. Otherwise adjust the density of the suspension by adding material from the 10% Formalin suspension or diluting with more water. If adjustment of the sediment was necessary or if the stool is very oily, repeat the wash procedure.
4. When using formalinized specimens, the specific gravity of the zinc sulfate solution must be adjusted to 1.2.^{1, 5} Fill the tube about half full with zinc sulfate solution and resuspend the sediment by mixing thoroughly with applicator sticks.
5. Add additional zinc sulfate solution to within a thumb-width of the top.
6. Centrifuge one minute at 1000-1200 xg.
7. Carefully remove the tube from the centrifuge and, avoiding agitation, place it upright in a test tube rack or other suitable holder.
8. Carefully fill the tube to the brim with zinc sulfate solution. Do not allow any overflow.
9. A clean coverslip may now be placed on top of the tube. If the coverslip does not contact the meniscus of the liquid, carefully add more zinc sulfate solution until it does.
10. Do not disturb the tube or coverslip for ten minutes.
11. With a quick, deft motion remove the coverslip straight upward so that a drop of liquid containing eggs and cysts adheres to the center of the coverslip.
12. The coverslip may now be placed on a clean glass slide. If an iodine mount is desired, place a small drop of iodine on the slide prior to adding the coverslip. Sealing the edge of the slide with Vaspar (Vaseline/paraffin (1:1) mixture) will prevent drying and distortion of the larger eggs.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

1. Visual inspection: Vials should contain approximately 15 mL of clear fluid and be free from gelling.
2. If gelled, the fixative may be cleared by placing in a 50 C water bath until clear and fluid. Fixative vials containing patient specimens may be cleared once or twice in this manner also, provided the product is used cooled.
3. On prolonged storage, a thin haze of sediment may develop in the bottom of the PVA vial. This will not affect adhesive, staining or fixative properties.
4. When a PVA-fixed film of stock trophozoite or human buffy coat is stained, the organisms or cells should appear well fixed and defined.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

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REFERENCES

1. Bartlett, Marilyn, et al. Comparative evaluation of a modified zinc sulfate floatation technique. J Clin Microbiol. 1978;7:524-528.
2. Brooke MM. and Goldman M. Polyvinyl alcohol fixative as a preservative and adhesive for protozoa in dysenteric stools and other liquid materials. J Lab and Clin Med. 1944;34:1554-1560.
3. Brooke MM. and Norman C. The effectiveness of the PVA fixative technique in revealing intestinal amoebae in diagnostics cultures. Am J Trop Med Hyd. 1955;4:479-482.
4. Brooke MM. PVA fixative technique in laboratory confirmation of amebiasis. Triangle 1960;4:326-335.
5. Brooke MM. Intestinal and urogenital protozoa. Manual of Clinical Microbiology, ASM, Washington, D.C. 2nd ed. 1974;582-601.
6. Burrows RB. Microscopic diagnosis of the parasites of man. Yale University Press, New Haven. 1965.
7. Erdman Dean. Clinical comparison of ethyl acetate and diethyl ether in the formalin-ether sedimentation technique. J Clin Microbiol. 1981;14:483-385.
8. Faust EC. et al. A critical study of clinical laboratory techniques for the diagnosis of protozoan cysts and helminth eggs in feces. Am J Trop Hyg Med. 1938;18:169-183.
9. Garcia LS and Ash Lawrence. Diagnostic Parasitology, C. V. Mosby Co., St. Louis, 1979;pp 9-24.
10. Harper K. et al. Advantages of the PVA fixative two-bottle stool collection technique in the detection and identification of intestinal parasites. Pub Health Lab 1957;15:96-108.
11. Markell EK and Quinn Patricia M. Comparison of immediate and polyvinyl alcohol (PVA) fixation with delayed Schaudin's fixation for the demonstration of protozoa in stool specimens. Am J Trop Med Hyg. 1977;26:1139-1142.
12. Melvin DM and Brooke MM. Laboratory procedures for the diagnosis of intestinal parasites, U. S. Department of Health, Education and Welfare Publication (CDC) – 1975;75-82821, Atlanta.
13. Ritchie LS. et al. A comparison of the zinc sulfate and the MGL (formalin-ether) techniques. J Parasit. 1952;38:16.
14. Scholten TH. and Yang J. Evaluation of unpreserved and preserved stools for the detection and identification of intestinal parasites. Am J Clin Path. 1974;62:563-567.
15. Swartzwalder J, Clyde Hunter GW, Frye WW. Manual of Tropical Medicine, W. B. Saunders Co., Philadelphia. 1966.
16. Simitch T. and Petrovich Z. Longevite de la forme vegetive de dysenterioe dans divers mileaux. Arch. Inst. Pasteur d'Algerie 1953;31:375-380.
17. Wheatley WB., A rapid staining procedure for intestinal amoebae and flagellates. Am J Clin Path. 1951;2:990-991.
18. Young Kirk H. et al. Ethyl acetate as a substitute for diethyl ether in the formalin-ether sedimentation technique. J Clin Microbiol. 1979;10:852-853.

SN10912

REV. 02/24

 Manufactured By	Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com <u>Contacts:</u> Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com
--	---

INTERNATIONAL SYMBOL CHART

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.