Para-Pak® TRICHOME STAIN

REF 400101, 400112

IVD

Rx Only

INTENDED USE

Para-Pak Trichrome Stain is a rapid staining procedure which provides excellent differentiation of internal structures of intestinal parasites as well as facilitating the separation of these organisms from background material and artifacts.

SUMMARY AND EXPLANATION OF THE TEST

The trichrome stain was first introduced by Masson in 1929 as an easy, rapid alternative to tissue stains then in use. It was particularly useful at the time for the study of muscles. A number of authors have modified the stain for suitability to the study of other tissues, but in 1949 Gomori developed a more rapid methodology which, after preliminary nuclear staining with hematoxylin, was broadly applicable to tissue sections and cytologic smears.^{1,3,8,10} In 1951 Wheatley noted the affinity of the chromatin material of *E. histolytica* for chromotrope 2R and found that with changes in fixation and dehydration a simple and rapid staining procedure for intestinal amoeba and flagellates was possible.¹¹ Para-Pak Trichrome Stain is suggested for use in the regimen of Wheatley for intestinal parasites which follows. It is also suitable for use in staining tissues with the Gomori technique. The contrast of the blue-green, purplish-tinged or red-stained organisms with the green background material is a considerable aid to locating organisms when compared with iron-hematoxylin.^{4,9}

REAGENTS/MATERIALS PROVIDED

1. Chromotrope 2R
2. Light green SF
3. Fast green FC
4. Phosphotungstic Acid
5. Acetic Acid
6. Distilled H₂O

MATERIALS NOT PROVIDED**Solution A – Schaudin's fixative:**

95% Ethanol 33.3 mL
Mercuric chloride 66.6 mL
(saturated aqueous solution about 14 g/100 mL H₂O). Just prior to use add 5 mL of glacial acetic acid/100 mL of solution.

Solution B – Iodine alcohol:

Prepare stock solution: To 70% ethanol add enough iodine crystals to make a dark, concentrated solution. Amount will vary from lot to lot of iodine.

Prepare working solution: Dilute stock solution to port wine, strong tea and amber colored solution by adding 70% ethanol. The exact concentration is unimportant.

Solution C – Acid ethanol:

To 100 mL 90% Ethanol add .5 mL glacial acetic acid.

Solution D – Carbol-xylene:

Melt sufficient crystals of carbolic acid (phenol) to produce 100 mL of liquid in a water bath and add to 300 mL xylene or add 100 gm phenol crystals to 300 mL xylene and dissolve by shaking.

PRECAUTIONS

1. All reagents are for *in vitro* diagnostic use only.
2. When utilizing PVA fixed specimens the slides must be absolutely dry. Accelerating the drying process by warming may distort some organisms.
3. When utilizing Schauden's fixed specimens the slides must not be allowed to dry at any time until finally mounted.
4. If necessary to interrupt the staining procedure the slides may be stored for long periods of time in the last ethanol bath prior to staining. (**Smear Preparation:** Step No. 3).
5. The purpose of the iodine alcohol (**Staining:** Step No. 1) is to remove the mercuric chloride in Schaudin's fixative and PVA fixative. Many problems experienced in the trichrome method are directly attributable to saturation of this solution; change it frequently.
6. The purpose of the first alcohol rinse after decolorizing is to halt decolorization. Since carryover of the acid alcohol will occur, change this solution frequently.

HAZARD and PRECAUTIONARY STATEMENTS

| | |
|--|--|
|  Trichrome Stain | Signal word Warning Hazard statements H313 - May be harmful in contact with skin Precautionary Statements - EU (§28, 1272/2008) P262 - Do not get in eyes, on skin, or on clothing P302 + P352 - IF ON SKIN: Wash with plenty of soap and water |
|--|--|

SHELF LIFE AND STORAGE

Shelf life of the Para-Pak Trichrome Stain is indicated on the outer package label. Store at 15-30 C. Excessive heat and cold should be avoided.

TEST PROCEDURE**Smear Preparation:**

For preparation of smears from PVA fixative or SAF fixative, please refer to the appropriate Meridian Bioscience package insert supplied with the collection kit or fixative. A complimentary insert will be provided on request.

1. Make a thin smear of the fecal material on a clean glass slide. The specimen may be emulsified in a drop of saline if necessary. Alternatively, a smear may be prepared in PVA fixative and allowed to dry overnight. In this case proceed directly to **Staining:** Step No. 1.
2. While the smear is still wet, immediately place the slide in a coplin jar containing Schaudin's fixative (Solution A below) for 5 minutes at 50 C or one hour at room temperature.
3. With the exception that the iodine may be omitted from the first alcohol solution for those slides prepared in fixatives not containing mercuric chloride (such as SAF) all specimens may now be processed identically.

Staining:

1. Iodine alcohol (Solution B below) (See Note 1) 10 minutes
 2. 70% ethanol (See Note 3) 3-5 minutes
 3. 70% ethanol (See Note 3) 3-5 minutes
 4. Para-Pak Trichrome Stain 6-8 minutes
 5. Acid ethanol (Solution C below) 5-10 seconds
- Observe carefully. When stain just begins to run from the smear, transfer to the next bath immediately. Two quick dips are usually sufficient.
6. 95% ethanol (See Notes 2 & 4) Dip twice
 7. 95% ethanol (See Note 4) 5 minutes
 8. 95% ethanol (See Note 4) 5 minutes
 9. 100% ethanol (absolute) or carbol-xylene (solution D) 3 minutes
 10. Xylene 3 minutes
- (or until refraction appears at the smear-xylene interface)
11. Mount with coverslip using the mounting medium of choice (e.g., Canada balsam in xylene, D.P.X., Permount, Crisalite or clear mount).

Note 1: See Precautions.

Note 2: See Precautions.

Note 3: "Reagent Alcohol", commercially available mixture of 95% ethanol, 5% methanol, may be substituted if diluted to 70% in water.

Note 4: "Reagent Alcohol", commercially available as a mixture of 95% ethanol, 5% methanol, may be substituted without dilution.

INTERPRETATION OF RESULTS

Background material and artifacts will stain green. Bacteria and red blood cells usually stain red, while yeasts, molds, etc., are usually green. The cytoplasm and cysts of *E. histolytica* will be blue-green; usually with a slight purplish tinge. *E. coli* cysts may be somewhat more purplish. Karyosomes of nuclei, chromatoid bodies and chromatin material will appear red to purplish as will helminth eggs and larvae.

TROUBLESHOOTING

The Wheatley trichrome method is a remarkably trouble-free procedure when used exactly as directed. Occasionally, as with any procedure, problems may arise. These may be divided into four categories:

1. Problems related to incomplete fixation of the smear resulting from failure to thoroughly emulsify the fecal specimen in the fixative.
 - a. Degenerate forms staining a pale green.
 - b. Unsatisfactory staining of cytoplasm and/or nucleus.

Be sure to thoroughly mix the specimen in the fixative.
2. Failure to observe the smear carefully while decolorizing (**Staining:** Step No. 5).
 - a. Lack of contrast between organisms and background material of cytoplasm and nucleus.

Remove the slide from the acid-ethanol and proceed to the ethanol rinse as soon as the color begins to run from the smear.
3. Saturation of the iodine-alcohol when using fixatives containing mercuric chloride:
 - a. Obscuration of the smear by dark crystals.

Prepare the iodine alcohol solution exactly as directed and change it frequently.
4. Carryover of staining/decolorizing solutions from one jar to the next:
 - a. Smear too green or lack of contrast.
 - b. Smear is "cloudy".

Lack of contrast usually results from carryover of the acid-alcohol to the ethanol rinses or failure to observe the slide while decolorizing. After repeated use the stain solution may be "watered down" by carryover of 70% ethanol. "Cloudy" smears are caused by incomplete dehydration after staining and decolorization. Carbol-xylene or absolute ethanol and xylene are dehydrating agents in this procedure.

Change all solutions periodically to avoid problems.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

1. A positive smear should be run through the series each time the solutions are changed and the results recorded.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

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REFERENCES

1. Belding DC, Textbook of clinical parasitology. New York; D. Appleton-Century Company: 1947.
2. Biological Stain Commission, Staining procedures. Baltimore; Williams and Wilkins: 1960
3. Gomori G, A rapid one step trichrome method. Am J Clin Pathol, 1950; 20:661-664.
4. Gurr E, Staining: practical and theoretical. Baltimore; Williams and Wilkins; 1962. p. 452-455.
5. U. S. Armed Forces Institute of Pathology, Laboratory manual of special staining technics. Washington (DC): 1953. p. 58-59.
6. Lillie RD, et. al. Ed. H. J. Conn's Biological stains. Baltimore; Williams and Wilkins; 1968. p. 76.
7. Markel AK, A comparison of three staining techniques for protozoan troponozites as applied to PVA-preserved fecal specimens. J Parasitol 1956; 42:47.
8. Masson P, Trichrome staining and their preliminary techniques. J Tech Meth 1929; 12:75-90.
9. Paik G and Morris T Suggs. Reagents, stains, and miscellaneous test procedures. Washington, DC: Chapter 96 in Manual of Clinical Microbiology, 2nd ed, A.S.M 1974.
10. Pollak OJ, 1944. A rapid trichrome stain. Arch. Path. Vol. 37:294.
11. Wheatley WB. A rapid staining procedure for intestinal amoeba and flagellates. Am J Clin Pathol 1951; 21:990-991.

SN10900

REV. 08/22

| | |
|---|--|
|  Manufactured By | <p>Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com</p> <p>Contacts: Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com</p> |
|---|--|

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product: Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simboloas, Zeichenerklärung)

| | | | |
|----------------|---|--|---|
| | Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis | CONTROL | Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle |
| LOT | Batch Code / Codice del lotto / Code du lot / Código de lote / chargebezeichnung | CONTROL | Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle |
| IVD | In vitro diagnostic medical device / Dispositivo medico-diagnóstico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum | EC REP | Authorized representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft |
| | Meridian products carrying the European Conformity (CE) mark fulfill the requirements of Directive 98/79/EC or the Regulation 2017/746 on in-vitro diagnostic medical devices / I prodotti Meridian recanti il marchio di Conformità Europea (CE) soddisfano i requisiti della Direttiva 98/79/CE o del Regolamento 2017/746 sui dispositivi medico-diagnostici in vitro / Les produits Meridian portant la marque de Conformité européenne (CE) sont conformes aux exigences de la Directive 98/79/CE ou du règlement 2017/746 sur les dispositifs médicaux de diagnostic in vitro / Los productos de Meridian que llevan la marca de conformidad europea (CE) cumplen los requisitos de la Directiva 98/79/CE o del Reglamento 2017/746 sobre los productos sanitarios para diagnóstico in vitro / Produkte von Meridian mit der CE-Kennzeichnung erfüllen die Anforderungen der EU-Richtlinie 98/79/EG bzw. der Verordnung (EU) 2017/746 über In-vitro-Diagnostika | SMP PREP DIL SPE | Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet |
| REF | Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer | | CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risques de danger / Precaución: Peligroso / WARNUNG: Risikogefahr |
| | Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten | BUF RXN | Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer |
| | Manufacturer / Fabbricante / Fabricant / Fabricante / Hersteller | | For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung |
| | Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" test / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen | SOLN STOP | Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung |
| | Temperature limitaion / Limiti di temperatura / Limites de température / Límite de tempperatura / Temperaturbegrenzung | CONJ ENZ | Enzyme Conjugate / Coniugato enzimatico / Conjugué enzymatique / Conjulado enzimático / enzymkonjugat |
| SN | Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer | CONTROL | Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest |
| TEST | Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / testgarät | REAG | Reagent / Reagente / Réactifs / Reactivos / Reagenzien |
| | Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum | BUF WASH | Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer |
| BUF | Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer | | Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise |
| CONJ | Conjugate / Coniugato / Conjugué / Conjulado / Konjugat | DIL SPE | Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer |
| SUBS | Substrate / Substrato / Substrat / Substrato / Substrat | BUF WASH 20X | Wash Buffer Concentration 20X / Soluzione dil lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat |
| Rx Only | Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig | DET REAG | Detection Reagent / Reagente Diretto / Réactif de Detection / Reactivo de Detección / Nachweis Reagenz |
| | Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / ne pas utiliser si le paquet est endommagé / No use si el paquete está dañado / Nicht verwenden, wenn die Verpackung beschädigt ist | TUBE | Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß |
| | Single Use Only / Prodotto Monouso / A usage unique / Para Un Solo Uso / nur für die einmalige Anwendung | | |

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.