For use with the Revogene®

FOR EXPORT USE ONLY; NOT FOR SALE IN THE USA

REF 410100

IVD For *in vitro* Diagnostic Use

### INTENDED USE

The Revogene® GBS DS assay performed on the Revogene® instrument is a qualitative *in vitro* diagnostic test designed to detect Group B *Streptococcus* (GBS) DNA from vaginal/rectal swabs obtained from pregnant women. The Revogene GBS DS assay utilizes automated sample processing and real-time polymerase chain reaction (PCR) to detect a *cfa* gene sequence specific to the *Streptococcus agalactiae* genome.

Revogene GBS DS assay is indicated for the identification of *intrapartum* GBS colonization and does not provide susceptibility results. It is not intended to diagnose or monitor treatment of GBS infection. Culture isolates are needed for performing susceptibility testing as recommended for penicillin-allergic women.

### SUMMARY AND EXPLANATION

Group B *Streptococcus* is a Gram-positive bacterium commonly found in the gastrointestinal, genital, and urinary tract of healthy adults. Approximately 10-30% of all pregnant women are colonized with GBS in the vagina or rectum.<sup>1</sup>

GBS colonized mothers typically show no symptoms or health effects. However, this bacterium can be transmitted to the newborn mainly during labor and birth (*intrapartum*)<sup>2</sup> and can be associated with serious illnesses in newborns such as sepsis, pneumonia, and meningitis. The fatality rate for infants with early on-set Group B Streptococcal disease (during the first week of life) is currently estimated at 4% to 6%.<sup>1</sup> Surviving infants may experience long-term disabilities including hearing loss, vision loss or mental retardation.<sup>2</sup>

The current guidelines for preventing neonatal GBS disease is screening pregnant women at 35-37 weeks of gestation (*antenpartum*) to determine their GBS colonization status. Most GBS testing is performed by culture and can take up to 48 hours for definitive identification of GBS following the initial 18 to 24 hours of incubation of vaginal/rectal swab in a selective broth medium.

If the colonization status is positive, an *intrapartum* antibiotic prophylaxis (IAP) is given to women in order to prevent perinatal Group B streptococcal disease. However, GBS colonization can be intermittent during pregnancy, resulting in an insufficient correlation between antenatal screening results and *intrapartum* GBS colonization.<sup>6</sup> Moreover, the *antenpartum* screening excludes women who have preterm delivery or women arriving at hospital without known status. In consequence, some women may not be appropriately treated for GBS colonization while other may be unnecessarily treated with antibiotics.<sup>7</sup>

The GBS DS assay can provide results from one (1) up to eight (8) samples as early as 40 minutes for positive specimens and in approximately 70 minutes for negative specimens, by directly using the swab specimen collected from woman in labor (*intrapartum*). The assay minimizes operator intervention from the time the single-use disposable microfluidic cartridge (named PIE hereinafter) containing the sample is placed into the Revogene carousel until results are available.

### PRINCIPLES OF THE PROCEDURE

The Revogene GBS DS assay is performed using the Revogene instrument which automates sample homogenization and dilution, cell lysis, DNA amplification and detection of the amplified PCR products. User intervention is only required to discharge the patient specimen into the Sample Buffer Tube (SBT), transfer the sample from SBT into the PIE and insert the PIEs into the Revogene carousel.

Each PIE is a fully integrated and closed device in which a sample is dispensed and processed through different microfluidic chambers and channels which allow for the sample processing (i.e. sample homogenization and dilution, cell lysis and DNA extraction) and subsequent real-time PCR steps (Figure 1). The liquid from a single sample is transferred by centrifugation from one (1) chamber to the next in sequence and all reagents specific for the PCR reaction are incorporated and dried within the PCR well. A Process Control (PrC) is incorporated into each PIE to verify sample processing and amplification steps including the verification of potential inhibitor substances as well as microfluidic, instrument or reagent failure. The amplified products are detected in real time using target-specific TaqMan® chemistry-based probes. No operator intervention is necessary once a PIE is loaded into the Revogene.


Figure 1. Top View of a PIE.

1: Sample Loading Chamber, 2: Homogenization Chamber containing PrC, 3: Overflow Chamber,  
4: Dilution/Lysis Chamber, 5: Three (3) PCR Wells (A to C from left to right) and one (1) Waste Chamber (at the right end).

The Revogene can process from one (1) up to eight (8) samples simultaneously in the same run. The carousel must contain eight (8) PIEs to maintain thermodynamic balance within the run. During the run and at run completion, the results are computed by the system from measured fluorescent signals and embedded calculation algorithms. Results that are displayed on the touchscreen may be printed, transferred and stored by the user using the USB port or the connectivity option.

An Early Positive Result Outcome (E-PRO) feature provides positive result if the signal from the target DNA reaches a predetermined threshold before the full PCR cycles have been completed. With E-PRO, a positive result could be obtained in approximately 40 minutes for highly positive samples. For negative samples, the time to result is approximately 70 minutes.

## REAGENTS AND MATERIALS

The Revogene GBS DS kit contains sufficient reagents and materials to process 24 samples. The kit contains the following materials:

1. 24 **GBS DS Sample Buffer Tubes** (SBT): Barcode-labeled tube containing TE 1X buffered solution (Tris-HCl pH 8.0/EDTA.Na<sub>2</sub>) as a dilution and preservation buffer for sample.
2. 24 **Disposable Transfer Tools** (DTT): Plastic pipette with minimal and maximal volume marks for transferring the sample from the SBT to the PIE.
3. 24 individual pouches containing one (1) **GBS DS Microfluidic cartridge** (PIE): Barcode-labeled integrated device, composed of dried reagents allowing sample process and real-time PCR steps for PrC DNA and GBS *cfb* gene DNA simultaneous amplification and detection. Each PIE contains PrC, PrC-specific primers and TaqMan® chemistry-based probe, GBS *cfb* gene-specific primers and TaqMan® chemistry-based probe, dNTPs, buffer and DNA polymerase.

## MATERIALS/EQUIPMENTS REQUIRED BUT NOT PROVIDED

1. Revogene® (cat# 610210)
2. Disposable powderless gloves
3. Vortex mixer with a maximal speed of ≥ 3000 rpm
4. Sample rack (cat# 132539; optional)
5. Copan Liquid Stuart (cat. # 141C)
6. MOCK PIE (cat# 610208; optional)
7. Clean scissors (optional)
8. Gauze (optional)

## WARNING AND PRECAUTIONS

1. The Revogene GBS DS assay can only be used on the Revogene.
2. Do not use the kit if the label that seals the outer box is broken upon arrival.
3. Do not use GBS DS PIEs if the protective pouches are open or broken upon arrival.
4. Do not interchange DTT, SBT and PIE between kit lots.
5. Each single-use GBS DS PIE and DTT are used to process one (1) sample. Do not reuse PIE or DTT.
6. Always handle specimens as if they are infectious and in accordance with Good Laboratory Practices such as those described in Biosafety in Microbiological and Biomedical Laboratories<sup>3</sup> and in CLSI Document M29-A4<sup>4</sup>.
7. Wear disposable powderless gloves while handling samples and thoroughly wash hands afterwards.
8. The GBS DS PIE contains dried reagents. The protective pouch should not be opened until ready to perform the test.
9. Dispose of unused materials, reagents and waste in accordance with country, federal, provincial, state and local regulations.
10. Do not open or break apart the PIE after use to avoid contamination with amplification products and/or infectious particles.
11. Do not use a PIE that has been dropped, shaken or inverted after the sample has been loaded, as this may cause invalid results.
12. The Revogene GBS DS assay does not provide susceptibility results. Additional time is required to culture isolates and perform susceptibility testing as recommended for penicillin-allergic women.
13. Do not use a kit that has passed its stated expiration date.
14. Do not refrigerate the loaded PIE.
15. Each run must be performed with eight (8) PIEs in the Revogene carousel to maintain thermodynamic and mechanical balance within the run.

## HAZARD AND PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

## STORAGE AND STABILITY

1. Vaginal/rectal swab specimens should be stored at 25 C if they are processed within two (2) days of collection. Otherwise, they should be stored at 2-8 C for up to seven (7) days.
2. Store the Revogene GBS DS kit at 2-25 C. The expiration date is indicated on the kit box label.
3. Do not open a pouch until ready to perform testing. Use the PIE within one (1) hour after opening the pouch.
4. Inoculated SBT can be stored at 25 C for up to two (2) days, or at 2-8 C for up to three (3) days.

## INSTRUCTION FOR USE

There are no known hazards associated with this product.

## SPECIMEN COLLECTION AND TRANSPORT

**Specimen type:** Vaginal/rectal swab taken from woman in labor (*intrapartum*).

Vaginal/rectal specimen collection from the same patient should be performed in accordance with published guidelines for collection of clinical specimens for culture of GBS.<sup>1</sup>

1. Collect the vaginal/rectal swab according to European guidelines<sup>7</sup>:
  - a. Swab the lower vagina (vaginal introitus), followed by the rectum (i.e. insert swab through the anal sphincter) using the same swab or two different swabs. Take care not to mix stool, water, urine or soap in the sample.
  - b. Place the swab(s) into a non-nutritive transport medium (Liquid Stuart) right after sample collection.
2. Label the transport device with the specimen or patient identification (ID) and transport to the laboratory for testing (refer to **Storage and Stability** section).

## SAMPLE PREPARATION AND HANDLING

**NOTE 1:** Start the test within one (1) hour after opening the pouch containing the PIE.


### SBT PREPARATION

1. For each specimen to be tested, obtain one (1) SBT from the kit box.
2. Identify (or label) the SBT with the appropriate specimen identification without obscuring or writing over the barcode. Place the SBT on the Sample rack, if used.
3. Unscrew the SBT cap. Remove the swab from its transport sheath and transfer it in the SBT. Using a gauze, break the exceeding part of the swab stem by holding the swab by the stem near the rim of the tube. Lift the swab a few millimeters (mm) from the bottom of the tube and bend the stem against the edge of the tube to break it. Alternatively, use clean scissors to cut the stem. Make sure that only one (1) SBT is open at once.
4. Replace the SBT cap, close it tightly, and replace the SBT on the Sample rack (if used).
5. Prepare any additional specimens for testing by repeating steps 1 to 4 then proceed to step 6.
6. When all samples are prepared, proceed to GBS DS PIE preparation (next section).

### PIE PREPARATION

**NOTE 1:** Process one (1) sample at a time.

7. Vortex the SBT for 30 seconds at maximal speed (≥ 3000 rpm) using a vortex mixer.
8. Unseal the PIE pouch and remove the PIE. Once the pouch is unsealed, the PIE must be used within one (1) hour.
9. Place the PIE on the Sample rack, if used or on a flat surface.
10. Obtain one DTT from the kit box. Unscrew the SBT cap, using the DTT, aspirate the inoculated Sample Buffer (SB) by squeezing the entire bulb. The liquid level into the DTT must be anywhere between the two (2) marks (**Figure 2**). If the liquid level is not between the two (2) marks, discharge completely the SB in the SBT and repeat this step.
11. Discharge completely the SB volume contained in the DTT into the Sample Loading Chamber of the PIE (**Figure 1**).
12. Close the cap of the PIE and the SBT cap tightly.
13. Repeat steps 7 to 12 for any additional prepared SBT then proceed to the Revogene operation section.

**Figure 2. Representation of an Appropriate Sample Buffer (SB) Level using the Disposable Transfer Tool (DTT).**

#### REVOGENE OPERATION

**NOTE 1:** Each run must be performed with eight (8) PIEs in the Revogene. When less than eight (8) samples are processed, fill the empty places with MOCK PIEs\*. **NOTE 2:** Refer to the Revogene Operator's Manual<sup>5</sup> for further information regarding instrument set-up and operation.

1. Power on the Revogene (if not already done). The software will launch automatically.
2. Log in by entering the <Username> and <Password> and tap <Login>. The main menu will appear automatically.
3. Tap <Setup Run>.
4. Enter the sample identification using either the barcode scanner or manual entry. Manual entry can be done by tapping the **pencil** icon of the <Scan or Enter Sample ID> line.
5. Enter the SBT and PIE barcodes using the Revogene barcode scanner by gently positioning the PIE vertically in front of the scanner. Alternatively, SBT and PIE barcodes may be entered manually (tap the **pencil** icon of their respective lines). Handle the PIE carefully without shaking or dropping it.
6. (Optional) Tap the **pencil** icon of the <Add Comments> line and type a comment.
7. Insert the PIE into the Revogene, at any position of the carousel. The software will automatically associate sample and SBT to the correct PIE.
8. Confirm that the PIE is inserted into instrument by tapping <OK> on the <insert PIE into instrument> line and repeat steps 4 to 8 for all samples. If less than eight (8) PIEs are being tested, load MOCK PIEs\* in the carousel remaining positions. No scan is required when inserting MOCK PIEs into the Revogene.
9. When all PIEs are inserted into the carousel, tap <Next>.
10. If you have entered MOCK PIEs into the carousel, follow the instructions on the screen.
11. Scan the retention ring and place it on the carousel. Close the instrument lid.
12. Initiate the test run by tapping <Start>. A timer on the screen and lights on the Revogene lid will show the progression of the test.
13. Store inoculated SBT in appropriate conditions (refer to **Storage and Stability** section) for further repeat testing, if needed.


\*If MOCK PIEs are not available, use unused assay PIEs filled with uninoculated SB (BLANK) or with External Controls.

#### VIEWING AND EXPORTING RESULTS

**NOTE 1:** Refer to the Revogene Operator's Manual<sup>5</sup> for further information regarding the acquisition of test results.

#### DURING THE RUN

1. If the lights on the Revogene lid start to blink, this is a notification for an Early Positive Result Outcome (E-PRO). An icon will also appear on the title bar of the screen.
- The E-PRO icon includes the "+" symbol and a number representing the quantity of positive results available at this time (Figure 3). The number will increase if additional positive results become available.

**Figure 3. Main Menu Display showing E-PRO icon.**

2. Enter <Username> and <password> and tap <login> if the user's session has logged-out.
3. Tap E-PRO icon.
4. Result from the current run are automatically listed on the screen (Figure 4). For each sample line, there will be either a Positive symbol or an In-progress symbol. The In-progress symbol, consisting of a rotating animation, is displayed until a positive result is obtained or the run is completed.
5. Tap on a Positive result symbol to open the interim report.
6. (Optional) Tap <Export> and save the interim report where appropriate (e.g. USB key or connectivity option).

Once the results screen has been consulted, the E-PRO notifications will disappear. New E-PRO notifications will appear upon the detection of an additional Early Positive result.

Any Sample result displayed in the interim report is final and will remain unchanged in the final assay result report.

	Sample ID	Patient ID	Assay	Result
1	1300006001		GBS DS	+ (green)
2	1300006002		GBS DS	(grey)
3 ✓	1235227901		GBS DS	+ (green)
4	1295227902		GBS DS	(grey)
5	1235227801		GBS DS	(grey)
6	Empty			
7	Empty			
8	Empty			

Started By: admin Complete At: 14:10 Report Run Report Export Edit Abort

Figure 4. In Progress Results Display.

#### AT THE END OF THE RUN

7. Once the run is completed, the lid opens automatically.
8. If the Revogene has logged-out, re-enter <Username> and <Password> and tap <Login>. The main menu will appear automatically.
9. Tap **Results** icon to access test results. The **Results** window shows the reported result(s) for each sample (Figure 5).
10. Tap <Last Run> to see the latest test results.

	Sample ID	Patient ID	Started	Assay	Result
	1235227901		08/11/2021 12:58	GBS DS	+ (green)
✓	1300006001		08/11/2021 12:58	GBS DS	+ (green)
	1235227801		08/11/2021 12:58	GBS DS	- (grey)
	1295227902		08/11/2021 12:58	GBS DS	- (grey)
	1300006002		08/11/2021 12:58	GBS DS	- (grey)

Last Run Search Report Run Report Export Export Data

Figure 5. Results window showing the reported result(s) for each sample.

11. From the <Last Run>, select samples for which results report(s) has (have) to be exported. All samples can be selected in one time by checking the box on the upper left corner of the screen.
12. Tap <Export> and save where appropriate (e.g. USB key or connectivity option).
13. (Optional) Tap <Search> to find a specific sample and its result.
14. Remove the retention ring and PIES from the Revogene. Used PIES should be discarded in appropriate waste containers according to the institution's standard practices.

#### REPEAT TESTING PROCEDURE

##### UNRESOLVED OR INDETERMINATE RESULT FOR A SPECIMEN

**NOTE 1:** Because of the *intrapartum* setting, repeat testing may not be feasible and will depend on practices and policies within each facility. Coordination between clinicians and the testing laboratory is important to not delay administration of antibiotics while results are pending.

When an Unresolved (UNR) or an Indeterminate (IND) result is obtained for a specimen, a repeat test from the corresponding inoculated SBT must be performed within the specified timeframe described in the **Storage and Stability** section. Only one (1) repeat testing from the SBT is allowed.

Vortex the SBT for a minimum of 30 seconds at maximal speed. Using a new pouch, follow steps 8 to 13 of the **Sample Preparation and Handling / PIE Preparation** section, then follow the **Revogene Operation** section.

##### UNRESOLVED, INDETERMINATE, FALSE NEGATIVE OR FALSE POSITIVE RESULT

##### FOR AN EXTERNAL CONTROL

When an Unresolved, an Indeterminate, a false negative or a false positive result is obtained for an External Control, the run is invalid. Samples included in the run should be repeated using the corresponding inoculated SBT, along with freshly prepared External Controls, within the specified timeframe described in the **Storage and Stability** section. Refer to the next **Quality Control** section for preparation of fresh External Controls.

For the repeat testing using the corresponding inoculated SBT, vortex the SBT for a minimum of 30 seconds at maximal speed. Using a new pouch, follow steps 8 to 13 of the **Sample Preparation and Handling / PIE Preparation** section, then follow the **Revogene Operation** section.

#### QUALITY CONTROL

Quality control procedures monitor the accuracy and precision of the analytical process. Each laboratory must establish the number, type and frequency of testing control materials per applicable regulations or accrediting agencies. The procedure described below may be employed, if appropriate, based on local policies and procedures.

#### PROCESS CONTROL

Each PIE contains a Process Control (PrC) that verifies the sample homogenization and dilution, cell lysis, inhibition of DNA amplification and assay reagents failure.

## EXTERNAL CONTROLS

NOTE 1: Separate DTT, SBT and PIE must be used for each External Control preparation.

- Good laboratory practice recommends the use of control materials. Users should follow the appropriate guidelines concerning the running of External Controls. It is recommended that one (1) Positive External Control and one (1) Negative External Control should be run at least on a daily basis until adequate process validation is achieved with the Revogene GBS DS assay on the Revogene in each laboratory setting.
- External control materials are not provided by Meridian Bioscience, Inc. External Controls are not used by the Revogene software for the purpose of sample test result interpretation. External Controls are treated as if they are specimens.
- Before reporting early positive sample results, ensure to have Valid Positive and Negative External Control results.
- Commercially available control material, previously characterized positive sample or negative sample spiked with well characterized organisms can be used as External Positive Control. Positive External Control should be used in accordance with the applicable regulations or accrediting agencies.
- Commercially available control material, previously characterized negative sample can be used as External Negative Control. Negative External Control should be used in accordance with the applicable regulations or accrediting agencies.

## RESULTS INTERPRETATION

The results are computed by the Revogene from measured fluorescent signals and embedded calculation algorithms and are available on the "Results" window of the Revogene. Possible results are:

Sample	Symbol	Result	Interpretation
Patient specimen		Positive	GBS target DNA detected.
		Negative	GBS target DNA not detected.
		Unresolved	Amplification/detection failure of the Process Control as well as for the GBS target DNA. Could be caused by inhibitory specimens, microfluidic or reagent failure. Repeat testing must be performed (refer to the <b>Repeat Testing Procedure</b> section for further guidance).
		Indeterminate	No reportable result due to possible Revogene detection error during the assay processing, the data analysis, or if the run is interrupted by the user. Repeat testing must be performed (refer to the <b>Repeat Testing Procedure</b> section for further guidance).
		In progress	Results not available yet.
Positive External Control		Positive	Valid Positive External Control result.
		Negative	An External Positive Control that yields a negative result is indicative of a sample handling/preparation problem. <b>The run is invalid.</b> Review the specimen handling/preparation technique. Repeat testing must be performed (refer to the <b>Repeat Testing Procedure</b> section for further guidance).
		Unresolved	Incorrect Positive External Control result. <b>The run is invalid.</b> Repeat testing must be performed (refer to the <b>Repeat Testing Procedure</b> section for further guidance).
		Indeterminate	Incorrect Positive External Control result. <b>The run is invalid.</b> Repeat testing must be performed (refer to the <b>Repeat Testing Procedure</b> section for further guidance).
		In progress	Results not available yet.
Negative External Control		Positive	An External Negative Control that yields a positive test result is indicative of a sample handling and/or contamination event. <b>The run is invalid.</b> Review the specimen handling technique. Repeat testing must be performed (refer to the <b>Repeat Testing Procedure</b> section for further guidance).
		Negative	Valid Negative External Control result.
		Unresolved	Incorrect Negative External Control result. <b>The run is invalid.</b> Repeat testing must be performed (refer to the <b>Repeat Testing Procedure</b> section for further guidance).
		Indeterminate	Incorrect Negative External Control result. <b>The run is invalid.</b> Repeat testing must be performed (refer to the <b>Repeat Testing Procedure</b> section for further guidance).
		In progress	Results not available yet.

## LIMITATIONS OF THE PROCEDURE

- The Revogene GBS DS assay must only be used with the Revogene by trained personnel.
- A positive test result does not necessarily indicate the presence of viable organisms. It is, however, indicative for the presence of GBS DNA.
- Use of the Revogene GBS DS assay only for clinical specimen collected with Copan Liquid Stuart, cat. # 141C. Use of the Revogene GBS DS assay for clinical specimen types other than those specified, has not been evaluated and performance characteristics are not established.
- Cervical, perianal, perirectal or perineal specimens are not acceptable, and a speculum should not be used for culture collection.
- Erroneous test results may occur from improper specimen collection, handling or storage, technical error, or sample mix-up. Careful compliance with the instructions of this insert, the Revogene Operator's Manual<sup>5</sup> and to established guidelines is necessary to avoid erroneous results.
- Contamination or false negative results may occur if a PIE cap is incorrectly closed.
- While there are no known strains/isolates of GBS lacking the *cfb* gene, the occurrence of such a strain would lead to an erroneous result (false negative) using the Revogene GBS DS assay.
- Mutations or polymorphisms in primer- or probe-binding regions may affect detection of GBS *cfb* gene variants, resulting in a false negative result with the Revogene GBS DS assay.
- Results from the Revogene GBS DS assay should be used as an adjunct to clinical observations and other information available to the physician.
- A negative result does not rule out the possibility of GBS colonization. False negative results may occur when the GBS concentration is below limit of detection of the assay.
- The test is not intended to differentiate carriers of GBS from those with streptococcal disease.
- Assay results may be affected by concurrent antimicrobial therapy as GBS DNA may continue to be detected.
- Moisturizing lotion (e.g. Aveeno<sup>®</sup>), whole blood or meconium may potentially interfere with the Revogene GBS DS assay when either of these substances is present in Sample Buffer at a concentration of > 6.50E-3% (w/v or v/v).
- Fungicide (e.g. Micatin<sup>®</sup>) or Amniotic fluid may potentially interfere with the Revogene GBS DS assay when either of these substances is present in sample buffer at a concentration of > 0.065% (w/v or v/v).
- Gastritis medication (e.g. Tums<sup>®</sup>) at concentration > 3.48E-6% (w/v), anti-diarrheal medications (e.g. Imodium<sup>®</sup>) at concentration > 1.50E-6% (w/v), anti-diarrheal medication (e.g. Pepto Bismol<sup>™</sup>) at concentration > 1.95E-3% (w/v) or laxative (e.g. Senokot<sup>®</sup>) at concentration > 3.19E-7% (w/v) may potentially interfere with the Revogene GBS DS assay when either of these substances is present in Sample Buffer.
- Body powder (Vagisil<sup>®</sup> deodorant powder) at the lower concentration tested, 6.50E-05 % (w/v), presented an inhibitory effect on the detection of GBS strains.

## EXPECTED VALUES

During clinical evaluation of the Revogene GBS DS assay, a total of 462 reportable results from specimens compliant at the reference method culture level were obtained. Overall, GBS prevalence rate was 26.4 % (112 of 462) with a 95%CI of 22.5 – 30.9 %.

## CLINICAL PERFORMANCE CHARACTERISTICS

A study was conducted at the Revogene facility to evaluate Revogene GBS DS assay clinical performance. A total of 462 *intrapartum* vaginal/rectal swabs were prospectively collected from pregnant women in two (2) Canadian sites. Only specimens that met the study inclusion criteria and did not meet any of the exclusion criteria were enrolled. Two (2) swabs per patient were sampled and tested. The first was tested with reference culture method while the second was discharged in SB and frozen at -70 C before testing with Revogene GBS DS assay on the Revogene. Of these 462 specimens, 422 specimens were tested with Revogene GBS DS assay at the Meridian Bioscience Canada, Inc. facility, the remaining 40 specimens were not available for testing.

The reference culture method was performed by one (1) of the Canadian sites. The reference culture method consisted in the transfer of a vaginal/rectal swab in LIM broth enrichment culture followed by a subculture on a non-selective blood agar plate. Colonies morphologically resembling GBS from blood agar plate were confirmed for their identity by Group B streptococcal grouping antisera and/or CAMP test.

The sensitivity, specificity, Positive Predictive Value (PPV) and Negative Predictive Value (NPV) values were calculated by comparing Revogene GBS DS assay results with the reference culture method. Discrepant analysis was performed by the site that has performed the reference culture method. A portion of specimens with discordant results between the Revogene GBS DS assay and the culture method were tested with a second Nucleic Acid Amplification Test (NAAT) method.

### Results

Overall performance characteristics of the Revogene GBS DS assay demonstrated 96.4 % sensitivity (107/111; 95%CI: 91.1 - 98.6 %) and 89.9 % specificity (277/308; 95%CI: 86.1 - 92.8 %). Performance results from this study are summarized in **Table 1**.

Of the 422 specimens tested with the Revogene GBS DS assay that were compliant at the specimen and PCR level, five (5) (1.2 %) were reported as Unresolved at initial testing. The unresolved rate after repeat testing was 0.7 % (3/422).

Out of the same 422 specimens tested with the Revogene GBS DS assay, one (1) (0.2 %) was initially reported as Indeterminate. No result remained Indeterminate upon repeat testing.

**Table 1.** Overall Performance Characteristics of the GBS DS Assay in Comparison to the Reference Method.

GBS DS assay		Reference Method		
		Positive	Negative	Total
		Positive	Negative	Total
Sensitivity		107	31 <sup>B</sup>	138
Specificity		4 <sup>A</sup>	277	281
Positive Predictive Value		77.5% (107/138; 95%CI: 69.9 - 83.7%)		
Negative Predictive Value		98.6% (277/281; 95%CI: 96.4 - 99.4%)		

<sup>A</sup> Two (2) specimens were tested using a second NAAT method; GBS DNA was detected in one (1) out of two (2) false negative specimens.

<sup>B</sup> Fifteen specimens were tested using a second NAAT method; GBS DNA was detected in 13 out 15 false positive specimens.

## ANALYTICAL PERFORMANCE CHARACTERISTICS

### ANALYTICAL SENSITIVITY

The analytical sensitivity (Limit of Detection or LoD) of the Revogene GBS DS assay was determined using simulated matrix, spiked with various concentrations of GBS bacterial suspensions. Three (3) strains of GBS (ATCC® 12403™, ATCC® 13813™ and ATCC® BAA-22™) were tested in 24 replicates per concentration. The LoD is defined as the lowest concentration at which 95% of all replicates tested positive.

In regard to the three (3) strains tested, the LoD of the GBS DS assay in simulated matrix ranged from 375 to 1 500 CFU/mL of sample buffer (SB). Results are summarized in **Table 2**. Equivalence of LoD was demonstrated between the clinical vaginal/rectal swab matrix previously tested negative for GBS and simulated matrix with strains ATCC® 12403™ and ATCC® 13813™.

**Table 2.** LoD of the Revogene GBS DS Assay

GBS strain (ATCC® number)	LoD in Simulated Matrix (CFU/mL of SB)
Serotype III (ATCC® 12403™)	750
Serotype III (ATCC® BAA-22™)	1 500
Non-hemolytic (ATCC® 13813™)	375

### INCLUSIVENESS

Inclusivity of the Revogene GBS DS assay was determined for 13 strains of GBS representing 11 known serotypes and one (1) non-hemolytic strain. Ten (10) strains were tested from a quantitated cell culture in presence of GBS-negative simulated matrix at a load of 1125 CFU/mL of SB, corresponding to three (3) times the LoD value of ATCC® 13813™ strain. Three (3) replicates per strain were tested using three (3) different GBS DS kit lots. Six (6) strains were detected with 100% positivity (in 3/3 replicates) by the Revogene GBD DS assay at a load of 1125 CFU/mL of SB. Two (2) strains (ATCC® 51487™ and ATCC® BAA-2669™) were detected with 100% positivity at a load of 1875 CFU/mL of SB. The strains ATCC® 27591™ and ATCC® 12973™ were detected with 100% positivity at a load of 2625 CFU/mL of SB and 15 000 CFU/mL of SB respectively. Strains tested are described in **Table 3**.

**Table 3.** GBS Strains Tested for Inclusivity with the Revogene GBS DS assay.

GBS strain	Serotype
ATCC® 12400™	Ia
ATCC® 51487™	Ib
ATCC® 27591™	Ic
ATCC® 12973™	II
ATCC® 12403™ <sup>1</sup>	III
ATCC® BAA-22™ <sup>1</sup>	III
ATCC® 49446™	IV
ATCC® BAA-611™	V
ATCC® BAA-2671™	VI
ATCC® BAA-2670™	VII
ATCC® BAA-2669™	VIII
ATCC® BAA-2668™	IX
ATCC® 13813™ <sup>1</sup>	Non-hemolytic

<sup>1</sup> Strains used for the LoD determination and therefore not retested

In addition, an *in silico* analysis was performed to assess inclusivity of the primers and probes of the Revogene GBS DS assay target regarding 36 additional GBS strains listed in the National Center for Biotechnology Information (NCBI) database of October 20<sup>th</sup>, 2016. The alignment results showed no mismatch with the selected sequences. The analysis predicted the detection of all these GBS strains.

## CROSS-REACTIVITY

The cross-reactivity of the Revogene GBS DS assay was assessed with high loads of organisms that are not targeted by the assay, or phylogenetically related to GBS, or present in the normal urogenital flora and intestinal tract. The study included 64 bacteria, four (4) yeasts, four (4) viruses, two (2) parasites and human DNA (**Table 4**). Bacteria and yeasts were tested at a load of  $\geq 10^6$  CFU/mL of SB. Nucleic acids from viruses, parasites and human DNA were tested at a load of  $\geq 10^5$  DNA or RNA copies/mL of SB. These organisms were tested using a quantitated cell cultures or nucleic acid solutions in presence of GBS-negative simulated matrix. Three (3) replicates per organism were tested using three (3) different GBS DS kit lots.

Under the conditions of the study, none of the organisms or nucleic acids tested were detected by the Revogene GBS DS assay.

**Table 4.** List of Organisms Tested for Cross-Reactivity with the Revogene GBS DS Assay

Bacteria	
<i>Acinetobacter baumannii</i>	<i>Mycoplasma genitalium</i> (gDNA)
<i>Aerococcus viridans</i>	<i>Mycoplasma hominis</i> (gDNA)
<i>Aeromonas hydrophila</i>	<i>Neisseria gonorrhoeae</i>
<i>Bacillus cereus</i>	<i>Peptostreptococcus anaerobius</i>
<i>Bacillus subtilis</i>	<i>Porphyromonas asaccharolytica</i>
<i>Bacteroides fragilis</i>	<i>Prevotella melaninogenica</i>
<i>Bifidobacterium adolescentis</i>	<i>Propionibacterium acnes</i>
<i>Bifidobacterium breve</i>	<i>Proteus mirabilis</i>
<i>Brevibacterium linens</i>	<i>Pseudomonas aeruginosa</i>
<i>Campylobacter jejuni</i>	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Dublin</i>
<i>Chlamydia trachomatis</i> (gDNA)	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Minneapolis</i>
<i>Citrobacter freundii</i>	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>typhimurium</i>
<i>Clostridium difficile</i>	<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar <i>Newport</i>
<i>Clostridium perfringens</i>	<i>Serratia marcescens</i>
<i>Corynebacterium diphtheriae</i>	<i>Shigella sonnei</i>
<i>Enterobacter aerogenes</i>	<i>Staphylococcus aureus</i>
<i>Enterobacter cloacae</i>	<i>Staphylococcus aureus</i> (Cowan)
<i>Enterococcus avium</i>	<i>Staphylococcus epidermidis</i>
<i>Enterococcus dispar</i>	<i>Staphylococcus saprophyticus</i>
<i>Enterococcus durans</i>	<i>Streptococcus anginosus</i>
<i>Enterococcus faecalis</i>	<i>Streptococcus bovis</i>
<i>Enterococcus faecium</i>	<i>Streptococcus dysgalactiae</i> subsp. <i>dysgalactiae</i>
<i>Escherichia coli</i>	<i>Streptococcus dysgalactiae</i> equisimilis
<i>Escherichia fergusonii</i>	<i>Streptococcus intermedius</i>
<i>Gardnerella vaginalis</i>	<i>Streptococcus oralis</i>
<i>Klebsiella oxytoca</i>	<i>Streptococcus pneumoniae</i>
<i>Klebsiella pneumoniae</i>	<i>Streptococcus pyogenes</i>
<i>Lactobacillus acidophilus</i>	<i>Streptococcus salivarius</i>
<i>Lactobacillus brevis</i>	<i>Streptococcus sanguinis</i>
<i>Lactobacillus casei</i>	<i>Streptococcus uberis</i>
<i>Lactobacillus delbreuckii</i>	<i>Ureaplasma urealyticum</i> (gDNA)
<i>Lactobacillus jensenii</i>	<i>Yersinia enterocolitica</i>
Yeasts	
<i>Candida albicans</i>	<i>Candida parapsilosis</i>
<i>Candida glabrata</i>	<i>Candida tropicalis</i>
Viruses	
<i>HerpesSimplexVirus-1</i> (gDNA)	<i>Norovirus GII</i> (synthetic RNA)
<i>HerpesSimplexVirus-2</i> (gDNA)	Human Papillomavirus (gDNA)
Parasites	
<i>Blastocystis hominis</i> (gDNA)	<i>Trichomonas vaginalis</i> (gDNA)
Human DNA	
Human DNA (gDNA)	-

In addition, an *in silico* analysis was performed to assess the cross-reactivity of the primers and probes of the Revogene GBS DS assay towards 62 additional organisms potentially found in urogenital and intestinal tract (**Table 5**). The analysis was performed on all sequences of these 62 organisms contained in the National Center for Biotechnology Information (NCBI) database on January 24<sup>th</sup>, 2017. The analysis suggested that none of these organisms should be reactive with the Revogene GBS DS assay.

**Table 5.** List of Organisms Verified by *in silico* Analysis for Cross-Reactivity with the Primers and Probes of the Revogene GBS DS Assay

Bacteria	
<i>Abiotrophia defectiva</i>	<i>Listeria monocytogenes</i>
<i>Achromobacter xylosoxidans</i>	<i>Micrococcus luteus</i>
<i>Alcaligenes faecalis</i>	<i>Mobiluncus mulieris</i>
<i>Campylobacter fetus</i>	<i>Moraxella catarrhalis</i>
<i>Chromobacterium violaceum</i>	<i>Moraxella lacunata</i>
<i>Clostridium bifermentans</i>	<i>Moraxella osloensis</i>
<i>Clostridium butyricum</i>	<i>Morganella morganii</i>
<i>Clostridium novyi</i>	<i>Mycobacterium smegmatis</i>
<i>Clostridium septicum</i>	<i>Neisseria flava</i>
<i>Clostridium sordelli</i>	<i>Neisseria flavescens</i>
<i>Clostridium sporogenes</i>	<i>Neisseria lactamica</i>
<i>Corynebacterium urealyticum</i>	<i>Neisseria meningitidis</i>
<i>Corynebacterium xerosis</i>	<i>Neisseria perflava</i>
<i>Deinococcus radiodurans</i>	<i>Plesiomonas shigelloides</i>
<i>Dexia gummosa</i>	<i>Proteus vulgaris</i>
<i>Eikenella corrodens</i>	<i>Providencia stuartii</i>
<i>Enterococcus gallinarum</i>	<i>Pseudomonas fluorescens</i>
<i>Gemella haemolysans</i>	<i>Pseudomonas putida</i>
<i>Haemophilus ducreyi</i>	<i>Rhodospirillum rubrum</i>
<i>Haemophilus influenzae</i>	<i>Serratia liquefaciens</i>
<i>Haemophilus influenzae type B</i>	<i>Shigella boydii</i>
<i>Helicobacter pylori</i>	<i>Shigella flexneri</i>
<i>Kingella denitrificans</i>	<i>Staphylococcus simulans</i>
<i>Kingella kingae</i>	<i>Stenotrophomonas maltophilia</i>
<i>Legionella pneumophila</i>	<i>Streptococcus gordonii</i>
<i>Streptococcus mitis</i>	<i>Vibrio parahaemolyticus</i>
<i>Streptococcus mutans</i>	-
Mold/Yeast	
<i>Cryptococcus neoformans</i>	<i>Saccharomyces cerevisiae</i>
Viruses	
<i>Adenovirus 40</i>	<i>Echovirus</i>
<i>Adenovirus 41</i>	<i>Epstein-Barr Virus</i>
<i>BK virus</i>	<i>Rotavirus</i>
<i>Coxsackievirus</i>	-

**INTERFERING ORGANISMS**

The potentially inhibitory effect of 29 non-targeted organisms that may be present in urogenital flora and intestinal tract was assessed using organisms selected from the cross-reactivity study (**Table 4**). Selection of organisms was based on their prevalence in the vaginal/rectal flora. Groups of one (1) to five (5) organisms were tested in triplicate with either 1125 CFU/mL of SB containing ATCC® 13813™ GBS strain, 4500 CFU/mL of SB containing GBS ATCC® BAA-22™ strain or GBS-negative samples in presence of GBS-negative simulated matrix to assess their potential interference on detection of GBS or PrC. Each organism within a group was diluted to reach a load of  $\geq 10^6$  CFU/mL or copies/mL of SB for bacteria and yeast, and  $\geq 10^5$  copies/mL of SB for virus and parasite. The 29 organisms included in the study are presented in **Table 6**.

**Table 6.** List of Organisms Tested for Interference with the Revogene GBS DS assay.

Group 1	
<i>Acinetobacter baumannii</i>	<i>Enterococcus faecalis</i>
<i>Bacillus subtilis</i>	<i>Enterococcus faecium</i>
<i>Bacteroides fragilis</i>	
Group 2	
<i>Campylobacter jejuni</i>	<i>Escherichia coli</i>
<i>Clostridium difficile</i>	<i>Gardnerella vaginalis</i>
<i>Corynebacterium genitalium</i>	
Group 3	
<i>Klebsiella oxytoca</i>	<i>Peptostreptococcus anaerobius</i>
<i>Lactobacillus acidophilus</i>	<i>Propionibacterium acnes</i>
<i>Neisseria gonorrhoeae</i>	
Group 4	
<i>Pseudomonas aeruginosa</i>	<i>Shigella sonnei</i>
<i>Salmonella enterica subsp. enterica serovar Dublin</i>	<i>Staphylococcus aureus</i>
<i>Serratia marcescens</i>	
Group 5	
<i>Staphylococcus epidermidis</i>	<i>Yersinia enterocolitica</i>
<i>Streptococcus pyogenes</i>	<i>Candida albicans</i>
Group 6	
<i>Chlamydia trachomatis</i>	<i>Human Papillomavirus (gDNA)</i>
<i>Herpes Simplex Virus-1 (gDNA)</i>	<i>Trichomonas vaginalis</i>
Group 7	
<i>Mycoplasma genitalium (gDNA)</i>	

None of the 29 organisms present at  $\geq 10^6$  CFU/mL or copies/mL of SB for bacteria and yeast and  $\geq 10^5$  copies/mL of SB for viruses and parasites interfered with detection of PrC and with the detection of the two (2) GBS strains (ATCC® 13813™ and ATCC® BAA-22™).

**INTERFERING SUBSTANCES**

The potentially inhibitory effect of 22 exogenous and eight (8) endogenous substances that may be present in urogenital flora and intestinal tract was assessed using GBS-negative samples and GBS strains ATCC® 13813™ and ATCC® BAA-22™ containing samples at three (3) times their respective LoD (1 125 CFU/mL of SB or 4 500 CFU/mL of SB respectively) in presence of GBS-negative simulated matrix. Substances were tested at a physiologically relevant concentration that could be found in a vaginal/rectal specimen or at a higher concentration. The results for the 30 substances are presented in **Table 7**.

Results demonstrated that fungicide (e.g. Micatin®), body powder (e.g. Vagisil ® deodorant powder), moisturizing lotion (e.g. Aveeno ® moisturizing lotion), gastritis medications (e.g. Tums®), anti-diarrheal medication (Imodium ®, Pepto Bismol™), laxatives (Senokot ®), whole blood, amniotic fluid and meconium showed a potentially inhibitory effect on the detection of PrC or GBS strains when any of these substances was present in SBT at the concentration above those presented in **Table 7**. When tested at concentration presented in **Table 7**, all these substances, except the body powder, showed no reportable interference with the Revogene GBS DS assay. The body powder at the lower concentration tested, 6.50E-05 % (w/v), presented an inhibitory effect on the detection of GBS strains.

**Table 7.** List of Exogenous and Endogenous Substances Tested with the Revogene GBS DS Assay.

Exogenous Substances		Concentration in SB without interference observed <sup>1</sup>
Substance (commercial name)		
Fungicide (Micatin®)		0.065 % w/v*
Hemorrhoid cooling gel (Preparation H ®)		6.50 % w/v
Lubricating gel (K-Y ® Personal lubricant)		0.65 % w/v
Body powder (Vagisil ® deodorant powder)		Concentration not determined*
Moisturizing lotion (Aveeno ® moisturizing lotion)		6.50E-03 % w/v*
Body oil (Neutrogena ® Body oil)		6.50 % v/v
Deodorant spray (Summer's Eve™ Spray)		6.50 % v/v
Enemas (Life BRAND™ Heavy Mineral Oil USP)		6.50 % v/v
Antimicrobials (Canesten ®)		6.50 % w/v
Enemas (Pentasa)		1.43E-03% w/v
Radiology Oral Compounds (Barium Sulfate)		7.50E-04 % w/v
Gastritis Medications (Nexium)		7.50E-04 % w/v
Antimicrobials (Flagyl)		4.64E-04 % w/v
Non-Steroidal Anti-Inflammatory (Aleve ®)		2.04E-03 % w/v
Antimicrobials (DIFLUCAN ® One)		5.57E-05 % w/v
Gastritis Medications (Tums ®)		3.48E-06 % w/v*
Anti-Diarrheal Medication (Imodium ®)		1.50E-06 % w/v*
Anti-Diarrheal Medication (Pepto Bismol™)		1.95E-03 % w/v*
Laxatives (Senokot ®)		3.19E-07 % w/v*
Spermicidal (Trojan ® with spermicidal Lubricant Condoms)		0.65 % v/v
Moist Towelettes (Equate™ Flushable Moist Wipes)		6.50 % v/v
Moist Towelettes (Wet Ones ®)		0.65 % v/v
Endogenous substances		Concentration in SB without interference observed <sup>1</sup>
Substance		
Whole blood		6.50E-03 % v/v*
Leukocytes		1.00E+06 cells/mL
Amniotic fluid		0.065 % v/v*
Seminal fluid		6.50 % v/v
Urine		6.50 % v/v
Feces		0.65 % v/v
Meconium		6.50E-03 % w/v*
Human DNA		1.55E+03 ng/mL

<sup>1</sup> % W/V: Weight/Volume (mg/mL); % V/V: Volume/Volume (µL/mL);

\*Concentrations above these numbers may cause interference (see limitations above).

**PRECISION/REPRODUCIBILITY**

Between-laboratory reproducibility study was performed in three (3) laboratories at Meridian Bioscience Canada, Inc. facility, by two (2) operators per laboratory, over five (5) distinct days using one (1) GBS DS kit lot (see **Table 8**). Results from laboratory 1 over five (5) days were used for the precision study (see **Table 8**).

Between-lot reproducibility study was performed on one (1) laboratory by two (2) operators, over five (5) days using three (3) GBS DS kit lots (see **Table 9**).

For the reproducibility and precision studies, a total of 60 replicates of negative samples and 90 replicates of each category of positive samples, all prepared in a GBS-negative simulated matrix, were tested. GBS strain ATCC® 13813™ (non-hemolytic) was used for positive samples.

Samples categories were described as follow:

- Low positive (LP): 1.95 x LoD or 731 CFU/mL of SB
- Moderate positive (MP): 3 x LoD or 1125 CFU/mL of SB
- True negative (TN): samples without GBS target

For the between-laboratory reproducibility, the overall percent agreement was 100% for TN, 98.8% for MP and 96.7% for LP of GBS ATCC® 13813™ strain (**Table 8**). For the precision study, the overall percent agreement was 100% for TN, 96.7% for MP and 96.7% for LP of GBS ATCC® 13813™ strain (**Table 8, Laboratory 1**). For the between-lot reproducibility, the overall percent agreement was 100% for TN, 97.8% for MP and 92.2% for LP of GBS ATCC® 13813™ strain (**Table 9**). Overall mean Cycle threshold (Ct) values with variance component (%CV) are shown in **Tables 8** and **9**. Summary of the overall reproducibility study results (between-laboratory and between-lot combined) are shown in **Table 10**.

**Table 8.** Precision and Between-Laboratory Reproducibility Study Results Using One (1) GBS DS Assay Kit Lot

Laboratory 1, instrument 1, lot 1				
Panel ID	Results/Total	% Agreement <sup>1</sup>	Ct Mean <sup>2</sup>	Ct %CV
LP	29/30	96.7%	38.2	5.6
MP	29/30	96.7%	37.3	4.2
TN	20/20	100%	30.2	3.1
Laboratory 2, instrument 2, lot 1				
Panel ID	Results/Total	% Agreement <sup>1</sup>	Ct Mean <sup>2</sup>	Ct %CV
LP	28/30	93.3%	36.7	2.9
MP	30/30	100%	36.9	4.7
TN	20/20	100%	32.5	3.1
Laboratory 3, instrument 3, lot 1				
Panel ID	Results/Total	% Agreement <sup>1</sup>	Ct Mean <sup>2</sup>	Ct %CV
LP	30/30	100%	37.9	5.5
MP	30/30	100%	36.4	3.7
TN	20/20	100%	30.1	1.8
Overall				
Panel ID	Results/Total	% Agreement <sup>1</sup>	Ct Mean <sup>2</sup>	Ct %CV
LP	87/90	96.7%	37.6	5.1
MP	89/90	98.8%	36.9	4.3
TN	60/60	100%	30.9	4.6

<sup>1</sup> For the TN category, percent agreement was calculated for negative results.<sup>2</sup> For the LP and MP categories, Ct values reported are for the GBS target. For the TN category, Ct values reported are for the PrC.

**Table 9.** Between-Lot Reproducibility Study Results Using Three (3) GBS DS Assay Kit Lots

Laboratory 1, instrument 1, lot 1				
Panel ID	Results/Total	% Agreement <sup>1</sup>	Ct Mean <sup>2</sup>	Ct %CV
LP	29/30	96.7%	38.2	5.6
MP	29/30	96.7%	37.3	4.2
TN	20/20	100%	30.2	3.1
Laboratory 1, instrument 1, lot 2				
Panel ID	Results/Total	% Agreement <sup>1</sup>	Ct Mean <sup>2</sup>	Ct %CV
LP	26/30	86.7%	38.6	6.1
MP	29/30	96.7%	37.0	4.8
TN	20/20	100%	29.9	2.6
Laboratory 1, instrument 1, lot 3				
Panel ID	Results/Total	% Agreement <sup>1</sup>	Ct Mean <sup>2</sup>	Ct %CV
LP	28/30	93.3%	37.1	3.5
MP	30/30	100%	38.0	5.7
TN	20/20	100%	29.9	2.6
Overall (all lot)				
Panel ID	Results/Total	% Agreement <sup>1</sup>	Ct Mean <sup>2</sup>	Ct %CV
LP	83/90	92.2%	37.9	5.4
MP	88/90	97.8%	37.4	5.0
TN	60/60	100%	30.0	2.8

<sup>1</sup> For the TN category, percent agreement was calculated for negative results.<sup>2</sup> For the LP and MP categories, Ct values reported are for the GBS target. For the TN category, Ct values reported are for the PrC.**Table 10.** Summary of the Reproducibility Study Results (Between-Laboratory and Between-Lot combined)

All laboratories and lots combined				
Panel ID	Results/Total	% Agreement <sup>1</sup>	Ct Mean <sup>2</sup>	Ct %CV
LP	141/150	94.0%	37.7	5.2
MP	148/150	98.7%	37.1	4.8
TN	100/100	100%	30.5	4.2

<sup>1</sup> For the TN category, percent agreement was calculated for negative results.<sup>2</sup> For the LP and MP categories, Ct values reported are for the GBS target. For the TN category, Ct values reported are for the PrC.**E-LABELING**

Documentation related to this product can be accessed online at [www.meridianbioscience.com/pi](http://www.meridianbioscience.com/pi). Additionally, paper copies are available upon request by contacting your local distributor or via the phone number listed on the kit box.

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The screenshot shows a software interface for Revogene. At the top, it says "revogene" and the date "08/11/2021 13:50". On the right, there's a user icon labeled "admin". Below is a table with columns: Sample ID, Patient ID, Assay, and Result. The table contains 8 rows, numbered 1 to 8. Row 3 is highlighted with a teal background and has a checkmark in the first column. A callout box points to this row with the text "Results GBS Positive". The bottom of the screen has navigation buttons: Home, Started By: admin, Complete At: 14:10, Report, Run Report, Export, Edit, and Abort.

	Sample ID	Patient ID	Assay	Result
1	1300006001		GBS DS	+
2	1300006002		GBS DS	-
3 ✓	1235227901		GBS DS	+
4	1235227902		GBS DS	-
5	1235227801		GBS DS	-
6	Empty			
7	Empty			
8	Empty			

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This screenshot shows the Revogene software after the run has completed. It displays the results for the last run. A callout box highlights the result for sample 1300006001, which is marked with a checkmark and shows a positive result for GBS. The bottom of the screen has buttons for Home, Last Run, Search, Report, Run Report, Export, and Export Data.

	Sample ID	Patient ID	Started	Assay	Result
	1235227901		08/11/2021 12:58	GBS DS	+
✓	1300006001		08/11/2021 12:58	GBS DS	+
	1235227801		08/11/2021 12:58	GBS DS	-
	1295227902		08/11/2021 12:58	GBS DS	-
	1300006002		08/11/2021 12:58	GBS DS	-

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**REFERENCES**

1. MMWR Morbidity and Mortality Weekly Report, Recommendations and Reports, Prevention of Perinatal Group B Streptococcal Disease: Revised Guidelines from CDC, November 19, 2010 / Vol. 59 / No. RR-10; 1-36
2. A Population-Based Comparison of Strategies to Prevent Early-Onset Group B Streptococcal Disease in Neonates. N Engl J Med 2002; 347:233-239
3. Biosafety in Microbiological and Biomedical Laboratories, 5<sup>th</sup> ed. U.S Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, National Institute of Health, HHS Publication No. (CDC) 21-1112, Revised December 2009
4. M29-A4 Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline-Fourth Edition, Clinical And Laboratory Standards Institute (CLSI), May 2014
5. SN134822 Revogene® Operator's Manual
6. Diagnostic Accuracy of a Rapid Real-Time Polymerase Chain Reaction Assay for Universal Intrapartum Group B Streptococcus Screening. Clin Infect Dis 2009; 49:417–23
7. Intrapartum GBS screening and antibiotic prophylaxis: a European consensus conference. J Matern Fetal Neonatal Med 2015; 28(7): 766–782

SN134753

REV. 05/23

 <b>Manufactured By</b>	<p><b>Meridian Bioscience, Inc.</b>  3471 River Hills Drive  Cincinnati, OHIO - 45244 USA  <a href="http://www.meridianbioscience.com">www.meridianbioscience.com</a></p> <p><b>Contacts:</b>  Main Telephone (+1) 513.271.3700  Customer Service/Orders 800.543.1980  Technical Support Center 800.343.3858  Information Fax: 513.272.5432  Ordering Fax: 513.271.0124  E-mail: <a href="mailto:info@meridianbioscience.com">info@meridianbioscience.com</a></p>
	<p><b>Meridian Bioscience Europe, SRL</b>  Via Dell'Industria 7, 20035 Villa Cortese  (Milano) ITALY  <a href="http://www.meridianbioscience.com">www.meridianbioscience.com</a></p> <p><b>Contacts:</b>  Main Telephone (+39) 0331.433636  E-mail: <a href="mailto:info@meridianbioscience.eu">info@meridianbioscience.eu</a>  Technical Support: <a href="mailto:MBE-TechService@meridianbioscience.eu">MBE-TechService@meridianbioscience.eu</a>  Customer Service/Orders:  <ul style="list-style-type: none"> <li>• For Italian Customers: <a href="mailto:ordini@meridianbioscience.com">ordini@meridianbioscience.com</a></li> <li>• For Distributors / International Customers: <a href="mailto:Export.CustomerService@meridianbioscience.eu">Export.CustomerService@meridianbioscience.eu</a></li> </ul> </p>
<b>AUSTRALIAN SPONSOR</b>	<p><b>Energo Australia</b>  Level 20, Tower II  Darling Park  201 Sussex Street  Sydney, NSW 2000  Australia</p>
	<p><b>MedEnvoy Switzerland</b>  Gotthardstrasse 28  6302 Zug  Switzerland</p>

## INTERNATIONAL SYMBOLS USAGE

You may see one or more of these symbols on the labeling/packaging of this product

**Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simblos, Zeichenerklärung)**

	Use-by date / Data di scadenza / Date de péremption / Fecha de caducidad / Verwendbar bis
<b>LOT</b>	Batch code / Codice di lotto / Code de lot / Código de lote / Chargennummer
<b>IVD</b>	In vitro diagnostic medical device / Dispositivo medico diagnostico <i>in vitro</i> / Instrument de test diagnostique <i>in vitro</i> / Producto sanitario para diagnóstico <i>In vitro</i> / <i>In-vitro</i> -Diagnostikum
<b>CE</b>	CE Mark / Marcatura CE / Symbole CE / Marcado CE / CE-Kennzeichen
<b>REF</b>	Catalog number / Numero di catalogo / Référence catalogue / Referencia / Bestellnummer
	Consult instructions for use / Consultare le istruzioni per l'uso / Consulter le mode d'emploi / Consultar las instrucciones de uso / Gebrauchsanweisung beachten
	Manufacturer / Produttore / Fabricant / Fabricante / Hersteller
	Contains sufficient for <n> tests / Contiene una quantità sufficiente per <n> test / Contient le matériel suffisant pour <n> tests / Contiene la cantidad suficiente para <n> ensayos / Inhalt ausreichend für <n> Tests
	Temperature limit / Limite di temperatura / Limite de température / Límite de temperatura / Temperaturgrenze
<b>EC REP</b>	Authorized representative in the European Community / Rappresentante autorizzato nella Comunità Europea / Représentant agréé dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter EU-Representant
	Do not reuse / Non riutilizzare / Ne pas réutiliser / No reutilizar / Nicht wiederverwenden
	Keep dry / Conservare all'asciutto / Conserver au sec / Mantener seco / Vor Feuchtigkeit schützen
	Contains # pouches: 1 Disposable Transfer Tool (DTT), 1 Sample Buffer Tube (SBT), 1 PIE / Contiene # buste: 1 PIE, 1 Sample Buffer Tube (SBT), 1 Disposable Transfer Tool (DTT) / Contient # sachets: Une cartouche PIE, Un tube échantillon, Un outil de transfert jetable (OTJ) / Incluye # bolsitas: 1 PIE, 1 Sample Buffer Tube (SBT), 1 Disposable Transfer Tool (DTT) / Enthält # Beutel: 1 PIE, 1 Sample Buffer Tube (SBT), 1 Disposable Transfer Tool (DTT)
	Humidity Limitation / Limitazione dell'umidità / Limite d'humidité / Limitación de humedad / Feuchtigkeitsbegrenzung
<b>Rx Only</b>	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Ne pas utiliser si l'emballage est endommagé / No utilizar si el envase está dañado / Nicht verwenden, wenn die Verpackung beschädigt ist
	Keep away from sunlight / Tenere lontano dalla luce del sole / Tenir à l'écart de la lumière du soleil / Mantener alejado de la luz solar / Vor Sonnenlicht schützen
	Contains 24 Disposable Transfer Loops (DTL) / Contiene 24 loop di trasferimento monouso (DTL) / Contient 24 boucles de transfert jetables (DTL) / Contiene 24 bucles de transferencia desechables (DTL) / Enthält 24 Einwegübertragungsschleifen (DTL)
<b>DTT</b>	Disposable Transfer Tool (DTT) / Disposable Transfer Tool 9DTT / Un outil de transfert jetable (OTJ) / Disposable Transfer Tool (DTT) / Dispaoable Transfer Tool (DTT)
<b>PIE</b>	Revogene Test Device / Dispositivo Test Revogene / Dispositif de test Revogene / Dispositivo para la Prueba Revogene / Revogene-Analysegefäß
<b>EUA</b>	For Emergency Use Authorization only / Solo per l'autorizzazione per l'uso di emergenza / pour autorisation d'utilisation d'urgence uniquement / para autorización de uso de emergencia solamente / nur für Notfallverwendungsautorisierung
<b>LOOP</b>	Disposable Transfer Loop / Loop di trasferimento monouso / Boucles de transfert jetables / Bucle de transferencia desechables / Einwegübertragungsschleifen
<b>SBT</b>	Sample Buffer Tube / Sample Buffer Tube / Un tube échantillon / Smample Buffer Tube / Sample Buffer Tube
	Revogene Test Device / Dispositivo Test Revogene / Dispositif de test Revogene / Dispositivo para la Prueba Revogene / Revogene-Analysegefäß
<b>CH REP</b>	Swiss Authorized Representative / Mandatario svizzero / Mandataire Suisse / Representante Autorizado Suizo / Schweizer Bevollmächtigter

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.

TaqMan is a registered trademark of Roche Molecular Systems, Inc.  
ATCC is a registered trademark of the American Type Culture Collection.  
Revogene and associated logos are trademarks of Meridian Bioscience, Inc.  
© 2020-12 Meridian Bioscience, Inc.  
Made in Canada.

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TaqMan è un marchio registrato di Roche Molecular Systems, Inc.
Revogene e i relativi loghi sono marchi di Meridian Bioscience, Inc.
TaqMan est une marque déposée de Roche Molecular Systems, Inc.
Revogene et les logos associés sont des marques de commerce de Meridian Bioscience, Inc.
TaqMan es una marca registrada de Roche Molecular Systems, Inc.
Revogene y los logotipos asociados son marcas comerciales de Meridian Bioscience, Inc.