INTENDED USE

The Revogene® Strep A assay, performed on the Revogene instrument, is an automated, qualitative *in vitro* diagnostic test that utilizes real-time polymerase chain reaction (PCR) for the direct detection of *Streptococcus pyogenes* (Group A β -hemolytic *Streptococcus*) nucleic acids from throat swab specimens obtained from patients with signs and symptoms of pharyngitis. The Revogene Strep A assay is intended for use as an aid in the diagnosis of Group A *Streptococcus* infection.

SUMMARY AND EXPLANATION

Streptococcus pyogenes, or Group A *Streptococcus* (GAS), is a facultative anaerobe, Gram-positive coccus which grows in chains and causes numerous infections in humans including pharyngitis, tonsillitis, scarlet fever, cellulitis, erysipelas, rheumatic fever, post-streptococcal glomerulonephritis, necrotizing fasciitis, myonecrosis and lymphangitis.

Group A *Streptococcus* (GAS) remains the most common bacterial etiology, and accounts for 15 to 30% of cases of acute pharyngitis in children and 5 to 20% in adults.¹ Because the signs and symptoms of GAS pharyngitis overlap extensively with other infectious causes (sore throat, fever with sudden onset), it makes difficult doing a diagnosis based solely on clinical findings. However, an accurate diagnosis is necessary to promptly and adequately treat the disease through antibiotherapy and avoid possible serious complications such as rheumatic fever.

Conventional methods for the detection of GAS infection include the use of Rapid Antigen Detection Tests (RADTs) or culturing pharyngeal swab samples followed by Lancefield group differentiation with latex agglutination. Although highly specific, the sensitivity of RADTs can be as low as 31 to 50%,² prompting the need for back-up cultures.^{3,4} On the other hand, culture can take up to 48 h, delaying appropriate antimicrobial treatment.⁵ The Revogene Strep A assay can provide results from one up to eight samples as early as 42 minutes for positive specimens and in approximately 70 minutes for negative specimens, without the need for culture confirmation of negative Revogene Strep A assay results. The assay minimizes operator intervention from the loading of the sample into the single use microfluidic cartridge (PIE) to final result.

PRINCIPLE OF THE PROCEDURE

The Revogene automates sample homogenization and dilution, cell lysis, DNA amplification and detection of the amplified PCR products. User intervention is only required to discharge the patient specimen into the Sample Buffer Tube (SBT), transfer the sample from SBT into the PIE and insert the PIEs into the Revogene carousel.

Each PIE is a fully integrated closed device in which a sample is dispensed and processed through different microfluidic chambers and channels which allow for the sample processing (i.e. sample homogenization, sample dilution, cell lysis and DNA extraction) and subsequent real-time PCR steps (Figure 1). The liquid from a single sample is transferred by centrifugation from one chamber to the next in sequence. All reagents specific for the PCR reaction are incorporated and dried within the PCR well. A Process Control (PrC) is incorporated into each PIE to verify sample processing and amplification steps including the verification of potential inhibitor substances as well as microfluidic, instrument or reagent failure. The amplified products are detected in real-time using target-specific TaqMan® chemistry-based probes. No operator intervention is necessary once a PIE is loaded into the Revogene.


Figure 1. Top View of PIE.
1: Sample Loading Chamber, 2: Overflow Chamber,
3: Homogenization Chamber Containing PrC, 4: Dilution/Lysis Chamber,
5: Three (3) PCR Wells (A to C from left to right) and one (1) Waste Chamber (at the right end).

The Revogene can process from one (1) up to eight (8) samples simultaneously in the same run. The carousel must contain eight (8) PIEs to maintain thermodynamic balance within the run. During the run and at run completion, the results are computed by the system from measured fluorescent signals and embedded calculation algorithms. Results that are displayed on the touchscreen may be printed, transferred and stored by the user using the USB port or the connectivity option.

An Early Positive Result Outcome (E-PRO) feature provides positive result if the signal from the target DNA reaches a predetermined threshold before the full PCR cycles have been completed. With E-PRO, a positive result could be obtained in approximately 42 minutes for strongly positive samples. For negative samples, the time to result is approximately 70 minutes.

REAGENTS AND MATERIALS PROVIDED

Revogene Strep A kit containing sufficient reagents and materials to process 24 samples. The kit contains the following materials:

1. **24 Strep A Sample Buffer Tubes (SBT):** Tube containing Tris-HCl pH 8.0/EDTA Na₂ (TE 1X) buffered solution as a dilution and preservation buffer for sample.
2. **24 Disposable Transfer Tools (DTT):** Plastic pipette with minimal and maximal volume marks for transferring the sample from the SBT to the PIE.
3. **24 individual pouches containing one (1) Strep A Microfluidic cartridge (PIE):** Integrated device, composed of dried reagents allowing sample preparation and real-time PCR steps for PrC DNA, GAS DNA amplification and detection. Each PIE contains PrC, PrC: specific primers and Taqman® chemistry-based probe, GAS gene-specific primers and Taqman® chemistry-based probes, dNTPs, buffer and DNA polymerase.

MATERIALS/EQUIPMENTS REQUIRED BUT NOT PROVIDED

1. Revogene® (cat# 610210)
2. Disposable powderless gloves
3. Rayon or polyester swabs with Liquid Stuart or Liquid Amies transport media compatible with throat specimen collection (Copan cat. # 139C, 138C, 140C and 141C; BD cat # 220146 and 220149).
4. Vortex mixer with a maximal speed of at least 3200 rpm
5. Gauzes
6. MOCK PIEs (cat# 610208; optional)
7. Sample Rack (cat# 132539; optional)

WARNING AND PRECAUTIONS

1. The Revogene Strep A assay can only be used on the Revogene.
2. Do not use the kit if the label that seals the outer box is broken upon arrival.
3. Do not use Strep A PIEs if the protective pouches are open or broken upon arrival.
4. Do not interchange DTT, SBT, and PIE between kit lots.
5. Each single-use DTT and Strep A PIE are used to process one (1) sample. Do not reuse DTT or PIE.
6. Do not use swab with gel's transport media.
7. Always handle samples as if they are infectious and in accordance with Good Laboratory Practices such as those described in Biosafety in Microbiological and Biomedical Laboratories⁶ and in CLSI Document M29-A4⁷.
8. Wear disposable powderless gloves while handling samples and thoroughly wash hands afterwards.
9. The Strep A PIE contains dried reagents. The protective pouch should not be opened until ready to perform the test.
10. Dispose of unused materials, reagents and waste in accordance with country, federal, provincial, state and local regulations.
11. Do not open or break apart the PIE after use to avoid contamination with amplification products and/or infectious particles.
12. Do not use a PIE that has been dropped, shaken or inverted after the sample has been loaded as this may cause invalid results.
13. Do not use a kit that has passed its stated expiration date.
14. Do not refrigerate the loaded PIE.
15. Each run must be performed with eight (8) PIEs in the Revogene carousel to maintain thermodynamic and mechanical balance within the run. Place MOCK PIEs or new assay PIEs loaded with Sample Buffer in empty positions if less than eight samples are tested.

HAZARD AND PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

STORAGE AND STABILITY

1. Throat swab specimens should be stored between 2 C and 25 C during transport.
2. Throat swab specimens can be stored at 25 C for up to two (2) days, or at 2-8 C for up to seven (7) days.
3. Inoculated SBT can be stored at 25 C for up to three (3) days, or at 2-8 C for up to seven (7) days.
4. Store the Revogene Strep A kit at 2-25 C. The expiration date is indicated on the kit box label.
5. Do not open a pouch until ready to perform testing. Use the PIE within one (1) hour after opening the pouch.

INSTRUCTION FOR USE

SPECIMEN COLLECTION AND TRANSPORT

Specimen type: Throat swab from patients with signs and symptoms of pharyngitis collected using rayon or polyester swab with Liquid Stuart or Liquid Amies transport media.

1. Swab the posterior pharynx, tonsils and other inflamed areas following standard clinical methods and secure the swab in its transport sheath.
2. Label the transport sheath with the specimen or patient identification (ID) and send to the laboratory for testing (refer to **Storage and Stability** section).

SBT PREPARATION

NOTE 1 : Make sure only one (1) SBT is open at any given moment and remove only one (1) swab from its transport sheath at a time.

Step 1

- a. Remove the appropriate number of pouches from the box.

Step 2

- a. For each specimen to be tested, obtain one SBT from the kit box.

NOTE : Unseal the pouch only when ready to load the PIE. Start the test within one (1) hour after opening the pouch containing the PIE.

Step 3

- a. Identify (or label) the SBT with the appropriate specimen identification without obscuring or writing over the barcode.
- b. Place the SBT on the Sample Rack, if used.

Step 4

- a. Unscrew the SBT cap.
- b. Remove the swab from its transport sheath and transfer it in the SBT.
- c. Break the exceeding part of the swab stem (**Figure 2**):
 - Hold the stem with a gauze.
 - Lift the stem a few millimeters (mm) from the bottom of the tube.
 - Bend the stem against the edge of the tube to break it.
 - Alternatively, use clean scissors to cut the stem. Disinfect scissors between each sample.
 - Dispose of the swab sheath and stem according to your institution's standard practices.
 - Close the SBT cap tightly.
- d. Repeat step 4 for each sample to be tested.

NOTE : A new gauze must be used for each sample to avoid cross-contamination between samples and aerosols.


Figure 2: Break the Swab in the SBT

Once broken, the swab stem should be shorter than the rim of the SBT to adequately discharge the sample in the SBT. The SBT cap should close without hindrance from the swab stem.

PIE PREPARATION

NOTE 1 : Process one (1) sample at a time.

NOTE 2 : Manipulate the DTT only by the bulb to avoid contamination of sample.

Step 5

- a. Mix the SBT for a minimum of 15 seconds at maximal speed (> 3200 rpm) using a vortex mixer.

Step 6

- a. Unseal the PIE pouch and remove the PIE.
- b. Place the Strep A PIE on the Sample Rack, if used or on a flat surface.

NOTE: Hold the PIE by the dilution chamber or by its sides to avoid smearing on PCR wells.

Step 7

- a. Obtain one DTT from the kit box and aspirate the inoculated Sample Buffer (SB) by squeezing the entire bulb, placing the pipette tip into the sample buffer and slowly releasing the pressure on the bulb. The liquid should rise between the two (2) marks on the DTT. (**Figure 3**). If the liquid level is not between the two (2) marks, discharge completely the SB in the SBT and repeat step B.


Figure 3. Representation of an appropriate Sample Buffer (SB) level using the Disposable Transfer Tool (DTT)

NOTE: Do not aspirate bubbles.

- b. Discharge completely the SB volume contained in the DTT into the Sample Loading Chamber of the PIE by squeezing the DTT bulb (**Figure 3**). Make sure not to touch the outer edges or the bottom of the Sample Loading Chamber with the DTT.
- c. Dispose of the DTT according to your institution's standard practices. Do not reuse DTT.
- d. Close the cap of the PIE tightly.
- e. Close the SBT cap tightly.
- f. Repeat steps 5 to 7 for any additional prepared SBT then proceed to step 8.

NOTE : Be cautious when manipulating a loaded PIE. Do not invert, shake or drop. Should it happen, discard and prepare a new PIE starting from step 5.

REVOGENE OPERATION

NOTE 1: Refer to the Revogene Operator's Manual[®] for further information regarding Revogene operation.

Step 8

- Log in by entering <Username> and <Password> and tap <Login> (Figure 4).
- Tap <Setup Run> on the main menu (it will automatically unlock the lid if it is closed) (Figure 5).
- Open the lid and remove the retention ring.


Figure 4. Login Screen

Figure 5. Main Menu Display (Hand on Setup Run Icon)

Step 9

- Enter the **Sample ID** using the barcode scanner or the on-screen keyboard:
 - To scan the SBT, hold it vertically in front of the scanner.
 - To use the on-screen keyboard, tap on the pencil icon. Manually entered Sample ID will need to be confirmed by entering it a second time.
- Scan the **SBT barcode** and the **PIE barcode**, holding the PIE vertically with the cap at the top (Figure 6).

Figure 6. Scanning PIE

NOTE: Do not invert, shake or drop the loaded PIE. Should it happen, discard and prepare a new PIE starting from step 5.

- (Optional) Tap the **pencil** icon of the <Patient ID> line and type the patient ID using the on-screen keyboard.
- (Optional) Tap the **pencil** icon of the <Add Comments> line and type a comment using the on-screen keyboard.
- Insert the PIE into the Revogene carousel while holding it by its sides or by the dilution chamber (Figure 7).
 - Any of the eight (8) positions of the carousel can be used. The software will automatically associate sample and SBT to the correct PIE.
- Tap <ok> on the <insert PIE into the instrument> line to confirm that the PIE is inserted.
- Repeat step 9 for all loaded PIE.

Figure 7. Inserting PIE

- Step 10**
- a. Tap **<Next>**.
 - b. If less than eight (8) PIEs are being tested, load MOCK PIEs in the carousel remaining positions (MOCK PIEs do not need to be scanned) (**Figure 8**). If MOCK PIEs are not available, use new assay PIEs filled with uninoculated SB (BLANK).
 - c. Tap the **Confirm checkbox**.
 - d. Scan the retention ring and place it on the carousel (**Figure 9**).
 - e. Close the instrument lid.
 - f. Tap **<Start>** to initiate the test (**Figure 10**).
 - A timer on the screen and lights on the Revogene lid will show the progression of the test.
 - g. Store the inoculated SBT in appropriate conditions (refer to **Storage and Stability** section) for further repeat testing, if needed.


Figure 8. Adding MOCK PIE

Figure 9. Placing Retention Ring


Figure 10. Confirm Run Start

VIEWING AND EXPORTING RESULTS

Step 11 - During the run

a. If the lights indicating run progression on the Revogene lid start to blink, this is a notification for an Early Positive Result Outcome (E-PRO). An icon will also appear on the title bar of the screen.

The **E-PRO** icon includes the "+" symbol and a number representing the number of Positive results available at this time (**Figure 11**). The number will increase if additional positive results become available.

Figure 11. E-PRO icon

b. Enter <Username> and <Password> and tap <Login> if the user's session has logged-out.

c. Tap the **E-PRO** icon.

d. Results from the current run are automatically listed on the screen (**Figure 12**). For each sample line, there will be either a Positive symbol or an In-Progress symbol. The In-Progress symbol, consisting of a rotating animation, is displayed until a positive result is obtained or the run is completed.

e. Tap on a positive result symbol to open the interim report.

f. (Optional) Tap <Export> and save the interim report where appropriate (e.g. USB key or connectivity option).


Figure 12. In progress Results Display

Once the results screen has been consulted, the E-PRO notifications will disappear. New E-PRO notifications will appear upon the detection of an additional Early Positive Result.

Any sample result displayed in the interim report is final and will remain unchanged in the final assay result report.

Step 12 - At the End of the Run

a. Once the run is completed, the lid opens automatically.

b. Enter <Username> and <Password> and tap <Login> if the user's session has logged-out.

c. Tap **Results** icon (**Figure 13**).

d. Results from the last run are automatically listed on the screen.

e. Select samples for which results report(s) has (have) to be exported:

- All samples can be selected simultaneously by checking the box on the upper left corner of the screen.

f. Tap <Export> and save the final assay result report where appropriate (e.g. USB key or via connectivity option) (**Figure 14**).

g. Alternatively:

- Tap <Search> to find a specific sample and its result.

Figure 13. Main Menu Display (Hand on Results Icon)


Figure 14. Final Results Display

Step 13

- a. Remove the retention ring.
- b. Remove the PIEs from the Revogene:
 - Dispose of the PIEs according to your institution's standard practices.

NOTE : Do not discard MOCK PIEs as they can be reused. Clean and decontaminate MOCK PIE according to the MOCK PIE Package Insert.

REPEAT TESTING PROCEDURE**UNRESOLVED OR INDETERMINATE RESULT FOR A SPECIMEN**

When an Unresolved (UNR) or Indeterminate (IND) result is obtained for a specimen, a repeat test from the corresponding inoculated SBT must be performed within the specified timeframe described in the **Storage and Stability** section. Only one (1) repeat testing from the SBT is allowed.

Use the corresponding inoculated SBT to load a new PIE from the same kit lot. Follow the **PIE Preparation** section from step 5, then follow the **Revogene Operation** section.

UNRESOLVED OR INDETERMINATE RESULT FOR AN EXTERNAL CONTROL

When an Unresolved (UNR) or Indeterminate (IND) result is obtained for an External Control, the run is invalid. A repeat test must be performed for both the Positive and Negative External controls from the corresponding inoculated SBT within the specified timeframe described in the **Storage and Stability** section. Specimens included in the run must be repeated using the corresponding inoculated SBT. Only one (1) repeat testing from the SBT is allowed.

Use the corresponding inoculated SBT to load a new PIE from the same kit lot. Follow the **PIE Preparation** section from step 5, then follow the **Revogene Operation** section.

FALSE NEGATIVE OR FALSE POSITIVE RESULT FOR AN EXTERNAL CONTROL

When a False Negative or a False Positive result is obtained for an External Control, the run is invalid. Both controls must be repeated using a new set of External Controls.

In addition, a repeat testing must be performed for all samples included in the run, using the corresponding inoculated SBT within the specified timeframe described in the **Storage and Stability** section with a new set of External Controls. Only one (1) repeat testing from the SBT is allowed.

Refer to the **Quality Control** section for preparation of a new set of External controls.

Using a new pouch, follow the **PIE Preparation** section from step 5, then follow the **Revogene Operation** section.

QUALITY CONTROL

Quality control procedures monitor the accuracy and precision of the analytical process. Each laboratory must establish the number, type and frequency of testing control materials per applicable regulations or accrediting agencies. The procedure described below may be employed, if appropriate, based on local policies and procedures.

PROCESS CONTROL

Each Strep A PIE contains a Process Control (PrC) that verifies for sample homogenization, sample dilution, cell lysis, inhibition of DNA amplification and assay reagents failure.

EXTERNAL CONTROLS

NOTE 1: Separate DTT, SBT and PIE must be used for each External Control preparation.

1. Good laboratory practices recommend the use of control materials. User should follow the appropriate guidelines concerning the running of External Controls. It is recommended that one (1) Positive External Control and one (1) Negative External Control should be run at least on a daily basis until adequate process validation is achieved with the Revogene Strep A assay on the Revogene in each laboratory setting.
2. External Control materials are not provided by Meridian Bioscience, Inc. External Controls are not used by the Revogene software for the purpose of sample test result interpretation. External Controls are treated as if they are specimens.
3. Before reporting early positive sample results, ensure to have Valid Positive and Negative External Control results.

Positive External Controls

- A cell suspension of a *Streptococcus pyogenes* strain, from commercially available strain (e.g. ATCC® 12344™), prepared at 0.5 ± 0.05 McFarland diluted 1/100 in saline (e.g., BD BBL™ Prepared Saline Solution, cat# 221819) is recommended as a Positive External Control.

Negative External Control

- A cell suspension of a *Streptococcus salivarius* strain, from commercially available strain (e.g. ATCC® 13419™), prepared at 0.5 ± 0.05 McFarland in saline (e.g., BD BBL™ Prepared Saline Solution, cat# 221819) is recommended for use as a Negative External Control.

1. Process and test External Control preparations according to the **Instruction for Use** section.
2. At step 4, remove an unused, compatible swab from its transport sheath and quickly dip it in the Positive or Negative External Control suspension (ensure to completely submerge the swab's tip in the liquid).
3. Transfer the swab in the SBT.
4. Continue with the instructions in the **Instruction for Use** section at step 4-C.

RESULTS INTERPRETATION

The results are computed by the Revogene from measured fluorescent signals and embedded calculation algorithms and are available on the "Results" window. Possible results are:

Sample	Symbol displayed on user screen	Reported Target Result	Interpretation
Patient Specimen		Positive	GAS target DNA detected.
		Negative	GAS target DNA not detected.
		Unresolved	Amplification/detection failure of the Process Control as well as for the GAS target DNA. Could be caused by inhibitory specimens, microfluidic or reagent failure. Repeat testing must be performed (refer to the Repeat Testing Procedure section for further guidance).
		Indeterminate	No reportable result due to possible Revogene detection error during the assay processing, the data analysis, or if the run is interrupted by the user. Repeat testing must be performed (refer to the Repeat Testing Procedure section for further guidance).
		In progress	Results not available yet
Positive External Control		Positive	Valid Positive External Control result.
		Negative	A Positive External Control that yields a negative result is indicative of a sample handling/preparation problem. The run is invalid. Review the sample handling/preparation technique.
		Unresolved	Incorrect Positive External Control result. The run is invalid. Repeat testing must be performed (refer to the Repeat Testing Procedure section for further guidance).
		Indeterminate	Incorrect Positive External Control result. The run is invalid. Repeat testing must be performed (refer to the Repeat Testing Procedure section for further guidance).
		In progress	Results not available yet
Negative External Control		Positive	A Negative External Control that yields a positive test result is indicative of a specimen handling and/or contamination event. The run is invalid. Review the specimen handling technique.
		Negative	Valid Negative External Control result.
		Unresolved	Incorrect Negative External Control result. The run is invalid. Repeat testing must be performed (refer to the Repeat Testing Procedure section for further guidance).
		Indeterminate	Incorrect Negative External Control result. The run is invalid. Repeat testing must be performed (refer to the Repeat Testing Procedure section for further guidance).
		In progress	Results not available yet

LIMITATIONS OF THE PROCEDURE

- The Revogene Strep A assay does not differentiate asymptomatic carriers of Group A *Streptococcus* from those exhibiting streptococcal infection.
- Additional follow-up testing using the culture method is required if the result is negative and clinical symptoms persist, or in the event of an acute rheumatic fever (ARF) outbreak.
- The Revogene Strep A assay must only be used with the Revogene by trained personnel.
- The Revogene Strep A assay does not provide susceptibility results. Culture isolates are needed for performing susceptibility testing.
- Clinical performance characteristics of the Revogene Strep A assay were established only with throat swab specimens collected with rayon swabs in Liquid Stuart or Liquid Amies transport media (Copan cat. # 138C and 139C). In addition, throat swab specimens collected with rayon or polyester swabs in Liquid Stuart or Liquid Amies transport media were evaluated analytically (Copan cat. # 140C and 141C; BD cat # 220146 and 220149). Use of BD CultureSwab™ MaxV rayon swab in Liquid Amies transport media is not recommended with the Revogene Strep A assay (BD cat.# 220232). Use of the Revogene Strep A assay with clinical specimen in collection/transport systems other than those specified has not been evaluated and performance characteristics are not established.
- Results from the Revogene Strep A assay should be used as an adjunct to clinical observations and other information available to the physician.
- Assay results may be affected by concurrent antimicrobial therapy as GAS DNA may continue to be detected.
- A positive assay result does not necessarily indicate the presence of viable organisms. It is, however, indicative for the presence of GAS DNA.
- Erroneous assay results may occur from improper specimen collection, handling or storage, technical error, sample mix-up, or because the sample does not contain a sufficient bacterial DNA load for detection by the assay.
- Careful compliance to the instructions of this insert, the Revogene Operator's Manual and to established guidelines is necessary to avoid erroneous results.
- Contamination or false negative results may occur if a PIE cap is incorrectly closed or if a droplet has been dropped on the edges of the Sample loading Chamber.
- Mutations or polymorphisms in primer- or probe-binding regions may affect detection of *S. pyogenes* gene variants, resulting in a false negative result with the Revogene Strep A assay.
- Analgesic / Antipyretic (e.g. Tylenol®), Nonsteroidal anti-inflammatory drug (Advil®), Bronchodilator (e.g. Albuterol sulfate) or Mucin may potentially interfere with the Revogene Strep A assay when either of these substances is present in Sample Buffer at a concentration of > 0.4% (w/v).
- Nonsteroidal anti-inflammatory drug (Aspirin®) or Whole Blood may potentially interfere with the Revogene Strep A assay when either of these substances is present in Sample Buffer at a concentration of > 0.1% (w/v or v/v).
- A combinatory effect of *Streptococcus mutans*, *Streptococcus oralis*, *Streptococcus salivarius* and *Moraxella catarrhalis*, at $\geq 10^5$ CFU/mL of Sample Buffer or the presence of *Streptococcus dysgalactiae* subsp. *dysgalactiae*, at $\geq 10^5$ CFU/mL of Sample Buffer, may have an inhibitory effect on the detection of *S. pyogenes*.

EXPECTED VALUES

In a prospective clinical trial, a total of 767 throat swab specimens were collected at seven (7) geographically diverse sites across the United States and Canada from February to July 2018. Of the 604 specimens meeting all inclusion criteria without meeting any of the exclusion criteria, 154 were positive based on the Reference Method for an observed prevalence of 25.5% [154/604; 95%CI 22.2% - 29.1%]. The percentage of positive results observed with the Revogene Strep A assay in the study population was 29.0% [175/604; 95%CI 25.5% - 32.7%].

CLINICAL PERFORMANCE CHARACTERISTICS

The performance of the Revogene Strep A assay was established in a prospective clinical trial conducted at eight (8) geographically diverse sites (seven (7) collecting sites and one (1) reference center) across the United States and Canada from February to July 2018.

Dual throat swab specimens were collected from patients with signs and symptoms of pharyngitis with the exception of patients under antibiotherapy that can affect the viability of the *Streptococcus* genus. One (1) swab was tested with the Revogene Strep A assay at one (1) of the seven (7) clinical sites while the second swab was sent to one (1) central laboratory per country (one (1) in Canada and one (1) in the United States) for Reference Method testing. Upon arrival at the central laboratory for Reference Method testing, a Blood Agar Plate (BAP) was inoculated with the second swab of the dual throat swab samples. Following incubation, colonies morphologically resembling *S. pyogenes* were isolated on a new BAP. Isolates were further characterized and identified by catalase test, Gram stain, Latex agglutination and MALDI-TOF.

The sensitivity, specificity, Positive Predictive Value (PPV) and Negative Predictive Value (NPV) values were calculated by comparing Revogene Strep A assay results with the Reference Method. Discrepant analysis was performed on specimens with discrepant results between the Revogene Strep A assay and Reference Method, using an alternative PCR for the amplification and detection of *S. pyogenes* followed by bi-directional sequencing of the amplified product.

RESULTS

Out of the 767 enrolled specimens, 163 specimens were excluded: three (3) failed to comply with the Revogene testing procedures, 13 exceeded the delay in shipment, 23 did not meet inclusion criteria, 27 did not have complete testing results, and 97 failed to comply with the reference culture protocol.

A total of 604 specimens were used to establish the performance of the Revogene Strep A assay in comparison with bacterial culture. The study population demographics are presented in **Table 1**.

Table 1. Study Population Demographics

Age group	All Subjects (N = 604)
> 24 months old – 12 years old	278 (46.03%)
13 – 21 years old	80 (13.25%)
22 – 64 years old	235 (38.91%)
≥ 65 years old	11 (1.82%)
Ethnicity	
White	186 (30.79%)
Black or African American	231 (38.25%)
Hispanic	97 (16.06%)
Asian	5 (0.83%)
Unknown	85 (14.07%)

COMPARISON WITH REFERENCE METHOD

The clinical performance of the Revogene Strep A assay compared to the Reference Method results is presented in **Table 2**.

The Revogene Strep A assay demonstrated 98.1% sensitivity (151/154; 95%CI: 94.4-99.3) and 94.7% specificity (426/450; 95%CI: 92.2-96.4) in comparison to the Reference Method.

Table 2. Overall Performance Characteristics of the Revogene Strep A Assay in Comparison to the Reference Method

Revogene Strep A	Reference Method		
	Positive	Negative	Total
	151	24 ^a	175
	3 ^b	426	429
	Total	450	604
Sensitivity	98.1% [95% CI: 94.4-99.3]		
Specificity	94.7% [95% CI: 92.2-96.4]		
Positive Predictive Value	86.3% [95% CI: 80.4-90.6]		
Negative Predictive Value	99.3% [95% CI: 98.0-99.8]		

^a Results from alternative PCR with bi-directional sequencing: 17 of 24 were Strep A Positive; six (6) of 24 were Strep A Negative and one (1) was inconclusive.

^b Results from alternative PCR with bi-directional sequencing: one (1) of three (3) was Strep A Positive and two (2) were inconclusive.

The overall initial unresolved rate was 0.5% (3/610) and 0.2% (1/610) after repeat testing. The overall initial indeterminate rate was 1.1% (7/610) and 0.8% (5/610) after repeat testing. The overall initial non-reportable rate was 1.6% (10/610) and 1.0% (6/610) after repeat testing.

ANALYTICAL PERFORMANCE CHARACTERISTICS

ANALYTICAL SENSITIVITY

The analytical sensitivity (Limit of Detection or LoD) of the Revogene Strep A assay was determined using clinical throat swab specimen matrix (previously tested negative for *S. pyogenes*) spiked with various concentrations of *S. pyogenes* bacterial suspensions. Three (3) strains of *S. pyogenes* (ATCC® 12344™, ATCC® 19615™, ATCC® 700942™) were tested in 24 replicates per concentration. The LoD is defined as the lowest concentration at which 95% or more of all replicates tested positive.

In regard to the three (3) strains tested, the LoD of the Revogene Strep A assay ranged from 333 to 1333 CFU/mL of SB (**Table 3**). The preliminary LoD concentrations were confirmed by testing twenty (20) replicates of each strain.

Table 3. LoD of the Revogene Strep A Assay

<i>S. pyogenes</i> strain (ATCC® number)	LoD (CFU/mL of SB)
ATCC® 12344™	333
ATCC® 19615™	1000
ATCC® 700942™	1333

INCLUSIVITY

Inclusivity of the Revogene Strep A assay was determined for 12 strains of *S. pyogenes* representing six (6) different *emm* sequence types from diverse geographic origins. Nine (9) strains were tested from a quantitated cell culture in presence of *S. pyogenes*-negative throat matrix at a load of 999 CFU/mL of SB, corresponding to three (3) times the LoD value of ATCC® 12344™ strain. Nine (9) replicates per strain were tested using three (3) different Revogene Strep A kit lots. Eight (8) strains were detected with 100% positivity (in 9/9 replicates) by the Revogene Strep A assay at a load of 999 CFU/mL of SB. The *S. pyogenes* ATCC® 49399™ strain was detected with 100% positivity at a load of 1667 CFU/mL of SB. Strains tested are described in **Table 4**.

Table 4. *S. pyogenes* Strains Tested for Inclusivity with the Revogene Strep A assay

<i>S. pyogenes</i> strains	Lancefield grouping	<i>emm</i> type
ATCC® 12384™	A	3
ATCC® 49399™	A	N/A
CCUG 33409	A	N/A
CCUG 39158	A	N/A
CCUG 53553	A	N/A
ATCC® 700952™	A	92
ATCC® 700294™	A	1
ATCC® 12357™	A	18
CCUG 65322	A	N/A
ATCC® 19615™ ¹	A	80
ATCC® 12344™ ¹	A	1
ATCC® 700942™ ¹	A	82

¹ Strains used for the LoD determination and therefore not retested

In addition, an *in silico* analysis was performed on August 8th, 2018 to assess inclusivity of the primers and probe of the Revogene Strep A assay target regarding 123 *S. pyogenes* strains listed in the National Center for Biotechnology Information (NCBI) database. This include sequences of *S. pyogenes* strains of significant *emm* types included in candidate vaccines against Group A Streptococci⁹. The alignment results showed no mismatch with 122 sequences. One (1) variant presented a mismatch in the region corresponding to the *S. pyogenes* probe, however, this variant is considered detectable by the Revogene Strep A assay. The analysis predicted the detection of all these *S. pyogenes* strains.

CROSS-REACTIVITY

The cross-reactivity of the Revogene Strep A assay was assessed with high loads of organisms that are not targeted by the assay, or phylogenetically related to *S. pyogenes*, or found in clinical throat swab specimens. The study included 42 bacteria, one (1) yeast, six (6) viruses and human DNA (Table 5). Bacteria and yeast were tested at a load of 10⁶ CFU/mL of SB. Nucleic acids from viruses and human DNA were tested at a load of 10⁵ DNA or RNA copies/mL of SB. These organisms were tested using a quantitated cell culture or nucleic acid solution in presence of *S. pyogenes*-negative throat matrix. Each organism was tested in PCR triplicates.

Under the conditions of the study, all organisms and nucleic acids were found non-reactive with the Revogene Strep A assay.

The cross-reactivity with primers and probes of the Revogene Strep A assay was evaluated by an *in silico* analysis performed on sequences contained in the NCBI database on August 8th, 2018. No other target than *S. pyogenes* was found to have significant level of homology with the *S. pyogenes* primers and probe.

Table 5. List of Organisms Tested for Cross-Reactivity with the Revogene Strep A Assay

Bacteria	
<i>Acinetobacter lwoffii</i>	<i>Staphylococcus aureus</i>
<i>Arcanobacterium haemolyticum</i>	<i>Staphylococcus epidermidis</i>
<i>Bacillus cereus</i>	<i>Stenotrophomonas maltophilia</i>
<i>Bordetella pertussis</i>	<i>Streptococcus agalactiae</i>
<i>Burkholderia cepacia</i>	<i>Streptococcus anginosus</i>
<i>Corynebacterium diphtheriae</i>	<i>Streptococcus anginosus</i> subsp. <i>whileyi</i>
<i>Enterococcus faecalis</i>	<i>Streptococcus bovis</i>
<i>Escherichia coli</i>	<i>Streptococcus canis</i>
<i>Fusobacterium necrophorum</i>	<i>Streptococcus constellatus</i> subsp. <i>viborgensis</i>
<i>Haemophilus influenza</i>	<i>Streptococcus dysgalactiae</i> subsp. <i>dysgalactiae</i>
<i>Klebsiella pneumoniae</i>	<i>Streptococcus dysgalactiae</i> subsp. <i>dysgalactiae</i> *
<i>Lactococcus lactis</i>	<i>Streptococcus gordonii</i>
<i>Legionella jordanis</i>	<i>Streptococcus intermedius</i>
<i>Legionella micdadei</i> (Tatlockia micdadei)	<i>Streptococcus mitis</i>
<i>Legionella pneumophila</i>	<i>Streptococcus mutans</i>
<i>Moraxella catarrhalis</i>	<i>Streptococcus oralis</i>
<i>Neisseria gonorrhoeae</i>	<i>Streptococcus pneumoniae</i>
<i>Neisseria subflava</i>	<i>Streptococcus salivarius</i>
<i>Parvimonas micra</i>	<i>Streptococcus sanguinis</i>
<i>Pseudomonas aeruginosa</i>	<i>Streptococcus suis</i>
<i>Serratia marcescens</i>	<i>Veillonella parvula</i>
Yeast	
<i>Candida albicans</i>	
Viruses	
Adenovirus Type 1	Adenovirus type 11 (Slobitski)
Influenza B/Taiwan/2/62	Influenza A/Aichi/2/68/H3N2
Rhinovirus Type 17	Parainfluenza virus 4b
Human DNA	
Human DNA	

*Different ATCC strain

INTERFERING ORGANISMS

The potentially inhibitory effect of 30 non-target organisms that may be present in the normal flora of the throat, mouth and/or respiratory tract was assessed using organisms selected from the cross-reactivity study (Table 6). Each organism category (i.e. bacteria, yeast, viruses) was represented with special attention to include the most frequent causative agents of respiratory tract infections. Groups of one (1) to four (4) organisms were tested in triplicate with either 999 CFU/mL of SB containing *S. pyogenes* ATCC® 12344™ strain, 3000 CFU/mL of SB containing *S. pyogenes* ATCC® 19615™ strain, or *S. pyogenes*-negative samples in presence of *S. pyogenes*-negative throat matrix to assess their potential interference on detection of *S. pyogenes* or PrC. Each organism within a group was diluted to reach a load of 10⁶ CFU/mL of SB for bacteria and yeast, and 10⁵ copies/mL of SB for virus. The 30 organisms included in the study are presented in Table 6.

Table 6. List of Organisms Tested for Interference with the Revogene Strep A Assay

Group 1	
<i>Acinetobacter lwoffii</i>	<i>Arcanobacterium haemolyticum</i>
<i>Haemophilus influenza</i>	<i>Bordetella pertussis</i>
Group 2	
<i>Klebsiella pneumoniae</i>	<i>Parvimonas micra</i>
<i>Staphylococcus aureus</i>	<i>Corynebacterium diphtheriae</i>
Group 3	
<i>Streptococcus anginosus</i>	<i>Streptococcus anginosus</i> subsp. <i>whileyi</i>
<i>Streptococcus constellatus</i> subsp. <i>viborgensis</i>	<i>Fusobacterium necrophorum</i>
Group 4	
<i>Streptococcus dysgalactiae</i> subsp. <i>dysgalactiae</i>	
Group 5	
<i>Streptococcus dysgalactiae</i> subsp. <i>dysgalactiae</i>	
Group 6	
<i>Streptococcus gordonii</i>	<i>Streptococcus intermedius</i>
<i>Streptococcus mitis</i>	<i>Legionella pneumophila</i>
Group 7	
<i>Streptococcus mutans</i>	<i>Streptococcus oralis</i>
<i>Streptococcus salivarius</i>	<i>Moraxella catarrhalis</i>
Group 8	
<i>Streptococcus sanguinis</i>	<i>Veillonella parvula</i>
<i>Candida albicans</i>	<i>Neisseria gonorrhoeae</i>
Group 9	
Adenovirus Type 1	Influenza B/Taiwan/2/62
Rhinovirus Type 17	Influenza A/Aichi/2/68/H3N2

None of the 30 organisms present at 10⁶ CFU/mL of SB for bacteria and yeast and 10⁵ copies/mL of SB for viruses interfered with detection of PrC or with *S. pyogenes* strain ATCC® 19615™.

Organisms included in group 5 and group 7 showed a potentially inhibitory effect on detection of *S. pyogenes* strain ATCC® 12344™. Nevertheless, the *S. dysgalactiae* subsp. *dysgalactiae* strain showed no interference at a load of 10⁵ CFU/mL of SB. When each organism from group 7 was tested individually at a load of 10⁵ CFU/mL of SB in presence of *S. pyogenes* strain ATCC® 12344™, no interference was observed.

INTERFERING SUBSTANCES

The potentially inhibitory effect of 16 exogenous and three (3) endogenous substances that may be present in a throat swab specimen was assessed using *S. pyogenes*-negative samples and *S. pyogenes* strains ATCC® 12344™ and ATCC® 19615™ containing samples at three (3) times the LoD (999 CFU/mL of SB and 3000 CFU/mL of SB, respectively) in presence of *S. pyogenes*-negative throat matrix. Substances were tested at their potentially highest concentration that could be found in a throat swab specimen. The results for the 19 substances are presented in **Table 7**.

Results demonstrated that Analgesic / Antipyretic (Tylenol®), Nonsteroidal anti-inflammatory drug (Aspirin®, Advil®), Bronchodilator (Albuterol sulfate), Whole Blood and Mucin showed a potentially inhibitory effect on the detection of PrC or *S. pyogenes* when any of these substances was present in SBT at a concentration of 4.3% (w/v or v/v). These substances showed no reportable interference with the Revogene Strep A assay when tested at 0.4% (w/v) for Tylenol®, Advil®, Albuterol sulfate and Mucin or at 0.1% (w/v) or 0.1% (v/v) for Aspirin® and Whole Blood, respectively.

Table 7. List of Exogenous and Endogenous Substances Tested with the Revogene Strep A Assay

Exogenous substances	
Substance (commercial name)	Concentration in SB without interference observed ¹
Analgesic/ Antipyretic (Tylenol®)	0.4 % (w/v)*
Nonsteroidal anti-inflammatory drug (Aspirin®)	0.1 % (w/v)*
Nonsteroidal anti-inflammatory drug (Advil®)	0.4 % (w/v)*
Bronchodilator (Albuterol sulfate)	0.4 % (w/v)*
Nasal decongestant (Zicam®)	4.3 % (w/v)
Cough lozenges (Halls Menthon-Lyptus)	4.3 % (w/v)
Cough lozenges (Cépacol® Extra Strength)	4.3 % (w/v)
Cough lozenges (Sucrets® Complete)	4.3 % (w/v)
Toothpaste (Crest® Pro-Health Advanced)	0.56 % (w/v)
Mouthwash (Crest® Scope® Classic)	4.3 % (v/v)
Mouthwash (Listerine® Ultraclean)	4.3 % (v/v)
Cough syrup (NyQuil™ Complete)	4.3 % (v/v)
Cough syrup (Dimetapp®)	4.3 % (v/v)
Cough syrup (Robitussin® Cough Control Extra Strength)	4.3 % (v/v)
Cough syrup (Robitussin Total® Cough, Cold & Flu Extra Strength Nighttime)	4.3 % (v/v)
Antihistamines (Benadryl® Allergy Elixir)	4.3 % (v/v)
Endogenous substances	
Substance	Concentration in SB without interference observed ¹
Human saliva	4.3 % (v/v)
Whole Blood	0.1 % (v/v)*
Mucin from bovine submaxillary glands	0.4 % (w/v)*

*w/v: Weight/Volume; v/v: Volume/Volume
*Concentrations above these numbers may cause interference (see limitations above)

CARRY-OVER AND CROSS-CONTAMINATION

The carry-over (between-run) and cross-contamination (within-run) were assessed using positive samples tested in presence of *S. pyogenes*-negative throat matrix. Positive samples were prepared with *S. pyogenes* ATCC® 12344™ strain at a final concentration of 10⁶ CFU/mL of SB. Negative samples were also tested in presence of *S. pyogenes*-negative throat matrix.

For the carry-over study, a total of ten (10) runs were performed by two (2) operators with the Revogene Strep A assay on one (1) Revogene. Four (4) high positive samples and four (4) negative samples were tested by alternating positive and negative samples in each run. For the cross-contamination study, a run of eight (8) replicates of high positive samples followed by a run of eight (8) replicates of negative samples were performed by two (2) operators, for a total of ten (10) runs on one (1) Revogene.

Absence of carry-over and cross-contamination was demonstrated.

REPRODUCIBILITY

Between-site reproducibility study was performed on three (3) sites, by two (2) operators per site, over five (5) distinct days using one (1) Revogene Strep A kit lot.

Between-lot reproducibility study was performed at one (1) site by two (2) operators, over 15 days using three (3) Revogene Strep A assay kit lots (five (5) days per kit lot). For each reproducibility study, a total of 120 replicates for negative samples and 90 replicates for each category of positive samples were tested in presence of *S. pyogenes*-negative throat matrix. Two (2) *S. pyogenes* strains were used for positive samples: ATCC® 12344™ and ATCC® 19615™.

Sample categories were described as follows:

- Low Positive (LP): *S. pyogenes* ATCC® 12344™ strain was tested at 650 CFU/mL of SB
- Moderate Positive (MP): *S. pyogenes* ATCC® 19615™ strain was tested at 3000 CFU/mL of SB
- Negative (NEG): samples without *S. pyogenes* strain

Results for each sample category tested during the between-site and the between-lot reproducibility studies are shown in **Tables 8 and 9** respectively. For the between-site reproducibility, the overall percent agreement for NEG, LP and MP was 100% (**Table 8**). For the between-lot reproducibility, the overall percent agreement for NEG, LP and MP was 100%, 97.8% and 100% respectively (**Table 9**). Overall mean Cycle Threshold (Ct) values with variance components (SD and %CV) are shown in **Tables 8 and 9**.

Table 8. Between-Site Reproducibility Study Results using One (1) Revogene Strep A Assay Kit Lot

Category	<i>S. pyogenes</i> strain (ATCC® number)	Site 1		Site 2		Site 3		Overall Results/ Total	Overall Percent Agreement ¹	Overall 95% CI	Ct Values ²		
		Results/ Total	Percent Agreement	Results/ Total	Percent Agreement	Results/ Total	Percent Agreement				Overall Mean	SD	%CV
LP	ATCC® 12344™	30/30	100%	30/30	100%	30/30	100%	90/90	100%	96.7%-100%	36.97	1.88	5.10
MP	ATCC® 19615™	30/30	100%	30/30	100%	30/30	100%	90/90	100%	96.7%-100%	34.98	0.75	2.15
NEG	N/A	40/40	100%	40/40	100%	40/40	100%	120/120	100%	97.5%-100%	28.52	0.69	2.42

¹ For the NEG category, percent agreement was calculated for negative results.
² For the LP and MP categories, Ct values reported are for the GAS target. For the NEG category, Ct values reported are for the PrC.

Table 9. Between-Lot Reproducibility Study Results on One (1) Site using Three (3) Revogene Strep A Assay Kit Lots

Category	<i>S. pyogenes</i> strain (ATCC® number)	Lot 1		Lot 2		Lot 3		Overall Results/ Total	Overall Percent Agreement ¹	Overall 95% CI	Ct Values ²		
		Results/ Total	Percent Agreement	Results/ Total	Percent Agreement	Results/ Total	Percent Agreement				Overall Mean	SD	%CV
LP	ATCC® 12344™	30/30	100%	29/30	96.6%	29/30	96.7%	88/90	97.8%	92.2%-99.7%	37.85	2.30	6.08
MP	ATCC® 19615™	30/30	100%	30/30	100%	30/30	100%	90/90	100%	96.7%-100%	35.49	1.61	4.53
NEG	N/A	40/40	100%	40/40	100%	40/40	100%	120/120	100%	97.5%-100%	28.55	0.74	2.61

¹ For the NEG category, percent agreement was calculated for negative results.
² For the LP and MP categories, Ct values reported are for the GAS target. For the NEG category, Ct values reported are for the PrC.

PRECISION

Precision study was performed on one (1) site, by two (2) operators, over 20 days using one (1) Revogene Strep A assay kit lot. For each sample category, two (2) runs per day with two (2) replicates per run were performed for a total of 80 replicates.

For each category, all samples were tested in presence of *S. pyogenes*-negative throat matrix. Two (2) *S. pyogenes* strains were used for positive samples: ATCC® 12344™ and ATCC® 19615™.

Sample categories were described as follows:

- Low Positive (LP): ATCC® 12344™ strain was tested at 650 CFU/mL of SB
- Moderate Positive (MP): ATCC® 19615™ strain was tested at 3000 CFU/mL of SB
- Negative (NEG): samples without *S. pyogenes* strain

Precision study results for NEG, LP and MP demonstrated 97.5%, 100% and 100% agreement, respectively (Table 10).

Table 10. Precision Study Percent Agreement on One (1) Site using One (1) Revogene Strep A Assay Kit Lot

Category	<i>S. pyogenes</i> strain ATCC® number	Overall Percent Agreement ¹	Overall 95% CI
LP	ATCC® 12344™	100%	96.3%-100%
MP	ATCC® 19615™	100%	96.3%-100%
NEG	N/A	97.5%	91.3%-99.7%

¹ For the NEG category, percent agreement was calculated for negative results.

E-LABELING

Documentation related to this product can be accessed online at www.meridianbioscience.com/pi. Additionally, paper copies are available upon request by contacting your local distributor or via the phone numbers listed on the kit box.

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REFERENCES

1. Bisno AL, Gerber MA, Gwaltney JM Jr, Kaplan EL, Schwartz RH, for the Infectious Diseases Society of America. Practice guidelines for the diagnosis and management of group A streptococcal pharyngitis. Clin Infect Dis. 2002;35(2):113–125.
2. Wegner DL, Witte DL, Schrantz RD. Insensitivity of rapid antigen detection methods and single blood agar plate culture for diagnosing streptococcal pharyngitis. JAMA 1992;267:695–697. 10.1001/jama.1992.03480050099033.
3. Shulman ST, Bisno AL, Clegg HW, Gerber MA, Kaplan EL, Lee G, Martin JM, Van Beneden C. Clinical practice guideline for the diagnosis and management of group A streptococcal pharyngitis: 2012 update by the Infectious Diseases Society of America. Clin Infect Dis 2012;55:e86–e102. 10.1093/cid/cis629
4. Gerber MA, Baltimore RS, Eaton CB, Gewitz M, Rowley AH, Shulman ST, Taubert KA. Prevention of rheumatic fever and diagnosis and treatment of acute streptococcal pharyngitis: a scientific statement from the American Heart Association Rheumatic Fever, Endocarditis, and Kawasaki Disease Committee of the Council on Cardiovascular Disease in the Young, the Interdisciplinary Council on Functional Genomics and Translational Biology, and the Interdisciplinary Council on Quality of Care and Outcomes Research: endorsed by the American Academy of Pediatrics. Circulation 2009;119:1541–1551. 10.1161/CIRCULATIONAHA.109.191959.
5. Snow V, Motur-Pilson C, Cooper RJ, Hoffman JR. 2001. Principles of appropriate antibiotic use for acute pharyngitis in adults. Ann Intern Med. 134:506–508. 10.7326/0003-4819-134-6-200103200-00018.
6. Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health. Biosafety in microbiological and biomedical laboratories 5th ed. HHS Publication No. (CDC) 21-1112. Revised December 2009
7. CLSI. Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline-4th ed. CLSI document M29-A4. Wayne, PA : Clinical and Laboratory Standards Institute; 2014.
8. SN134822 Revogene® Operator's Manual
9. Dale JB, Penfound TA, Chiang EY, Walton WJ. New 30-valent M protein-based vaccine evokes cross-opsonic antibodies against non-vaccine serotypes of group A streptococci. Vaccine. Oct 2011;26(29(46):8175-8. Page 22

Manufactured By

Meridian Bioscience, Inc.
Corporate Office
3471 River Hills Drive
Cincinnati, Ohio 45244 USA
Telephone: 513.271.3700
Orders/Customer Service:
800.543.1980
Technical Support Center:
800.343.3858
Information Fax: 513.272.5432
Ordering Fax: 513.271.0124

Authorized Representative

Meridian Bioscience Europe S.r.L.
Via dell'Industria, 7
20035 Villa Cortese, Milano
ITALY
Tel: +39 0331 43 36 36
Fax: +39 0331 43 36 16
Email: info@meridianbioscience.eu
WEB: www.meridianbioscience.com/eu

Meridian Bioscience Europe s.a./n.v.
2 Avenue du Japon - 1420 Braine l'Alleud
BELGIUM
Tel: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
Email: info.bnl@meridianbioscience.eu

Meridian Bioscience Europe France
34 rue de Ponthieu - 75008 Paris
FRANCE
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 9 70 06 62 10
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe B.V.
Zekeringstraat 17 A
1014BM Amsterdam
NETHERLANDS
Tel: +31 (0) 411 62 11 66
Fax: +31 (0) 411 62 48 41
Email: Info.bnl@meridianbioscience.eu

INTERNATIONAL SYMBOLS USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

	Use-by date / Data di scadenza / Date de péremption / Fecha de caducidad / Verwendbar bis
	Batch code / Codice di lotto / Code de lot / Código de lote / Chargennummer
	In vitro diagnostic medical device / Dispositivo medico diagnostico <i>in vitro</i> / Instrument de test diagnostique <i>in vitro</i> / Producto sanitario para diagnóstico <i>In vitro</i> / <i>In-vitro</i> -Diagnostikum
	CE Mark / Marcatura CE / Symbole CE / Marcado CE / CE-Kennzeichen
	Catalog number / Numero di catalogo / Référence catalogue / Referencia / Bestellnummer
	Consult instructions for use / Consultare le istruzioni per l'uso / Consulter le mode d'emploi / Consultar las instrucciones de uso / Gebrauchsanweisung beachten
	Manufacturer / Produttore / Fabricant / Fabricante / Hersteller
	Contains sufficient for <n> tests / Contiene una quantità sufficiente per <n> test / Contient le matériel suffisant pour <n> tests / Contiene la cantidad suficiente para <n> ensayos / Inhalt ausreichend für <n> Tests
	Temperature limit / Limite di temperatura / Limite de température / Limite de temperatura / Temperaturgrenze
	Authorized representative in the European Community / Rappresentante autorizzato nella Comunità Europea / Représentant agréé dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter EU-Repräsentant
	Do not reuse / Non riutilizzare / Ne pas réutiliser / No reutilizar / Nicht wiederverwenden
	Keep dry / Conservare all'asciutto / Conserver au sec / Mantener seco / Vor Feuchtigkeit schützen
	Contains # pouches: 1 Disposable Transfer Tool (DTT), 1 Sample Buffer Tube (SBT), 1 PIE / Contiene # buste: 1 PIE, 1 Sample Buffer Tube (SBT), 1 Disposable Transfer Tool (DTT) / Contient # sachets: Une cartouche PIE, Un tube échantillon, Un outil de transfert jetable (OTJ) / Incluye # bolsitas: 1 PIE, 1 Sample Buffer Tube (SBT), 1 Disposable Transfer Tool (DTT) / Enthält # Beutel: 1 PIE, 1 Sample Buffer Tube (SBT), 1 Disposable Transfer Tool (DTT)
	Humidity Limitation / Limitazione dell'umidità / Limite d'humidité / Limitación de humedad / Feuchtigkeitsbegrenzung
	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Ne pas utiliser si l'emballage est endommagé / No utilizar si el envase está dañado / Nicht verwenden, wenn die Verpackung beschädigt ist
	Keep away from sunlight / Tenere lontano dalla luce del sole / Tenir à l'écart de la lumière du soleil / Mantener alejado de la luz solar / Vor Sonnenlicht schützen
	Contains 24 Disposable Transfer Loops (DTL) / Contiene 24 loop di trasferimento monouso (DTL) / Contient 24 boucles de transfert jetables (DTL) / Contiene 24 bucles de transferencia desechables (DTL) / Enthält 24 Einwegübertragungsschleifen (DTL)
	Disposable Transfer Tool (DTT) / Disposable Transfer Tool 9DTT / Un outil de transfert jetable (OTJ) / Disposable Transfer Tool (DTT) / Dispoable Transfer Tool (DTT)
	Revogene Test Device / Dispositivo Test Revogene / Dispositif de test Revogene / Dispositivo para la Prueba Revogene / Revogene-Analysegefaß
	For Emergency Use Authorization only / Solo per l'autorizzazione per l'uso di emergenza / pour autorisation d'utilisation d'urgence uniquement / para autorización de uso de emergencia solamente / nur für Notfallverwendungsautorisierung
	Disposable Transfer Loop / Loop di trasferimento monouso / Boucles de transfert jetables / Bucles de transferencia desechables / Einwegübertragungsschleifen
	Sample Buffer Tube / Sample Buffer Tube / Un tube échantillon / Smample Buffer Tube / Sample Buffer Tube
	Revogene Test Device / Dispositivo Test Revogene / Dispositif de test Revogene / Dispositivo para la Prueba Revogene / Revogene-Analysegefaß

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order call Customer Service Department at 800-543-1980.

TaqMan is a registered trademark of Roche Molecular Systems, Inc. ATCC is a trademark of the American Type Culture Collection.
BBL and CultureSwab are trademarks of Becton, Dickinson and Company. Revogene and associated logos are trademarks of Meridian Bioscience, Inc.
© 2020-12 Meridian Bioscience, Inc. Made in Canada

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Revogene e i loghi associate sono marchi commerciali di Meridian Bioscience, Inc.
© 2020-12 Meridian Bioscience, Inc. Prodotto in Canada.
Revogene et les logos connexes sont des marques déposées de Meridian Bioscience, Inc.
© 2020-12 Meridian Bioscience, Inc. Fabriqué au Canada.
Revogene así como los logotipos asociados son marcas comerciales de Meridian Bioscience, Inc.
©2020-12 Meridian Bioscience, Inc. Hecho en Canadá