alethia[®]

CMV

External Control Kit

External Control Materials for use with
Alethia CMV DNA Amplification Assay

REF 479880

IVD

Rx Only

INTENDED USE

The Alethia CMV Assay Test System includes separately provided test kits for the Alethia CMV DNA Amplification Assay and the Alethia CMV External Control Reagents.

The Alethia CMV External Control Reagents are used as part of a routine quality control program to aid the user in detection of unexpected conditions that may lead to test errors. The external controls are intended for use with the Alethia CMV DNA Amplification Assay; the controls are not intended for use with other assays or systems.

SUMMARY AND EXPLANATION OF THE TEST

External controls are used as part of a routine quality control program to aid the user in detection of unexpected conditions that may lead to test errors.

Quality control testing is performed to detect factors such as reagent deterioration, adverse environmental or test conditions, or variance in operator performance that can cause test errors. External control reagents, such as the Alethia CMV External Positive Control and Alethia CMV External Negative Control, are reagents that are not built into the test system but are amplified and detected in the same manner as patient specimens.

BIOLOGICAL PRINCIPLES

The Alethia CMV DNA Amplification test kit should be tested with known positive and negative external control reagents on a periodic basis to verify that the kits are performing correctly. The frequency with which a laboratory performs external controls will be affected by:

1. The proficiency level of the laboratory,
2. The laboratory's own internal requirements,
3. The requirements of the laboratory's accrediting agencies,
4. The number of new and different operators performing the test,
5. Whether a new kit lot is being added to testing,
6. Whether the kit lot is from a different shipment, and
7. Whether deviations from the manufacturer's stated storage or handling conditions have occurred.

NOTE: When unacceptable quality control test results are obtained, all test results should be considered invalid. QC test failures are an indication that either reagents, test environment or operator performance have changed.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. Alethia CMV Positive Control: Tris-buffered solution with plasmid containing DNA inserts (Cytomegalovirus and human mitochondrial DNA inserts) and azide (0.09%) as a preservative.
2. Alethia CMV Negative Control: Tris-buffered solution with plasmid containing human mitochondrial DNA insert and azide (0.09%) as a preservative.

MATERIALS PROVIDED SEPARATELY

1. Alethia CMV DNA Amplification Test Kit, Meridian Bioscience, Inc. Catalog Number: 481325

MATERIALS NOT PROVIDED

1. Disposable latex gloves, powder free
2. DNase/RNase free, aerosol resistant pipette tips

EQUIPMENT NOT PROVIDED

1. Interval timer
2. Vortex Mixer
3. Micropipette(s) capable of dispensing 50 µL
4. Micropipette(s) capable of dispensing 100 µL
5. Alethia[®] Instrument, Meridian Bioscience, Inc. Catalog Number: 610189

- PRECAUTIONS**
1. All reagents are for in vitro diagnostic use only.
 2. This is a quality control reagent and is used to evaluate the performance of the Alethia CMV DNA Amplification Assay. It is not directly used to test patient samples.
 3. Transport and store these reagents according to the kit label. Do not freeze.
 4. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
 5. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
 6. Quality Control Programs for Molecular Testing Laboratories should be employed.
 7. Alethia Test Devices contain lyophilized reagents. The protective pouch should not be opened until ready to perform the assay.
 8. Alethia Test Devices include a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
 9. Dispose of used Alethia CMV Test Devices and tubes immediately after processing. Leave the Test Device latch securely in place. Do NOT open the Test Device after processing. Opening the device after amplification may result in contamination of the test area with amplification product.

HAZARD AND PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-30 C.

REAGENT PREPARATION

1. Bring Alethia CMV External Controls and all Alethia CMV kit components to room temperature (19-30 C) before use. Incorrect results may be obtained if Control Materials and components are not brought to room temperature prior to use.
2. Use one Alethia CMV Test Device for each Positive Control and Negative Control to be tested.

QUALITY CONTROL SAMPLE PREPARATION

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

NOTE: Ensure that the Alethia instrument is powered on and the required performance verifications have been completed prior to initiation of Quality Control Sample Preparation. Refer to the Alethia instrument Operator's Manual for further information regarding instrument set-up and operation.

NOTE: The Positive and Negative Controls require separate CMV Buffer I & II Tubes

NOTE: When vortexing solutions, care should be taken to avoid foaming to prevent splash, touch and aerosol contamination of target viral nucleic acids. Good laboratory practices should be adhered to when performing the alethia CMV assay.

1. Vortex the Alethia CMV Negative Control for approximately 10 seconds to mix.
2. Remove the cap from the Alethia CMV Buffer I Tube. Transfer 100 µL of Alethia CMV Negative Control to the vial and recap. Vortex tube for approximately 10 seconds to mix.
3. Incubate at 19-30 C for 2 minutes.
4. Remove the cap from the Alethia CMV Buffer II Tube. Transfer 50 µL of the processed control to the Alethia CMV Buffer II tube. Replace the Alethia CMV Buffer II cap.
5. Vortex the Alethia CMV Buffer II tube for approximately 10 seconds.
6. Repeat Steps 1-5 for Alethia CMV Positive Control.

QUALITY CONTROL TEST PROCEDURE

A maximum of 10 samples can be processed in a single Alethia instrument run.

1. Remove 1 Alethia CMV Test Device from its protective pouch for the Negative Control. Carefully open the device by lifting the latch, holding the chambers such that the lyophilized reagent will not fall out upon opening. Place device on a flat surface or in a rack that can accommodate the device.
2. Using a micropipette, transfer 50 µL of the prepared Negative Control to both the TEST (White Bead) and CONTROL (Yellow Bead) chambers of the Alethia Test Device. Take care not to introduce air to the reaction mixture. Close the Alethia test device and fasten the latch securely.
3. Repeat Test Procedure Steps 1-2 for the prepared Positive Control Sample to be tested.
4. Tap device on the bench top or mix to remove air bubbles. Carefully examine the Alethia Test Devices for rehydration of the Control/Test Bead, for air bubbles left in the chamber and liquid in the top of the devices. If undissolved beads, air bubbles or liquid in the top of the devices are noted, tap the devices on the bench top and repeat visual inspection. Amplification and detection should be initiated within 10 minutes.
5. Insert the Alethia Test Devices into the Alethia instrument and initiate amplification reaction and detection. Results will be displayed at the conclusion of the run.

INTERPRETATION OF RESULTS

Sample ID	Reported Result	Interpretation
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, Alethia instrument performing correctly.
	NEGATIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, Alethia instrument performing correctly.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No Alethia Test Device in the Alethia instrument Well. OR The Alethia Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. Repeat the test using original sample.

LIMITATIONS OF THE PROCEDURE

1. The Positive Control is manufactured in an aqueous solution matrix. Although specimen matrix interference has not been observed with this assay, the aqueous matrix of the controls may not adequately control for specimen matrix effects. If the user wishes to supply controls in the sample matrix, the user is referred to the Clinical and Laboratory Standards Institute guideline EP14-A3, Evaluation of Commutability of Processed Samples: Approved Guideline – Third Edition, August 2014.

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