alethia®

Group B Streptococcus

External Control Kit

External Control Materials for use with
Alethia Group B *Streptococcus* (GBS) DNA Amplification Assay

REF 479900

IVD

R_x Only

INTENDED USE

The Alethia GBS External Control Kit contains Positive and Negative Control Reagents for use with the Alethia Group B *Streptococcus* test kit. External controls are used as part of a routine quality control program.

SUMMARY AND EXPLANATION OF THE TEST

Quality control testing is performed to detect factors such as reagent deterioration, adverse environmental or test conditions, or variance in operator performance that can cause test errors. External control reagents, such as Alethia GBS External Positive Control and Negative Control, are reagents that are not built into the test system but are tested in the same manner as patient specimens.

BIOLOGICAL PRINCIPLES

To verify that Alethia GBS DNA Amplification test kits are performing correctly; they should be tested with known positive and negative external control reagents on a periodic basis. The frequency with which a laboratory performs external controls will be affected by:

1. The proficiency level of the laboratory,
2. The laboratory's own internal requirements,
3. The requirements of the laboratory's accrediting agencies,
4. The number of new and different operators performing the test,
5. Whether a new kit lot is being added to testing,
6. Whether the kit lot is from a different shipment, and
7. Whether deviations from the manufacturer's stated storage or handling conditions have occurred.

When unacceptable quality control test results are obtained, all test results should be considered invalid. QC test failures are an indication that either reagents, test environment or operator performance have changed.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. **Alethia GBS Positive Control:** Tris-buffered solution containing non-infectious Plasmid DNA (*S. aureus* and *S. agalactiae* inserts) with azide (0.09%) as a preservative.
2. **Alethia Negative Control:** Tris-buffered solution containing non-infectious Plasmid DNA (*S. aureus* insert) with azide (0.09%) as a preservative.

MATERIALS NOT PROVIDED

1. Alethia Group B *Streptococcus* (GBS) DNA Amplification Test Kit (Product Code 480350)
2. Disposable latex gloves, powder free
3. DNase/RNase free, aerosol resistant pipette tips

EQUIPMENT NOT PROVIDED

1. Vortex Mixer
2. Interval Timer
3. Micropipette capable of dispensing 50 µL
4. Alethia Incubator/Reader

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. This is a quality control reagent and is used to evaluate the performance of the Alethia GBS DNA Amplification Assay. It is not directly used to test patient samples.
3. Transport and store these reagents at 2-8 C when not in use. Do not freeze.
4. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
5. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
6. CDC/NIH manual "Biosafety in Microbiology and Biomedical Laboratories", 2009
7. Quality Control Programs for Molecular Testing Laboratories should be employed.
8. The Alethia GBS Test Device contains lyophilized reagents. The protective pouch should not be open until ready to perform the assay.
9. The Alethia GBS Test Device includes a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
10. Dispose of used Alethia Test Devices immediately after processing, leaving the device latch securely in place. Opening the device after amplification may result in contamination of the test area with amplification product.

HAZARD and PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-8 C.

QUALITY CONTROL TEST PROCEDURE

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

1. Bring Alethia External Controls and all kit components to room temperature (21-27 C) before testing. Incorrect results may be obtained if Control Materials are not brought to room temperature prior to use.
2. Use 1 Alethia GBS Test Device for each Positive Control and Negative Control reagent to be tested.
3. Transfer 50 µL of Negative Control to an appropriately labeled Alethia Reaction Buffer tube.
4. Vortex the Reaction Buffer Tube containing Negative Control for approximately 10 seconds.
5. Remove 1 Alethia GBS Test Device from its protective pouch for the Negative Control. Carefully open the device by lifting the latch, holding the tubes such that the lyophilized reagent will not fall out upon opening. Place device on a flat surface or in a rack that can accommodate the device.
6. Using a micropipette, first transfer 50 µL from the Reaction Buffer tube containing **Negative Control** to the TEST chamber (left/White Bead), and then transfer 50 µL from the Reaction Buffer tube containing Negative Control to the CONTROL chamber (right/Yellow Bead) of the Alethia Test Device. Take care not to introduce air to the reaction mixture. Do not mix reactions with pipette. Close the Alethia Test Device and fasten the latch securely.
7. Transfer 50 µL of Positive Control to an appropriately labeled Alethia Reaction Buffer tube.
8. Vortex the Reaction Buffer Tube containing Positive Control for approximately 10 seconds.
9. Remove 1 Alethia GBS Test Device from its protective pouch for the Positive Control. Carefully open the device by lifting the latch, holding the tubes such that the lyophilized reagent will not fall out upon opening. Place device on a flat surface or in a rack that can accommodate the device.
10. Using a micropipette, first transfer 50 µL from the Reaction Buffer tube containing **Positive Control** to the TEST chamber (left/White Bead), and then transfer 50 µL from the Reaction Buffer tube containing Positive Control to the CONTROL chamber (right/Yellow Bead) of the Alethia Test Device. Take care not to introduce air to the reaction mixture. Do not mix reactions with pipette. Close the Alethia Test Device and fasten the latch securely.
11. Gently mix each device by tapping or flicking, being careful to not introduce air bubbles or get the reaction mixture into the device lid. Carefully examine the test device for air bubbles or liquid in the lid. If air bubbles or liquid is noted, tap the device on the bench top and repeat the visual inspection.
12. Insert each Alethia Test Device into the Alethia Incubator/Reader and initiate amplification reaction and detection. Results will be displayed at the conclusion of the run.

INTERPRETATION OF RESULTS

Sample ID	Reported Result	Interpretation
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, Alethia Incubator/Reader performing correctly.
	NEGATIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, Alethia Incubator/Reader performing correctly.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No Alethia Test Device in the Alethia Incubator/Reader. OR The Alethia Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. Repeat the test using original sample.

LIMITATIONS OF THE PROCEDURE

The Positive Control is manufactured in an aqueous solution matrix. Although specimen matrix interference has not been observed with this assay, the aqueous matrix of the controls may not adequately control for specimen matrix effects. If the user wishes to supply controls in the sample matrix, the user is referred to the Clinical and Laboratory Standards Institute guideline EP14-A2, Evaluation of Matrix Effects: Approved Guideline – second edition, January 2005.

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