Group A Streptococcus External Control Kit

External Control Materials for use with
Alethia Group A Streptococcus (Group A Strep) DNA Amplification Assay

REF 479910

IVD

Rx Only

INTENDED USE
The Alethia Group A Streptococcus (Group A Strep) External Control Kit contains Positive Reagent for use with the Alethia Group A Streptococcus DNA Amplification Assay. The External Positive Control Reagent is used in conjunction with the Alethia Sample Preparation Apparatus II/Negative Control III reagent provided with the Group A Strep Kit as part of routine Quality Control Testing. External controls are used to aid the user in detecting unexpected conditions that may lead to test errors.

SUMMARY AND EXPLANATION OF THE TEST

Quality control testing is performed to detect factors such as reagent deterioration, adverse environmental or test conditions, or variance in operator performance that can cause test errors. External control reagents, such as the Alethia Group A Strep External Positive Control, are reagents that are not built into the test system but are amplified and detected in the same manner as patient specimens.

BIOLOGICAL PRINCIPLES

To verify that Alethia Group A Streptococcus DNA Amplification test kits are performing correctly; they should be tested with known positive and negative external control reagents on a periodic basis. The frequency with which a laboratory performs external controls will be affected by:

1. The proficiency level of the laboratory,
2. The laboratory's own internal requirements,
3. The requirements of the laboratory's accrediting agencies,
4. The number of new and different operators performing the test,
5. Whether a new kit lot is being added to testing,
6. Whether the kit lot is from a different shipment, and
7. Whether deviations from the manufacturer's stated storage or handling conditions have occurred.

When unacceptable quality control test results are obtained, all test results should be considered invalid. QC test failures are an indication that either reagents, test environment or operator performance have changed.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. Alethia Group A Streptococcus Positive Control: Tris-buffered solution containing non-infectious plasmid DNA (*S. aureus* and *S. pyogenes* target DNA inserts) with azide (0.09%) as a preservative.

MATERIALS NOT PROVIDED

1. Alethia Group A Streptococcus (Group A Strep) DNA Amplification Assay (Product Code 480150)
2. Disposable latex gloves, powder free
3. DNase/RNase free, aerosol resistant pipette tips

EQUIPMENT NOT PROVIDED

1. Dry-bath with 12 mm heat block capable of 95 C
2. Digital Thermometer with Max/Min Temperature Memory (e.g. Traceable® Lollipop™ Waterproof/Shockproof Thermometer)
3. Vortex Mixer
4. Interval Timer
5. Micropipette capable of dispensing 50 μ L
6. Alethia reader Catalog Number: 610189

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. This is a quality control reagent and is used to evaluate the performance of the Alethia Group A Strep DNA Amplification Assay. It is not directly used to test patient samples.
3. Transport and store these reagents at 2-8 C when not in use. Do not freeze.
4. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
5. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
6. Quality Control Programs for Molecular Testing Laboratories should be employed.
7. The Alethia Group A Strep Test Device contains lyophilized reagents. The protective pouch should not be open until ready to perform the assay.
8. The Alethia Group A Strep Test Device includes a latch feature that is designed to prevent contamination of the test area with amplification product. DO NOT use Test Devices with broken latches.
9. Dispose of used Alethia Test Devices immediately after processing, leaving the device latch securely in place. Opening the device after amplification may result in contamination of the test area with amplification product.

HAZARDS and PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-8 C.

REAGENT PREPARATION

1. Bring Alethia Group A Strep Positive Control and all kit components to room temperature (21-27 C) before use. Incorrect results may be obtained if Control Materials and components are not brought to room temperature prior to use. External Controls brought to room temperature are stable for up to 8 hours.
2. Use one Alethia Group A Strep Test Device for each Positive Control and Negative Control to be tested.

QUALITY CONTROL SAMPLE PREPARATION

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

NOTE: Ensure that the Alethia reader is powered on and required performance verifications have been completed prior to initiation of Quality Control Sample Preparation. Refer to the Alethia incubator/reader Operator's Manual for further information regarding instrument set-up and operation.

1. Label 1 Alethia Heat Treatment Tube for each Control to be tested.
2. **For Negative Control:** Vortex the Sample Preparation Apparatus II/Negative Control III (**SMP PREP**) for a minimum of 10 seconds. Remove the tip cap from the Alethia **SMP PREP** and squeeze 5 to 10 drops of sample into a clean Alethia Heat Treatment Tube.
3. **For Positive Control:**
 - Add 50 μ L of Alethia Group A Strep Positive Control to the **Alethia SMP PREP**. Replace and secure the **Alethia SMP PREP** cap. Vortex the Sample Preparation Apparatus II/Negative Control III with Positive Control for a minimum of 10 seconds.
 - Remove the tip cap from the **Alethia SMP PREP**, and squeeze five to ten drops of sample into a clean Alethia Heat Treatment Tube.
4. Heat each Alethia Heat Treatment Tube in a dry-bath/heat block at 95 \pm 5 C for 10 \pm 2 minutes. Monitor heat-treatment step with digital thermometer and interval timer.
5. Remove each Heat Treatment Tube from the dry-bath/heat block and proceed to Test Procedure.

QUALITY CONTROL TEST PROCEDURE

A maximum of 10 samples can be processed in a single Alethia reader run.

1. Vortex heat-treated samples for approximately 10 seconds.
2. Remove 1 Alethia Group A Strep Test Device from its protective pouch for the Negative Control. Carefully open the device by lifting the latch, holding the tubes such that the lyophilized reagent will not fall out upon opening. Place device on a flat surface or in a rack that can accommodate the device.
3. Using a micropipette, first transfer 50 μ L of the heat-treated Negative Control sample to the TEST chamber (left/White Bead), and then transfer 50 μ L of the heat treated Negative Control to the CONTROL chamber (right/Yellow Bead) of the Alethia Test Device. Take care not to introduce air to the reaction mixture. Do not mix reactions with pipette.
4. Remove 1 Alethia Group A Strep Test Device from its protective pouch for the Positive Control. Carefully open the device by lifting the latch, holding the tubes such that the lyophilized reagent will not fall out upon opening. Place device on a flat surface or in a rack that can accommodate the device.
5. Using a micropipette, first transfer 50 μ L of the heat-treated Positive Control to the TEST chamber (left/White Bead), and then transfer 50 μ L of the heat treated Positive Control to the CONTROL chamber (right/Yellow Bead) of the Alethia Test Device. Take care not to introduce air to the reaction mixture. Do not mix reactions with pipette. Close the Alethia Test Device and fasten the latch securely.
6. Tap device on the bench top or mix to remove air bubbles. Carefully examine the test device for dissolution of the Control/Test Bead, for air bubbles left in the tube and liquid in the top of the device. If undissolved beads, air bubbles or liquid is noted, tap the device on the bench top and repeat the visual inspection.
7. Insert each Alethia Test Device into the Alethia reader and initiate amplification reaction and detection. Results will be displayed at the conclusion of the run.

INTERPRETATION OF RESULTS

Sample ID	Reported Result	Interpretation
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, Alethia reader performing correctly.
	NEGATIVE	Incorrect control result. Repeat control testing as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services at 1-800-343-3858 (US) or your local distributor.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	Incorrect control result. Repeat control testing as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, Alethia reader performing correctly.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No Alethia Test Device in the Alethia reader. OR The Alethia Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. Repeat the test using original sample.

LIMITATIONS OF THE PROCEDURE

The Positive Control is manufactured in an aqueous solution matrix. Although specimen matrix interference has not been observed with this assay, the aqueous matrix of the controls may not adequately control for specimen matrix effects. If the user wishes to supply controls in the sample matrix, the user is referred to the Clinical and Laboratory Standards Institute guideline EP14-A2, Evaluation of Matrix Effects: Approved Guideline – second edition, January 2005.

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Meridian Bioscience Europe B.V.