alethia[®]

HSV 1&2

External Control Kit

External Control Materials for use with
Alethia HSV 1&2 DNA Amplification Assay

REF 479960

IVD

R_x Only

INTENDED USE

The Alethia HSV 1&2 External Control Kit contains Positive and Negative Control Reagents for use with the Alethia HSV 1&2 DNA Amplification Assay. External controls are used as part of a routine quality control program to aid the user in detection of unexpected conditions that may lead to test errors.

SUMMARY AND EXPLANATION OF THE TEST

Quality control testing is performed to detect factors such as reagent deterioration, adverse environmental or test conditions, or variance in operator performance that can cause test errors. External control reagents, such as the Alethia HSV 1&2 External Positive Control and Alethia HSV 1&2 External Negative Control, are reagents that are not built into the test system but are amplified and detected in the same manner as patient specimens.

BIOLOGICAL PRINCIPLES

The Alethia HSV 1&2 DNA Amplification test kit should be tested with known positive and negative external control reagents on a periodic basis to verify that the kits are performing correctly. The frequency with which a laboratory performs external controls will be affected by:

1. The proficiency level of the laboratory,
2. The laboratory's own internal requirements,
3. The requirements of the laboratory's accrediting agencies,
4. The number of new and different operators performing the test,
5. Whether a new kit lot is being added to testing,
6. Whether the kit lot is from a different shipment, and
7. Whether deviations from the manufacturer's stated storage or handling conditions have occurred.

When unacceptable quality control test results are obtained, all test results should be considered invalid. QC test failures are an indication that either reagents, test environment or operator performance have changed.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. **Alethia HSV 1&2 Positive Control:** Tris-buffered solution with plasmid DNA containing target inserts for herpes simplex virus types 1 & 2 and sodium azide 0.09% as a preservative.
2. **Alethia Negative Control IV:** Tris-buffered solution with sodium azide 0.09% as a preservative.

MATERIALS PROVIDED SEPARATELY

1. Alethia HSV 1&2 DNA Amplification Test Kit, Meridian Bioscience, Inc. (Catalog Number 480650)

MATERIALS NOT PROVIDED

1. Disposable latex gloves, powder free
2. DNase/RNase-free, aerosol resistant pipette tips

EQUIPMENT NOT PROVIDED

1. Vortex mixer
2. Dry-bath with 12 mm heat block capable of 95 C
3. Digital thermometer with Max/Min temperature memory (eg, Traceable® Lollipop™ Waterproof/Shockproof Thermometer)
4. Interval timer
5. Micropipette capable of dispensing 50 µL
6. Alethia Reader Meridian Bioscience, Inc. Catalog Number: 610189

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. This is a quality control reagent and is used to evaluate the performance of the Alethia HSV 1&2 DNA Amplification Assay. It is not directly used to test patient samples.
3. Transport and store these reagents at 2-8 C when not in use. Do not freeze.
4. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
5. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
6. Quality Control Programs for Molecular Testing Laboratories should be employed.
7. The Alethia HSV 1&2 Test Devices contain lyophilized reagents. The protective pouch should not be opened until ready to perform the assay.
8. The Alethia HSV 1&2 Test Devices include a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
9. Dispose of used Alethia Test Devices immediately after processing, leaving the device latch securely in place. Opening the device after amplification may result in contamination of the test area with amplification product.

HAZARD AND PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-8 C.

REAGENT PREPARATION

1. Bring Alethia HSV 1&2 External Controls and all Alethia kit components to room temperature (19-30 C) before use. Incorrect results may be obtained if Control Materials and components are not brought to room temperature prior to use.
2. Use one Alethia HSV 1 and one HSV 2 Test Device for each Positive Control and Negative Control to be tested.

QUALITY CONTROL SAMPLE PREPARATION

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

NOTE: Ensure that the Alethia reader is powered on and required performance verifications have been completed prior to initiation of QUALITY CONTROL SAMPLE PREPARATION. Refer to the Alethia reader Operator's Manual for further information regarding instrument set-up and operation.

For Negative Control:

1. Vortex the Alethia Negative Control IV for approximately 10 seconds.
2. Using a micropipette, transfer 50 µL of Alethia Negative Control IV into one Alethia Sample Preparation Apparatus III (SMP PREP III).
3. Replace SMP PREP III cap and vortex for approximately 10 seconds.
4. Remove the tip cap from the SMP PREP III and squeeze the entire contents of the SMP PREP III into a Heat Treatment Tube.
5. Heat each Heat Treatment Tube containing Negative Control in a dry-bath/heat block at 95 ± 5 C for 10 ± 2 minutes. Monitor heat-treatment step with digital thermometer and interval timer.
6. Remove each Heat Treatment Tube from the dry-bath/heat. Heat-treated samples may be held at room temperature (19-30 C) for up to 15 minutes prior to testing.

For Positive Control:

1. Vortex the Alethia HSV 1&2 Positive Control for approximately 10 seconds.
2. Using a micropipette, transfer 50 µL of Alethia HSV 1&2 Positive Control into one SMP PREP III.
3. Replace the SMP PREP III cap and vortex for approximately 10 seconds.
4. Remove the tip cap from the SMP PREP III and squeeze the entire contents of the SMP PREP III into a Heat Treatment Tube.
5. Heat each Heat Treatment Tube containing Positive Control in a dry-bath/heat block at 95 ± 5 C for 10 ± 2 minutes. Monitor heat-treatment step with digital thermometer and interval timer.
6. Remove each Heat Treatment Tube from the dry-bath/heat. Heat-treated samples may be held at room temperature (19-30 C) for up to 15 minutes prior to testing.

QUALITY CONTROL TEST PROCEDURE

NOTE: A maximum of 5 samples can be processed for both HSV 1 and HSV 2 in a single Alethia reader run. **HSV 2 Test Devices are identified by an orange band on the Test Device closure tab.**

- Vortex heat-treated samples for approximately 10 seconds.
- Remove 1 Alethia HSV 1 or 1 Alethia HSV 2 Test Device from its protective pouch for the **Negative Control**. Carefully open the device by lifting the latch, holding the device upright so that the lyophilized reagents do not fall out upon opening. Place device on a flat surface or in a rack that can accommodate the device.
- Using a micropipette, transfer 50 µL of the prepared Negative Control sample to both the TEST (Left/White Bead) and CONTROL (Right/Yellow Bead) chambers of the Alethia Test Device. Do not introduce extraneous air to the reaction mixture. Do not mix reactions with pipette.
- Add 1 drop of Mineral Oil to both the TEST and CONTROL chambers of the Alethia Test Device. Close the Alethia Test Device and fasten the latch securely.
- Repeat steps 2 through 4 for remaining **Negative Control** samples to be tested.
- Remove 1 Alethia HSV 1 or 1 Alethia HSV 2 Test Device from its protective pouch for the **Positive Control**. Carefully open the device by lifting the latch, holding the device upright so that the lyophilized reagents do not fall out upon opening. Place device on a flat surface or in a rack that can accommodate the device.
- Using a micropipette, transfer 50 µL of the prepared Positive Control sample to both the TEST (Left/White Bead) and CONTROL (Right/Yellow Bead) chambers of the Alethia Test Device. Do not introduce extraneous air to the reaction mixture. Do not mix reactions with pipette.
- Add 1 drop of Alethia Mineral Oil to both the TEST and CONTROL chambers of the Alethia Test Device. Close the Alethia Test Device and fasten the latch securely.
- Repeat steps 6 through 9 for remaining **Positive Control** samples to be tested.
- Tap device on the bench top or mix to remove air bubbles. Carefully examine the Test Device for rehydration of the Control/Test Bead, for air bubbles left in the chamber and liquid in the top of the device. If undissolved beads, air bubbles or liquid in the top of the device are noted, tap the device on the bench top and repeat visual inspection. Amplification and detection should be initiated within 15 minutes.
- Insert each Alethia Test Device into the Alethia reader and initiate run. Results will be displayed at the conclusion of the run.

INTERPRETATION OF RESULTS

Sample ID	Reported Result	Interpretation
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, Alethia reader performing correctly.
	NEGATIVE	Incorrect control result. Repeat control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	Incorrect control result. Repeat control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, Alethia reader performing correctly.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No Alethia Test Device in the Alethia reader Well. OR The Alethia Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. Repeat the test using original sample.

LIMITATIONS OF THE PROCEDURE

The Positive Control is manufactured in an aqueous solution matrix. Although specimen matrix interference has not been observed with this assay, the aqueous matrix of the controls may not adequately control for specimen matrix effects. If the user wishes to supply controls in the sample matrix, the user is referred to the Clinical and Laboratory Standards Institute guideline EP14-A2, Evaluation of Matrix Effects: Approved Guideline – second edition, January 2005.

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