alethia[®]

Malaria

External Control Kit

External Control Material for use with
Alethia Malaria DNA Amplification Assays

REF 479970

IVD

R_x Only

INTENDED USE

The Alethia Malaria External Control Kit contains Positive and Negative Control Reagents for use with the Alethia Malaria DNA Amplification Assays. External controls are used as part of a routine quality control program to aid the user in detection of unexpected conditions that may lead to test errors.

SUMMARY AND EXPLANATION OF THE TEST

Quality control testing is performed to detect factors such as reagent deterioration, adverse environmental or test conditions, or variance in operator performance that can cause test errors. External control materials are not built into the test system but are amplified and detected in the same manner as patient specimens.

BIOLOGICAL PRINCIPLES

The Alethia Malaria DNA Amplification test kits should be tested with known positive and negative external control reagents on a periodic basis to verify that the kits are performing correctly. The frequency with which a laboratory performs external controls will be affected by:

1. The proficiency level of the laboratory,
2. The laboratory's own internal requirements,
3. The requirements of the laboratory's accrediting agencies,
4. The number of new and different operators performing the test,
5. Whether a new kit lot is being added to testing,
6. Whether the kit lot is from a different shipment, and
7. Whether deviations from the manufacturer's stated storage or handling conditions have occurred.

When unacceptable quality control test results are obtained, all test results should be considered invalid. QC test failures are an indication that either reagents, test environment or operator performance have changed.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. **Alethia Malaria Positive Control:** Tris-buffered solution with plasmid containing DNA inserts (*Plasmodium sp.* and human mitochondrial DNA inserts) and azide (0.09%) as a preservative.
2. **Alethia Malaria Negative Control:** Tris-buffered solution with plasmid containing human mitochondrial DNA inserts) and azide (0.09%) as a preservative.

MATERIALS PROVIDED SEPARATELY

1. Alethia Malaria DNA Amplification Test Kit, Meridian Bioscience, Inc. Catalog Number: 480925

MATERIALS NOT PROVIDED

1. Disposable latex gloves, powder free
2. DNase/RNase free, aerosol resistant pipette tips

EQUIPMENT NOT PROVIDED

1. Interval timer
2. Vortex Mixer (optional)
3. Micropipette capable of dispensing 50 µL
4. Alethia Incubator/Reader, Meridian Bioscience, Inc. Catalog Number: 610189

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. This is a quality control reagent and is used to evaluate the performance of the Alethia Malaria DNA Amplification Assays. It is not directly used to test patient samples.
3. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
4. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
5. Quality Control Programs for Molecular Testing Laboratories should be employed.
6. The Alethia Malaria Test Devices contain lyophilized reagents. The protective pouch should not be opened until ready to perform the assay.
7. The Alethia Malaria Test Devices include a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
8. Dispose of used Alethia Test Devices and tubes immediately after processing, leaving the device latch securely in place. Opening the device after amplification may result in contamination of the test area with amplification product.

HAZARD AND PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-30 C.

REAGENT PREPARATION

1. Bring Alethia Malaria External Positive and Negative Controls, and all Alethia Malaria kit components to room temperature (19-30 C) before use. Incorrect results may be obtained if control materials and components are not brought to room temperature prior to use.
2. Use one Alethia Malaria Test Device for each Positive Control and Negative Control to be tested.

QUALITY CONTROL SAMPLE PREPARATION

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

NOTE: Ensure that the Alethia Incubator/Reader is powered and the required performance verifications have been completed prior to initiation of Quality Control Sample Preparation. Refer to the Alethia Incubator/Reader Operator's Manual for further information regarding instrument set-up and operation.

1. **Preparation of controls using Alethia Malaria**
Note: The Positive and Negative Control require separate Buffer tubes and SMP PREP IV.
 - a. Mix Positive Control by inversion 5 times or by vortexing 10 seconds. Using a micropipette, add 50 µL of the appropriate Positive or Negative Control to one Buffer I tube. Mix by inversion 5 times or by vortexing for approximately 10 seconds. Incubate the sample for 2 minutes.
 - b. Mix by inversion 5 times or by vortexing for approximately 10 seconds. Using a micropipette, immediately transfer 50 µL of prepared sample into SMP PREP IV. Mix by inversion 5 times or vortex for 10 seconds.
 - c. Gently squeeze the SMP PREP IV and slowly collect 5 to 10 drops into a clean Tube I. Proceed to Quality Control Test Procedure.

QUALITY CONTROL TEST PROCEDURE

NOTE: A maximum of 10 samples can be processed in a single Alethia Incubator/Reader run.

1. Remove 1 Alethia Malaria Test Device from its protective pouch for the Negative Control. Carefully open the device by lifting the latch, holding the chambers such that the lyophilized reagents will not fall out upon opening. Place the device on a flat surface or in a rack that can accommodate the device.
2. Using a micropipette, transfer 50 µL of the prepared Negative Control to both the TEST (Left/White Bead) and CONTROL (Right/Yellow Bead) chambers of the Alethia Malaria Test. Take care to not introduce air to the reaction mixture. Do not mix reactions with pipette.
3. Close the Alethia Test Device and fasten the latches securely.
4. Repeat Test Procedure Steps 1-3 for the prepared Positive Control sample to be tested.
5. Tap device(s) on the bench top or mix to remove air bubbles. Carefully examine the Test Device(s) for rehydration of the Control/Test Bead, for air bubbles left in the chamber and liquid in the top of the device. If undissolved beads, air bubbles or liquid in the top of the device are noted, tap the device on the bench top and repeat visual inspection. Amplification and detection should be initiated within 15 minutes.
6. Insert the Alethia Test Devices into the Alethia Incubator/Reader and initiate run **using the Malaria Program**. Results will be displayed at the conclusion of the run.

INTERPRETATION OF RESULTS

Sample ID	Reported Result	Interpretation
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, Alethia Incubator/Reader performing correctly.
	NEGATIVE	Incorrect control result. Repeat control testing as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	Incorrect control result. Repeat control testing as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, Alethia Incubator/Reader performing correctly.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No Alethia Test Device in the Alethia Incubator/Reader Well. OR The Alethia Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. Repeat the test using original sample.

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