# alethia™

## Mycoplasma Direct DNA Amplification Assay

### DNA Amplification Assay for the Detection of *Mycoplasma pneumoniae* in human throat swab specimens

REF 480250

IVD

Rx Only

#### INTENDED USE

The Alethia Mycoplasma Direct DNA amplification assay, performed on the Alethia Reader, is a qualitative in vitro diagnostic test for the direct detection of DNA from *Mycoplasma pneumoniae* in human throat swabs obtained from patients suspected of having *Mycoplasma pneumoniae* infection.

The Alethia Mycoplasma Direct assay utilizes loop-mediated isothermal DNA amplification (LAMP) technology to detect *Mycoplasma pneumoniae* by targeting a segment of the *Mycoplasma pneumoniae* genome.

Results from the Alethia Mycoplasma Direct DNA amplification assay should be used in conjunction with clinical presentation, other laboratory findings, and epidemiological risk factors as an aid in the diagnosis of *Mycoplasma* infection and should not be used as the sole basis for treatment or other patient management. Positive results do not rule out co-infection with other organisms and negative results in persons with respiratory tract infections may be due to pathogens not detected by this assay. Lower respiratory tract infections due to *M. pneumoniae* may not be detected by this assay. If lower respiratory tract infection due to *M. pneumoniae* is suspected, additional laboratory testing using methods other than the Alethia Mycoplasma Direct DNA Amplification Assay may be necessary.

#### SUMMARY AND EXPLANATION OF THE TEST

The Alethia Mycoplasma Direct DNA Amplification Assay is based on loop-mediated amplification (LAMP) technology.<sup>1,2</sup> The assay targets a 208 base pair (bp) sequence of the *Mycoplasma pneumoniae* genome. The target DNA sequence is found in the intracellular protease-like protein gene and is found in all sequences comprising twenty genomes of *Mycoplasma pneumoniae*.

Loop-mediated amplification uses specially designed primers to provide for specific and continuous isothermal DNA amplification. A by-product of amplification is the formation of magnesium pyrophosphate, which forms a white precipitate leading to a turbid reaction solution. Reaction solution absorbance characteristics are monitored by the Meridian Alethia Incubator/Reader. Changes in reaction solution absorbance characteristics created by precipitation of magnesium pyrophosphate indicate the presence of target DNA. The absence of target DNA results in no significant change in sample absorbance.

The Alethia Mycoplasma Direct kit includes Alethia Sample Preparation Apparatus II/Negative Control III (SMP PREP), Alethia Mycoplasma Direct Test Devices, and Heat Treatment Tubes. The Alethia Sample Preparation Apparatus II, used for specimen dilution and preparation, is filled with a buffered solution containing formalin treated *E. coli* harboring *Staphylococcus aureus* DNA. The Alethia Mycoplasma Direct Test Device contains one lyophilized amplification reagent bead in each of two chambers: a TEST chamber with *Mycoplasma*-specific primers and a CONTROL chamber with *S. aureus*-specific primers. The *S. aureus* DNA from the Control Reagent and the *S. aureus*-specific primers in the CONTROL chamber function as the Internal Control for the assay. During specimen preparation, each patient specimen is added to the Sample Preparation Apparatus and combined with the *S. aureus* DNA prior to amplification. Addition of *S. aureus* DNA to the patient specimen allows for parallel processing of target DNA and Control DNA through sample preparation, amplification and detection. The Internal Control monitors amplification inhibition, assay reagent performance and sample processing effectiveness. The control *S. aureus* target must be amplified and detected in the final reaction or the test is considered invalid and patient results are not reported.

The Alethia Reader monitors changes in absorbance characteristics by measuring transmission of light through the Test and Control reaction solutions. Light transmission is checked at the assay Run Start (Signal<sub>initial</sub>, S) and at the assay Run End (Signal<sub>final</sub>, S). The Alethia Reader calculates the change in light transmission between Run End and Run Start (S<sub>r</sub>:S<sub>i</sub>) and compares the ratio to a fixed cut-off value.

Fixed cut-off values for the TEST chamber are used to report sample results. TEST chamber S<sub>r</sub>:S<sub>i</sub> ratios less than 82% are reported as 'POSITIVE'; TEST chamber S<sub>r</sub>:S<sub>i</sub> ratios greater than or equal to 82% are reported as 'NEGATIVE'. Numerical values are not reported.

Fixed cut-off values for the CONTROL chamber are used to determine validity. CONTROL chamber S<sub>r</sub>:S<sub>i</sub> ratios less than 90% are considered valid and allow for reporting of TEST chamber results (POSITIVE, NEGATIVE). CONTROL chamber S<sub>r</sub>:S<sub>i</sub> ratios greater than or equal to 90% are considered invalid and prevent reporting of TEST chamber results. Invalid CONTROL chamber reactions are reported as 'INVALID'. Numerical values are not reported.

More stringent cut-off criteria are applied to the CONTROL chamber reaction to ensure amplification is not inhibited, reagents are performing as intended and that sample processing was performed appropriately.

#### BIOLOGICAL PRINCIPLES

*Mycoplasma pneumoniae* is a common cause of human upper and lower respiratory infections, including pharyngitis, acute bronchitis and pneumonia.<sup>3</sup> *Mycoplasma* became recognized as a human pathogen in the 1960s, with *Mycoplasma pneumoniae* the best known and most studied.<sup>4</sup> *M. pneumoniae* contains a triple-layered cell membrane rather than a structured cell wall,<sup>4</sup> which prevents identification by traditional gram stain methods and visualization by light microscopy.

*M. pneumoniae* has been associated with up to 40%<sup>3-4</sup> of community-acquired pneumonia cases. Infection occurs in both children and adults without geographical, gender or climate-related restrictions.<sup>4,7</sup> *M. pneumoniae* is most often associated with atypical pneumonia, presenting with symptoms that include headache, malaise, myalgia, fever, and sore throat accompanied by dry, paroxysmal cough.<sup>4,5,7</sup> While the clinical course of mycoplasmal pneumonia is generally mild and self-limiting, it has been associated with a mortality rate of approximately 1.4%.<sup>5</sup> An estimated 2 million cases of *M. pneumoniae* infection occur annually with approximately 100,000 pneumonia-related hospitalizations in the United States each year.<sup>9</sup> It is estimated as many as 18% of children requiring hospitalization with community-acquired pneumonia are caused by *M. pneumoniae*.<sup>4</sup>

Transmission of *M. pneumoniae* is generally from person to person by aerosol with a reported incubation period of one to four weeks.<sup>6</sup> As patients with active infection carry *Mycoplasma* in the nose, throat, trachea, and sputum, the spread of the disease is facilitated by its accompanying cough.<sup>4</sup> This mode of transmission lends itself to outbreaks in close personal contact settings such as schools, military barracks, businesses, summer camps, colleges or institutions. Similarly, the disease is commonly spread among family members within a household.<sup>4-7</sup>

Diagnosis of acute *M. pneumoniae* infection is difficult. Basic diagnostic strategies in clinical practice include culture and serology. *M. pneumoniae* culture is often impractical for patient management as the organism may take as many as six weeks to culture.<sup>3</sup> Serology products are widely available and offer practical diagnostic solutions, however, may be limited in their ability to detect acute infection. Serological assays are often directed at specific IgM and IgG antibodies. IgM antibodies are generally not detectable within the first seven days of symptom onset and may persist for months after active infection.<sup>4</sup> Studies indicate that less than 50% of patients with acute infection demonstrate a positive IgG response.<sup>4</sup>

New diagnostic techniques, including DNA-amplification methods, may enable more rapid diagnosis and earlier identification of outbreaks for prevention of secondary cases through implementation of control measures.<sup>6</sup>

#### REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. **Alethia Sample Preparation Apparatus II/Negative Control III (SMP PREP):** Tris-buffered solution containing formalin-treated *E. coli* harboring *S. aureus* DNA and sodium azide (0.09%) as a preservative.
2. **Alethia Mycoplasma Test Device:** Two-chambered device containing lyophilized amplification reagents (DNA polymerase, deoxynucleotide triphosphates) and either *Mycoplasma*-specific primers (TEST Chamber) or control primers (CONTROL Chamber).
3. **Alethia Heat Treatment Tubes**

#### MATERIALS PROVIDED SEPARATELY

1. Alethia Mycoplasma Direct External Control Kit, Catalog Number: 479890

#### MATERIALS NOT PROVIDED

1. Disposable latex gloves, powder free
2. DNase/RNase-free, aerosol resistant pipette tips
3. Specimen collection and transport materials:
  - **Swabs with a breakable, plastic shaft (single or double swab):** Polyester, Flocked Nylon, or Rayon
  - **Non-nutritive Transport Medium on Pledge/Sponge** (Maximum Volume: 1.2 mL): Liquid Amies (without charcoal) or Liquid Stuart

#### EQUIPMENT NOT PROVIDED

1. Dry-bath with 12 mm heat block capable of 95 C
2. Digital thermometer with max/min temperature memory (e.g., Traceable® Lollipop™ waterproof/shockproof thermometer)
3. Vortex mixer
4. Interval timer
5. Micropipette(s) capable of dispensing 50 µL
6. Alethia Reader, Meridian Bioscience, Inc. Catalog Number: 610189

#### PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Mini-tip swabs (absorbance capacity less than 60 µL) are not acceptable for Alethia Mycoplasma Direct specimen collection.
3. Do not interchange the Sample Preparation Apparatus or Test Devices between lots. Heat treatment tubes are interchangeable provided they are within assigned expiration dates when used.
4. Follow Biosafety Level 2 and Good Laboratory practices during testing.<sup>8</sup> Treat all specimens and used Test Devices as capable of transmitting infectious agents. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
5. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards. Change gloves often.
6. Quality Control Programs for Molecular Testing Laboratories should be employed.<sup>9</sup> To avoid risk of false positives due to organism/amplicon contamination, employ prevention measures (eg, clean work area with suitable cleaner that destroys nucleic acid, use only aerosol resistant pipette tips, handle one specimen at a time) as required by Good Molecular Practices.
7. The Alethia Mycoplasma Test Device contains lyophilized reagents. The protective pouch should not be opened until ready to perform the assay.
8. The Alethia Mycoplasma Test Device includes a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
9. Dispose of used Alethia Test Devices immediately after processing, leaving the device latch securely in place. Do NOT open the Test Device after processing. Opening the device after amplification may result in contamination of the test area with amplification product.

#### HAZARD AND PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

#### SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit components at the temperature indicated on the label.

#### REAGENT PREPARATION

Ensure kit reagents are at room temperature (19-29 C) before use. Incorrect results may be obtained if reagents are not brought to room temperature prior to use.

#### SPECIMEN COLLECTION AND PREPARATION

**Sample type:** Throat swab specimens.

**Sample Collection:** Sample collection should be performed in accordance with institutional guidelines for collection of throat swab specimens. Samples should be collected by vigorously swabbing the tonsils and the posterior pharynx. Swab sample(s) should be collected with suitable swab types (e.g. Rayon, Flocked Nylon, or Polyester). **NOTE: Mini-tip swabs (absorbance capacity less than 60 µL) are not acceptable for specimen collection.**

Place swab(s) in a non-nutritive transport medium (e.g. Liquid Amies, without charcoal; Liquid Stuart) containing pledget and transport to the laboratory. Samples should be held at 2-29 C during transport.

Samples may be held at room temperature (19-29 C) for up to 24 hours prior to testing. When testing will not be initiated within this time, the sample may be stored at 2-8 C for up to 14 days prior to testing.

**SPECIMEN PREPARATION**

**NOTE:** Ensure that the Alethia Reader is powered on and required performance verifications have been completed prior to initiation of SPECIMEN PREPARATION. Refer to the Alethia Operator’s Manual for further information regarding instrument set-up and operation.

1. Add single throat swab to the Alethia SMP PREP. Break off the handle of the swab. Replace and secure SMP PREP cap. Vortex the SMP PREP for a minimum of 10 seconds.
2. Remove the tip cap from the SMP PREP and squeeze 5 - 10 drops of sample into a clean Alethia Heat Treatment Tube.
3. Repeat Sample Preparation Steps for all samples to be processed.
4. Heat each Sample/Control mixture in a dry-bath/heat block at 95 ± 5 C for 10 ± 2 minutes. Monitor heat-treatment step with digital thermometer and interval timer.
5. Remove each Heat Treatment Tube from the dry-bath/heat block. Heat-treated samples may be held at 19-29 C for up to 3 hours for repeat testing. **NOTE: Heat treatment must be re-performed (95 ± 5 C for 10 ± 2 minutes) according to Step 4 if heat-treated samples are stored prior to testing.**

**TEST PROCEDURE**

**NOTE:** A maximum of 10 samples can be processed in a single Alethia Reader run. Refer to the Alethia Operator’s Manual for further information regarding instrument set-up and operation.

1. Vortex heat-treated samples for approximately 10 seconds.
2. Remove 1 Alethia Mycoplasma Test Device from its protective pouch per sample. Carefully open the device, holding the chambers such that the lyophilized reagent will not fall out upon opening. Place device on a flat surface or in a rack that can accommodate the device.
3. Using a micropipette, first transfer 50 µL of the heat-treated sample to the TEST chamber (left/White Bead), and then transfer 50 µL of the heat treated sample to the CONTROL chamber (right/Yellow Bead) of the Alethia Test Device. Take care not to introduce air to the reaction mixture. Do not mix reactions with pipette. Close the Alethia Test Device and fasten the latch securely.
4. Tap device on the bench top or mix to remove air bubbles. Carefully examine the test device for dissolution of the Control/Test Bead, for air bubbles left in the tube and liquid in the top of the device. If undissolved beads, air bubbles or liquid is noted, tap the device on the bench top and repeat the visual inspection. Amplification and detection should be initiated within 15 minutes.
5. Insert the Alethia Test Device into the Alethia Reader and initiate amplification reaction and detection. Results will be displayed at the conclusion of the run.

**INTERPRETATION OF RESULTS**

Sample ID	Reported Result	Interpretation
Patient Specimen	POSITIVE	Sample contains <i>Mycoplasma pneumoniae</i> target DNA.
	NEGATIVE	No <i>Mycoplasma pneumoniae</i> DNA detected.
	INVALID	<b>No reportable result. Repeat the test using the original sample.</b> Inhibitory patient specimen, improper sample preparation, reagent failure, instrument failure or internal control failure.
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, Alethia Reader performing correctly.
	NEGATIVE	<b>Incorrect control result. Patient results not reportable.</b> Repeat the control test as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian’s Technical Services at 1-800-343-3858 (US) or your local distributor.
	INVALID	<b>No reportable result. Repeat entire assay run using original samples.</b> Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	<b>Incorrect control result. Patient results not reportable.</b> Repeat the control test as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian’s Technical Services at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, Alethia Reader performing correctly.
	INVALID	<b>No reportable result. Repeat entire assay run using original samples.</b> Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No Alethia Test Device in the Alethia Reader Well. <b>OR</b> The Alethia Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. <b>Repeat the test using original sample.</b>

**QUALITY CONTROL**

**This test should be performed per applicable local, state, or federal regulations or accrediting agencies.**

1. Each device contains an internal control chamber that controls for amplification inhibition, assay reagents and sample processing effectiveness.
2. The heat-treatment step is monitored with an external thermometer and interval timer. Use the max/min temperature memory of the thermometer to ensure that a temperature of 95 ± 5 C is maintained. Use the interval timer to ensure that heat-treatment duration is 10 ± 2 minutes.

- Good laboratory practice recommends the use of control materials. Users should follow the appropriate federal, state and local guidelines concerning the running of external quality controls.
- Alethia Mycoplasma Direct External Positive Control Reagent is supplied separately (Catalog 479890). The Sample Preparation Apparatus II/Negative Control III reagent supplied within the Alethia Mycoplasma Direct kit (Catalog 480250) serves as the External Negative Control. It is recommended that the reactivity of each new lot and each new shipment of Alethia Mycoplasma Direct be verified on receipt or before use. External control tests should be performed thereafter in accordance with appropriate federal, state and local guidelines. The Alethia Mycoplasma Direct test kit should not be used in patient testing if the external controls do not produce the correct results.
- A separate Test Device must be used for each external control reagent.

**EXPECTED VALUES**

The incidence of *Mycoplasma pneumoniae* was established during clinical studies conducted in 2015-2016. The overall prevalence of *M. pneumoniae* in prospective throat swab specimens was 7.5% (34/456). Prevalence by patient age is provided below.

Patient Age	Total Positive	Total Samples	Prevalence (%)
0 to 1 Month	0	7	0.0%
2 Months to 2 Years	5	157	3.2%
3 Years to 12 Years	22	147	15.0%
13 Years to 21 Years	5	38	13.2%
22 Years to 65 Years	2	86	2.3%
>65 Years	0	21	0.0%

**LIMITATIONS OF THE PROCEDURE**

- This product can be used only with the Alethia instrument.
- The Alethia Mycoplasma Direct DNA assay is a qualitative assay and does not provide quantitative values or information about organism load.
- This device has not been evaluated for monitoring treatment of *Mycoplasma pneumoniae* infections.
- Respiratory infections can be caused by *Mycoplasma pneumoniae* as well as other pathogens. Positive results do not preclude coinfection with other respiratory pathogens.
- Performance of Alethia Mycoplasma Direct has been established with prospective human throat swab specimens only. Performance has not been established for other specimen types (nasopharyngeal specimens, samples from immunocompromised individuals, or from asymptomatic patients.)
- The detection of nucleic acid is dependent upon proper specimen collection, handling, transportation, storage and preparation. Failure to observe proper procedure in any one of these steps can lead to incorrect results.
- Prevalence of *M. pneumoniae* infection will affect the test's predictive value.
- Organism nucleic acids may persist *in vivo* independent of organism viability. Detection of analyte DNA target may not imply causative agents for clinical symptoms. The Alethia Mycoplasma Direct assay does not differentiate between carriers and infected individuals.
- Results should be interpreted together with clinical history, epidemiological data, and other information available to the physician. As with all molecular based tests, (A) False negative results may occur from the presence of inhibitors, technical error, sample mix-up or low numbers of organisms in the clinical specimen, (B) False positive results may occur from the presence of cross-contamination by target organisms, their nucleic acids or amplified product, and from non-specific signals.
- Nasal decongestants containing Phenylephrine HCl at concentrations greater than 0.595 mg/mL may produce false negative results with low positive specimens.
- Whole blood at concentrations greater than 2% may interfere with the Alethia Mycoplasma Direct assay.
- It has been observed that Starplex collection devices can generate a higher rate of invalid results when used in conjunction with the Alethia Mycoplasma Direct assay.

**SPECIFIC PERFORMANCE CHARACTERISTICS**

The performance characteristics of the Alethia Mycoplasma Direct DNA Amplification Assay were established in clinical studies conducted in 2015-2016 at independent clinical test sites representing three geographically distinct regions throughout the United States. Performance characteristics of the assay were compared to a reference molecular in vitro diagnostic method, the Alethia Mycoplasma DNA Amplification Assay (Meridian Bioscience, Inc., Cincinnati, OH). A total of 458 prospective, deidentified, throat swab specimens collected under informed consent from symptomatic male and female patients were evaluated. Two samples (0.4%) were excluded from the clinical sample population due to an instrument error during testing (1) or an unacceptable sample type (1) for a total of 456 eligible samples for analysis. There were no invalid results generated by the Alethia Mycoplasma Direct assay during the clinical study (0.0%; 95% CI: 0.0 – 0.8%). Performance characteristics of the Alethia Mycoplasma Direct assay are summarized below.

Alethia Mycoplasma Direct	Alethia Mycoplasma			Total
	Positive	Negative	Invalid	
Positive	24	10 <sup>a</sup>	0	34
Negative	1 <sup>b</sup>	421	0	422
Invalid	0	0	0	0
<b>Total</b>	<b>25</b>	<b>431</b>	<b>0</b>	<b>456</b>
				<b>95% CI</b>
<b>Positive Percent Agreement</b>	24/25	96.0%		80.5 - 99.3%
<b>Negative Percent Agreement</b>	421/431	97.7%		95.8 - 98.7%
<b>Overall Percent Agreement</b>	445/456	97.6%		95.7 - 98.6%
<b>Invalid Rate</b>	0/456	0.0%		0.0 - 0.8%

<sup>a</sup> 4/10 samples were identified positive by Alethia Mycoplasma after testing with an additional frozen sample.

<sup>b</sup> Repeat testing by Alethia Mycoplasma with the original patient sample and an additional frozen sample produced negative results.

The eligible study population included throat swab specimens from pediatric, adult, and geriatric patients ranging in age from 3 weeks to 97 years. There was no difference in assay performance based on patient age. The study population included 230 (50.4%) female and 226 (49.6%) male specimens. There was no difference in assay performance based on gender.

**ANALYTICAL SENSITIVITY**

The analytical sensitivity or Limit of Detection for the Alethia Mycoplasma Direct assay was determined for two *M. pneumoniae* strains (FH and M129). Limit of Detect was determined using a minimum of 20 replicates for each measurand and a stated probability (e.g., 95% where 19/20 replicates are positive) of obtaining positive responses. Analytical sensitivity testing is summarized below:

<i>Mycoplasma pneumoniae</i> Strain	CFU/mL
FH (ATCC 15531)	2350
M129 (ATCC 29342)	200

**ASSAY REACTIVITY**

The following *M. pneumoniae* strains were tested with Alethia Mycoplasma and produced positive reactions at or below stated assay limit of detect of 2350 CFU/mL: PI 1428 (ATCC 29085), MAC (ATCC 15492), M52 (ATCC 15293), Bru (ATCC 15377), M129-B170 (ATCC 29343), Mutant 22 (ATCC 39505) UAB 55612, UAB 56317, UMTB-10G (ATCC 49899). Testing demonstrated that both Type 1 and Type 2 strains are detected with the assay. The same reactivity is expected with the Alethia Mycoplasma Direct assay.

**REPRODUCIBILITY**

Blind-coded panels of 10 samples were supplied to three independent laboratories for reproducibility studies. Samples were randomly sorted within each panel to mask sample identities. The panels included contrived manufactured as low positive samples (near the limit of detect, n=3), moderate positive samples (n=3), and high negative samples (n=3). The panel also included one natural negative sample. Testing was performed by different operators at each site on the same day (intra-assay variability) for five days (inter-assay variability). Three lots of Alethia Mycoplasma Direct and eight Alethia instruments were used in the study. Positive and Negative Controls were tested with each panel. The results are provided in the table below:

Sample Type	Site 1		Site 2		Site 3		Total	
	Percent Agreement		Percent Agreement		Percent Agreement		Percent Agreement	
High Negative	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
Low Positive	30/30	100.0%	30/30	100.0%	29/30	96.7%	89/90	98.9%
Moderate Positive	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
Negative	9/10	90.0%	10/10	100.0%	10/10	100.0%	29/30	96.7%
Negative Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
Positive Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%

**CROSSREACTIVITY**

Crossreactivity studies were previously performed with the Alethia Mycoplasma DNA Amplification assay. Positive and negative specimens were inoculated with potentially interfering bacterial or fungal organisms to minimum final concentrations of 1.0 x 10<sup>8</sup> CFU/mL. Viruses were tested at concentrations greater than 1.0 x 10<sup>5</sup> TCID<sub>50</sub>/mL or 1.0 x 10<sup>6</sup> copies/mL. Human DNA was tested at 2.0 ng/test with no crossreactivity observed. Positive specimens contained *M. pneumoniae* concentrations near the limit of detection. None of the following organisms or materials were identified as crossreactive or interfering in the Alethia Mycoplasma assay:

*Acinetobacter baumannii*, *Acinetobacter calcoaceticus*, *Actinomyces odontolyticus*, *Bacillus subtilis*, *Bacteroides fragilis*, *Bordetella parapertussis*, *Bordetella pertussis*, *Burkholderia cepacia*, *Candida albicans*, *Candida glabrata*, *Candida parapsilosis*, *Chlamydia pneumoniae*, *Citrobacter freundii*, *Clostridium difficile*, *Corynebacterium diphtheriae*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Escherichia coli* (ESBL), *Fusobacterium nucleatum*, *Haemophilus ducreyi*, *Haemophilus influenzae*, *Haemophilus parainfluenzae*, *Helicobacter pylori*, *Klebsiella pneumoniae*, *Klebsiella pneumoniae* (KPC), *Legionella pneumophila*, *Listeria monocytogenes*, *Moraxella catarrhalis*, *Mycoplasma genitalium*, *Mycoplasma hominis*, *Neisseria cinerea*, *Neisseria gonorrhoeae*, *Neisseria meningitidis*, *Nocardia asteroides*, *Peptostreptococcus anaerobius*, *Proteus mirabilis*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Salmonella paratyphi* (Group A), *Salmonella typhimurium* (Group B), *Serratia liquefaciens*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Streptococcus agalactiae* (Group B), *Streptococcus anginosus* (Group F), *Streptococcus bovis* (Group D), *Streptococcus canis* (Group G), *Streptococcus equisimilis*, *Streptococcus mitis*, *Streptococcus pneumoniae*, *Streptococcus pyogenes*, *Streptococcus salivarius*, *Ureaplasma urealyticum*, Adenovirus, Coronavirus, Coxsackievirus, Cytomegalovirus, Epstein Barr virus, Herpes simplex virus 1, Herpes simplex virus 2, Human metapneumovirus, Influenza A, Influenza B, Parainfluenza virus 1, Parainfluenza virus 2, Parainfluenza virus 3, Respiratory syncytial virus A, Respiratory syncytial virus B, Rhinovirus, Human DNA. No crossreactivity is expected with the Alethia Mycoplasma Direct assay.

**TESTS FOR INTERFERING SUBSTANCES**

Potentially interfering substances were tested with simulated negative and contrived low positive (*M. pneumoniae* strains M129 and FH) samples. The following substances, at the saturated solvent/diluent concentrations indicated, do not interfere with Alethia Mycoplasma Direct assay: Acetaminophen (18.1 mg/mL), Albuterol Sulfate (20 mg/mL), Aspirin (9.1 mg/mL), Azithromycin dehydrate (2.0 mg/mL), Cepacol® Mouthwash [Alcohol (1.4% v/v), Cetylpyridinium chloride (0.005% v/v)], Contac® Cold + Flu Tablets [Acetaminophen (14.8 mg/mL), Chlorpheniramine maleate (0.06 mg/mL), Phenylephrine HCl (0.15 mg/mL)], Diphenhydramine HCl (2.6 mg/mL), Erythromycin (20.0 mg/mL), HALLS® Cough Drops [Menthol (0.06 mg/mL)], Ibuprofen (12.7 mg/mL), Phenylephrine HCl (0.595 mg/mL), Prednisone (20.0 mg/mL), "Robitussin® Cough+Chest Congestion Cough Syrup [Dextromethorphan HBr (0.20 mg/mL), Guaifenesin (2.0 mg/mL)], Saine Nasal Spray [Sodium chloride (0.65 mg/mL)], Mucus (5.0mg/mL), and White Blood Cells (0.5% v/v).

Whole blood was found to produce invalid results at concentrations greater than 2%.

Phenylephrine HCl at concentrations greater than 0.595 mg/mL may produce false negative results with low positive samples.

\* *Robitussin® Cough+Congestion Cough Syrup initially produced 1/3 invalid results with Positive Sample P2 that produced acceptable results during repeat testing (1/1 positive). Robitussin® Cough+Congestion Cough Syrup is not considered an interferent.*

**ITALIANO**

# alethia™

## Mycoplasma Direct DNA Amplification Assay

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## REFERENCES

1. Nagamine K, Hase T, Notoni, T. Accelerated reaction by loop-mediated isothermal amplifications using loop primers. Mol Cell Probes 2002;16:223-29.
2. Mori, Y., Kitao, M., Tomita, N., Notoni, T. Real-time turbidimetry of LAMP reaction for quantitating template DNA. J Biochem Biophys 2004;59:145-47.
3. Waites KB, et al. New insights into the pathogenesis and detection of *Mycoplasma pneumoniae* infections. Future Microbiol. 2008 December; 3(6):635-648.
4. Waites KB, Talkington DF. *Mycoplasma pneumoniae* and Its role as a human pathogen. Clin Microbiol Rev. 2004 October; 17(4):697-728.

5. Thibodeau KP, Viera, AJ. Atypical pathogens and challenges in community-acquired pneumonia. *American Family Physician*. 2004 April; Vol 66, Number 7.
6. US Department of Health and Human Services, Centers for Disease Control and Prevention. *Mycoplasma pneumoniae* Infection; <http://www.cdc.gov/pneumonia/atypical/mycoplasma/index.html>
7. American Academy of Pediatrics. *Red Book: 2012 Mycoplasma pneumoniae and Other Mycoplasma Species Infections*. Pickering LK, ed. 29th ed. Elk Grove Village, IL: American Academy of Pediatrics; 2012.
8. US Department of Health and Human Services PHS/CDC/NIH. Biosafety in microbiology and biomedical laboratories, Washington DC: US Government Printing Office, 2007.
9. CLSI: MM3-A2 Molecular diagnostic methods for infectious disease; approved guideline, 2nd ed. Wayne PA: Clinical Laboratory Standards Institute. 2006.

SN11039

REV. 10/18

 Manufactured By	<p><b>Meridian Bioscience, Inc.</b>  <b>Corporate Office</b>          3471 River Hills Drive          Cincinnati, Ohio 45244 USA  <b>Telephone: 513.271.3700</b>  <b>Orders/Customer Service:</b>          800.543.1980  <b>Technical Support Center:</b>          800.343.3858  <b>Information Fax: 513.272.5432</b>  <b>Ordering Fax: 513.271.0124</b></p>
<div style="border: 1px solid black; padding: 5px; display: inline-block;"> <b>EC REP</b> </div> Authorized Representative	<p><b>Meridian Bioscience Europe S. r. L</b>          Via dell' Industria, 7          20020 Villa Cortese, Milano          ITALY  <b>Tel: +39 0331 43 36 36</b>  <b>Fax: +39 0331 43 36 16</b>  <b>Email: <a href="mailto:info@meridianbioscience.eu">info@meridianbioscience.eu</a></b>  <b>WEB: <a href="http://www.meridianbioscience.com/eu">www.meridianbioscience.com/eu</a></b></p>

**Meridian Bioscience Europe s.a./n.v.**  
 2 Avenue du Japon - 1420 Braine l'Alleud  
 BELGIUM  
 Tel: +32 (0) 67 89 59 59  
 Fax: +32 (0) 67 89 59 58  
 Email: [info.bnl@meridianbioscience.eu](mailto:info.bnl@meridianbioscience.eu)


**Meridian Bioscience Europe France**  
 34 rue de Ponthieu - 75008 Paris  
 FRANCE  
 Tel: +33 (0) 1 42 56 04 40  
 Fax: +33 (0) 9 70 06 62 10  
 Email: [info.fr@meridianbioscience.eu](mailto:info.fr@meridianbioscience.eu)


**Meridian Bioscience Europe b.v.**  
 Postbus 301 - 5460 AH Veghel  
 NETHERLANDS  
 Tel: +31 (0) 411 62 11 66  
 Fax: +31 (0) 411 62 48 41  
 Email: [info.bnl@meridianbioscience.eu](mailto:info.bnl@meridianbioscience.eu)

## INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

**Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)**

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	<b>CONTROL</b> +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
<b>LOT</b>	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	<b>CONTROL</b> -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
<b>IVD</b>	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	<b>EC REP</b>	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostika 98/79/EG.	<b>SMP PREP DIL SPE</b>	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet
		<b>STERILE R</b>	Sterilization by gamma irradiation / Sterilizzazione con raggi gamma / Sterilisation par irradiation aux rayons gamma / Esterilizado por irradiación gamma / Sterilisation durch Gammastrahlen
<b>REF</b>	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	<b>RoHS</b>	Restriction of Hazardous Substances / Restrizione all'uso di sostanze pericolose / Limitation de substances dangereuses / Restricción de Substancias Nocivas / Beschränkung der Verwendung bestimmter gefährlicher Stoffe
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		Caution, consult accompanying documents / Attenzione, vedere le istruzioni per l'uso / Attention voir notice d'instructions / Atención, ver instrucciones de uso / Achtung, Begleitdokumente beachten
<b>SN</b>	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	<b>STERILE EO</b>	Sterilization by ethylene oxide / Sterilizzazione con ossido di etilene / Stérilisation par oxyde d'éthylène / Esterilizado por óxido de etileno / Sterilisation durch Ethylenoxid
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	<b>BUF RXN</b>	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung		ETL Registered Mark Certified / Marchio di certificazione registrato a livello nazionale / Certifié Conforme ETL / Marca de Certificación Registrada Nacional / ETL Konform beglaubigt
	Female / Femmineile / De sexe féminin / Hembra / Frau		Recycle - do not dispose of as general waste / Riciclare - non eliminare come rifiuto generico / Recycler - ne pas jeter dans une poubelle / Recicle - no desecho como basura general / Recycling- dieses Produkt nicht über den Hausmüll entsorgen
	Single Use Only / Prodotto Monouso / A usage unique / Para Un Solo Uso / nur für die einmalige Anwendung		

	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät		Heat Treatment Tube / Provetta per il Trattamento termico / Tube pour le traitement thermique / Tubo de tratamiento de calor / Röhrchen zur Hitzebearbeitung
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	LASER RADIATION: Avoid Exposure to Beam / RADIAZIONE LASER: Evitare l'esposizione al raggio / RAYONNEMENT LASER: Eviter toute exposition au faisceau / Radiación Laser: Evite Exposición a los Rayos / LASERSTRAHLUNG: Direkten Kontakt mit dem Strahl vermeiden		HOT SURFACE: Keep hands Away from Hot Surfaces / Superficie calda: tenere le mani lontane dalle superfici calde / SURFACES CHAUDES: Ne pas toucher les surfaces chaudes / Superficie Caliente: Mantenga las manos alejadas de la superficie caliente / Heiße Oberfläche: Kontakt mit heißen Oberflächen vermeiden
	CAUTION: Laser Radiation / ATTENZIONE: Radiazione Laser / AVERTISSEMENT: Rayonnement laser / Precaución: Radiación Laser / WARNUNG: Lasersstrahlung		CAUTION: Protect from water / ATTENZIONE: Proteggere dall'acqua / AVERTISSEMENT: Protéger de l'humidité / Precaución: Proteja del agua / WARNUNG: Vor Feuchtigkeit schützen
	CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risques de danger / Precaución: Peligroso / WARNUNG: Riskogefahr		Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
	Buffer / Soluzione tampne / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
	Mineral Oil / Olio Minerale / Huile Minérale / Aceite Mineral / Mineralöl		Media / Terreno di trasporto / Milieux / Medio / Medium
	Screw Top Tube / Provetta con tappo a vite / Tube à bouchon vissé / Tubo con tapa de rosca / Röhrchen mit Schnappverschluss		Sample Preparation Column / Colonna di preparazione dei campione / Colonne pour la préparation de l'échantillon / Columna de preparación de muestra / Säule zur Probenaufarbeitung
	Sample Buffer / Soluzione tampone per il campione / Tampon de l'échantillon / Tampón de muestra / Probenpuffer		Pretreatment Reagent / Reagente di Prettrattamento / Réactif de prétraitement / Reactivo de pretratamiento / Reagenz für die Vorbehandlung
	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig		Sample Preparation / Preparazione del campione / Préparation de l'échantillon / Preparación de Muestra / Probenvorbereitung
	Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß		

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.