alethia[™]

Pertussis

DNA Amplification Assay

DNA Amplification Assay for the Detection of *Bordetella pertussis* in nasopharyngeal swab samples

REF 480750

IVD

R_x Only

INTENDED USE

The Alethia Pertussis DNA Amplification Assay, performed on the Alethia Reader, is a qualitative in vitro diagnostic test for the direct detection of *Bordetella pertussis* in human nasopharyngeal swab samples taken from patients suspected of having respiratory tract infection attributable to *Bordetella pertussis*.

The Alethia Pertussis assay utilizes loop-mediated isothermal DNA amplification (LAMP) technology to detect *B. pertussis* by targeting the IS481 insertional element of the *B. pertussis* genome. The IS481 insertional element can also be found in *B. holmesii* and some *B. bronchiseptica* strains. Respiratory infections with *B. pertussis*, *B. holmesii* or *B. bronchiseptica* may yield positive test results in IS481 assays. *B. holmesii* infection may cause clinical illness similar to *B. pertussis*, and mixed outbreaks involving both *B. pertussis* and *B. holmesii* infection have been reported. Additional testing should be performed if necessary to differentiate *B. holmesii* and *B. pertussis*. *B. bronchiseptica* is a rare cause of infection in humans. When clinical factors suggest that *B. pertussis* may not be the cause of respiratory infection, other clinically appropriate investigation(s) should be carried out in accordance with published guidelines.

Negative results for the Alethia Pertussis DNA Amplification Assay do not preclude *Bordetella pertussis* infection and positive results do not rule out co-infection with other respiratory pathogens. Results from the Alethia Pertussis assay should be used in conjunction with information obtained during the patient's clinical evaluation as an aid in diagnosis of *B. pertussis* infection and should not be used as the sole basis for treatment or other patient management decisions.

Alethia Pertussis is intended for use in hospital, reference or state laboratory settings. The device is not intended for point-of-care use.

SUMMARY AND EXPLANATION OF THE TEST

The Alethia Pertussis DNA Amplification Assay is based on loop-mediated amplification (LAMP)¹⁻² technology. The assay targets a 198 base pair sequence of the *Bordetella pertussis* genome residing in a region of the IS481 insertional element sequence.

Loop-mediated amplification uses specially designed primers to provide for specific and continuous isothermal DNA amplification. A by-product of this amplification is the formation of magnesium pyrophosphate, which forms a white precipitate leading to a turbid reaction solution. Reaction solution absorbance characteristics are monitored by the Meridian Alethia Incubator/Reader. Changes in reaction solution absorbance characteristics created by precipitation of magnesium pyrophosphate indicate the presence of target DNA. The absence of target DNA results in no significant change in sample absorbance.

The Alethia Pertussis kit includes Alethia Pertussis Assay Control/Negative Control Reagent, Alethia Pertussis Test Devices, Alethia Pertussis Sample Buffer, and Mineral Oil. The Alethia Assay Control/Negative Control, used for specimen dilution and preparation, is a Tris-buffered solution containing formalin-treated *E. coli* harboring *Staphylococcus aureus* DNA. The Alethia Pertussis Test Device contains one lyophilized amplification bead in each of two chambers: a TEST chamber with IS481-specific primers and a CONTROL chamber with *S. aureus*-specific primers. The *S. aureus* DNA in the Assay Control/Negative Control Reagent and the *S. aureus*-specific primers in the CONTROL chamber function as the Internal Control for the assay. During specimen preparation, a patient specimen is added to the Assay Control/Negative Control Reagent. Adding the specimen to the Assay Control/Negative Control reagent combines specimen DNA with the Control *S. aureus* DNA and allows for parallel processing of target DNA and Control DNA through amplification and detection. The Internal Control monitors amplification inhibition, assay reagent performance and sample processing effectiveness. The Control *S. aureus* target must be amplified and detected in the final reaction or the test is considered invalid and patient results are not reported.

The Alethia Reader monitors changes in absorbance characteristics by measuring transmission of light through the Test and Control reaction solutions. Light transmission is checked at the assay Run Start (Signal_{Initial}, S_i) and at the assay Run End (Signal_{Final}, S_f). The Alethia Reader calculates the change in light transmission between Run End and Run Start (S_f:S_i) and compares the ratio to a fixed cut-off value.

Fixed cut-off values for the TEST chamber are used to report sample results. TEST chamber S_f:S_i ratios less than 82% are reported as 'POSITIVE'; TEST chamber S_f:S_i ratios greater than or equal to 82% are reported as 'NEGATIVE'. Numerical values are not reported.

Fixed cut-off values for the CONTROL chamber are used to determine validity. CONTROL chamber S_f:S_i ratios less than 90% are considered valid and allow for reporting of TEST chamber results (POSITIVE, NEGATIVE). CONTROL chamber S_f:S_i ratios greater than or equal to 90% are considered invalid and prevent reporting of TEST chamber results. Invalid CONTROL chamber reactions are reported as 'INVALID'. Numerical values are not reported.

More stringent cut-off criteria are applied to the CONTROL chamber reaction to ensure amplification is not inhibited, reagents are performing as intended and that sample processing was performed appropriately.

BIOLOGICAL PRINCIPLES

B. pertussis is a human pathogen exclusively responsible for the endemic respiratory disease whooping cough (pertussis); a publically reportable disease that impacted over 200,000 people worldwide in 2012. Whooping cough is currently the only disease in the US that displays an increasing trend of reportable cases and substantial morbidity despite the availability of a vaccine.³ The US Centers for Disease Control (CDC) reported approximately 48,300 cases of pertussis in 2012, an average of 15.4 cases per 100,000 people in the United States.⁴ Recent pertussis outbreaks have indicated waning immunity in children after receiving the final booster acellular vaccine dose, allowing for susceptibility of pertussis infection in both vaccinated and unvaccinated populations of children ages 7-11.^{5,6}

Pertussis outbreaks are particularly challenging to control as the early stages of disease resemble other respiratory infections and patients can remain highly infectious for up to 5 weeks after onset of symptoms.³ In addition, pertussis is not officially diagnosed until the presence of the characteristic "whooping" cough that is displayed almost 2 weeks post-symptom onset.³ Clinical diagnosis of pertussis can be performed using culture, serology, or nucleic acid amplification tests (NAAT). While culture is highly specific, sensitivity is low, results can take up to 7 days, and the viable cells required for culture decrease with disease progression (after 2 weeks of symptom onset).⁷ Unlike culture, NAAT assays can provide rapid test results and do not require viable bacteria. Nucleic acid amplification testing for *B. pertussis* is most sensitive during the first three weeks of cough as bacterial DNA is present in the nasopharynx. The amount of bacteria present in the nasopharynx dramatically decreases after the fourth week of cough, making the likelihood of obtaining false-negative test results higher during the later stages of infection.⁸ In addition, only symptomatic patients with cough should be tested by NAAT in order to prevent false positive results by asymptomatic close contacts.⁹ Serology can only be used for diagnosis in late stages of disease, approximately 2-8 weeks after cough onset.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. **Alethia Pertussis Assay Control/Negative Control Reagent:** Tris-buffered solution containing formalin-treated *E. coli* harboring plasmid containing a segment of the *S. aureus* genome and sodium azide (0.09%) as a preservative. The Assay Control/Negative Control Reagent functions as the Assay Control during patient testing and as the External Negative Control during routine Quality Control testing.

2. **Alethia Pertussis Test Device:** Two-chambered device containing lyophilized amplification reagents (DNA polymerase, deoxynucleotide triphosphates), and either IS481-specific primers (TEST Chamber) or control primers (CONTROL Chamber).
3. **Alethia Pertussis Sample Buffer:** Tris-EDTA solution containing sodium azide (0.09%) as a preservative.
4. **Mineral Oil** (bottle with dropper tip)
5. **Alethia Pretreatment Reagent (for optional sample treatment procedure):** Lyophilized reagent containing neutralizing agent and sodium azide (0.1%) as a preservative. Each vial contains enough reagent to treat 50 specimens.

MATERIALS PROVIDED SEPARATELY

Alethia Pertussis External Control Kit, Catalog Number: 479930

MATERIALS NOT PROVIDED

1. Disposable latex gloves, powder free
2. DNase/RNase-free, aerosol resistant pipette tips
3. Specimen collection and transport system
- Nasopharyngeal Swabs:** Polyester (minimum capacity 18 µL, eg, Puritan Medical Products catalog 25-801D 50); Rayon (minimum capacity 31 µL, eg, Puritan Medical Products catalog 25-801R 50); or Flocked Nylon (Minimum Capacity 69 µL, eg, Copan catalog 503CS01 or the swab of Catalog 482C).
- Nonnutritive Transport Medium:** Liquid Amies without charcoal; or Liquid Stuart (Maximum Volume: 1.2 mL in pledget/sponge)
- Collection System:** Copan Liquid Amies Elution Swab (ESwab[™]) Collection and Transport System (Copan catalog 482C)
4. Deionized water
5. Microcentrifuge tubes

EQUIPMENT NOT PROVIDED

1. Dry-bath with 12 mm heat block capable of 95 C
2. Digital thermometer with max/min temperature memory (eg, Traceable® Lollipop[™] Waterproof/Shockproof Thermometer)
3. Vortex mixer
4. Interval timer
5. Micropipette capable of dispensing 50 µL
6. Pipette capable of delivering 2.5 mL
7. Cutting equipment (eg, scissors, shears, safety snips)
8. Alethia Reader, Meridian Bioscience, Inc. Catalog Number: 610189

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Do not interchange Assay Control/Negative Control Reagent or Test Devices between lots. Sample Buffer, Pretreatment Reagent and Mineral Oil are interchangeable provided they are within assigned expiration dates when used.
3. Follow Biosafety Level 2 and Good Laboratory practices during testing.⁹ Treat all specimens and used Test Devices as capable of transmitting infectious agents. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
4. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
5. Quality Control Programs for Molecular Testing Laboratories, including proper use and care of equipment, should be employed.¹⁰
6. The Alethia Pertussis Test Device contains lyophilized reagents. The protective pouch should not be opened until ready to perform the assay.
7. The Alethia Pertussis Test Device includes a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
8. Dispose of used Alethia Test Devices immediately after processing, leaving the device latch securely in place. Do NOT open the Test Device after processing. Opening the device after amplification may result in contamination of the test area with amplification product.
9. Invalid test results may result due to improper sample preparation, reagent failure, instrument failure or internal control failure. In these cases, the invalid result can be corrected by repeating the test. The nasopharynx of some patients contains a biological inhibitor that will also cause an invalid result. In these cases, the inhibitor will cause repeated invalid results unless the sample is treated to remove the interference. (See section on OPTIONAL SAMPLE TREATMENT/PRETREATMENT for the treatment method.)
10. Local, state, and federal regulations for notification of reportable disease are continually updated and include a number of organisms for surveillance and outbreak investigations.^{a,b} Additionally, the Centers for Disease Control (CDC) recommends that when pathogens from reportable diseases are detected by a culture independent diagnostic test, the laboratory should facilitate obtaining the isolate or clinical materials for submission to the appropriate public health laboratory to aid in outbreak detection and epidemiological investigations. Laboratories are responsible for following their state and/or local regulations and should consult their local and/or state public health laboratories for isolate and/or clinical sample submission guidelines.
 - a. Summary of Notifiable Diseases. MMWR <http://www.cdc.gov/mmwr/PDF/wk/mm6153.pdf>
 - b. CIFOR Analysis of State Legal Authorities. <http://www.cifor.us/documents/CIFORAnalysisStateLegalAuthorities.pdf>

HAZARD and PRECAUTIONARY STATEMENTS

Refer to the SDS, available at www.meridianbioscience.com for Hazards and Precautionary Statements.

SHELF LIFE and STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-30 C.

REAGENT PREPARATION

Ensure kit reagents are at room temperature (21-30 C) before use. Incorrect results may be obtained if reagents are not brought to room temperature prior to use.

SPECIMEN COLLECTION and PREPARATION

Sample type: Nasopharyngeal swabs

Sample Collection: Nasopharyngeal swab specimen collection should be performed in accordance with institutional procedures for collection of clinical specimens for *B. pertussis* infection. Nasopharyngeal swab samples should be collected with suitable swab types (eg, Polyester, Flocked Nylon or Rayon).

Place swab(s) in nonnutritive transport medium (eg, Liquid Amies, without charcoal, Liquid Stuart or Copan ESwab medium) or store unpreserved in a sterile tube without medium.

Unpreserved samples or samples stored in transport media should be tested as soon as possible, but may be held at room temperature (21–30 C) for up to 5 days or refrigerated (2-8 C) for up to 7 days prior to testing. Do not freeze samples.

SPECIMEN PREPARATION:

NOTE: Ensure that the Alethia Reader is powered on and required performance verifications have been completed prior to initiation of SPECIMEN PREPARATION. Refer to the Alethia Operator's Manual for further information regarding instrument set-up and operation.

NOTE: Laboratory equipment used for cutting swabs (eg, scissors, shears, safety snips), should be treated with a molecular grade cleaning agent (eg, 10% bleach), prior to each use. Equipment should be completely dry before cutting swabs.

1. Specimen Preparation
 - a. Swabs (rayon, nylon or flocked):
 - i. Place the swab specimen in a labeled Sample Buffer tube. Cut the swab shaft to ensure sample fits in the tube. Elute the sample by vortexing for 45-60 seconds. **Eluted samples may be held at room temperature (21-30 C) for up to 48 hours or refrigerated (2-8 C) for up to 7 days prior to testing.**
 - ii. Add 50 µL of eluted sample to a labeled Assay Control/Negative Control tube and cap.
 - iii. Repeat SPECIMEN PREPARATION Steps for all swab samples to be processed.

- b. ESwab Medium: The medium within the ESwab System collection tube contains eluted patient sample.
 - i. Vortex the ESwab tube containing the swab sample for a minimum of 10 seconds to elute the sample.
 - ii. Transfer 25 μ L of the medium to a labeled Assay Control/Negative Control tube and cap. (Do not transfer this material to the Sample Buffer tube. It goes directly to the Assay Control/Negative Control Tube)
 - iii. Repeat SPECIMEN PREPARATION steps for all ESwab medium samples to be processed.
2. Vortex each Assay Control/Negative Control tube containing eluted sample for approximately 10 seconds.
3. Heat each Sample/Control tube in a dry-bath at 95 \pm 5 C for 10 \pm 2 minutes. Monitor the heat-treatment step with a digital thermometer and interval timer.
4. Remove each Sample/Control tube from the dry-bath. Heat-treated samples may be held at room temperature (21-30 C) for up to 15 minutes prior to testing.
5. Vortex for approximately 10 seconds.

TEST PROCEDURE

NOTE: A maximum of 10 samples can be processed in a single Alethia Reader run.

1. Remove 1 Alethia Pertussis Test Device from its protective pouch per sample. Carefully open the device, holding the chambers such that the lyophilized reagent will not fall out upon opening. Place device on a flat surface or in a rack that can accommodate the device.
2. Using a micropipette, first transfer 50 μ L of the heat-treated sample to the TEST chamber (left/White Bead) and then transfer 50 μ L of the heat treated sample to the CONTROL chamber (right/Yellow Bead) of the Alethia Test Device. Take care not to introduce air to the reaction mixture. Do not mix reactions with pipette.
3. Add 1 drop of Mineral Oil to both the TEST chamber and CONTROL chamber. Close the Alethia Test Device and fasten the latch securely.
4. Tap device on the bench top or mix to remove air bubbles. Carefully examine the Test Device for rehydration of the Control/Test Bead, for air bubbles left in the chamber and liquid in the top of the device. If undissolved beads, air bubbles or liquid in the top of the device are noted, tap the device on the bench top and repeat visual inspection. Amplification and detection should be initiated within 15 minutes.
5. Insert the Alethia Test Device into the Alethia Reader and initiate amplification reaction and detection. Results will be displayed at the conclusion of the run.

INTERPRETATION OF RESULTS

Sample ID	Reported Result	Interpretation
Patient Specimen	POSITIVE	Untreated or Pretreated Specimen: Sample contains <i>B. pertussis</i> IS481 target DNA. ^a
	NEGATIVE	Untreated or Pretreated Specimen: No <i>B. pertussis</i> DNA detected.
	INVALID	Untreated Specimen: No reportable result. Repeat the test using a Pretreatment Reagent-treated/eluted sample. Inhibitory patient specimen, improper sample preparation, reagent failure, instrument failure or internal control failure. Pretreated Specimen: No reportable result. DO NOT retreat specimens already treated with Pretreatment Reagent. Pretreated INVALID samples can be retested, if desired, and should be reported as INVALID if retesting generates the same result. Samples that produce INVALID results before and following pretreatment should be reported as INVALID.
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, Alethia Reader performing correctly.
	NEGATIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	INVALID	No reportable result. Repeat entire assay run using original samples. If the repeated control result is INVALID, please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor. Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, Alethia Reader performing correctly.
	INVALID	No reportable result. Repeat entire assay run using original samples. If the repeated control result is INVALID, please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor. Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No Alethia Test Device in the Alethia Reader Well. OR The Alethia Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. Repeat the test using original sample.

^aIS481 is found in multiple copies in *B. pertussis* (50 to 238 copies/genome), in *B. holmesii* (8 to 10 copies/genome) and less frequently in *B. bronchiseptica*.

OPTIONAL SAMPLE PRETREATMENT/TREATMENT TO ELIMINATE REACTION INHIBITORS

The following method can be used to remove the inhibitory activity of biological substances in a patient's sample.

Preparation of Pretreatment Reagent

1. Reconstitute lyophilized Pretreatment Reagent by adding 2.5 mL deionized water to the vial. Recap the vial.
2. Allow the reconstituted Pretreatment Reagent to stand at 21-30 C for 30 minutes, then mix gently to avoid foaming.
3. Store the reconstituted Pretreatment Reagent in one of two ways:
 - a. Store reconstituted Pretreatment Reagent in its original vial at 2-8 C for up to one month. Label the vial with its new expiration date, or
 - b. Immediately transfer 0.1 mL or larger aliquots of the reconstituted reagent to labeled microcentrifuge tubes, cap securely and freeze the aliquots in a nonfreezing freezer at \leq -20 C. Frozen aliquots are stable to the original expiration of the lyophilized material. Thaw each frozen aliquot no more than once before use. Thawed aliquots can be stored at 2-8 C for up to one week.
4. Warm reconstituted Pretreatment Reagent to 21-30 C before use. Mix gently to avoid foaming. Do not use the Pretreatment Reagent solution if it appears cloudy or turbid.

Treatment of Previously Prepared Samples for Retesting:

- 1a. Swabs (rayon, nylon or flocced):
 - i. Add 50 μ L of reconstituted Pretreatment Reagent to the Sample Buffer tube containing the patient's sample (remaining from SPECIMEN PREPARATION step 1.a.i.) and vortex for 10-15 seconds to mix the contents. The treated sample should not be stored beyond the original eluted specimen hold time described in step 1 of the SPECIMEN PREPARATION section above.
 - ii. Transfer 50 μ L of the medium to a labeled Assay Control/Negative Control tube and cap.
- 1b. ESwab medium:
 - i. Transfer 500 μ L of ESwab medium inoculated with the patient's sample (from the ESwab collection tube) to an empty tube.
 - ii. Add 50 μ L of reconstituted Pretreatment Reagent to the tube and vortex for 10-15 seconds to mix the contents. The treated sample should not be stored beyond the original eluted specimen hold time described in step 1 of the SPECIMEN PREPARATION section above.
 - iii. Transfer 25 μ L of the medium to a labeled Assay Control/Negative Control tube and cap.
2. Complete Steps 2-5 of the SPECIMEN PREPARATION section and proceed to testing.

Pretreatment of Untested Samples:

- 1a. Swabs (rayon, nylon or flocced):
 - i. Add 50 μ L of reconstituted Pretreatment Reagent to a labeled Sample Buffer tube.
 - ii. Place the swab specimen in the Sample Buffer tube. Cut the swab shaft to ensure the sample fits in the tube. Elute the sample by vortexing for 45-60 seconds. The treated, eluted sample may be held at room temperature (21-30 C) for up to 48 hours or refrigerated (2-8 C) for up to 7 days prior to testing.
 - iii. Transfer 50 μ L of the medium to a labeled Assay Control/Negative Control tube and cap.
- 1b. ESwab medium:
 - i. Vortex the ESwab tube containing the swab sample for a minimum of 10 seconds to elute the sample.
 - ii. Transfer 500 μ L of ESwab medium inoculated with the patient's sample to an empty tube.
 - iii. Add 50 μ L of reconstituted Pretreatment Reagent to the tube and vortex for 10-15 seconds to mix the contents. The treated sample should not be stored beyond the original eluted specimen hold time described in step 1 of the SPECIMEN PREPARATION section above.
 - iv. Transfer 25 μ L of the medium to a labeled Assay Control/Negative Control tube and cap.
2. Complete Steps 2-5 of the SPECIMEN PREPARATION section and proceed to testing.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

1. Each device contains an internal control chamber that controls for amplification inhibition, assay reagents and sample processing effectiveness.
2. The heat-treatment step is monitored with an external thermometer and interval timer. Use the max/min temperature memory of the thermometer to ensure that a temperature of 95 \pm 5 C is maintained. Use the interval timer to ensure that heat-treatment duration is 10 \pm 2 minutes.
3. Good laboratory practice recommends the use of control materials. Users should follow the appropriate federal, state and local guidelines concerning the running of external quality controls.
4. Alethia Pertussis External Control Reagents are supplied separately (Catalog 479930). It is recommended that reactivity of each new lot and each new shipment of Alethia Pertussis be verified on receipt or before use. External control tests should be performed thereafter in accordance with appropriate federal, state and local guidelines. The Alethia Pertussis test kit should not be used in patient testing if the external controls do not produce the correct results.
5. A separate device must be used for each external control reagent.

EXPECTED VALUES

Overall incidence of *B. pertussis* as detected by the Alethia Pertussis Assay in prospectively and retrospectively collected, non-selected specimens (all comers) during the period of this study was 8.2% (57/692).

LIMITATIONS OF THE PROCEDURE

1. The Alethia Pertussis assay targets the IS481 insertional element of the Bordetella genome. The IS481 insertional element is present in *B. pertussis*, *B. holmesii* and some strains of *B. bronchiseptica*.
2. This product can be used only with the Alethia instrument.
3. The Alethia Pertussis DNA assay is a qualitative assay and does not provide quantitative values or information about organism load.
4. This device has not been evaluated for monitoring treatment of *B. pertussis* infections.
5. This test has not been evaluated for specimens other than nasopharyngeal swab specimens, for immunocompromised individuals or from patients not suspected of infection with *B. pertussis*.
6. Results from this test must be correlated with the clinical history, epidemiological data, and any other data available to the clinician.
7. Prevalence of *B. pertussis* will affect the positive and negative predictive values for the assay.
8. *B. parapertussis* which causes a pertussis-like illness is not detected by the Alethia Pertussis DNA assay. Illness caused by *B. parapertussis* is generally milder than illness caused by *B. pertussis* because the bacteria do not produce pertussis toxin.
9. Respiratory infections can be caused by *B. pertussis* as well as other pathogens. Positive results do not preclude coinfection with other respiratory pathogens. False-negative *B. pertussis* results are more likely if patients are tested later in the disease course (more than two weeks after symptom onset), due to declining Bordetella DNA. False-negative results may also be increased in patients treated with antibiotic therapy.
10. Environmental contamination of an exam room from a prior patient or a recent pertussis vaccination administration may result in false-positive test results.
11. The detection of nucleic acid is dependent upon proper specimen collection, handling, transportation, storage and preparation. Failure to observe proper procedure in any one of these steps can lead to incorrect results.
12. Organism nucleic acids may persist *in vivo*, independent of organism viability. The Alethia Pertussis assay does not distinguish between viable and nonviable organisms.
13. As with all molecular based diagnostic tests, (A) False negative results may occur from the presence of inhibitors, technical error, sample mix-up or low numbers of organisms in the clinical specimen; (B) False positive results may occur from the presence of cross-contamination by target organisms, their nucleic acids or amplified product, and from non-specific signals.
14. Acetyl salicylic acid, as found in aspirin, produced invalid results when tested at concentrations above 5 mg/mL during *B. pertussis* strain BAA-589 Limit of Detection replicate testing.

SPECIFIC PERFORMANCE CHARACTERISTICS

The Alethia Pertussis DNA Amplification Assay was evaluated from December 2012 to July 2013 by independent clinical test sites representing geographically distinct regions throughout the United States. A total of 729 qualified nasopharyngeal (NP) swab specimens collected from patients suspected of respiratory infections due to *B. pertussis* were evaluated with the test device to establish performance characteristics. Specimens were leftover, de-identified specimens that had previously been submitted for routine *B. pertussis* testing. Specimens included in performance evaluation were prospective (never frozen) and retrospective (frozen prior to Alethia testing). The retrospective population included non-selected (all-comers) and selected specimens.

The performance of Alethia Pertussis was compared to a Composite Reference Method that included two manufacturer validated, IS481-targeted real-time PCR assays (PCR1 and PCR2) followed by bi-directional sequencing of amplicon from PCR positive specimens.

The PCR1 and PCR2 assay protocols included 40 amplification cycles. Bi-directional sequencing was performed for all specimens producing amplicon prior to the end of the 40-cycle amplification. Specimens were considered positive when bi-directional sequencing results from either comparator PCR assay confirmed the presence of *B. pertussis* amplicon. Specimens were considered negative when neither comparator PCR assay produced amplicon at the end of the 40-cycle amplification.

A total of 508 (69.7%) prospective and 221 (30.3%) retrospective specimens were tested. Invalid results were obtained for 13 specimens (1.8%); two specimens remained invalid after repeat testing. Two prospective and two retrospective specimens produced indeterminate comparator PCR results. Table 1 summarizes Alethia performance characteristics.

Table 1. Alethia Assay Performance

Specimen Description	Positive Specimens			Negative Specimens			Invalid Results ^b
	Alethia vs. Comparator	PPA ^a	95% CI	Alethia vs. Comparator	NPA	95% CI	
Composite Method Comparator, All Comers							
Prospective	39/45	86.7%	73.8 – 93.7%	447/459	97.4%	95.5 – 98.5%	2 (13)
Retrospective	4/4	100.0%	51.0 – 100.0%	176/178	98.9%	96.0 – 99.7%	0
Total:	43/49	87.8%	75.8 – 94.3%	623/637	97.8%	96.3 – 98.7%	2 (13)
Composite Method Comparator, Selected Specimens							
Retrospective	19/21	90.5%	71.1 – 97.3%	14/16	87.5%	64.0 – 96.5%	0

^a Eight specimens produced false-negative Alethia results when compared to the Composite Comparator Method. Six of the eight specimens produced detectable levels of DNA between 35 and 40 comparator assay amplification cycles and were confirmed positive by bi-directional sequencing. Three of these six specimens gave positive results in only one of the two comparator PCR/Sequencing assays.

^b 11/13 initial invalid specimens produced valid results upon repeat testing.

False-negative Alethia results were individually evaluated at the conclusion of clinical testing. The Cycle threshold (Ct) values produced during comparator assay testing were above 35 for one or both PCR/Bi-directional sequencing assays for six of the eight specimens evaluated. High Ct values in PCR assays may indicate that low levels of DNA are present.⁵ False-negative Alethia results and corresponding Composite Comparator data are shown in Table 2.

Table 2. False-Negative Alethia Specimens, Comparator Assay Results

Specimen Designation	Specimen Status	PCR1		PCR2	
		Ct Value	Bi-Directional Sequencing Result	Ct Value	Bi-Directional Sequencing Result
1-19	Prospective	34.90	+	33.85	+
1-29	Prospective	Negative	N/A	37.69	+
1-259	Prospective	34.07	+	35.09	+
1-269	Prospective	35.72	+	Negative	N/A
1-275	Prospective	Negative	N/A	39.13	+
3-33	Prospective	39.28	+	35.80	+
4-710	Retrospective	36.87	+	36.06	+
4-712	Retrospective	37.44	+	38.41	+

Performance characteristics for the Alethia assay were evaluated according to Clinical Test Site. Table 3 shows performance characteristics by Clinical Test Site.

Table 3. Alethia Assay Performance by Clinical Test Site

Specimen Description	Positive Specimens			Negative Specimens			Invalid Results ^b
	Alethia vs. Comparator	PPA	95% CI	Alethia vs. Comparator	NPA	95% CI	
Composite Method Comparator, All Comers							
Site 1	35/40	87.5%	73.9 – 94.5%	440/450	97.8%	96.0 – 98.8%	2 (13)
Site 2	4/4	100.0%	51.0 – 100.0%	67/69	97.1%	90.0 – 99.2%	0 (2)
Site 3	0/1	0.0%	0.0 – 79.3%	7/7	100.0%	64.6 – 100.0%	0 (0)
Site 4	4/4	100.0%	51.0 – 100.0%	109/111	98.2%	93.7 – 99.5%	0 (0)
Composite Method Comparator, Selected Specimens							
Site 1	15/15	100.0%	79.6% - 100.0%	6/8	75.0%	40.9 – 92.9%	0 (0)
Site 4	4/6	66.7%	30.0 – 90.3%	8/8	100.0%	67.6 – 100.0%	0 (0)

Clinical studies were conducted with multiple nasopharyngeal swab and sample elution buffer types. Sample buffers tested during clinical studies included 0.85% Saline (n=30 or 4.1%), Tris EDTA (n=8 or 1.1%) and Molecular Grade Water (n= 687 or 94.2%). All sample buffers were used in 0.5 mL volumes. Analytical studies were performed with 0.85% Saline, Tris EDTA, Phosphate Buffered Saline (PBS) and Molecular Grade Water. Analytical studies established equivalence between all sample elution buffer types.

Age information was known for 723 (99.2%) of the patients from whom samples were tested. Patient age ranged from 1 month to 88 years. Thirty-eight (5.2%) patients were less than 1 year of age; 13 (1.8%) were between 1 and 2 years old; 296 (40.6%) were between 2 and up to 12 years, 157 (21.5%) were between 12 and up to 21 years, 190 (26.0%) were above 21 but below 65 years, and the remaining 29 (4.0%) patients were above 65 years of age.

The study population included 413 (56.7%) female and 308 (42.2%) male patients. Gender was unknown for 8 (1.1%) patients included in the study. There is no expectation that the Alethia Pertussis assay performance characteristics are influenced by patient gender.

Studies using Pretreatment Reagent

Use of the Pretreatment Reagent was evaluated in external studies conducted from April to June 2015 at four selected test sites with prospectively collected specimens from symptomatic patients. Of the 164 samples collected for this study, 145 met study criteria. The majority of samples submitted to the test sites were collected with rayon tipped swabs; one sample was collected with a flocced nylon swab and three samples with a polyester swab. The majority of the samples were tested on the same day of collection. Collected specimens were tested immediately after collection, first using the predicate testing swab method, with retesting of all invalid results. Following testing with the predicate method the same specimens were then treated with the Pretreatment Reagent and tested after pretreatment with any invalid result retested. All predicate, invalid and pretreatment testing was completed within a 24 hour period. Of the 19 samples producing invalid test results by the predicate method, 15 were repeatedly invalid on retesting; three were negative and one positive on retesting. Following addition of Pretreatment Reagent the 19 samples that initially gave invalid results produced 17 negative and two positive results. Six samples gave a result with the predicate method, but produced an invalid result upon treatment with the Pre-treatment Reagent. The data is presented in Table 4 below.

Table 4. Method Comparison: Pretreatment Reagent versus Predicate Method (without Pretreatment Reagent)

	Test Sites Percent Agreement (Samples in Agreement/Total Tested)				
	Site 1	Site 2*	Site 3	Site 4	Total
PPA	87.5% (6/7) (CI: 48.7-97.4%)	0% (0/0)	100% (2/2) (CI: 34.2-100%)	100% (6/6) (CI: 61.8-100%)	93.3% (14/15) (CI: 70.2-98.8%)
NPA	94.3% (60/53) (CI: 84.6-98.1%)	93.3% (14/15) (CI: 70.2-98.8%)	100% (15/15) (CI: 79.6-100%)	100% (26/26) (CI: 87.1-100%)	96.3% (105/109) (CI: 90.9-98.6%)
Initial Invalid Rate without Pretreatment Reagent	1.6% (1/63)	0% (0/16)	37.0% (10/27)	20.5% (8/39)	13.1% (19/145)
Initial invalid rate with Pretreatment Reagent	4.8% (3/63)	6.25% (1/16)	3.7% (1/27)	2.6% (1/39)	4.1% (6/145)

* Site 2 data represents first pass data only as this site did not retest INVALID untreated samples, as required by the package insert, and before treating INVALID samples with Pretreatment Reagent. All other sites repeated INVALID test samples before pretreating samples.

ANALYTICAL SENSITIVITY

The analytical sensitivity or Limit of Detect for the Alethia Pertussis assay was determined for *B. pertussis* strain ATCC BAA-589, Tahoma 1.

Limit of Detection was determined using 60 replicates *B. pertussis* strain ATCC BAA-589 and a stated probability (eg. 95%, where 57/60 replicates are positive) of obtaining positive responses. Analytical sensitivity testing is summarized below.

<i>B. pertussis</i> Strain Description	CFU/mL	CFU/Test
ATCC BAA-589 (Tahoma 1)	3265	1.48

ASSAY REACTIVITY

The following *B. pertussis* strains were tested and produced positive reactions at 3265 CFU/mL or 1.48 CFU/Test with Alethia Pertussis: ATCC 12743; ATCC 8478; ATCC 8467; ATCC 9797; ATCC 53894; ATCC 10380; ATCC 12742; and A639. The following *B. pertussis* strains were tested and produced positive reactions at 3500 CFU/mL or 1.59 CFU/test with Alethia Pertussis: ATCC 51445 and ATCC BAA-1335.

REPRODUCIBILITY

Reproducibility studies were carried out by three of the four participating Clinical Sites. Blind-coded panels of 10 samples were supplied to participating laboratories. Samples were randomly sorted within each panel to mask sample identities. The panels included contrived samples manufactured as moderate positive samples (1.31 x 10⁴ CFU/mL or 6 CFU/test), low positive samples (4.89 x 10³ CFU/mL or 2 CFU/test), and high negative samples (8.2 CFU/mL or 0.004 CFU/test). The panel also included one negative sample, positive control and negative control. Testing was performed by different operators at each site on the same day (intra-assay variability) for five days (inter-assay variability). Three lots of Alethia Pertussis and six Alethia instruments were used in this study. Positive and Negative Controls were tested each day of testing. The results are provided in the table below:

Sample Type	Site 1		Site 2		Site 4		Total	
	Percent Agreement	Percent Agreement	Percent Agreement	Percent Agreement	Percent Agreement	Percent Agreement	Percent Agreement	
Moderate Positive	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
Low Positive	27/30	90.0%	29/30	96.7%	30/30	100.0%	86/90	95.6%
High Negative	26/30	86.7%	23/30	76.7%	29/30	96.7%	78/90	86.7%
Negative	10/10	100.0%	9/10	90.0%	10/10	100.0%	29/30	96.7%
Negative Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
Positive Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%

CROSSREACTIVITY

Crossreactivity studies were performed with positive and negative nasal wash specimens inoculated with bacterial or fungal organisms to a final concentration of 1.0 x 10⁶ CFU/mL or virus at a minimum of 1.0 x 10⁶ TCID₅₀/mL. None of the following organisms reacted with Alethia Pertussis: *Acinetobacter baumannii*, *Acinetobacter calcoaceticus*, *Acinetobacter lwoffii*, *Actinomyces odontolyticus*, *Arcanobacterium haemolyticum*, *Bacillus subtilis*, *Bacteroides fragilis*, *Bordetella avium*, *Bordetella hinzii*, *Bordetella parapertussis*, *Bordetella pertussis*, *Bordetella trematum*, *Burkholderia cepacia*, *Candida albicans*, *Candida glabrata*, *Chlamydia pneumoniae*, *Chlamydia trachomatis*, *Citrobacter freundii*, *Clostridium difficile*, *Corynebacterium diphtheriae*, *Corynebacterium pseudodiphtheriticum*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Escherichia coli*, *Escherichia coli* (ESBL), *Fusobacterium nucleatum*, *Haemophilus influenzae*, *Haemophilus parainfluenzae*, *Klebsiella oxytoca*, *Klebsiella pneumoniae* (KPC), *Lactobacillus acidophilus*, *Lactobacillus plantarum*, *Legionella jordanis*, *Legionella longbeachae*, *Legionella micdadei*, *Legionella pneumophila*, *Listeria monocytogenes*, *Moraxella catarrhalis*, *Mycobacterium tuberculosis*, *Mycoplasma genitalium*, *Mycoplasma hominis*, *Mycoplasma pneumoniae*, *Neisseria cinerea*, *Neisseria elongata*, *Neisseria gonorrhoeae*, *Neisseria meningitidis*, *Nocardia asteroides*, *Peptostreptococcus anaerobius*, *Proteus mirabilis*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Pseudomonas fluorescens*, *Serratia liquefaciens*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Stenotrophomonas maltophilia*, *Streptococcus anginosus* (Group F), *Streptococcus bovis* (Group D), *Streptococcus canis* (Group G), *Streptococcus dysgalactiae* ssp. *dysgalactiae*, *Streptococcus dysgalactiae* ssp. *equisimilis* *Streptococcus intermedius*, *Streptococcus mitis*, *Streptococcus mutans*, *Streptococcus pneumoniae*, *Streptococcus pyogenes*, *Streptococcus salivarius*, *Streptococcus suis*, *Ureaplasma urealyticum*, Adenovirus, Coronavirus, Coxsackievirus, Cytomegalovirus, Epstein Barr Virus, Herpes Simplex Virus 1, Herpes Simplex Virus 2, Human Metapneumovirus, Influenza A, Influenza B, Measles virus, Mumps virus, Parainfluenza virus 1, Parainfluenza virus 2, Parainfluenza virus 3, Respiratory syncytial virus A, Respiratory syncytial virus B, Rhinovirus.

B. bronchiseptica Strain 4617 and *B. holmesii* were tested at 1.0 x 10⁶ CFU/mL and were found to react with the Alethia Pertussis assay.

Unexpected results were observed during original testing of specimens containing *B. hinzii*, *H. parainfluenzae* and *M. genitalium*. One of three negative sample replicates containing *B. hinzii* produced a false-positive result that was not confirmed with further testing (20/20 replicates). Three of three negative sample replicates containing *H. parainfluenzae* produced false-positive results that were not confirmed with further testing (20/20 replicates). Three of three *B. pertussis* positive sample replicates produced invalid results that were not confirmed with further testing (10/10 replicates). Three of three negative sample replicates containing *M. genitalium* produced invalid results that were not confirmed with further testing (10/10 replicates). As repeat testing using a heightened number of replicates did not confirm original results, *B. hinzii*, *H. parainfluenzae* and *M. genitalium* are not considered cross-reactive or interferents in Alethia Pertussis testing.

TESTS FOR INTERFERING SUBSTANCES

The following chemical substances, at the specified saturated solvent/diluent concentrations, do not interfere with test results: Acetaminophen (10 mg/mL), Advil® [ibuprofen (10 mg/mL)], Afrin® Decongestant Nasal Spray [Oxymetazoline hydrochloride (0.0005% w/v)], Albuterol Sulfate [salbutamol sulfate (1% w/v)], Aspirin (5 mg/mL), Coricidin® HBP Cold/Flu Tablets [Acetaminophen (3.26 mg/mL), Chlorpheniramine maleate (0.02 mg/mL), Diphenhydramine HCl (0.25 mg/mL), Erythromycin (2% w/v), Mupirocin (2% w/v), Petroleum Jelly [white petrolatum (1% w/v)], Robitussin® Cough+Chest Congestion DM Cough Syrup [dextromethorphan HBr (0.1 mg/mL), Guaifenesin (1.0 mg/mL)], Sudhrine PE [phenylephrine HCl (0.3 mg/mL)], Saline Nasal Spray [sodium chloride (0.0065% w/v)], Smokeless Tobacco (snuff) (1% w/v), Tobramycin (0.6 mg/mL), Vicks® VaporRub® [camphor (0.48% w/v), eucalyptus oil (0.12% w/v), menthol (0.26% w/v)].

Ibuprofen (10 mg/mL) produced invalid results (3/6 replicates) during original testing of contrived *B. pertussis* specimens. All repeat testing produced positive results (10/10 replicates). As the original results were not confirmed, Ibuprofen (10 mg/mL) is not considered an interfering substance.

Aspirin was found to interfere with Alethia Pertussis testing at concentrations greater than 5 mg/mL.

The following biological substances, at the specified saturated solvent/diluent concentrations, do not interfere with test results: Human DNA (200 ng/μL), Mucin [bovine submaxillary gland type I-S (1% w/v)], Whole blood (1% v/v).

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REFERENCES

- Nagamine K, Hase T, Notoni T. Accelerated reaction by loop-mediated isothermal amplifications using loop primers. *Mol Cell Probes* 2002; 16:223-29.
- Mori Y, Kitao M, Tomita N, Notoni T. Real-time turbidimetry of LAMP reaction for quantitating template DNA. *J Biochem Biophys* 2004; 59:145-47.
- Centers for Disease Control and Prevention. 2005 CDC Guidelines: Recommended Antimicrobial Agents for the Treatment and Postexposure Prophylaxis of Pertussis. 2005; 54(RR14); 1-16.
- Centers for Disease Control and Prevention. 2012 Final Pertussis Surveillance Report. 2013; 61(52).
- Klein NP, et al. Waning Protection after Fifth Dose of Acellular Pertussis Vaccine in Children. *N Engl J Med* 2012; 367: 1012-1019.
- Misegades LK, et al. Association of Childhood Pertussis with Receipt of 5 Doses of Pertussis Vaccine by Time Since Last Vaccine Dose, California, 2010. *JAMA* 2012; 308(20): 2126-2132.
- Andre P, et al. Comparison of Serological and Real-Time PCR Assays to Diagnose *Bordetella pertussis* Infection in 2007. *J Clin Micro* 2008; 46(5): 1672-1677.
- Centers for Disease Control and Prevention. Best Practices for Health Care Professionals on the use of Polymerase Chain Reaction (PCR) for Diagnosing Pertussis.
- US Department of Health and Human Services PHS/CDC/NIH. Biosafety in microbiology and biomedical laboratories, Washington DC: US Government Printing Office, 2007.
- CLSI: MM3-A2 Molecular diagnostic methods for infectious disease; approved guideline, 2nd ed. Wayne PA: Clinical Laboratory Standards Institute, 2006.
- Tatti, K., Sparks, K., Boney, K., Tondella, M. Novel Multitarget Real-Time PCR Assay for Rapid Detection of *Bordetella* Species in Clinical Specimens. *J Clin Micro* 2011; 49:4059-4066.

SN11025

REV. 10/18

Manufactured By

Meridian Bioscience, Inc.
Corporate Office
3471 River Hills Drive
Cincinnati, Ohio 45244 USA
Telephone: 513.271.3700
Orders/Customer Service:
800.543.1980
Technical Support Center:
800.343.3858
Information Fax: 513.272.5432
Ordering Fax: 513.271.0124

Authorized Representative

Meridian Bioscience Europe S. r. l.
Via dell'Industria, 7
20020 Villa Cortese, Milano
ITALY
Tel: +39 0331 43 36 36
Fax: +39 0331 43 36 16
Email: info@meridianbioscience.eu
WEB: www.meridianbioscience.com/eur

Meridian Bioscience Europe s.a./n.v.
2 Avenue du Japon - 1420 Braine l'Alleud
BELGIUM
Tel: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
Email: info.bn@meridianbioscience.eu






Meridian Bioscience Europe France
34 rue de Ponthieu - 75008 Paris
FRANCE
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 9 70 06 62 10
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
Postbus 301 - 5460 AH Veghel
NETHERLANDS
Tel: +31 (0) 411 62 11 66
Fax: +31 (0) 411 62 48 41
Email: info.bn@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product.

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 90/79/EEC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 90/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 90/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 90/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostika 90/79/EEG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdüpfungspuffer befindet
		STERILE R	Sterilization by gamma irradiation / Sterilizzazione con raggi gamma / Sterilisation par irradiation aux rayons gamma / Esterilizado por irradiación gamma / Sterilisation durch Gammastrahlen
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	RoHS	Restriction of Hazardous Substances / Restrizione all'uso di sostanze pericolose / Limitation de substances dangereuses / Restricción de Substancias Nocivas / Beschränkung der Verwendung bestimmter gefährlicher Stoffe
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		Caution, consult accompanying documents / Attention, vedere le istruzioni per l'uso / Attention voir notice d'instructions / Atención, ver instrucciones de uso / Achtung, Begleitdokumente beachten
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	STERILE EO	Sterilization by ethylene oxide / Sterilizzazione con ossido di etilene / Sterilisation par oxyde d'éthylène / Esterilizado por oxidido de etileno / Sterilisation durch Ethylenoxid
	Contains sufficient for n tests / Contenuto sufficiente per n saggi / Contenu suffisant pour n tests / Contenido suficiente para n ensayos / Inhalt ausreichend für n Prüfungen	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Temperature limitation / Limiti di temperatura / Limites de température / Límite de temperatura / Temperaturbegrenzung		ETL Registered Mark Certified / Marchio di certificazione registrato a livello nazionale / Certifié Conforme ETL / Marca de Certificación Registrada Nacional / ETL Konform beglaubigt
	Female / Femminile / De sexe féminin / Hembra / Frau		Recycle - do not dispose of as general waste / Riciclare - non eliminare come rifiuto generico / Recycler - ne pas jeter dans une poubelle / Recycle - no deshecho como basura general / Recycling- dieses Produkt nicht über den Hausmüll entsorgen
	Single Use Only / Prodotto Monouso / A usage unique / Para Un Solo Uso / nur für die einmalige Anwendung		Heat Treatment Tube / Provetta per il Trattamento Termico / Tube pour le traitement thermique / Tubo de tratamiento de calor / Röhren zur Hitzebearbeitung
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät		For IVD Performance Evaluation Only / Solamente per valutazione delle prestazioni / Réactif IVD réservé à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum		HOT SURFACE: Keep hands away from Hot Surfaces / Superficie calda: tenere le mani lontane dalla superficie calda / SURFACES CHAUDES: Ne pas toucher les surfaces chaudes / Superficie Caliente: Mantenga las manos alejadas de la superficie caliente / Heiße Oberfläche: Kontakt mit heißen Oberflächen vermeiden
	CAUTION: Laser Radiation / ATTENZIONE: Radiazione Laser / AVERTISSEMENT: Rayonnement laser / Precaución: Radiación Laser / WARNUNG: Laserstrahlung	IPX-0	CAUTION: Protect from water / ATTENZIONE: Proteggere dall'acqua / AVERTISSEMENT: Protéger de l'humidité / Precaución: Proteja del agua / WARNUNG: Vor Feuchtigkeit schützen
	CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risques de danger / Precaución: Peligroso / WARNUNG: Risikogefahr	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
MIN OIL	Mineral Oil / Olio Minerale / Huile Minérale / Aceite Mineral / Mineralöl	MEDIA	Mediu / Terreno di trasporto / Milieux / Medio / Medium
ST TUBE	Screw Top Tube / Provetta con tappo a vite / Tube à bouchon vissé / Tubo con tapa de rosca / Röhren mit Schnappverschluss	COL	Sample Preparation Column / Colonna di preparazione del campione / Colonne pour la préparation de l'échantillon / Columna de preparación de muestra / Säule zur Probenaufarbeitung
BUF SMP	Sample Buffer / Soluzione tampone per il campione / Tampón de échantillon / Tampón de muestra / Probenpuffer	PRE REAG	Pretreatment Reagent / Reagente di Pretattamento / Réactif de prétraitement / Reactivo de pretreatment / Reagenz für die Vorbehandlung
Rx Only	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig	SMP PREP	Sample Preparation / Preparazione del campione / Préparation de l'échantillon / Preparación de Muestra / Probenvorbereitung
TUBE	Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß		

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.