# alethia®

## Malaria

### DNA Amplification Assay

#### DNA Amplification Assays for the Detection of *Plasmodium* sp.

REF 480925

IVD

Rx Only

#### INTENDED USE

The Alethia Malaria DNA amplification assays, performed on the Alethia Incubator/Reader, are qualitative in vitro diagnostic tests for the direct detection of *Plasmodium* sp. DNA in human venous EDTA whole blood specimens from individuals with signs and symptoms of malarial infection. Results from Alethia Malaria assays are intended to be used as an aid in the diagnosis of human malaria infection.

Alethia Malaria assays utilize loop-mediated isothermal DNA amplification (LAMP) technology to detect *Plasmodium* sp. DNA by targeting segments of the *Plasmodium* genome. Alethia Malaria assays do not distinguish between *Plasmodium* species.

Alethia Malaria is intended for use in hospital, reference or state laboratory settings. The devices are not intended for nonlaboratory point-of-care use.

#### SUMMARY AND EXPLANATION OF THE TEST

The Alethia Malaria molecular assays are based on loop-mediated amplification (LAMP)<sup>1,2</sup> technology. The assays target a region of the *Plasmodium* genome that is conserved across *Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium ovale*, *Plasmodium malariae*, and *Plasmodium knowlesi*. The Alethia Malaria assays target a 214 base pair (bp) sequence of a *Plasmodium* sp. mitochondrial DNA noncoding region.

Loop-mediated amplification uses specially designed primers to provide for specific and continuous isothermal DNA amplification. A by-product of amplification is magnesium pyrophosphate, which forms a white precipitate leading to a turbid reaction solution. Reaction solution absorbance characteristics are monitored by the Meridian Alethia Incubator/Reader Incubator/Reader. Changes in reaction solution absorbance characteristics created by precipitation of magnesium pyrophosphate indicate the presence of target DNA. The absence of target DNA results in no significant change in sample absorbance.

The Alethia Malaria kit includes Alethia Malaria Test Devices, Sample Preparation Apparatus IV (SMP PREP IV), Alethia Buffer I, and tubes. Whole blood specimens are first treated with Buffer I; cells are disrupted and nucleic acids are released. The lysate is added to SMP PREP IV and filtered into a clean Tube I by gently squeezing. Sample effluent collected after filtration contains nucleic acids for use with the Alethia Malaria Test Device.

The Alethia Malaria Test Device contains one lyophilized amplification reagent bead in each of two chambers: a TEST chamber with *Plasmodium* sp.-specific primers and a CONTROL chamber with human mitochondrial DNA-specific primers. Human mitochondrial DNA in whole blood samples and the human mitochondrial DNA-specific primers in the Test Device CONTROL chambers function as the Internal Control for the assay. During specimen preparation, human mitochondrial DNA is extracted with the *Plasmodium* sp. mitochondrial DNA to allow for parallel processing of target DNA and Control DNA through amplification and detection. The Internal Control monitors DNA extraction, amplification inhibition, assay reagent performance and sample processing effectiveness. The Control target must be amplified and detected in the final reaction or the test is considered invalid and results are not reported.

The Alethia Incubator/Reader monitors changes in absorbance characteristics by measuring transmission of light through the Test and Control reaction solutions. Light transmission is checked at the assay Run Start (Signal<sub>initial</sub>, S<sub>i</sub>) and at the assay Run End (Signal<sub>final</sub>, S<sub>f</sub>). The Alethia Incubator/Reader calculates the change in light transmission between Run End and Run Start (S<sub>f</sub>:S<sub>i</sub>) and compares the ratio to a fixed cut-off value.

Fixed cut-off values for the TEST chamber are used to report sample results. TEST chamber S<sub>f</sub>:S<sub>i</sub> ratios less than 70% are reported as 'POSITIVE'; TEST chamber S<sub>f</sub>:S<sub>i</sub> ratios greater than or equal to 70% are reported as 'NEGATIVE'. Numerical values are not reported.

Fixed cut-off values for the CONTROL chamber are used to determine validity. CONTROL chamber S<sub>f</sub>:S<sub>i</sub> ratios less than 85% are considered valid and allow for reporting of TEST chamber results (POSITIVE, NEGATIVE). CONTROL chamber S<sub>f</sub>:S<sub>i</sub> ratios greater than or equal to 85% are considered invalid and prevent reporting of TEST chamber results. Invalid CONTROL chamber reactions are reported as 'INVALID'. Numerical values are not reported.

More stringent cut-off criteria are applied to the CONTROL chamber reaction to ensure amplification is not inhibited, reagents are performing as intended and that sample processing was performed appropriately.

#### BIOLOGICAL PRINCIPLES

Malaria is a public health challenge world-wide, with an estimated 3.2 billion people at risk for infection in 97 countries.<sup>3</sup> The World Health Organization (WHO) estimates 198 million cases occurred in 2013 that led to approximately 584,000 deaths.<sup>3</sup> Malaria is caused by the intracellular parasite *Plasmodium* and is transmitted to humans by female *Anopheles* mosquitoes.<sup>3</sup> Once humans are infected, *Plasmodium* exhibit multiple life stage morphologies for initial growth in liver cells, followed by invasion of red blood cells.<sup>4</sup> *Plasmodium* replicate within red blood cells, causing cell rupture and clinical symptoms of the disease.<sup>4</sup>

There are five *Plasmodium* species known to cause malaria in humans (*Plasmodium falciparum*, *Plasmodium vivax*, *Plasmodium malariae*, *Plasmodium ovale*, and *Plasmodium knowlesi*), with the majority of infections caused by *P. falciparum* and *P. vivax*. Malaria treatment is dependent on the *Plasmodium* species, clinical status of the patient, and antimicrobial susceptibility.<sup>5</sup> Infections caused by *P. falciparum* and *P. knowlesi* are typically more severe and require rapid treatment.<sup>5</sup>

The WHO recommends diagnostic testing for all individuals with suspected malaria.<sup>3</sup> Thin and thick smear microscopy is currently the gold standard diagnostic method. However, accuracy and precision of microscopy is highly dependent on the technician skill (limit of detection ranging from 4-100 parasites/ $\mu$ L or higher), with false negative and false positive rates between 10-25% and approximately 2-3 fold variability in quantitation.<sup>6,7,8</sup> In addition, errors in the *Plasmodium* species identification are common in routine microscopy examination.<sup>9</sup> Immunoassay-based rapid diagnostic tests (RDTs) are also an available diagnostic method. Quality assessment testing performed by the WHO demonstrated that at low parasite densities (200 parasites/ $\mu$ L), only 76% and 42% of available RDTs could detect *P. falciparum* and *P. vivax*, respectively.<sup>3</sup>

The Alethia Malaria assays provide rapid results with higher sensitivity than microscopy and RDTs.

#### REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. **Alethia Malaria Test Device:** Two-chambered device containing lyophilized amplification reagents (DNA polymerase, deoxynucleotide triphosphates) and either *Plasmodium* sp.-specific primers (TEST Chamber) or human mitochondrial DNA-specific primers (CONTROL Chamber).
2. **Alethia Buffer I:** Lysis solution containing 0.2N sodium hydroxide.
3. **Alethia Sample Preparation Apparatus IV (SMP PREP IV):** Tris-buffer solution containing 0.09% azide as a preservative.
4. **Alethia Tube I:** 1.5 mL tubes

#### MATERIALS PROVIDED SEPARATELY

1. Alethia Malaria External Control Kit, Meridian Bioscience, Inc. Catalog Number: 479970

#### MATERIALS NOT PROVIDED

1. Disposable latex gloves, powder free
2. DNase/RNase-free, aerosol resistant pipette tips
3. Blood collection tubes with EDTA anticoagulant

#### EQUIPMENT NOT PROVIDED

1. Interval Timer
2. Vortex mixer (optional)
3. Micropipette capable of dispensing 50  $\mu$ L
4. Alethia Incubator/Reader, Meridian Bioscience, Inc. Catalog Number: 610189

#### PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Do not interchange kit reagents and Test Devices between lots. Individual lots of Tube I and ST Tubes are interchangeable within each kit type provided they are within the assigned expiration date when used.
3. The Alethia Incubator/Reader may produce incorrect results if the Malaria assay program is not used for testing.
4. Follow Biosafety Level 2 and Good Laboratory practices during testing.<sup>9</sup> Treat all specimens and used Test Devices as capable of transmitting infectious agents. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
5. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
6. Quality Control Programs for Molecular Testing Laboratories, including proper use and care of equipment, should be employed.<sup>10</sup>
7. The Alethia Malaria Test Devices contain lyophilized reagents. The protective pouches should not be opened until ready to perform the assay.
8. The Alethia Malaria Test Devices include a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
9. Dispose of used Alethia Malaria Test Devices and tubes immediately after processing. Leave the Test Device latch securely in place. Do NOT open the Test Device after processing. Opening the device after amplification may result in contamination of the test area with amplification product.

## HAZARD AND PRECAUTIONARY STATEMENTS

 Buffer I	<b>Signal Word</b> Warning <b>Hazard Statements</b> H315 – Causes skin irritation H319 – Causes serious eye irritation <b>Precautionary Statements - EU (§28, 1272/2008)</b> P280 – Wear protective gloves/protective clothing/eye protection/face protection
--	---

## SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-30 C.

## REAGENT PREPARATION

Ensure kit reagents are at room temperature (19-30 C) before use. Ensure Buffer I has been warmed to room temperature completely prior to use and no precipitate is visible. Incorrect results may be obtained if reagents are not brought to room temperature prior to use.

## SPECIMEN COLLECTION AND PREPARATION

**Sample type:** Human venous whole blood samples with EDTA as a preservative.

**Sample Collection:** Venous whole blood samples should be collected into a specimen tube containing EDTA according with institutional guidelines for collection of clinical specimens for malaria infection. Ensure the blood collection tube is inverted immediately after collection at least 8 times, or according to the manufacturer's instructions.

EDTA venous whole blood samples may be stored at 2-30 C after collection and during transportation to the laboratory. Samples should be tested as soon as possible, but may be stored for up to 7 days at room temperature (19-30 C) or up to 14 days refrigerated (2-8 C) prior to testing. Samples that will not be tested within this time frame should be frozen immediately at -20 C for up to 30 days until tested. Samples may be frozen and thawed 2 times after storage at -20 C prior to testing with the Alethia Malaria assays.

## SPECIMEN PREPARATION

**NOTE:** Ensure that the Alethia Incubator/Reader is powered on and required performance verifications have been completed prior to initiation of SPECIMEN PREPARATION. Refer to the Alethia Incubator/Reader Operator's Manual for further information regarding instrument set-up and operation.

**NOTE:** Ensure specimens are at room temperature (19-30 C) before specimen preparation.

### 1. Alethia Malaria Specimen Preparation

- Invert EDTA whole blood specimen 2-3 times to mix.
- Add 50 µL of the collected venous whole blood sample (with EDTA) to one tube of Buffer I. Mix by inversion 5 times or by vortexing for approximately 10 seconds. Hold the sample for 2 minutes.
- Mix by inversion 5 times or by vortexing for approximately 10 seconds and immediately transfer 50 µL of lysate into SMP PREP IV. Mix by inversion 5 times or by vortexing for 10 seconds.
- Gently squeeze the SMP PREP IV and slowly collect 5 to 10 drops into a clean Tube I. Visually verify that the eluate is tinted red to reddish brown. Label the tube with the specimen identification and proceed to the Test Procedure.

## TEST PROCEDURE

**NOTE:** A maximum of 10 samples can be processed in a single Alethia Incubator/Reader run.

- Remove 1 Alethia Malaria Test Device from its protective pouch per sample. Carefully open the device by lifting the latch, holding the chambers such that the lyophilized reagents will not fall out upon opening. Place the device on a flat surface or in a rack that can accommodate the device.
- Using a micropipette, transfer 50 µL of the sample to both the TEST (Left/White Bead) and CONTROL (Right/Yellow Bead) chambers of the Alethia Malaria Test. Take care to not introduce air to the reaction mixture. Do not mix reactions with pipette.
- Close the Alethia Test Device and fasten the latches securely.
- Tap device(s) on the bench top or mix to remove air bubbles. Carefully examine the Test Device(s) for rehydration of the Control/Test Bead, for air bubbles left in the chamber and liquid in the top of the device. If undissolved beads, air bubbles or liquid in the top of the device are noted, tap the device on the bench top and repeat visual inspection. Amplification and detection should be initiated within 15 minutes.
- Repeat Test Procedure Steps for all samples to be tested.
- Insert the Alethia Test Devices into the Alethia Incubator/Reader and initiate run **using the Malaria Program**. Results will be displayed at the conclusion of the run.

## INTERPRETATION OF RESULTS

Sample ID	Reported Result	Interpretation
Patient Specimen	POSITIVE	Sample contains <i>Plasmodium sp.</i> target DNA.
	NEGATIVE	No <i>Plasmodium sp.</i> DNA detected.
	INVALID	<b>No reportable result. Repeat the test using the original sample.</b> Inhibitory patient specimen, improper sample preparation, reagent failure, instrument failure or internal control failure.
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, Alethia Incubator/Reader performing correctly.
	NEGATIVE	<b>Incorrect control result.</b> Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	INVALID	<b>No reportable result. Repeat entire assay run using original samples.</b> Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	<b>Incorrect control result.</b> Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, Alethia Incubator/Reader performing correctly.
	INVALID	<b>No reportable result. Repeat entire assay run using original samples.</b> Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No Alethia Test Device in the Alethia Incubator/Reader Well. <b>OR</b> The Alethia Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. <b>Repeat the test using original sample.</b>

## QUALITY CONTROL

**This test should be performed per applicable local, state, or federal regulations or accrediting agencies.**

- Each device contains an internal control that controls for amplification inhibition, assay reagents, DNA preparation, and sample processing effectiveness. Human mitochondrial DNA, which serves as the internal control DNA, is isolated from the whole blood sample and processed through all steps of the procedure. Primers for amplification of internal control DNA are present in the Control Chamber of the Alethia Test Device.
- Good laboratory practice recommends the use of control materials. Users should follow the appropriate federal, state and local guidelines concerning the running of external quality controls.
- Alethia Malaria External Positive Control Reagent is supplied separately (Catalog 479970). Alternatively, previously characterized clinical or contrived *Plasmodium sp.* positive blood samples can be used as an external positive control. A qualified negative human whole blood sample with EDTA may be used as an external negative control. It is recommended that the reactivity of each new lot and each new shipment of Alethia Malaria be verified on receipt and before use. External control tests should be performed thereafter in accordance with appropriate federal, state and local guidelines. The Alethia Malaria test kit should not be used in patient testing if the external controls do not produce the correct results.
- A separate Test Device must be used for each external control reagent.

## EXPECTED VALUES

The incidence of *Plasmodium sp.* as detected by the Alethia Malaria assay during the 2015 clinical study was 68.1% (147/216) using the Alethia Malaria sample preparation method. The overall incidence of each *Plasmodium* species is provided below.

Species	Prevalence by <i>Plasmodium</i> species	
	Alethia Malaria (n=216)	
	Total Positive	Prevalence
<i>P. falciparum</i>	137	63.4%
<i>P. vivax</i>	0	0%
<i>P. ovale</i>	1	0.5%
<i>P. malariae</i>	1	0.5%
<i>P. knowlesi</i>	0	0%
Unknown*	8	3.7%

\*Samples were negative by microscopy and therefore species could not be determined.

## LIMITATIONS OF THE PROCEDURE

- This product can only be used with the Alethia Incubator/Reader instrument.
- Alethia Malaria do not distinguish between *Plasmodium* species.
- Hematocrit (>57.5%) and hemoglobin (>30 g/dL) above normal physiological levels may produce invalid results with the Alethia Malaria assays.
- Performance of the Alethia Malaria assays has not been established for packed red blood cell specimens.
- Performance of *Plasmodium knowlesi* with the Alethia Malaria assay was established using purified genomic DNA only; whole organism testing was not performed.
- Alethia Malaria are qualitative assays and do not provide quantitative values or information about organism load.
- The detection of nucleic acids is dependent upon proper specimen collection, handling, transportation, storage and preparation. Failure to observe proper procedure in any one of these steps can lead to incorrect results.
- Organism nucleic acid may persist in vivo, independent of organism viability. Alethia Malaria do not distinguish between viable and nonviable organisms.
- As with all molecular based diagnostic tests, (A) False negative results may occur from the presence of inhibitors, technical error, sample mix-up or low numbers of organisms in the clinical specimen; (B) False positive results may occur from the presence of cross-contamination by target organisms, their nucleic acids or amplified product, and from non-specific signals.

### SPECIFIC PERFORMANCE CHARACTERISTICS

The performance characteristics of the Alethia Malaria DNA Amplification Assays were established in clinical studies conducted in November 2015 in Senegal, Africa. Performance characteristics of the assay were compared to the reference method, thin and thick film microscopy. *Plasmodium* species identification for all positive samples was determined by microscopy.

Samples included prospective and retrospective venous whole blood specimens with EDTA from patients with signs and symptoms of *Plasmodium* infection. A total of 216 eligible, de-identified whole blood specimens were evaluated. 66 specimens were tested prospectively and 150 specimens were prospectively-collected and stored frozen prior to testing by the Alethia Malaria assay (retrospective). Each specimen was prepared by Alethia Malaria sample preparation methods. All samples were tested prospectively by microscopy.

Analysis of assay performance with prospective and retrospective specimens indicated there is no difference in performance for fresh and frozen specimens with the Alethia Malaria assays. The tables below summarize the performance of the Alethia Malaria DNA Amplification Assays for both prospective and retrospective samples combined.

#### Alethia Malaria Performance Compared to Microscopy

	Microscopy			Alethia	Performance				
	Pos	Neg	Total	INV <sup>a</sup>				95% CI	
Alethia Malaria	Pos	139	8	147	0(1)	Sensitivity	100.0%	139/139	97.3-100.0%
	Neg	0	67	67	2(2)	Specificity	89.3%	67/75	80.3-94.5%
	Total	139	75	214	2				

<sup>a</sup>Initial invalid results are reported within the parentheses. The final number of invalid samples remaining after repeat testing is shown before the parenthesis.

Samples were collected from patients ranging in age from 4 to 77 years. 15 were between 1 – 12 years of age; 75 were between 13 – 21 years of age; 125 were ≥22 years old; the age of one patient was not defined. There was no noted performance difference based on age. The study population included 84 (38.9%) female and 131 (60.6%) male patients; the gender of one patient was not defined. There is no expectation that assay performance is influenced by gender.

### ANALYTICAL SENSITIVITY

Limit of Detection (LoD) was determined using *Plasmodium falciparum* and *Plasmodium vivax*. The LoD was confirmed using 20 replicates and a stated probability (eg, 95%, where 19/20 replicates are positive) of obtaining positive responses.

	Alethia Malaria
<b>Plasmodium Species</b>	<b>Parasites/μL</b>
<i>P. falciparum</i> (3D7)	2
<i>P. vivax</i> (India VII)	0.125

### ASSAY REACTIVITY

Quantified stocks of *P. malariae*, *P. ovale*, and *P. knowlesi* DNA were diluted in negative whole blood to approximately three times the limit of detect for Alethia Malaria sample preparation method (approximately 72 copies/μL or 6 parasites/μL). All species reacted at the concentrations tested.

### REPRODUCIBILITY

Reproducibility studies were carried out by three internal laboratories. Blind-coded panels of 10 samples were supplied to participating laboratories. The panels included contrived *Plasmodium falciparum* (strain 3D7) samples manufactured as moderate positive samples, low positive samples, or high negative samples. The panel also included one negative whole blood sample. Testing was performed by at least 2 different operators at each laboratory on the same day (intra-assay variability) for five days (inter-assay variability). Three lots each of Alethia Malaria and 8 Alethia Incubator/Reader instruments were used in this study. External Positive and Negative Controls were tested with each panel; a qualified negative blood sample was used as the Negative Control.

Reproducibility Study Summary: Alethia Malaria Sample Preparation Method								
Sample Type	Site 1		Site 2		Site 3		Total	
	Percent Agreement		Percent Agreement		Percent Agreement		Percent Agreement	
Moderate Positive (8 parasites/μL)	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
Low Positive (2 parasites/μL)	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
High Negative (0.458 parasites/μL)	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
Negative	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
Negative Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
Positive Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%

### CROSSREACTIVITY STUDIES

Crossreactivity studies employed positive (*P. falciparum* 3D7) and negative whole blood specimens inoculated with bacterial or fungal organisms to a minimum concentration of  $1.0 \times 10^5$  CFU/mL, virus at a minimum of  $1.0 \times 10^5$  TCID<sub>50</sub>/mL, or protozoans to a minimum concentration of  $1.0 \times 10^5$  organisms/mL. Where whole organisms were not available,  $1.0 \times 10^5$  copies/mL for genomic DNA was tested. None of the following organisms or their genetic material reacted with the Alethia Malaria assays:

*Babesia microti*, *Borrelia burgdorferi*, *Haemophilus influenzae*, *Klebsiella pneumoniae*, *Leishmania donovani*, *Leishmania infantum*, *Leptospira interrogans*, *Mycobacterium smegmatis*, *Salmonella typhi*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Streptococcus pneumoniae*, *Treponema pallidum*, *Trypanosoma cruzi*, *Trypanosoma rangeli*, *Vibrio parahaemolyticus*, Cytomegalovirus (CMV), Epstein-Barr virus (EBV), Hepatitis B virus, Hepatitis C virus, Herpes simplex virus 1 (HSV 1), HIV-1, Human Papilloma Virus (HPV), and Rubella virus.

The following organisms or their DNA could not be obtained for testing and were evaluated through in silico analysis. None of the following organisms are expected to react with the Alethia Malaria assays based on this analysis: *Anaplasma phagocytophilum*, *Clostridium botulinum*, *Orientia tsutsugamushi*, *Rickettsia conorii*, *Rickettsia rickettsii*, *Rickettsia typhi*, Chikungunya virus, Dengue virus (types 1-4), West Nile virus, and Yellow Fever virus.

Human genomic DNA was nonreactive at  $1.0 \times 10^6$  copies/mL.

### TESTS FOR INTERFERING SUBSTANCES

Interference testing was performed in the presence of worst-case concentrations of chemical and biological substances introduced directly into contrived low positive (*P. falciparum* 3D7) and negative whole blood samples. No interference was observed with the following substances for either Alethia Malaria sample preparation methods:

Acetaminophen (1 mg/mL), Amoxicillin (0.1 mg/mL), Artesunate (1 mg/mL), Aspirin (1 mg/mL), Atovaquone (0.1 mg/mL), Cephalexin (0.1 mg/mL), Chloroquine (1 mg/mL), Ciprofloxacin (0.1 mg/mL), Clindamycin (1 mg/mL), Doxycycline hydrochloride (1 mg/mL), Erythromycin (0.1 mg/mL), Hydroxychloroquine sulfate (1 mg/mL), Ibuprofen (1 mg/mL), Lumefantrine (1 mg/mL), Mefloquine (1 mg/mL), Primaquine phosphate (1 mg/mL), Proguanil hydrochloride (1 mg/mL), Pyrimethamine (1 mg/mL), Quinine sulfate (1 mg/mL), Sodium Citrate (0.11 M), elevated Bilirubin (>0.15 mg/dL), elevated leukocytes (buffy coat) >10% v/v, serum albumin (>0.03 g/dL), and triglycerides (>9.9 mg/dL).

Hematocrit (>57.5%) and hemoglobin (>30 g/dL) above normal physiological levels may produce invalid results with the Alethia Malaria assays.

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**REFERENCES**

1. Nagamine K, Hase T, Notoni T. Accelerated reaction by loop-mediated isothermal amplifications using loop primers. Mol Cell Probes 2002;16:223-29.
2. Mori Y, Kitao M, Tomita N, Notoni T. Real-time turbidimetry of LAMP reaction for quantitating template DNA. J Biochem Biophys 2004;59:145-47.
3. World Malaria Report 2014. Global Malaria Programme, World Health Organization. [http://www.who.int/malaria/publications/world\\_malaria\\_report/en/](http://www.who.int/malaria/publications/world_malaria_report/en/).
4. Cowman AF and Crabb BS. Invasion of Red Blood Cells by Malaria Parasites. Cell 2006; 124: 755-766.
5. Treatment of Malaria (Guidelines for Physicians). Centers for Disease Control and Prevention, July 2013. [www.cdc.gov/malaria/resources/pdf/clinicalguidance.pdf](http://www.cdc.gov/malaria/resources/pdf/clinicalguidance.pdf).
6. O'Meara WP, Barcus M, Wongsrichanalai C, Muth S, Maguire JD, Jordan RG, Prescott WR, and McKenzie FE. Reader Technique as a Source of Variability in Determining Malaria Parasite Density by Microscopy. Malaria J 2006;5:118.
7. McKenzie FE, Sirichaisinthop J, Miller RS, Gasser RA, and Wongsrichanalai C. Dependence of Malaria Detection and Species Diagnosis by Microscopy on Parasite Density. Am J Trop Med Hyg 2003;69(4):372-376.
8. Wongsrichanalai C, Barcus MJ, Muth S, Sutamihardja A, and Wernsdorfer WH. A Review of Malaria Diagnostic Tools: Microscopy and Rapid Diagnostic Test (RDT). Am J Trop Med Hyg 2007;77(6):119-127.
9. US Department of Health and Human Services PHS/CDC/NIH. Biosafety in microbiology and biomedical laboratories, Washington DC: US Government Printing Office, 2007.
10. CLSI: MM3-A2 Molecular diagnostic methods for infectious disease; approved guideline, 2nd ed. Wayne PA: Clinical Laboratory Standards Institute. 2006.

SN11035

REV. 04/23

 <p>Manufactured By</p>	<p><b>Meridian Bioscience, Inc.</b>          3471 River Hills Drive          Cincinnati, OHIO - 45244 USA  <a href="http://www.meridianbioscience.com">www.meridianbioscience.com</a></p> <p><u>Contacts:</u>          Main Telephone (+1) 513.271.3700          Customer Service/Orders 800.543.1980          Technical Support Center 800.343.3858          Information Fax: 513.272.5432          Ordering Fax: 513.271.0124          E-mail: <a href="mailto:info@meridianbioscience.com">info@meridianbioscience.com</a></p>
 <p>Authorized Representative</p>	<p><b>Meridian Bioscience Europe, SRL</b>          Via Dell'Industria 7, 20035 Villa Cortese          (Milano) ITALY  <a href="http://www.meridianbioscience.com">www.meridianbioscience.com</a></p> <p><u>Contacts:</u>          Main Telephone (+39) 0331.433636          E-mail: <a href="mailto:info@meridianbioscience.eu">info@meridianbioscience.eu</a>          Technical Support: <a href="mailto:MBE-TechService@meridianbioscience.eu">MBE-TechService@meridianbioscience.eu</a>          Customer Service/Orders:</p> <ul style="list-style-type: none"> <li>• For Italian Customers:  <a href="mailto:ordini@meridianbioscience.com">ordini@meridianbioscience.com</a></li> <li>• For Distributors / International Customers:  <a href="mailto:Export.CustomerService@meridianbioscience.eu">Export.CustomerService@meridianbioscience.eu</a></li> </ul>
<p><b>AUSTRALIAN SPONSOR</b></p>	<p><b>Emergo Australia</b>          Level 20, Tower II          Darling Park          201 Sussex Street          Sydney, NSW 2000          Australia</p>
 <p>Authorized Representative</p>	<p><b>MedEnvoy Switzerland</b>          Gotthardstrasse 28          6302 Zug          Switzerland</p>

**INTERNATIONAL SYMBOL USAGE**

You may see one or more of these symbols on the labeling/packaging of this product:

**Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)**

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	<b>CONTROL +</b>	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
<b>LOT</b>	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	<b>CONTROL -</b>	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
<b>IVD</b>	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	<b>EC REP</b>	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	Meridian products carrying the European Conformity (CE) mark fulfill the requirements of Directive 98/79/EC or the Regulation 2017/746 on in-vitro diagnostic medical devices (I prodotti Meridian recanti il marchio di Conformità Europea (CE) soddisfano i requisiti della Direttiva 98/79/CE o del Regolamento 2017/746 sui dispositivi medico-diagnostici in vitro / Les produits Meridian portant la marque de Conformité européenne (CE) sont conformes aux exigences de la Directive 98/79/CE ou du Règlement 2017/746 relatifs aux dispositifs médicaux de diagnostic in vitro / Los productos de Meridian que llevan la marca de conformidad europea (CE) cumplen los requisitos de la Directiva 98/79/CE o del Reglamento 2017/746 sobre los productos sanitarios para diagnóstico in vitro / Produkte von Meridian mit der CE-Kennzeichnung erfüllen die Anforderungen der EU-Richtlinie 98/79/EG bzw. der Verordnung (EU) 2017/746 über In-vitro-Diagnostika	<b>SMP PREP DIL SPE</b>	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenzubereitung, in dem sich Probenverdünnungspuffer befindet
		<b>STERILE R</b>	Sterilization by gamma irradiation / Sterilizzazione con raggi gamma / Stérilisation par irradiation aux rayons gamma / Esterilizado por irradiación gamma / Sterilisation durch Gammastrahlen
<b>REF</b>	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	<b>RoHS</b>	Restriction of Hazardous Substances / Restrizione all'uso di sostanze pericolose / Limitation de substances dangereuses / Restricción de Substancias Nocivas / Beschränkung der Verwendung bestimmter gefährlicher Stoffe
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		Caution, consult accompanying documents / Attenzione, vedere le istruzioni per l'uso / Attention voir notice d'instructions / Atención, ver instrucciones de uso / Achtung, Begleitdokumente beachten
<b>SN</b>	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	<b>STERILE EO</b>	Sterilization by ethylene oxide / Sterilizzazione con ossido di etilene / Stérilisation par oxyde d'éthylène / Esterilizado por oxidode etileno / Sterilisation durch Ethylenoxid
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	<b>BUF RXN</b>	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung		ETL Registered Mark Certified / Marchio di certificazione registrato a livello nazionale / Certifé Conforme ETL / Marca de Certificación Registrada Nacional / ETL Konform beglaubigt
	Female / Femmine / De sexe féminin / Hembra / Frau		Recycle - do not dispose of as general waste / Riciclare - non eliminare come rifiuto generico / Recycler - ne pas jeter dans une poubelle / Recicle - no desecho como basura general / Recycling- dieses Produkt nicht über den Hausmüll entsorgen
	Single Use Only / Prodotto Monouso / A usage unique / Para Un Solo Uso / nur für die einmalige Anwendung		Heat Treatment Tube / Provetta per il Trattamento termico / Tube pour le traitement thermique / Tubo de tratamiento de calor / Röhrchen zur Hitzebearbeitung
<b>TEST</b>	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum		HOT SURFACE: Keep hands Away from Hot Surfaces / Superficie calda: tenere le mani lontane dalle superfici calde / SURFACES CHAUDES: Ne pas toucher les surfaces chaudes / Superficie Caliente: Mantenga las manos alejadas de la superficie caliente / Heiße Oberfläche: Kontakt mit heißen Oberflächen vermeiden
	CAUTION: Laser Radiation / ATTENZIONE: Radiazione Laser / AVERTISSEMENT: Rayonnement laser / Precaución: Radiación Laser / WARNUNG: Laserstrahlung	<b>IPX-0</b>	CAUTION: Protect from water / ATTENZIONE: Proteggere dall'acqua / AVERTISSEMENT: Protéger de l'humidité / Precaución: Proteja del agua / WARNUNG: Vor Feuchtigkeit schützen
	CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risques de danger / Precaución: Peligroso / WARNUNG: Risikogefahr	<b>CONTROL</b>	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
<b>BUF</b>	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
<b>MIN OIL</b>	Mineral Oil / Olio Minerale / Huile Minérale / Aceite Mineral / Mineralöl	<b>MEDIA</b>	Media / Terreno di trasporto / Milieux / Medio / Medium
<b>ST TUBE</b>	Screw Top Tube / Provetta con tappo a vite / Tube à bouchon vissé / Tubo con tapa de rosca / Röhrchen mit Schnappverschluss	<b>COL</b>	Sample Preparation Column / Colonna di preparazione del campione / Colonne pour la préparation de l'échantillon / Columna de preparación de muestra / Säule zur Probenaufarbeitung
<b>BUF SMP</b>	Sample Buffer / Soluzione tampone per il campione / Tampon de l'échantillon / Tampón de muestra / Probenpuffer	<b>PRE REAG</b>	Pretreatment Reagent / Reagente di Pretattamento / Réactif de prétraitement / Reactivo de pretratamiento / Reagenz für die Vorbehandlung
<b>Rx Only</b>	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig	<b>SMP PREP</b>	Sample Preparation / Preparazione del campione / Préparation de l'échantillon / Preparación de Muestra / Probenzubereitung
<b>TUBE</b>	Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß	<b>CH REP</b>	Swiss Authorized Representative / Mandatario svizzero / Mandataire Suisse / Representante Autorizado Suizo / Schweizer Bevollmächtigter

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.