merifluor®

VZV

REF 500113

IVD

Rx Only

INTENDED USE

Merifluor VZV is intended for direct detection of varicella-zoster virus (VZV), in vesicle smears and biopsy specimens, by the direct fluorescent antibody technique, in addition, the reagent can be used to confirm the presence of VZV in cell culture.

SUMMARY AND EXPLANATION OF THE TEST^{2, 3, 5}

VZV is a DNA virus of the *herpesvirus* group which causes chicken pox (varicella) and shingles (zoster). Varicella is the primary infection while zoster is caused by reactivation of latent virus. Although VZV is commonly identified by cell culture isolation, Schmidt, et. al., have demonstrated that cell culture detects only 23% of positive cases. In contrast, a direct fluorescent assay detected all positive cases tested.

BIOLOGICAL PRINCIPLES^{1, 3}

Merifluor VZV provides a direct fluorescent antibody stain (FA) for VZV identification. A monoclonal antibody (3B3), specific for VZV glycoprotein and conjugated to fluorescein isothiocyanate (FITC), is reacted with cells suspected of VZV infection. This results in a direct fluorescent label of VZV-infected cells which can be visualized using a fluorescence microscope.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. 3 mL prediluted FITC-conjugated monoclonal anti-VZV (3B3) antibody in 0.01M phosphate-buffered saline (PBS) containing 0.1% sodium azide, 1% carrier protein and Evans blue counterstain.
2. 5 mL FA Mounting Medium containing Tris-buffered glycerol with photobleach inhibitor.

MATERIALS NOT PROVIDED

The supplied staining reagent will permit identification of VZV. Additional materials and equipment needed for identification from direct smear specimens, cell culture or from biopsy specimens are as follows:

1. Acetone (reagent grade)
2. Phosphate-buffered saline (PBS), pH 7.4 ± 0.2, calcium and magnesium free
3. Clean wash bottle
4. Microscope slides with 6-8 mm wells
5. Coverslips, No. 1 thickness
6. Humidity chamber
7. Sodium hypochlorite, 0.5% (1:10 dilution of household bleach)
8. Fluorescence microscope equipped with a filter system suitable for fluorescein isothiocyanate (maximum excitation wavelength = 490 nm, mean emission wavelength = 520 nm). Lenses must be of high quality and have a final power of 100X to 400X.

NOTE: A well-functioning fluorescence microscope is essential to proper test interpretation. Variations in bulb wattage and alignment, illumination type (epi-or incident) and intensity, and filter type may affect performance. Use the positive control to verify adequate performance of the microscope.

9. Timer

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Sodium azide (NaN₃) is used in the reagent as a preservative and may react with lead or copper in drain lines to form explosive metal azides. Upon disposal, flush with a large volume of water to prevent azide build-up in drains.
3. Use a safety pipeting device for all pipeting. Never pipet by mouth.
4. All specimens, pipets and inoculated tissue culture tubes are to be considered potentially infectious and should be disposed of in 0.5% sodium hypochlorite or placed in biohazard bags and autoclaved.
5. Avoid microbial contamination of the reagents.

HAZARD and PRECAUTIONARY STATEMENTS

 Detection Reagent	<p>Signal Word Danger Hazard Statements H301 – Toxic if swallowed. H310 – Fatal in contact with skin H411 – Toxic to aquatic life with long lasting effects Precautionary Statements - EU (§28, 1272/2008) P301 + P310 – IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician P321 – Specific treatment (see .? on this label) P280 – Wear eye protection/ face protection P361 – Remove/Take off immediately all contaminated clothing P322 – Specific measures (see .? on this label)</p>
--	---

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2 to 8 C. Do not freeze or store in strong light. Reagents should not be warmed in a 37 C water bath. Bring to room temperature 18 to 25 C before use.

TEST PROCEDURE^{2, 5}

NOTE

1. Do not use reagents beyond expiration date.
2. Reagents are supplied ready to use. Do not dilute.
3. **Do not let slides dry during the staining procedure.**
4. Do not store the reagent or perform the staining procedure in strong light, such as direct sunlight.

Preparation of Smear Specimens

Incise vesicle peripherally with a sterile scalpel. Lift back the top of the lesion and remove excess fluid by gentle blotting with a swab. This swab can be used as a culture specimen. Scrape the base of the lesion thoroughly with a scalpel blade. Gross bleeding should be avoided. Spread cellular material collected on the edge of the blade *thinly* over the 6-8 mm diameter well area on a clean microscope slide and air dry.

Fix in fresh, cold acetone 2 to 8 C for 10 minutes at room temperature 18 to 25 C and air dry, or transport to laboratory for prompt fixation. Stain within two hours or store the fixed slide *desiccated* at 2 to 8 C.

Preparation of Tissue Culture Specimens

For staining tissue culture isolations, remove cells demonstrating typical VZV cytopathic effect (CPE) per routine tissue culture procedures. Drop cells onto a slide, preferably within a well size of 6-8 mm diameter. Air dry. Fix in fresh, cold acetone 2 to 8 C for 10 minutes at room temperature 18 to 25 C and air dry. Stain within two hours or store *desiccated* at 2 to 8 C.

Warning: Isolation of virus by tissue culture should be attempted only by those laboratories experienced in tissue culture procedures.

Preparation of Biopsy Specimens

Collect biopsy specimens using a punch biopsy or standard technique. Biopsy specimens may provide more reliable results than smear specimens.⁵ Biopsy specimens should be frozen, sectioned and air dried. Fix in fresh, cold acetone 2 to 8 C for 10 minutes at room temperature 18 to 25 C and air dry. Stain within two hours or store *desiccated* at 2 to 8 C.

Staining

The reagents are supplied prediluted. Do not dilute. Allow the reagents to warm to room temperature 18 to 25 C before use.

1. Add one to three drops of FITZ-conjugated monoclonal anti-VZV to well.
2. Incubate in humidity chamber at room temperature 18 to 25 C for 30 minutes.
3. Rinse each slide carefully using a squeeze bottle of PBS. BE SURE NOT TO DISRUPT CELLS WITH BUFFER STREAM.
4. Air dry (a small fan or blow dryer on cool setting can be used).
5. Use one drop of FA Mounting Medium. Use a coverslip for mounting slides.
6. Evaluate each slide with a fluorescence microscope at 100X to 400X magnification. Read slides immediately for best results. If necessary, slides may be stored in the dark at 2 to 8 C for up to 24 hours without significant fading of the fluorescent stain.

INTERPRETATION OF RESULTS

Direct Smears

Specimens exhibiting typical apple-green fluorescence in two or more cells in the entire smear are considered positive for VZV.

Specimens with fewer than one cell per high power field (400X), and not exhibiting positive fluorescence, should be reported as inadequate for testing.

Specimens with an adequate number of cells, but not exhibiting specific fluorescence in two or more cells are considered negative.

Trapping of conjugate under cells, or by debris and mucus, should not be misinterpreted. It is highly recommended that the user gain experience by extensive comparison of direct FA results with cell culture results (see LIMITATIONS OF THE PROCEDURE #1).

Cell Culture Isolation

Both positive and negative cells should be observed. Positive cells exhibit apple-green cytoplasmic fluorescence. Negative cells stain a faint red due to the counterstain. Culture tubes exhibiting few (less than five) foci of VZV-infected cells may not yield positive results due to inadequate sampling to prepare smear. Only samples with greater than five infected foci which demonstrate no specific staining should be reported negative.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations ro accrediting agencies.

As with all direct staining tests, specificity controls cannot be performed on the same smear. However, the reagent should be tested concurrently by staining known positive and negative cell culture monolayers. In addition, when confirming a virus isolate, uninfected cells from the same lot as those used to isolate the virus should be stained in the same manner as the patient sample.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

LIMITATIONS OF THE PROCEDURE

1. Because cell preparations or cytology smears are not always uniform monolayers, nonspecific reagent trapping can occur.
2. Bacterial or fungal contamination may nonspecifically bind the antibody and interfere with interpretation.
3. Isolation of varicella-zoster is dependent on the presence of live virus in the specimen. Proper collection, transportation and culture techniques are necessary for successful isolation.

SPECIFIC PERFORMANCE CHARACTERISTICS^{1, 3, 4, 6}

Vesicle smears

Vesicle smear specimens from 26 patients with chicken pox or zoster diagnosed by virus isolation or by serologic studies were stained using the monoclonal 3B3. All samples were positive for varicella-zoster antigen. Vesicle smears from four patients with HSV infections were positive when stained with anti-HSV antisera but negative when stained using 3B3.

Tissue Sections

Vesicle biopsy specimens were obtained from a patient with zoster four days after onset. Sections were frozen, sectioned and stained using 3B3 monoclonal antibody. Specific staining was observed in the basal and malpighian layers of the epidermis. No varicella-zoster antigens were detected in the outer keratin layer or in the underlying dermis.

Lung sections from a patient with interstitial pneumonitis caused by CMV were stained using the VZV monoclonal. No staining was seen. The same technique was positive when antibody to CMV was utilized.

Cell Culture

Human fibroblast cells infected with VZV-32, VZV-OKa, or five different clinical isolates, were stained using 3B3 in an indirect fluorescent antibody procedure. All strains of VZV exhibited positive results. Monolayers infected with HSV 1 (McIntyre), HSV 2 (MS) or CMV (AD-169) were stained with 3B3 and exhibited negative results. The following cell lines demonstrated no cross reactivity:

Human embryonic lung fibroblast	Transformed guinea pig (ATCC 104C1)
Human melanoma cells (Mewo strain)	VERO

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REFERENCES

1. Grose C, Edwards DP, Friedrichs WE, Weigle KA and McGuire WL. Monoclonal antibodies against three major glycoproteins of varicella-zoster virus. Infect Immun 1983 40:381-388.
2. Schmidt NJ, Gallo D, Devlin V, Woodie JD and Emmons RW. Direct immunofluorescence staining for detection of herpes simplex and varicella-zoster virus antigens in vesicular lesions and certain tissue specimens. J Clin Microbiol 1980 12 (5):651-655.
3. Weigle KA and Grose C. Common expression of varicella-zoster viral glycoprotein antigens in vitro and in chickenpox and zoster vesicles. J Infect Dis 1983 148(4):630-638.
4. Weigle KA and Grose C. Varicella pneumonia: immunodiagnosis with a monoclonal antibody. J Pediatrics 1984 August 105(2):265-9. 1979
5. Weller TH. 1979. Varicella and Herpes Zoster, in *Diagnostic procedures for viral, rickettsial, and chlamydial infections*, Chapter 12, Lennette EH and Schmidt NJ. (eds.), p. 375-398. American Public Health Association, Washington, D.C.
6. Unpublished study. Data on file, Meridian Bioscience, Inc.

SN11068

REV. 11/23

 Manufactured By	<p>Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com</p> <p><u>Contacts:</u> Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com</p>
	<p>Meridian Bioscience Europe, SRL Via Dell'Industria 7, 20035 Villa Cortese (Milano) ITALY www.meridianbioscience.com</p> <p><u>Contacts:</u> Main Telephone (+39) 0331.433636 E-mail: info@meridianbioscience.eu Technical Support: MBE-TechService@meridianbioscience.eu Customer Service/Orders: <ul style="list-style-type: none"> • For Italian Customers: ordini@meridianbioscience.com • For Distributors / International Customers: Export.CustomerService@meridianbioscience.eu </p>

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simblos, Zeichenerklärung)

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For technical assistance, call Technical Support Services at (800) 343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at (800) 543-1980.