


premier®


Cryptococcal Antigen

TEST PROCEDURE


How to perform the test




1. Detach microwells needed, place in holder. Add 1 drop of Positive Control, 50 μ L of Sample Diluent (Negative Control) and 50 μ L of patient sample to designated wells. Mix by gently shaking. Incubate at room temperature for 10 minutes.




2. Wash 4 times with 1X Wash Buffer. See package insert for proper washing procedure.




3. Add 1 drop of Enzyme Conjugate to all wells, and shake gently. Incubate at room temperature for 10 minutes.



4. Wash 4 times with 1X Wash Buffer. See package insert for proper washing procedure.



5. Add 2 drops of Substrate I to all wells. Mix by gently shaking for 15-20 seconds. Incubate at room temperature for 10 minutes.



6. Add 2 drops of Stop Solution I to all wells. Mix by gently shaking and wait 2 minutes before reading.

Semi-quantitative Assay

- Label tubes #1 through #5. In tube #1 combine 100 μ L specimen and 100 μ L of Sample Diluent.
- Place 200 μ L of Sample Diluent in tubes #2-#5. Transfer 50 μ L from tube #1 to tube #2 and continue this through tube #5.
- For increased precision, perform the procedure using a pipette in place of the dropper bottle tips.

Interpretation of Results:

Visual

Negative= colorless
Positive= definite yellow color

Spectrophotometric Single Wavelength (450 nm)

Negative= $OD_{450} < 0.100$
Indeterminant= $OD_{450} \geq 0.100$ and < 0.150
Positive= $OD_{450} \geq 0.150$

Dual Wavelength (450/630 nm)

Negative= $OD_{450/630} < 0.070$
Indeterminant= $OD_{450/630} \geq 0.070$ and < 0.100
Positive= $OD_{450/630} \geq 0.100$

EIA Titer =

Absorbance Value x Multiplication Factor

Example: The patient specimen has an absorbance of 1.2 at a 1:50 dilution.

$$1.2 \times 500 = 600$$

Which is reported as an EIA Titer of 1:600

Calculation Multiplication Factors =

Tube	1	2	3	4	5
Dilution	1:2	1:10	1:50	1:250	1:1250
Multiplication Factor	20	100	500	2500	12500



This illustration is representative of the current Package Insert at the time of publication. Please refer to the most current version of the Package Insert for complete instructions.

USA/CORPORATE OFFICE

3471 River Hills Drive
Cincinnati, Ohio 45244 USA
Telephone: 513-271-3700

MERIDIAN BIOSCIENCE EUROPE

Belgium/Luxembourg Tel: +32 (0)6789 5959
France Tel: +33 (0)1 4256 0440
Germany Tel: +49 (0)3371 60 222 31
Italy Tel: +39 0331 43 3636
Netherlands Tel: +31 (0)411 62 11 66
United Kingdom Tel: +44 (0)20 8453 7970

Orders/Customer Service: 1-800-543-1980
Technical Support: 1-800-343-3858
Information Fax: 513-272-5432
Ordering Fax: 513-271-0124
meridianbioscience.com

meridian BIOSCIENCE™