PREMIER® CRYPTOCOCCAL ANTIGEN

EIA for the Detection of Cryptococcal Antigen in Serum and CSF

REF 602096

IVD

Rx Only

INTENDED USE

The PREMIER Cryptococcal Antigen enzyme immunoassay (EIA) is a screening or a semi-quantitative test system for the detection of capsular polysaccharide antigens of *Cryptococcus neoformans* in serum and cerebrospinal fluid (CSF).

SUMMARY AND EXPLANATION OF THE TEST

Cryptococcosis is typically initiated as a subclinical pulmonary infection caused by the encapsulated soil yeast *C. neoformans*.¹ Although *C. neoformans* can infect any organ, the most serious complications occur when the organism invades the central nervous system of the immunocompromised patient.^{1,3} Four major serotypes (A-D) exist within the capsular polysaccharide antigen, with 80-85% of isolates being either A or B.⁴

Detection of cryptococcal antigens in serum or CSF is used successfully in the diagnosis of meningitis associated with infection by *C. neoformans*.^{1,3,7,9} When both serum and CSF are evaluated the sensitivity of antigen detection is greater than 90%.^{1,2,9} Antigen detection is more sensitive than the India ink mount^{8,9} and is a valuable aid in establishing a diagnosis when the culture is negative.⁹

Determining the relative changes in antigen levels has prognostic value and is used to monitor the efficacy of drug therapy.^{1,2,5,6}

The PREMIER Cryptococcal Antigen EIA detects and provides a reliable semi-quantitative determination of all cryptococcal antigen serotypes present in specimens.

BIOLOGICAL PRINCIPLES

The PREMIER Cryptococcal Antigen EIA utilizes anticytotoxic polyclonal antibodies adsorbed to microwells and a detection system based upon a monoclonal-peroxidase conjugate. If cryptococcal antigens are present in the sample, a complex is formed between the antigens, enzyme conjugate, and the adsorbed antibody. After washing to remove unbound conjugate, a substrate solution is added. Color develops in the presence of bound enzyme. The amount of color is related to the quantity of cryptococcal antigen present.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this kit is listed on the outer box.

- PREMIER Enzyme Conjugate** – Cryptococcal antigen specific monoclonal antibody conjugated to horseradish peroxidase in buffered protein solution containing 0.02% Thimerosal as a preservative.
- PREMIER 20X Wash Buffer II** – Concentrated wash buffer containing 0.2% Thimerosal as a preservative.
- PREMIER Substrate I** – Buffered solution containing peroxidase and tetramethylbenzidine.
- PREMIER Stop Solution I** – 1M Phosphoric Acid. **CAUTION:** Avoid contact with skin. Flush with water if contact occurs.
- PREMIER Positive Control** – Cryptococcal antigen in buffered protein solution containing 0.01% Thimerosal and 0.16% Gentamicin as preservatives.
- PREMIER Sample Diluent** – Buffered protein solution containing 0.02% Thimerosal as a preservative.
- PREMIER Antibody Coated Microwells** – Breakaway plastic microwell strips coated with polyclonal antibodies specific for cryptococcal antigen.
- Microwell Strip Holder**

MATERIALS NOT PROVIDED

- Pipettes capable of delivering 50 µL, 100 µL and 200 µL.
- Test tubes (12 x 75 mm) for dilution of sample.
- Distilled or deionized water.
- EIA plate reader capable of reading absorbance at 450 nm or 450/630.
- Squirt bottle or automatic EIA plate washer.
- Timer.
- Graduated cylinder for making 1X Wash Buffer.

PRECAUTIONS

- All reagents are for in vitro diagnostic use only.
- Patient specimens may contain infectious agents and should be handled and disposed of as potential biohazards. Wear disposable gloves while handling patient specimens and performing the test procedure.
- All reagents should be gently mixed and at room temperature before use.
- Do not interchange reagents from different kit lot numbers or use expired reagents.
- Hold reagent vials vertically to ensure proper drop size and delivery.
- Return caps to correct vials.
- The Positive Control contains inactivated cryptococcal antigen. This reagent, used wash buffer, used microwells, etc., should be disposed of as potentially hazardous material. The microwell plate holder should be disinfected and retained for future assays.
- Avoid skin contact with Stop Solution I (1M Phosphoric Acid). Flush with water immediately if contact occurs.
- Unused microwells must be placed back inside the foil ziplock pouch and sealed immediately. It is important to protect strips from moisture.**
- Avoid splashing when dispensing diluted sample into microwells by placing the pipette tip about halfway down the microwell and dispensing slowly down the side of the microwell.
- Perform microwell washing as directed in the assay procedure. **Inadequate washing may cause elevated background.**
- Incubation times and reagent concentrations have been optimized for sensitivity and specificity. Avoid any deviations from set incubation times or dilutions.
- Do not store specimens in a frost-free type freezer. Freezing and thawing of a specimen repeatedly can affect test results.
- Do not use microwells if the foil pouch appears damaged (ie, has holes, punctures).
- Do not use microwells where the desiccant indicator has changed from blue to pink.

HAZARD and PRECAUTIONARY STATEMENTS

Refer to the SDS, available at www.meridianbioscience.com for Hazards and Precautionary Statements

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-8 C and return the kit promptly to the refrigerator after each use.

REAGENT PREPARATION

- Bring the entire kit including microwell pouch to room temperature before use.
- Prepare enough 1X Wash Buffer for use by diluting 1 part 20X Wash Buffer with 19 parts purified water. The 1X Wash Buffer can be stored at room temperature for up to one month.

SPECIMEN COLLECTION AND PREPARATION

- Cerebrospinal Fluid**
 - Collect specimen aseptically following accepted procedures. (See PRECAUTION #2).
 - Centrifuge at 1000 xg for 15 minutes to ensure the removal of all white cells and particulate matter.
 - Carefully aspirate the CSF into a sterile container and seal.
 - Specimen may be tested immediately, stored up to 3 days at 2-8 C, or preserved by freezing. Do not store in a frost-free freezer.
- Serum**
 - Collect whole blood aseptically according to in-house procedures. (See PRECAUTION #2). The specimen must not contain anticoagulants as this may invalidate the test.

- Permit blood to clot for 10 minutes or more at room temperature in a collection tube.
- Centrifuge at 1000 xg for 10 minutes.
- Carefully aspirate the serum into a sterile container and seal.
- Specimen may be tested immediately, stored up to 3 days at 2-8 C, or preserved by freezing. Do not store in a frost-free freezer.

NOTE: Pretreatment of specimens is not recommended for the PREMIER Cryptococcal Antigen EIA. Studies indicate that pronase or heat treatment of specimens (56 C, 15 minutes) does not appear to adversely effect screening assay results. Occasionally, pronase or heat treatment alters EIA Titers. Therefore, if pronase treatment is used, it is essential that pronase be used for all serial specimens from an individual patient which are utilized to monitor therapy.

Specimens to be used for future comparisons should be aliquotted and frozen, as storage longer than 72 hours at 2-8 C may affect titers.

SCREENING ASSAY PROCEDURE

The screening assay procedure can be performed directly with serum or CSF prepared as above. The screening assay is intended solely for differentiation of positive and negative specimens. Values obtained with the screening assay cannot be used for calculation of EIA Titers. For accurate determination, the semi-quantitative assay should be performed.

- Snap off a sufficient number of microwells for samples and two controls. Insert microwells into the holder and record sample positions. (See PRECAUTION #9).
- Add 50 µL of each specimen, Sample Diluent (negative control), and 1 drop of Positive Control to the bottom of separate microwells. Mix by gently shaking. **Caution should be taken not to cross contaminate samples, as this could cause erroneous results.**
- Incubate at room temperature (20-30 C) for 10 minutes.
- Remove samples from each microwell by aspiration or by inverting the microwells and firmly tapping on an absorbent pad. Fill each microwell with 1X Wash Buffer (approximately 300 µL). Invert the microwells and firmly tap on an absorbent pad. Repeat this wash procedure three times for a total of four washes. (See PRECAUTION #11).
- Add 1 drop of Enzyme Conjugate to each microwell. Mix by gently shaking. Incubate at room temperature (20-30 C) for 10 minutes.
- Remove solution from each microwell by aspiration or by inverting the microwells and firmly tapping on an absorbent pad. Fill each microwell with 1X Wash Buffer (approximately 300 µL). Invert the microwells and firmly tap on an absorbent pad. Repeat this wash procedure three times for a total of four washes. (See PRECAUTION #11).
- Add 2 drops of Substrate Solution I to each microwell. Mix by gently shaking 15-20 seconds.
- Incubate at room temperature (20-30 C) for 10 minutes.
- Add 2 drops of Stop Solution I to each microwell. Mix by gently shaking and wait 2 minutes before reading.
- Observe Reactions:
 - Visual determination – Read within 15 minutes after adding Stop Solution.
 - Spectrophotometric Determination – Zero EIA reader on air. Wipe the underside of microwells with a lint free tissue. Read absorbance at 450 nm or 450/630 nm within 15 minutes after adding Stop Solution.
- Disinfect and retain microwell holder. Discard the used assay materials as biohazardous waste.

INTERPRETATION OF SCREENING ASSAY

- Visual Reading
 - Negative = Colorless
 - Positive = Definite yellow color

Use an EIA reader if there is uncertainty in determining a definite yellow color.
- Spectrophotometric Single Wavelength (450 nm)
 - Negative = OD₄₅₀ < 0.100
 - Indeterminant = OD₄₅₀ ≥ 0.100 and < 0.150
 - Positive = OD₄₅₀ ≥ 0.150
- Spectrophotometric Dual Wavelength (450/630 nm)
 - Negative = OD_{450/630} < 0.070
 - Indeterminant = OD_{450/630} ≥ 0.070 and < 0.100
 - Positive = OD_{450/630} ≥ 0.100

Indeterminant results should be repeated. If the repeated results are still indeterminate, a second specimen should be obtained and tested. Extremely strong reactions may yield a purple precipitate within a few minutes of stopping; these reactions should be reported as positive.

SEMI-QUANTITATIVE ASSAY PROCEDURE

The semi-quantitative assay procedure can be performed on any positive CSF or serum specimen. This procedure is intended to yield values necessary for calculation of an accurate EIA Titer.

If drug therapy is to be monitored using sequential specimens, the preferred method is to repeat each specimen at the same time as the new specimen is titrated. It is not recommended that EIA titers from different assay runs be compared directly to monitor therapy. To prevent multiple freeze-thawing of specimens, several aliquots of each specimen should be frozen.

A. Dilutions

The following dilution series is recommended when performing the semi-quantitative procedure. Combine 100 µL of specimen and 100 µL of Sample Diluent in a test tube, labeled #1, and vortex 1-3 seconds. From this mixture prepare a series of five-fold dilutions as follows:

- Place 200 µL of Sample Diluent in each of 4 test tubes labeled 2-5.
- Using a clean pipet, transfer 50 µL of the diluted patient specimen from tube #1 into tube #2 and mix well.
- Transfer 50 µL from tube #2 into tube #3 and mix well. Continue this dilution procedure through tube #5. Further 1:5 dilutions may be made, if necessary. Final dilutions are: 1:2, 1:10, 1:50, 1:250, 1:1250

B. Assay Procedure

For increased precision, all reagents must be delivered using a pipetter. Use the dropper bottle to deliver the required number of drops (see below) into a clean test tube. Transfer the required volumes of reagents to the microwells using a pipetter.

NOTE: Do not remove dropper tips from reagent bottles; contamination may compromise activity.

Number of Drops Required

	Enzyme Conjugate	Substrate I	Stop Solution I
1 Titer (7 wells)	7	22	12
2 Titers (12 wells)	11	35	19
3 Titers (17 wells)	15	49	27
4 Titers (22 wells)	19	62	30
10 titers (52 wells)	42	143	79

- Snap off a sufficient number of microwells for diluted samples and two controls. Insert into the microwell holder and record sample positions. (See PRECAUTION #9).
- Pipet 50 µL of each diluted sample, Positive Control and Sample Diluent (negative control) into the bottom of separate microwells. **Caution should be taken not to cross contaminate samples, as this could cause erroneous results.**
- Incubate at room temperature (20-30 C) for 10 minutes.
- Remove solution from each microwell by aspiration or by inverting microwells and firmly tapping on an absorbent pad. Fill each microwell with 1X Wash Buffer (approximately 300 µL). Invert microwells and firmly tap on an absorbent pad. Repeat this wash procedure three times for a total of four washes. (See PRECAUTION #11).
- Pipet 50 µL of Enzyme Conjugate into each microwell, including control microwells. Mix by gently shaking. Incubate at room temperature (20-30 C) for 10 minutes.
- Remove solution from each microwell by aspiration or by inverting the microwells and firmly tapping on an absorbent pad. Fill each microwell with 1X Wash Buffer (approximately 300 µL). Invert the microwells and firmly tap on an absorbent pad. Repeat this wash procedure three times for a total of four washes. (See PRECAUTION #11).
- Pipet 100 µL of Substrate Solution I into each microwell. Mix by gently shaking 15-30 seconds.
- Incubate at room temperature (20-30 C) for 10 minutes.

9. Pipet 100 µL of Stop Solution I into each microwell. Mix by gently shaking and wait 2 minutes before reading.
10. Spectrophotometric Determination – Zero EIA reader on air. Wipe the underside of wells with a lint free tissue. Read absorbance at 450 nm or 450/630 nm within 15 minutes after adding Stop Solution.
11. Disinfect and retain microwell holder. Discard the used assay materials as biohazardous waste.

INTERPRETATION OF SEMI-QUANTITATIVE ASSAY

NOTE: Quality Control procedures listed after the Screening Assay should be followed when performing the semi-quantitative assay.

When using a single wavelength reader to obtain results for the semi-quantitative assay, the absorbance of the negative control (Sample Diluent) must be subtracted from all other absorbance values before calculation of the EIA Titer.

Dual wavelength absorbance values need not be adjusted.

A. Selections of Absorbance Values for Calculations

1. Use only absorbance values which are from the semi-quantitative assay to calculate the EIA titer.
 2. Select the highest absorbance value within the acceptable range of 0.100 and 1.500 to calculate the EIA titer. If the next greater dilution also yields an absorbance within the acceptable range, this value should also be selected.
 3. If the 1:2 dilution of the specimen yields an absorbance < 0.100, repeat the assay on the undiluted specimen.
 4. Using the selected absorbance, calculate the EIA titer as outlined in B below. If two absorbances are selected, calculate each titer separately and average the results.
- NOTE:** Appearance of a purple precipitate does not invalidate a positive result, but readings from such wells cannot be used for calculation of EIA Titer.

B. Calculation

1. Calculate the EIA titer for a particular specimen with the following formula:

$$\text{Absorbance Value} \times \text{Multiplication Factor} = \text{EIA titer}$$

	Multiplication Factors*				
Tube/Well Number	1	2	3	4	5
Dilution	1:2	1:10	1:50	1:250	1:1250
Multiplication Factor	20	100	500	2,500	12,500

*Use a Multiplication Factor of 10 to calculate an EIA titer for undiluted specimens.

2. Sample Calculations:

For example, if a patient specimen had an absorbance of 1.2 at a 1:50 dilution, the EIA Titer would be: 1.2 x 500 = 600 which is reported as 1:600

C. Reporting EIA Titers

Numbers derived from the EIA calculations may be reported to physicians as the EIA Titer. A higher titer reflects a greater concentration of antigen. However, do not attempt to compare EIA Titer with latex titer. Due to serotype variations and differing test technologies, EIA Titers are not numerically equivalent to latex titers. Decreases in titers of less than four-fold are not considered clinically significant.⁶

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

At the time of each use, kit components should be visually examined for obvious signs of microbial contamination, freezing, or leakage.

The positive and negative controls must be used with each batch of specimens to provide quality assurance of the reagents.

The negative control (Sample Diluent) absorbance should be > 0.000 and < 0.100 at 450 nm or > 0.000 and < 0.070 at 450/630 nm. The Positive Control absorbance should be > 0.700 at either 450 nm or 450/630 nm.

If absorbance values do not correspond to visual readings, wipe the underside of the microwell, check the alignment of the microwells within the plate holder as well as the alignment of the plate holder, and reread.


It is suggested that the results of each quality control check be recorded in an appropriate log book to facilitate monitoring testing procedures and to comply with regulatory agencies.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

EXPECTED VALUES

Cryptococcal antigen in CSF or serum of the untreated patient indicates an active infection.¹⁻³ The PREMIER Cryptococcal Antigen EIA should have both diagnostic and prognostic value since antigen concentrations rise during disease progression and fall during successful treatment.^{1,2,5,6} Failure of titers to decline may indicate inadequate therapy. However, in some treated patients, antigen levels may remain stationary for extended periods of time when no viable organisms can be demonstrated.⁸

Comparison of EIA and Latex Titers on an Individual Patient

The primary use for semi-quantitative values derived with the latex or the EIA is to follow the course of disease and treatment within an individual patient. These data suggest that the PREMIER Cryptococcal Antigen EIA will be a valuable method for monitoring the efficacy of drug therapies.

AIDS patients may have very high EIA Titers. The dilution series recommended in the procedure may not be adequate to calculate an EIA Titer. Further 1:5 dilutions may be made as needed.

LIMITATIONS OF THE PROCEDURE

The PREMIER Cryptococcal Antigen EIA is intended for use with CSF or serum specimens. Use of this assay is not recommended with other specimens.

A negative result does not preclude diagnosis of cryptococcosis, particularly if only a single specimen has been tested and the patient shows symptoms consistent with cryptococcosis.

A positive result implies the presence of cryptococcal antigen; however, all test results should be reviewed in light of other clinical data by the physician. *Trichosporon beigelli* yielded a false positive in a single case with latex⁹ and a similar cross-reaction is possible with the EIA.

SPECIFIC PERFORMANCE CHARACTERISTICS

The PREMIER Cryptococcal Antigen EIA was evaluated at two major clinical laboratories (279 specimens) in the United States and at Meridian Bioscience, Inc., (88 specimens, including one indeterminate specimen. Insufficient quantity precluded retesting of the indeterminate, and this specimen is not listed below). The CALAS[®] Latex Agglutination System was used as the reference assay. The combined test results are summarized in the table below.

		Latex	
		+	-
PREMIER	+	108	6*
EIA	-	0	252
Sensitivity		100%	
Specificity		98%	
Agreement		98%	

*In four of six discrepancies listed above, previous or subsequent specimens confirmed that the patients had cryptococcosis.

The EIA appeared to be more sensitive by detecting low concentrations of cryptococcal antigen at early or late stages of disease.

Patients with AIDS have emerged as a major group susceptible to cryptococcosis. Cryptococcosis is the fourth most common infection complicating AIDS.² A comparison of clinical diagnosis and PREMIER Cryptococcal Antigen results in an AIDS patient population is listed below.

Clinical Diagnosis Cryptococcosis

		Latex	
		+	-
PREMIER	+	44	2
EIA	-	0	25

ANALYTICAL SENSITIVITY

Limits of detection for serotypes A and B have been demonstrated to be 0.63 ng/mL for EIA and 7.6 ng/mL for the latex assay.¹⁰

REPRODUCIBILITY

Intra-assay Variability – 16 replicates of 3 known positive sera were tested in one assay to determine the intra-assay reproducibility.

Sample	Mean Absorbance	% C.V.
1	0.67	2.1
2	1.74	2.2
3	2.70	3.0

Coefficients of Variation for the same 3 samples using the screening assay with dropper-tips (16 replicates each) were 12, 3.5, and 3.0% respectively.

Inter-assay Variability – 16 replicates of 3 known positive sera were tested in multiple assays by more than one operator.

Sample	Mean Absorbance	% C.V.
4	0.27	11.8
5	0.73	11.5
6	1.76	7.2

Coefficients of Variation for the same 3 samples using the screening assay with dropper-tips (16 replicates each) were 30, 5.8, and 18.9% respectively.

CROSSREACTIVITY

The following list represents organisms that were negative when tested in the EIA:

Bacteria	Fungal Antigen Preparation	Virus	Yeast
<i>Escherichia coli</i> K1	<i>Aspergillus flavis</i>	Coxsackie A-16, B-5	<i>Candida albicans</i>
<i>Haemophilus influenzae</i>	<i>Aspergillus fumigatus</i>	Echovirus II	<i>Rhodotorula rubra</i>
<i>Neisseria meningitidis</i> 301, 362	<i>Aspergillus niger</i>	Poliovirus I, II, III	<i>Torulopsis glabrata</i>
<i>Pseudomonas aeruginosa</i>	<i>Blastomyces dermatitidis</i>		
<i>Staphylococcus aureus</i>	<i>Coccidioides immitis</i>		
<i>Staphylococcus aureus</i> Cowan strain 1	<i>Histoplasma capsulatum</i>		
<i>Staphylococcus epidermidis</i>			
<i>Streptococcus agalactiae</i>			
<i>Streptococcus pneumoniae</i>			

The only organism found to crossreact in the EIA was *T. beigelli*. This organism is not likely to be found in typical specimens being tested for cryptococcal antigen.⁹

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REFERENCES

1. Perfect JR. Cryptococcosis. Infectious Disease Clinics of North America. 1989; 3:77-102.
2. Patterson TF and Andriole VT. Current concepts in cryptococcosis. Eur. J Clin Microbiol Infect Dis 1989; 8:457-465.
3. Rinaldi MG. Cryptococcosis. Lab. Medicine 1982; 13:11-19.
4. Murphy JM. Cryptococcosis, in immunology of the fungal diseases. 1989; Chapter 4, pp. 98-138. CRC Press, Boca Raton, FL.
5. Bennett JE, Hasenclever HF and Tynes BS. Detection of cryptococcal polysaccharide in serum and spinal fluid: value in diagnosis and prognosis. Trans Assoc Am Physicians 1964; 77:145-150.
6. Diamond D and Bennett E. Prognostic factors in cryptococcal meningitis. Ann Int med. 1974; 80:176-181.
7. Goodman JS, Kaufman L and Koenig MG. Diagnosis of cryptococcal meningitis: Value of immunologic detection of cryptococcal antigen. New Eng J of Med 1971; 285:434-436.
8. Kaufman L and Blumer S. Cryptococcosis: The Awakening Giant. Proc. Of the Fourth International Conference on the Mycoses. Pan American Health Organization Scientific Publication #356.1977; pp. 176-182.
9. McManus EJ and Jones JM. Detection of a trichosporon beigelii antigen cross-reactive with cryptococcus neoformans capsular polysaccharide in serum from a patient with disseminated trichosporon infection. J Clin Microbiol 1985; 21:681-685.
10. Data on file (presented at ICAAC, October 1990, reprints available. Manuscript to be submitted to *Journal of Clinical Microbiology*, 1991).

SN11125

REV. 08/17

<p>Manufactured By</p>	<p>Meridian Bioscience, Inc. Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 USA Telephone: 513.271.3700 Orders/Customer Service: 800.543.1980 Technical Support Center: 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124</p>
<p>Authorized Representative</p>	<p>Meridian Bioscience Europe S. r. l. Via dell'Industria, 7 20020 Villa Cortese, Milano ITALY Tel: +39 0331 43 36 36 Fax: +39 0331 43 36 16 Email: info@meridianbioscience.eu WEB: www.meridianbioscience.eu</p>

Meridian Bioscience Europe s.a./n.v.
2 Avenue du Japon - 1420 Braine l'Alleud
BELGIUM
Tel: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
Email: info.bn@meridianbioscience.eu

Meridian Bioscience Europe France
34 rue de Ponthieu - 75008 Paris
FRANCE
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 9 70 06 62 10
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
Postbus 301 - 5460 AH Veghel
NETHERLANDS
Tel: +31 (0) 411 62 11 66
Fax: +31 (0) 411 62 48 41
Email: info.bn@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo medico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenzubereitung, in dem sich Probenverdünnungspuffer befindet
REF	Catalogue number / Numero de catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactif IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur für IVD Leistungsbewertung
	Contains sufficient for n tests / Contenuto sufficiente per n saggi / Contenu suffisant pour n tests / Contenido suficiente para n ensayos / Inhalt ausreichend für n Proben	SOLN STOP	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Límites de temperatura / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Coniugato enzimatico / Conjugat enzymatique / Conjugado enzimático / Enzymkonjugat
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	REAG	Reagent / Reagents / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Miso En Garde / Advertencia / Warnhinweise
CONJ	Conjugate / Coniugato / Conjugat / Conjugado / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
SUBS	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
R. Only	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig	DET REAG	Detection Reagent / Reagente Diretto / Réactif de Detection / Reactivo de Detección / Nachweis Reagenz
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Ne pas utiliser si le paquet est endommagé / No use si el paquete esta dañado / Nicht verwenden, wenn die Verpackung beschädigt ist	TUBE	Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß

For technical assistance, call Technical Support Services at (800) 343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at (800) 543-1980.