PREMIER® COCCIDIOIDES

An EIA for the Qualitative Detection of IgM and IgG Antibodies Directed Against the TP and CF Antigens of *Coccidioides immitis* in Serum and Cerebrospinal Fluid (CSF)

REF 603096

IVD

Rx Only

INTENDED USE

The PREMIER *Coccidioides* enzyme immunoassay (EIA) is for the qualitative detection of IgM and IgG antibodies directed against the TP and CF antigens of *Coccidioides immitis* in serum and cerebrospinal fluid (CSF). It is recommended that EIA positive results be confirmed by immunodiffusion assay.

SUMMARY AND EXPLANATION OF THE TEST

C. immitis is a major causative agent of deep-seated mycoses. Though endemic in the southwestern United States and Mexico, increased travel to the endemic areas has also increased the incidence in nonendemic areas.^{1,2}

Coccidioidomycosis presents a diagnostic challenge to the physician and laboratorian. Lesions produced by *C. immitis* may be difficult to distinguish from other fungal diseases, or from tuberculous lesions and neoplasms (benign and malignant). Symptoms are often unremarkable and can mimic various pneumonias, sarcoidosis, cancer and other maladies.^{1,2} Culturally and histologically the organisms can be difficult to demonstrate, even after repeated attempts.^{1,2}

Frequently serology offers the only evidence available to form a preliminary diagnosis. Serological evidence can then guide the selection of more definitive diagnostic techniques such as culture or biopsy. Complement fixation, immunodiffusion (ID), and latex agglutination are the primary serologic methods.¹ The tube precipitin (TP) antigen, is a 120 kDa glycoprotein with carbohydrate determinants.^{3,4} Antibodies to the TP antigen can be detected with either the tube precipitin assay or an immunodiffusion system.^{5,8} The presence of antibodies to the TP antigen is interpreted as an indication of acute coccidioidal disease and is primarily an IgM response. The complement fixation (CF) antigen is a heat labile protein. Antibodies to the CF antigen can be detected by either ID or complement fixation systems during the later stages of disease.^{1,2,7,9} The complement fixation assay is used to quantify the IgG response and monitor therapy.^{1,2,7} Latex agglutination offers a sensitive and rapid, negative screening capability, but lacks the specificity to be used without confirmation. No single commercial assay system has provided the possibility of screening for both IgM and IgG responses to both TP and CF antigens.

The PREMIER *Coccidioides* EIA is a rapid test for the detection of IgM and IgG to both TP and CF coccidioidal antigens. When used in conjunction with other laboratory and clinical evidence, the PREMIER *Coccidioides* EIA can be a useful tool to aid in the diagnosis of infection due to *C. immitis*.

BIOLOGICAL PRINCIPLES

The PREMIER *Coccidioides* EIA utilizes a mixture of purified TP and CF antigens adsorbed to microwells to capture relevant IgM and IgG antibodies. Diluted specimen is added to each of two duplicate microwells. If antibodies to the coccidioidal antigens are present, they become bound to the adsorbed antigens. After washing to remove unbound specimen components, an anti-IgM enzyme conjugate is added to one microwell and an anti-IgG enzyme conjugate is added to the other. If patient antibodies are bound, then a sandwich is formed between the adsorbed antigens, patient antibodies, and one or both conjugates. After washing to remove unbound conjugate, a substrate solution is added. Color develops in the presence of bound enzyme conjugate.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

- PREMIER *Coccidioides* Antigen Coated Microwells – Breakaway plastic microwells, each coated with a mixture of TP and CF antigens.
- PREMIER *Coccidioides* Positive Control – Prediluted positive human serum containing 0.10% Sodium Azide as a preservative. Do not dilute further.
- PREMIER *Coccidioides* Sample Diluent – Buffered protein solution containing 0.10% Sodium Azide as a preservative.
- PREMIER 20X Wash Buffer II – Concentrated wash buffer containing 0.2% Thimerosal as a preservative.
- PREMIER *Coccidioides* IgM Enzyme Conjugate – Affinity purified goat anti-human IgM antibodies conjugated to horseradish peroxidase in buffered protein solution containing 0.02% Thimerosal as a preservative.
- PREMIER *Coccidioides* IgG Enzyme Conjugate – Affinity purified goat anti-human IgG antibodies conjugated to horseradish peroxidase in buffered protein solution containing 0.02% Thimerosal as a preservative.
- PREMIER Substrate I – Buffered solution containing peroxide and tetramethylbenzidine.
- PREMIER Stop Solution I – 1M Phosphoric Acid. CAUTION: Avoid contact with skin. Flush with water if contact occurs.
- Microwell Strip Holder

MATERIALS NOT PROVIDED

- Pipettes capable of delivering 10 µL, 20 µL, 50 µL, 100 µL and 200 µL.
- Test tubes (12 x 75 mm) for dilution of sample
- Distilled or deionized water
- EIA plate reader capable of reading absorbance at 450 nm. A dual wavelength reader is preferred, using a second filter of 630 nm to correct for light scatter
- Squirt bottle or automatic EIA plate washer
- Timer
- Graduated cylinder for making 1X Wash Buffer

PRECAUTIONS

- All reagents are for in vitro diagnostic use only.
- Reagent concentration, incubation times, and temperatures have been optimized for sensitivity and specificity. Best results are obtained by adhering to these specifications.
- Patient specimen and the Positive Control may contain infectious agents. The Positive Control contains human sera which were screened for HBsAg and antibody to HIV-1 and found to be negative. However, no test can offer complete assurance that human blood will not transmit HIV, hepatitis or other infectious agents. Patient specimens, Positive Control and all materials which they contact should be handled at Biosafety Level 2 as recommended in the CDC/NIH manual "Biosafety in Microbiology and Biomedical Laboratories". All reagents should be gently mixed and at room temperature before use.
- Do not interchange Microwells, IgM Enzyme Conjugate, IgG Enzyme Conjugate or Positive Control Reagent between different kits. 20X Wash Buffer II, Sample Diluent, Substrate I and Stop Solution I can be interchanged provided they are within their expiration dates at the time of testing.
- Hold reagent vials vertically to insure proper drop size and delivery.
- Replace colored caps on correct vials.
- Avoid skin contact with Stop Solution I (1M Phosphoric acid). Flush with water immediately if contact occurs.
- Avoid splashing when dispensing diluted sample into microwell by placing pipette tip about halfway down the well and dispensing slowly down side of well.
- It is important that microwell washing be performed precisely as directed in assay procedure. **Inadequate washing may be the cause of elevated background in any EIA protocol.**
- Do not use microwells if the foil pouch appears damaged (ie, has holes, punctures).
- Do not use microwells where the desiccant indicator has changed from blue to pink.

HAZARD and PRECAUTIONARY STATEMENTS

Refer to the SDS, available at www.meridianbioscience.com for Hazards and Precautionary Statements.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-8 C and return the kit promptly to the refrigerator after each use. **Unused microwells must be removed from the microwell holder and placed back inside the resealable foil pouch and sealed immediately after opening pouch. It is important to protect strips from moisture.**

REAGENT PREPARATION

Bring the entire kit, including the microwell pouch, to 22-25 C before and during use. **Warming requires at least one hour.** Prepare sufficient 1X Wash Buffer for use by measuring one part 20X Wash Buffer II and diluting with 19 parts purified water. The 1X Wash Buffer can be stored at 22-25 C for up to one month. Discard if the 1X Wash Buffer becomes contaminated. For performing washes, the 1X Wash Buffer can be transferred to a wash or "squirt" bottle. Alternatively, an automated plate washer can be used. NOTE: Automated plate washers tend to waste buffer and may deplete the 20X Wash Buffer II prematurely. Call Technical Services (1-800-543-1980) to order additional PREMIER 20X Wash Buffer II if needed.

SPECIMEN COLLECTION AND PREPARATION

Serum from clotted blood or CSF may be tested with this kit. Heat-inactivated specimens (56 C for 30 minutes) are also acceptable. Do not use bloody CSF. Specimens should be tested as soon as possible, but may be stored up to five days at 2-8 C prior to testing. If longer storage is required, it is recommended that several aliquots of each specimen be frozen (-20 to -80 C) to avoid multiple freeze-thaw cycles. Do not store in a frost-free freezer.

Dilute serum and CSF with Sample Diluent as described below. Performing the dilutions by other schemes may exhaust the supply of Sample Diluent prematurely.

Serum

Dilute serum 1:441 with Sample Diluent as follows:

- Obtain two test tubes for each serum specimen. Transfer 200 µL of Sample Diluent into the first tube and 400 µL into the second tube.
- Transfer 10 µL of serum to the first tube and mix thoroughly.
- Transfer 20 µL of the first dilution into the second tube and mix thoroughly.

CSF

Dilute CSF 1:21 with Sample Diluent as follows:

- Obtain a test tube for each CSF specimen and transfer 400 µL of Sample Diluent into each tube.
- Transfer 20 µL of CSF into the dilution tube and mix thoroughly.

TEST PROCEDURE

NOTE: Each microwell is coated with both TP and CF antigens. It is recommended that both IgM and IgG assays be performed simultaneously for all sera.

- Snap off a sufficient number of microwells for samples and controls and insert them into the microwell holder. Record sample positions, as in the example below. **Avoid testing more specimens than can be loaded in 6 minutes in step 3.**

	1	2	3	4	5	6	
A	PC	SD	Pt. #1	Pt. #2	Pt. #3	Pt. #4	IgM
B	PC	SD	Pt. #1	Pt. #2	Pt. #3	Pt. #4	IgG

- Add 100 µL of Positive Control (PC) to each of 2 separate microwells (A1 and B1) and 100 µL of Sample Diluent (SD) to each of 2 separate microwells (A2 and B2). **The latter wells are considered Negative Controls.**
- Add 100 µL of each specimen (diluted as described in the previous section) to 2 separate microwells (in rows A and B, as needed). Caution should be taken not to cross contaminate samples, as this could cause erroneous results. Mix by gently shaking 10-15 seconds on the countertop.
- Incubate at 22-25 C for 30 minutes after the addition of all specimens and controls.
- Hold the plate firmly from the bottom and gently squeeze it.
 - Dump plate contents firmly into a biohazard receptacle and strike the inverted plate firmly on a clean stack of paper towels or other absorbent material.
 - Fill all wells with 1X Wash Buffer (use wash bottle or other device) directing stream of buffer to the sides of the wells to avoid foaming. Dump the plate contents and strike plate on towels.
 - Repeat wash cycle (step 5b) 2 times for a total of 3 washes. After the final wash, strike plate on fresh towels hard enough to remove as much excess wash buffer as possible but do not allow wells to completely dry at any time.
- Add 2 drops of IgM Enzyme Conjugate to the microwells of row A, or as required according to the number of specimens.
- Add 2 drops of IgG Enzyme Conjugate to the microwells of row B, or as required according to the number of specimens. Mix by gently shaking 10-15 seconds on the countertop.
- Incubate at 22-25 C for 30 minutes after the Enzyme Conjugates have been added to all microwells.
- Remove Enzyme Conjugates from each microwell and wash 3 times as in step 5.
- Add 2 drops of Substrate I to each microwell. Start timer with addition of Substrate to first well. Mix by gently shaking 10-15 seconds on the counter top.
- Incubate at 22-25 C for 5 minutes.
- Add 2 drops of Stop Solution I to each microwell in the same order as step 10. Mix by gently shaking 10-15 seconds on the counter top and wait 2 minutes before reading. Readings should be made within 15 minutes.
- Carefully, wipe the underside of the microwells with a lint free tissue and measure the absorbance.
 - A dual wavelength reader is preferred, with absorbance read at 450 nm and 630 nm. BLANK ON AIR.
 - If a single wavelength reader is used, read the absorbance at 450 nm as follows:
 - Blank on air.
 - Read the Negative Control microwell(s), value(s) should be < 0.100.
 - Then reblank the reader on the Negative Control microwell(s) to compensate background.
- Disinfect and retain microwell holder. Discard used assay materials as biohazardous waste.

INTERPRETATION OF RESULTS

Results

Negative = Absorbance Value < 0.150
Indeterminate = Absorbance Value ≥ 0.150 but ≤ 0.199
Positive = Absorbance Value ≥ 0.200

A negative result with both IgM and IgG indicates that antibody to *C. immitis* antigens is either absent, below the level of detection of the assay or the specimen was obtained too early in the response. A positive result with either IgM or IgG implies the presence of antibody to *C. immitis* antigens and a positive response with either conjugate should be reported. An early acute phase patient may only present an IgM response, while the chronic or convalescent patient may only present an IgG response.

Retest specimens which yield an indeterminate result and report the result of the retest. If the retest result is still indeterminate, a second specimen should be obtained. Extremely strong positive reactions may yield a purple precipitate within a few minutes of stopping the reaction. Absorbances obtained with such reactions may be lower than expected but will still be positive.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

- The Positive and Negative Controls must be assayed with each batch of patient specimens to provide quality assurance of the reagents.
- The Negative Control (Sample Diluent) microwells should yield an A₄₅₀ value < 0.100 with both IgM and IgG Conjugates when blanked on air or an A_{450,630} value < 0.050. Negative values ≤ -0.005 may indicate a problem with EIA plate reader.
- The Positive Control microwells should have a definite yellow color and yield an A₄₅₀ or A_{450,630} value ≥ 0.500 and ≤ 2.500 with both IgM and IgG Conjugates.
- At the time of each use, kit components should be visually examined for obvious signs of microbial contamination, freezing, or leakage. A flaky precipitate may form in the Sample Diluent which will not affect its performance.
- It is suggested that the results of each quality control check be recorded in an appropriate log book to maintain high quality testing procedures and compliance with regulatory agencies.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

EXPECTED VALUES

Expected values for a given population should be determined for each laboratory. The rate of positivity may vary depending on geographic location, method of specimen collection, handling and transportation, test employed and general health environment of patient population under study.

The PREMIER *Coccidioides* EIA test has shown positive rates of 42.5%, 4.0%, and 2.5% from three prospective clinical studies. The frequency of coccidioid infection was higher from the endemic region than the nonendemic region.

LIMITATIONS OF THE PROCEDURE

- The PREMIER *Coccidioides* EIA is intended for use with CSF or serum specimens. Use of this assay is not recommended with other specimens.
- A negative result with both IgM and IgG does not preclude diagnosis of coccidioidomycosis, particularly if only a single specimen has been tested and the patient shows symptoms consistent with a positive diagnosis.
- A positive result with either IgM or IgG implies the presence of antibody to *C. immitis*. Positive EIA results should be confirmed by immunodiffusion assay. However, because of the relative insensitivity of the ID procedures, and ID negative does not preclude the possibility of coccidioidomycosis.¹¹ Diagnosis of coccidioidomycosis is based on laboratory and clinical findings as well as the presence of antibody.
- Positive results with either IgM or IgG (but not both) also suggest coccidioid infection, but in different disease states. An early acute phase patient may only present an IgM response, while the chronic or convalescent patient may only present an IgG response. Such EIA results should be compared with patient symptoms to determine whether there is a logical correlation. Because the IgM and IgG results do not confirm each other in these cases, additional confirmation may be prudent, particularly if symptoms do not correlate with the single EIA result. Positive results with both IgM and IgG strongly suggest active disease with a high degree of confidence since the assays confirm each other.

SPECIFIC PERFORMANCE CHARACTERISTICS

The PREMIER *Coccidioides* EIA was evaluated at two Medical Centers in an endemic area, one central reference laboratory and Meridian Bioscience, Inc. Of the 1,315 samples submitted for *Coccidioides* testing, 1,242 were sera and 73 were CSF. A comparison of the results obtained with the EIA and complement fixation methods are presented in the tables below. When the EIA and complement fixation results did not agree, immunodiffusion (ID) and/or tube precipitin (TP) results were used to resolve the discrepancies. The EIA results for IgM and IgG are presented separately and combined. For the combined comparison, a specimen was considered positive if either IgM, IgG or both results were positive. The results indicate that the greatest sensitivity is obtained when both IgM and IgG EIA assays are performed. This observation is consistent with the notion that patients in the early acute phase may present with only an IgM response while the chronic or convalescent patient may present with only an IgG response.

PREMIER EIA IgM Only	Complement Fixation		Resolved	
	Pos	Neg	Pos	Neg
Pos	160	67	189	38
Neg	75	971	65	981
Ind	10	32	12	30

n=1,315

	vs. Comp. Fix.	vs. Resolved
Relative Sensitivity	68.1%	74.4%
Relative Specificity	93.5%	96.3%
Overall Agreement	86.0%	89.0%

When the discrepancies between the EIA IgM results and complement fixation were resolved, 29 of the 67 specimens which were positive by EIA but negative by complement fixation were confirmed to be positive by ID assays. Ten of the 75 specimens negative by EIA and positive by complement fixation were ID negative.

PREMIER EIA IgG Only	Complement Fixation		Resolved	
	Pos	Neg	Pos	Neg
Pos	214	69	248	35
Neg	24	987	22	989
Ind	7	14	8	13

n=1,315

	vs. Comp. Fix.	vs. Resolved
Relative Sensitivity	89.9%	91.9%
Relative Specificity	93.5%	96.6%
Overall Agreement	91.3%	94.1%

When the discrepancies between the EIA IgG results and complement fixation were resolved, 34 of the 69 specimens which were positive by EIA but negative by complement fixation were confirmed to be positive by ID assays. Two of the 24 specimens negative by EIA and positive by complement fixation were ID negative.

PREMIER EIA IgM and IgG	Complement Fixation		Resolved	
	Pos	Neg	Pos	Neg
Pos	228	99	267	60
Neg	14	934	7	941
Ind	3	37	3	37

n=1,315

	vs. Comp. Fix.	vs. Resolved
Relative Sensitivity	94.2%	97.4%
Relative Specificity	90.4%	94.0%
Overall Agreement	91.1%	94.7%

When the discrepancies between the EIA IgG and IgM results and complement fixation were resolved, 39 of the 99 specimens which were positive by EIA but negative by complement fixation were confirmed to be positive by ID assays. Seven of the 14 specimens negative by EIA and positive by complement fixation were ID negative. At one site in the endemic area, three patients were determined to be EIA positive, while each was negative by all of the classical serologic methods. Retesting of a subsequent specimen from 2-14 days later, confirmed that these patients had acute disease and the EIA had correctly demonstrated antibody at an earlier stage than the classical methods.

REPRODUCIBILITY

Intra-assay Variability – Twelve replicates of three different known positive sera were tested for both IgM and IgG in one assay to determine the intra-assay reproducibility. All 12 replicates of each specimen were performed using a single 1:441 dilution.

Sample	Intra-assay Variability		
	Mean IgM Abs.	S.D.	% C.V.
1	0.478	0.026	5.4
2	0.477	0.033	6.9
3	0.928	0.045	4.8
Sample	Intra-assay Variability		
	Mean IgG Abs.	S.D.	% C.V.
1	2.179	0.072	3.3
4	1.132	0.069	6.1
5	0.511	0.048	9.4

Inter-assay Variability – Three known positive specimens were tested for IgM and IgG in multiple assay runs by two or more operators. Each operator performed the 1:441 dilutions and then performed the assay protocol. Each run consisted of triplicates of each specimen run by a single operator.

Sample	Inter-assay Variability		
	Mean IgM Abs	S.D.	% C.V.
6	0.532	0.033	6.1
7	0.562	0.064	11.4
8	1.020	0.071	7.0
Sample	Inter-assay Variability		
	Mean IgG Abs.	S.D.	% C.V.
6	2.082	0.170	8.2
9	1.147	0.080	7.0
10	0.520	0.069	13.4

CROSSREACTIVITY

Numerous reports have noted a significant cross-reaction problem with complement fixation, latex, and other EIA assays for coccidioid antibodies,¹⁻⁴ particularly with patients suffering from other fungal diseases such as Histoplasmosis, Blastomycosis, and Aspergillosis. The data below demonstrate that when specimens serologically positive for antibodies to various organisms as well as specimens containing rheumatoid factor, antinuclear antibodies of Cryptococcal antigen were tested with PREMIER *Coccidioides*, the overall specificity was 95.7%.

	Mean IgM	Mean IgG	Pos./Total	Specificity
<i>Histoplasma</i>	0.043	0.041	4/79	95%
<i>Blastomyces</i>	0.026	0.025	0/9	100%
<i>Aspergillus</i>	0.062	0.039	2/20	90%
Rheumatoid Factor	0.060	0.015	0/10	100%
Antinuclear Abs	0.026	0.010	0/7	100%
<i>Mycoplasma</i>	0.026	0.008	0/11	100%
<i>Cryptococcus</i>	0.007	0.008	0/8	100%
Influenza	0.016	0.009	0/5	100%
HIV	0.008	0.031	1/15	93%

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REFERENCES

- Pappagianis D and BL Zimmer. Serology of coccidioidomycosis. Clin Microbiol Reviews 1990;3:247-268.
- Bronnmann DA, JN Galgiani. Coccidioidomycosis Eur J Clin Microbiol Infect Dis 1989;8:466-473.
- Cole GT, D Kruse, S Zhu, KR Seshan, and RW Wheat. Composition, serologic reactivity, and immunolocalization of a 120-kilodalton tube precipitin antigen of *Coccidioides immitis*. Infect Immun 1990;58:179-188.
- Cole GT, D Kruse, and KR Seshan. Antigen complex of *Coccidioides immitis* which elicits a precipitin antibody response in patients. Infect Immun 1991;59:2434-2446.
- Huppert M and JW Bailey. The use of immunodiffusion tests in coccidioidomycosis. Amer J Clin Path 1965;44:364-368.
- Kaufman L. Value of immunodiffusion tests in the diagnosis of systemic mycotic diseases. Ann Clin Lab Sci 1973;3:141-146.
- Smith CE, MT Saito, and SA Simons. Pattern of 39,500 serologic tests in coccidioidomycosis. J Amer Med Assoc 1956;160:546-552.
- Cox, RA and LA Britt. Isolation of a coccidioidin component that reacts with immunoglobulin M precipitin antibody. Infect Immun 1986;53:449-453.
- Zimmer BL and D Pappagianis. Characterization of a soluble protein of *Coccidioides immitis* with activity as an immunodiffusion-complement fixation antigen. J Clin Microbiol 1988;26:2250-2256.
- Data on file (submitted to ASM for May 1992 meeting. Manuscript to be submitted to J Clin Microbiol, 1992).
- Kaufman L and MJ Clark. Value of the concomitant use of complement fixation and immunodiffusion tests in the diagnosis of coccidioidomycosis. Applied Microbiol 1974;28:641-643.

SN11130

REV. 08/17

 Manufactured By	Meridian Bioscience, Inc. Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 USA Telephone: 513.271.3700 Orders/Customer Service: 800.543.1980 Technical Support Center: 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124
 Authorized Representative	Meridian Bioscience Europe S. r. L. Via dell'Industria, 7 20020 Villa Cortese, Milano ITALY Tel: +39 0331 43 36 36 Fax: +39 0331 43 36 16 Email: info@meridianbioscience.eu WEB: www.meridianbioscience.eu

Meridian Bioscience Europe s.a./n.v.
 2 Avenue du Japon - 1420 Braine l'Alleux
 BELGIUM
 Tel: +32 (0) 67 89 59 59
 Fax: +32 (0) 67 89 59 58
 Email: info.bn@meridianbioscience.eu

Meridian Bioscience Europe France
 34 rue de Ponthieu - 75008 Paris
 FRANCE
 Tel: +33 (0) 1 42 56 04 40
 Fax: +33 (0) 9 70 06 62 10
 Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
 Postbus 301 - 5460 AH Veghel
 NETHERLANDS
 Tel: +31 (0) 411 62 11 66
 Fax: +31 (0) 411 62 48 41
 Email: info.bn@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo medico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostika 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenzubereitung, in dem sich Probenverdünnungspuffer befindet
REF	Catalogue number / Numero de catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consultar las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Reagents IVD reserves à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur für IVD Leistungsbewertung
	Contains sufficient for n tests / Contiene sufficiente per n test / Contenu suffisant pour n tests / Contenido suficiente para n ensayos / Inhalt ausreichend für n Prüfungen	SOLN STOP	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Límites de temperatura / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Coniugato enzimatico / Conjugat enzymatique / Conjugado enzimático / Enzymkonjugat
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	REAG	Reagent / Reagents / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Miso En Garde / Advertencia / Warnhinweise
CONJ	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
SUBS	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
R. Only	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig	DET REAG	Detection Reagent / Reagente Diretto / Réactif de Détection / Reactivo de Detección / Nachweis Reagenz
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Ne pas utiliser si le paquet est endommagé / No use si el paquete esta dañado / Nicht verwenden, wenn die Verpackung beschädigt ist	TUBE	Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.