EIA for the Detection and Verification of IgG Antibody to
Helicobacter pylori in Human Serum and Plasma

REF 606096

IVD

Rx Only

INTENDED USE

Premier *H. pylori* is an enzyme immunoassay (EIA) for the *in vitro* qualitative detection of IgG antibodies to *Helicobacter pylori* in human serum and plasma. Test results are intended to aid in the diagnosis of *H. pylori* infection.

SUMMARY AND EXPLANATION OF THE TEST

The importance of *Helicobacter pylori* in gastro-intestinal diseases has increased greatly since Marshall and Warren described the presence of campylobacter-like organisms in the antral mucosa of patients with histological evidence of antrum gastritis and peptic ulcers especially duodenal ulcers.^{3, 15} Though a causative relationship has not yet been fully established, the strong correlation between the presence of *H. pylori* and histologically confirmed gastritis, peptic ulcer disease and gastric carcinoma can no longer be debated.^{3, 11, 20}

The ecological niche in humans appears to be restricted to the stomach and duodenum. Patients who harbor the organism are divided into two basic groups. The first are those who are said to be "colonized". These patients have the organism, yet have no signs of gastrointestinal disease. Those with symptoms and *H. pylori* are considered to be infected. The process by which a colonized individual becomes infected remains unclear. The process by which patients become colonized is also still under investigation.^{3, 7, 11, 13, 14}

The diagnostic strategies for the determination of *H. pylori* have been developed along two lines: 1) direct detection of the organism, and 2) detection of antibodies made against *H. pylori*.^{8, 9, 22}

Direct detection by invasive methods requires that a biopsy be taken from the upper gastrointestinal tract. The presence of *H. pylori* is then confirmed by direct microscopic examination, rapid urease testing or culturing of the organism from the biopsy material. This strategy has the advantage of being able to detect active infections while being highly specific with a very high positive predictive value. The difficulties associated with this approach are risk and discomfort to the patient. In addition, *H. pylori* tends to colonize in patches and may be missed totally by the biopsy. Culturing of *H. pylori* from biopsy material tends to be difficult and time consuming. These technical difficulties may lead to false negative results.^{1, 2, 3, 19}

The urea breath test (using ¹⁴C or ¹³C -urea) is a noninvasive method that detects *H. pylori* by exploiting its highly active urease.^{3, 8, 10} Though highly sensitive and specific, the test has limitations. It is expensive, time consuming, and requires ingestion of isotopically labeled urea as well as specialized instrumentation for the detection of ¹⁴C or ¹³C. The most common noninvasive approach to the detection of *H. pylori* is the serological identification of specific antibodies in infected patients. Though an indirect approach, the correlation between histological gastritis, the presence of *H. pylori* and seropositivity is strong. The disadvantages of these tests are that many have a lower specificity due to cross reactions from other organisms, or antibody may be present as a "serological" scar after eradication of the organism. Their major advantage is that they are rapid and offer a high sensitivity.^{5, 18, 21}

Premier *H. pylori* is a micowell-based enzyme immunoassay that detects IgG antibody by using a sonicated bacterial cell lysate as the capture antigen. No calculations are required and the visual color change makes the interpretation of results objective and simple.

BIOLOGICAL PRINCIPLES

Premier *H. pylori* utilizes a sonicated *H. pylori* bacterial cell lysate adsorbed to plastic microwells. Diluted patient serum is incubated in the wells at 19-27 C. After washing, a peroxidase conjugated monoclonal antibody to human IgG is added and the plate reincubated. The unbound conjugate is washed away and Premier Substrate I added. Color develops in the presence of bound enzyme. Premier Stop Solution I is added and the results are interpreted visually or spectrophotometrically.

REAGENTS/MATERIALS PROVIDED

The maximum number of test obtained from this test kit is listed on the outer box.

1. Antigen Coated Microwells - Breakaway plastic wells coated with *H. pylori* antigen.
2. Positive Control - Human serum containing IgG antibodies to *H. pylori* in a buffer (pH 7.3) with 0.09% sodium azide.
3. Negative Control - Nonreactive human serum in a buffer (pH 7.3) with 0.09% sodium azide.
4. Sample Diluent - Protein solution with a preservative (0.1% sodium azide).
5. Premier 20X Wash Buffer I - Concentrated wash buffer with a preservative (0.2% thimerosal).
6. Enzyme Conjugate - Monoclonal anti-human IgG labeled with peroxidase in a buffer with preservative (0.02% thimerosal).
7. Premier Substrate I - Buffered solution containing peroxide and tetramethylbenzidine.
8. Premier Stop Solution I - 1M Phosphoric Acid. **CAUTION:** Avoid skin contact. Flush with water if contact occurs.
9. Micowell strip holder
10. Micowell strip sealer

MATERIALS NOT PROVIDED

1. Test tubes (12 x 75 mm)
2. Distilled or deionized water
3. Precision pipettes to deliver 10, 100, and 500 µL
4. Squirt bottle or automatic EIA plate washer
5. Timer
6. Graduated cylinder
7. **OPTIONAL:** EIA plate reader capable of reading absorbance of 450 nm or 450/630 nm

PRECAUTIONS

1. All reagents are for *in vitro* diagnostic use only.
2. Patient specimens may contain infectious agents and should be handled and disposed of as potential biohazards. Wear disposable gloves while handling specimens and performing the test procedure.
3. **CAUTION:** Source material from which the controls were derived was negative for HBsAg/anti-HCV and HIV nonreactive. However, no test method can offer complete assurance that human blood will not transmit HIV, hepatitis or other potentially infectious agents. These reagents should be handled at Biosafety Level 2 as recommended in the CDC/NIH manual Biosafety in Microbiology and Biomedical Laboratories, 2007.
4. All reagents should be mixed gently before use.
5. Reproducible results depend on careful pipetting, correct incubation periods, and temperature, as well as washing the test strips and thorough mixing of all prepared solutions.
6. All reagents except the Premier 20X Wash Buffer I and Controls are provided already diluted to the proper concentration.
7. Do not scratch coated well during washing and aspiration. Wash and fill all microwells without interruption. While washing, check that microwells are filled evenly with washing solution, and that there are no residues in the microwells. Don't allow wells to dry out. **Inadequate washing may cause elevated background.**
8. Instructions for using appropriate photometers are to be observed; check adjustment of proper wave length (450 nm) and reference wave length (630 nm) respectively.
9. If comparisons with other methods are required, always perform both tests simultaneously.
10. Premier Stop Solution I should be handled carefully because it can cause burns or irritation to the skin.
11. All reagents must be stored at 2-8 C and warmed to 19-27 C prior to running the assay.
12. Replace colored caps on correct vials.

13. Use only plastic disposable tips to pipet reagents.
14. Store the unused microwell strips in the original pouch at 2-8 C. It is important to protect strips from moisture.
15. Do not mix reagents from other kit lots, or use kit components beyond labeled expiration date.
16. Use separate pipette tips for each sample, control and reagent. Cross contamination of samples could cause false results.
17. Do not reuse microwells.
18. Dispose of all used test materials in appropriate container. Treat as potential biohazard.
19. Hold reagent vials vertically to insure proper drop size and delivery.
20. Do not use microwells where the desiccant indicator has changed from blue to pink.
21. Do not use microwells with pouches that have been damaged (ie, show holes or punctures).

HAZARD and PRECAUTIONARY STATEMENTS

 Sample Diluent	Signal Word Danger Hazard Statements H300 – Fatal if swallowed H310 – Fatal in contact with skin H410 – Very toxic to aquatic life with long lasting effects Precautionary Statements - EU (§28, 1272/2008) P301 + P310 – IF SWALLOWED: Immediately call a POISON CENTER or doctor/ physician P280 – Wear eye protection/ face protection P361 – Remove/Take off immediately all contaminated clothing
 20X Wash Buffer I	Signal Word Warning Hazard Statements H312 - Harmful in contact with skin Precautionary Statements - EU (§28, 1272/2008) P280 - Wear protective gloves/protective clothing/eye protection/face protection. P321 - Specific treatment (see supplemental instructions on the administration of antidotes on this label).
 Stop Solution I	Signal Word Danger Hazard Statements H412 - Harmful to aquatic life with long lasting effects H314 - Causes severe skin burns and eye damage Precautionary Statements - EU (§28, 1272/2008) P260 - Do not breathe dust/fume/gas/mist/vapors/spray P264 - Wash face, hands and any exposed skin thoroughly after handling P280 - Wear protective gloves/ protective clothing/ eye protection/ face protection P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing P310 - Immediately call a POISON CENTER or doctor/ physician P303 + P361 + P353 - IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/shower P363 - Wash contaminated clothing before reuse P304 + P340 - IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing P310 - Immediately call a POISON CENTER or doctor/ physician P301 + P330 + P331 - IF SWALLOWED: Rinse mouth. Do NOT induce vomiting P405 - Store locked up P501 - Dispose of contents/ container to an approved waste disposal plant.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-8 C and return the kit promptly to the refrigerator after each use.

SPECIMEN COLLECTION AND PREPARATION

1. Serum or plasma can be used in this assay. Sodium citrate and EDTA are acceptable anticoagulants. The performance of **Premier H. pylori** from plasma obtained using other anticoagulants is unknown. Specimens must be collected using established laboratory procedures.
2. Serum and plasma can be stored in the refrigerator at 2-8 C for up to five days. If testing is to be performed at a later date, the specimen should be frozen. Avoid repeated freezing and thawing. A cloudy serum/plasma specimen or one containing particulate matter, may be centrifuged prior to testing.

REAGENT PREPARATION

1. Bring entire kit, including microwell pouch, to 19-27 C before use.
2. Prepare 1X Wash Buffer as needed. For example; 4.0 mL of Premier 20X Wash Buffer I + 76.0 mL of distilled or deionized water is sufficient to wash one strip. Place in a clean squirt bottle. The 1X Wash Buffer can be stored at room temperature for up to three months.
3. The Controls (Negative and Positive) must be diluted in Sample Diluent prior to use in assay.

TEST PROCEDURE

1. Use 1 microwell for each sample and control to be tested:
 - a. Negative Control microwell.
 - b. Positive Control microwell.
 - c. Patient serum microwells.
2. Record location of each microwell in the holder.
3. Prepare 1:50 dilution of each patient specimen, Positive Control and Negative Control in Sample Diluent. Prepare dilution by mixing 10 μ L of specimen, Positive or Negative Control with 0.5 mL Sample Diluent.
4. Pipet 100 μ L of each diluted Control (Negative and Positive) to the appropriate microwell. Pipet 100 μ L of diluted patient specimen into corresponding microwells. (Load time should not exceed 5 minutes).
5. Incubate at 19-27 °C for 20 minutes. Protect plate from drafts and dust which may cause evaporation and contamination.
6. Remove solution from each microwell by inverting the microwells and tapping on an absorbent pad, firmly compressing the holder to prevent wells from falling out. Fill each microwell with 1X Wash Buffer (approximately 350 μ L). Invert the microwells and firmly tap on an absorbent pad. Repeat this procedure 2 times for a total of 3 washes.
7. Add 2 drops of Conjugate to all wells.
8. Incubate at 19-27 °C for 20 minutes.
9. Wash the microwells 3 times with 350 μ L 1X Wash Buffer per well as in Step #6.
10. Add 2 drops of Premier Substrate Solution I to all wells. Agitate plate gently to mix.
11. Incubate at 19-27 °C for 10 minutes.
12. Add 2 drops of Premier Stop Solution I to each well. Agitate plate gently for 30 seconds to mix. **Note:** Initial color of positive reaction is blue, which changes to yellow upon addition of Premier Stop Solution I.
13. Observe Reactions:
 - a. Visual Determination - Read within 10 minutes after adding Premier Stop Solution I.
 - b. Spectrophotometric Determination - Zero EIA reader on air. Wipe underside of wells with a lint free tissue. Read absorbance at 450 nm or 450/630 nm within 10 minutes of adding Premier Stop Solution I.

INTERPRETATION OF RESULTS

1. **Visual Reading**
Negative = colorless
Positive = definite yellow color
2. **Spectrophotometric Single Wavelength (450 nm)**
Negative = $OD_{450} < 0.120$
Positive = $OD_{450} \geq 0.120$
3. **Spectrophotometric Dual Wavelength (450/630 nm)**
Negative = $OD_{450/630} < 0.070$
Positive = $OD_{450/630} \geq 0.070$

A positive test result indicates the presence of IgG antibody to *H. pylori*. A negative test result indicates the absence of IgG antibody to *H. pylori* or that the IgG level is below that which can be detected by this test.

Extremely strong positive reactions may yield a purple precipitate within a few minutes of stopping the reaction.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

The Positive and Negative Controls must be used with each batch of specimens to provide quality assurance of the reagents.

The Negative Control should read < 0.120 at 450 nm and < 0.070 at 450/630 nm but greater than 0.00. If the Negative Control is < 0.00, re-blank the plate reader to air and re-read the plate. The Negative Control should be colorless when read visually. The Positive Control must read ≥ 0.400 at both 450 nm and 450/630 nm.

Any positive well without visible color should be repositioned, wiped on the underside of the well and reread.

At the time of each use, kit components should be visually examined for obvious signs of microbial contamination, freezing, or leakage.

Repeat the control tests as the first step in determining the root cause of the failure. **If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.**

EXPECTED VALUES

Studies on the epidemiology of *H. pylori* have shown that this organism is present worldwide.^{12, 16} Gastritis caused by *H. pylori* has been shown to correlate with age, ethnic background, family size and socioeconomic class.^{6, 10} While the epidemiological studies on *H. pylori* infection or colonization are ever expanding, the route by which individuals acquire this organism is, at present, unknown.^{3, 7, 11, 13, 14}

Studies indicate that the incidence of infection in the United States may increase 1-2% annually.¹⁰ It is not uncommon to see positivity rates of 50% in those who are 60 years or older. The frequency of *H. pylori* infection in patients diagnosed with duodenal ulcers however, has been shown for patients in every age group to be approximately 80%.¹² The positivity rates seen with *H. pylori* clearly depend on the sample population and how that population is defined.

The **Premier *H. pylori*** test detects the presence of IgG antibodies to *H. pylori* in serum and plasma. Expected values for a given population should be determined for each laboratory. The rate of positivity may vary depending on geographic location, method of specimen collection, handling and transportation, test employed and general health environment of patient population under study.

LIMITATIONS OF THE PROCEDURE

1. The test should be used only for individuals with symptoms suggestive of gastrointestinal disease. **Premier *H. pylori*** is not intended for use with asymptomatic patients.
2. The test is qualitative and no quantitative interpretation should be made with respect to the values.
3. Samples obtained too early during infection may not contain detectable antibodies. If *H. pylori* infection is suspected and a negative result is found, a second sample should be obtained 2 to 7 weeks later and tested in parallel with the first sample.
4. A positive test result (i.e., presence of antibody to *H. pylori*) in symptomatic individuals can only infer active infection and, whenever possible, should be confirmed by bacterial isolation or other diagnostic testing. The test cannot distinguish between active and inactive infection.
5. Test results should be used in conjunction with information available from the patient clinical evaluation and other diagnostic procedures.
6. **Premier *H. pylori*** is a qualitative test that only indicates the presence of IgG antibody to *H. pylori*. This test does not indicate that gastrointestinal disease is present.
7. A negative test result indicates the absence of IgG antibody to *H. pylori* or that the IgG level is below that which can be detected by this test.

SPECIFIC PERFORMANCE CHARACTERISTICS

The Premier *H. pylori* test was evaluated at three hospitals. A total of 462 serum specimens were obtained from patients with gastrointestinal complaints. All patients underwent endoscopy. The Premier *H. pylori* test was compared to results obtained from the biopsied material using rapid biopsy urease tests. Discrepant results were resolved by using stained tissue sections and/or clinical diagnosis.⁴

Premier <i>H. pylori</i>	Biopsy Urease Test		Resolved by Tissue and Diagnosis	
	Positive	Negative	Positive	Negative
	Positive	Negative	Positive	Negative
Premier <i>H. pylori</i>	166	102	259	9
Positive	3	191	3	191

n = 462

Relative Sensitivity = 98.2% (99.4-94.9)*(166/169) 98.8% (99.6-96.7) (259/262)
 Relative Specificity = 65.2% (70.4-59.6) (191/293) 95.5% (97.6-91.7) (191/200)

Relative Agreement = 77.3% (80.8-73.2) (357/462) 97.4% (98.5-95.5) (450/462)

*95% confidence intervals calculated by the normal method are indicated in parenthesis.

Ninety-three biopsy urease test negative/Premier *H. pylori* positive specimens were positive by either tissue or clinical diagnosis. Nine biopsy urease test negative/Premier *H. pylori* positive specimens had no indication of *H. pylori* by either resolving method. The three biopsy urease positive/Premier *H. pylori* negative samples were not considered for resolution. Premier *H. pylori* was compared to the HM-CAP micowell based EIA also for the detection of IgG antibody to *H. pylori*. Discrepant data was resolved by the biopsy urease test and/or the determination for the presence of *H. pylori* by histological staining methods.

Premier <i>H. pylori</i>	Competitor EIA			Resolved by Biopsy Urease Test and/or Histology	
	Positive	Negative	Ind.	Positive	Negative
	Positive	Negative	Ind.	Positive	Negative
Premier <i>H. pylori</i>	256	8	4	260	8
Positive	17	172	5	2	192

n=462

Relative Sensitivity = 93.8% (96.1-90.2)*(256/273) 99.2% (99.8-97.2) (260/262)
 Relative Specificity = 95.5% (97.7-91.5) (172/180) 96.0% (98.0-92.3) (192/200)
 Relative Agreement = 94.5% (96.2-92.0) (428/453) 97.8% (98.8-96.1) (452/462)

Fifteen EIA positive/Premier *H. pylori* negative and all five EIA indeterminate/Premier *H. pylori* negative were negative by both histology and biopsy urease tests. Two EIA positive/Premier *H. pylori* negative were either biopsy urease positive or histology positive. Three *H. pylori* positive/EIA negative were positive by the biopsy tests. One EIA indeterminate/Premier *H. pylori* positive was biopsy urease test positive. Five EIA negative/Premier *H. pylori* positive and three EIA indeterminate/Premier *H. pylori* positive were negative by both histology and biopsy urease tests.

REPRODUCIBILITY

Each of three individual serum specimens were tested independently 10 times during one day. This process was repeated for a total of three days. All assays were performed at a clinical trial site using a negative, borderline positive and a high positive serum sample. The statistical analysis is found in the following tables.⁴

Within Run Assay Variance

SERUM TYPE	GRAND MEAN	SD	CV
Negative	0.013	0.003	23.1%
Borderline Positive	0.201	0.024	11.9%
High Positive	2.518	0.146	5.8%

Between Run Assay Variance

SERUM TYPE	GRAND MEAN	SD	CV
Negative	0.013	0.002	23.1%
Borderline Positive	0.201	0.014	15.9%
High Positive	2.518	0.302	12.0%

Total Run Assay Variance

SERUM TYPE	GRAND MEAN	SD	CV
Negative	0.013	0.004	31.0%
Borderline Positive	0.201	0.039	19.4%
High Positive	2.518	0.332	13.2%

ASSAY SPECIFICITY

A. Bacterial Studies

The specificity of Premier *H. pylori* was tested by performing adsorption studies with the following clinical isolates (CI) or ATCC strains. A seropositive or seronegative serum was reacted with 3×10^9 organisms and the adsorbed serum tested by Premier *H. pylori*. The seropositive serum remained positive after adsorption with every organism except the *H. pylori* strain. Additionally, all other organisms when reacted with a seronegative serum did not cause false positive reactions.⁴

Campylobacter coli CI	Klebsiella pneumoniae ATCC 13883
Campylobacter fetus CI	Proteus vulgaris ATCC 8427
Campylobacter hyoilealis ATCC 35217	Pseudomonas aeruginosa ATCC 9027
Campylobacter jejuni ATCC 33292	Pseudomonas fluorescens CI
Campylobacter lari CI	Salmonella sp CI
Candida albicans ATCC 10231	Serratia liquefaciens ATCC 35551
Citrobacter freundii ATCC 8090	Shigella flexneri CI
Enterobacter aerogenes ATCC 13048	Shigella sonnei CI
Enterobacter cloacae ATCC 13407	Staphylococcus aureus (Cowan I) ATCC 12598
Escherichia coli ATCC 8739	Streptococcus faecalis CI
Helicobacter cinaedi ATCC 35683	Yersinia enterolitica CI

B. Serum Studies

Sera from patients who were culture positive or had serologic evidence of the following diseases did not react when tested in Premier *H. pylori*.⁴

Anti-Nuclear Antibody (10)	Rotavirus (4)	Rheumatoid Factor (4)	Cytomegalovirus (3)
Lupus Erythematosus (3)	Rubella (5)	Influenza A (4)	Coccidioides (3)
Histoplasma (1)	Mononucleosis (16)	Herpes Simplex Virus (4)	Cryptococcus (1)

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REFERENCES

1. Alpert LC, Graham DY, Evans DJ, Yoshimura, HH Jr, Hazell SL, Evans DG, and Klein PD. 1989. Diagnostic possibilities for *Campylobacter pylori* infection. Eur. J Gastroenterol Hepatol. 1:17-26.
2. Barthel JS, and Everett ED. 1990. Diagnosis of *Campylobacter pylori* infections: The "gold standard" and the alternatives. Rev. Infect. Dis. 12:S107-S114.
3. Buck GE. *Campylobacter pylori* and gastroduodenal disease. Clin Microbiol 1990;Rev. 3:1-12.
4. Data on file, Meridian Bioscience, Inc.
5. Evans DJ Jr, Evans DG, Graham DY, and Klein PD.. A sensitive and specific serologic test for detection of *Campylobacter pylori* infection. Gastroenterology 1989;96:1004-1008.
6. Fiedorek SC, Malaty HM, Evans DL, Pumphrey CL, Casteel HB, Evans DJ, Jr, and Graham DY. Factors influencing the epidemiology of *Helicobacter pylori* infection in children. Pediatrics. 1991;88:578-582.
7. Graham DY. I. *Helicobacter pylori*: Its epidemiology and its role in duodenal ulcer disease. J Gastroenterol Hepatol 1991;6:105-113.
8. Graham DY, Klein PD, Evans DJ Jr, Evans DG, Alpert LC, Opekun AR, and Boutton TW.. *Campylobacter pylori* detected noninvasively by the ¹³C-urea breath test. Lancet i:1174-1177. 1987
9. Graham DY, Klein PD, Opekun AR, and Boutton TW. Effect of age on the frequency of active *Campylobacter pylori* infection diagnosed by the [¹³C] Urea breath test in normal subjects and patients with peptic ulcer disease. J Infect Dis 1988;157:777-780.
10. Graham DY, Malaty HM, Evans DG, Evans DJ Jr, Klein PD, and Adam E. Epidemiology of *helicobacter pylori* in an asymptomatic population in the United States. Effect of age, race, and socioeconomic status. gastroenterology. 1991;100:1495-1501.
11. Gregson DB and Simor AE. 1991. *Helicobacter pylori* and inflammatory disease of the stomach and duodenum. Stomach and bowel, Sept. 1993. pp.29-30,35-37.
12. Loeffeld RJLF, Stobberingh E, Van Spreeuwel JP, Flendrig JA and Arends JW. The prevalence of anti-*helicobacter* (*campylobacter*) *pylori* antibodies in patients and healthy blood donors. J Med Microbiol 1991;32:105-109.
13. Malaty HM, Graham DI, Klein PD, Evans DG, Adam E, and Evans DJ.. Transmission of *helicobacter pylori* infection. Studies in families of healthy individuals. Scand J Gastroenterol 1991;26:927-932.
14. Marshall BJ, Armstrong JA, McGechie DB and Glancy RJ. Attempt to fulfill Koch's postulates for pyloric campylobacter. Med J Aust 1985;142:436- 439.
15. Marshall BJ and Warren JR. 1984. Unidentified curved bacilli in the stomach of patients with gastritis and peptic ulceration. Lancet i:1311-1314.
16. Morris A, Nicholson G, Lloyd G, Haines D, Rogers A, and Taylor D. Seroepidemiology of *campylobacter pyloridis*. N Z Med J 1986;99:657-659.
17. Murray DM. Clinical relevance of infection by *Helicobacter pylori*. Clin Microbiol Lett. 1993;15:33-37.
18. Nerdeff DG.. Identification of the outer membrane proteins of *campylobacter pyloridis* and antigenic cross-reactivity between *C. pyloridis* and *C. jejuni*. J Gen Microbiol 1987;133:163-170.
19. Nichols L, Sugihay M, DeGirolami PC, Balogh K, Pleskow D, Eichelberger K, and Santos M. Evaluation of diagnostic methods for *helicobacter pylori* gastritis. Am J Clin Path 1991;95:769-773.
20. Parsonnet J, Friedman GD, Vandersteen DP, Chanu Y, Vogelman JH, Orentreich N, and Sibley RK. *Helicobacter pylori* infection and the risk of gastric carcinoma. N Engl J Med 1991;325:1127-1131.
21. Talley NJ, Newell DG, Orman JE, Carpenter HA, Wilson WR, Zinsmeister AR, Perez-Perez GI, and Blaser MJ. Serodiagnosis of *helicobacter pylori*: Comparison of enzyme-linked immunosorbent assays. J Clin Microbiol 1991;29:1635-1639.
22. Westblom TU, Madan E, Kemp J, and Subik MA. Evaluation of a rapid urease test to detect *Campylobacter pylori* infection. J Clin Microbiol 1988;26:1393-1394.

SN11140

REV. 12/23

 Manufactured By	<p>Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com</p> <p>Contacts: Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com</p>
	<p>Meridian Bioscience Europe, SRL Via Dell'Industria 7, 20035 Villa Cortese (Milano) ITALY www.meridianbioscience.com</p> <p>Contacts: Main Telephone (+39) 0331.433636 E-mail: info@meridianbioscience.eu Technical Support: MBE-TechService@meridianbioscience.eu Customer Service/Orders:<ul style="list-style-type: none">• For Italian Customers: ordini@meridianbioscience.com• For Distributors / International Customers: Export.CustomerService@meridianbioscience.eu</p>

INTERNATIONAL SYMBOL USAGE CHART

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simblos, Zeichenerklärung)

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For technical assistance, call Technical Support Services at (800) 343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at (800) 543-1980.