revogene® check

For use with the Revogene®

REF 610220

IVD For in vitro Diagnostic Use

R_x Only

INTENDED USE

The Revogene® Check is a single-use, qualitative test utilizing automated real-time Polymerase Chain Reaction (PCR) to confirm that the optical and thermal parameters of the Revogene® instrument are within specifications. It is intended for use by qualified laboratory personnel who have been trained and are proficient in performing testing on the Revogene instrument, or by Meridian Bioscience, Inc. service personnel.

Package Contents:

Each box contains two (2) Revogene Check pouches, each pouch containing one (1) Sample Buffer Tube (SBT), one (1) Disposable Transfer Tool (DTT), and one (1) Revogene Check microfluidic cartridge (PIE). Two (2) Quick Start Guides with their distinct Master barcode are also included in each box.

PRINCIPLE OF THE PROCEDURE

The Revogene instrument requires periodic verification of optical and thermal parameters, through a Revogene Check run, to confirm they are within specifications. The Revogene Check is designed to detect issues on specific instrument modules and extend the Revogene calibration deadline upon verification success by using the internal process control (PrC) of DNA-based Revogene assays.

The Revogene automates PrC homogenization, dilution, lysis, DNA amplification and detection of the amplified PCR product. User intervention is only required to transfer sample buffer from the Sample Buffer Tube (SBT) into the PIE and insert the PIE into the Revogene carousel. Each Revogene Check PIE is a fully integrated closed device into which sample buffer is dispensed and processed through different microfluidic chambers and channels, which allows for PrC processing (i.e. homogenization, dilution and lysis) and subsequent real-time PCR steps (Figure 1). The sample buffer is transferred by centrifugation from one chamber to the next in sequence.

Revogene Check PIE contains all reagents specific for the thermal and optical verification of the Revogene. The thermal verification is performed through a PCR reaction within one (1) of the PCR wells into which dried reagents are incorporated. The amplified product is detected in real-time using a target-specific TaqMan® chemistry-based probe. The optical verification is done via constant signal chemistry-based probes. No operator intervention is necessary once a PIE is loaded into the Revogene.


Figure 1. Top View of a PIE
1: Sample Loading Chamber, 2: Overflow Chamber, 3: Homogenization Chamber, 4: Dilution/Lysis Chamber, 5: Three (3) PCR wells (A to C from left to right) and one (1) Waste Chamber (at the right end).

The Revogene carousel must contain eight (8) PIEs to maintain thermodynamic balance within the run. At run completion, the results are computed by the system from measured fluorescent signals and embedded calculation algorithms. Results are displayed on the touchscreen and may be printed, transferred to a USB flash drive, or transmitted to the laboratory information system.

REAGENTS AND MATERIALS PROVIDED

The Revogene Check kit contains sufficient reagents and materials to execute two (2) tests. The kit contains two (2) individual pouches and each pouch contains the following materials:

1. One (1) **Revogene Check Sample Buffer Tube (SBT)**: Tube containing Tris-HCl pH 8.0/EDTA.Na₂ (TE 1X) buffered solution as a dilution buffer for the PrC.
2. One (1) **Disposable Transfer Tool (DTT)**: Plastic pipette with minimal and maximal volume marks for transferring sample buffer from the SBT into the PIE.
3. One (1) **Revogene Check microfluidic cartridge (PIE)**: Integrated disposable device, which comprises dried reagents allowing PrC processing and real-time PCR steps for amplification and detection of PrC DNA. Each PIE contains PrC-specific primers and probe, constant signal chemistry-based probes, dNTPs, buffer and DNA polymerase.

Each pouch is provided with one (1) **Quick Start Guide** including one (1) specific **Master barcode**. This barcoded-label is required to start a Revogene Check run. Each Master barcode contains the settings used for verification of the Revogene performance.

MATERIALS REQUIRED BUT NOT PROVIDED

1. Revogene® (cat# 610210)
2. Disposable powderless gloves
3. Vortex mixer
4. MOCK PIEs (cat# 610208)

WARNING AND PRECAUTIONS

1. The Revogene Check can only be used on the Revogene instrument.
2. Do not use the kit if the label that seals the outer box is broken upon arrival.
3. Do not use PIEs if the protective pouches are open or broken upon arrival.
4. Do not interchange DTT, SBT, PIE and Master barcode between kit lots.
5. Each single-use DTT and PIE is used to execute one (1) test. Do not reuse DTT or PIE.
6. Do not add sample or specimen into the Revogene Check PIE.
7. Wear disposable powderless gloves while handling PIEs and thoroughly wash hands afterwards.
8. The PIE contains dried reagents. The protective pouch should not be opened until ready to perform the test.
9. Dispose of unused materials, reagents and waste in accordance with country, federal, provincial, state and local regulations.
10. Do not open or break apart the PIE after use. The cap and the seal in the PIE prevent contamination with amplification product.
11. Do not use a PIE that has been dropped, shaken or inverted after the sample has been loaded as this may cause erroneous results.
12. Do not refrigerate the loaded PIE.
13. Each run must be performed with one (1) Revogene Check PIE and seven (7) MOCK PIEs in the Revogene carousel to maintain thermodynamic and mechanical balance within the run.

Any serious incident that has occurred in relation to the device should be reported to Meridian Bioscience, Inc., 3471 River Hills Drive, Cincinnati, Ohio 45244 USA or Technical Support Center 800.343.3858 and competent authority of the EU Member State in which the clinician and/or patient is established.

HAZARD and PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

STORAGE AND STABILITY

1. Store the Revogene Check kit at 2-25 C. The expiration date is indicated on the kit box label and on each pouch; do NOT use the Revogene Check kit beyond the expiration date indicated.
2. Do not open a pouch until ready to perform testing. Use the PIE within one (1) hour after opening the pouch.

INSTRUCTION FOR USE

NOTE 1: Open the kit box and remove one (1) Quick Start Guide holding one (1) pouch by pulling the Quick Start Guide tab.

NOTE 2: Do not discard the Quick Start Guide as it contains a Master barcode required to initiate the Revogene Check run.

SBT PREPARATION

NOTE 1: The contents of one (1) pouch and one (1) Master barcode are required to perform a Revogene Check test.

NOTE 2: The remaining pouch and Master barcode is to be used only in the event a repeat test is necessary.

1. Unseal the right side of the pouch (when facing label) containing the DTT and SBT and remove only the SBT from the pouch.
2. Identify (or label) the SBT with the appropriate identification without obscuring or writing over the barcode. Place the SBT on a sample rack, if used.

PIE PREPARATION

3. Mix the SBT for a minimum of 15 seconds at maximal speed using a vortex mixer.
4. Unseal the left side of the pouch (when facing label) containing the PIE, removing it from the pouch. Once unsealed, the PIE must be used within one (1) hour.
5. Remove the DTT from the right side of the pouch and use it to aspirate the sample buffer (SB) by squeezing the entire bulb. The liquid level in the DTT must be anywhere between the two (2) marks (Figure 2). If the liquid level is not between the two (2) marks, discharge the SB completely into the SBT and repeat this step until the liquid level is between the two (2) marks as indicated in **Figure 2**.
6. Discharge completely the SB volume into the sample loading chamber of the PIE (Figure 1). Make sure not to touch the outer edges or the bottom of the sample loading chamber with the DTT.
7. Close the cap of the PIE tightly. Do not refrigerate the loaded PIE.

NOTE 1: Hold the PIE by the dilution/lysis chamber or by its sides to avoid smearing on PCR wells.

NOTE 2: Be cautious when manipulating a loaded PIE. Do not invert, shake or drop.


Figure 2. Representation of an Appropriate Sample Buffer (SB) Level Using the Disposable Transfer Tool (DTT).

REVOGENE SYSTEM OPERATION

NOTE 1: Each run must be performed with one (1) Revogene Check PIE in the Revogene. Insert MOCK PIEs into the remaining seven (7) positions of the carousel. All eight (8) positions of the carousel must be filled.

NOTE 2: Refer to the Revogene Operator's Manual¹ for further information regarding Revogene set-up and operation.

1. Ensure that the Revogene is powered on.
2. Log in by entering <Username> and <Password> and tap <Login>. The main menu will appear automatically.
 - In order to start a Revogene Check run, tap <Tools> from the main menu, then tap <Revogene Check>.
 - If the calibration has reached its due date and a Revogene Check run needs to be performed, a visual warning will be shown on the "Setup Run" button displayed on the main menu.
3. Scan the single-use Master barcode provided on the Quick Start Guide by holding it vertically in front of the Revogene barcode scanner.
4. Enter the SBT and PIE barcodes using the Revogene barcode scanner. Gently position the PIE almost vertically in front of the scanner. Alternatively, SBT and PIE barcodes may be entered manually (tap the **pencil** icon of their respective lines). Handle the PIE carefully without dropping, shaking or inverting it.
5. Insert the PIE into the Revogene at any position of the carousel.
6. Confirm that the PIE is inserted by tapping <OK> at the <Insert PIE into the instrument> step.
7. Confirm that the run setup is completed by tapping <Next>.
8. Since one (1) PIE is used for the Revogene Check run, load seven (7) MOCK PIEs in the carousel remaining positions. No scan is required when inserting MOCK PIEs into the Revogene. Confirm that they were inserted by tapping the checkbox.
9. Scan the retention ring and place it on the carousel.
10. Gently rotate carousel around to check all PIEs are loaded correctly with no friction or resistance to rotation.
11. Close the instrument lid with both hands and hold on until you hear the lock activate. Then, gently pull up on the lid to guarantee that the lock is engaged.
12. Initiate the test run by tapping <Start>.

VIEWING AND EXPORTING RESULTS

NOTE 1: Refer to the Revogene Operator's Manual¹ for further information regarding the acquisition of test results.

- Once the run is completed, the lid opens automatically.
- Enter <Username> and <Password> and tap <Login> if the user's session has logged-out.
- The Revogene Check completion message (see **Results Interpretation** section) is automatically shown on the screen.
- Tap the <Show Revogene Check Menu> button of the Revogene Check completion message to access the Revogene Check Report or tap the <Close> button to return to the main menu.
- If the <Show Revogene Check Menu> is chosen, tap the <Show reports> button (**Figure 3**) to consult the pdf report file.
 - Use the right side scroller to see the detailed result by module.
 - Tap <Export> to export the report and save where appropriate (e.g. USB flash drive or via connectivity option).

Note that Revogene Check data can be exported as part of the "troubleshooting package" as needed. Consult the "Tools-Support-Troubleshooting" section in the Operator's Manual¹ for further information.

- Remove the retention ring and the PIE/MOCK PIEs from the Revogene:
 - Dispose of the PIE in accordance with your institution's standard practices.
 - Do not discard MOCK PIEs as they can be reused.


Figure 3. Revogene Check Settings

RESULTS INTERPRETATION

Results are computed by the Revogene from measured fluorescent signals and embedded calculation algorithms and are displayed on the "Results" window. Possible results are listed below.

i. Revogene Check run succeeded.

Revogene verification completed. The instrument calibration deadline is extended for an additional 12 months. Resume testing of patient specimens.




ii. Revogene Check needs to be repeated.

Revogene verification indicates that calibration is required (see an example of Overall Result under the screenshot). At this point, the instrument is locked preventing its use with Revogene assays. Repeat the process, as recommended by the on-screen Revogene Check completion message, with the second PIE in the Revogene Check kit box.

iii. Instrument needs attention, please contact technical support.

Revogene verification indicates that calibration is required (see an example of Overall Result under the screenshot). At this point, the instrument is locked preventing its use with Revogene assays.

If the prompted window shown below appears upon completion of Revogene Check, contact Technical Support at: 1-800-343-3858. Do not attempt to run the Revogene Check more than twice.

One (1) example of possible Overall Result in the Revogene Check Report is presented below.

One (1) example of possible Overall Result in the Revogene Check Report is presented below.

Overall Result

PASS

Detailed Result By Module

Performance Module Tested	Result
PCR System	PASS
Optical Detection	PASS
Interface Communications	PASS

Overall Result

Calibration Required

Detailed Result By Module

Performance Module Tested	Result
PCR System	Calibration Required
Optical Detection	Calibration Required
Interface Communications	PASS

Overall Result


Calibration Required

Detailed Result By Module

Performance Module Tested	Result
PCR System	Calibration Required
Optical Detection	Calibration Required
Interface Communications	PASS

❖ **Revogene Check run has encountered an error.**

If an error is encountered by the instrument, Revogene verification cannot be completed. No Revogene Check result is available at this point. The user will be required to perform a self-test before being allowed to execute another Revogene Check run. Refer to the Revogene Operator's Manual¹ for further information regarding the self-test function.

Overall Result

ERROR

Detailed Result By Module

Performance Module Tested	Result
PCR System	ERROR
Optical Detection	ERROR
Interface Communications	ERROR

LIMITATIONS OF THE PROCEDURE

1. The Revogene Check must only be used with the Revogene instrument by trained personnel.
2. Only Revogene Check PIEs and MOCK PIEs can be used during the run.
3. Performing the Revogene Check does not replace use of external controls to verify Revogene assay reagent performance when used on the Revogene instrument.

E-LABELING


Documentation related to this product can be accessed online at www.meridianbioscience.com/pi. Additionally, paper copies are available upon request by contacting your local distributor or via the phone number listed on the kit box.

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Overall Result	PASS
Detailed Result By Module	
Performance Module Tested	Result
PCR System	PASS
Optical Detection	PASS
Interface Communications	PASS

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Overall Result	Calibration Required
Detailed Result By Module	
Performance Module Tested	Result
PCR System	Calibration Required
Optical Detection	Calibration Required
Interface Communications	PASS

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INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

	Use-by date / Data di scadenza / Date de péremption / Fecha de caducidad / Verwendbar bis
	Batch code / Codice di lotto / Code de lot / Código de lote / Chargennummer
	In vitro diagnostic medical device / Dispositivo medico diagnostico <i>in vitro</i> / Instrument de test diagnostique <i>in vitro</i> / Producto sanitario para diagnóstico <i>in vitro</i> / In-vitro-Diagnostikum
	CE Mark / Marcatura CE / Symbole CE / Marcado CE / CE-Kennzeichen
	Catalog number / Numero di catalogo / Référence catalogue / Referencia / Bestellnummer
	Consult instructions for use / Consultare le istruzioni per l'uso / Consulter le mode d'emploi / Consultar las instrucciones de uso / Gebrauchsanweisung beachten
	Manufacturer / Produttore / Fabricant / Fabricante / Hersteller
	Contains sufficient for <n> tests / Contiene una quantità sufficiente per <n> test / Contient le matériel suffisant pour <n> tests / Contiene la cantidad suficiente para <n> ensayos / Inhalt ausreichend für <n> Tests
	Temperature limit / Limite di temperatura / Limite de température / Limite de temperatura / Temperaturgrenze
	Authorized representative in the European Community / Rappresentante autorizzato nella Comunità Europea / Représentant agréé dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter EU-Repräsentant
	Do not reuse / Non riutilizzare / Ne pas réutiliser / No reutilizar / Nicht wiederverwenden
	Keep dry / Conservare all'asciutto / Conserver au sec / Mantener seco / Vor Feuchtigkeit schützen
	Contains # pouches: 1 Disposable Transfer Tool (DTT), 1 Sample Buffer Tube (SBT), 1 PIE / Contiene # buste: 1 PIE, 1 Sample Buffer Tube (SBT), 1 Disposable Transfer Tool (DTT) / Contient # sachets: Une cartouche PIE, Un tube échantillon, Un outil de transfert jetable (DTJ) / Incluye # bolsitas: 1 PIE, 1 Sample Buffer Tube (SBT), 1 Disposable Transfer Tool (DTT) / Enthält # Beutel: 1 PIE, 1 Sample Buffer Tube (SBT), 1 Disposable Transfer Tool (DTT)
	Humidity Limitation / Limitazione dell'umidità / Limite d'humidité / Limitación de humedad / Feuchtigkeitsbegrenzung
	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Ne pas utiliser si l'emballage est endommagé / No utilizar si el envase está dañado / Nicht verwenden, wenn die Verpackung beschädigt ist
	Keep away from sunlight / Tenere lontano dalla luce del sole / Tenir à l'écart de la lumière du soleil / Mantener alejado de la luz solar / Vor Sonnenlicht schützen
	Contains 24 Disposable Transfer Loops (DTL) / Contiene 24 loop di trasferimento monouso (DTL) / Contient 24 boucles de transfert jetables (DTL) / Contiene 24 bucles de transferencia desechables (DTL) / Enthält 24 Einwegübertragungsschleifen (DTL)
	Disposable Transfer Tool (DTT) / Disposable Transfer Tool 9DTT / Un outil de transfert jetable (DTJ) / Disposable Transfer Tool (DTT) / Disposable Transfer Tool (DTT)
	Revogene Test Device / Dispositivo Test Revogene / Dispositif de test Revogene / Dispositivo para la Prueba Revogene / Revogene-Analysegefäß
	For Emergency Use Authorization only / Solo per l'autorizzazione per l'uso di emergenza / pour autorisation d'utilisation d'urgence uniquement / para autorización de uso de emergencia solamente / nur für Notfallverwendungsautorisierung
	Disposable Transfer Loop / Loop di trasferimento monouso / Boucles de transfert jetables / Bucle de transferencia desechables / Einwegübertragungsschleifen
	Sample Buffer Tube / Sample Buffer Tube / Un tube échantillon / Sample Buffer Tube / Sample Buffer Tube
	Revogene Test Device / Dispositivo Test Revogene / Dispositif de test Revogene / Dispositivo para la Prueba Revogene / Revogene-Analysegefäß
	Swiss Authorized Representative / Mandatario svizzero / Mandataire Suisse / Representante Autorizado Suizo / Schweizer Bevollmächtigter

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order call Customer Service Department at 800-543-1980.

TaqMan is a registered trademark of Roche Molecular Systems, Inc.
Revogene and associated logos are registered trademarks of Meridian Bioscience, Inc.
© 2020-12 Meridian Bioscience, Inc.
Made in Canada

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TaqMan è un marchio registrato di Roche Molecular Systems, Inc.
Revogene e i relativi loghi sono marchi registrati di Meridian Bioscience, Inc.
TaqMan est une marque commerciale déposée de Roche Molecular Systems, Inc.
Revogene et les logos connexes sont des marques déposées de Meridian Bioscience, Inc.
© 2020-12 Meridian Bioscience, Inc. Fabriqué au Canada
TaqMan es una marca registrada de Roche Molecular Systems, Inc.