



PREMIER
C. difficile GDH TEST PROCEDURE *How to perform the test*


SPECIMEN
 PROCESSING



1. Measure 200 μ L of Sample Diluent/ Negative Control into a test tube.




2. Mix stool thoroughly. Using a transfer pipette add 50 μ L (1st mark from tip) to the tube.




3. Vortex for 15 seconds.


SAMPLE / ENZYME
 CONJUGATE INCUBATION



1. Detach microwells needed, place in holder. Using transfer pipette, add 100 μ L (2nd mark from tip) of the diluted stool to the appropriate wells.




2. Add 2 drops of Positive Control and 100 μ L (2nd mark from tip) of Sample Diluent/Negative Control to the appropriate wells.




3. Add 1 drop of Enzyme Conjugate to all wells. Shake firmly for 30 seconds. Seal the plate and incubate at 35-39 C for 50 minutes or at 37 C for 20 minutes while rotating the plate at 1000 rpm (using the Stat-Fax rotator ONLY).

SUBSTRATE
 INCUBATION



1. Wash 5-7 times with 1X Wash Buffer. (See package insert for proper wash procedure.)



2. Add 2 drops of Substrate I to all wells. Shake firmly for 30 seconds and incubate for 10 minutes at 21-27 C.



3. Add 2 drops of Stop Solution I to all wells and shake firmly for 30 seconds.

INTERPRETATION
 OF RESULTS

Spectrophotometric Single Wavelength (450 nm)
 (read within 30 minutes)
 Negative = OD 450 < 0.200
 Positive = OD 450 \geq 0.200

Spectrophotometric Dual Wavelength (450/630 nm)
 (read within 30 minutes)
 Negative = OD 450/630 < 0.150
 Positive = OD 450/630 \geq 0.150



Please refer to the most current version of the package insert for complete instructions.



© 2010 SN11083 REV.03/11