PREMIER®

C. difficile GDH

Enzyme immunoassay for the detection of *Clostridium difficile* GDH in stool specimens

REF 611096

IVD In vitro diagnostic medical device

INTENDED USE

Premier *C. difficile* GDH is a qualitative enzyme immunoassay screening test to detect *Clostridium difficile* antigen, glutamate dehydrogenase, in fecal specimens from symptomatic persons suspected of having *C. difficile* infection (CDI). This test does not distinguish between toxigenic and non-toxicogenic strains of *C. difficile*. Samples from symptomatic patients that produce positive results with this test must be further tested with an assay designed to detect toxigenic *C. difficile* strains and assist with the diagnosis of CDI.

SUMMARY AND EXPLANATION OF THE TEST

Clostridium difficile has become the most important etiologic agent in nosocomial diarrhea worldwide. *C. difficile* infection (CDI) can range from mild diarrhea to life-threatening toxic megacolon. Classic risk factors for nosocomial CDI include antibiotic use and old age. Recently however CDI has been seen more frequently outside the hospital environment without the typical risk factors historically associated with *C. difficile* infection.¹ *C. difficile* produces two toxins (toxin A and toxin B) that have both been associated with disease caused by *C. difficile*. Strains of *C. difficile* that do not harbor the toxin A and toxin B genes are considered to be non-pathogenic. All strains of *C. difficile* produce glutamate dehydrogenase (GDH), also known as common antigen.^{2,3} Non-molecular laboratory approaches to *C. difficile* detection can be divided into two general groups: toxin-based and non-toxin based. Toxin-based assays include enzyme-immunoassay (EIA) for toxin A or toxins A and B and cytotoxin, which primarily detects toxin B. Non-toxin based assays, such as EIA for *C. difficile* GDH, while inherently more sensitive, must be followed up with testing for toxigenic *C. difficile*.

BIOLOGICAL PRINCIPLES

Premier *C. difficile* GDH is a qualitative enzyme immunoassay screening test to detect *Clostridium difficile* antigen, glutamate dehydrogenase, in fecal specimens from symptomatic persons suspected of having *C. difficile* infection (CDI). Breakaway microwells are coated with rabbit anti-GDH antibodies. Diluted patient specimen and a horseradish peroxidase conjugated anti-GDH antibody are added to the microwells and incubated. Upon completion of the incubation, a wash step is performed to remove unbound material. If *C. difficile* GDH is present, an antibody-enzyme complex is formed. A chromagenic substrate is added to the microwells and incubated. A blue color develops in the presence of bound enzyme. Stop solution is added, changing the initial blue reaction to yellow. Test results are interpreted spectrophotometrically.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. Premier *C. difficile* GDH Microwells: Polyclonal antibody-coated microwells. The antibodies are specific to *C. difficile* GDH.
2. Premier *C. difficile* GDH Enzyme Conjugate: HRP-conjugated polyclonal antibodies specific to *C. difficile* GDH in a buffered protein solution containing 0.1% ProClin® and 0.03% gentamicin as preservatives.
3. Premier 20X Wash Buffer II: Concentrated wash buffer containing 0.2% thimerosal as a preservative. Reagent is diluted before use.
4. Premier Substrate I: Buffered solution containing urea peroxide and tetramethylbenzidine at pH 5.0.
5. Premier Stop Solution I: 1M Phosphoric acid.
6. Premier *C. difficile* GDH Sample Diluent/Negative Control: Buffered protein solution containing 0.1% ProClin® and 0.03% gentamicin as preservatives.
7. Premier *C. difficile* GDH Positive Control: *C. difficile* GDH in a buffered protein solution containing 0.1% ProClin® and 0.03% gentamicin as preservatives.
8. Transfer pipettes
9. Microwell strip holder
10. Microwell plate sealers

MATERIALS NOT PROVIDED

1. Disposable latex gloves
 2. Test tubes (eg, 10 x 75 mm or 12 x 75 mm) for dilution of sample
 3. Distilled or deionized water
 4. Squirt bottle
 5. Graduated cylinder for making 1X Wash Buffer
 6. Absorbent paper
 7. Waste container with disinfectant and/or autoclavable biohazard bags
 8. Vortex mixer
 9. Interval timer
 10. EIA microplate reader capable of reading absorbance at 450 nm or 450/630 nm*
 11. StatFax™-2200 Incubator/Shaker (StatFax™ is a trademark of Awareness Technology, Inc.) (optional)*
 12. Semiautomated microplate washer (optional)*
- * Note: It is the operator's responsibility to validate StatFax™, semiautomated plate washers and readers prior to their use with this product.

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Do not interchange the Microwells, Enzyme Conjugate, Premier Substrate I, or Positive Control reagents between lots. (The Sample Diluent/Negative Control, Premier 20X Wash Buffer II and Premier Stop Solution I are interchangeable provided the reagents are within their assigned expiration dates when used.)
3. Do not use kit components beyond labeled expiration date.
4. Do not use vials that lack a label, a lot number, or an expiration date.
5. Do not use any reagent if it is discolored or turbid. Discoloration or turbidity may be a sign of microbial contamination.
6. Allow reagents to warm to 21–27 C before use.
7. All reagents should be mixed gently before use.
8. Hold reagent vials vertically at suitable distance above the well to insure proper drop size and delivery.
9. Replace colored caps on correct vials.
10. Do not reuse microwells.
11. Unused microwells must be placed back inside resealable pouch. It is important to protect strips from moisture.
12. The transfer pipettes provided with this kit must be used for specimen preparation and transfer. Use one per specimen.
13. Avoid splashing when dispensing diluted stool into microwells by placing the transfer pipette tip about halfway down the well and dispensing slowly down the side of well.
14. Microwell washing is to be performed precisely as directed in the assay procedure. **Inadequate washing may be the cause of elevated background in any EIA protocol.**
15. All reagents except the Premier 20X Wash Buffer II are provided already diluted to the proper concentration.
16. Any deviation below or above set incubation times may affect sensitivity and specificity and should be avoided.
17. Stool must be mixed thoroughly (regardless of consistency) to ensure a representative sample prior to pipetting.

WARNINGS

1. Patient specimens may contain infectious agents and should be handled and disposed of as potential biohazards.
2. Dispose of used wash buffer and all test materials in an appropriate container. Treat waste as a potential biohazard.
3. The Positive Control reagent contains *C. difficile* GDH. It should be handled as a potential biohazard.
4. Avoid skin contact with Premier Stop Solution I. Flush with water immediately if contact occurs.

HAZARD and PRECAUTIONARY STATEMENTS

Refer to the SDS, available at www.meridianbioscience.com for Hazards and Precautionary Statements.

SHELF LIFE AND STORAGE

The kit expiration date is indicated on the kit label. Store the kit at 2–8 C and return the kit promptly to the refrigerator after each use.

PROCEDURAL NOTES

The Premier *C. difficile* GDH transfer pipette is diagrammed below.

REAGENT PREPARATION

1. Bring the entire kit, including microwell pouch, to 21–27 C before use.
2. Prepare 1X Wash Buffer as needed. For example: 10.0 mL of Premier 20X Wash Buffer II + 190.0 mL distilled or deionized water is sufficient to wash one strip. 1X Wash Buffer can be stored at 21–27 C for up to three months.

SPECIMEN COLLECTION AND PREPARATION

Human stool samples, unpreserved: Samples should be received in an airtight transport container and stored at 2–8 C prior to testing. Samples should be tested as soon as possible, but may be held up to five days at 2–8 C. (See SPECIMEN PREPARATION section for instructions on diluting samples.) Samples that will not be tested within five days should be frozen immediately upon receipt and stored at ≤ –20 C until tested. Specimens may be held frozen for up to two months and may be frozen and thawed twice.

SPECIMEN PREPARATION

Mix stool as thoroughly as possible prior to pipetting.

1. **Human stool samples, unpreserved:**
 - a. With the dropper assembly (or equivalent), add 200 µL of Sample Diluent/Negative Control to a small test tube.
 - b. **Liquid/Semi-Solid stools:**
 - i. Using a transfer pipette provided with the kit, add 50 µL (first mark from pipette tip) of thoroughly mixed stool to the Sample Diluent/Negative Control test tube.
 - ii. Vortex the mixture for a minimum of 15 seconds.
 - iii. Save the transfer pipette in the test tube for future use.
 - c. **Stool diluted in Sample Diluent:** If needed, stool diluted in Sample Diluent can be held at 2–27 C for up to 24 hours before testing provided the tube is sealed.
 - d. Stool specimens may be centrifuged after dilutions. Centrifuge at approximately 2750 x G for 3 minutes or until solid matter separates from liquid. Proceed with the assay after recovering the supernatant.

TEST PROCEDURE

1. After the pouch has reached temperature (21–27 C), break off the required number of microwells (1 well for each specimen, plus 1 positive and 1 negative control well per batch). Place the microwells in the microwell strip holder and record the location of all wells. Unused microwells must be resealed in the pouch immediately.
2. Using the specimen transfer pipette, add 100 µL of diluted stool (second mark from the tip of the pipette) to the appropriate well. (Place the pipette tip halfway into well and let the sample slowly run down side of well.)
3. Add 2 free falling drops of Positive Control. With a transfer pipette, add 100 µL (second mark from the tip of the pipette) of Sample Diluent/Negative Control to the appropriate wells.
4. Add 1 free falling drop (approximately 50 µL) of Enzyme Conjugate to each well. Firmly shake/swirl the plate for 30 seconds.
5. Cut the plate sealer to size and press it firmly onto the top of the microwells to seal. Incubate the plate for 50 minutes at 35–39 C. **Alternatively, laboratories equipped with a heated plate shaker (StatFax™-2200) can incubate and rotate the plate for 20 minutes at 37 C at 1000 rpm (setting #5).**
6. Carefully remove the plate sealer and wash the wells:
 - a. **Manual method:**
 - i. Dump the plate contents firmly into a biohazard receptacle.
 - ii. Invert and bang the plate on a clean stack of paper towels.
 - iii. Fill all wells with 1X Wash Buffer, directing a stream of buffer to the sides of the wells to avoid foaming.
 - iv. Repeat the wash cycle (dump, bang on fresh towels, fill) 4 to 6 more times for a total of 5–7 wash cycles. After the last fill, dump and bang plate on fresh towels hard enough to remove as much excess wash buffer as possible, but do not allow wells to completely dry at any time.
 - b. **Semi-automated wash method:**
 - i. Dump the plate contents firmly into a biohazard receptacle.
 - ii. Invert and bang the plate on a clean stack of paper towels, then place empty plate on washer device.
 - iii. Fill the wells to the top (approx. 300–350 µL/well) with 1X Wash Buffer and then aspirate. The washer manifold should be adjusted to ensure no foaming occurs during the filling of the wells and that the wells are thoroughly aspirated after each wash.
 - iv. Repeat step iii 4 more times. Following the last wash, test wells should be thoroughly aspirated to remove as much moisture as possible, but do not allow wells to completely dry at any time.
7. Add 2 free-falling drops (approximately 100 µL) of Premier Substrate I to each well. Firmly shake/swirl the plate for 30 seconds. Incubate the plate for 10 minutes at 21–27 C.
8. Add 2 free-falling drops (approximately 100 µL) of Premier Stop Solution I to each well. Firmly shake/swirl the plate for 30 seconds. **Note:** The initial color of a positive reaction is blue, which changes to yellow upon addition of Premier Stop Solution I.
9. Clean the underside of all wells with a lint-free tissue.
10. Inspect and record reactions. Test results should be read using a spectrophotometric reader. Zero EIA reader on air. Read absorbance at 450 nm or 450/630 nm within 30 minutes of adding Premier Stop Solution I.
11. Suspect failure of the washing method or device if the Negative and/or Positive Controls consistently produce out-of-specification results. Increasing the number of washes, washing more vigorously, decanting more thoroughly or recalibrating washing devices should correct the problem. If problems persist, contact Meridian's Technical Support Services at 800-343-3858 (USA) or your authorized Meridian distributor for assistance.

INTERPRETATION OF RESULTS

Spectrophotometric Single Wavelength (450 nm)

Negative: < 0.200
Positive: ≥ 0.200
Negative Control: < 0.150
Positive Control: ≥ 0.600

Spectrophotometric Dual Wavelength (450/630 nm)

Negative: < 0.150
Positive: ≥ 0.150
Negative Control: < 0.100
Positive Control: ≥ 0.600

If a Negative Control is < 0.000, reblank the plate reader to air and reread the plate.

A positive result indicates the presence of *C. difficile* GDH. A negative result indicates the absence of *C. difficile* GDH, or that the level of GDH is below what can be detected by the assay. The magnitude of the OD above the cutoff is neither indicative of the severity or extent of *C. difficile* infection, nor can it be correlated to an endpoint titer. Extremely strong positive samples may yield either an intense yellow color or a purple precipitate within a few minutes of stopping the reaction. In this case, the spectrophotometer may yield an "out" reading. This reading is considered a positive.

If the frequency of low positive results (OD between 0.200 and 0.250 for single wavelength and between 0.150 and 0.200 with dual wavelength) is greater than 5% of specimens tested, this may indicate insufficient washing. More vigorous washing or increasing the washes to seven washes in step 6 of the procedure is recommended.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

- At the time of each use, kit components should be visually examined for obvious signs of microbial contamination, freezing or leakage. Do not use contaminated or suspect reagents.
- The performance of Premier C. difficile GDH should be verified using the Positive Control and Sample Diluent/Negative Control with each test batch. See the INTERPRETATION OF RESULTS section above for a description of the expected results for control reagents. Tests should be considered invalid when either control reagent does not produce the expected results. In such cases, repeat tests and controls. If, on repeat testing, the expected reactions are still not observed with in-date reagents, call Meridian's Technical Support Services at 800-343-3858 (USA) or your authorized Meridian distributor for assistance.
- The controls are used to monitor reagent reactivity. Failure of the controls to produce the expected results can mean that one or more of the reagents are no longer reactive at the time of use, the test was not performed correctly, or that reagents or samples were not added. **If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.**
- Specimen matrix interference has not been observed in this assay as samples are significantly diluted before testing in Sample Diluent. For this reason, the Positive and Negative Controls supplied as part of this assay are prepared in matrices similar to the Sample Diluent. If control materials that are identical in composition to test samples are preferred, the user can prepare these by diluting known positive and negative specimens in Sample Diluent according to the SPECIMEN PREPARATION section of this insert. Add 100 µL of user prepared controls to test wells.

EXPECTED VALUES

The frequency of antibiotic-associated diarrhea caused by C. difficile is dependent on several factors including: patient population, type of institution and epidemiology. The reported incidence of C. difficile infection in patients suspected of having antibiotic-associated diarrhea is 15-25%⁴ although different facilities may find positivity rates outside this range.

LIMITATIONS OF THE PROCEDURE

- For in vitro diagnostic use only.
- The test is qualitative and no quantitative interpretation should be made with respect to the values.
- A positive result does not define the existence of C. difficile disease or the presence of toxigenic C. difficile, but only indicates the detection of the organism. All Premier C. difficile GDH positive stools should be tested to verify the presence of toxigenic C. difficile.
- Overincubation of the tests may lead to an increase in false-positive test results. Conversely, incubation for periods less than those defined in this insert can result in an increase in false-negative tests. Follow incubation times defined in this insert.
- Staphylococcus aureus (Cowan strain I) and Clostridium sporogenes were found to be crossreactive with Premier C. difficile GDH.
- This test is for use with unfrozen unpreserved stool samples only. Performance characteristics of other clinical specimen types have not been established.

SPECIFIC PERFORMANCE CHARACTERISTICS

Premier C. difficile GDH was evaluated by five independent test sites located in the Southwestern, Southeastern, and Midwestern regions of the United States. A total of 733 qualified patient samples were evaluated; all specimens were prospectively collected samples. Samples were collected from males (46%) and females (53%). In the case of 1% of the patients, the gender was not known. The age groups of the patients ranged from 22 days to 99 years of age. No differences in test performance were observed based on patient age or gender. The following tables show the assay performance by clinical site and patient age.

Table 1. Performance characteristics by clinical site

Site	Positive Samples			Negative Samples		
	Premier GDH/ Culture	Sensitivity %	95% CI	Premier GDH/ Culture	Specificity %	95% CI
Total Sites	108/117	92.3%	86.0 – 95.9%	590/616	95.8%	93.9 – 97.1%
Site 1	13/16	81.3%	57.0 – 93.4%	84/87	96.6%	90.3 – 98.8%
Site 2	28/30	93.3%	78.7 – 98.2%	132/140	94.3%	89.1 – 97.1%
Site 3	44/46	95.7%	85.5 – 98.8%	147/153	96.1%	91.7 – 98.2%
Site 4	15/15	100.0%	79.6 – 100.0%	169/175	96.6%	92.7 – 98.4%
Site 5	8/10	80.0%	49.0 – 94.3%	58/61	95.1%	86.5 – 98.4%

- Discrepant specimens were evaluated using an FDA-cleared ELISA test for the detection of C. difficile GDH.
- Sixteen of the 26 false positive specimens were positive when tested with another FDA cleared GDH assay.
- Eight of the nine false negative specimens were negative when tested with another FDA cleared GDH assay.

Table 2. Performance characteristics by patient age

Patient Age	Positive Samples			Negative Samples		
	Premier GDH/ Culture	Sensitivity %	95% CI	Premier GDH/ Culture	Specificity %	95% CI
0-28 days	0/0	N/A	N/A	1/1	100.0%	20.7 – 100.0%
29 days to 2 years	22/24	91.7%	74.2 – 97.7%	56/60	93.3%	84.1 – 97.4%
> 2 years to < 12 years	24/25	96.0%	80.5 – 99.3%	106/109	97.2%	92.2 – 99.1%
12 years to < 18 years	9/11	81.8%	52.3 – 94.9%	59/61	96.7%	88.8 – 99.1%
18 years to 21 years	3/4	75.0%	30.1 – 95.4%	22/23	95.7%	79.0 – 99.2%
> 21 years	50/53	94.3%	84.6 – 98.1%	345/361	95.6%	92.9 – 97.3%
Not Defined	0/0	N/A	N/A	1/1	100.0%	20.7 – 100.0%

REPRODUCIBILITY

Assay precision, intra-assay variability and inter-assay variability were assessed with a reference panel prepared from pools of negative samples spiked with C. difficile GDH antigen. The reproducibility panel consisted moderately positive (n=3), low positive (n=3), high negative (n=3), and negative (n=1) specimens. The low positive and high negative specimens were prepared near the assay limit of detection. Each sample was evaluated twice per day for five days by three different laboratories. Reproducibility was 100% with no intra-assay and inter-assay variability for samples prepared above or below the limit of detection for the assay.

CROSSREACTIVITY STUDIES

Crossreactivity studies were performed with positive and negative stool specimens inoculated with bacterial or fungal organisms to a final concentration of 1.2 x 10⁶ CFU/mL or virus concentration greater than 1 x 10⁵ TCID₅₀/mL. None of the following organisms in stool reacted with Premier C. difficile GDH:

Aeromonas hydrophila, Bacillus cereus, Bacillus subtilis, Bacteroides fragilis, Campylobacter coli, Campylobacter fetus, Campylobacter jejuni, Candida albicans, Citrobacter freundii, Clostridium butyricum, Clostridium bifermens, Clostridium histolyticum, Clostridium novyi, Clostridium perfringens, Clostridium septicum, Clostridium sordellii, Clostridium tetani, Enterobacter aerogenes, Enterobacter cloacae, Enterococcus faecalis, Escherichia coli, Escherichia coli O157:H7, Escherichia hermannii, Escherichia fergusonii, Helicobacter pylori, Klebsiella pneumoniae, Lactococcus lactis, Listeria monocytogenes, Peptostreptococcus anaerobius, Pleisiomonas shigelloides, Porphyromonas asaccharolytica, Proteus vulgaris, Pseudomonas aeruginosa, Pseudomonas fluorescens, Salmonella Group B, Salmonella Group C, Salmonella Group D, Salmonella Group E, Serratia liquefaciens, Serratia marcescens, Shigella boydii, Shigella flexneri, Shigella sonnei, Staphylococcus aureus, Staphylococcus epidermidis, Vibrio parahaemolyticus, Yersinia enterocolitica, Adenovirus Type 40, Adenovirus Type 41, Coxsackievirus Strain B4, Echovirus Strain 11, Rotavirus Strain WA

Stools spiked with Staphylococcus aureus (Cowan strain I) and Clostridium sporogenes, were found to be crossreactive with Premier C. difficile GDH.

ANALYTICAL SENSITIVITY

The analytical sensitivity of the Premier C. difficile GDH test is 8 ng/mL when spiked in stool.

ASSAY REACTIVITY

The following C. difficile stock cultures from different sources were tested and produced positive reactions at 5.7 x 10⁷ CFU/mL with Premier C. difficile GDH:
Toxigenic C. difficile strains: 8864, 10463, 43598, 2004052, 2004111, 2004118, 2004205, 2004206, 2005070, 2005257, 2005325, 2006240, 2007431, 2007435, 2007858, 2008016, 2008029, 2008162, 2008188, 2008341, 2008351, 2009018, 2009065, 2009066, 2009099, 2009132, 2009155, 2009277, B1, B17, B18, BK6, CF1, G1, J7, K12, X151

Non-toxigenic C. difficile strains: 11186, 234, 586, 611, 620, 2C62, 2C165, C122, UNC19904, X15076

TESTS FOR INTERFERING SUBSTANCES

The following substances, at the specified saturated solvent/diluent concentrations, do not interfere with test results in the final concentrations listed:

Barium sulfate at 5 mg/mL, Stearic acid/Palmitic acid (fecal fat) at 2.65mg/mL/1.3 mg/mL, Hemoglobin at 3.2 mg/mL, Imodium AD® (Loperamide HCl) at 6.67x10⁻³ mg/mL, Kaopectate® (Bismuth subsalicylate) at 0.87 mg/mL, Metronidazole at 12.5 mg/mL, Mucin at 3.33 mg/mL, Mylanta® (Ammonium hydroxide w/ magnesium hydroxide) at 4.2 mg/mL, Pepto Bismol® (Bismuth subsalicylate) at 0.87 mg/mL, Polyethylene glycol 3350 at 79.05 mg/mL, Prilosec® (Omeprazole) at 0.5 mg/mL, Simethicone at 0.625 mg/mL, Tagamet® (Cimetidine) at 0.5 mg/mL, Tums® (Calcium carbonate) at 0.5 mg/mL, Vancomycin hydrochloride at 2.5 mg/mL, White blood cells at 5%, Whole blood at 25%

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REFERENCES

- Dawson LF, Valente E, Wren B. *Clostridium difficile* - A continually evolving and problematic pathogen. Infect Genet Evol 2009;9:1410-1417.
- Lyerty DM, Barroso LA, Wilkins TD. Identification of the latex-reactive protein of *Clostridium difficile* as glutamate dehydrogenase. J Clin Micro 1991;29:2639-2642.
- Willis, DH, Kraft JA. Confirmation that the latex-reactive protein of *Clostridium difficile* is a glutamate dehydrogenase. J Clin Micro 1992;30:1363-1364.
- Bartlett JG, Gerding D. Clinical Recognition and Diagnosis of *Clostridium difficile* infection. Clin Inf Dis 2008;46:S12-18.

SN11082

REV. 12/14

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de símbolos, Erläuterung der graphischen symbole)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenzubereitung, in dem sich Probenverdünnung befindet
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Reactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	SOLN STOP	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Conjugato enzimatico / Conjugué enzymatique / Conjugado enzimático / Enzymkonjugat
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	CONTROL	Assay Control / Controllo dei test / Test de contrôle / Control de Ensayo / Kontrolltest
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	REAG	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweis
CONJ	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
SUBS	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
		DET REAG	Detection Reagent / Reagente Diretto / Réactif de Détection / Reactivo de Detección / Nachweis Reagenz

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.

Manufactured by

Meridian Bioscience, Inc.
 USA/Corporate Office
 3471 River Hills Drive
 Cincinnati, Ohio 45244
 Telephone: 513.271.3700
 Orders/Customer Service:
 800.543.1980
 Technical Support Center:
 800.343.3858
 Information Fax: 513.272.5432
 Ordering Fax: 513.271.0124

Authorized Representative

Meridian Bioscience Europe S. r. L.
 Via dell'Industria, 7
 20020 Villa Cortese, Milano
 ITALY
 Tel: +39 0331 43 36 36
 Fax: +39 0331 43 36 16
 Email: info@meridianbioscience.eu
 WEB: www.meridianbioscience.eu

Meridian Bioscience Europe s.a./n.v.
 2 Avenue du Japon - 1420 Braine l'Alleud
 BELGIUM
 Tel: +32 (0) 67 89 59 59
 Fax: +32 (0) 67 89 59 58
 Email: info.bn@meridianbioscience.eu

Meridian Bioscience Europe France
 34 rue de Ponthieu - 75008 Paris
 FRANCE
 Tel: +33 (0) 1 42 56 04 40
 Fax: +33 (0) 9 70 06 62 10
 Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
 Postbus 301 - 5460 AH Veghel
 NETHERLANDS
 Tel: +31 (0) 411 62 11 66
 Fax: +31 (0) 411 62 48 41
 Email: info.bn@meridianbioscience.eu