# PREMIER® TOXINS A & B

## Enzyme Immunoassay for the Detection of *Clostridium difficile* Toxin A and Toxin B in Stool Specimens

REF 616096

IVD In vitro diagnostic medical device

### INTENDED USE

Premier Toxins A&B is a qualitative enzyme immunoassay for the detection of *Clostridium difficile* toxin A and toxin B in stool from patients with antibiotic-associated diarrhea. Premier Toxins A&B is intended for use as an aid in the diagnosis of *C. difficile* associated disease (CDAD).

### SUMMARY AND EXPLANATION OF THE TEST

Toxigenic *Clostridium difficile* is a major cause of antibiotic-associated diarrhea and colitis<sup>1</sup> and is the causative agent for virtually all cases of pseudomembranous colitis.<sup>2,3</sup> Although about 2% of normal healthy adults are colonized with *C. difficile*,<sup>4</sup> many patients acquire this organism through nosocomial infection.<sup>5</sup> Exposure to most antibiotics is thought to allow proliferation of toxigenic *C. difficile* by disrupting the normal intestinal flora.<sup>6</sup> Two toxins, toxin A and toxin B, are associated with disease caused by *C. difficile*.<sup>7</sup> These toxins are immunochemically and biologically distinct. Antiserum prepared against purified toxin A or toxin B does not cross-react with the other toxin.<sup>8</sup> Toxin A has been described as an enterotoxin and causes an increase in intestinal permeability with subsequent enteric fluid accumulation and diarrhea.<sup>9</sup> Toxin B is a potent cytotoxin which causes rounding of cells in culture.<sup>7,10</sup> In hamsters, toxin B is lethal when administered intravenously by itself or intragastrically if combined with sublethal doses of toxin A.<sup>11</sup> The contribution of toxin B to the development of disease in the gut is not understood. It has, however, been hypothesized that the two proteins may act synergistically *in vivo*.<sup>11,12</sup>

The most frequently used test in the diagnosis of *C. difficile* colitis is determination of toxin B in patients' stool by cell culture with neutralization of the toxin by specific antiserum. Although this assay is extremely sensitive, and the presence of toxin B has a positive correlation to 90-100% of patients with severe disease,<sup>13</sup> it requires cell culture capability and up to 72 hours of incubation. In addition, the cytotoxin assay is not standardized and the procedures and cell lines used vary considerably.<sup>9</sup>

### BIOLOGICAL PRINCIPLES

Premier Toxins A&B is an enzyme immunoassay for the direct detection of *Clostridium difficile* toxin A and toxin B in stool samples. Breakaway microwells are coated with toxin specific monoclonal and polyclonal antibodies. Diluted patient specimens and HRP-conjugated anti-toxin A and B polyclonal antibodies are added to microwells. If either toxin is present in the diluted patient samples, HRP-conjugated toxin polyclonal antibodies (specific for both toxins) complexes are formed which remain in the microwells after washing. After a final washing step, a substrate/chromagen (peroxide and tetramethylbenzidine) is added to the wells. Any bound conjugate converts the substrate/chromagen to a blue color. Addition of acid (Stop Solution I) converts the blue to a yellow color.

### REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. Premier Toxins A&B Antibody-coated Microwells – Breakaway plastic wells coated with mouse monoclonal anti-toxin A antibody and polyclonal goat anti-toxin B antibody.
2. Premier Toxins A&B Positive Control – Inactivated toxin A and B in buffered protein solution containing gentamicin and thimerosal (0.02%) as preservatives.
3. Premier Toxins A&B Sample Diluent/Negative Control – Buffered protein solution containing gentamicin and thimerosal (0.02%) as preservatives.
4. Premier 20X Wash Buffer II – Concentrated wash buffer containing thimerosal (0.2%) as a preservative.
5. Premier Toxins A&B Enzyme Conjugate - Polyclonal goat anti-toxin A and anti-toxin B antibodies conjugated to horseradish peroxidase in buffered protein solution containing gentamicin and thimerosal (0.02%) as preservatives.
6. Premier Substrate I - Buffered solution containing peroxide and tetramethyl benzidine.
7. Premier Stop Solution I – 1M Phosphoric acid. CAUTION: Avoid contact with skin. Flush with water if contact occurs.
8. Sample Transfer pipettes
9. Microwell plate sealers

### MATERIALS NOT PROVIDED

1. Test tubes for dilution of sample
2. Distilled or deionized water
3. OPTIONAL: EIA microwell reader capable of reading absorbance at 450 nm or 450/630 nm\*
4. Squirt bottle
5. Pipettes and graduated cylinder for making 1X Wash Buffer
6. Wooden applicator sticks
7. Absorbent paper
8. Incubator 37 ± 2 C
9. Timer
10. Vortex mixer
11. Waste container with disinfectant (ie, 10% solution of household bleach) and/or autoclavable biohazard bags
12. Disposable gloves
13. OPTIONAL: Centrifuge
14. OPTIONAL: Stat Fax™ - 2200 Incubator/Shaker for use in 20 minute incubation procedure\*
15. OPTIONAL: Semiautomated microplate washer\*

\* It is the operator's responsibility to validate the semiautomated microplate washer, incubators and readers prior to their use with this product.

Stat Fax™ is a trademark of Awareness Technology, Inc.

### PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Kit reagents should be warmed to room temperature (21-27 C) and gently mixed before use.
3. Do not mouth pipet samples or reagents. Avoid contact with skin or mucous membranes.
4. Do not smoke, drink or eat in areas where specimens or kit reagents are handled.
5. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
6. Do not interchange the microwells, Enzyme Conjugate, Substrate I, and Positive Control reagents between lots. (The Sample Diluent/Negative Control, 20X Wash Buffer II and Stop Solution I are interchangeable provided the reagents are within their assigned expiration dates.)
7. Do not use kit beyond the expiration date given on the kit label.
8. Stop Solution I contains 1M phosphoric acid. If contact with skin or mucous membranes occurs, flush with water immediately.
9. Patient specimens may contain infectious agents and should be handled at Biosafety Level 2 as recommended in the CDC/NIH manual "Biosafety in Microbiology and Biomedical Laboratories".
10. Avoid splashing or forming of aerosols when handling, diluting or transferring specimens.
11. Avoid microbial contamination of reagents or incorrect results may occur.
12. Transfer pipettes provided must be used for specimen preparation and transfer. Use 1 per specimen. Cross contamination of samples or reagents may cause incorrect results.
13. All stool samples must be mixed thoroughly, regardless of consistency, to insure a representative sample prior to pipetting.
14. Reagent concentration, incubation times and temperatures have been optimized for sensitivity and specificity. Best results are obtained by adhering to these specifications.
15. Positive Control contains inactivated toxin A and toxin B. However, it should be handled as a potential biohazard.
16. Hold all vials vertically to insure proper drop size and delivery.

17. Do not allow the tips of any vial to touch the microwells.
18. Replace colored caps on correct vials.
19. All reagents are provided already diluted to the proper concentration (except the Premier 20X Wash Buffer II). Do not dilute further.
20. Any deviation below or above set incubation times may affect sensitivity and specificity and should be avoided.
21. Do not use microwells if the foil pouch appears damaged (ie, has holes, punctures).

### HAZARD and PRECAUTIONARY STATEMENTS

Refer to the SDS, available at [www.meridianbioscience.com](http://www.meridianbioscience.com) for Hazards and Precautionary Statements.

### SHELF LIFE AND STORAGE

1. Store all kit reagents at 2-8 C. Return kit to the refrigerator promptly after use.
2. Keep microwells in their pouch until the pouch reaches room temperature to avoid condensation. Return all unused microwells and the desiccant to the ziplock pouch and seal tightly.
3. Diluted Wash Buffer II may be kept at room temperature and used for up to three months.


### REAGENT PREPARATION

1. Bring the entire kit, including microwell pouch, to room temperature (21-27 C) before use. Return to 2-8 C immediately after use.
2. Prepare decontamination vessel for discarding reagents and materials.
3. Do not allow microwells to dry out between steps.
4. Reproducibility in any EIA is largely dependent on the consistency and thoroughness with which the microwells are washed. Carefully follow the recommended wash procedure as outlined in the EIA test procedure. An automated washer may be used.
5. Prepare 1X Wash Buffer as needed. For example: 5.0 mL of 20X Wash Buffer II + 95.0 mL of distilled or deionized wash is sufficient to wash 1 strip. Place in clean squirt bottle. The 1X Wash Buffer can be stored at 21-27 C for up to 3 months.
6. One Positive and 1 Negative Control well must be included with each run of specimens. Use the Positive Control as provided. **DO NOT DILUTE.**
7. Use EIA plate sealers to cover the assay during incubation steps. Cut to size, then remove paper backing before use.

### SPECIMEN COLLECTION AND PREPARATION (1:5 Sample Dilution)

Collect stool specimens into a clean, airtight container with no preservative. All stool specimens should be stored at 2-8 C and tested as soon as possible. Ideally, stool specimens should be tested within 24 hours but specimens may be held at 2-8 C for up to 72 hours prior to testing. If specimens cannot be tested within 72 hours, they should be frozen immediately upon receipt at -20 C or colder. A single freeze thaw cycle should not affect results.<sup>13</sup> Repeated freezing and thawing of samples should be avoided. Use only the Sample Diluent provided with this kit for diluting specimens. Specimen can be stored diluted in Sample Diluent in a sealed tube for up to 8 hours at 2-8 C.

1. Add 200 µL of Sample Diluent to a clean test tube with dropper assembly (or equivalent).
2. Mix stool as thoroughly as possible prior to pipetting.
  - a. For liquid or semi-solid stool: Using the disposable transfer pipettes and measuring to the 1st calibration mark (50 µL), transfer 50 µL of stool into the dilution tube containing Sample Diluent. Rinse the transfer pipette by drawing the stool suspension in and out of the pipette several times. Remove pipette and vortex the suspension thoroughly (15 seconds), and replace pipette in tube.
  - b. Solid Stools: Using a wooden applicator stick, transfer a small 3-4 mm diameter portion of thoroughly mixed stool into Sample Diluent. Emulsify stool using the wooden applicator stick then vortex 15 seconds. Place a transfer pipette in the tube. Promptly proceed to assay.
3. Stool specimens may be centrifuged after dilutions. Centrifuge at approximately 2750 x G for 5 minutes or until solid matter separates from liquid. Proceed with the assay after recovering the supernate.

### TEST PROCEDURE

**NOTE:** With large numbers of specimens, repetitive or multichannel pipettes may be used for dispensing the reagents.

1. After the pouch has reached 21-27 C, break off the required number of microwells (1 well for each specimen plus 1 positive and 1 negative control well per batch). Place the microwells in the microwell strip holder and record the location of all wells. Unused microwells must be resealed in the pouch immediately (see SHELF LIFE AND STORAGE).
  2. Using the original transfer pipette, draw up diluted stool to the 100 µL calibration point (2nd mark from tip of the pipette) and add to the appropriate well (place pipette tip halfway into well and let sample slowly run down side of well).
  3. Add 2 free-falling drops of Positive Control to the appropriate well. Add 100 µL (2nd mark on transfer pipette) of Negative Control (Sample Diluent) to the appropriate well.
  4. Add 1 free-falling drop of Enzyme Conjugate (50 µL) to all wells. Mix wells by firmly shaking/swirling the plate for 30 seconds.
  5. Cut plate sealer to size and press firmly onto top of microwells to seal. Incubate the plate for 50 minutes at 35-39 C. Alternatively, laboratories equipped with a heated plate shaker (Stat Fax™-2200) can incubate and rotate the plate for 20 minutes at 37 C at 1000 rpm (setting #5).
  6. Carefully, remove the plate sealer and wash wells (see REAGENT PREPARATION):
    - a. Dump plate contents firmly into a biohazard receptacle.
    - b. Bang the inverted plate on a clean stack of paper towels
    - c. Fill all wells with 1X Wash Buffer directing stream of buffer to the sides of the wells to avoid foaming. Immediately proceed to Step 6.d.

**OPTIONAL: Alternatively, a semi-automated plate washer may be used. Select the appropriate setting according to manufacture instructions. After automated washing is complete, proceed to Step 7.**

    - d. Repeat wash cycle (dump, bang on fresh towels, fill) 4 to 6 times (for a total of 5-7 washes). After the last fill, dump and bang plates on fresh towels hard enough to remove as much excess wash buffer as possible but do not allow wells to completely dry at any time.
  7. Clean the underside of all wells with a lint free tissue.
  8. Add 2 free-falling drops of Substrate I (100 µL) to each well.
  9. Firmly shake/swirl the plate for 10-15 seconds then incubate for 10 minutes at 21-27 C.
  10. Add 2 free-falling drops of Stop Solution I (100 µL) to all wells. Shake/swirl the plate firmly for 30 seconds to assure complete mixing. After addition of Stop Solution I, wait 2 minutes before reading (Step 11).
- NOTE:** initial color of positive reaction is blue, which changes to yellow upon addition of Stop Solution I.
11. Observe reactions:
    - a. Visual Determination – Read within 15 minutes after adding Stop Solution I.
    - b. Spectrophotometric Determination – Zero EIA reader on air. Wipe the underside of wells with a lint free tissue. Read absorbance at 450 nm or 450/630 nm within 30 minutes of adding Stop Solution I.

### INTERPRETATION OF RESULTS

1. **Visual Reading**  
Negative = colorless to faint (barely visible) yellow  
Positive = definite yellow color
2. **Spectrophotometric Single Wavelength (450 nm)**  
Negative = OD<sub>450</sub> < 0.150  
Positive = OD<sub>450</sub> ≥ 0.150
3. **Spectrophotometric Dual Wavelength (450/630 nm)**  
Negative = OD<sub>450/630</sub> < 0.100  
Positive = OD<sub>450/630</sub> ≥ 0.100

Extremely strong positive reaction may yield a purple precipitate within a few minutes of stopping the reaction.

A positive result indicates the presence of *C. difficile* toxin A and/or B. A negative result indicates the absence of toxins A and B, or that the level of toxin is below that which can be detected by the assay. The magnitude of the OD, above the cut-off, is not indicative of the severity or extent of the *C. difficile* infection.

#### QUALITY CONTROL

**This test should be performed per applicable local, state, or federal regulations or accrediting agencies.**

- The Positive and Negative Controls must be used with each batch of specimens to provide quality assurance of the reagents and the procedure. It is suggested that the results of each quality control check be recorded in an appropriate log book to maintain high quality testing procedures and compliance with regulatory agencies.
- The Negative Control should read < 0.150 at 450 nm and < 0.100 at 450/630 nm but greater than 0.00. If control is < 0.00, re-blank the plate reader to air and reread the plate.
- The Negative Control should be colorless to faint (barely visible) yellow when read visually.
- The Positive Control should read < 2.999 but > 0.600 at either 450 nm or 450/630 nm. The Positive Control should have a definite yellow color when read visually.
- If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
- Any positive well without visible color should be repositioned, wiped on the underside of the well and reread.
- At the time of each use, kit components should be visually examined for obvious signs of microbial contamination, freezing, or leakage.

#### EXPECTED VALUES

The frequency of antibiotic-associated diarrhea caused by *C. difficile* is dependent on several factors including: patient population, type of institution and epidemiology. The reported incidence of *C. difficile*-associated disease in patients suspected of having antibiotic-associated diarrhea is 15-25%<sup>1</sup> although different facilities may find positivity rates outside this range.

#### LIMITATIONS OF THE PROCEDURE

- Premier Toxins A&B detects the presence of *Clostridium difficile* toxin A and toxin B in stool. Failure to detect toxin A or toxin B in stool from patients with suspected *C. difficile*-associated disease may not preclude actual disease but may be due to such factors as improper collection, handling and storage of the specimen or toxin concentrations in stool below the detection limit of this kit. The Premier Toxins A&B kit will detect toxin A at levels  $\geq 1.4$  ng/mL stool and toxin B at levels  $\geq 2.4$  ng/mL stool. Reactivity of positive samples as determined by Premier Toxins A&B may decrease with time due to degradation of the toxins. As with any diagnostic test, the results of Premier Toxins A&B should be interpreted with respect to patient history and other clinical and laboratory findings.
- The performance of this assay has not been evaluated in a pediatric population.<sup>14</sup>
- Certain strains of *C. sordellii* produce toxins which are immunologically crossreactive with *C. difficile* toxins A and B. *C. sordellii* has not been isolated, however, from stool obtained from patients with antibiotic-associated diarrhea or pseudomembranous colitis while *C. difficile* was found to be present.<sup>15,16</sup>
- While relatively stable at 2-8 C, *C. difficile* toxins – particularly at low concentrations – easily degrade at room temperatures. The rate at which the toxins degrade differs from patient sample to patient sample and therefore the rate cannot be predicted. For this reason, best practice requires that samples be refrigerated or frozen within two hours of collection and the samples tested within the timeframes recommended in this insert. Do not accept samples that have not been collected, handled or transported properly.

#### SPECIFIC PERFORMANCE CHARACTERISTICS

The performance of Premier Toxins A&B was evaluated in a clinical study performed at two sites in the United States. The kit was compared to the cellular cytotoxicity assay and discrepant results were resolved using toxigenic culture and a competitor's EIA for toxins A and B.

Premier Toxins A&B Results	Cytotoxin Result: Site 1		Cytotoxin Result: Site 2		Cytotoxin Result: All Sites	
	Pos	Neg	Pos	Neg	Pos	Neg
Pos	55	7	35	6	90	13
Neg	3	257	2	208	5	465
Performance Statistic	Value	95% CI	Value	95% CI	Value	95% CI
Sensitivity	94.8%	85.6-98.9%	94.6%	81.8-99.3%	94.7%	88.1-98.3%
Specificity	97.3%	94.6-98.9%	97.2%	94.0-99.0%	97.3%	95.4-98.5%
Positive Predictive Value	88.7%	78.1-95.3%	85.4%	70.8-94.4%	87.4%	81.0-93.8%
Negative Predictive Value	98.8%	96.7-99.8%	99.0%	96.6-99.9%	98.9%	97.5-99.7%
Correlation	96.9%	94.4-98.5%	96.8%	93.8-98.6%	96.9%	95.1-98.1%

Performance of the Premier Toxins A&B EIA was equivalent at each of the clinical trial sites. Overall sensitivity and specificity compared to the reference cytotoxin method were 94.7% and 97.3% respectively. Examination of the 13 false-positive specimens revealed that one was positive by toxigenic culture (ie, yielded a *C. difficile* isolate producing both toxins A and B). Four other specimens were positive by a competitor's toxin A+B EIA. Thus, 5/13 (38%) of the false-positive specimens had additional independent evidence supporting the presence of toxigenic *C. difficile*. Toxigenic culture was negative on 3/5 of the Premier Toxins A&B negative, cytotoxin positive specimens. In addition, all five were negative by a competitor's toxin A+B EIA.

A competitor's A+B EIA (Wampole™ Laboratories) was also run on at each trial site. Results are compared to cytotoxin and Premier Toxins A&B below:

Premier Toxins A&B Results	Wampole A/B EIA		Wampole Toxins A/B Results	Cytotoxin Result	
	Pos	Neg		Pos	Neg
Pos	92	11	88	13	465
Neg	9	461	7	465	95% CI
Relative Performance	Value	95% CI	Value	95% CI	
Co-Positivity	91.1%	93.8-95.8%	92.6%	84.5-97.0%	
Co-Negativity	97.7%	95.6-98.7%	97.3%	95.4-98.5%	
Agreement	96.5%	94.5-97.7%	98.5%	97.0-99.4%	
			96.5%	94.7-97.9%	

Agreement between the two EIA's for toxins A&B was 96.5%. Two of the 11 specimens positive by Premier Toxins A&B and negative by the competitor's EIA were positive by cytotoxin and another specimen was positive by toxigenic culture. None of the 9 specimens positive by the competitor's EIA and negative by Premier Toxins A&B were positive by cytotoxin.

Wampole™ is a trademark of Wampole Laboratories

#### REPRODUCIBILITY

The reproducibility of the Premier Toxins A&B test was established using negative, low, medium and strong positive samples and controls tested in triplicate in three different batches at three sites. Inter- and intra-assay variances were calculated and are given below.

Source of Variance	Positive Control	Negative Control	High Positive	Medium Positive	Low Positive	Negative
Mean Absorbance	2.010	0.013	2.250	1.146	0.280	0.009
Within Run CV	4.1%	24.5%	7.3%	6.9%	15.9%	28.9%
Between Run CV	7.0%	16.2%	6.2%	13.9%	14.6%	31.7%

#### CROSSREACTIVITY

The specificity of Premier Toxins A&B was tested by following the following bacterial or viral strains. Positive and negative stools were spiked with  $\geq 1 \times 10^8$  bacteria/mL and tested by Premier Toxins A&B. The only non-*C. difficile* microorganisms reactive with the Premier Toxins A&B were two strains of *Clostridium sordellii*, ATCC 9714 and VPI 9048. Both of these strains produce toxin A and B homologues HT and LT respectively. All other organisms were found to be negative when spiked into the negative stool. In addition, they did not interfere with the positive specimen.<sup>13</sup>

#### Microorganism or virus (# strains tested)

Adenovirus 40 & 41 (2)	<i>Aeromonas hydrophila</i> (1)	<i>Bacillus cereus</i> (1)
<i>Bacillus subtilis</i> (1)	<i>Bacteroides fragilis</i> (1)	<i>Campylobacter coli</i> (1)
<i>Campylobacter jejuni</i> (1)	<i>Candida albicans</i> (1)	<i>Clostridium butyricum</i> (1)
<i>Clostridium difficile</i> non-tox.(7)	<i>Clostridium haemolyticum</i> (1)	<i>Clostridium histolyticum</i> (1)
<i>Clostridium novyi</i> (1)	<i>Clostridium perfringens</i> (1)	<i>Clostridium septicum</i> (1)
<i>Clostridium sporogenes</i> (1)	<i>Clostridium tetani</i> (1)	<i>Enterobacter aerogenes</i> (1)
<i>Enterobacter cloacae</i> (1)	<i>Escherichia coli</i> (3)	<i>Helicobacter pylori</i> (1)
<i>Klebsiella pneumoniae</i> (1)	<i>Peptostreptococcus anaerobius</i> (1)	<i>Porphyromonas asaccharolytica</i> (1)
<i>Proteus vulgaris</i> (1)	<i>Pseudomonas aeruginosa</i> (1)	Rotavirus (1)
<i>Salmonella typhimurium</i> (1)	<i>Serratia liquefaciens</i> (1)	<i>Shigella dysenteriae</i> (1)
<i>Shigella flexneri</i> (1)	<i>Shigella sonnei</i> (1)	<i>Staphylococcus aureus</i> (1)
<i>Staphylococcus aureus</i> (Cowan I) (1)	<i>Staphylococcus epidermidis</i> (1)	<i>Streptococcus faecalis</i> (1)
<i>Vibrio cholerae</i> (1)	<i>Vibrio parahaemolyticus</i> (1)	<i>Yersinia enterocolitica</i> (1)

#### ASSAY REACTIVITY

The Premier Toxins A&B test was also evaluated using several reference strains of *C. difficile*. Strains were grown in BHI broth for 48 hours and tested for reactivity with the Premier Toxins A&B test. Results summarized below, indicated that the test correctly identified all toxigenic strains of *C. difficile*, even those producing only toxin B. No cross-reactivity was observed with non-toxigenic strains of *C. difficile*.

<i>C. difficile</i> Type	# of Premier Toxins A&B Pos / Total (% Correct)
A+ / B+	25/25 (100%)
A- / B+	3/3 (100%)
A- / B-	0/14 (100%)

#### TESTS FOR INTERFERING SUBSTANCES

The results of Premier Toxins A&B testing are not affected by blood, barium sulfate, metronidazole or vancomycin present in stool specimens.<sup>13</sup>

The following substances were found to have no effect on results when present in stool samples at the concentrations indicated: Barium Sulfate (stool sample diluted 1:2 with 10% Barium Sulfate), Metronidazole (2.8 µg/well), Vancomycine (2.8 µg/well), and Whole blood (50%).

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## REFERENCES

- Bartlett JG. *Clostridium difficile*: clinical considerations. Rev Infect Dis. 1990;12 S243-S251.
- Bartlett JG, TW Chang, M Gurwith, SL Gorbach and AB Onderdonk. Antibiotic associated pseudomembranous colitis due to toxin producing clostridia. N Engl J Med 1978;298: 531-534.
- George WL, VL Sutter, EJC Goldstein, SL Ludwig and SM Finegold. Aetiology of anti-microbial agent associated colitis. Lancet. 1978;1: 802-803.
- Rolle RD. Asymptomatic intestinal colonization by *Clostridium difficile*. In: RD Rolfe and SM Finegold (eds.). *Clostridium difficile*: its role in intestinal disease. Academic Press, Inc. New York. 1988.
- McFarland LV, ME Mulligan, RYY Kwok and WE Stamm. Nosocomial acquisition of *Clostridium difficile* infection. New Engl J Med. 1989;320:204-210.
- Johnson S, A Adelmann, CR Clabots, LR Peterson, DN Gerding. Recurrences of *Clostridium difficile* diarrhea not caused by the original infecting organism. J Infect Dis. 1989;159 (2):340-343.
- Sullivan NM, S Pellett, TD Wilkins. Purification and characterization of Toxins A and B of *Clostridium difficile*. Infect Immun. 1982;35(3): 1032-1040.
- Lyerly DM, CJ Phelps, J Toth and TD Wilkins. Characterization of toxins A and B of *Clostridium difficile* with monoclonal antibodies. Infect Immun. 1986;54(1):70-76.
- Lyerly DM, CJ Phelps, and TM Wilkins. Monoclonal and specific polyclonal antibodies for immunoassay of *Clostridium difficile* Toxin A. J Clin Microbiol. 1985;21(1): 12-14.
- Chang TW, M Lauermann and JG Bartlett. Cytotoxicity assay in antibiotic-associated colitis. J Infect Dis. 1980;140:765-770.
- Lyerly DM, KE Saum, DK McDonald, and TD Wilkins. Effects of *Clostridium difficile* toxins given intragastrically to animals. Infect Immun. 1985;47:349-352.
- Meador J, and RK Tweten. Purification and characterization of toxin B from *Clostridium difficile*. Infect Immun. 1988;56(7):1708-1714.
- Data on file, Meridian Bioscience, Inc.
- Cooperstock, M. *Clostridium difficile* in infants and children. In: RD Rolfe and SM Finegold (eds.). *Clostridium difficile*: its role in intestinal disease. Academic Press, Inc. New York. 1988.
- Lyerly DM and TD Wilkins. Purification and properties of toxins A and B of *Clostridium difficile*. In: RD Rolfe and SM Finegold (eds.). *Clostridium difficile*: its role in intestinal disease. Academic Press, Inc. New York. 1988.
- Willey SH and JG Bartlett. Culture for *Clostridium difficile* in stools containing a cytotoxin neutralized by *Clostridium sordellii* antitoxin. J Clin Micro. 1979;10: 860-864.

SN11161

REV. 01/15

 Manufactured By	Meridian Bioscience, Inc. USA/Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 Telephone: 513.271.3700 Orders/Customer Service: 800.543.1980 Technical Support Center: 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124
 Authorized Representative	Meridian Bioscience Europe S. r. l. Via dell'Industria, 7 20020 Villa Cortese, Milano ITALY Tel: +39 0331 43 36 36 Fax: +39 0331 43 36 16 Email: info@meridianbioscience.eu WEB: www.meridianbioscience.eu

Meridian Bioscience Europe s.a./n.v.  
 2 Avenue du Japon - 1420 Braine l'Alleud  
 BELGIUM  
 Tel: +32 (0) 67 89 59 59  
 Fax: +32 (0) 67 89 59 58  
 Email: info.bn@meridianbioscience.eu

Meridian Bioscience Europe France  
 34 rue de Ponthieu - 75008 Paris  
 FRANCE  
 Tel: +33 (0) 1 42 56 04 40  
 Fax: +33 (0) 9 70 06 62 10  
 Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.  
 Postbus 301 - 5460 AH Veghel  
 NETHERLANDS  
 Tel: +31 (0) 411 62 11 66  
 Fax: +31 (0) 411 62 48 41  
 Email: info.bn@meridianbioscience.eu

## INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de símbolos, Erläuterung der graphischen symbole)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	<b>CONTROL +</b>	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
<b>LOT</b>	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	<b>CONTROL -</b>	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
<b>IVD</b>	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo medico para diagnóstico in vitro / In-Vitro-Diagnostikum	<b>EC REP</b>	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In-Vitro-Diagnostika	<b>SMP   PREP   DIL   SPE</b>	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenzubereitung, in dem sich Probenverdünnungspuffer befindet
<b>REF</b>	Catalogue number / Numero di catalogo / Référence ou catalogue / Numero de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	<b>BUF RXN</b>	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	<b>SOLN   STOP</b>	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	<b>CONJ   ENZ</b>	Enzyme Conjugate / Conjugato enzimatico / Conjugué enzymatique / Conjugado enzimático / Enzymkonjugat
<b>SN</b>	Serial number / Numero di serie / Numero de série / Numero de serie / Seriennummer	<b>CONTROL</b>	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
<b>TEST</b>	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	<b>REAG</b>	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	<b>BUF   WASH</b>	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
<b>BUF</b>	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
<b>CONJ</b>	Conjugate / Conjugato / Conjugué / Conjugado / Konjugat	<b>DIL   SPE</b>	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
<b>SUBS</b>	Substrate / Substrato / Substrat / Substrato / Substrat	<b>BUF   WASH   20X</b>	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
		<b>DET   REAG</b>	Detection Reagent / Reagente Diretto / Réactif de Détection / Reactivo de Detección / Nachweis Reagenz

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.