# PREMIER® CAMPY

## Enzyme Immunoassay for the Detection of *Campylobacter* Antigens in Stool Specimens

REF

618906

IVD

In vitro diagnostic medical device

### INTENDED USE

Premier CAMPY enzyme immunoassay (EIA) is an in vitro qualitative procedure for the detection of specific *Campylobacter* antigens in stool samples from patients with signs and symptoms of gastroenteritis. Premier CAMPY detects *C. jejuni* and *C. coli* in human stool that may be either unpreserved or preserved in Cary-Blair-based transport media. Test results are to be used in conjunction with information obtained from the patient's clinical evaluation and other diagnostic procedures.

Premier CAMPY is not intended for point-of-care use. The device is intended for use in hospital, reference, regional, private or state laboratory settings.

### SUMMARY AND EXPLANATION OF THE TEST

*Campylobacter* is a gram-negative, microaerophilic bacterium. Virtually all human illness is caused by one or two species. These two species are *C. coli* and *C. jejuni*. The Centers for Disease Control reports that 99% of illness is caused by *C. jejuni* while other studies have shown that globally, more than 90% of *Campylobacter* infections are caused by *C. jejuni*, followed by *C. coli* with 5-10%.<sup>2</sup> The disease caused by the genus *Campylobacter* is referred to as *Campylobacteriosis*. Most people with *Campylobacteriosis* get diarrhea, cramping, abdominal pain, and fever within two to five days after exposure to the organism. The diarrhea may be bloody and can be accompanied by nausea and vomiting. The illness typically lasts one week. Some persons who are infected with *Campylobacter* do not have any symptoms at all. In persons with compromised immune systems, *Campylobacter* occasionally spreads to the bloodstream and causes a serious life-threatening infection.

*Campylobacter* is one of the most common bacterial causes of diarrheal illness in the United States. Virtually all cases occur as isolated, sporadic events, not as a part of large outbreaks. Active surveillance through FoodNet indicates about 15 cases are diagnosed each year for each 100,000 persons in the population.<sup>1, 2</sup> Many more cases go undiagnosed or unreported, and *Campylobacteriosis* is estimated to affect over one million persons every year, or 0.5% of the general population. *Campylobacteriosis* occurs much more frequently in the summer months than in the winter. The organism is isolated from infants and young adults more frequently than from other age groups and from males more frequently than females. Although *Campylobacter* doesn't commonly cause death, it has been estimated that approximately 100 persons with *Campylobacter* infections may die each year.<sup>2</sup>

### BIOLOGICAL PRINCIPLES

Premier CAMPY is an enzyme immunoassay for the direct detection of *Campylobacter* antigens in human stool samples. Breakaway microwells are coated with *Campylobacter*-specific monoclonal antibodies. Diluted patient specimen is added to the microwells and incubated. Upon completion of the incubation, a wash step is performed to remove unbound material and a Horseradish Peroxidase (HRP)-Anti-*Campylobacter* conjugate is added to the washed microwells. If *Campylobacter* antigens are present, an antibody-enzyme complex is formed. A second wash step is performed to remove unbound materials and a chromagen substrate is added to the microwells. A blue color develops in the presence of bound enzyme. Premier Stop Solution I is added, changing the initial blue reaction to yellow. Test results are interpreted visually or spectrophotometrically.

### REAGENTS AND MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

- Premier CAMPY Microwells:** Monoclonal antibody-coated microwells. The antibodies are specific to *C. jejuni* and *C. coli*.
- Premier CAMPY Enzyme Conjugate:** HRP-conjugated monoclonal antibodies specific to *C. jejuni* and *C. coli* in a buffered protein solution containing 0.1% ProClin® and 0.03% gentamicin as preservatives.
- Premier 20X Wash Buffer III:** Concentrated wash buffer containing 0.1% ProClin® and 0.3% gentamicin as preservatives.
- Premier Substrate I:** Buffered solution containing urea peroxide and tetramethylbenzidine at pH 5.0.
- Premier Stop Solution I:** 1M Phosphoric acid.
- Premier CAMPY Sample Diluent/Negative Control:** Buffered protein solution containing 0.1% ProClin® and 0.03% gentamicin as preservatives.
- Premier CAMPY Positive Control:** Inactivated *C. jejuni* in a buffered protein solution containing 0.094% sodium azide and 0.03% gentamicin as preservatives.
- Transfer pipettes
- Microwell strip holder
- Microwell plate sealers

### MATERIALS NOT PROVIDED

- Disposable latex gloves
- Test tubes (eg, 10 x 75 mm or 12 x 75 mm) for dilution of sample
- Distilled or deionized water
- Squirt bottle
- Graduated cylinder for making 1X Wash Buffer
- Absorbent paper
- Wooden applicator sticks
- Waste container with disinfectant and/or autoclavable biohazard bags
- Vortex mixer
- Interval timer
- EIA microplate reader capable of reading absorbance at 450 nm or 450/630 nm\* (optional)
- StatFax™-2200 Incubator/Shaker (StatFax™ is a trademark of Awareness Technology, Inc.)\* (optional)
- Semiautomated microplate washer\* (optional)

\* **Note:** It is the operator's responsibility to validate StatFax™, semiautomated plate washers and readers prior to their use with this product.

### PRECAUTIONS

- All reagents are for in vitro diagnostic use only.
- Directions should be read and followed carefully.
- Do not interchange the Microwells, Enzyme Conjugate, Premier Substrate I or Positive Control reagents between lots. (The Sample Diluent/Negative Control, Premier 20X Wash Buffer III and Premier Stop Solution I are interchangeable provided the reagents are within their assigned expiration dates when used.)
- Do not use kit components beyond labeled expiration date.
- Do not use vials that lack a label, a lot number, or an expiration date.
- Do not use any reagent if it is discolored or turbid. Discoloration or turbidity may be a sign of microbial contamination.
- Allow reagents to warm to 21–27 C before use.
- All reagents should be mixed gently before use.
- Hold reagent vials vertically at suitable distance above the well to insure proper drop size and delivery.
- Replace colored caps on correct vials.
- Do not reuse microwells.
- Unused microwells must be placed back inside resealable pouch. It is important to protect strips from moisture.
- The transfer pipettes provided with this kit must be used for specimen preparation and transfer. Use one per specimen.
- Avoid splashing when dispensing diluted stool into microwells by placing the transfer pipette tip about halfway down the well and dispensing slowly down the side of well.

- Microwell washing is to be performed precisely as directed in the assay procedure. **Inadequate washing may be the cause of elevated background in any EIA protocol.**
- All reagents except the Premier 20X Wash Buffer III are provided already diluted to the proper concentration.
- Any deviation below or above set incubation times may affect sensitivity and specificity and should be avoided.
- Stool must be mixed thoroughly (regardless of consistency) to ensure a representative sample prior to pipetting.

### WARNINGS

- Patient specimens may contain infectious agents and should be handled and disposed of as potential biohazards.
- Dispose of used wash buffer and all test materials in an appropriate container. Treat waste as a potential biohazard.
- The Positive Control reagent contains inactivated *C. jejuni*. It should be handled, however, as a potential biohazard.
- Avoid skin contact with Premier Stop Solution I. Flush with water immediately if contact occurs.

### HAZARD and PRECAUTIONARY STATEMENTS


Refer to the SDS, available at [www.meridianbioscience.com](http://www.meridianbioscience.com) for Hazards and Precautionary Statements.

### SHELF LIFE and STORAGE

The kit expiration date is indicated on the kit label. Store the kit at 2–8 C and return the kit promptly to the refrigerator after each use.

### PROCEDURAL NOTES

The Premier CAMPY transfer pipette is diagrammed below.

### REAGENT PREPARATION

- Bring the entire kit, including microwell pouch, to 21–27 C before use.
- Prepare 1X Wash Buffer as needed. For example: 10.0 mL of Premier 20X Wash Buffer III + 190.0 mL distilled or deionized water is sufficient to wash one strip. 1X Wash Buffer can be stored at 2–27 C for up to one month. Bring to 21–27 C prior to use.

### SPECIMEN COLLECTION

**Human stool samples, unpreserved:** Samples should be received in an airtight transport container and stored at 2–8 C prior to testing. Samples should be tested as soon as possible, but may be held up to 96 hours at 2–8 C. (See SPECIMEN PREPARATION section for instructions on diluting samples.) Samples that will not be tested within 96 hours should be frozen immediately upon receipt and stored at ≤ –20 C until tested. Specimens may be frozen and thawed twice.

**Human stool samples, preserved in Cary-Blair-based media:** Samples should be stored at 2–8 C prior to testing. Samples should be tested as soon as possible, but may be held up to 96 hours at 2–8 C. (See SPECIMEN PREPARATION section for instructions on diluting samples.) Samples that will not be tested within 96 hours should be frozen immediately upon receipt and stored at ≤ –20 C until tested. Specimens may be frozen and thawed twice.

### SPECIMEN PREPARATION

Mix stool as thoroughly as possible prior to pipetting.

- Human stool samples, unpreserved:**
  - With the dropper assembly (or equivalent), add 200 µL of Sample Diluent/Negative Control to a small test tube.
  - Formed/Solid stools:**
    - Use a wooden applicator stick (or equivalent) to add a small portion (3–4 mm diameter) of thoroughly mixed stool to the Sample Diluent/Negative Control test tube.
    - Emulsify the stool using the wooden applicator stick.
    - Vortex the emulsion for a minimum of 15 seconds.
  - Liquid/Semi-Solid stools:**
    - Using a transfer pipette provided with the kit, add 50 µL (first mark from pipette tip) of thoroughly mixed stool to the Sample Diluent/Negative Control test tube.
    - Vortex the mixture for a minimum of 15 seconds.
    - Save the transfer pipette in the test tube for future use.
- Human stool samples, preserved in Cary-Blair-based media:**
  - With the dropper assembly (or equivalent), add 200 µL of Sample Diluent/Negative Control to a small test tube.
  - Using a transfer pipette provided with the kit, add 200 µL (third mark from pipette tip) of thoroughly mixed stool to the Sample Diluent/Negative Control test tube.
  - Vortex the mixture for a minimum of 15 seconds.
  - Save the transfer pipette in the test tube for future use.

### TEST PROCEDURE

- After the pouch has reached temperature (21–27 C), break off the required number of microwells (1 well for each specimen, plus 1 positive and 1 negative control well per batch). Place the microwells in the microwell strip holder and record the location of all wells. Unused microwells must be resealed in the pouch immediately.
- Using the specimen transfer pipette, add 100 µL of diluted stool (second mark from the tip of the pipette) to the appropriate well. (Place the pipette tip halfway into well and let the sample slowly run down side of well.)
- Add 2 free falling drops of Positive Control. With a transfer pipette, add 100 µL (second mark from the tip of the pipette) of Sample Diluent/Negative Control to the appropriate wells. Cut the plate sealer to size and press it firmly onto the top of the microwells to seal.
- Incubate the plate for 1 hour at 21–27 C. Alternatively, laboratories equipped with a shaker/incubator (StatFax™-2200) can incubate and rotate the plate for 30 minutes at 24 C and at Setting 5 (1000 rpm).
- Carefully remove the plate sealer and wash the wells:
  - Manual method:**
    - Dump the plate contents firmly into a biohazard receptacle.
    - Invert and bang the plate on a clean stack of paper towels.
    - Fill all wells with 1X Wash Buffer, directing a stream of buffer to the sides of the wells to avoid foaming.
    - Repeat the wash cycle (dump, bang on fresh towels, fill) 4 more times for a total of 5 wash cycles. After the last fill, dump and bang plate on fresh towels hard enough to remove as much excess wash buffer as possible, but do not allow wells to completely dry at any time.
  - Semi-automated wash method:**
    - Dump the plate contents firmly into a biohazard receptacle.
    - Invert and bang the plate on a clean stack of paper towels, then place empty plate on washer device.
    - Fill the wells to the top (approx. 300–350 µL/well) with 1X Wash Buffer and then aspirate. The washer manifold should be adjusted to ensure no foaming occurs during the filling of the wells and that the wells are thoroughly aspirated after each wash.
    - Repeat step iii 4 more times. Following the last wash, test wells should be thoroughly aspirated to remove as much moisture as possible, but do not allow wells to completely dry at any time.
- Add 2 free falling drops (approximately 100 µL) of Enzyme Conjugate to each well. Firmly shake/swirl the plate for 30 seconds.
- Cut the plate sealer to size and press it firmly onto the top of the microwells to seal. Incubate the plate for 30 minutes at 21–27 C. Alternatively, laboratories equipped with a shaker/incubator (StatFax™-2200) can incubate and rotate the plate for 15 minutes at 24 C and at Setting 5 (1000 rpm).
- Carefully remove the plate sealer and wash the wells as described in Step 5.

- Add 2 free-falling drops (approximately 100 µL) of Premier Substrate I Solution to each well. Firmly shake/swirl the plate for 30 seconds. Incubate the plate for 10 minutes at 21–27 °C.
- Add 2 free-falling drops (approximately 100 µL) of Premier Stop Solution I to each well. Firmly shake/swirl the plate for 30 seconds. **Note:** The initial color of a positive reaction is blue, which changes to yellow upon addition of Premier Stop Solution I.
- Clean the underside of all wells with a lint-free tissue
- Inspect and record reactions. Test results can be read visually or using a spectrophotometric reader.
  - Visual Determination – Read within 15 minutes after adding Premier Stop Solution I.
  - Spectrophotometric Determination – Zero EIA reader on air. Read absorbance at 450 nm or 450/630 nm within 15 minutes of adding Premier Stop Solution I.

## INTERPRETATION OF RESULTS

### Visual Reading

Negative = Colorless to very faint yellow

Positive = Definite yellow color

A very faint yellow color must be evaluated by a spectrophotometric reading.

### Spectrophotometric Single Wavelength (450 nm)

Negative: < 0.150

Positive: ≥ 0.150

Negative Control: < 0.150

Positive Control: ≥ 0.600

### Spectrophotometric Dual Wavelength (450/630 nm)

Negative: < 0.100

Positive: ≥ 0.100

Negative Control: < 0.100

Positive Control: ≥ 0.600

If a Negative Control is < 0.000, reblank the plate reader to air and reread the plate.

A positive result indicates the presence of *Campylobacter* antigens. A negative result indicates the absence of *Campylobacter* antigens, or that the level of antigens is below what can be detected by the assay. The magnitude of the OD above the cutoff is neither indicative of the severity or extent of *Campylobacter* infection, nor can it be correlated to an endpoint titer. Extremely strong positive samples may yield either an intense yellow color or a purple precipitate within a few minutes of stopping the reaction. In this case, the spectrophotometer may yield an "out" reading. This reading is considered a positive.

If the frequency of low positive results (OD between 0.150 and 0.200 for single wavelength and between 0.100 and 0.150 with dual wavelength) is greater than 5% of specimens tested, this may indicate insufficient washing. More vigorous washing or increasing the washes to seven washes in step 5 of the procedure is recommended.

## QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

- At the time of each use, kit components should be visually examined for obvious signs of microbial contamination, freezing or leakage. Do not use contaminated or suspect reagents.
- The performance of Premier CAMPY should be verified using the Positive Control and Sample Diluent/Negative Control with each test batch. See the INTERPRETATION OF RESULTS section above for a description of the expected results for control reagents. Tests should be considered invalid when either control reagent does not produce the expected results. In such cases, repeat tests and controls. If, on repeat testing, the expected reactions are still not observed with in-date reagents, call Meridian's Technical Support Services at 800-343-3858 (USA) or your authorized Meridian distributor for assistance.
- The controls are used to monitor reagent reactivity. Failure of the controls to produce the expected results can mean that one or more of the reagents are no longer reactive at the time of use, that the test was not performed correctly, or that reagents or samples were not added. **Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.**
- Specimen matrix interference has not been observed in this assay as samples are significantly diluted before testing in Sample Diluent. For this reason, the Positive and Negative Controls supplied as part of this assay are prepared in matrices similar to the Sample Diluent. If control materials that are identical in composition to test samples are preferred, the user can prepare these by diluting known positive and negative specimens in Sample Diluent according to the SPECIMEN PREPARATION section of this insert. Add 100 µL of user prepared controls to test wells.

## EXPECTED VALUES

The performance of Premier CAMPY was evaluated during 2008 in several geographic regions of the United States. The incidence of positive samples (either as *C. jejuni* or *C. coli*) during this period was approximately 1.5%. The expected frequency for an individual laboratory may differ from this number since it is dependent on factors such as locale, the time of the year, and whether an outbreak has occurred.

## LIMITATIONS OF THE PROCEDURE

- The test is qualitative and no quantitative interpretation should be made with respect to the values.
- Premier CAMPY cannot be used as the sole means of determining Campylobacteriosis. Test results must be used in conjunction with information available from the patient clinical evaluation and other diagnostic procedures.
- Overincubation of the tests may lead to an increase in false-positive test results. Conversely, incubation for periods less than those defined in this insert can result in an increase in false-negative tests. Follow incubation times defined in this insert.
- Suspect failure of the washing method or device if the Negative and/or Positive Controls consistently produce out-of-specification results. Increasing the number of washes, washing more vigorously, decanting more thoroughly or recalibrating washing devices should correct the problem. If problems persist, contact Meridian's Technical Support Services at 800-343-3858 (USA) or your authorized Meridian distributor for assistance.

## SPECIFIC PERFORMANCE CHARACTERISTICS

Premier CAMPY was evaluated in 2008 by four independent test sites located in the Western, Midwestern and Southeastern regions of the United States. A total of 2073 qualified patient samples were evaluated; of these, 166 were retrospective frozen samples. The majority (1862/2073) were collected in a Cary-Blair-based transport and preservative medium. The remaining 211 samples were tested in the unpreserved state. Samples were collected from males (41%) and females (57%). In the case of 2% of the patients, the sex was not known. The age groups of the patients ranged from less than one month of age to 97 years. No differences in test performance were observed based on patient age or sex. The following tables show the assay performance by clinical site, patient age and sample type.

Table 1. Performance characteristics by clinical site

Site	Positive Samples			Negative Samples		
	Premier/Culture	Sensitivity %	95% CI	Premier/Culture	Specificity %	95% CI
Site 1	22/23	95.7%	79.0 – 99.2%	189/193	97.9%	94.8 – 99.2%
Site 2	0/0	N/A	N/A	51/51	100%	93.0 – 100%
Site 3	26/27	96.3%	81.7 – 99.3%	1429/1511	94.6%	93.3 – 95.6%
Site 4	11/11	100%	74.1 – 100%	255/257	99.2%	97.2 – 99.8%
Total Sites	59/61	96.7%	88.8 – 99.1%	1924/2012	95.6%	94.6 – 96.4%

Table 2. Performance characteristics by patient age

Patient Age	Positive Samples			Negative Samples		
	Premier/Culture	Sensitivity %	95% CI	Premier/Culture	Specificity %	95% CI
Birth to 1 month	0/0	N/A	N/A	12/13	92.3%	66.7 – 98.6%
> 1 month to 2 years	3/3	100%	43.9 – 100%	338/348	97.1%	94.8 – 98.4%
> 2 years to 12 years	5/5	100%	56.6 – 100%	374/388	96.4%	94.0 – 97.8%
> 12 years to 21 years	2/2	100%	34.2 – 100%	139/145	95.9%	91.3 – 98.1%
> 21 years	32/34	94.1%	80.9 – 98.4%	1052/1109	94.9%	93.4 – 96.5%
Not Defined	17/17	100%	81.6 – 100%	9/9	100%	70.1 – 100%

Table 3. Performance characteristics by sample type (preserved vs unpreserved)

Specimen Type	Positive Samples			Negative Samples		
	Premier/Culture	Sensitivity %	95% CI	Premier/Culture	Specificity %	95% CI
Preserved	42/43	97.7%	87.9 – 99.6%	1733/1819	95.3%	94.2 – 96.2%
Unpreserved	17/18	94.4%	74.2 – 99.0%	191/193	99.0%	96.3 – 99.7%

Table 4. Performance characteristics of fresh vs frozen samples

Fresh/Frozen	Positive Samples			Negative Samples		
	Premier/Culture	Sensitivity %	95% CI	Premier/Culture	Specificity %	95% CI
Fresh	17/18	94.4%	74.2 – 99.0%	1810/1889	95.8%	94.8 – 96.6%
Frozen	42/43	97.7%	87.9 – 99.6%	114/123	92.7%	86.7 – 96.1%

## ANALYTICAL SENSITIVITY

The analytical sensitivity of this assay for *C. jejuni* and *C. coli* was based on 45 tests for each measurand and with a stated probability (e.g. 95%) of obtaining positive responses at the following levels of the measurands: *C. jejuni* 1.2 x 10<sup>5</sup> cells/mL; *C. coli* 8.0 x 10<sup>6</sup> cells/mL.

## ASSAY REACTIVITY

The following *Campylobacter* stock cultures from different sources were tested and produced positive reactions at 1.1 x 10<sup>8</sup> CFU/mL with Premier CAMPY: *C. coli* strains 10956, 17755, 36994 and 53138 and *C. jejuni* strains 6951, 10940, 12081, 29411 and 38106.

## CROSSREACTIVITY STUDIES

Crossreactivity studies were performed with positive and negative stool specimens inoculated with bacterial or fungal organisms to a final concentration of 1.1 x 10<sup>8</sup> CFU/mL or virus ranging from 1.3 x 10<sup>4</sup> to 3.1 x 10<sup>5</sup> TCID<sub>50</sub>/mL. None of the following organisms in stool reacted with Premier CAMPY: *Aeromonas hydrophila*, *Bacteroides fragilis*, *Campylobacter fetus*, *Campylobacter lari*, *Candida albicans*, *Clostridium difficile*, *Clostridium perfringens*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Escherichia coli*, *Escherichia coli* O157:H7, *Escherichia fergusonii*, *Escherichia hermannii*, *Helicobacter pylori*, *Klebsiella pneumoniae*, *Lactococcus lactis*, *Listeria monocytogenes*, *Peptostreptococcus anaerobius*, *Plesiomonas shigelloides*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Pseudomonas fluorescens*, *Salmonella* Groups B-E, *Serratia liquefaciens*, *Serratia marcescens*, *Shigella boydii*, *Shigella dysenteriae*, *Shigella flexneri*, *Shigella sonnei*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Vibrio parahaemolyticus*, *Yersinia enterocolitica*, Adenovirus Types 40 and 41, Coxsackievirus, Echovirus, Rotavirus.

## TESTS FOR INTERFERING SUBSTANCES

The following substances, at the specified saturated solvent/diluent concentrations, do not interfere with test results in the final concentrations listed: Barium sulfate (5 mg/mL), fecal fat (equivalent to 2.65 mg stearic plus 1.3 mg palmitic acids per mL), hemoglobin (as methemoglobin) (3.2 mg/mL), Imodium AD® (0.00667 mg/mL), Kaopectate® (0.87 mg/mL), mucin (3.33 mg/mL), Mylanta® (4.2 mg/mL), Pepto-Bismol® (0.87 mg/mL), Prilosec® (0.5 mg/mL), Tagamet® (0.5 mg/mL), TUMS® (0.5 mg/mL), whole blood (5% v/v).

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### REFERENCES

- Coker AO, Isokpehi RD, Thomas BN, Amisu KO, Obi CL. Human Campylobacteriosis in developing countries. EID. 2002 Mar 2;8(3):237-43.
- Friedman CR, Hoekstra RM, Samuel M, Marcus R, Bender J, Shiferaw B et al. Risk factors for sporadic *Campylobacter* infection in the United States: A case-control study in FoodNet sites. CID 2004;38 (Suppl 3):S265-96.

SN11174

REV. 01/15

Manufactured By

Meridian Bioscience, Inc.  
 USA/Corporate Office  
 3471 River Hills Drive  
 Cincinnati, Ohio 45244  
 Telephone: 513.271.3700  
 Orders/Customer Service:  
 800.543.1980  
 Technical Support Center:  
 800.343.3858  
 Information Fax: 513.272.5432  
 Ordering Fax: 513.271.0124

Meridian Bioscience Europe S. r. L.  
 Via dell'Industria, 7  
 20020 Villa Cortese, Milano  
 ITALY  
 Tel: +39 0331 43 36 36  
 Fax: +39 0331 43 36 16  
 Email: info@meridianbioscience.eu  
 WEB: www.meridianbioscience.eu

Authorized Representative

Meridian Bioscience Europe s.a.n.v.  
 2 Avenue du Japon - 1420 Braine l'Alleud  
 BELGIUM  
 Tel: +32 (0) 67 89 59 59  
 Fax: +32 (0) 67 89 59 58  
 Email: info.bn@meridianbioscience.eu

Meridian Bioscience Europe France  
 34 rue de Ponthieu - 75008 Paris  
 FRANCE  
 Tel: +33 (0) 1 42 56 04 40  
 Fax: +33 (0) 9 70 06 62 10  
 Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.  
 Postbus 301 - 5460 AH Veghel  
 NETHERLANDS  
 Tel: +31 (0) 411 62 11 66  
 Fax: +31 (0) 411 62 48 41  
 Email: info.bn@meridianbioscience.eu

## INTERNATIONAL SYMBOL USAGE CHART

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de símbolos, Erläuterung der graphischen symbole)

	Use By / Utilizzare entro / Utiliser jusqu'à / Fecha de caducidad / Verwendbar bis	<b>CONTROL +</b>	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
<b>LOT</b>	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	<b>CONTROL -</b>	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
<b>IVD</b>	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	<b>EC REP</b>	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostico in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	<b>SMP PREP DIL SPE</b>	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet
<b>REF</b>	Catalogue number / Numero de catalogo / Référence au catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	<b>BUF RXN</b>	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	<b>SOLN STOP</b>	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	<b>CONJ ENZ</b>	Enzyme Conjugate / Coniugato enzimatico / Conjugué enzymatique / Conjugado enzimático / Enzymkonjugat
<b>SN</b>	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	<b>CONTROL</b>	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
<b>TEST</b>	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	<b>REAG</b>	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	<b>BUF WASH</b>	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
<b>BUF</b>	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenza / Mise En Garde / Advertencia / Warnhinweis
<b>CONJ</b>	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	<b>DIL SPE</b>	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
<b>SUBS</b>	Substrate / Substrato / Substrat / Substrato / Substrat	<b>BUF WASH 20X</b>	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
		<b>DET REAG</b>	Detection Reagent / Reagente Diretto / Réactif de Détection / Reactivo de Detección / Nachweis Reagenz

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.