PREMIER® ADENOCLONE® Type 40/41

EIA for the detection of Adenovirus Serotypes 40 and 41 in Human Fecal Samples

REF 696006

IVD In vitro diagnostic medical device

INTENDED USE

Premier Adenoclone-Type 40/41 is an enzyme immunoassay (EIA) for the qualitative detection of enteric adenovirus serotypes 40 and 41 in human fecal specimens.

SUMMARY AND EXPLANATION OF THE TEST

Human adenoviruses consist of 41 known types, classified primarily by serology and DNA electrophoretotyping. The viruses are nonenveloped, DNA-containing icosahedrons composed of 252 capsomers. The virion contains a number of proteins associated with genus and type specificity. Those associated with genus (group) specificity are the alpha portion of the hexon, the penton base, the delta portion of the fiber and polypeptide IIIa. Those associated with species (type) specificity are the epsilon portion of the hexon and the gamma portion of the fiber. Those associated with both genus and type specificity are the major core protein and polypeptide IX.¹¹ The monoclonal antibodies in Premier Adenoclone-Type 40/41 specific EIA are directed against the type-specific portion of the hexon of adenovirus types 40 and 41. The type-specific monoclonal antibodies in the EIA will not react with any of the other adenovirus serotypes.¹

The enteric adenovirus types 40 and 41 are the only adenoviruses which have been consistently associated with gastroenteritis in infants and young children, and may be the second major cause of gastroenteritis after rotaviruses.^{2,4, 8, 13-17} Until recently, these viruses could be detected only by electron microscopy (EM). These adenovirus types can now be detected by cultivation in Graham 293 cells, a human embryonic kidney cell line transformed by adenovirus type 5 DNA, but cannot readily be cultivated in all lines traditionally used for adenovirus isolation, e.g. HEP-2 cells.³ Identification of adenovirus in tissue culture isolation is followed by the appearance of characteristic cytopathic effect (CPE). Generally CPE occurs within 3 to 7 days but may require up to 28 days.⁷⁻¹¹ Once the virus is isolated, it must be further characterized by analysis of restriction endonuclease digests of the virus DNA followed by gel electrophoresis.⁹

The Premier Adenoclone-Type 40/41 test kit for direct stool specimens is a rapid and specific test for adenovirus 40 and 41 that eliminates the need for tissue culture and DNA electrophoretotyping.

BIOLOGICAL PRINCIPLES

Premier Adenoclone-Type 40/41 utilizes monoclonal antibodies in a solid phase sandwich type assay. Breakaway microwells are coated with a monoclonal antibody directed against the group specific antigen for all known human adenoviruses. An aliquot of fecal suspension is added to the microwell and incubated simultaneously with anti-adenovirus types 40 and 41 monoclonal antibodies, conjugated to horseradish peroxidase, resulting in the adenovirus antigen being sandwiched between the solid phase and enzyme conjugate. After a 60 minute incubation at room temperature, the sample well is washed with distilled or deionized water to remove unbound specimen and enzyme labelled antibodies. Enzyme substrate A (urea peroxide) and substrate B (tetramethylbenzidine) are added to the wells and incubated for 10 minutes at room temperature. Any bound enzyme-conjugate in the wells converts the colorless substrate to a blue color.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

- Microwells** - Coated with Group Specific Anti-Adenovirus Murine Monoclonal Antibody.
- Sample Diluent/Negative Control** - Buffered saline with 0.02% thimerosal as preservative.
- Positive Control** - Inactivated Adenovirus (AD41) in buffer with 0.02% thimerosal as preservative.
- Enzyme Conjugate** - Horseradish Peroxidase conjugated to anti-adenovirus types 40 and 41 murine monoclonal antibodies in a buffered protein solution with gentamicin and 0.02% thimerosal as preservatives.
- Substrate A** - Substrate buffer contains urea peroxide.
- Substrate B** - Substrate buffer contains tetramethylbenzidine.
- Stop Solution** - Contains 1 N sulfuric acid. **Caution:** Avoid contact with skin. If skin contact occurs, flush with water.
- Microwell Holder**
- Sample Transfer Pipettes**

MATERIALS NOT PROVIDED

- 12 x 75 mm test tubes, test tube rack
- Deionized or distilled water
- Absorbent paper
- Precision micropipette tips to deliver 100 µL and 1000 µL (optional)
- Waste container with a 1:10 dilution of household bleach. For autoclaving, use an iodophor disinfectant
- Microwell plate reader capable of reading absorbance at 450 nm (optional)
- Device for dispensing wash solution such as multi-channel pipette, syringe with manifold, wash bottle, etc
- Timer (minimum 1 hour)

PRECAUTIONS

- All reagents are for in vitro diagnostic use only.
- Do not mouth pipette samples or reagents. Avoid contact with broken skin or mucous membranes.
- Do not smoke, eat or drink in areas where specimens or kit reagents are handled.
- Wear disposable gloves while handling samples and wash hands after assay is complete. Patient specimens, assay controls and all materials coming into contact with them should be handled at Biosafety Level 2 as recommended in the CDC/NIH manual "Biosafety in Microbiology and Biomedical Laboratories".
- Avoid skin contact with stop solution (1N sulfuric acid). It may cause irritation and burns. Flush with water if contact occurs.
- Dispose of all materials used to perform the test by autoclaving at 121 C for at least one hour. Liquid waste may be disposed by mixing with a 1:10 dilution of household bleach for a minimum of 30 minutes. **CAUTION:** Liquid waste containing stop solution must be neutralized before addition of household bleach. Avoid splashing or generation of aerosols.
- Do not use Premier Adenoclone-Type 40/41 reagents beyond the kit expiration date. Each reagent has been optimized for maximum performance. Dilution or adulteration of these reagents may result in a loss of sensitivity. Incubation time and temperatures other than those specified may give erroneous results.
- Do not interchange or mix different lots of Premier Adenoclone-Type 40/41 reagents.
- Avoid microbial contamination of reagents or incorrect results may occur. Contamination of samples could cause erroneous results.
- Use separate pipettes or pipette tips for each sample, control and reagent.
- DO NOT REUSE MICROWELLS.

HAZARD and PRECAUTIONARY STATEMENTS

Refer to the SDS, available at www.meridianbioscience.com for Hazards and Precautionary Statements.

SHelf LIFE AND STORAGE

Store kit reagents at 2-8 C. Bring kit reagents to room temperature (20-30 C) before use and promptly return to 2-8 C after use. Keep microwells in pouch until the pouch reaches room temperature. Return all unused microwells to their original foil pouch. The microwells are stable for the shelf life of the kit after opening.

INDICATION OF INSTABILITY OR DETERIORATION

The following conditions may indicate reagent deterioration:

- Any evidence of microbial contamination or heavy precipitation.
- Any blue color in the substrate solutions before addition to microwells.
- A negative control value greater than 0.150 absorbance units at 450 nm may indicate deterioration of reagents.

- A positive control value of less than 0.3 absorbance units may indicate deterioration of reagents. If any of the above conditions are observed, contact our Technical Services Department at 1-800-343-3858 or your local distributor.

REAGENT PREPARATION

- Bring all reagents to room temperature (20-30 C) before use.
- Return all reagents to 2-8 C immediately after use.
- Do not allow microwells to dry between steps.
- Prepare decontaminating vessel for discarding reagents and materials.
- Reproducibility in any EIA assay is largely dependent upon the consistency with which the microwells are washed. Carefully follow the recommended washing sequence as outlined in the EIA test procedure.


SPECIMEN COLLECTION AND PREPARATION

Add 1 mL of Sample Diluent/Negative Control to properly marked tube, using a transfer pipette or precision pipette. Add samples by one of the following:

- Stool specimens or rectal swabs should be collected in containers that do not contain media, preservatives, animal serum or detergent as any of these additives may interfere with the Premier Adenoclone-Type 40/41 assay. Caution should be taken when using rectal swabs to insure that sufficient fecal sample (30 to 50 mg of raw stool) is obtained.
- Solid stool — press sample into transfer pipette to first mark.
- Liquid stool — aspirate sample into transfer pipette to first mark.
- Rectal swab — swirl swab in 1 mL of sample diluent to release fecal material. Firmly press swab against side of tube to remove liquid. Mix thoroughly. *Do not dilute further.*
- Diluted samples may be stored at 2-8 C for three days without interfering with the assay performance. For long term storage of undiluted specimens, -20 C or colder is recommended. Repeated freezing and thawing of samples is not recommended and may cause erroneous results. Do not store in self-defrosting freezers.

Use the sample transfer pipette to add two drops (100 µL) of diluted stool specimen to each well.

Sample Transfer Pipette

TEST PROCEDURE

- Snap off a sufficient number of wells for samples and the controls and insert into the microwell holder. Record sample position.
- Add 2 drops (100 µL) of diluted fecal sample, negative control (sample diluent) or positive control to the bottom of separate microwells.
- Add 2 drops (100 µL) of enzyme conjugate to each microwell. Mix by gently swirling on table top.
- Incubate at room temperature for 60 ± 5 minutes.
- Pour the liquid out of the wells into a discard vessel. Tap the microwell holder upside down vigorously against absorbent paper to ensure complete removal of liquid from the wells.
- Fill all the wells to overflowing with deionized water and pour the liquid out as in Step 5.
- Repeat the washing procedure (Steps 5 & 6) four more times for a total of 5 washes.
- Add 2 drops (100 µL) of Substrate A (urea peroxide) solution, to each microwell.
- Add 2 drops (100 µL) of Substrate B (TMB) solution to each microwell.
- Incubate at room temperature for 10 minutes.
- Discard the used wells and absorbent paper into biohazardous waste after reading and recording results. Visual determinations can be made after incubating for 10 minutes in Step 10, and should be read within 10 minutes. Samples with blue color greater than the negative control are positive. Samples showing equal or less color than negative control are negative.

Optional - Spectrophotometric Procedure

Spectrophotometric determinations can be made by adding 2 drops (100 µL) of stop solution (1 N sulfuric acid) to each microwell after the incubation in step 10. Read the absorbance of each microwell at 450 nm using a >600 nm reference filter (optional) against an air blank within 60 minutes.

INTERPRETATION OF RESULTS

- Visual Determination:
Any sample with blue color more intense than that of the negative control is considered positive. Any sample with color equal to or less intense than the negative control is considered negative.
- Spectrophotometric Determination:
Specimens with absorbance units (A_{450}) equal to or greater than 0.150 are considered positive. Specimens with absorbance less than 0.150 are considered negative. Occasionally, a discrepant result may occur between visual and spectrophotometric determinations with samples containing low amounts of antigen. Spectrophotometric determination, being an objective method, is slightly more accurate.

NOTE: A precipitate may form in high positive samples. This will not affect the results.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

Premier Adenoclone-Type 40/41 is supplied with a positive and negative (sample diluent) control which are to be used every time an assay is performed. These reagents are used to insure that the kit reagents are working properly. The positive control well should be clearly blue with an absorbance at 450 nm of greater than 0.30. The negative control must have no visible color, and an absorbance at 450 nm of less than 0.150.

If the positive and negative controls do not perform properly, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

EXPECTED VALUES

Rate of positivity may vary depending on geographic location, method of specimen collection, handling and transportation, test employed and general health environment of patient population under study.

The frequency of adenovirus infection will vary with the clinical syndrome and the age of the individual. It has been reported that approximately 5 to 10% of viral gastroenteritis in infants and children under two years is caused by adenovirus types 40 and 41.^{2,4, 6, 8, 13}

Performance characteristics have not yet been fully determined for direct testing of neonatal stool samples. Lack of false positive results should be determined in each individual laboratory before neonatal stools are tested directly.

LIMITATIONS OF THE PROCEDURE

Premier Adenoclone-Type 40/41 is highly specific and sensitive for adenovirus types 40 and 41 antigen. Failure to detect adenovirus 40/41 may be a result of factors such as collection of specimen at an improper time in the disease, improper sampling or handling of the specimen.

A negative result for a stool sample may not preclude the presence of an infection with non-enteric adenovirus in other body fluids. If respiratory or ophthalmic infection is suspected specimens from the site of infection should be cultured and confirmed with Meridian Bioscience's Premier Adenoclone (Catalog #696007).

All positive results must be interpreted with caution since adenovirus is capable of latency and recrudescence.¹² Asymptomatic shedding may occur for months after infection.¹³ Enteric adenoviruses may be found in the stools of asymptomatic children.¹³ Test results should be interpreted in conjunction with information available from epidemiological studies or clinical evaluation of patient or other diagnostic procedures.¹²

While the relationship between adenovirus 40/41 and gastroenteritis is well established, co-infection with other viral and bacterial pathogens (including non-enteric adeno-viruses) is possible. Therefore, bacteriological tests should be performed in parallel with this test to rule out bacteriological etiology. Meridian Bioscience's Premier Rotaclo® rotavirus diagnostic kit (Catalog #696004) may be performed in parallel to rule out rotaviral gastroenteritis. False positives could occur with high levels of *S. aureus* possessing protein A such as Cowan strain. *Staphylococcal enterocolitis* is an uncommon disease generally described in adult patients and is extremely rare in infants and children.^{4, 10}

SPECIFIC PERFORMANCE CHARACTERISTICS

Premier Adenoclone-Type 40/41 for enteric adenovirus antigen test was tested in collaboration with institutions in Arizona, Massachusetts, Ontario and Manitoba. The results were compared with initial Electron Microscopy followed by DNA restriction endonuclease analysis of adenovirus positive samples. All of the samples were tested by Premier Adenoclone-Type 40/41 EIA and scored visually then read spectrophotometrically. There was 100% correlation between visual and spectrophotometric determinations.

EM/DNA Restriction Analysis

	+	-
Premier Adenoclone-Type 40/41 EIA	58	3
	1	169
Sensitivity	98%	
Specificity	98%	
Total Agreement	98%	

ANALYTICAL SENSITIVITY

Virus particle counts from known adenovirus 40 and 41 stool samples were determined by electron microscopy. Serial dilutions were tested by Premier Adenoclone-Type 40/41. These results show that adenovirus particle counts as low as 10⁴ per mL can be detected by Premier Adenoclone-Type 40/41.

Virus Particle Count	Ad 40 (A ₄₅₀)	Ad 41 (A ₄₅₀)
2 x 10 ⁵	> 2.0	—
1 x 10 ⁶	> 2.0	> 2.0
5 x 10 ⁵	> 2.0	> 2.0
2.5 x 10 ⁵	> 2.0	> 2.0
1.25 x 10 ⁵	1.510	> 2.0
6.25 x 10 ⁴	1.001	> 2.0
3.1 x 10 ⁴	0.595	1.445
1.6 x 10 ⁴	0.307	0.781
8 x 10 ³	0.167	0.425
4 x 10 ³	0.108	0.219
2 x 10 ³	—	0.142

CROSSREACTIVITY

- Non-enteric Adenovirus Subgroups: Adenoviral 3 + CPE culture supernatants were tested by Premier Adenoclone-Type 40/41. The data shows no cross reactivity with the Premier Adenoclone-Type 40/41.
PLEASE NOTE: CPE culture supernates are not approved samples for use with this product. The information derived from tests with culture filtrates that is cited here is provided for information purposes only. Do not use culture supernates to verify specificity or sensitivity.

Adenoviruses		Premier Adenoclone-Type 40/41	Premier Adenoclone-Group specific*
Subgroup	Type	(A ₄₅₀)	(A ₄₅₀)
A	12	0.043	> 3.0
A	18	0.042	> 3.0
A	31	0.037	> 3.0
B	3	0.050	> 3.0
B	7	0.061	> 3.0
B	21	0.027	> 3.0
C	1	0.042	> 3.0
C	2	0.041	> 3.0
C	6	0.039	> 3.0
D	22	0.030	> 3.0
D	37	0.033	> 3.0
D	39	0.035	> 3.0
E	4	0.047	> 3.0
Negative control		0.049	0.060
Positive control		1.071	0.580

* The Adenovirus subgroups were tested by Meridian Bioscience's Premier Adenoclone (Catalog #696007) group specific EIA, which detects all 41 types of human adenovirus to demonstrate the presence of adenovirus.

- Microorganisms and Viruses: Common intestinal pathogens and other organisms occasionally present in feces were tested in Premier Adenoclone-Type 40/41 and showed no cross reactivity in the test.

Microorganisms

<i>Acholeplasma laidlawii</i>	<i>Escherichia coli</i>	<i>Proteus mirabilis</i>
<i>Campylobacter lariidis</i>	<i>Klebsiella pneumoniae</i>	<i>Pseudomonas aeruginosa</i>
<i>Campylobacter fetus</i>	<i>Mycoplasma arginini</i>	<i>Salmonella dysenteriae</i>
<i>Candida albicans</i>	<i>Mycoplasma hominis</i>	<i>Serratia marcescens</i>
<i>Clostridium difficile</i>	<i>Mycoplasma hyorhinis</i>	<i>Shigella sonnei</i>
<i>Enterobacter aerogenes</i>	<i>Mycoplasma orale</i>	<i>Staphylococcus aureus</i> * (Cowan)
<i>Enterobacter cloacae</i>	<i>Mycoplasma pneumoniae</i>	

*False positives could occur with high levels of *S. aureus*, possessing protein A such as Cowan strain. *Staphylococcal enterocolitis* is an uncommon disease generally described in adult patients and is extremely rare in infants and children.^{9,10}

Viruses

Astrovirus	Echovirus 22	Reovirus 1
Calicivirus	Human Coronavirus	Reovirus 2
Chlamydia trachomatis (serotype A)	Influenza B	Reovirus 3
Coxsackievirus A-9	Measles	Respiratory Syncytial virus
Coxsackievirus B-3	Myxo/Paramyxoviruses SV-5	Rhinovirus 3
Coxsackievirus B-4	Norwalk virus	Rhinovirus 17
Coxsackievirus B-5	Papovavirus SV-40	Rhinovirus 25
Cytomegalovirus	Parainfluenza-1	Rhinovirus 43
Echovirus 4	Parainfluenza-2	Rhinovirus 62
Echovirus 8	Parainfluenza-3	Rotavirus
Echovirus 11	Poliovirus I	Simian foamyvirus (Types 2 & 8)
Echovirus 16	Poliovirus II	Small Round virus (SRV)
Echovirus 20	Poliovirus III	Varicella zoster virus

A urine sample was tested in the Premier Adenoclone-Type 40/41 and found to be negative.

REPRODUCIBILITY

Inter-assay variation

The inter-assay precision was determined by diluting known adenovirus 40 and adenovirus 41 positive samples to obtain high, medium, and low absorbance values. These diluted samples were tested in 21 different assays on different days by several technicians and the mean and coefficient of variation determined.

Adenovirus Type 40	High	Medium	Low	Negative
Mean	1.456	0.856	0.262	0.042
S.D.	0.174	0.088	0.025	0.017
CV	12.0%	10.3%	9.5%	40.5%

Adenovirus Type 41	High	Medium	Low
Mean	0.892	0.565	0.193
S.D.	0.072	0.068	0.019
CV	8.1%	12.0%	9.8%

Intra-assay variation

The intra-assay precision was determined by diluting known adenovirus 40 and adenovirus 41 positive samples to obtain high, medium, and low absorbance values. These diluted samples were tested 21 times within a single assay and the mean and coefficient of variation determined.

Adenovirus Type 40	High	Medium	Low	Negative
Mean	1.522	0.894	0.193	0.064
S.D.	0.038	0.038	0.011	0.022
CV	3.8%	4.3%	5.7%	34.3%

Adenovirus Type 41	High	Medium	Low
Mean	0.924	0.647	0.156
S.D.	0.086	0.039	0.012
CV	9.3%	6.0%	7.7%

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REFERENCES

- Herrmann, JE, Perron-Henry, DM, Stobbs-Walro, D, and Blacklow, NR. Preparation and characterization of monoclonal antibodies to enteric adenovirus types 40 and 41. Arch. Virol. 1987;94:259-265.
- Uhnnoo, L, Wadell, G., Svensson, L., and Johansson, M. Two new serotypes of adenovirus from human infantile diarrhea. Dev. Biol. Stand. 1983;55:311-318.
- Brandt, CD, Kim, HW, Yolken, RH, Kapikian, AZ, Arrobio, JO, Rodriguez, WJ, Wyatt, RG, Chanock, RM, and Parrott, RH Comparative epidemiology of two rotavirus serotypes and other viral agents associated with pediatric gastroenteritis. Amer. J. Epidemiol. 1979;110: 243-254.
- Uhnnoo, L, Wadell, G, Svensson, L, Johansson, ME. Importance of enteric adenoviruses 40 and 41 in acute gastroenteritis in infants and young children. J. Clin. Microbiol. 1984;20:365-372.
- Takiff, HF, Straus SE, and Garon, CF. Propagation and *in vitro* studies of previously non-cultivable enteric adenoviruses in 293 cells. Lancet ii: 1981:832-834.
- Wadell, G. Molecular epidemiology of human adenoviruses. In: Doerfler W., ed. The molecular biology of adenoviruses 2. 30 years of adenovirus research 1953-1983. Current topics in microbiology and immunology. Berlin: Springer-Verlag, 1984;110:191-220.
- Herrmann, JE, Personal communication, University of Massachusetts Medical School, Worcester, Massachusetts, 1987.
- Brandt, CD, Kim, HW, Rodriguez, WJ, Arrobio, JO, Jeffries, BC, Stallings, EP, Lewis, C, Miles, AJ, Chanock, RM, Kapikian, AK, Parrott, RH. Pediatric viral gastroenteritis during eight years of study. J. Clin. Microbiol. 1983;18:71-78.
- Gutman, LT, Idriss, ZH, Gehlbach, S, Blackman, L. Neonatal staphylococcal enterocolitis: association with indwelling feeding catheters and *S. aureus* colonization. J. Pediatr. 1976;88:836-839.
- Hay, P, and McKenzie, P. Side effects of oxytetracycline therapy. Lancet i: 1954:945-950.
- Horwitz, MS, Adenoviruses and Their Replication. In B.N. Fields ed. Virology, Ravew Press, New York, pp. 433-438.
- Horwitz, MS, Adenoviruses and Their Replication. In B.N. Fields ed. Virology, Ravew Press, New York, pp. 480.
- Cukor, G, and Blacklow, NR, Human viral gastroenteritis. Microbiological Reviews. 1984;48:157-179.
- Chiba, S, Wakamura, I, Urasawa, S, Nakata, S, Taniguchi, K, Fujinaga, K, Nakao, T. Outbreak of infantile gastroenteritis due to type 40 adenovirus. Lancet ii: 1983:954-957.
- Johansson, ME, Uhnnoo, I, Kidd, AH, Madeley, CR, Wadell, G. Direct identification of enteric adenovirus, a candidate new serotype, associated with infantile gastroenteritis. J. Clin. Micro., 1980;12:95-100.
- Wigand, R, Baumeister, HG, Maass, G, Kuhn, J, Nammer, HJ. Isolation and identification of enteric adenoviruses. J. Med. Virol. 1983;11:233-240.
- de Jong, JC, Wigand, R, Kidd, AH, Wadell, G, Kapsenberg, JG, Muzerie, CJ, Wermenbol, AG, Firtzloff, RG. Candidate adenoviruses 40 and 41: fastidious adenoviruses from human infant stool. J. Med. Virol. 1983;11:215-231.

SN11244

REV. 12/14

 Manufactured By

Meridian Bioscience, Inc.
 USA/Corporate Office
 3471 River Hills Drive
 Cincinnati, Ohio 45244
 Telephone: 513.271.3700
 Orders/Customer Service:
 800.543.1980
 Technical Support Center:
 800.343.3858
 Information Fax: 513.272.5432
 Ordering Fax: 513.271.0124

Meridian Bioscience Europe S. r. l.
 Via dell'Industria, 7
 20020 Villa Cortese, Milano
 ITALY
 Tel: +39 0331 43 36 36
 Fax: +39 0331 43 36 16
 Email: info@meridianbioscience.eu
 WEB: www.meridianbioscience.eu

 Authorized Representative

Meridian Bioscience Europe s.a./n.v.
 2 Avenue du Japon - 1420 Braine l'Alleud
 BELGIUM
 Tel: +32 (0) 67 89 59 59
 Fax: +32 (0) 67 89 59 58
 Email: info.bn@meridianbioscience.eu








Meridian Bioscience Europe France
 34 rue de Ponthieu - 75008 Paris
 FRANCE
 Tel: +33 (0) 1 42 56 04 40
 Fax: +33 (0) 9 70 06 62 10
 Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
 Postbus 301 - 5460 AH Veghel
 NETHERLANDS
 Tel: +31 (0) 411 62 11 66
 Fax: +31 (0) 411 62 48 41
 Email: info.bn@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simbolos, Erläuterung der graphischen symbole)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis		Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung		Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum		Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.		Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon: diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet
	Catalogue number / Numero de catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten		Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen		Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung		Enzyme Conjugate / Coniugato enzimatico / Conjugé enzymatique / Conjugado enzimático / Enzymkonjugat
	Serial number / Numero de serie / Numéro de série / Número de serie / Seriennummer		Assay Control / Controllo dei test / Test de contrôle / Control de Ensayo / Kontrolltest
	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät		Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum		Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
	Buffer / Solution tamponne / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat		Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillon / Diluyente de muestra / Probenverdünnungspuffer
	Substrate / Substrato / Substrat / Substrato / Substrat		Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
			Detection Reagent / Reagente Diretto / Réactif de Détection / Reactivo de Detección / Nachweis Reagenz

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.