ImmunoCard[®] *H. pylori*

A Rapid Enzyme Immunoassay for the Detection of IgG to *Helicobacter pylori* in Human Serum, Plasma and Whole Blood

REF 710030

IVD

Rx Only

INTENDED USE

The ImmunoCard *H. pylori* enzyme immunoassay (EIA) is an in vitro qualitative procedure for the detection of IgG to *Helicobacter pylori* in human serum, plasma, or whole blood. Test results are intended to aid in the diagnosis of *H. pylori* infection. The ImmunoCard *H. pylori* assay may be performed in clinical laboratories and physician's offices.

SUMMARY AND EXPLANATION OF THE TEST

The importance of *Helicobacter pylori* in gastrointestinal diseases has increased greatly since Marshall and Warren described the presence of campylobacter-like organisms in the antral mucosa of patients with histological evidence of antrum gastritis and peptic ulcers especially duodenal ulcers.^{3,15} Though a causative relationship has not yet been fully established, the strong correlation between the presence of *H. pylori* and histologically confirmed gastritis, peptic ulcer disease and gastric carcinoma can no longer be debated.^{3,11,20}

The ecological niche in humans appears to be restricted to the stomach and duodenum. Patients who harbor the organism are divided into two basic groups. The first are those who are said to be "colonized". These patients have the organism, yet have no signs of gastrointestinal disease. Those with symptoms and *H. pylori* are considered to be infected. The process by which a colonized individual becomes infected remains unclear. The process by which patients become colonized is also still under investigation.^{3,7,11,13,14}

The diagnostics strategies for the determination of *H. pylori* have been developed among two lines. The first involves direct detection of the organism and, to date, is invasive. The second, a noninvasive approach, involves the detection of antibodies made against *H. pylori* or the detection of the organism's urease.^{8,9,22}

The invasive route requires that a biopsy be taken from the upper gastrointestinal tract. The presence of *H. pylori* is then confirmed by direct microscopic examination or culturing of the organism from the biopsy material. This strategy has the advantage of being able to detect active infections while being highly specific with a very high positive predictive value. The difficulties associated with this approach is that there is a risk and discomfort to the patient. In addition, *H. pylori* tends to colonize in patches and may be missed totally by the biopsy. Culturing of *H. pylori* from biopsy material tends to be difficult and time consuming. These technical difficulties may lead to false negative results.^{1,3,19}

The urea breath test is a noninvasive method that detects *H. pylori* by exploiting its highly active urease.^{3,8,10} The most common noninvasive approach to the detection of *H. pylori* is the serological identification of specific antibodies in infected patients. Through an indirect approach, the correlation between histological gastritis, the presence of *H. pylori*, and seropositivity is extremely strong. The disadvantages of these tests are that most require instruction to be interpreted and many have a lower specificity due to cross reactions from other organisms. Their major advantage is that they are rapid and offer a high sensitivity.^{5,18,21}

ImmunoCard *H. pylori* is an enzyme immunoassay in a self contained device which detects IgG antibodies to *H. pylori*. Results are interpreted by observing for a visual color change within seven minutes after addition of specimen.

BIOLOGICAL PRINCIPLES

The test detects IgG antibodies to *H. pylori* through the use of an enzyme immunoassay. Each test card consists of a plastic housing and a fiber membrane which provides the solid support for the EIA. There are two large sample ports and two small reaction ports. The left reaction port serves as a **CONTROL** and has human IgG antibodies immobilized on the fiber membrane. The right reaction port is the **TEST** port and contains immobilized *H. pylori* antigens.

Patient serum, plasma, or whole blood is added to each of the two sample ports. After allowing the sample to enter the device and migrate along the membrane to the reaction ports, anti-human IgG alkaline phosphatase conjugate is added to the reaction ports and allowed to enter the device. Wash Buffer and Substrate are then added to each of the reaction ports. Reaction ports are observed for the development of blue color.

The development of blue color in the **TEST** reaction port indicates a positive test result. The absence of blue color in the **TEST** reaction port denotes a negative test result. The **CONTROL** reaction port serves as the procedural control. Color development confirms that the reagents are performing and that the procedure was done correctly.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

- Test Cards** — Individually foil pouched cards containing immobilized *H. pylori* antigens (**TEST** reaction port) and human IgG (**CONTROL** reaction port).
- Positive Control** — Human Serum containing IgG antibodies to *H. pylori* in a buffer (pH 7.3) with 0.1% sodium azide.
- Negative Control** — Nonreactive Human Serum in a buffer (pH 7.3) with 0.1% sodium azide.
- Enzyme Conjugate** — Monoclonal anti-human IgG labeled with alkaline phosphatase in a buffer (pH 7.25) containing 0.1% sodium azide.
- Wash Buffer** — Buffer (pH 7.2) containing 9.5% (weight/vol.) guanidine hydrochloride.
- Substrate Reagent** — Buffer (pH 9.8) containing 5-bromo-4 chloro-3 indolyl phosphate and 0.1% sodium azide.
- Whole Blood Sample Diluent** — Buffer (pH 7.0) containing sodium dodecyl sulfate and 0.1% sodium azide.
- Transfer Pipettes**

MATERIALS NOT PROVIDED

- Timer
- Test tubes (for whole blood only)
- Capillary tubes (for whole blood only)

PRECAUTIONS

- All reagents are for in vitro diagnostic use only.
- Reagent concentrations, incubation times and temperature (22-27 C) have been optimized for sensitivity and specificity. Best results are obtained by adhering to these specifications. Once the assay has been started, complete all subsequent steps without interruption.**
- Patient specimens and controls may contain infectious agents. The positive and negative controls contain human sera which were screened for HBsAg and antibody to HIV and found to be negative. However, no test can offer complete assurance that human blood will not transmit HIV, hepatitis or other infectious agents. Therefore, patient specimens, controls, used Test Cards and other potentially hazardous material should be handled at Biosafety Level 2 as recommended in the CDC/NIH manual, "Biosafety in Microbiology and Biomedical Laboratories".
- All reagents and patient specimens should be at room temperature (22-27 C) before use. Additionally, all reagents should be removed from the foam insert and mixed prior to use.**
- Do not interchange reagents from different kit lot numbers or use expired reagents.
- Hold reagent vials vertically to insure proper drop size and delivery. Do not allow the tips of bottles or pipets to touch either the sample or reaction ports.
- Replace caps on correct vials.
- The Substrate Reagent may be light sensitive and should not be exposed to excessive illumination. If this reagent exhibits a blue color, it should be discarded.
- Use only one transfer pipet per specimen. Discard after use. Do not attempt to reuse.
- Discard any color reactions in the sample (lower) ports. Results are determined by color development in the reaction (upper) ports only.

WARNING: Some reagents in this kit contain sodium azide, which is a skin irritant. Avoid skin contact with the kit reagents. Disposal of reagents containing sodium azide into lead or copper plumbing can result in the formation of explosive metal azides. This can be avoided by flushing with a large volume of water during such disposal.

HAZARD and PRECAUTIONARY STATEMENTS


There are no known hazards associated with this product.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store at 2-8 C and return the kit promptly to the refrigerator after each use.

PROCEDURAL NOTES

- The ImmunoCard *H. pylori* Test Card format is diagrammed below:

- Batch processing of samples or controls is possible provided all procedural steps and incubation times are maintained for each Test Card.
- The **CONTROL** (left) side of each card provides a procedural control which detects problematic specimens and tests the reagent performance.

SPECIMEN COLLECTION AND PREPARATION

- Serum, plasma, or whole blood can be used in this assay. Sodium heparin, lithium heparin, sodium citrate and EDTA are acceptable anticoagulants. The performance of the ImmunoCard *H. pylori* assay with plasma obtained using other anticoagulants is unknown. Specimens must be collected using established laboratory procedures.
- Serum and plasma can be stored in the refrigerator at 2-8 C for up to five days. If testing is to be performed at a later date, the specimen should be frozen. Avoid repeated freezing and thawing. A cloudy serum/plasma specimen, or one containing particulate matter, may be centrifuged prior to testing.
- Whole blood with anticoagulant can be used without centrifugation. With lithium heparin, EDTA, or sodium citrate, whole blood specimens can be stored at 2-8 C for up to 72 hours.
- For fingerstick whole blood, completely fill capillary tube and test immediately.

TEST PROCEDURE


PLASMA, SERUM AND CONTROLS

- Remove the appropriate number of Test Cards from their envelopes. Label with appropriate identification. Use one card for each control or specimen to be tested. Test Cards must be maintained in a horizontal position during the performance of the assay procedure.
- Using a transfer pipet, dispense two drops of serum or plasma to both lower sample ports.
Optional Controls: The Positive Control is run by directly adding two drops to both lower sample ports of a card. The Negative Control can be run in the same manner on a separate card.
- Wait one minute after the sample or control has been added to both lower ports.
- Add two drops of Enzyme Conjugate to both upper ports.
- As soon as the Enzyme Conjugate is fully absorbed, add two drops of Wash Buffer to both upper ports.
- As soon as the Wash Buffer is fully absorbed, add two drops of Substrate Solution to both upper ports.
- Wait five minutes.**
- Without delay, visually read upper wells only for the presence or absence of blue color in the same order that the Substrate Solution was added.**

WHOLE BLOOD AND CONTROLS

- Remove the appropriate number of Test Cards from their envelopes. Label with appropriate identification. Use one card for each specimen to be tested. Test Cards must be maintained in a horizontal position during the performance of the assay procedure.
- Add four drops of Whole Blood Sample Diluent to labeled glass test tube.
- Using a transfer pipet, add two drops of whole blood to the test tube.
NOTE: For fingerstick procedure, completely fill a capillary tube with whole blood. Add the entire contents to a test tube that contains four drops of Whole Blood Sample Diluent.
Optional Controls: The Positive and Negative Controls are tested by adding two drops to the test tube which contains four drops of Whole Blood Sample Diluent.
- Using a new transfer pipet, dispense two drops of the diluted whole blood or controls to both lower sample ports.
- Wait three minutes after the sample or control has been added to both lower ports.
CAUTION: Some whole blood specimens may clog the membrane. If the sample does not completely wet the test and control reaction ports, the test is invalid.
- Add two drops of Enzyme Conjugate to both upper ports.
- As soon as the Enzyme Conjugate is fully absorbed, add two drops of Wash Buffer to both upper ports.
- As soon as the Wash Buffer is fully absorbed, add two drops of Substrate Solution to both upper ports.
- Wait five minutes.
- Without delay, visually read upper wells only for the presence or absence of blue color in the same order that the Substrate Solution was added.

INTERPRETATION OF RESULTS

READ ONLY UPPER REACTION PORTS TO INTERPRET RESULTS

Positive Test Result: Visually detectable blue color in **BOTH** reaction ports. Occasionally a positive test reaction may show evidence of a gradation of blue color from the bottom of the port to the top. This should be considered a positive result. A positive test result indicates the presence of IgG to *H. pylori*.

Negative Test Result: Visually detectable blue color in **CONTROL** reaction port (upper left) only. The **TEST** reaction port (upper right) should be colorless to faint grey. Occasionally, the **TEST** reaction port (upper right) may show evidence of a hint of blue color in the right or left side of the port, with the rest of the port remaining colorless. This should be considered a negative test result. A negative test result indicates that IgG to *H. pylori* is absent or below the level of detection of the assay.

Invalid Test Result: No detectable color in **CONTROL** reaction port. Invalid test results may be due to an ImmunoCard reagent problem, a procedural error, or clogging of the membrane by specimen: repeat test.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies. The Positive and Negative Controls should be assayed with each new kit lot or shipment and yield the following results:

1. Positive Control: Visually detectable blue color in both reaction (upper) ports (Figure 1).
2. Negative Control: Visually detectable blue color in CONTROL (upper left) reaction port only. The TEST (upper right) reaction port should be colorless (Figure 2).

The Procedural Control present in the upper left port of each Test Card tests the individual specimen and reagent performance. Failure of the Procedural Control to yield a blue color with any specimen or control reagent indicates an invalid test and the test should be repeated.

At the time of each use, kit components should be visually examined for obvious signs of microbial contamination, freezing, or leakage. It is suggested that the results of each quality control check be recorded in an appropriate log book to maintain high quality testing records. The Positive Control stresses the test at the assay cut-off. The Procedural Control present in the upper left portion of each Test Card verifies the reactivity of each reagent.

If the Procedural, Positive or Negative Control does not generate the appropriate and expected results, the test is invalid and the patient results should not be reported. If the expected reactions are not observed and the reagents are still within their expiration date, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

EXPECTED VALUES

Studies on the epidemiology of *H. pylori* have shown that this organism is present worldwide.^{12, 16} Gastritis caused by *H. pylori* has been shown to correlate with age, ethnic background, family size and socioeconomic class.^{6, 10} While the epidemiological studies on *H. pylori* infection or colonization are over expanding, the route by which individuals acquire this organism is, at present, unknown.^{3, 7, 11, 13, 14}

Studies indicate that the incidence of infection in the United States may increase 1-2% annually.¹⁰ It is not uncommon to see positivity rates of 50% in those who are 60 years or older. The frequency of *H. pylori* infection in patients diagnosed with duodenal ulcers, however, has been shown for patients in every age group to be approximately 80%.¹² The positivity rates seen with *H. pylori* clearly depend on the sample population and how that population is defined.

The ImmunoCard *H. pylori* test detects the presence of IgG antibodies to *H. pylori* in serum, plasma, and whole blood. Sera from 43.2% of a general symptomatic population were antibody positive by ImmunoCard *H. pylori*. There was a 99.1% correlation of whole blood and serum results in a general population of symptomatic and asymptomatic individuals. Expected values for a given population should be determined for each laboratory. The rate of positivity may vary depending on geographic location, method of specimen collection, handling and transportation, test employed, and general health environment of the patient population under study.

LIMITATIONS OF THE PROCEDURE

1. ImmunoCard *H. pylori* test results should be used in conjunction with information from the patient clinical evaluation and other diagnostic procedures.
2. This test should only be used on patients who present symptoms suggestive of gastrointestinal disease. ImmunoCard *H. pylori* is not intended for use with asymptomatic patients.
3. A Positive test result does not distinguish between an active infection and colonization by *H. pylori*.
4. ImmunoCard *H. pylori* is a qualitative test that only indicates the presence of IgG antibody to *H. pylori*. This test does not indicate that gastrointestinal disease is present.
5. A Negative test result indicates the absence of IgG antibody to *H. pylori* or that the IgG level is below that which can be detected by this test.
6. ImmunoCard *H. pylori* has not been tested in pediatric populations. Therefore, no performance characteristics have been established for this group.
7. ImmunoCard *H. pylori* was evaluated for possible interferences from high levels of hemoglobin, lipid, and bilirubin. No interference was observed in serum specimens that were either positive or negative for antibodies to *H. pylori*. Whole blood specimen results were not evaluated for interferences by these compounds.
8. Whole blood specimens with hematocrit levels > 61 volume % have not been evaluated and may clog the membrane.

SPECIFIC PERFORMANCE CHARACTERISTICS

The ImmunoCard *H. pylori* test was evaluated at two hospitals and one gastroenterology laboratory.⁴ A total of 382 serum specimens were obtained from patients with gastrointestinal complaints. All patients underwent endoscopy. The ImmunoCard *H. pylori* test was compared to results obtained from the biopsied material using rapid biopsy urease tests.

ImmunoCard <i>H. pylori</i>	Biopsy Urease Test	
	Positive	Negative
Positive	115	50
Negative	2	215

n= 382
Sensitivity = 98.2% (94.0-99.5)*
Specificity = 81.1% (76.0-85.4)
Agreement = 86.4% (82.5-89.5)

*95% confidence intervals calculated by the normal method are indicated in parenthesis

Forty-seven biopsy urease test negative/ImmunoCard *H. pylori* positive specimens were positive by a competitor's microtiter well EIA. Forty-six of these specimens were positive by histological staining and/or the pathology diagnosis was associated with *H. pylori* disease. The additional information indicates that 47 of the 50 relative false positive results were, in fact, true antibody positives that were not positive by biopsy urease test. The two urease test positive/ImmunoCard *H. pylori* negative specimens were antibody negative by a competitor's microtiter well EIA, indicating that these relative false negative results were true antibody negative results.

ImmunoCard *H. pylori* was compared to a competitor's microwell based EIA also for the detection of IgG to *H. pylori*.

ImmunoCard <i>H. pylori</i>	Competitor EIA		
	Positive	Negative	I
Positive	161	4	0
Negative	5	211	1

n= 382
Relative Sensitivity = 97.0% (93.1-98.7)*
Relative Specificity = 98.1% (95.3-99.3)
Relative Agreement = 97.6% (95.6-98.7)

One EIA negative/ImmunoCard *H. pylori* positive specimen was positive by the biopsy urease test. Three EIA negative/ImmunoCard *H. pylori* positive specimens were negative by the biopsy urease test. All five ImmunoCard *H. pylori* negative/EIA positive specimens were negative by the biopsy urease test. The one EIA indeterminate/ImmunoCard *H. pylori* negative specimen was biopsy urease test negative.

A specimen matrices study evaluated whole blood collected in lithium heparin, EDTA and sodium citrate and serum collected from 108 symptomatic and asymptomatic individuals. The relative comparison of whole blood to serum was not traceable to biopsy results. There was a 99.1% (95% CI = ± 1.8%) agreement between whole blood results, regardless of anticoagulant (30 positives) and serum results (31 positives).

A comparison of the results of capillary whole blood obtained by finger stick and serum from 77 individuals indicated 100% agreement; 15 positives and 62 negatives. The specimen results were not traceable to biopsy results.

REPRODUCIBILITY

Two clinical test sites ran a negative, low positive and high positive serum specimen 10 times a day each of three days. The Positive and Negative Controls were run three times a day over this testing period. Three clinical test sites ran the Positive and Negative Controls and four negative and four positive serum specimens three times a day for three days. A total of 450 assays were performed at all the testing sites. The inter and intra-assay reproducibility of ImmunoCard *H. pylori* with all specimens was 100%.⁴

ASSAY SPECIFICITY

A. Bacterial Studies

Absorption studies demonstrated that ImmunoCard *H. pylori* does not suffer cross reactivity problems with the following clinical isolates (CI) or ATCC strains.⁴


<i>Campylobacter coli</i>	CI	<i>Klebsiella pneumoniae</i>	AMC 13883
<i>Campylobacter fetus</i>	CI	<i>Proteus vulgaris</i>	ATCC 8427
<i>Campylobacter hyointestinalis</i>	ATCC 29543	<i>Pseudomonas aeruginosa</i>	ATCC 9027
<i>Campylobacter jejuni</i>	ATCC 33292	<i>Pseudomonas fluorescens</i>	CI
<i>Campylobacter lari</i>	CI	<i>Salmonella sp</i>	CI
<i>Campylobacter rectus</i>	ATCC 35217	<i>Serratia liquefaciens</i>	ATCC 3551
<i>Candida albicans</i>	ATCC 10231	<i>Shigella flexneri</i>	CI
<i>Citrobacter freundii</i>	ATCC 8090	<i>Shigella sonnei</i>	CI
<i>Enterobacter aerogenes</i>	ATCC 13048	<i>Staphylococcus aureus</i> (Cowan I)	ATCC 12598
<i>Enterobacter cloacae</i>	ATCC 13407	<i>Streptococcus faecalis</i>	CI
<i>Escherichia coli</i>	ATCC 8739	<i>Wolfinella succinogenes</i>	ATCC 29543
<i>Helicobacter cinaedi</i>	ATCC 35683	<i>Yersinia enterocolitica</i>	CI

B. Serum Studies

Sera from patients who were culture positive or had serologic evidence of the following diseases did not react when tested in ImmunoCard *H. pylori*.⁴

Anti-Nuclear Antibody	Cytomegalovirus	Rheumatoid Factor	Rubella
Lupus Erythematosus	Toxoplasma	Influenza A	Coccidioides
Histoplasma	Mononucleosis	Herpes Simplex Virus	Cryptococcus

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REFERENCES

- Alpert LC, Graham DY, Jr Evans DJ, Yoshimura HH, Hazell SaL, Evans DG, et al. Diagnostic possibilities for *Campylobacter pylori* infection. Eur J Gastroenterol. Hepatol, 1989; 1:17-26.
- Barthel JS, Everett ED. Diagnosis of *Campylobacter pylori* infections: The "gold standard" and the alternatives. Rev Infect Dis. 1990; 12:S107-S114.
- Buck GE. *Campylobacter pylori* and gastroduodenal disease. Clin Microbial Rev. 1990; 3:1-12.
- Data on file, Meridian Bioscience, Inc.
- Jr Evans DJ, Evans DG, Graham DY, Klein PD. A sensitive and specific serologic test for detection of *Campylobacter pylori* infection. Gastroenterology. 1989; 96:1004-1008.
- Fiedorek SC, Malaty HM, Evans DL, Pumphrey CL, Casteel HB, Jr Evans DJ, et al. Factors influencing the epidemiology of *Helicobacter pylori* infection in children. Pediatrics. 1991; 88:578-582.
- Graham DY. I. *Helicobacter pylori*: Its epidemiology and its role in duodenal ulcer disease. J Gastroenterol Hepatol. 1991; 6:105-113.
- Graham DY, Klein PD, Jr Evans DJ, Evans DG, Alpert LC, Opekun AR, et al. *Campylobacter pylori* detected noninvasively by the ¹³C-urea breath test. Lancet i: 1987; 1174-1177.
- Graham DY, Klein PD, Opekun AR, Boutton TW. Effect of age on the frequency of active *Campylobacter pylori* infection diagnosed by the [¹³C] Urea breath test in normal subjects and patients with peptic ulcer disease. J Infect Dis. 1988; 157:777-780.
- Graham DY, Malaty HM, Evans DG, Jr Evans DJ, Klein PD, Adam E. Epidemiology of *Helicobacter pylori* in an asymptomatic population in the United States. Effect of age, race, and socioeconomic status. Gastroenterology. 1991; 100:1495-1501.
- Gregson DB, Sidor AE. *Helicobacter pylori* and inflammatory disease of the stomach and duodenum. Stomach and Bowel. Sept. 1991; 29-30, 35-37.
- Loffeld RJLF, Stobberingh E, Van Spreuwel JP, Flendrig JA, Arends JW. The prevalence of anti-*Helicobacter (Campylobacter) pylori* antibodies in patients and healthy blood donors. J Med Microbiol. 1991; 32:105-109.
- Malaty HM, Graham DI, Klein PD, Evans DG, Adam E, Evans DJ. Transmission of *Helicobacter pylori* infection. Studies in families of healthy individuals. Scand J Gastroenterol. 1991; 26:927-932.
- Marshall BJ, Armstrong JA, McGeachie DB, Glancy RJ. Attempt to fulfill Koch's postulates for pyloric campylobacter. Med J Aust. 1985; 142:436-439.
- Marshall BJ, Warren JR. Unidentified curved bacilli in the stomach of patients with gastritis and peptic ulceration. Lancet i: 1984; 1311-1314.
- Morris A, Nicholson G, Lloyd G, Haines D, Rogers A, Taylor D. Seroprevalence of *Campylobacter pyloridis*. N Z Med J. 1986; 99:657-659.
- Murray DM. Clinical relevance of infection by *Helicobacter pylori*. Clin Microbial Lett. 1993; 15:33-37.
- Newell DG. Identification of the outer membrane proteins of *Campylobacter pyloridis* and antigenic cross-reactivity between C. pyloridis and C. jejuni. J Gen Microbiol. 1987; 133:163-170.
- Nichols L, Sughayer M, DeGirolami PC, Balogh K, Pleskow D, Eichelberger K, et al. Evaluation of diagnostic methods for *Helicobacter pylori* gastritis. Am J Clin Path. 1991; 95:769-773.
- Parsonnet J, Friedman GD, Vandersteen DP, Chang Y, Vogelman JH, Orentreich N, et al. *Helicobacter pylori* infection and the risk of gastric carcinoma. N Engl J Med. 1991; 325:1127-1131.
- Talley NJ, Newell DG, Ormand JE, Carpenter HA, Wilson WR, Zinsmeister AR, et al. Serodiagnosis of *Helicobacter pylori*. Comparison of enzyme-linked immunosorbent assays. J Clin Microbiol. 1991; 29:1635-1639.
- Westblom TU, Madan E, Kemp J, Subik MA. Evaluation of a rapid urease test to detect *Campylobacter pylori* infection. J Clin Microbiol 1988; 26:1393-1394.

INTERNATIONAL SYMBOL CHART

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo medico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
CE	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In-Vitro Diagnostica 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon échant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenverdünnung, in dem sich Probenverdünnungspuffer befindet
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Solamente per valutazione delle prestazioni / Reactifs IVD réservés à l'évaluation des performances / Solo para evaluación de funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contiene sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	SOLN STOP	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Límites de temperatura / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Coniugato enzimatico / Conjugué enzymatique / Conjugado enzimático / Enzymkonjugat
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	REAG	Reagent / Reagents / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
CONJ	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
SUBS	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
R. Only	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Per Receta / verschreibungspflichtig	DET REAG	Detection Reagent / Reagente Diretto / Réactif de Détection / Reactivo de Detección / Nachweis Reagent

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.

SN11233

REV. 12/16

 Manufactured By	Meridian Bioscience, Inc. USA/Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 Telephone: 513.271.3700 Orders/Customer Service: 800.543.1980 Technical Support Center: 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124
 Authorized Representative	Meridian Bioscience Europe S. r. L. Via dell'Industria, 7 20020 Villa Cortese, Milano ITALY Tel: +39 0331 43 36 36 Fax: +39 0331 43 36 16 Email: info@meridianbioscience.eu WEB: www.meridianbioscience.eu

Meridian Bioscience Europe s.a./n.v.
 2 Avenue du Japon - 1420 Braine l'Alleud
 BELGIUM
 Tel: +32 (0) 67 89 59 59
 Fax: +32 (0) 67 89 59 58
 Email: info.bn@meridianbioscience.eu

Meridian Bioscience Europe France
 34 rue de Ponthieu - 75008 Paris
 FRANCE
 Tel: +33 (0) 1 42 56 04 40
 Fax: +33 (0) 9 70 06 62 10
 Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
 Postbus 301 - 5460 AH Veghel
 NETHERLANDS
 Tel: +31 (0) 411 62 11 66
 Fax: +31 (0) 411 62 48 41
 Email: info.bn@meridianbioscience.eu