ImmunoCard[®] TOXINS A & B

A Rapid Enzyme Immunoassay for the Detection of *Clostridium difficile* Toxins A and B in Stool Specimens

REF 712050

IVD

Rx Only

INTENDED USE

ImmunoCard Toxins A & B is a rapid, qualitative, horizontal-flow enzyme immunoassay (EIA) for detecting *Clostridium difficile* Toxins A and B in human stool. This assay is used as an aid in the diagnosis of *C. difficile*-associated diseases.

In the US, ImmunoCard Toxins A & B is not intended for point-of-care use. The device is intended for moderately complex laboratories. In Canada, this device is not intended for point-of-care use.

SUMMARY AND EXPLANATION OF THE TEST

Toxigenic *Clostridium difficile* is the leading cause of nosocomial infectious diarrhea in developed countries. An estimated 300,000 cases of *C. difficile* associated disease (CDAD) are seen per year in U.S. hospitals alone.^{1,2} *Clostridium difficile* is the etiologic agent in approximately 25% of all cases of antibiotic-associated diarrhea. Virtually any antibiotic can predispose a patient to CDAD. The clinical presentation for CDAD ranges from asymptomatic colonization to life-threatening pseudomembranous colitis and toxic megacolon.² Most pathogenic strains of *C. difficile* produce two biologically and immunologically distinct toxins: toxin A (enterotoxin) and toxin B (cytotoxin). Toxin A was once thought to be responsible for most of the pathology seen in human CDAD until reports of clinically relevant disease caused by strains of *C. difficile* that produce only toxin B began to appear in the late 1990's.¹ The most accurate assay overall for the detection of *C. difficile* is the cytotoxin assay, however the method requires tissue culture facilities, 48-72 hour incubation and is not standardized.² The use of a rapid test such as ImmunoCard Toxins A & B enables the physician to verify infection quickly, begin proper treatment and to initiate enteric isolation precautions in a hospital setting.²

BIOLOGICAL PRINCIPLES

ImmunoCard Toxins A & B consists of a membrane held in a plastic frame with two sample ports and two reaction ports. The membrane carries immobilized antibodies to toxins A and B. The Enzyme Conjugate Reagent consists of antibodies to toxins A and B coupled to horseradish peroxidase. To perform the test, patient stool sample is diluted with Sample Diluent and Enzyme Conjugate and the mixture is incubated for five minutes to 24 hours. During the incubation, molecules of toxin, if present, are bound to the anti-toxin antibodies of the conjugate. Following incubation, an aliquot of the mixture is added to each of the two sample ports and the test is incubated for an additional five minutes at 20-26 C. During the second incubation the toxin-conjugate complex is separated from particulate matter as the fluid portion of the sample flows through the membrane to the TEST and CONTROL reaction ports. The toxin-conjugate complexes are then captured at the TEST reaction port by immobilized antitoxin in the reaction membrane. (The second of the two reaction ports serves as an internal control.) Both reaction ports are subsequently washed with ImmunoCard Wash Buffer I to reduce interference by contaminating proteins before ImmunoCard Substrate I is added. The reaction ports are incubated for an additional five minutes during which time the Enzyme Conjugate modifies the ImmunoCard Substrate I. The result is the appearance of a blue color. Reactions are read visually. Development of a blue color in the TEST reaction port indicates a positive test. In the CONTROL port, the anti-toxin antibodies of the conjugate bind directly to the immobilized toxin. The appearance of blue in the CONTROL reaction port indicates that sample was added, that reagents were active at the time of use and that proper sample migration occurred.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this kit is listed on the outer box.

- ImmunoCard Toxins A & B Test Card:** A membrane pad housed in a plastic frame and enclosed in a foil pouch with a desiccant. The pad carries immobilized monoclonal anti-toxin A and goat polyclonal anti-toxin B at the TEST reaction port and *C. difficile* toxin at the CONTROL reaction port. Store the test cards at 2-8 C when not in use.
- Sample Diluent*:** A buffered protein solution containing gentamicin and thimerosal (0.02%) as a preservative. Store at 2-8 C when not in use.
- Positive Control*:** Inactivated crude *C. difficile* toxin in a buffered solution containing thimerosal (0.02%) as a preservative. The reagent is supplied ready for use. Store at 2-8 C when not in use.
- Enzyme Conjugate*:** A blend of goat polyclonal antibodies to toxins A and B conjugated to horseradish peroxidase and suspended in a buffered protein solution containing gentamicin and thimerosal (0.02%). Store at 2-8 C when not in use.
- ImmunoCard Wash Buffer I*:** A buffered solution containing thimerosal (0.01%) as a preservative. Store at 2-8 C when not in use.
- ImmunoCard Substrate I*:** A buffered solution containing tetramethyl-benzidine and peroxide. Store at 2-8 C when not in use.
- Plastic transfer pipettes with measuring marks for 25 μ L and 150 μ L.
* In a plastic dropper vial. Use as supplied.

MATERIALS NOT PROVIDED

- Disposable latex gloves that should be used during the handling of the fecal samples as they are considered potentially biologically hazardous material
- Vortex for suspending of the stool sample in the Sample Diluent
- Interval timer
- Applicator sticks
- Small test tubes, such as 10 x 75 mm or 12 x 75 mm

PRECAUTIONS

- All reagents are for in vitro diagnostic use only.
- Do not store at temperatures above 2-8 C. Do not freeze.
- Do not deviate from the established insert method or falsely positive or falsely negative results may occur.
- Some patient samples contain infectious agents; therefore all patient samples should be handled and disposed of as if they are biologically hazardous.
- DO NOT interchange the Test Card, Enzyme Conjugate and Positive Control reagents from different kit lot numbers. The other reagents (Sample Diluent, ImmunoCard Wash Buffer I, ImmunoCard Substrate I) can be interchanged between kits providing the reagents are within their expiration periods. Do not use any reagent beyond its labeled expiration date.
- All stool samples must be mixed thoroughly before testing, regardless of consistency, to ensure a representative sample prior to testing.
- Failure to bring samples and reagents to room temperature (20-26 C) before testing may decrease assay sensitivity.
- Inspect Test Cards before removing the foil pouch. Do not use Test Cards that have holes in the foil pouch or where the pouch has not been completely sealed. Do not use Test Cards where the desiccant indicator has changed from blue to pink. False-negative reactions may result due to deterioration of the improperly stored Test Cards.
- Do not use the Sample Diluent or Positive Control if they are discolored or turbid. Discoloration or turbidity may be a sign of microbial contamination. It is normal for the ImmunoCard Wash Buffer I to appear slightly turbid.
- The ImmunoCard Substrate I is light sensitive and should be stored in the dark. The appearance of the reagent should be clear and colorless. If the reagent exhibits any color or is turbid it should be discarded.
- Hold reagent vials vertically when dispensing drops to ensure consistent drop size and delivery.

HAZARD and PRECAUTIONARY STATEMENTS


There are no known hazards associated with this product

SAMPLE COLLECTION AND PREPARATION

The stool sample should be transported in an airtight container and stored at 2-8 C until tested. The sample should be tested as soon as possible, but may be held up to four days at 2-8 C prior to testing. If testing cannot be performed within this time frame, samples should be frozen immediately on receipt and stored frozen (≤ -20 C) until tested. Samples that are frozen and thawed up to three times are acceptable for use.

PROCEDURAL NOTES


- The ImmunoCard test card format is diagrammed below:

- Batch processing of any number of samples or controls is possible provided that, for each card, the appropriate steps, sequence of reagent addition, incubation (wait) and result reading times are maintained. Each procedural step should be completed with each sample before the next step is started.
- The Control (left) side of each card provides an internal control for each sample. This tests the individual sample for sufficient dilution and reagent performance.


TEST PROCEDURE

- Bring all test cards, reagents and samples to room temperature (20-26 C) before testing. Remove the reagents from the kit box to warm. **Reagents may take up to 60 minutes to warm following refrigeration. RETURN TO 2-8 C IMMEDIATELY AFTER USE.**
- Label 1 test tube for each patient sample to be tested.
- Add 200 μ L of Sample Diluent to a test tube using the dropper in the bottle.
- Add 3 drops of Enzyme Conjugate (150 μ L) to each tube.
- Immediately add stool sample or control as follows:
 - Liquid/Semi-solid stools** – Mix the stool sample thoroughly. Using a clean transfer pipette (included with the kit), draw the stool to the first mark from the end of the tip (25 μ L). Dispense the stool into the Diluent/Conjugate mixture. Use the same pipette to mix the stool by gently aspirating, then dispensing the suspension several times. Vortex the final mixture for 10 seconds. Leave the pipette in the tube for use later in the test.*
 - Formed/Solid stools** – Mix the stool sample thoroughly. Using a wooden applicator stick, transfer a small portion (approximately 2 mm diameter) of the stool into the Diluent/Conjugate mixture. Emulsify the sample using the applicator stick, and then thoroughly mix with a vortex for 10 seconds. Add a pipette to the tube and leave for later use in testing.*
 - Optional control:** Using the container droppers, add 1 drop of Positive Control or Negative Control (Sample Diluent) (provided with the kit) to the Diluent/Conjugate mixture. Vortex 10 seconds. Add a pipette to the tube and leave for later testing.

- *NOTE: The transfer of too little sample, or failure to mix and suspend the sample in the Diluent/Conjugate mixture completely may result in a false-negative test result. Addition of too much stool sample may cause invalid results due to restricted sample flow.
- Let each diluted sample stand for 5 minutes to 24 hours at 20-26 C. Seal tube containing diluted sample or control if incubating beyond 6 minutes.
 - Use 1 ImmunoCard Test Card for each sample or control. When ready to perform testing, remove the ImmunoCard Test Card from its foil pouch. Discard the pouch and desiccant. Label the device with the name of the patient or the control.
 - Vortex each sample or control for 10 seconds before use.
 - Add 150 μ L of each diluted sample or control to each of the 2 sample ports on 1 test card (150 μ L represents the second mark on the transfer pipette).
 - Incubate for 5 \pm 1 minutes at 20-26 C. NOTE: At the end of incubation both reaction ports must appear completely wet. If either reaction port is not completely wet, discard the Test Card and repeat the procedure. Proceed with testing if sample remains in the sample ports yet both reaction ports are completely wetted.
 - Hold the ImmunoCard Wash Buffer I vial vertically and dispense exactly 3 drops to each of the reaction ports.
 - When the ImmunoCard Wash Buffer I has been completely absorbed, hold the vial vertically and add exactly 3 drops of ImmunoCard Substrate I.
 - Incubate the test card at 20-26 C for 5 minutes.
 - Visually read the results of each card within 30 seconds of the end of incubation.**

INTERPRETATION OF RESULTS

NOTE: Disregard color development in the sample ports. The color of the sample ports is not evaluated during testing.

Positive test result: Blue color in the TEST (upper right) and CONTROL (upper left) reaction ports. The intensity of the blue color of the TEST port may vary from the bottom to the top of port. A positive test result indicates that toxin A and/or B is/are present in the sample. PLEASE NOTE: The TEST port may appear less blue than the CONTROL port.

Negative test result: Blue color in the CONTROL reaction (upper left) port only. The TEST reaction (upper right) port should be colorless. (The wetting of the membrane may make the TEST port appear to be slightly grey.) A negative test result indicates that *C. difficile* toxins are absent or below the limit of detection of the assay. Occasionally the TEST PORT (upper right) may show evidence of a hint of blue color in the left or right side of the port, with the rest of the port colorless. This should be considered a negative test result.

Invalid test results:

- No detectable blue color in the CONTROL reaction port. Failure of the CONTROL invalidates any test result. The invalid test may be due to the failure of a reagent or the Test Card to perform properly, failure to add sample, failure of the sample to migrate, or failure to dilute the sample correctly leading to over inoculation of the Test Card. Samples with extremely high levels of toxins may also present with a positive TEST yet negative CONTROL port.

- A blue ring on the plastic frame surrounding the TEST (upper right) port during the test procedure. When invalid results are obtained, redilute the sample and repeat the test on a new card.

ITALIANO

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies. Each Test Card contains an internal control. The appearance of a blue color in the CONTROL reaction port verifies the assay was active at the time of use, that sample was added and that there was adequate migration of the sample.

The reactivity of each new lot and each new shipment of the ImmunoCard Toxin A & B should be verified on receipt or before use. Positive Control reagent is supplied with the kit. The control is used to verify the reactivity of the other reagents associated with the assay, and is not intended to ensure precision at the analytical assay cut-off. Sample Diluent is used for the Negative Control. Additional tests can be performed with the controls to meet the requirements of local, state and/or federal regulations and/or accrediting organizations. The results expected with the controls are described in the section on INTERPRETATION OF RESULTS. The Test Card should not be used if control tests do not produce the correct results. Proper results obtained with the Control Port, Positive Control and Negative Control (Sample Diluent) serve as indicators that the test was performed correctly, that the antibodies embedded in the membrane and the Enzyme Conjugate are active at the time of testing, and that the membrane supports proper sample flow. Failure of the internal and external control to produce the expected results suggests the test was not performed correctly (ie, incorrect volume of reagents added; incorrect incubation temperature or times used or that reagents were not brought to room temperature prior to testing). If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

EXPECTED VALUES

The frequency of antibiotic-associated diarrhea caused by *C. difficile* is dependent on several factors, including patient population, type of institution and epidemiology. The reported incidence of *C. difficile*-associated disease in patients suspected of having antibiotic-associated diseases is 15-20%, although different facilities may find positive rates above or below this range.^{2,3,5} A positive result in an uninfected patient is not expected.

LIMITATIONS OF THE PROCEDURE

- The test is qualitative and no quantitative interpretation should be made with respect to the intensity of the color when reporting the result.
- A positive test confirms the presence of toxins A and/or B only. A physician must use the test results in conjunction with other diagnostic procedures and the patient's clinical condition to establish a diagnosis of *C. difficile*-associated disease. The levels of toxin do not always correlate directly with the presence or severity of disease.¹
- Isolates of *C. sordellii* may react with this test due to the immunological identity of the *C. sordellii* toxins.³
- The performance of samples from pediatric patients has not been evaluated.
- Two distinct groups have been identified that can harbor *C. difficile* asymptomatically at very high rates. Colonization rates of up to 50% and higher have been reported in infants² and rates of up to 32% in cystic fibrosis patients.⁴
- While relatively stable at 2-8 C, *C. difficile* toxins -- particularly at low concentrations -- easily degrade at room temperatures. The rate at which the toxins degrade differs from patient sample to patient sample and therefore the rate cannot be predicted. For this reason, best practice requires that sample be refrigerated or frozen within two hours of collection and the samples tested within the timeframes recommended in this insert. Do not accept samples that have not been collected, handled or transported properly.

SPECIFIC PERFORMANCE CHARACTERISTICS

Comparative studies

Three independent laboratories tested samples with ImmunoCard Toxins A & B and the standard cell cytotoxicity assay with neutralization. The results of the parallel tests are given below.

Prospective Samples	ICTAB			95% CI
	Pos	Neg	Total	
Cytotoxin Pos (Std)	67	5	72	
Cytotoxin Neg (Std)	2	176	178	
Total	69	181	250	
Clinical sensitivity	67/72 (93.1%)			87.1-98.9%
Clinical specificity	176/178 (98.9%)			97.6-100%
Predictive value positive test	67/69 (97.1%)			92.9-100%
Predictive value negative test	176/181 (97.2%)			94.5-99.5%
Retrospective Samples	ICTAB			95% CI
	Pos	Neg	Total	
Cytotoxin Pos (Std)	33	0	33	
Cytotoxin Neg (Std)	5	303	308	
Total	38	303	341	
Clinical sensitivity	33/33 (100%)			N/A
Clinical specificity	303/308 (98.4%)			96.4-99.6%
Predictive value positive test	33/38 (86.8%)			76.2-97.8%
Predictive value negative test	303/303 (100%)			N/A

ANALYTICAL SENSITIVITY

The lower limit of detection of this assay is approximately 3 ng of toxin A and 3 ng of toxin B per mL of stool. This limit does not vary from solid to liquid/semi-solid stool.

REPRODUCIBILITY

The reproducibility of ImmunoCard Toxins A & B was determined with known negative (n=2) and positive (n=6) samples that were coded and randomly sorted to prevent their identification during testing. Two of the positive samples were near the limit of detection for the assay. Three laboratories tested the samples on three consecutive days. The samples produced the expected results 100% of the time.

ANALYTICAL SPECIFICITY

Stool samples inoculated with the following microbial agents (to a final sample concentration of $\geq 1 \times 10^8$ organisms/mL) do not react with ImmunoCard Toxins A & B: Adenovirus 40, Adenovirus 41, *Aeromonas hydrophila*, *Bacillus cereus*, *Bacillus subtilis*, *Bacteroides fragilis*, *Campylobacter coli*, *Campylobacter jejuni*, *Candida albicans*, *Clostridium butyricum*, nontoxigenic *C. difficile* (10 strains), *Clostridium perfringens*, *Clostridium septicum*, *Clostridium sordellii* strain 9714, *Clostridium sporogenes*, *Enterobacter cloacae*, *Escherichia coli* (four strains including two O157:H7 strains), *Helicobacter pylori*, *Klebsiella pneumoniae*, *Peptostreptococcus anaerobius*, *Porphyromonas asaccharolytica*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, Rotavirus, *Salmonella typhimurium*, *Serratia liquefaciens*, *Shigella dysenteriae*, *Shigella flexneri*, *Shigella sonnei*, *Staphylococcus aureus*, *Staphylococcus aureus* Cowan I, *Staphylococcus epidermidis*, *Streptococcus faecalis*, *Vibrio parahaemolyticus*, *Yersinia enterocolitica*. The only non-*C. difficile* microorganism reactive with ImmunoCard Toxins A & B was *C. sordellii* VPI 9048. This strain produces toxin A and B homologs HT and LT.

TESTS FOR INTERFERING SUBSTANCES

The following substances had no effect on test results when present in stool in the concentrations indicated: fecal fat (4.8% w/v), whole blood (40% v/v), mucin (3.5% w/v), Metronidazole (0.25% w/v), Vancomycin (0.25% w/v), barium sulfate (5% w/v), Imodium AD® (5% w/v), Kaopectate® Caplets (5 mg/mL), Pepto Bismol® (5% v/v).

Additional information regarding the performance characteristics of this product can be obtained by contacting Meridian's Technical Support Department at (800) 343-3858 for assistance (USA) or your local distributor.

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REFERENCES

- Wilkins T, Lyerly D. *Clostridium difficile* testing: After 20 years, still challenging. J Clin Microbiol 2003;41: 531-534.
- Mylonakis E, Ryan E, Calderwood S. *Clostridium difficile*-associated diarrhea. Arch Int Med 2001;161: 525-533.
- Lyerly DM, Krivan, HC, Wilkins, TD. *Clostridium difficile*. Its disease and toxins. Clin Microbiol Rev 1988;1:1-18.
- Pokorny CS, Bye PTP, MacLeod C, Selby WS. Antibiotic-associated colitis and cystic fibrosis. Dig Diseases Sci 1992;37:1464-1468.
- Bartlett JG. *Clostridium difficile*: Clinical consideration. Rev Infect Dis 1990;12:S243-51.

SN11223

REV. 12/16

<p>Manufactured By</p>	<p>Meridian Bioscience, Inc. USA/Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 Telephone: 513.271.3700 Orders/Customer Service: 800.543.1980 Technical Support Center: 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124</p>
	<p>Meridian Bioscience Europe S. r. l. Via dell'Industria, 7 20020 Villa Cortese, Milano ITALY Tel: +39 0331 43 36 36 Fax: +39 0331 43 36 16 Email: info@meridianbioscience.eu WEB: www.meridianbioscience.eu</p>

Meridian Bioscience Europe s.a./n.v.
2 Avenue du Japon - 1420 Braine l'Alleud
BELGIUM
Tel: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
Email: info.bn@meridianbioscience.eu

Meridian Bioscience Europe France
34 rue de Ponthieu - 75008 Paris
FRANCE
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 9 70 06 62 10
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
Postbus 301 - 5460 AH Veghel
NETHERLANDS
Tel: +31 (0) 411 62 11 66
Fax: +31 (0) 411 62 48 41
Email: info.bn@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de símbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo medico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenverdünnung, in dem sich Probenverdünnungspuffer befindet
	Catalogue number / Numero di catalogo / Référence ou catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Solamente per valutazione delle prestazioni / Reagents IVD reserves à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per <n> saggi / Contenu suffisant pour <n> tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	SOLN STOP	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stoplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Coniugato enzimatico / Coniugado enzimático / Enzymkonjugat
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	REAG	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenza / Mise En Garde / Advertencia / Warnhinweise
CONJ	Conjugate / Coniugato / Coniugado / Conjugado / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
SUBS	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentrierter Waschpufferkonzentrat
R. Only	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig	DET REAG	Detection Reagent / Reagente Diretto / Réactif de Détection / Reactivo de Detección / Nachweis Reagent

For technical assistance, please call our Technical Support Center at 800-343-3858 (USA) between the hours of 8AM and 6PM, Eastern Standard Time. To place an order, please call our Customer Service Department at 800-543-1980 (USA).