## A Rapid Enzyme Immunoassay for the Detection of *Clostridium difficile* Antigens in Stool

REF 716050

IVD

Rx Only

### INTENDED USE

ImmunoCard *C. difficile* GDH is a rapid qualitative enzyme immunoassay screening test to detect *Clostridium difficile* antigen, glutamate dehydrogenase, in fecal specimens from persons suspected of having *C. difficile* infection (CDI). This test does not distinguish between toxigenic and non-toxicogenic strains of *C. difficile*. Samples from patients that produce positive results with this test must be further tested with an assay designed to detect toxigenic *C. difficile* strains and assist with the diagnosis of CDI.

### SUMMARY AND EXPLANATION OF THE TEST

*Clostridium difficile* has become the most important etiologic agent in nosocomial diarrhea worldwide. *C. difficile* infection (CDI) can range from mild diarrhea to life-threatening toxic megacolon. Classic risk factors for nosocomial CDI include antibiotic use and old age. Recently however CDI has been seen more frequently outside the hospital environment without the typical risk factors historically associated with *C. difficile* infection.<sup>1</sup> All strains of *C. difficile* produce glutamate dehydrogenase (GDH), also known as common antigen.<sup>2,3</sup> Some strains of *C. difficile* produce two toxins (toxin A and toxin B) that have both been associated with disease. Strains of *C. difficile* that do not harbor the toxin A and toxin B genes are considered to be non-pathogenic. Non-molecular laboratory approaches to *C. difficile* detection can be divided into two general groups: toxin-based and non-toxin based. Toxin-based assays include enzyme immunoassay (EIA) for toxin A or toxins A and B and cytotoxin, which primarily detects toxin B. Non-toxin based assays such as EIA for *C. difficile* GDH, while inherently more sensitive, must be followed up with testing for toxigenic *C. difficile*.

### BIOLOGICAL PRINCIPLES

ImmunoCard *C. difficile* GDH consists of a membrane held in a plastic frame with two sample ports and two reaction ports. The membrane carries immobilized antibodies to glutamate dehydrogenase. The Enzyme Conjugate Reagent consists of antibodies to glutamate dehydrogenase coupled to horseradish peroxidase. To perform the test, patient stool sample is diluted with Sample Diluent and Enzyme Conjugate and the mixture is incubated for 15 minutes. During the incubation, GDH antigen, if present, is bound to the anti-GDH antibodies of the conjugate. Following incubation, an aliquot of the mixture is added to each of the two sample ports and the test is incubated for an additional five minutes at 20-26 C. During the second incubation, the GDH-conjugate complex is separated from particulate matter as the fluid portion of the sample flows through the membrane to the TEST and CONTROL reaction ports. The GDH-conjugate complexes are then captured at the TEST reaction port by immobilized antibodies in the reaction membrane. The second of the two reaction ports serves as an internal control. Both reaction ports are subsequently washed with Wash Reagent to reduce interference by contaminating proteins before Substrate Reagent is added. The reaction ports are incubated for an additional five minutes during which time the Enzyme Conjugate modifies the Substrate Reagent. The result is the appearance of a blue color. Reactions are read visually. Development of a blue color in the TEST reaction port indicates a positive test. In the CONTROL port, the anti-GDH antibodies of the conjugate bind directly to the immobilized antibodies. The appearance of blue in the CONTROL reaction port indicates that sample was added, that reagents were active at the time of use and that proper sample migration occurred.

### REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this kit is listed on the outer box.

- ImmunoCard *C. difficile* GDH Test Cards:** A membrane pad housed in a plastic frame and enclosed in a foil pouch with a desiccant. The pad carries immobilized polyclonal anti-GDH antibodies at the TEST reaction port and Goat anti-mouse antibodies at the CONTROL reaction port. Store the test cards at 2-8 C when not in use.
- ImmunoCard *C. difficile* GDH Positive Control\*:** *C. difficile* GDH in a buffered protein solution containing 0.1% ProClin® and 0.03% gentamicin as preservatives. Store at 2-8 C when not in use.
- ImmunoCard *C. difficile* GDH Sample Diluent/Negative Control\*:** A buffered protein solution containing 0.03% gentamicin and 0.02% thimerosal as preservatives. Store at 2-8 C when not in use.
- ImmunoCard *C. difficile* GDH Enzyme Conjugate\*:** Monoclonal anti-GDH antibody conjugated to horseradish peroxidase and suspended in a buffered protein solution containing gentamicin and thimerosal (0.02%). Store at 2-8 C when not in use.
- ImmunoCard Wash Buffer I\*:** A buffered solution containing thimerosal (0.01%) as a preservative. Store at 2-8 C when not in use.
- ImmunoCard Substrate I\*:** A buffered solution containing tetramethylbenzidine and peroxide. Store at 2-8 C when not in use.
- Plastic transfer pipettes with measuring marks for 25 µL and 150 µL.**

\* In a plastic dropper vial. Use as supplied.

### MATERIALS NOT PROVIDED

- Disposable latex gloves should be used during the handling of the fecal samples as they are considered potentially biologically hazardous material.
- Vortex for suspending of the stool sample in the Sample Diluent.
- Interval timer
- Wooden applicator sticks
- Small test tubes, such as 10 x 75 mm or 12 x 75 mm

### PRECAUTIONS

- All reagents are for in vitro diagnostic use.
- Do not store at temperatures above 2-8 C. Do not freeze.
- Do not deviate from the established insert method or falsely positive or falsely negative results may occur.
- Some patient samples contain infectious agents; therefore all patient samples should be handled and disposed of as if they are biologically hazardous.
- DO NOT interchange the Test Card, Enzyme Conjugate and Positive Control reagents from different kit lot numbers. The other reagents (Sample Diluent/Negative Control, ImmunoCard Wash Buffer I, ImmunoCard Substrate I) can be interchanged between kits providing the reagents are within their expiration periods. Do not use any reagent beyond its labeled expiration date.
- All stool samples must be mixed thoroughly before testing, regardless of consistency, to ensure a representative sample prior to testing.
- Failure to bring samples and reagents to room temperature (20-26 C) before testing may decrease assay sensitivity.
- Inspect Test Cards before removing the foil pouch. Do not use Test Cards that have holes in the foil pouch or where the pouch has not been completely sealed. Do not use Test Cards where the desiccant indicator has changed from blue to pink. False-negative reactions may result due to deterioration of the improperly stored Test Cards.
- Do not use the Sample Diluent/Negative Control or Positive Control if they are discolored or turbid. Discoloration or turbidity may be a sign of microbial contamination. It is normal for the ImmunoCard Wash Buffer I to appear slightly turbid.
- ImmunoCard Substrate I is light sensitive and should be stored in the dark. The appearance of the reagent should be clear and colorless. If the reagent exhibits any color or is turbid it should be discarded.
- Hold reagent vials vertically when dispensing drops to ensure consistent drop size and delivery.

### WARNINGS

- Patient specimens may contain infectious agents and should be handled and disposed of as potential biohazards.
- Dispose of all used test materials in an appropriate container. Treat waste as a potential biohazard.

### HAZARD and PRECAUTIONARY STATEMENTS


Refer to the SDS, available at [www.meridianbioscience.com](http://www.meridianbioscience.com) for Hazards and Precautionary Statements.

### SHELF LIFE AND STORAGE

The kit expiration date is indicated on the kit label. Store the kit at 2-8 C and return the kit promptly to the refrigerator after each use.

### PROCEDURAL NOTES

- The ImmunoCard test card and pipette are diagrammed below:

- Batch processing of any number of samples or controls is possible provided that, for each card, the appropriate steps, sequence of reagent addition, incubation (wait) and result reading times are maintained. Each procedural step should be completed with each sample before the next step is started.
- The Control (left) side of each card provides an internal control for each sample. This tests the individual sample for sufficient dilution and reagent performance.

### REAGENT PREPARATION

Bring the entire kit, including test card pouch, to 20-26 C before use. No additional reagent preparation is required.


### SPECIMEN COLLECTION AND PREPARATION

The stool sample should be transported in an airtight container and stored at 2-8 C until tested. The sample should be tested as soon as possible, but may be held up to five days at 2-8 C prior to testing. Samples that will not be tested within five days should be frozen immediately on receipt and stored frozen ( $\leq -20$  C) until tested. Sample may be held for up to 30 days and may be frozen and thawed twice.

### TEST PROCEDURE

- Bring all test cards, reagents and samples to room temperature (20-26 C) before testing. Remove the reagents from the kit box to warm. Reagents may take up to 60 minutes to warm following refrigeration. RETURN TO 2-8 C IMMEDIATELY AFTER USE.
- Label 1 test tube for each patient sample to be tested.
- Add 200 µL of Sample Diluent/Negative Control to a test tube using the dropper in the bottle.
- Immediately add stool sample or control as follows:
  - Liquid/Semi-solid stools** – Mix the stool sample thoroughly. Using a clean transfer pipette, (included with the kit), draw the stool to the first mark from the end of the tip (25 µL). Dispense the stool into the Diluent. Use the same pipette to mix the stool by gently aspirating, then dispensing the suspension several times. Vortex the mixture for 10 seconds. Leave the pipette in the tube for use later in the test.\*
  - Formed/Solid stools** – Mix the stool sample thoroughly. Using a wooden applicator stick, transfer a small portion (approximately 2 mm diameter) of the stool into the Diluent. Emulsify the sample using the applicator stick, and then thoroughly mix with a vortex for 10 seconds. Add a pipette to the tube and leave for later using in testing.\*
  - Optional control:** Using the container droppers, add 1 drop of Positive Control or Sample Diluent/Negative Control to the Diluent. Vortex for 10 seconds. Add a pipette to the tube and leave for later testing.
- \*NOTE: The transfer of too little sample, or failure to mix and suspend the sample in the Diluent completely may result in a false-negative test result. Addition of too much stool sample may cause invalid results due to restricted sample flow.
- Add 3 drops of ImmunoCard *C. difficile* GDH Enzyme Conjugate (150 µL) to each tube. Vortex the mixture for 10 seconds. Leave the pipette in the tube for use later in the test.
- Let each diluted sample stand for 15 minutes at 20-26 C.
- Use 1 ImmunoCard Test Card for each sample or control. When ready to perform testing, remove the ImmunoCard Test Card from its foil pouch. Discard the pouch and desiccant. Label the device with the name of the patient or the control.
- Vortex each sample or control for 10 seconds before use.
- Add 150 µL of each diluted sample or control to each of the 2 sample ports on 1 test card (150 µL represents the second mark on the transfer pipette).
- Incubate for 5 minutes at 20-26 C.
- NOTE: At the end of incubation both reaction ports must appear completely wet. If either reaction port is not completely wet, discard the Test Card and repeat the procedure. Proceed with testing if sample remains in the sample ports yet both reaction ports are completely wetted.
- Hold the ImmunoCard Wash Buffer I vial vertically and dispense exactly 3 drops to each of the reaction ports.
- When the ImmunoCard Wash Buffer I has completely absorbed, hold the vial vertically and add exactly 3 drops of ImmunoCard Substrate I reagent.
- Incubate the test card at 20-26 C for 5 minutes.
- Visually read the results of each card within 30 seconds of the end of incubation.

### INTERPRETATION OF RESULTS

READ ONLY UPPER REACTION PORTS TO INTERPRET RESULTS

NOTE: Disregard color development in the sample ports. The color of the sample ports is not evaluated during testing.

**Positive test result:** Blue color in the TEST (upper right) and CONTROL (upper left) reaction ports. The intensity of the blue color of the TEST port may vary from the bottom to the top of port. A positive test result indicates that *C. difficile* GDH is present in the sample. PLEASE NOTE: The TEST port may appear less blue than the CONTROL port.

**Negative test result:** Blue color in the CONTROL reaction (upper left) port only. The TEST reaction (upper right) port should be colorless. The wetting of the membrane may make the TEST port appear to be slightly grey. A negative test result indicates that *C. difficile* GDH is absent or below the limit of detection of the assay. Occasionally the TEST PORT (upper right) may show evidence of a hint of blue color in the left or right side of the port, with the rest of the port colorless. This should be considered a negative test result.

### Invalid test results:

- No detectable blue color in the CONTROL reaction port. Failure of the CONTROL invalidates any test result. The invalid test may be due to the failure of a reagent or the Test Card to perform properly, failure to add sample, failure of the sample to migrate, or failure to dilute the sample correctly leading to over inoculation of the Test Card. Samples with extremely high levels of GDH may also present with a positive TEST yet negative CONTROL port.
- A blue ring on the plastic frame surrounding the TEST (upper right) port during the test procedure. When invalid results are obtained, redilute the sample and repeat the test on a new card.

### QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

- At the time of each use, kit components should be visually examined for obvious signs of microbial contamination, freezing or leakage. Do not use contaminated or suspect reagents.
- The performance of each new kit lot or each new shipment of the ImmunoCard C. difficile GDH kit should be verified on receipt or before use. See the INTERPRETATION OF RESULTS section above for a description of the expected results for control reagents. Tests should be considered invalid when either control reagent does not produce the expected results. In such cases, repeat tests and controls. If, on repeat testing, the expected reactions are still not observed with in-date reagents, call Meridian's Technical Support Services at 800-343-3858 (USA) or your authorized Meridian distributor for assistance.
- The controls are used to monitor reagent reactivity. Failure of the controls to produce the expected results can mean that one or more of the reagents are no longer reactive at the time of use, the test was not performed correctly, or that reagents or samples were not added. If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated, please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
- Each Test Card contains an internal control. The appearance of a blue color in the CONTROL reaction port verifies the assay was active at the time of use, that sample was added and that there was adequate migration of the sample.
- Specimen matrix interference has not been observed in this assay as samples are significantly diluted before testing in Sample Diluent. For this reason, the Positive and Negative Controls supplied as part of this assay are prepared in matrices similar to the Sample Diluent. If control materials that are identical in composition to test samples are preferred, the user can prepare these by diluting known positive and negative specimens in Sample Diluent according to the SPECIMEN PREPARATION section of this insert.

### EXPECTED VALUES

The frequency of antibiotic-associated diarrhea caused by *C. difficile* is dependent on several factors including: patient population, type of institution and epidemiology. The reported incidence of *C. difficile* infection in patients suspected of having antibiotic-associated diarrhea is 15-25%<sup>4</sup> although different facilities may find positivity rates outside this range. In clinical trials conducted for this assay, samples were collected from 975 male (45.7%) and female (54.3%) patients. Pediatric patients ranged from 14 days to 21 years. The prevalence in the populations in this clinical trial was 13.6%.

### LIMITATIONS OF THE PROCEDURE

- The test is qualitative and no quantitative interpretation should be made with respect to the intensity of the color when reporting the result.
- A positive result does not define the existence of *C. difficile* disease or the presence of toxigenic *C. difficile*, but only indicates the detection of the organism. All ImmunoCard C. difficile GDH positive stools should be tested by another method to verify the presence of toxigenic *C. difficile*.
- Staphylococcus aureus* (Cowan strain I) was found to be crossreactive with ImmunoCard C. difficile GDH.
- Samples containing whole blood at concentrations greater than 40% have the potential to produce false negative results.
- SHEA/IDSA Guideline defines the best strategy for *C. difficile* testing as testing performed only on diarrheal (unformed) stool, unless ileus due to *C. difficile* is expected. The proper specimen for the diagnosis of *C. difficile* infection is a watery, loose or unformed stool. Because 10% or more of hospitalized patients may be colonized with *C. difficile*, evaluating a formed stool can decrease the specificity of diagnosis of CDAD.<sup>5</sup>
- There is insufficient data to evaluate the performance of this device with specimens from patients 0-28 days old and 18-21 years old.

### SPECIFIC PERFORMANCE CHARACTERISTICS

ImmunoCard C. difficile GDH was evaluated by five independent test sites located in the Western, Southwestern, Southeastern, and Midwestern regions of the United States. A total of 975 qualified patient samples were evaluated; all specimens were prospectively collected samples. Samples were collected from males (46%) and females (54%) and categorized as solid (17.1%), semi-solid (50.8%), bloody (0.5%) and watery (31.5%). The age groups of the patients ranged from 14 days to 111 years of age. No differences in test performance were observed based on patient age or gender. The following tables show the assay performance by clinical site and patient age.

Table 1. Performance characteristics by clinical site

Site	Positive Samples			Negative Samples		
	ImmunoCard C. difficile GDH/Culture	Sensitivity %	95% CI	ImmunoCard C. difficile GDH/Culture	Specificity % <sup>1</sup>	95% CI
Total Sites	124/127	97.6% <sup>1</sup>	93.3 – 99.2%	738/848	87.0% <sup>2</sup>	84.6 – 90.1%
Site 1	32/32	100.0%	89.3 - 100.0%	165/205	80.5%	74.5 – 85.3%
Site 2	39/41	95.1%	83.9 – 98.7%	257/290	88.6%	84.4 – 91.8%
Site 3	35/36	97.2%	85.8 – 99.5%	140/168	83.3%	77.0 – 88.2%
Site 4	8/8	100.0%	67.6 – 100.0%	76/80	95.0%	87.8 – 98.0%
Site 5	10/10	100.0%	72.2 – 100.0%	100/105	95.2%	89.3 – 97.9%

<sup>1</sup> Three of 3 false negatives tested negative by another FDA cleared GDH assay.

<sup>2</sup> Seventy-four of the 110 false positives tested positive by another FDA cleared GDH assay.

<sup>3</sup> Refer to Table 2 footnote.

Table 2. Performance characteristics by patient age

Patient Age	Positive Samples			Negative Samples		
	ImmunoCard C. difficile GDH/Culture	Sensitivity %	95% CI	ImmunoCard C. difficile GDH/Culture	Specificity % <sup>1</sup>	95% CI
0-28 days	0/0	N/A	N/A	1/1	100.0%	20.7 – 100.0%
29 days to 2 years	22/22	100.0%	85.1 - 100.0%	72/98	73.5%	64.0 - 81.2%
> 2 years to < 12 years	30/30	100.0%	88.6 - 100.0%	130/160	81.3%	74.5 - 86.5%
12 years to < 18 years	14/15	93.3%	70.2 - 98.8%	80/88	90.9%	83.1 - 95.3%
18 years to 21 years	5/5	100.0%	56.6 - 100.0%	25/29	86.2%	69.4 - 94.5%
> 21 years	53/55	96.4%	87.7 - 99.0%	430/472	91.1%	88.2 - 93.3%

The ImmunoCard C. difficile GDH assay was designed to be a highly sensitive tier-one test with a Negative Predictive Value (NPV) close to 100% to increase the likelihood that the assay would identify most potentially positive samples. Because of the assay design, it is expected that the lower bounds of the confidence intervals could be below 90%. The ImmunoCard C. difficile GDH assay has a NPV of 99.6% (CI = 98.8 - 99.9%), which demonstrates that the assay is sufficiently sensitive to identify most potentially positive samples, therefore increasing the likelihood that all toxigenic samples would be identified by the subsequent testing.

### ANALYTICAL SENSITIVITY

The analytical sensitivity of the ImmunoCard C. difficile GDH test is 10 ng/mL in stool and 1.0 ng/mL when diluted in a non-matrix system.

### REPRODUCIBILITY

Assay precision, intra-assay variability and inter-assay variability were assessed with a reference panel prepared from pools of negative samples spiked with *C. difficile* GDH antigen. The reproducibility panel consisted moderately positive (n=3), low positive (n=3), high negative (n=3), and negative (n=1) specimens. The low positive and high negative specimens were prepared near the assay limit of detection. Each sample was evaluated twice per day for five days by three different laboratories. Reproducibility was 99.7% with no intra-assay and inter-assay variability for samples prepared above or below the limit of detection for the assay.

### CROSSREACTIVITY

Crossreactivity studies were performed with positive and negative stool specimens inoculated with bacterial or fungal organisms to a final concentration of 1.2 x 10<sup>8</sup> CFU/mL or virus concentration greater than 1 x 10<sup>5</sup> TCID<sub>50</sub>/mL. None of the following organisms in stool reacted with ImmunoCard C. difficile GDH:

*Aeromonas hydrophila*, *Bacillus cereus*, *Bacillus subtilis*, *Bacteroides fragilis*, *Campylobacter coli*, *Campylobacter fetus*, *Campylobacter jejuni*, *Candida albicans*, *Citrobacter freundii*, *Clostridium bifermentans*, *Clostridium butyricum*, *Clostridium haemolyticum*, *Clostridium histolyticum*, *Clostridium novyi*, *Clostridium perfringens*, *Clostridium septicum*, *Clostridium sordellii*, *Clostridium sporogenes*, *Clostridium tetani*, *Enterobacter aerogenes*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Escherichia coli*, *Escherichia coli* O157:H7, *Escherichia hermannii*, *Escherichia fergusonii*, *Helicobacter pylori*, *Klebsiella pneumoniae*, *Lactococcus lactis*, *Listeria monocytogenes*, *Peptostreptococcus anaerobius*, *Plesiomonas shigelloides*, *Porphyromonas asaccharolytica*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Pseudomonas fluorescens*, *Salmonella* Group B, *Salmonella* Group C, *Salmonella* Group D, *Salmonella* Group E, *Serratia liquefaciens*, *Serratia marcescens*, *Shigella boydii*, *Shigella flexneri*, *Shigella sonnei*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Vibrio parahaemolyticus*, *Yersinia enterocolitica*, Adenovirus Type 40, Adenovirus Type 41, Coxsackievirus Strain B4, Echovirus Strain 30, Rotavirus Strain WA.

Stools spiked with *Staphylococcus aureus* (Cowan strain I) were found to be crossreactive with ImmunoCard C. difficile GDH.

### ASSAY SPECIFICITY

The following *C. difficile* stock cultures from different sources were tested and produced positive reactions at a concentration of 1.2 x 10<sup>7</sup> CFU/mL with the ImmunoCard C. difficile GDH assay:

Toxigenic *C. difficile* strains: 8864, 10463, 43598, 2004052, 2004111, 2004118, 2004205, 2004206, 2005070, 2005257, 2005325, 2006240, 2007431, 2007435, 2007858, 2008016, 2008029, 2008162, 2008188, 2008341, 2008351, 2009018, 2009065, 2009066, 2009099, 2009132, 2009155, 2009277, B1, B17, B18, BK6, CF1, G1, J7, K12, Y1

Non-toxigenic *C. difficile* strains: 11186, 234, 586, 611, 620, 2C62, 2C165, 1C272, UNC19904, X15076

### TESTS FOR INTERFERING SUBSTANCES

The following substances, at the specified saturated solvent/diluent concentrations, do not interfere with test results in the final concentrations listed:

Barium sulfate (5 mg/mL), Fecal fat (2.65 mg stearic acid and 1.3 mg palmitic acid/mL), Hemoglobin (3.2 mg/mL), Iodium AD<sup>®</sup> (Loperamide HCl) (6.67 x 10<sup>-3</sup> mg/mL), Kaopectate<sup>®</sup> (Bismuth subsalicylate) (0.87 mg/mL), Metronidazole (12.5 mg/mL), Mucin (3.33 mg/mL), Mylanta<sup>®</sup> (Aluminum hydroxide w/ magnesium hydroxide) (4.2 mg/mL), Pepto-Bismol<sup>®</sup> (Bismuth subsalicylate) (0.87 mg/mL), Polyethylene glycol (79.05 mg/mL), Prilosec<sup>®</sup> (Omeprazole) (0.5 mg/mL), Simethicone (0.625 mg/mL), Tagamet<sup>®</sup> (Cimetidine) (0.5 mg/mL), Tums<sup>®</sup> (Calcium carbonate) (5.0 mg/mL), Vancomycin HCl (2.5 mg/mL), Whole blood (40%), White blood cells (5%).

Samples containing whole blood at concentrations greater than 40% have the potential to produce false negative results.

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**REFERENCES**

- Dawson LF, Valiente E and Wren B *Clostridium difficile* - A continually evolving and problematic pathogen. Infection, Genetics and Evolution. 2009;9:1410-1417.
- Lyerly DM, Barroso LA and Wilkins TD Identification of the latex-reactive protein of *Clostridium difficile* as glutamate dehydrogenase. J. Clin. Micro 1991;29:2639-2642.
- Willis DH and Kraft JA Confirmation that the latex-reactive protein of *Clostridium difficile* is a glutamate dehydrogenase. J Clin Micro 1992;30:1363-1364.
- Barlett JG, Gerding D. Clinical recognition and diagnosis of *Clostridium difficile* infection. Clin Inf Dis 2008;46:S12-18.
- Cohen SH et al. SHEA/IDSA Guideline. Clinical practice guidelines for *Clostridium difficile* infection in adults: 2010 update by the Society for Healthcare Epidemiology of America (SHEA) and the Infectious Diseases Society of America (IDSA) 2010;31:431-455.

SN11084

REV. 12/16

**Manufactured By**

**Meridian Bioscience, Inc.**  
 USA/Corporate Office  
 3471 River Hills Drive  
 Cincinnati, Ohio 45244  
 Telephone: 513.271.3700  
 Orders/Customer Service:  
 800.543.1980  
 Technical Support Center:  
 800.343.3858  
 Information Fax: 513.272.5432  
 Ordering Fax: 513.271.0124

**Meridian Bioscience Europe s.a./n.v.**  
 2 Avenue du Japon - 1420 Braine l'Alleud  
 BELGIUM  
 Tel: +32 (0) 67 89 59 58  
 Fax: +32 (0) 67 89 59 59  
 Email: info.bn@meridianbioscience.eu

---

**Authorized Representative**

**Meridian Bioscience Europe S. r. l.**  
 Via dell' Industria, 7  
 20020 Villa Cortese, Milano  
 ITALY  
 Tel: +39 0331 43 36 36  
 Fax: +39 0331 43 36 16  
 Email: info@meridianbioscience.eu  
 WEB: www.meridianbioscience.eu

**Meridian Bioscience Europe b.v.**  
 Postbus 301 - 5460 AH Veghel  
 NETHERLANDS  
 Tel: +31 (0) 411 62 11 66  
 Fax: +31 (0) 411 62 48 41  
 Email: info.bn@meridianbioscience.eu

**INTERNATIONAL SYMBOL USAGE**

You may see one or more of these symbols on the labeling/packaging of this product:  
**Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)**

	Use By / Utilizzare entro / Utiliser jusqu' / Fecha de caducidad / Verwendbar bis	<b>CONTROL +</b>	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
<b>LOT</b>	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	<b>CONTROL -</b>	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
<b>IVD</b>	In vitro diagnostic medical device / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	<b>EC REP</b>	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostika 98/79/EG.	<b>SMP   PREP   DIL   SPE</b>	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenzubereitung, in dem sich Probenverdünnungspuffer befindet
<b>REF</b>	Catalogue number / Numero di catalogo / Référence ou catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	<b>BUF   RXN</b>	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Reagents IVD reserves a Révaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contiene sufficiente per "n" tests / Contenu suffisant pour "n" tests / Inhalt ausreichend für <n> Prüfungen	<b>SOLN   STOP</b>	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	<b>CONJ   ENZ</b>	Enzyme Conjugate / Coniugato enzimatico / Conjugé enzymatique / Conjugado enzimático / Enzymkonjugat
<b>SN</b>	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	<b>CONTROL</b>	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
<b>TEST</b>	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	<b>REAG</b>	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	<b>BUF   WASH</b>	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Lavage / Tampón de lavado / Waschlösung
<b>BUF</b>	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenza / Mise En Garde / Advertencia / Warnhinweise
<b>CONJ</b>	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	<b>DIL   SPE</b>	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillon / Diluyente de muestra / Probenverdünnungspuffer
<b>SUBS</b>	Substrate / Substrato / Substrat / Substrato / Substrat	<b>BUF   WASH   20X</b>	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentrierter Waschlösung
<b>R. Only</b>	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Sólo Para Uso Por Receta / verschreibungspflichtig	<b>DET   REAG</b>	Detection Reagent / Reagente Diretto / Réactif de Détection / Reactivo de Detección / Nachweis Reagenz

For technical assistance, please call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, Eastern Standard Time. To place an order, please call Customer Service Department at 800-543-1980.