## A Rapid Immunoassay for the Detection of *Helicobacter pylori* Antigens in Stool Specimens

REF

750020

IVD

### INTENDED USE

The ImmunoCard STAT! HpSA HD is a rapid in vitro qualitative procedure for the detection of *Helicobacter pylori* antigens in human stool. The stool antigen detection is intended to aid in the diagnosis of *H. pylori* infection.

### SUMMARY AND EXPLANATION OF THE TEST

Since its discovery over 20 years ago by Marshall and Warren,<sup>1</sup> *Helicobacter pylori* is now recognized as one of the most common and medically important pathogens worldwide. *Helicobacter pylori* has been firmly established as an etiologic agent in chronic gastritis and peptic ulcer disease, and has been associated with mucosa-associated lymphoid tissue lymphoma and gastric adenocarcinoma.<sup>2-6</sup>

The ecological niche in humans appears to be restricted to the stomach and the duodenum. Patients who harbor the organism are divided into two basic groups. The first group shows no signs or symptoms of gastrointestinal disease and is considered as "colonized". The second group shows gastrointestinal signs and symptoms and is considered as "infected". The process by which an individual becomes colonized or infected is still under investigation.<sup>2, 3, 7-9</sup> Many possible routes of transmission of *Helicobacter pylori* to humans such as animals, contaminated water and oral reservoirs have been suggested.<sup>10</sup>

Diagnostic tests for *H. pylori* can be categorized as invasive (endoscopy, biopsy) or noninvasive (serology, urea breath test and stool antigen test). In invasive testing, a biopsy is taken from the upper gastrointestinal tract and examined microscopically. The tissue is also cultured for *H. pylori* or evaluated in the rapid urease test. This strategy offers the advantages of detecting an active infection and has high specificity and a high positive predictive value. The disadvantages of invasive testing include risk and discomfort to the patient and colonization in patches that might be missed by biopsy. Culture of biopsy material is time consuming and can yield false-negative results due to inherent technical difficulties.<sup>2, 11-14</sup>

The urea breath test (UBT) is a type of noninvasive test that detects the highly active urease of *H. pylori*. Although UBT is highly sensitive and specific, it has a number of significant drawbacks. UBT is time consuming, requires specialized detection equipment and involves the ingestion of isotopically labeled urea by the patient.<sup>2, 15, 16</sup> Serological tests, also noninvasive, based on the detection of IgG against *H. pylori* are useful for primary screening of patients that present with uncomplicated infections, yet they do not distinguish between past exposure and active infection.<sup>5, 10, 17</sup> The stool antigen test has been evaluated extensively and has been accepted as an accurate non-invasive test both before and after treatment.<sup>18-20</sup> The recent Maastricht 4 Consensus Report recommends the use of the stool antigen and UBT tests as an aid in the diagnosis of *H. pylori* disease in the primary care setting.<sup>21</sup>

ImmunoCard STAT! HpSA HD is designed to detect *H. pylori* antigen in human stool.

## BIOLOGICAL PRINCIPLES

ImmunoCard STAT! HpSA HD is a rapid lateral flow immunoassay that utilizes a monoclonal anti-*H. pylori* antibody as the capture and detector antibodies. A diluted patient stool sample is dispensed into the sample port of the test device and the appearance of a pink-red line in the reading window next to the letter T within 5-15 minutes of incubation at room temperature indicates a positive result.

## REAGENTS/MATERIALS PROVIDED

**The maximum number of tests obtained from this test kit is listed on the outer box.**

1. **ImmunoCard STAT! HpSA HD Test Device:** A chromatography strip housed in a plastic frame and enclosed in a foil pouch with a desiccant. The strip carries monoclonal anti-*H. pylori* capture antibody for the test and an animal protein for a control. The strips also contain red-latex conjugated anti-*H. pylori* and blue latex-conjugated anti-protein as the detector antibodies for tests and controls, respectively. Store the devices at 2-8 C when not in use. Do not freeze.
2. **Sample Diluent:** A buffered Tris solution containing 1% bovine albumin and detergents. Sodium azide (0.095%) is added as a preservative. The Diluent is supplied in a red-capped plastic dropper vial with an applicator tip. Use as supplied. Store at 2-8 C when not in use.
3. **Positive Control:** A suspension of inactivated *H. pylori* in a balanced salt solution containing sodium azide (0.095%) as a preservative. The reagent is supplied ready for use in a plastic dropper vial. Use as supplied. Store at 2-8 C when not in use.
4. 100 µL transfer pipettes

## MATERIALS NOT PROVIDED

1. Disposable latex gloves that should be used during the handling of the fecal samples as they are considered potentially hazardous material
2. Vortex for suspending the stool specimen in the Sample Diluent (Optional)
3. Interval timer

## PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Patient specimens may contain infectious agents and should be handled and disposed of as potentially biohazardous.
3. Do not interchange reagents from different kit lot numbers. Do not use kit components beyond the labeled expiration date of the kit.
4. Inspect Test Devices before removing the foil pouch. Do not use Test Devices that have holes in the foil pouch or where the pouch has not been completely sealed. False-negative reactions may result due to deterioration of the improperly stored Test Device.
5. Do not use the Sample Diluent or Positive Control if turbid. Turbidity may be a sign of microbial contamination.
6. Handle the Positive Control as if it is potentially infectious, even though it contains inactivated *H. pylori*.
7. Do not deviate from the method described here or falsely positive or falsely negative results may occur.
8. Test instructions should be thoroughly read before performing any testing.
9. Do not use a device if its pouch was punctured prior to use.

## HAZARDS AND PRECAUTION PHRASES

There are no known hazards associated with this product.

## SHELF LIFE AND STORAGE

Store the kit at 2-8 C when not in use. The shelf life (expiry) for this product is listed on the kit box label.

## PROCEDURAL NOTES

1. Allow kit components and specimens to reach the room temperature (20-26 C) before performing a test, as cold reagents and/or specimens may decrease assay sensitivity. Reagents may take 20-30 minutes to warm following refrigeration.
2. Stool samples must be mixed thoroughly (regardless of consistency) to ensure a representative sample prior to sampling.
3. Hold reagent vials vertically when dispensing drops to ensure consistent drop size and delivery.
4. On occasion, particulate matter may interfere with sample flow. In cases where the Test Device does not readily absorb the diluted specimen, gently touch the bottom of the sample port with an applicator stick, moving the stool solid particle that might prevent the absorption. Alternatively, a new aliquot of the sample can be withdrawn from the Diluent and retested. Diluted samples containing a heavy concentration of particulate matter may be centrifuged (1-5 minutes at 700 x G) or allowed to stand for 3-5 minutes before proceeding in order to allow particulate matter to settle to the bottom of the test tube.

## REAGENT PREPARATION

Reagents are supplied ready for use. No preparation is needed.

## SPECIMEN COLLECTION AND PREPARATION

**DO NOT USE stool in transport media, on swabs, or mixed with preservatives.** The specimen should be transported in an airtight container and stored at 2-8 C until tested. The specimen should be tested as soon as possible, but may be held up to 72 hours at 2-8 C prior to testing. If testing cannot be performed within this time frame, specimens should be frozen immediately on receipt and stored frozen ( $\leq -20$  C) until tested. Specimens may be frozen and thawed twice. Mix stool samples thoroughly (regardless of consistency) before testing.

1. **Liquid or Semi-solid stools** – Unscrew the red cap from the Sample Diluent vial (red capped vial). Use a clean calibrated transfer pipette (supplied with the kit) to draw the mixed sample to the second mark from the pipette tip (100  $\mu$ L). Dispense the sample into the Sample Diluent vial. Use the same transfer pipette to mix the diluted sample thoroughly, but gently, by squeezing the pipette bulb three times. Recap the vial tightly and mix thoroughly but gently by swirling the contents for 15 seconds. Alternatively, mix for 15 seconds using a vortex mixer. **NOTE: Care should be taken when pipetting semisolid stool. The addition of less than 100  $\mu$ L of stool may cause a false-negative test. The addition of more than 100  $\mu$ L of stool may cause invalid results due to restricted sample flow.**


Figure of 100  $\mu$ L pipette

2. **Formed/Solid stools** – Unscrew the red cap of the Sample Diluent vial (red capped vial). Use the white plastic applicator stick in the red cap to collect a small portion of stool (5-6 mm pellet). Transfer the pellet to the Sample Diluent vial. Recap the vial tightly and mix thoroughly but gently by swirling the contents of the vial for 15 seconds. Alternatively, mix for 15 seconds using a vortex mixer. Wooden applicator sticks may also be used to transfer solid stool to the Sample Diluent. **NOTE: The transfer of too little stool, or failure to mix and suspend the stool in Sample Diluent completely may result in a false-negative test results. Care should be taken to transfer no less and no more than the amount indicated. The addition of more than 100  $\mu$ L of stool may cause invalid results due to restricted sample flow.**

## TEST PROCEDURE

### A. Test

1. Bring all test devices, reagents and samples to room temperature (20-26 C) before testing.
2. Use 1 ImmunoCard STAT! HpSA HD Test Device for each patient sample.
3. Remove the ImmunoCard STAT! HpSA HD Test Device from its foil pouch. The Test Device is marked to indicate where test and control lines will appear. The round window marked with an arrow is the test window where sample is added.
4. Label the device with the patient's name. Prepare the specimen according to the instructions in the SPECIMEN COLLECTION AND PREPARATION section above.
5. Hold the diluted specimen vial upright and tap the bottom gently on the countertop before proceeding.
6. Cover the top of the diluted sample vial with absorbent paper to avoid splatter.
7. Break off the red tip on the outside of the red cap. (Do not break off the white applicator stick on the inside of the cap.)
8. Hold the vial upside down and dispense 3 drops of diluted sample into the round window (at arrow) of the Test Device. Do not touch the tip of the vial to the Test Device.
9. Set a timer and incubate the test at 20-26 C for 5 minutes.
10. At the end of 5 minutes, read the results. Negative tests may be incubated for an additional 10 minutes and then reread. See the INTERPRETATION OF RESULTS section below for a description of positive and negative test results.

### B. Controls

Positive and Negative Controls are designed to show all reagents are reactive, specific, and capable of producing the expected results.

1. Bring all control reagents to 20-26 C before testing.
2. Use 1 ImmunoCard STAT! HpSA HD Test Device each for a Positive and Negative Control. Label each device with the control to be tested.
3. Hold reagent vials upside down to dispense reagents.
4. Add 3 drops of the Positive Control to the test window (at arrow) of 1 device. **Do not allow the tip of the vial to touch the sample port.**
5. Break off the red tip on the outside of the red cap of an unused vial of Sample Diluent.
6. Dispense 3 drops of the Sample Diluent to the test window (at arrow) of another Test Device.
7. Set a timer and incubate the tests at 20-26 C for 5 minutes.
8. After 5 minutes, read the results within 1 minute of test completion.

## INTERPRETATION OF RESULTS

**Negative test result:** Only one BLUE colored band (Control Line) appears across the central window of the device close to the letter "C". (*H. pylori* antigens are absent or below the level of detection.) No other bands should be seen. The background should not interfere with reading the test.

**Positive test result:** In addition to the BLUE band (Control Line), a distinguishable PINK-RED band (Test Line) also appears across the central window of the device close to the letter "T". The intensity of the band will vary depending on the antigen concentration in the specimen. Any pink-red line, even very weak, must be considered as a positive result. (A positive test line indicates that *H. pylori* antigens are in the specimen.) The background should not interfere with reading the test.

### Invalid test results:

1. The BLUE band (Control Line) is absent, with or without a visually detectable PINK-RED band (Test Line).
2. A PINK-RED band appears at the letter "T" in the window after 15 minutes, or there is a line at this position of another color other than pink-red,

3. No Control Line band appears close to the letter "C". (The test is invalid since a shift in or absence of the control line indicates that the test procedure was performed improperly or that deterioration of the reagents has occurred.)

**If any test is difficult to interpret, the test should be repeated with the same sample to eliminate the potential for error. Obtain a new sample and retest when the original sample repeatedly produces unreadable results.**

#### **QUALITY CONTROL**

***This test should be performed per applicable local, state, or federal regulations or accrediting agencies.***

The reactivity of ImmunoCard STAT! HpSA HD Test Devices should be verified on receipt using the external Positive and Negative Control reagents provided in the kit. The number of additional tests performed with the external controls will be determined by the requirements of local, state or national regulations or accrediting agencies.

**Internal controls: Internal controls are contained within the Test Device and therefore are evaluated with each test.**

1. A colored band appearing at the control line serves as a positive control and indicates the test has been performed correctly, that sample was added, that it flowed properly, and that the test reagents were active at the time of use.
2. A clear background around the control or test lines serves as a negative control. A background that obscures the reading of results invalidates the test and is an indication of reagent deterioration, inappropriate sample or improper test performance.

The results expected with the Controls are described in the section on INTERPRETATION OF RESULTS. The Test Devices should not be used if control tests do not produce the correct results. Failure to achieve the expected results indicates either that the Test Devices are defective or that the test was not performed correctly. **If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.**

The Positive and Negative Control reagents are manufactured in an aqueous solution matrix. Although specimen matrix interference has not been observed with this assay, the aqueous matrix of the controls may not adequately control for specimen matrix effects.

#### **EXPECTED VALUES**

Studies on the epidemiology of *H. pylori* have shown that this organism is present worldwide.<sup>17, 22, 23</sup> Gastritis caused by *H. pylori* has been shown to correlate with age, ethnic background, family size and socioeconomic class.<sup>24, 25</sup> The prevalence of *H. pylori* infection in a given population can vary from 20% to 90%. In patients diagnosed with duodenal ulcers, however, it has been shown in every age group to be approximately 80%.<sup>17</sup> Currently recommended eradication treatments have an efficacy rate between 75% and 90%.

The ImmunoCard STAT! HpSA HD test detects the presence of *H. pylori* antigens in human stool. Expected values for a given population should be determined for each laboratory. The rate of positives may vary depending on geographic location, method of specimen collection, handling and transportation, test employed and general health environment of patient population under study.

#### **LIMITATIONS OF THE PROCEDURE**

1. The test is qualitative and no quantitative interpretation should be made with respect to the intensity of the positive line when reporting the result.
2. Test results are to be used in conjunction with information available from the patient clinical evaluation and other diagnostic procedures.

- Antimicrobials, proton pump inhibitors and bismuth preparations are known to suppress *H. pylori*, and ingestion of these prior to *H. pylori* testing (culture, histology, rapid urease, UBT, antigen) may cause false-negative results. If a patient has ingested these compounds within two weeks prior to performing the ImmunoCard STAT! HpSA HD test, a false-negative result may occur. In such cases, the test should be repeated on a new specimen obtained two weeks after discontinuing treatment. A positive result for a patient ingesting these compounds within two weeks prior to performing the ImmunoCard STAT! HpSA HD test, should be considered accurate.
- Failure to add sufficient stool to the Specimen Diluent may result in a falsely negative test result. Addition of too much stool may result in invalid test results due to the inhibition of proper sample flow.
- Overincubation of tests may lead to false-positive test results. Incubating tests at reduced temperatures or times may lead to falsely negative results.

### SPECIFIC PERFORMANCE CHARACTERISTICS

Two independent laboratories compared the performance of ImmunoCard STAT! HpSA HD in parallel with the reference EIA method, Premier Platinum HpSA PLUS (Meridian Bioscience, Inc, Cincinnati, OH). Sixty one prospective (fresh) and 114 retrospective (frozen) samples were included in the evaluation. Table 1 provides the performance data for ImmunoCard STAT! HpSA HD. Table 2 lists the OD readings for ImmunoCard STAT! HpSA HD false-negative samples.

Table 1. Performance data for ImmunoCard STAT! HpSA HD

Premier Platinum HpSA PLUS	ImmunoCard STAT! HpSA HD		Total
	Positive	Negative	
Positive	50	3	53
Negative	1	121	122
Total	51	124	175
Positive Agreement (95% CI)	94.3% (84.6 to 98.1%)		
Negative Agreement (95% CI)	99.2% (95.5 to 99.9%)		

Table 2. EIA OD readings for False-Negative ImmunoCard STAT! HpSA HD samples

Sample No.	OD reading (cut off 0.100)	Interpretation
64	0.117	EIA weak pos; just above negative cut off
96	0.924	EIA pos
118	0.193	EIA weak pos; just above negative cut off

### ANALYTICAL SENSITIVITY

The lower limit of detection for this assay is 16 ng/mL in tests with sonicated antigen prepared from *H. pylori* strain TV1970. The limit does not vary from formed (solid) to semisolid stool samples.

### REPRODUCIBILITY

The repeatability of this assay was determined using (1) a natural negative and (2) contrived positive samples prepared from varying concentrations of the *H. pylori* standard. Ten replicates of each sample were tested by one operator at one time. Concordance in results was 100%.

Inter-day precision: The inter-day precision was determined by preparing four serial dilution curves (sensitivity curves) of the *H. pylori* standard. The dilutions were tested by the same operator on four different days with a concordance in results of 100%.

Inter-operator precision: Inter-operator precision was determined using three serial dilutions (sensitivity curves) of the *H. pylori* standard. The dilutions were then analyzed by three operators on the same day with a concordance in results of 100%.

Inter-batch precision: Inter-batch precision was determined by analyzing the standard dilutions with three different lots of product. Tests were performed by one operator on one day. Concordance in results was 100%.

#### CROSSREACTIVITY/ASSAY SPECIFICITY

Crossreactivity was evaluated utilizing the following bacterial and yeast strains. Positive and negative stools were spiked with  $\geq 1 \times 10^6$  bacteria/ fungus /yeast. None of the microorganisms tested yielded a positive result in negative stool or interfered with detection of a positive stool. Both negative and positive stool were positive when spiked with *H. pylori* strain 43504.

*Campylobacter coli*, *Campylobacter jejuni*, *Clostridium difficile*, *Helicobacter felis*, *Klebsiella pneumoniae*, *Pseudomonas aeruginosa*, *Salmonella ssp.*, *Staphylococcus aureus* (Cowan I), *Shigella ssp*

#### TESTS FOR INTERFERING SUBSTANCES

The following substances were found to have no effect on results when present in stool at the concentrations indicated.

Tums® Antacid (5 mg/mL), Tagamet® (5 mg/mL), Prilosec® (5 mg/mL), Mylanta® Antacid (1:20), Pepto-Bismol® (1:20), Barium sulfate (5%), Whole Blood (50%), Mucin (3.4%)

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## REFERENCES

1. Marshall BJ, Warren JR. Unidentified curved bacilli in the stomach of patients with gastritis and peptic ulceration. *Lancet* 1984; i:1311-1314.
2. Dunn BE, Cohen H, Blaser MJ. *Helicobacter pylori*. *Clin Microbiol Rev* 1997;10:720-41.
3. Gregson DB, Simor AE. *Helicobacter pylori* and inflammatory disease of the stomach and duodenum. *Stomach and Bowel*, 1991; Sep:29-30, 35-7.
4. Parsonnet J, Friedman GD, Vandersteen DP, et al. *Helicobacter pylori* infection and the risk of gastric carcinoma. *N Engl J Med* 1991; 325:1127-31.
5. Sanders MK, Peura DA. *Helicobacter pylori*-associated diseases. *Curr Gastro Reports* 2002; 4:448-54.
6. Uemura N, Okamoto S, Yamamoto S, et al. *Helicobacter pylori* infection and the development of gastric cancer. *N Engl J Med* 2001; 345:784-9.
7. Graham DYI. *Helicobacter pylori*: Its epidemiology and its role in duodenal ulcer disease. *J Gastroenterol Hepatol* 1991; 6:105-13.
8. Malaty HM, Graham DY, Klein PD, et al. Transmission of *Helicobacter pylori* infection. Studies in families of healthy individuals. *Scand J Gastroenterol* 1991; 26:927-32.

9. Marshall BJ, Armstrong JA, McGeachie DB, Glancy RJ. Attempt to fulfill Koch's postulates for pyloric campylobacter. *Med J Austral* 1985; 142:436-439.
10. Veralovic J, Fox JG. *Helicobacter*. In: Murray PR et al eds. *Manual of clinical microbiology*, 8<sup>th</sup> ed. Washington DC:ASM Press, 2003; 915-28.
11. Alpert LC, Graham DY, Evans Jr DJ, et al. Diagnostic possibilities for *Campylobacter pylori* infection. *Eur J Gastroenterol Hepatol* 1989; 1:17-26.
12. Barthel JS, Everett ED. Diagnosis of *Campylobacter pylori* infections: The "gold standard" and the alternatives. *Rev Infect Dis* 1990; 12:S107-14.
13. Nichols L, Sughayer M, DeGirolami PC, et al. Evaluation of diagnostic methods for *Helicobacter pylori* gastritis. *Am J Clin Pathol* 1991; 95:769-73.
14. Vaira D, Vakil N. Blood, urine, stool, breath, money and *Helicobacter pylori*. *Gut* 2001; 48:287-9.
15. Graham DY, Klein PD, Evans Jr DJ, et al. *Campylobacter pylori* detected noninvasively by the 13C-urea breath test. *Lancet* 1987; i:1174-7.
16. Graham DY, Malaty HM, Evans DG, et al. Epidemiology of *Helicobacter pylori* in an asymptomatic population in the United States. Effect of age, race, and socioeconomic status. *Gastroenterology* 1991; 100:1495-1501.
17. Talley NJ, Newell DG, Ormand JE, et al. Serodiagnosis of *Helicobacter pylori* Comparison of enzyme-linked immunosorbent assays. *J Clin Microbiol* 1991; 29:1635-9.
18. Vaira D, Malfertheiner P, Mégraud F, et al. Diagnosis of *Helicobacter pylori* infection with a new non-invasive antigen-based assay. *Lancet* 1999; 354:30-3.
19. Vaira D, Malfertheiner P, Mégraud F, et al. Non invasive antigen-based assay for assessing *Helicobacter pylori* eradication: a European Multicentre Study. *Am J Gastroenterol* 2002; 95:925-9.
20. Montiero L, de Mascarel A, Sarrasqueta A, et al. Diagnosis of *Helicobacter pylori* Infection: Noninvasive Methods Compared to Invasive Methods and Evaluation of Two New Tests. *Am J Gastroenterol* 2001; 96:353-22.
21. Malfertheiner P, Megraud F, Morain O, et al. Current concepts in the management of *Helicobacter pylori* infection – The Maastricht 2-2000 Consensus Report. *Aliment Pharmacol Ther* 2002; 16:167-80.
22. Murray DM. Clinical relevance of infection by *Helicobacter pylori*. *Clin Microbiol Newlett* 1993; 15:33-7.
23. Loffeld RJLF, Stobberingh E, Van Spreeuwel JP, et al. The prevalence of anti-*Helicobacter (Campylobacter) pylori* antibodies in patients and healthy blood donors. *J Med Microbiol* 1991; 32:105-9.
24. Morris A, Nicholson G, Lloyd G, et al. Seroepidemiology of *Campylobacter pyloridis*. *N Z Med J* 1986; 99:657-9.
25. Fiedorek SC, Malaty HM, Evans DL, et al. Factors influencing the epidemiology of *Helicobacter pylori* infection in children. *Pediatrics* 1991; 88:578-582.

REV. 04/18

Manufactured By

**Meridian Bioscience, Inc.**  
Corporate Office  
3471 River Hills Drive  
Cincinnati, Ohio 45244 USA  
Telephone: 513.271.3700  
Orders/Customer Service:  
800.543.1980  
Technical Support Center:  
800.343.3858  
Information Fax: 513.272.5432  
Ordering Fax: 513.271.0124

Authorized Representative

**Meridian Bioscience Europe S. r. L**  
Via dell'Industria, 7  
20020 Villa Cortese, Milano  
ITALY  
Tel: +39 0331 43 36 36  
Fax: +39 0331 43 36 16  
Email: [info@meridianbioscience.eu](mailto:info@meridianbioscience.eu)  
WEB: [www.meridianbioscience.eu](http://www.meridianbioscience.eu)

**Meridian Bioscience Europe s.a./n.v.**  
2 Avenue du Japon - 1420 Braine l'Alleud  
BELGIUM  
Tel: +32 (0) 67 89 59 59  
Fax: +32 (0) 67 89 59 58  
Email: [info.bnl@meridianbioscience.eu](mailto:info.bnl@meridianbioscience.eu)

**Meridian Bioscience Europe France**  
34 rue de Ponthieu - 75008 Paris  
FRANCE  
Tel: +33 (0) 1 42 56 04 40  
Fax: +33 (0) 9 70 06 62 10  
Email: [info.fr@meridianbioscience.eu](mailto:info.fr@meridianbioscience.eu)

**Meridian Bioscience Europe b.v.**  
Postbus 301 - 5460 AH Veghel  
NETHERLANDS  
Tel: +31 (0) 411 62 11 66  
Fax: +31 (0) 411 62 48 41  
Email: [info.bnl@meridianbioscience.eu](mailto:info.bnl@meridianbioscience.eu)

## INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	<b>CONTROL +</b>	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
<b>LOT</b>	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	<b>CONTROL -</b>	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
<b>IVD</b>	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	<b>EC REP</b>	Authorized representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	<b>SMP PREP DIL SPE</b>	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenzubereitung, in dem sich Probenverdünnungspuffer befindet
<b>REF</b>	Catalogue number / Numero di catalogo / Référence du catalogue / Numero de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einlagieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	<b>BUF RXN</b>	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Solo per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contiene sufficiente per "n" saggi / Contenu suffisant pour "n" test / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	<b>SOLN STOP</b>	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	<b>CONJ ENZ</b>	Enzyme Conjugate / Coniugato enzimatico / Conjugado enzimático / Conjugado enzimático / enzymkonjugat
<b>SN</b>	Serial number / Numero di serie / Numéro de série / Numero de serie / Seriennummer	<b>CONTROL</b>	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
<b>TEST</b>	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / testgerät	<b>REAG</b>	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	<b>BUF WASH</b>	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
<b>BUF</b>	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
<b>CONJ</b>	Conjugate / Coniugato / Conjugado / Conjugado / Konjugat	<b>DIL SPE</b>	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
<b>SUBS</b>	Substrate / Substrato / Substrat / Substrato / Substrat	<b>BUF WASH 20X</b>	Wash Buffer Concentration 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
<b>R. Only</b>	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig	<b>DET REAG</b>	Detection Reagent / Reagente Diretto / Réactif de Détection / Reactivo de Detección / Nachweis Reagens
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / ne pas utiliser si le paquet est endommagé / No use si el paquete esta dañado / Nicht verwenden, wenn die Verpackung beschädigt ist	<b>TUBE</b>	Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service department at 800-543-1980.