# ImmunoCard STAT! ROTAVIRUS

A RAPID IMMUNOASSAY FOR THE DETECTION OF ROTAVIRUS ANTIGEN IN HUMAN STOOL

REF 750030

IVD

R<sub>x</sub> Only

**INTENDED USE**

The ImmunoCard STAT! Rotavirus Immunoassay is a rapid in vitro qualitative procedure for the detection of rotavirus antigen in human stool. The test can be used to aid in the diagnosis of rotavirus associated gastroenteritis.

**SUMMARY AND EXPLANATION OF THE TEST**

Rotavirus is a major cause of acute gastroenteritis, especially in children 6 to 24 months in age.<sup>1-3</sup> In addition, rotavirus infections can produce severe illness as well as asymptomatic infection in adults.<sup>1,2,4,5</sup> The incubation period of rotavirus infection is usually one to three days, followed by gastroenteritis with an average duration of five to eight days. Virus titers in stool reach a maximum shortly after the onset of illness, then decline.<sup>2,6,7</sup> Due to inadequacies in existing culture methods, human rotavirus is not routinely isolated from rotavirus-containing specimens.<sup>2,8</sup> For many years, electron microscopy has been the standard method for rotavirus detection. However, newly introduced enzyme immunoassays and latex agglutination assays with increased sensitivities and specificities are now the methods of choice.<sup>9-15</sup> The ImmunoCard STAT! Rotavirus test offers a simple, rapid method for detecting rotavirus antigen in patient stool.

**BIOLOGICAL PRINCIPLES**

The ImmunoCard STAT! Rotavirus assay detects the presence of rotavirus antigen in stool. Patient specimen is diluted 1:15 in Sample Diluent. The suspension is mixed and 150 µL is added to the sample port of the device. The sample mobilizes gold particles coated with monoclonal antibody to rotavirus and migrates along the membrane through the **Test** (polyclonal anti-rotavirus antibody) and **Control** zones. After ten minutes, the **Test** and **Control** zones are observed for the presence of red/purple lines across the membrane surface. If rotavirus is present in the sample, a complex is formed between the capture antibody and the monoclonal antibody-gold conjugate which can be seen visually as a red/purple line in the **Test** zone. No red/purple line in the **Test** zone indicates a negative result. The **Control** line serves as a procedural control to assure that the sample has migrated the appropriate distance along the membrane.

**REAGENTS/MATERIALS PROVIDED**

The maximum number of tests obtained from this kit is listed on the outer box.

1. **ImmunoCard STAT! Rotavirus Devices** - Individually foil pouched devices containing immobilized rabbit anti-rotavirus antibody (**Test** zone), goat anti-mouse IgG antibody (**Control** zone) and embedded monoclonal antibody to rotavirus VP-6 conjugated to gold particles.
2. **Positive Control** - Inactivated rotavirus (SA-11) in a buffer containing 0.02% thimerosal as a preservative.
3. **Sample Diluent** - Buffer containing 0.1% sodium azide as a preservative.
4. **Transfer Pipettes**

**MATERIALS NOT PROVIDED**

1. 12 x75 mm test tubes
2. Applicator sticks (recommended for solid stools)
3. Timer

**PRECAUTIONS**

1. All reagents are for in vitro diagnostic use only.
2. Reagent concentration, incubation times and temperatures (21-25 C) have been optimized for sensitivity and specificity. Best results are obtained by adhering to these specifications. Once the assay has been started, complete all subsequent steps without interruption.
3. Patient specimens and used ImmunoCard STAT! Rotavirus devices may contain infectious agents and should be handled at Biosafety Level 2 as recommended in the CDC/NIH manual "Biosafety in Microbiology and Biomedical Laboratories".
4. Positive Control reagent contains inactivated rotavirus. However, it should be handled as a potential biohazard.
5. All reagents should be gently mixed and at 21-25 C before use.
6. Do not interchange reagents from different kit lot numbers, or use expired reagents.
7. Hold Positive Control vial vertically to insure proper drop size and delivery. Do not allow the tips of the vial or pipette to touch the sample port.
8. Use only one transfer pipette per control or specimen. Discard after use. Do not attempt to reuse.
9. Stool must be mixed thoroughly (regardless of consistency) to insure a representative sample prior to pipetting. Do not use stools that have dried out.
10. Dilution of stool as described under **Sample Dilution** is important. Over-inoculation of stool into the Sample Diluent may restrict movement within the ImmunoCard STAT! Rotavirus device so as to produce an invalid result.

**WARNING:** Some reagents in this kit contain sodium azide, which is a skin irritant. Avoid skin contact with reagents. Disposal of reagents containing sodium azide into lead or copper plumbing can result in the formation of explosive metal azides. This can be avoided by flushing with a large volume of water during such disposal.

**HAZARD and PRECAUTIONARY STATEMENTS**


There are no known hazards associated with this product.

**SHELF LIFE AND STORAGE**

The expiration date is indicated on the kit label. Store the kit at 2-8 C.

**PROCEDURAL NOTES**

1. The ImmunoCard STAT! Rotavirus format is diagrammed below.

2. Batch processing of samples or controls is possible provided correct incubation time is maintained for each device.
3. The **Control** zone of each device is a procedural control to assure that the sample has migrated sufficiently in the device to permit a valid test result to be read.

**REAGENT PREPARATION**

1. All reagents come ready to use (no further dilution is required).
2. Allow kit components to reach 21-25 C prior to use.
3. Gently mix liquid reagents prior to use.

**SPECIMEN COLLECTION AND PREPARATION**

For the best results, specimens should be collected after onset of symptoms. Several authors have reported declining numbers of rotavirus particles after day eight or nine, with peak counts occurring on days three through five. Samples collected after day eight or nine may be less reactive than those collected earlier in the course of the disease.<sup>2,6,7</sup> Specimens containing high levels of blood may fail to flow in the ImmunoCard STAT! Rotavirus device, resulting in an invalid test result. Testing of an additional specimen is recommended under such circumstances.

Stool specimens may be collected by the method routinely utilized by the laboratory, providing no dilution of the specimen occurs. A sample may be collected into a clean, dry container free of detergent residue, or obtained from a diaper.

The specimen should be tested as soon as possible, but may be stored up to 72 hours at 2-8 C prior to testing. If testing cannot be performed within this time frame, specimen should be frozen in a non-defrosting freezer (-20 C or lower).

**TEST PROCEDURE**

**Note: All components must be 21-25 C prior to use.**

1. **Sample Dilution**
  - a. Add 350 µL Sample Diluent to one 12 x 75 mm test tube for each specimen to be tested.
  - b. Mix stool thoroughly, regardless of consistency.
  - c. Liquid or semi-solid stool: Using a transfer pipette, draw stool to the 25 µL calibration point (first mark from tip of the pipette). Dispense the stool into the Sample Diluent in appropriate 12 x 75 mm tube. Using the same pipette, gently withdraw and expel the stool suspension several times, then vortex ten seconds. Leave transfer pipette in tube for further use. **Note:** Do not pipette more than 25 µL of stool. Over-inoculation with stool may produce invalid results. Proceed to Step 2 within 30 minutes.
  - d. Solid stool (nonpipetable stool from a diaper): Using a wooden applicator stick, transfer a 2 mm diameter portion of stool into the Sample Diluent in the appropriate 12 x 75 mm tube. Emulsify the stool thoroughly using the applicator stick, then vortex ten seconds. Place transfer pipette in the tube. Proceed to Step 2 within 30 minutes.
2. Remove appropriate number of ImmunoCard STAT! Rotavirus devices from their pouches. Label appropriately. Use one device per control or sample.
3. Vortex each diluted specimen for ten seconds. Using the original specimen transfer pipette, draw diluted sample to the 150 µL calibration point (second mark from tip of the pipette) and add to **Sample** port.
4. Incubate ten minutes at 21-25 C. **NOTE:** During the ten minute incubation, diluted specimen must move past the **Control** zone.
5. Visually read **Control** and **Test** zones for the presence or absence of a red/purple line at the end of the incubation period.

**INTERPRETATION OF RESULTS**

**Positive Test Result:** Visually detectable red/purple **Test** and **Control** lines. A positive result indicates the presence of rotavirus antigen.<sup>1</sup>

**Negative Test Results:** Visually detectable red/purple **Control** line. No red/purple **Test** line present. A negative result indicates that rotavirus antigen is absent or below the level of detection.

**Invalid Test Result:** No visually detectable red/purple **Control** line, with or without a visually detectable red/purple **Test** line.<sup>2</sup>

<sup>1</sup>Levels of the detection and capture antibodies were adjusted to achieve performance equivalence with the Premier Rotacloone EIA.

<sup>2</sup>Invalid test results may be due to a Reagent/Device problem, a procedural error, or over-inoculation of stool into Sample Diluent during specimen dilution. Re-dilute stool and repeat test. Stools containing high levels of blood may fail to flow properly, resulting in an invalid result. Testing with an additional specimen is recommended. On occasion, a stool may have very high levels of rotavirus antigen and will yield a visible **Test** line and no visible **Control** line. In such cases, the specimen may be diluted twofold or greater, beyond original 1:15 dilution (ex. 25 µL stool + 750 µL Sample Diluent) and retested.

**QUALITY CONTROL**

**This test should be performed per applicable local, state or federal regulations or accrediting agencies.**

The positive and negative controls should be assayed with each new kit lot or new shipment. Thereafter, Meridian recommends that the procedural control, performed with each test, is sufficient to verify reagent integrity and assay performance. The Negative and Positive Controls are intended to monitor for reagent failure, but will not ensure precision at the analytical assay cut-off. Additional controls may be tested in accordance with guidelines or requirements of local, state, and/or federal regulations and accrediting organizations.

Add three drops of Positive Control or, using a transfer pipette, add 150 µL Sample Diluent (Negative Control) directly to **Sample** port of appropriate device (do not dilute Positive Control).

1. The Positive Control should yield visually detectable red/purple **Test** and **Control** lines.
2. The Negative Control should yield a visually detectable red/purple **Control** line. No **Test** line should be present.


**If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.**

**EXPECTED VALUES**

The ImmunoCard STAT! Rotavirus test detects the presence of rotavirus antigen in stool. Expected values for a given population should be determined for each laboratory. The rate of positivity may vary depending on patient age, geographic location, season, method of specimen collection, handling and transportation, and a general health environment of the patient population under study.<sup>2,16</sup>

It has been reported that in neonates, when rotavirus was present, the disease was mild or totally asymptomatic.<sup>2,6,17</sup> However, during cooler months, rotavirus may account for approximately 50% or more of the gastroenteritis found in hospitalized children.<sup>3,16</sup> In adults, the incidence of serious gastroenteritis caused by the virus is relatively low and when infected, adults tend to be asymptomatic.<sup>2,4</sup> Studies from nursing homes and hospital geriatric wards show that this population is at an increased risk and susceptible to rotavirus associated disease.<sup>2,17</sup>

Clinical studies with ImmunoCard STAT! Rotavirus conducted at three different hospitals (two Children's Hospitals, Midwest and West Coast, and one General Hospital, East Coast) examined the incidence of rotavirus in different age groups (Figure 1). Values over the frequency bars indicate the number of ImmunoCard STAT! Rotavirus positive specimens.

## LIMITATIONS OF THE PROCEDURE

- The ImmunoCard STAT! Rotavirus test does not define the presence of rotavirus associated gastroenteritis, but only demonstrates the presence of the antigen in stool. As with all *in vitro* diagnostic procedures, test results should be interpreted by a physician in conjunction with other clinical information.
- Limit of detection in stool specimens is 1.8-3.7x10<sup>9</sup> rotavirus particles (SA-11) per test volume.
- The use of meconium stools in this assay is not recommended as their performance characteristics have not been evaluated.
- A positive result does not preclude the presence of other infective organisms.

## SPECIFIC PERFORMANCE CHARACTERISTICS

The ImmunoCard STAT! Rotavirus test was evaluated on stools from pediatric patients submitted for rotavirus testing at three sites. All specimens were tested by ImmunoCard STAT! Rotavirus, a competitor's flow-through membrane EIA, Premier Rotaclo<sup>®</sup> microwell EIA and electron microscopy (EM). In the table below, test results are compared to electron microscopy.<sup>18</sup>

Comparison of ImmunoCard STAT! Rotavirus, Competitor's Membrane EIA and Premier Rotaclo<sup>®</sup> EIA to Electron Microscopy Results for Rotavirus (n=249)

Electron Microscopy Results	ImmunoCard STAT! Results		Competitor's Membrane EIA Results		Premier Rotaclo <sup>®</sup> Results	
	Pos	Neg	Pos	Neg	Pos	Neg
Pos	121	9	122	8	122	8
Neg	5	115	11	109	7	113
Performance Statistic	Value	± 95% Conf.	Value	± 95% Conf.	Value	± 95% Conf.
Sensitivity	93.1%	4.4%	93.8%	4.1%	93.8%	4.1%
Specificity	95.8%	3.6%	90.8%	5.2%	94.2%	4.2%
Predictive Value Pos.	96.0%	3.4%	91.7%	4.7%	94.6%	3.9%
Predictive Value Neg.	92.7%	4.6%	93.2%	4.6%	93.4%	4.4%
Correlation	94.4%	2.9%	92.4%	3.3%	94.0%	2.9%

All five of the ImmunoCard STAT! Rotavirus positive, electron microscopy negative specimens were positive by Premier Rotaclo<sup>®</sup> and the competitor's membrane EIA. Six of nine ImmunoCard STAT! Rotavirus negative, electron microscopy positive specimens were negative by all three immunoassays. Furthermore, quantitative estimates by the electron microscopist indicated that eight of these specimens had few rotavirus particles.

ImmunoCard STAT! Rotavirus Compared to Premier Rotaclo<sup>®</sup> Results (n=250)

ImmunoCard STAT! Results	Premier Rotaclo <sup>®</sup> Results	
	Positive	Negative
Positive	126	0
Negative	3	122
Statistic	Value	± 95% Conf.
Relative Sensitivity	97.7%	2.6%
Relative Specificity	100.0%	0.0%
Relative Correlation	98.8%	1.3%

Two of three ImmunoCard STAT! Rotavirus negative, Premier Rotaclo<sup>®</sup> positive specimens were negative by electron microscopy. One of the three was electron microscopy positive.

## ASSAY SPECIFICITY

The ImmunoCard STAT! Rotavirus test was evaluated for specificity using the reference bacterial and viral strains listed below. A stool positive for rotavirus and a stool negative for rotavirus were spiked with bacteria (≥ 6 X 10<sup>8</sup> cfu/mL stool) or viruses (≥ 1 X 10<sup>2.7</sup> TCID<sub>50</sub>/mL).

## Microorganism or Virus (# strains tested)

<i>Aeromonas hydrophila</i> (1)	<i>Escherichia coli</i> (O157:H7) (1)	<i>Staphylococcus aureus</i> (1)
<i>Bacillus cereus</i> (1)	<i>Escherichia coli</i> (non-O157:H7, SLT+) (1)	<i>Staphylococcus aureus</i> (Cowan 1) (1)
<i>Bacillus subtilis</i> (1)	<i>Helicobacter pylori</i> (1)	<i>Streptococcus faecalis</i> (1)
<i>Bacteroides fragilis</i> (1)	<i>Klebsiella pneumoniae</i> (1)	<i>Vibrio cholerae</i> (1)
<i>Campylobacter coli</i> (1)	<i>Proteus mirabilis</i> (1)	<i>Vibrio parahaemolyticus</i> (1)
<i>Campylobacter fetus</i> (1)	<i>Pseudomonas aeruginosa</i> (1)	Adenovirus, types 2, 40, 41
<i>Campylobacter jejuni</i> (1)	<i>Salmonella</i> (1)	Coronavirus
<i>Citrobacter freundii</i> (1)	<i>Serratia liquefaciens</i> (1)	Coxsackie types A9, B1, B6
<i>Clostridium difficile</i> (1)	<i>Shigella dysenteriae</i> (1)	Echovirus, types 22, 32
<i>Clostridium sordellii</i> (1)	<i>Shigella flexneri</i> (1)	Enterovirus, type 69
<i>Enterobacter aerogenes</i> (1)	<i>Shigella sonnei</i> (1)	Poliiovirus, type 1
<i>Escherichia coli</i> (1)		

All organisms were found to be negative when spiked into the negative stool. In addition, they did not interfere with either the positive specimen or the procedural control.

Studies with different serotypes of group A rotavirus (the most common form) showed that the ImmunoCard STAT! Rotavirus test detects serotypes 1, 2, 3, and 4. Reactivity with serotypes 8 and 9 has not been evaluated.

## ASSAY PRECISION

Three physician office laboratories tested six specimens and the two controls, in triplicate, on each of three different days. The Test line of the ImmunoCard STAT! Rotavirus device showed 100% reproducibility with negative and medium positive stools, and 96% reproducibility on low positive stools. Control line reproducibility was 100% with all specimens.

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#### REFERENCES

- Burke B, Desselberger U. Minireview: Rotavirus Pathogenicity. *Virology*. 1996; 218: 299-305.
- Christensen ML. Human viral gastroenteritis. *Clin Microbiol*. 1989; Rev. 2: 51-89.
- Truant AL, Chonmaitree T. Incidence of rotavirus infection in different age groups of pediatric patients with gastroenteritis. *J Clin Microbiol*. 1982; 16: 568-569.
- Echeverria P, Blacklow NR, Cukor G, Vibulbandhikit S, Changchawalit S, Boonthai P. Rotavirus as a cause of severe gastroenteritis in adults. *J Clin Microbiol*. 1983; 18: 663-667.
- Marrie TJ, Lee SHS, Faulkner RS, Ethier J, Young CJ. Rotavirus infection in a geriatric population. *Arch Intern Med*. 1982; 142: 313-316.
- Coulson BS, Holmes IH. An improved enzyme linked immunosorbant assay for the detection of rotavirus in faeces of neonates. *J Virol Meth*. 1984; 8: 165-179.
- Nagayoshi S, Yamaguchi Y, Ichikawa T, Miyazu M, Morshima T, Ozaki T, et al. Changes of the rotavirus concentration in faeces during the course of acute gastroenteritis as determined by the immune adherence hemagglutination test. *Eur J Pediatr*. 1980; 134: 99-102.
- Sato K, Inaba I, Shinozaki T, Fujii R, Matsumoto M. Isolation of human rotaviruses in cell cultures. *Arch Virol*. 1981; 69: 155-160.
- Denney PH, Gauntlett DR, Spangenberg SE. Choice of reference assay for the detection of rotavirus in fecal specimens: Electron Microscopy versus Enzyme Immunoassay. *J Clin Microbiol*. 1990; 28: 1280-1283.
- Denney PH, Gauntlett DR, Tente WE. Comparison of nine commercial immunoassays for the detection of rotavirus in fecal specimens. *J Clin Microbiol*. 1988; 26: 1630-1634.
- Gilchrist MJR, Bretl TS, Moulthey K, Knowlton DR, Ward RL. Comparison of seven kits for detection of rotavirus in fecal specimens with a sensitive, specific enzyme immunoassay. *Diagn Microbiol Infect Dis*. 1987; 8: 221-228.
- Lipson SM, Zelinsky-Papez KA. Comparison of four latex agglutination (LA) and three enzyme-linked immunosorbant assays (ELISA) for the detection of rotavirus in fecal specimens. *Am J Clin Path*. 1989; 92: 637-643.
- Mathewson JJ, Jr Winsor DK, DuPont RL, Secor SL. Evaluation of assay systems for the detection of rotavirus in stool specimens. *Diagn Microbiol Infect Dis*. 1989; 12: 139-141.
- Rubenstein AS, Miller MF. Comparison of an enzyme immunoassay with electron microscopic procedures for detecting rotavirus. *J Clin Microbiol*. 1982; 15: 938-944.
- Thomas EE, Puterman ML, Kawano E, Curran M. Evaluation of seven immunoassays for detection of rotavirus in pediatric stool samples. *J Clin Microbiol*. 1988; 26: 1189-1193.
- Brandt CD, Kim HW, Rodriguez WJ, Arrobio JO, Jeffries BC, Stallings EP, et al. Pediatric gastroenteritis during eight years of study. *J Clin Microbiol*. 1983; 18: 71-78.
- Champsaur J, Questiaux E, Prevot J, Henry-Amar M, Goldszmidt D, Bourjouane M, et al. Rotavirus carriage, asymptomatic infection, and disease in the first two years of life. I. Virus shedding. *J Infect Dis*. 1984; 149: 667-674.
- Data on file, Meridian Bioscience, Inc.

SN11240

REV. 12/16

**CE**

Manufactured By

**EC REP**

Authorized Representative

**Meridian Bioscience Europe S. r. l.**  
Via dell'Industria, 7  
20020 Villa Cortese, Milano  
ITALY  
Tel: +39 0331 43 36 36  
Fax: +39 0331 43 36 16  
Email: info@meridianbioscience.eu  
WEB: www.meridianbioscience.eu

**Meridian Bioscience Europe s.a./n.v.**  
2 Avenue du Japon - 1420 Braine l'Alleux  
BELGIUM  
Tel: +32 (0) 67 89 59 59  
Fax: +32 (0) 67 89 59 58  
Email: info.bn@meridianbioscience.eu

**Meridian Bioscience Europe France**  
34 rue de Ponthieu - 75008 Paris  
FRANCE  
Tel: +33 (0) 1 42 56 04 40  
Fax: +33 (0) 9 70 06 62 10  
Email: info.fr@meridianbioscience.eu

**Meridian Bioscience Europe b.v.**  
Postbus 301 - 5460 AH Veghel  
NETHERLANDS  
Tel: +31 (0) 411 62 11 66  
Fax: +31 (0) 411 62 48 41  
Email: info.bn@meridianbioscience.eu

## INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	<b>CONTROL +</b>	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
<b>LOT</b>	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	<b>CONTROL -</b>	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
<b>IVD</b>	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	<b>EC REP</b>	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostico in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	<b>SMP PREP DIL SPE</b>	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet
<b>REF</b>	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le Istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	<b>BUF RXN</b>	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contiene sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contiene suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	<b>SOLN STOP</b>	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stoplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	<b>CONJ ENZ</b>	Enzyme Conjugate / Coniugato enzimatico / Conjugué enzymatique / Conjugado enzimático / Enzymkonjugat
<b>SN</b>	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	<b>CONTROL</b>	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
<b>TEST</b>	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	<b>REAG</b>	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	<b>BUF WASH</b>	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
<b>BUF</b>	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
<b>CONJ</b>	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	<b>DIL SPE</b>	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
<b>SUBS</b>	Substrate / Substrato / Substrat / Substrato / Substrat	<b>BUF WASH 20X</b>	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
<b>R. Only</b>	Prescription Use Only / Per uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig	<b>DET REAG</b>	Detection Reagent / Reagente Diretto / Réactif de Détection / Reactivo de Detección / Nachweis Reagenz

For technical assistance, call Technical Support Services at (800) 343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at (800) 543-1980.