immunocard STAT![®]

C. difficile GDH-AB

Rapid one-step Immunoassay for the Simultaneous Detection of *Clostridium difficile* Glutamate Dehydrogenase Antigen (GDH) and Toxins A and B in human Stool

REF 750520

IVD

INTENDED USE

Immunocard STAT! C. difficile GDH-AB is a rapid, qualitative, immunochromatographic assay for the simultaneous detection, in human stool, of *Clostridium difficile* Glutamate Dehydrogenase antigen (GDH, also called "common antigen") and *Clostridium difficile* Toxins A and B. This assay is used as an aid in the diagnosis of C. difficile-associated disease and results should be used by the clinician in conjunction with clinical picture, other laboratory findings and epidemiological risk factors.

SUMMARY AND EXPLANATION OF THE TEST

Clostridium difficile is a spore-forming, gram-positive, anaerobic bacterium that can be present asymptotically in up to 5% of the healthy population¹. Toxicogenic strains of *Clostridium difficile* represent the leading cause of nosocomial infectious diarrhea in developed countries, also representing the etiologic agent in approximately 25% of all cases of antibiotic-associated diarrhea². An estimated 300,000 cases of C. difficile infections (CDAD) are seen per year in U.S. hospitals alone.^{1, 2} Virtually any antibiotic can predispose a patient to CDAD. The clinical picture for CDAD ranges from asymptomatic colonization to life-threatening pseudomembranous colitis and toxic megacolon.² The pathogenic strains of C. difficile produce at least one of two biologically and immunologically distinct toxins³: toxin A (enterotoxin) and toxin B (cytotoxin), which are the main virulence factors responsible for the clinical signs of the disease.

Glutamate dehydrogenase (GDH), is an enzyme produced in large quantities by toxigenic and non-toxigenic strains of C. difficile, thus making it an excellent marker for determining the presence of this microorganism itself⁴.

The Society for Healthcare Epidemiology of America (SHEA), the Infectious Diseases Society of America (IDSA) and the European Society of Clinical Microbiology and Infectious Diseases (ESCMID) recommend different testing approaches for the laboratory detection of toxigenic C. difficile⁵, including different combinations of GDH, toxins and NAAT assays. One of the suggested testing algorithms consists in the use of a two-step protocol, involving an initial screening of the sample with a GDH test and, if the result is positive, use of a second test to confirm the toxigenicity of the strain, as only the toxigenic and pathologic strains of this microorganism produce them. Immunocard STAT! C. difficile GDH-AB enables a lab to adopt this approach, obtaining the two results simultaneously with a single-step assay.

A positive result in the test for the glutamate dehydrogenase of C. difficile confirms the presence of this organism in a fecal specimen; a negative result indicates the absence of the organism. A positive result in the test for toxins A and/or B confirms the presence of toxigenic C. difficile.

BIOLOGICAL PRINCIPLES

The test Immunocard STAT! C. difficile GDH-AB contains two strips placed in a double cassette, one for GDH and one for Toxins A & B:

1. GDH strip uses a combination of:
 - a. Pink-Red latex particles conjugated to specific antibodies against C. difficile GDH and GDH specific antibodies immobilized in the Test Position, below the Control Band. During incubation, if GDH antigen is present, is bound to both the antibodies, developing a pink-red line at Test Position.
 - b. Blue latex particles conjugated to an antigen recognized by a specific antibody for this antigen, bound to the membrane, defining the Control Band of the test.
2. TOXINS strip uses a combination of:
 - a. Pink-Red latex particles conjugated to specific antibodies against toxin B and antibodies specific for toxin B immobilized in B Position, below the Control Band. During incubation, if Toxin B is present, is bound to both the antibodies, developing a pink-red line at B Position.
 - b. Pink-Red latex particles conjugated to specific antibodies against toxin A and antibodies specific for toxin A immobilized in A position, below the toxin B band. During incubation, if Toxin A is present, is bound to both the antibodies, developing a pink-red line at A Position.
 - c. Blue latex particles conjugated to an antigen recognized by a specific antibody for this antigen, bound to the membrane, defining the Control band of the test.

A diluted patient stool sample is dispensed into both sample ports of the cassette and migrates along the membranes through the Test and Control zones. After 15 minutes of incubation at room temperature, the appearance of a specific colored line in the reading window next to the corresponding letter indicates a positive result in the presence of the Control line (see Fig. 1).

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this kit is listed on the outer box.

1. Immunocard STAT! C. difficile GDH-AB device: Two reactive strips housed in a plastic frame and enclosed in a foil pouch with a desiccant.
2. Sample Diluent vials: Buffered solution containing 0.095% of sodium azide as preservative. The diluent is supplied in a single use plastic dropper vial with an applicator stick. Use as supplied.
3. Disposable Transfer pipettes

MATERIALS NOT PROVIDED

1. External Control Set with Positive and Negative Controls. Meridian Cat. No 750501
2. Vortex
3. Interval Timer

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only. For professional use only.
2. Do not deviate from the method described here or falsely positive or falsely negative results may occur. Once the assay has started, complete all subsequent steps without interruption.

3. Patient specimen and used Immunocard STAT! C. *difficile* GDH-AB device may contain infectious agents and should be handled at Biosafety Level 2 as recommended in the CDC/NIH manual "Biosafety in Microbiology and Biomedical Laboratories".
4. Do not interchange reagents from different kit lot numbers and do not use expired reagents.
5. Do not use the sample dilution buffer if there is evidence of contamination or precipitation.
6. The sample dilution buffer contains sodium azide which is a skin irritant. Avoid skin contact with reagents. Disposal of reagents containing sodium azide into lead or copper plumbing can result in the formation of explosive metal azides. This can be avoided by flushing with a large volume of water during such disposal.
7. Correct stool storage and correct stool dilution are essential to ensure correct results. Over-inoculation of stool into the Sample Diluent may restrict the flow within the Immunocard STAT! C. *difficile* GDH-AB device, determining invalid results. Incorrect stool storage or under-inoculation into the Sample Diluent may determine false negative results.
8. In case the primary packaging is damaged (foil pouch or diluent buffer vial) the product should be discarded and not used.
9. Do not use this product if a colored line appears in the result area of any strip before you start to use it.
10. If the test is stored refrigerated, allow all the kit components and faecal samples reach room temperature, because cold reagents and/or samples can reduce test functionality.
11. Do not discard the kit box until all the content has been used. This box contains essential information about the CE mark of the product and the batch number.
12. The used product should be discarded in compliance with current local laws

HAZARD AND PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

SHELF LIFE AND STORAGE

Store the Immunocard STAT! C. *difficile* GDH-AB kit at 2-30 C when not in use. The kit expiration date is indicated on the kit label.

PROCEDURAL NOTES

1. Allow kit components and specimens to reach room temperature (19-27 C) before performing a test, as cold reagents and/or specimens may decrease assay sensitivity. Reagents may take up to 60 minutes to warm up following refrigeration.
2. Stool samples must be mixed thoroughly (regardless of consistency) prior to sampling to ensure a representative sample.
3. Hold reagent vials vertically when dispensing drops to ensure consistent drop size and delivery.
4. On occasion, particulate matter may interfere with sample flow. In cases when the Test Device does not readily absorb the diluted specimen, gently touch the bottom of the sample port with an applicator stick, moving the stool solid particle that might prevent the absorption. Alternatively, a new aliquot of the sample can be withdrawn from the Diluent and retested. Diluted samples containing a heavy concentration of particulate matter may be centrifuged (1-5 minutes at 700 x G) or allowed to stand for 3-5 minutes before proceeding.
5. Ensure to add the appropriate amount of sample (about 110 mg of stool) to the Sample Dilution Vial. If the sample is semi-liquid (difficult to be dispensed with a pipette) take a small portion of stool of about 5 mm diameter (quantity adequate for covering the grooves of the stick in the Sample Diluent cap). If the sample is liquid, add 110 µL to the Sample Diluent vial (4 drops if the disposable pipettes provided with the kit are used). Do not test any kind of specimen without the proper dilution in the Sample Diluent.
6. Ensure to add the correct volume of the extracted sample to the two sample ports marked with an arrow in the plastic device. If the volume is lower than indicated, the flow may not occur due to insufficient sample reaching the reaction strips; if the volume is higher, brown lines may appear instead of the expected ones (see Fig. 1).

REAGENT PREPARATION

Reagents are supplied ready to use. Allow kit components and specimens to reach room temperature (19-27 C) prior to use. Gently mix liquid reagent prior to use. Open the device pouch only when ready to run the assay.

SPECIMEN COLLECTION AND STORAGE

The test is validated for fresh untreated samples. Do not use samples collected in transport media, or those with added preservative agents (such as formalin, SAF, PVA or similar) or enrichment media, as their presence could interfere with correct performance of the test.

Stool samples should be transported in an airtight container. The sample should be tested as soon as possible but may be held up to 2 days at 2-8 C prior to testing. If testing cannot be performed within this timeframe, samples should be frozen immediately upon receipt and stored frozen (-20 C) until tested. Samples may be frozen and thawed twice. Ensure the samples have reached room temperature before testing.

TEST PROCEDURE

Bring all test cards, sample diluents and samples to room temperature (19-27 C) before testing. Remove the reagents from the kit box to warm. Reagents may take up to 60 minutes to warm following refrigeration.

1. Label one Sample Diluent Vial for each patient sample to be tested.
2. Unscrew the cap from the vial with caution in order not to spill the sample diluent buffer.
3. Immediately add stool sample or controls as follows:
 - a. Formed/Solid stools – Mix the stool sample thoroughly. Using the stick attached to the vial cap, collect a small portion of 5 mm diameter.
 - b. Semi-solid stools - Mix the stool sample thoroughly. Using the stick attached to the vial cap, collect a sample amount that completely covers the grooves of the stick.
 - c. Liquid Stools - Using the disposable pipettes provided in the kit, mix the stool sample thoroughly, collect and dispense in the vial 4 drops of stool (corresponding to a volume of a 110 µL).
 - d. External Positive or Negative Control: refer to EXTERNAL CONTROL SET Instructions for Use.
4. Carefully add the sample into the corresponding vial containing the dilution buffer. Screw the cap firmly and shake it vigorously to ensure mixture homogenization.
5. Use 1 Immunocard STAT! Test Card for each Sample or Control. When ready to perform testing, remove the Test Card from its foil pouch. Discard the pouch and desiccant. Label the device with the name of the patient or the control.
6. Break the top of the vial cap using a piece of paper to prevent leakage.
7. Invert the vial and, keeping it in a vertical position, add 4 drops of diluted sample in the Sample Port of each strip (rectangular windows marked with an arrow).
8. Incubate the test card at 19-27 C for 15 minutes.
9. Visually read the results of each card within 30 seconds at the end of incubation.

INTERPRETATION OF RESULTS

The cassette contains two strips: on the left side the strip for GDH, on the right side the strip for toxins A and B. The possible outcomes are shown in Fig. 1.


Fig. 1. Patterns of possible results

NEGATIVE RESULTS: A **BLUE** band at the Control Line position (C). No other bands are present in the two strips of the cassette.

A negative result in the GDH antigen strip indicates *C. difficile* antigen either is absent in the specimen or is below the detection limit of the assay. A negative result in the Toxins strip indicates either *C. difficile* toxins are absent in the specimen or they are below the detection limit of the assay.

POSITIVE RESULTS FOR GDH: A **PINK-RED** band at the TEST (T) Line Position, in the presence of a **BLUE** band at the Control Line Position (C) in the LEFT STRIP of the Cassette. A positive result indicates the presence of *C. difficile* in the sample.

POSITIVE RESULTS FOR TOXINS: A **PINK-RED** band at the TOXIN A (A) and/or the TOXIN B (B) Line Position, in the presence of a **BLUE** band at the Control Line Position (C) in the RIGHT STRIP of the Cassette. A positive result indicates the presence of *C. difficile* toxin(s).

Read immediately at the end of the Incubation period. Any **PINK-RED** band appearing at either the GDH, toxin A or toxin B test line position of the device after extended incubation times shall not be reported as positive.

INVALID RESULTS:

The test shall be considered invalid, and the results cannot be reported in case of:

- ABSENCE OF THE CONTROL LINES (**BLUE**) at the designated position for the Control Line (C). The test is invalid since the absence of a control band indicates that the test procedure was performed improperly, or deterioration of reagents has occurred.
- Presence of at least one line of the **WRONG COLOR** (for example: dark red, brown-red, and grey lines instead of **PINK-RED**). Bands with colors other than **PINK-RED** may indicate reagent deterioration, procedural mistakes, or the presence of interfering substances.

If any result is INVALID or difficult to interpret, the test shall be repeated on a new preparation of the with the same specimen (use a new Sample Diluent vial). If the original sample repeatedly produces unreadable results, obtain a new sample and retest.

IMPORTANT ADVICE:

Low percentage of specimens may result **negative for the GDH antigen but positive for the toxin(s)**. In such events, the test shall be considered "**not reportable**", and repeated after a new preparation of the original stool specimen (see above).

If the result does not change, the test must be repeated on a newly obtained stool specimen or it must be considered the use of an alternative additional diagnostic test, such as a NAAT assay for Tox A/Tox B genes detection on the original sample.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

At the time of each use, kit components should be visually examined for obvious signs of microbial contamination, freezing, leakage or any damage.

INTERNAL CONTROLS: Internal controls are contained within the test strip and therefore are evaluated with each test.

- The **BLUE** bands appearing at the Control Line Positions serve as a procedural control and indicate that the test has been performed correctly, proper flow occurred, and the test reagents were active at the time of use.
- A clean background around the Control and Test Lines also serves as a procedural control. Control or test lines that are obscured by heavy background color may invalidate the test and may be an indication of reagent deterioration, use of an inappropriate sample or improper test performance.

EXTERNAL CONTROLS: External Positive and Negative Controls should be assayed with each new kit lot or new shipment. The number of tests performed with the external controls will be determined by the requirements of local, state, federal regulations or accrediting agencies. External controls are used to monitor reagent reactivity and test performance. Failure of the controls to produce the expected results can mean that one of the reagents or components is no longer reactive at the time of use, the test was not performed correctly or that reagents or samples were not added.

The results expected with the controls are described in the section INTERPRETATION OF RESULTS. If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. The kit should not be used if control tests do not produce the correct results.

EXPECTED VALUES

The frequency of antibiotic-associated diarrhea caused by *C. difficile* is dependent on several factors, including patient population, type of institution and epidemiology. The reported incidence of *C. difficile*-associated disease in patients suspected of having antibiotic-associated diseases is 15-20%, although different facilities may find positive rates above or below this range.

LIMITATIONS OF THE PROCEDURE

1. The assay is intended to confirm the presence of the GDH common antigen and toxins in patient's stool. A Negative Test result does not completely preclude an infection of a toxigenic *C. difficile* strain. Assay results are to be used in conjunction with information available from the patient clinical evaluation and other diagnostic procedures.
2. The test is qualitative and no quantitative interpretation should be made with respect to the intensity of the positive line when reporting the result.
3. The test is intended to be used on unpreserved human stool; any other kind of sample has not been validated.
4. Correct stool preservation is essential for reliable results. Sample degradation may determine False Negative results. While relatively stable at 2-8 C, *C. difficile* toxins - particularly at low concentrations - easily degrade at room temperatures. The rate at which the toxins degrade differs from patient sample to patient sample and therefore the rate cannot be predicted. For this reason, best practice requires that samples be refrigerated or frozen within two hours of collection and the samples tested within the timeframes recommended in this insert. Do not accept samples that have not been collected, handled or transported properly.
5. Two distinct groups have been identified that can harbor *C. difficile* asymptotically at very high rates. Colonization rates of up to 50% and higher have been reported in infants and rates of up to 32% in cystic fibrosis patients.
6. Failure to add the correct amount of stool to the Sample Diluent Vial can lead to falsely negative or falsely positive results.
7. Over incubation of tests may lead to false-positive test results. Incubating tests at reduced temperatures or times may lead to falsely negative results.
8. Highly hemorrhagic samples may interfere with the assay by determining a false negative result or the appearance of aspecific bands. This event is often accompanied by the alteration of the color of the bands. Do not report any result if the bands are not of the correct color (Control Lines must be **BLUE**, Test Lines must be **PINK-RED**) – Refer to Fig. 1.
9. International Guidelines recommend avoiding repeated testing and test of cure with any kind of *C. difficile* assay, because patients can shed spores, DNA and/or toxins of *C. difficile* for a prolonged time after resolution of diarrhoea or clinical recovery. For this reason, the performances of the test cannot be established to verify the efficacy of the eradication (test of cure).

PERFORMANCE CHARACTERISTICS

Clinical performances of Immunocard STAT! *C. difficile* GDH-AB have been evaluated on a total of 304 samples for GDH and 180 samples for Toxins A and B. Specimen have been collected retrospectively and stored frozen. Two independent evaluations have been performed: one in the Manufacturer's Laboratories, one by an external independent clinical laboratory.

The samples have been characterized using commercial immunoassays for GDH and for TOXINS A&B. Discrepant results have been resolved with additional testing, including the gold standard "cell cytotoxicity & neutralization assay".

The Positive Agreement is used to assess the Sensitivity of the Immunocard STAT! *C.difficile* GDH-AB assay, while the Negative Agreement is used to assess its Specificity.

Comparative results below:

TABLE 1 – Manufacturer's Internal Evaluation

| COMPARATOR: ELISA C. DIFF CHECK – 60 | | | |
|---|-------|-------------|-------------|
| Immunocard STAT! <i>C.difficile</i> GDH-AB GDH assay | POS | NEG | TOTAL |
| POSITIVE | 53 | 0 | 53 |
| NEGATIVE | 0 | 87 | 87 |
| Total | 53 | 87 | 140 |
| | | % | CI 95% |
| Positive Agreement (Clinical sensitivity) | 53/53 | 100 | 93.2 – 100 |
| Negative Agreement (Clinical specificity) | 87/87 | 100 | 95.7 – 100 |
| COMPARATOR: ELISA PROSPECT | | | |
| Immunocard STAT! <i>C.difficile</i> GDH-AB TOXINS assay | POS | NEG | TOTAL |
| POSITIVE | 45 | 0 | 45 |
| NEGATIVE | 6* | 85 | 91 |
| Total | 51 | 85 | 136 |
| | | % | CI 95% |
| Positive Agreement (Clinical sensitivity) | 45/51 | 88.2 | 76.6 – 94.5 |
| Negative Agreement (Clinical specificity) | 85/85 | 100 | 95.6 – 100 |

*POSITIVE RESULTS confirmed by CECTNA (cell cytotoxicity & neutralization assay)

TABLE 2 – Independent Laboratory Evaluation

| Immunocard STAT1 C.difficile GDH-AB GDH assay | COMPARATOR: C. DIFF QUIK CHECK® | | |
|--|---|-------------|-------------|
| | POS | NEG | TOTAL |
| POSITIVE | 42 | 1** | 43 |
| NEGATIVE | 4* | 117 | 121 |
| Total | 46 | 118 | 164 |
| | | % | CI 95% |
| Positive Agreement (Clinical sensitivity) | 42/46 | 91.3 | 79.7 – 96.6 |
| Negative Agreement (Clinical specificity) | 117/118 | 99.2 | 95.4 – 99.8 |
| Immunocard STAT1 C.difficile GDH-AB TOXINS assay | COMPARATOR: C. DIFF QUIK CHECK® COMPLETE | | |
| | POS | NEG | TOTAL |
| POSITIVE | 19 | 0 | 19 |
| NEGATIVE | 5*** | 20 | 25 |
| Total | 24 | 20 | 44 |
| | | % | CI 95% |
| Positive Agreement (Clinical sensitivity) | 19/24 | 79.2 | 59.5 – 90.7 |
| Negative Agreement (Clinical specificity) | 20/20 | 100 | 83.9 – 100 |

RESOLUTION OF DISCREPANT RESULTS

The discrepant samples were analysed using the Techlab ELISA C. DIFF CHECKTM 60.

*Of the 4 samples Immunocard STAT GDH NEG/COMPARATOR POS, 2 were resolved as NEGATIVE, 1 was confirmed as POSITIVE.

** The Immunocard STAT GDH POS/ COMPARATOR NEG was resolved as NEGATIVE.

*** Of the 5 Immunocard STAT TOXINS NEG/COMPARATOR POS, 2 were resolved as NEGATIVE with the third assay (ELISA), 2 showed a Toxin concentration next to the LoD of the assay, 1 was POSITIVE.

ANALYTICAL SENSITIVITY

Limit of Detection of the assay (LoD) has been determined as follows:

- GDH assay: serial dilutions of *C. difficile* glutamate dehydrogenase from TCG BIOMICs was used to determine the LoD value. A mean value of **0.195 ng/mL** was obtained.
- TOX A&B assay: serial dilutions of toxins A and B from TCG BIOMICs were used to determine the LoD value. A mean value of **1.56 ng/mL** for toxin A and a mean value of **0.75 ng/mL** for toxin B were obtained.

REPRODUCIBILITY

INTER-DAY

Using the same batch of the Simple GDH-Toxins test, a sensitivity curve was measured, and real samples established as NC ("negative sample"), LPC ("weak positive sample") and PC ("positive sample") for each analyte, in duplicate, over a five-day period spaced outover time. The results were very reproducible obtaining maximum differences of ½ dilution.

INTER-OPERATOR PRECISION

In triplicate, three operators measured a sensitivity curve for each analyte and real samples established as NC ("negative sample"), LPC ("weak positive sample") and PC ("positive sample") for each analyte. The observed differences were not over ½ dilution with the standard in any of the cases, and the same results were obtained with the real samples.

INTER-BATCH PRECISION

Three different batches of the Simple GDH-Toxins test were evaluated in parallel by analysing the sensitivity curves and real samples established as NC ("negative sample"), LPC ("weak positive sample") and PC ("positivesample") for each analyte. The analysis was performed by one person on the same day. The maximum differences were less than or equal to ½ dilution with the standard, and the same results were obtained with the real samples, providing evidence of the test's high inter-batch precision.

The differences found in the "Reproducibility" sections are acceptable for a qualitative immunochromatographic technique with its inherent variability.

PROZONE / HOOK EFFECT

Very high concentrations of the three analytes detected by the assay have been tested without observing any decrease in the intensity of the positive signals. These concentration values (higher than the maximum values that can be found among the population) are the following ones:

- GDH: 6000 ng/mL, about 30000-fold the assay LoD.
- Toxin A: 3000 ng/mL, about 1900-fold the assay LoD.
- Toxin B: 3000 ng/mL, about 4000-fold the assay LoD.

INTERFERING SUBSTANCES

The substances indicated in the below at the concentration specified did not interfere with the results of the test when added to stool samples (positive and negative ones):

Whole blood 20% (v/v)
Vancomycin 1.17 mg/mL
Metronidazole 1.17 mg/mL
Pepto-Bismol 1% (v/v)
Omeprazole 11.7 µg/mL
Cimetidine 7.38 mg/mL
Paracetamol 2.34 mg/mL

Ibuprofen 1.41 mg/mL
Loperamide 9.38 µg/mL
Acetylsalicylic acid 2.34 mg/mL
Barium sulphate 4.4 mg/mL
Mucin 5% (w/w)
Stearic-palmitic acid (1:1)

CROSS REACTIVITY

Cross reactivity of the assay was evaluated for different bacteria that can be present in the intestinal tract. Two independent evaluations have been performed: one in the Manufacturer's Laboratories, one by an external independent clinical laboratory.

Manufacturer's Internal Evaluation

The test did not present cross-reactivity with any of the bacterial suspensions listed below:

| | |
|---------------------------------------|----------------------------------|
| <i>Escherichia coli</i> | 0.25 x 10 ⁸ (CFU/mL) |
| <i>Salmonella</i> spp. | 0.25 x 10 ⁸ (CFU/mL) |
| <i>Campylobacter coli</i> | 0.25 x 10 ⁸ (CFU/mL) |
| <i>Campylobacter jejuni</i> | 0.25 x 10 ⁸ (CFU/mL) |
| <i>Yersinia enterocolitica</i> | 0.25 x 10 ⁸ (CFU/mL) |
| <i>Entamoeba histolytica</i> | 1.56 x 10 ⁵ (CFU/mL) |
| <i>Giardia lamblia</i> (trophozoites) | 5.5 x 10 ⁶ (cells/mL) |
| <i>Giardia lamblia</i> (cysts) | 1.7 x 10 ⁶ (cells/mL) |

The test was evaluated against stools strongly positive for the following microorganisms without any cross-reactivity observed: Adenovirus, Rotavirus, Norovirus GII, Astrovirus, *Campylobacter* spp, *Helicobacter pylori*, *Giardia lamblia*, *Cryptosporidium parvum*.

A cross-reaction with a stool sample strongly positive for *Entamoeba histolytica* has been registered (a weak positive result for Toxin B). However, further additional testing did not confirm the cross reaction.

Independent External Laboratory Evaluation

Stool samples inoculated with the following microbial agents (to a final sample concentration of ~ 1 x 10⁸ organisms/mL)

Aeromonas caviae, *Aeromonas hydrophila*, *Bacillus cereus*, *Bacillus* spp, *Bacteroides nordii*, *Campylobacter jejuni*, *Candida albicans*, *Clostridium butyricum*, *Clostridium cadaveris*, *Clostridium perfringens*, *Clostridium sordellii*, *Enterobacter aerogenes*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Enterococcus faecium*, *Escherichia coli* 1, *Escherichia coli* 2, *Klebsiella pneumoniae*, *Lactobacillus gossleri*, *Listeria monocitogenes*, *Plesiomonas shigelloides*, *Proteus mirabilis*, *Pseudomonas aeruginosa*, *Salmonella enteritidis*, *Salmonella typhi*, *Salmonella typhimurium*, *Serratia marcescens*, *Shigella dysenteriae*, *Shigella flexnerii*, *Shigella sonnei*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Vibrio cholerae*, *Vibrio parahaemolyticus*, *Yersinia enterocolitica*.

NONE of the listed organisms cross reacted with the assay.

INFORMATION & SUPPORT

For any technical inquiry and complaint, contact Meridian Bioscience Technical Support at +390331433636 – email: TS@meridianbioscience.com or call your Country Distributor (www.meridianbioscience.com/diagnostics/distributors/)

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REFERENCES

1. Rupnik M, Dupuy B, Fairweather NF, Gerding DN, Johnson S, Just I, Lyerly DM, Popoff MR, Rood JI, Sonenshein AL, Thelestam M, Wren BW, Wilkins TD and von Eichel-Streiber C. Revised nomenclature of Clostridium difficile toxins and associated genes. J Med Microbiol 2005 Feb 54(2):113-117.
2. Vonberg RP, Reichardt C, Behnke M, Schwab F, Zindler S and Gastmeier P. Costs of nosocomial Clostridium difficile associated diarrhoea. J Hosp Infect 2008 Sep 70(1):15-20.
3. Kuehne SA, Cartman ST, Heap JT, Kelly ML, Cockayne A and Minton NP. The role of toxin A and toxin B in Clostridium difficile infection. Nature 2010 Oct 467(7316):711-713.
4. Shetty N, Wren MW, Coen PG. The role of glutamate dehydrogenase for the detection of Clostridium difficile in faecal samples: a meta-analysis. J Hosp Infect 2011 Jan 77(1):1-6.
5. Cohen SH, Gerding DN, Johnson S, Kelly CP, Loo VG, McDonald LC, Pepin J, Wilcox MH; Society for Healthcare Epidemiology of America; Infectious Diseases Society of America. Clinical practice guidelines for Clostridium difficile infection in adults: 2010 update by the society for healthcare epidemiology of America (SHEA) and the infectious diseases society of America (IDSA). Infect. Control Hosp. Epidemiol. 2010 May 31(5):431-55.
6. Wilcox MH, Planche T, Fang FC and Gilligan P. What is the current role of algorithmic approaches for diagnosis of Clostridium difficile infection. J Clin Microbiol 2010 Dec 48(12):4347-53.
7. Wren MW, Sivapalan M, Kinson R and Shetty NR. Laboratory diagnosis of Clostridium difficile infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory. Br J Biomed Sci 2009;66 (1):1-5.
8. Crobach MJ, Planche T, Eckert C, Barbut F, Terveer EM, Dekkers OM, Wilcox MH, Kuijper EJ European Society of Clinical Microbiology and Infectious Diseases: update of the diagnostic guidance document for Clostridium difficile infection. Clin Microbiol Infect 2016 Aug 22(Suppl 4):S63-81.

SN750520IFU

REV. 06/23

| | |
|---|--|
|  | <p>Meridian Bioscience Europe SRL Via Dell'industria 7 – 20035 Villa Cortese (Milano) – Italy Telephone: (+39) 0331433636 - E-mail: info@meridianbioscience.eu www.meridianbioscience.com</p> |
| | <p>Product Support: mbe-techservice@meridianbioscience.eu For IFUs & MSDS: www.meridianbioscience.com/diagnostics</p> |
|  <p>Authorized Representative in Switzerland</p> | <p>MedEnvoy Switzerland Gotthardstrasse 28 6302 Zug Switzerland</p> |

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simbolos, Zeichenerklärung)

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