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LITERATURA

- Rupnik M, Dupuy B, Fairweather NF, Gerding DN, Johnson S, Just I, Lyster DM, Popoff MR, Rood JI, Sonenshein AL, Thelestam M, Wren BW, Wilkins TD and von Eichel-Streiber C. Revised nomenclature of *Clostridium difficile* toxins and associated genes. J Med Microbiol 2005 Feb 54(2):113-117.
- Vonberg RP, Reichardt C, Behnke M, Schwab F, Zindler S and Gastmeier P. Costs of nosocomial *Clostridium difficile* associated diarrhoea. J Hosp Infect 2008 Sep 70(1):15-20.
- Kuehne SA, Cartman ST, Heap JT, Kelly ML, Cockayne A and Minton NP. The role of toxin A and toxin B in *Clostridium difficile* infection. Nature 2010 Oct 467(7316):711-713.
- Shetty N, Wren MW, Coen PG. The role of glutamate dehydrogenase for the detection of *Clostridium difficile* in faecal samples: a meta-analysis. J Hosp Infect 2011 Jan 77(1):1-6.
- Cohen SH, Gerding DN, Johnson S, Kelly CP, Loo VG, McDonald LC, Pepin J, Wilcox MH; Society for Healthcare Epidemiology of America; Infectious Diseases Society of America. Clinical practice guidelines for *Clostridium difficile* infection in adults: 2010 update by the society for healthcare epidemiology of America (SHEA) and the infectious diseases society of America (IDSA). Infect. Control Hosp. Epidemiol. 2010 May 31(5):431-55.
- Wilcox MH, Planche T, Fang FC and Gilligan P. What is the current role of algorithmic approaches for diagnosis of *Clostridium difficile* infection. J Clin Microbiol 2010 Dec 48(12):4347-53.
- Wren MW, Sivapalan M, Kinson R and Shetty NR. Laboratory diagnosis of *Clostridium difficile* infection. An evaluation of tests for faecal toxin, glutamate dehydrogenase, lactoferrin and toxigenic culture in the diagnostic laboratory. Br J Biomed Sci 2009;66 (1):1-5.
- Crobach MJ, Planche T, Eckert C, Barbut F, Terveer EM, Dekkers OM, Wilcox MH, Kuijper EJ European Society of Clinical Microbiology and Infectious Diseases: update of the diagnostic guidance document for *Clostridium difficile* infection. Clin Microbiol Infect 2016 Aug 22(Suppl 4):S63-81.

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