FLU/RSV POSITIVE CONTROL

External Positive Control Reagent for TRU FLU and TRU RSV

REF **751110**

IVD

R_x Only

INTENDED USE

FLU/RSV Positive Control is an external control reagent to be used with TRU FLU and TRU RSV Test Kits as part of a routine quality control program.

SUMMARY AND EXPLANATION OF THE TEST

Quality control testing is performed to detect factors such as reagent deterioration, adverse environmental or test conditions, or variance in operator performance that can cause test errors. External control reagents, such as FLU/RSV Positive Control, are reagents that are not built into the test system, but are sampled and tested in the same manner as patient specimens.

BIOLOGICAL PRINCIPLES

To verify that TRU FLU and/or TRU RSV test kits are performing correctly, they should be tested with known positive and negative external control reagents on a periodic basis. FLU/RSV Positive Control is used in conjunction with Sample Diluent provided with each test kit to create a known positive and negative control panel. Alternatively, a confirmed negative patient sample can be used as an external negative control. The frequency with which a laboratory performs external controls will be affected by:

- The proficiency level of the laboratory,
- The laboratory’s own internal requirements,
- The requirements of the laboratory’s accrediting agencies,
- The number of new and different operators performing the test,
- Whether a new kit lot is being added to testing,
- Whether the kit lot is from a different shipment, and
- Whether deviations from the manufacturer’s stated storage or handling conditions have occurred.

When unacceptable quality control test results are obtained, all test results should be considered invalid. QC test failures are an indication that either reagents, test environment or operator performance have changed.

REAGENT/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

- FLU/RSV Positive Control: Inactivated influenza A, influenza B and RSV antigens in a buffered diluent containing sodium azide (0.095%) as a preservative. Use as supplied. Store at 2-8 C when not in use.
- FLU/RSV Positive Control in plastic dropper vial. Each vial is sufficient to perform 8 TRU QC tests.

MATERIALS NOT PROVIDED

- TRU FLU (Product code 751230) and/or TRU RSV (product code 751330).
- Disposable latex gloves (Respiratory samples are considered potentially biologically hazardous material.)

PRECAUTIONS

- All reagents are for in vitro diagnostic use only.
- This is a quality control reagent and is used to evaluate the performance of TRU FLU and TRU RSV. It is not directly used to test patient samples.
- Transport and store this reagent at 2-8 C when not in use. Do not freeze.
- This reagent can be used with any TRU FLU or TRU RSV kit lot numbers.
- Do not use if discolored or turbid. Discoloration or turbidity may be a sign of microbial contamination.
- Hold reagent vials vertically when dispensing drops to ensure consistent drop size and delivery.
- FLU/RSV Positive Control contains inactivated influenza A, influenza B and RSV virus antigens and should be handled as if it were potentially infectious.

HAZARD and PRECAUTIONARY STATEMENTS

There are no known hazards associated with this product.

QUALITY CONTROL TEST PROCEDURE

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

- Bring FLU/RSV Positive Control and all test kit components (test devices or strips, reagents and samples) to room temperature (20-25 C) before testing.
- Use one Test Strip (TRU FLU, TRU RSV) for each positive and negative control reagent to be tested. Remove each Test Strip from its foil pouch. Label each device or strip with the name of the control to be tested.
- Perform QC tests according to the external QC instructions described in the TRU FLU or TRU RSV insert.
- At the end of testing, read reactions as described in each product’s insert.

LIMITATIONS OF THE PROCEDURE

The Positive Control is manufactured in an aqueous solution matrix. Although specimen matrix interference has not been observed with this assay, the aqueous matrix of the controls may not adequately control for specimen matrix effects. If the user wishes to supply controls in the sample matrix, the user is referred to the Clinical and Laboratory Standards Institute guideline EP14-A2, Evaluation of Matrix Effects: Approved Guideline.

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