Rapid Strip ROTA-ADENO

Rapid Immunoassay for the Detection of
of Rotavirus and Adenovirus Antigens in Stool Specimens

REF

IVD

751120

For in vitro diagnostic use only

ENGLISH p. 2

Rapid Immunochromatographic assay for the Detection of Rotavirus and
Adenovirus Antigens in Stool Specimens

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meridian BIOSCIENCE®
LIFE DISCOVERED. LIFE DIAGNOSED.

ENGLISH

Rapid Strip Rota-Adeno

REF 751120 – 20 tests

A Rapid Immunochromatographic assay for the Detection of Rotavirus and Adenovirus Antigens in Stool Specimens

INTENDED USE

The Rapid Strip Rota -Adeno is an *in vitro* qualitative procedure for the detection of Rotavirus and Adenovirus antigens in human stool.

EXPLANATION

Rotavirus is the primary causative agent of acute gastro-enteritis, especially in children less than 2 years old. Its discovery in 1973 and its association with infantile gastro-enteritis represented a very important advance in the study of gastro-enteritis not caused by acute bacterial infection. Rotavirus is transmitted by oral-faecal contact with an incubation period of 1-3 days. Although sample collections taken within the second and fifth day of the illness are ideal for antigen detection, the rotavirus may still be found whilst diarrhoea continues.

Because rotavirus is extremely difficult to culture, it is unusual to use isolation of the virus in diagnosing an infection. Instead, a variety of techniques have been developed to direct detect rotavirus in faeces. The Rapid Strip Rota-Adeno provides a simple, direct detection assay of rotavirus in faeces, which is highly sensitive, specific and rapid.

Rotavirus has been described as the most common causing agent of diarrhoea in the world. It has been estimated that in the developing countries the infection by rotavirus is causing one million deaths per annum. It is also responsible for 20 to 25% of the deaths by diarrhoea, and it is causing 6% of the deaths in children under 5 years (Cook et al, Bull WHO, 1990;68:171-7).

The Adenovirus is the second most common cause of viral gastro-enteritis in children (10-15%). This virus may also cause respiratory diseases and, depending on the serotype, also diarrhoea, conjunctivitis, cystitis, etc. At least 47 serotypes of adenovirus have been described, all sharing a common hexon antigen. Serotypes 40 and 41 are the ones associated with gastro-enteritis, whose main symptom is diarrhoea that may last between 9 and 12 days associated with fever and vomit.

The diagnosis of these two viruses is difficult by culture, therefore the detection of viral antigen particles in faeces samples is becoming a common diagnostic technique. The **Rapid Strip Rota-Adeno** is a rapid 5-minute immunoassay, based on a lateral flow chromatography technique that detects *Rota* and *Adenovirus* antigens present in human stool.

BIOLOGICAL PRINCIPLES

The rapid strip Rota-Adeno employs a combination of 1) pink-red latex-conjugated monoclonal antibodies against antigen VP6 of group A of rotavirus, and solid-phase specific rotavirus antibodies. 2) blue latex-conjugated monoclonal antibodies against the adenovirus hexon antigen, (which is present in all the human subtypes of adenovirus) and solid-phase specific adenovirus antibodies.

In this test the specimen is first treated with a sample diluent to extract rotavirus and adenovirus antigens from the faeces. Following extraction, the only step required is to add the extract and the strip in a tube.

As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured particles. Different coloured lines will be visible, depending upon the virus content of the sample. These lines, after 5 minutes of incubation at room temperature, are used to interpret the result.

MATERIALS PROVIDED

The kit contains all the necessary material to perform 20 tests:

- Reaction strips (20) individually foil pouched, containing immobilized anti-*Rota* and anti-*Adeno* antibody (test zone).
- One Sample diluent Bottle (30 ml)

MATERIALS NOT PROVIDED:

1. Disposable latex gloves, that should be used during the handling of the faecal samples as they are considered potential hazard material.
2. Swabs/applicators for stool sample collection
3. Vortex for suspension of the stool specimen in the extraction buffer.
4. Test tubes for the extraction
5. Transfer pipettes
6. Timer

PRECAUTIONS

1. All the reagents are for "In Vitro" diagnostic use only.
2. Do not use the test beyond its expiry date.
3. Do not interchange reagents from different kit lot numbers.
4. Allow kit components and specimens to reach the room temperature before use, as cold reagents and/or specimens may decrease assay performances. 20-30 minutes are recommended.
5. Do not use the test if any coloured line can be seen in the test zone before performing the test.
6. The test will not perform properly if the amount of sample is less or more than required.
7. Stool must be mixed thoroughly (regardless of consistency) to insure a representative sample prior to sampling.
8. All the specimens should be handled as if they could transmit infectious agents.
9. Dispose of all used materials in a container suitable for biohazardous waste.

The sample diluent contains Sodium Azide (0,095 %) as preservative. Avoid contact with the skin.

SHELF LIFE AND STORAGE

The expiry date is printed on the pouch of each reaction strip.

The product can be stored at any temperature between 2 and 30° C. Do not freeze.

SPECIMEN HANDLING

Stool samples must be collected as soon as the symptoms appear. Viral particles decrease in number after 8 days, making the diagnosis more difficult.

The specimen should be received in an airtight transport container and stored at 2°-8°C until tested. The samples can be stored in the refrigerator for 2 days. For longer storage they must be kept frozen at -20°C. In this case, the sample should be totally thawed, and brought to room temperature and homogenised before testing.

NOTE: Stool in transport media, swabs, or preservatives are inappropriate for testing

PROCEDURE

1. Transfer 1.0 ml of sample diluent in a test tube or vial
2. Add a sample portion of approximately 5-6 mm size (30-50mg), with a swab, a wooden applicator, or a bacteriology loop. For liquid or semi-solid stools add 100 microliters of stool using an appropriate pipette. Shake gently to resuspend it into the buffer.
3. Vortex for 15 seconds
4. Wait at least 3 minutes until the solid particles settle or centrifuge for 1 minute at 700 xg and transfer with a pipette 500 microliters of supernatant to another test tube,
5. Dip the reaction strip in the second test tube with the arrow pointing to the bottom.
IMPORTANT: the liquid must not reach the blue area above the arrowheads. If needed, use a larger tube or reduce the amount of sample.
6. Incubate the test at 19-27 C for 10 minutes and then read the results as described under INTERPRETATION OF RESULTS section below. Any line that may appear after 10 minutes does not have diagnostic value.

Note: the strip can also be introduced into the first tube or vial, provided it is large enough to avoid the liquid to reach the blue area above the arrowheads. However, in rare cases, when using a too concentrated diluted specimen, the absorption might be difficult. Alternatively, the strip can be dipped for 10 seconds in the first tube, avoiding exceeding the arrowheads, and then transferred to the bench top for the incubation.

INTERPRETATION OF RESULTS

Negative test result: only one **GREEN** band (Control Line) appears across the white central area of the reaction strip. *Rota* and *Adeno* antigens are absent or below the level of detection.

Positive test result: In addition to the **GREEN** band (Control Line), a distinguishable **PINK-RED** band (indicating that there are detectable *Rotavirus* antigens in the specimen) and/or a distinguishable **BLUE** band (indicating that there are detectable *Adenovirus* antigens in the specimen) also appears across the white central zone of the reaction strip.

The intensity of the band will be variable depending on the antigen concentration in the specimen. Any pink-red or blue line, even very weak, must be considered as a positive result. Any line or colour appearing after 10 minutes has no diagnostic value.


Invalid test result: The result window remains white, there are no lines visible or there is one blue line only or one red line only. The specimen should be re-tested using a new reaction device.

Interpretation code for the images on the strip.

R: Rotavirus

A: Adenovirus

N: Negative

Final diagnosis should not be based exclusively on a single test result. It should be based on the consistency of test results with other applicable data and the clinical symptoms.

QUALITY CONTROL

If no Control Line band appears, the test is invalid, since improper test procedure was carried out or deterioration of the reagents has occurred.

WARNING: The inclusion of a control with known results is recommended to ensure that the data obtained are correct.

LIMITATIONS OF THE PROCEDURE

1. The test should be used only for the detection of rotavirus or adenovirus antigen in faecal samples.
2. The test is qualitative, and no quantitative interpretation should be made with respect to the intensity of the positive line, when reporting the result.
3. More than 200 samples were evaluated to assure the correct performance of the test. The correlation of the results with other techniques (ELISA) was excellent. However, interferences in the performance of the tests should not be excluded.
4. An excess of faecal sample may cause to appear **brown lines** instead of red and blue lines. These brown lines have no diagnostic value. In this case repeat the test with a smaller amount of sample or dilute the extract already made.
5. No cross-reactions with other viruses or substances were observed during the evaluation of the test. A negative result does not totally exclude a possible rotavirus or adenovirus infection. The significance of the results must be evaluated in relation to the patient's clinical symptoms.
6. The analysis of some samples could show lines with indefinite colour, caused mostly by negative samples. Therefore, when lines with indefinite coloration appear, the test should be repeated. If the same result is obtained, it is suggested to use another analysis method.

PERFORMANCE CHARACTERISTICS

ANALYTICAL SENSITIVITY

The lower detection limit for the Rota-Adeno strip is 31 ng/mL for Rotavirus and Adenovirus antigens.

DIAGNOSTIC SENSITIVITY AND SPECIFICITY

ASSAY 1

A CE-marked ELISA product is used as reference product. Rapid Strip Rota-Adeno is evaluated against an ELISA for the detection of Rotavirus antigens and against a second ELISA for the detection of Adenovirus antigens.

The results obtained were:

ROTAVIRUS.

		Rapid Strip	
		+	-
ELISA	+	12	0
	-	2	88

Concordance percentage: $100\%[(12+88)/102] = 98,04\%$.

Positive concordance: $100\%[12 / (12+0)] = >99,9\%$.

Negative concordance: $100\%[88 / (88+2)] = 97,8\%$.

ADENOVIRUS.

		Rapid Strip	
		+	-
ELISA	+	6	0
	-	0	96

Concordance percentage: $100\%[(6+96)/102] = >99,9\%$.

Sensitivity Concordance estimate: $100\%[6 / (6+0)] = >99,9\%$.

Specificity Concordance estimate: $100\%[96 / (0+96)] = >99,9\%$.

ASSAY 2

A second reference result is used, two ELISAs for the detection of each antigen as above, approved in the Japanese market, to corroborate the results obtained.

ROTAVIRUS

		Rapid Strip	
		+	-
ELISA	+	80	1
	-	26	307

Concordance percentage: $100\%[(80+307)/414] = 93,5\%$.

Sensitivity Concordance estimate: $100\%[80 / (80+1)] = 98,8\%$.

Specificity Concordance estimate: $100\%[307 / (307+1)] = 92,2\%$.

ADENOVIRUS.

		Rapid Strip	
		+	-
ELISA	+	40	1
	-	13	360

Concordance percentage: $100\%[(40+360)/414] = 96,6\%$.

Sensitivity Concordance estimate: $100\%[40 / (40+1)] = 97,6\%$.

Specificity Concordance estimate: $100\%[360 / (360+13)] = 96,5\%$.

INTRA-ASSAY PRECISION - REPETEABILITY

Three duplicates of each concentration from the sensitivity curve are tested using three different lots and the same results are obtained.

PRECISION – REPRODUCIBILITY

INTER-DAY PRECISION

Using 1 lot of the product, ten duplicates of the sensitivity curve are performed throughout ten consecutive days. Only a difference of a half dilution is observed, acceptable and tolerable for the assay.

INTER-LAB PRECISION

Three different laboratories/operators test those same samples, presenting high precision and concordance. Only a difference of a half dilution is observed, acceptable and tolerable for the assay.

INTER-LOT PRECISION

Using 3 different lots of the product, a sensitivity curve is constructed in duplicate. The analysis is performed by the same person in the same day. Only a difference of a half dilution is observed, acceptable and tolerable for the assay.

HOOK EFFECT

The amount tested was 2,000 ng/ml Rotavirus + 2,000 ng/ml Adenovirus. The limit of the test is 31 ng/ml Rotavirus + 31 ng/ml Adenovirus. This amount is 65 times the limit of the test. One lot is tested in duplicate. No hook effect has been observed at a concentration equivalent to 65 times the limit of the test.

INTERFERING SUBSTANCES

The substances listed in the table, at the indicated concentration, didn't interfere with the results. Three lots were used for the study.

Bilirubin F: 0.9 – 1.1 mg/ml	Atropine (40 mg/dl)
Bilirubin C: 0.9 – 1.1 mg/ml	Caffeine (40 ng/dl)
Hemoglobin: 22 – 27 mg/ml	Gentisic acid (40 mg/dl)
Acetaminophen (20 mg/dl)	Glucose (2,000 mg/dl)
Acetylsalicylic acid (20 mg/dl).	Urea (4,000 mg/dl)
Ampicillin (40 mg/dl).	Uric acid (10 mg/dl)
Ascorbic acid (100 mg/dl)	

INTERFERING MICROORGANISMS

The microorganisms listed, at the indicated concentration, did not interfere with the results.

Bacteria. CONCENTRATION (McF:5(1x10E8/ml)

<i>E.coli</i>	Group G Strep	<i>H.influenzae</i>
<i>K.pneumoniae</i>	<i>S.agalactiae</i>	<i>C.albicans</i>
<i>E.faecalis</i>	<i>S.mutans</i>	<i>B.catarrhalis</i>
<i>S.aureus</i>	<i>S.pneumoniae</i>	<i>B.pertussis</i>
<i>S.epidermidis</i>	<i>S.marcescens</i>	<i>P.vulgaris</i>
<i>S.haemolyticus</i>	<i>P.aeruginosa</i>	<i>S.saprophyticus</i>
Group B Strep	<i>N.gonorrhoeae</i>	<i>S.milleri</i>
Group C Strep		<i>E.faecium</i>

Virus-1

Virus	Subtype	Final concentration (TCID ₅₀ /ml)
Influenza type A	H1N1	1.4×10 ⁶
Influenza type A	H3N2	3.0×10 ⁶
Influenza type B	Victoria	1.6×10 ⁶
Influenza type B	Yamagata	1.5×10 ⁷

Virus-2

Virus	Subtype	Final concentration (TCID ₅₀ /ml)
Echovirus	6	7.2×10 ⁵
Poliovirus	2	2.2×10 ⁷
Parainfluenza	2	2.2×10 ⁷

For Support, please contact Meridian Bioscience Technical Support at +390331433636 – email: TS@meridianbioscience.com or call your Country Distributor.

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REFERENCES / BIBLIOGRAPHIE / BIBLIOGRAFIA / BIBLIOGRAFÍA / LITERATURANGABEN

- F. Bon et al. *Prevalence of a group A rotavirus, human calicivirus, astrovirus, and adenovirus type 40 and 41 infections among children with acute gastroenteritis in Dijon, France* Journal of Clinical Microbiology, Sept. 1999, p. **3055-3058**
- Bodo R. Eing et al. *Evaluation of two enzyme immunoassays for detection of human rotaviruses in fecal specimens*, Journal of Clinical Microbiology, Dec. 2001, p.**4532-4534**
- Umesh D. Parashar et al. *Global illness and deaths caused by rotavirus disease in children*, Emerging Infectious Diseases, vol. 9, No.5, May 2003, p. **565-572**

REV. 0623

Manufactured by

Meridian Bioscience Europe S.R.L.

Via Dell'industria 7 – 20035 Villa Cortese (Milano) – Italy
Telephone: (+39) 0331433636
E-mail: info@meridianbioscience.eu

Product Support: mbe-techservice@meridianbioscience.eu

For IFUs & MSDS www.meridianbioscience.com

For technical assistance and ordering, contact Meridian Bioscience Europe or your Local Distributor (www.meridianbioscience.com/diagnostics/distributors/)

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simblos, Zeichenerklärung)

	Use By / Utilizzabile entro / Utiliser jusqu'à / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controlla positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controlla negativo / Contrôle négatif / Control negativo / Negativer Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico diagnostico in vitro / Dispositivo medico para diagnóstico in vitro / In-Vitro-Diagnosegerät	EC REP	Authorized representative in the European Community / Rappresentante autorizzato nella Comunità Europea / Représentant autorisé dans la Communauté Européenne / Autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medici diagnostici in vitro / Ce produit répond aux exigences de la Directrice 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Esta producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnósticos in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In-Vitro-Diagnosegeräte 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus contains Sample Dilution / Dispositivo per la preparazione di campioni contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente / System für Probenaufbereitung, in dem sich Probenverdünnergriffel befindet
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht eingefrieren
	Consort Instructions for Use / Consulter les instructions pour l'usage / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampon de réaction / Solution de réaction chimique / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluations Only / Soltanto para evaluaciones de las prestaciones / Réservé à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur zur IVD Leistungsevaluierung
	Contains sufficient for "n" tests / Contiene suficiente para "n" test / Contient suffisante para "n" ensayos / Inhalt ausreichend für "n" Prüfungen	SOLN STOP	Stopping Solution / Solution de Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limite di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Conjugato enzimatico / Conjugué enzymatique / Conjunto enzimático / Enzymsconjuguat
SN	Serial number / Número de serie / Numéro de série / Número de serie / Seriennummer	CONTROL	Assay Control / Controlla del test / Test de contrôle / Control de ensayo / Kontrolltest
TEST	Test Device / Dispositivo de prueba / Dispositif de test / Dispositivo de Prueba / Testgerät	REAG	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Témpon de lavado / Waschpuffer
BUF	Buffer / Soluzione tamponi / Solution tamponnée / Tampon / Puffer		Warning / Avvertenza / Mise En Garde / Advertencia / Warnhinweise
CONJ	Conjugate / Conjugato / Conjugué / Conjugado / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Complemento / Diluant échantillons / Diluyente de muestras / Probenniedersetzungs puffer
SUBS	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentration 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución de lavado 20X / 20fach konzentriertes Waschpufferkonzentrat
R. Only	Prescription Use Only / Pour l'uso su prescrizione medica / Usuaggio sur prescription / Solo Para Use Por Receta / verschreibungspflichtig	DET REAG	Detection Reagent / Reagente Diretto / Réactif de détection / Reactivo de Detección / Nachweis Reagens
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Ne pas utiliser si le paquet est endommagé / No usar si el paquete está dañado / Nicht verwenden, wenn die Verpackung beschädigt ist	TUBE	Empty Tube / Provaia vacua / Tube vide / Tubo vacío / Leerres Gefäß