Rapid Strip ROTA-ADENO

REF 751120

For in vitro diagnostic
use only

20 tests

Meridian Bioscience Europe Srl

Via dell'Industria,7
I-20020 Villa Cortese ,Milano-Italy
Tel:+39 (0331) 433636
Fax:+39 (0331) 433616
E-mail:info@mdeur.com

Meridian Bioscience Europe salnv

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ENGLISH p. 2

Rapid Immunochromatographic assay for the Detection of Rotavirus and Adenovirus Antigens in Stool Specimens

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Rapid Strip Rota-Adeno

Catalog # 751120 – 20 tests

A Rapid Immunochromatographic assay for the Detection of Rotavirus and Adenovirus Antigens in Stool Specimens

INTENDED USE

The Rapid Strip Rota -Adeno is an *in vitro* qualitative procedure for the detection of Rotavirus and Adenovirus antigens in human stool

EXPLANATION

Rotavirus is the primary causative agent of acute gastro-enteritis, especially in children less than 2 years old. Its discovery in 1973 and its association with infantile gastro-enteritis represented a very important advance in the study of gastro-enteritis not caused by acute bacterial infection. Rotavirus is transmitted by oral-faecal contact with an incubation period of 1-3 days. Although sample collections taken within the second and fifth day of the illness are ideal for antigen detection, the rotavirus may still be found whilst diarrhoea continues.

Because rotavirus is extremely difficult to culture, it is unusual to use isolation of the virus in diagnosing an infection. Instead, a variety of techniques have been developed to direct detect rotavirus in faeces. The Rapid Strip Rota-Adeno provides a simple, direct detection assay of rotavirus in faeces, which is highly sensitive, specific and rapid.

Rotavirus has been described as the most common causing agent of diarrhoea in the world. It has been estimated that in the developing countries the infection by rotavirus is causing one million deaths per annum. It is also responsible for 20 to 25% of the deaths by diarrhoea, and it is causing 6% of the deaths in children under 5 years (Cook et al, Bull WHO, 1990;68: 171-7).

The Adenovirus is the second most common cause of viral gastro-enteritis in children (10-15%). This virus may also cause respiratory diseases and, depending on the serotype, also diarrhoea, conjunctivitis, cystitis, etc. At least 47 serotypes of adenovirus have been described, all sharing a common hexon antigen. Serotypes 40 and 41 are the ones associated with gastro-enteritis, whose main symptom is diarrhoea that may last between 9 and 12 days associated with temperature and vomits.

The diagnosis of these two viruses is difficult by culture, therefore the detection of viral antigen particles in faeces samples is becoming a common diagnostic technique. The **Rapid Strip Rota-Adeno** is a rapid 5 minutes immunoassay, based on a lateral flow chromatography technique that detects *Rota* and *Adenovirus* antigens present in human stool.

BIOLOGICAL PRINCIPLES

The rapid strip Rota-Adeno employs a combination of: 1) red latex-conjugated monoclonal antibodies against antigen VP6 of group A of rotavirus, and solid-phase specific rotavirus antibodies; 2) blue latex-conjugated monoclonal antibodies against the adenovirus hexon antigen, (which is present in all the human subtypes of adenovirus) and solid-phase specific adenovirus antibodies.

In this test the specimen is first treated with a sample diluent to extract rotavirus and adenovirus antigens from the faeces. Following extraction, the only step required is to add the extract and the strip in a tube.

As the sample flows through the test membrane, the coloured particles migrate. In the case of a positive result the specific antibodies present on the membrane will capture the coloured particles. Different coloured lines will be visible, depending upon the virus content of the sample. These lines, after 5 minutes of incubation at room temperature, are used to interpret the result.

MATERIALS PROVIDED

The kit contains all the necessary material to perform 20 tests:

- Reaction strips (20) individually foil pouched, containing immobilized anti-Rota and anti-Adeno antibody (test zone).
- One Sample diluent Bottle (30 ml)

MATERIALS NOT PROVIDED:

1. Disposable latex gloves, that should be used during the handling of the fecal samples as they are considered potential hazard material.
2. Swabs/applicators for stool sample collection
3. Vortex for suspension of the stool specimen in the extraction buffer.
4. Test tubes for the extraction
5. Transfer pipets
6. Timer

PRECAUTIONS

1. All the reagents are for "In Vitro" diagnostic use only.
2. Do not use the test beyond its expiry date.
3. Do not interchange reagents from different kit lot numbers
4. Allow kit components and specimens to reach the room temperature before use, as cold reagents and/or specimens may decrease assay performances. 20-30 minutes are recommended
5. Do not use the test if any coloured line can be seen in the test zone before performing the test.
6. The test will not perform properly if the amount of sample is too large or too small.
7. Stool must be mixed thoroughly (regardless of consistency) to insure a representative sample prior to sampling
8. All the specimens should be handled as if they were capable of transmitting infectious agents.
9. Dispose of all used materials in a container suitable for biohazardous waste.

The sample diluent contains Sodium Azide (0,095 %) as preservative. Avoid contact with the skin

SHELF LIFE AND STORAGE

The expiry date is printed on the pouch of each reaction strip.

Stored the kit between 2 and 30° C. Do not freeze.

SPECIMEN HANDLING

Stool samples must be taken as soon as the symptoms appear. Viral particles decrease in number after 8 days, making the diagnosis more difficult.

The specimen should be received in an airtight transport container and stored at 2°-8°C until tested. The samples can be stored in the refrigerator for 1 to 2 days. For longer storage they must be kept frozen at -20°C. In this case, the sample should be totally thawed, and brought to room temperature and homogenised before testing.

NOTE: Stool in transport media, swabs, or preservatives are inappropriate for testing

PROCEDURE

1. Transfer 1.0 ml of sample diluent in a test tube or vial
2. Add a sample portion of approximately 5-6 mm size (30-50mg), with a swab, a wooden applicator or a bacteriology loop. For liquid or semi-solid stools add 100 microliters of stool using an appropriate pipette. Shake gently in order to resuspend it into the buffer.
3. Vortex for 15 seconds
4. Wait at least 3 minutes until the solid particles settle or centrifuge for 1 minute at 700 xg and transfer with a pipette 500 microliters of supernatant to another test tube.
5. Dip the reaction strip in the second test tube with the arrow pointing to the bottom. **IMPORTANT:** the liquid must not reach the blue area above the arrowheads. If needed, use a larger tube or reduce the amount of sample
6. Incubate the test at room temperature for no less than 5 but no more than 10 minutes. Read the test result in the white area.

Note: the strip can also be introduced into the first tube or vial, provided it is large enough to avoid the liquid to reach the blue area above the arrowheads. However, in rare cases, when using a too concentrated diluted specimen, the absorption might be difficult.

Alternatively, the strip can be dipped for 10 seconds in the first tube, avoiding to exceed the arrowheads, and then transferred to the bench top for the incubation.

INTERPRETATION OF RESULTS

Negative test result: only one GREEN band (Control Line) appears across the white central area of the reaction strip. *Rota* and *Adeno* antigens are absent or below the level of detection.

Positive test result: In addition to the GREEN band (Control Line), a distinguishable RED band (indicating that there are detectable *Rotavirus* antigens in the specimen) and/or a distinguishable BLUE band (indicating that there are detectable *Adenovirus* antigens in the specimen) also appears across the white central zone of the reaction strip.

The intensity of the band will be variable depending on the antigen concentration in the specimen. Any red or blue line, even very weak, must be considered as a positive result. Any line or colour appearing after 10 minutes has no diagnostic value.

Invalid test result: The result window remains white, there are no lines visible or there is one blue line only or one red line only. The specimen should be re-tested using a new reaction device

QUALITY CONTROL

If no Control Line band appears, the test is invalid, since improper test procedure was carried out or deterioration of the reagents has occurred.

LIMITATIONS OF THE PROCEDURE

The test should be used only for the detection of rotavirus or adenovirus antigen in faecal samples.

The test is qualitative and no quantitative interpretation should be made with respect to the intensity of the positive line, when reporting the result.

More than 200 samples were evaluated to assure the correct performance of the test. The correlation of the results with other techniques (ELISA) was excellent. However, interferences in the performance of the tests should not be excluded.

An excess of faecal sample may cause to appear brown lines instead of red and blue lines. These brown lines have no diagnostic value. In this case repeat the test with a smaller amount of sample or dilute the extract already made.

No cross-reactions with other viruses or substances were observed during the evaluation of the test. A negative result does not totally exclude a possible rotavirus or adenovirus infection. The significance of the results must be evaluated in relation to the patient's clinical symptoms.

The analysis of some samples could show lines with indefinite colour, caused mostly by negative samples. Therefore, when lines with indefinite coloration appear, the test should be repeated. If the same result is obtained, it is suggested to use another analysis method

PERFORMANCE CHARACTERISTICS

Sensitivity

After running a big amount of tests with different people reading in several conditions and with diverse lots, we have established a typical sensitivity of 8 ng/ml for rotavirus and 4 ng/ml for adenovirus at 10 min.

The diagnosis utility for this product has been evaluated from stool samples provided to the Microbiology Dpt. (virology section) of the Valladolid Clinico Universitario Hospital (mostly pediatric patient's samples)

Rotavirus: Comparison against the better analysis system, ELISA Ridascreen. Adenovirus: Comparison against a set of the used systems in the analysis; positive samples are those which has been detected at least for 2 of the employed tests.

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REFERENCES / BIBLIOGRAPHIE / BIBLIOGRAFIA / BIBLIOGRAFIA / LITERATURANGABEN

- F. Bon et al. *Prevalence of a group A rotavirus, human calicivirus, astrovirus, and adenovirus type 40 and 41 infections among children with acute gastroenteritis in Dijon, France* Journal of Clinical Microbiology, Sept. 1999, p. 3055-3058
- Bodo R. Eing et al. *Evaluation of two enzyme immunoassays for detection of human rotaviruses in fecal specimens*, Journal of Clinical Microbiology, Dec. 2001, p. 4532-4534
- Umesh D. Parashar et al. *Global illness and deaths caused by rotavirus disease in children*, Emerging Infectious Diseases, vol. 9, No.5, May 2003, p. 565-572

Expiry date / Date de péremption/ Scadenza/ Fecha de caducidad/ Halbtbarkeitsdatum

Lot number / Numéro de lot/ Numero di lotto/ Número de lote / Chargenbezeichnung

For in vitro diagnostic use / Usage in vitro/ Per uso diagnostico in vitro Uso diagnostico in vitro / in vitro Diagnosticum

This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro/ Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro/ Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro/ Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG

Catalogue number / Référence article/Codice dell'articolo/ Número de catálogo / Artikelnummer

Please read pack insert / Lire attentivement le mode d'emploi / Leggere il foglietto informativo/ Leer instrucciones de uso / Bitte Packungsbeilage beachten

Manufactured by / Fabriqué par/ Prodotto da / Fabricado por / Hergestellt von

Contains sufficient for <n> tests / Contenu suffisant pour <n> tests/ Contenu sufficiente per <n> prove/ Contenido suficiente para <n> ensayos / Enthält genügendes für <n> Tests

Conservare al Store at / Conserver à/ Conservar al/ Lagening bei