tru

Strep Pneumo™

A Rapid Immunoassay for the Detection of *Streptococcus pneumoniae* in Human Urine and Cerebrospinal Fluid Specimens

REF 752330

IVD

INTENDED USE

The Tru Strep Pneumo assay is an in vitro, rapid, lateral-flow immunoassay for the qualitative detection of *Streptococcus pneumoniae* in human urine or cerebrospinal fluid (CSF) specimens. The test is performed directly on urine or CSF samples from patients suspected of *Streptococcus pneumoniae* infection.

Test results are to be used as an aid in diagnosis of *Streptococcus pneumoniae* infection. Test results are to be used in conjunction with information obtained from the patient's clinical evaluation and other diagnostic procedures.

The Tru Strep Pneumo assay is intended for use in hospital, reference, or state laboratory settings. The device is not intended for point-of-care use.

SUMMARY AND EXPLANATION OF THE TEST

Streptococcus pneumoniae is the leading cause of community-acquired pneumonia (CAP) worldwide and is responsible for approximately 30-50% of adult CAP cases requiring hospitalization in the US and Europe.¹ *S. pneumoniae* additionally causes acute otitis media and a broad range of invasive pneumococcal disease (IPD), such as bacteremia and meningitis. The eCDC estimates that the incidence of invasive pneumococcal disease within 26 EU/EEA countries is approximately 5.59 per 100,000 people; the majority of these cases in infants < 1 year and adults age 65 and older.² *S. pneumoniae* is the most common cause of bacterial meningitis in adults.³ Pneumococcal meningitis is associated with a 10-30% mortality rate and a higher risk of both death and permanent disability than other bacterial etiologies of meningitis.³

Young children are known asymptomatic carriers of *S. pneumoniae*, including serotypes associated with IPD.⁴ While vaccines have significantly reduced infection and transmission of IPD, current vaccines are only effective against a limited number of known infectious serotypes and emergence of non-vaccine serotypes is continually increasing.² In addition, empiric treatment and antibiotic exposure has contributed to the emergence of multi-drug resistant serotypes that continue to challenge treatment of *S. pneumoniae* infection.² IPD and pneumonia associated with *S. pneumoniae* cannot be directly identified from other microorganisms and therefore clinicians rely on diagnostic testing to determine the causative agent. The current gold standard in laboratory testing is blood or sputum culture. While culture is highly specific, it is also very insensitive with sensitivity reported to be less than 30% and 57% for blood and sputum cultures, respectively.⁵ Culture also typically requires at least 24 hours, which could delay appropriate antimicrobial therapy that is critical in decreasing the high morbidity associated with invasive pneumococcal infections.⁵

Urine and CSF antigen tests provide rapid test results with high performance to direct appropriate antibiotic therapy.

BIOLOGICAL PRINCIPLES

Tru Strep Pneumo is a single use capture immunoassay to detect *Streptococcus pneumoniae* in human urine and CSF specimens. The test consists of a Conjugate Tube, a Test Strip, and Sample Diluent/Negative Control. The Conjugate Tube contains a lyophilized bead of colloidal gold-linked monoclonal antibody to *Streptococcus pneumoniae*. The Test Strip carries a nitrocellulose membrane with dried capture antibodies placed at a designated Test Line for *S. pneumoniae*. The Test Strip holder caps the Conjugate Tube during testing and subsequent disposal to reduce exposure to potential pathogens.

The conjugate bead is first rehydrated in the Conjugate Tube with the Sample Diluent/Negative Control. Patient sample is then added, the contents mixed and the Test Strip added. *Streptococcus pneumoniae* antigens, if present, bind to the antibody-colloidal gold conjugate. When the sample migrates up the Test Strip to the Test Line, the antigen-conjugate complex is bound to the capture antibody, yielding a pink-red line. When no antigen is present, no complexes are formed and no pink-red line appears at the Test Line. An internal control line helps determine whether the test has been executed properly, if kit reagents are performing appropriately and that adequate flow has occurred through the Test Strip during a test run. A visible pink-red line at the Control position of the Test Strip should be present each time a specimen or control is tested. If no pink-red Control Line is seen, the test is considered invalid.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

- Test Strip:** A Test Strip attached to a plastic holder enclosed in a foil pouch with desiccant. The Test Strip carries anti-*Streptococcus pneumoniae* capture antibody for the Test Line. The holder is used to stopper the Conjugate Tube. The strip frame portion of the holder indicates where Test and Control Lines should appear. Store the pouch at 2-27 C when not in use. Do not use a device if the desiccant indicator (line in center of desiccant) has changed from blue to pink.
- Conjugate Tube:** A capped plastic tube containing a conjugate bead. The tube is enclosed in a foil pouch with desiccant. The conjugate consists of gold-conjugated anti-*Streptococcus pneumoniae*, which serves as the detector antibody. Store the foil pouch at 2-27 C when not in use. Do not store in the freezer. Do not remove the cap until just before use.
- Sample Diluent/Negative Control:** A buffered protein solution provided in a plastic vial. Sodium azide (0.094%) and Gentamicin (0.032%) added as preservatives. Use as supplied. Store at 2-27 C when not in use.
- Positive Control:** Inactivated *Streptococcus pneumoniae* in a buffered diluent containing sodium azide (0.094%) as a preservative. The reagent is supplied ready to use. Store at 2-27 C when not in use.
- Plastic transfer pipettes** (with 100, 200 and 300 µL volume marks).
- Tru Strep Pneumo Conjugate Tube labels (to differentiate Tru Strep Pneumo Conjugate Tubes from Conjugate Tubes of other Tru assays).

MATERIALS NOT PROVIDED

- Disposable latex gloves
- Vortex (optional)
- Interval timer
- Marking pen
- Optional: Urine preservative containing boric acid, such as the BD Vacutainer® C&S Preservative Plus Urine Tube

PRECAUTIONS

- All reagents are for in vitro diagnostic use only.
- DO NOT interchange the Test Strips, Conjugate Tubes and Positive Control reagents from different kit lot numbers. The Sample Diluent/Negative Control can be interchanged between kits providing the reagents are within their expiration periods. Do not use any reagent beyond its labeled expiration date.
- Do not use components that lack a label, lot number, or an expiration date.
- Do not use any reagent if it is discolored or turbid. Discoloration or turbidity may be a sign of microbial contamination.
- Allow reagents to warm to 19-27 C before use.
- Dilution errors may affect test performance.
- All reagents should be mixed gently before use.
- The Positive Control reagent vial should be held vertically when dispensing drops to ensure consistent drop size and delivery.
- The transfer pipettes provided with the kit must be used for specimen preparation and transfer. Use one per specimen.
- Any deviation below or above set incubation times may affect sensitivity and specificity and should be avoided.
- Test Strips and Conjugate Tubes are packaged in foil pouches that exclude moisture during storage. Inspect each foil pouch before use. Do not use Test Strips or Conjugate Tubes from pouches that have holes or where the pouch has not been completely sealed. Do not use the Test Strip if the desiccant indicator has changed from blue to pink. The change in the desiccant color is an indicator the Test Device has been exposed to moisture. False-negative reactions may result if Test Strips or Conjugate Tubes are exposed to moisture.
- All samples must be mixed thoroughly before testing to ensure a representative sample prior to pipetting.
- Sample Diluent/Negative Control must be added to the Conjugate Tube within one minute after removing the cap from the tube.

WARNINGS

- Patient specimens may contain infectious agents and should be handled and disposed of as potential biohazards.
- Dispose of all used test materials in an appropriate container. Treat waste as a potential biohazard.
- Positive Control contains inactivated *Streptococcus pneumoniae* but should be handled as if it were potentially infectious.

HAZARDS and PRECAUTIONARY STATEMENTS


There are no known hazards associated with this product.

SHELF LIFE AND STORAGE


The expiration date is indicated on the kit label. Store the kit at 2-27 C.

PROCEDURAL NOTES

The Tru Strep Pneumo transfer pipette is diagrammed below.

The Tru Strep Pneumo Test Strip and Conjugate Tube are diagrammed below.

SPECIMEN COLLECTION AND PREPARATION

Human urine samples, unpreserved: Samples should be received in an airtight transport container and may be stored at 2-8 C prior to testing. Samples should be tested as soon as possible, but may be held up to seven days at 2-8 C. Samples that will not be tested within seven days should be frozen immediately upon receipt and stored at ≤ -20 C for up to 30 days until tested. Samples may be frozen and thawed once after storage at ≤ -20 C prior to testing.

Human urine samples, boric acid preserved: Samples should be received in an airtight Boric Acid Urine Tube and stored at 2-8 C prior to testing. Samples should be tested as soon as possible, but may be held up to seven days at 2-8 C. Samples that will not be tested within seven days should be frozen immediately upon receipt and stored at ≤ -20 C for up to 30 days until tested. Samples may be frozen and thawed once after storage at ≤ -20 C prior to testing.

Human cerebrospinal fluid (CSF) samples: Samples should be tested as soon as possible, but may be held up to seven days at 2-8 C. Samples that will not be tested within seven days should be frozen immediately upon receipt and stored at ≤ -20 C for up to 30 days until tested. Samples may be frozen and thawed once after storage at ≤ -20 C prior to testing. **NOTE: Samples that will be used for culture should be cultured first or split the CSF specimen as the Tru Strep Pneumo test materials are not sterile.**

SPECIMEN PREPARATION

Mix sample thoroughly prior to pipetting. Bring specimens and reagents to room temperature (19-27 C) before testing.

1. Human urine samples, unpreserved and preserved, and CSF samples:

- Remove one Conjugate Tube from its foil pouch and discard the pouch. Label the tube with the patient's name. Apply a Tru Strep Pneumo label to the tube.
- Remove the cap from the Conjugate Tube and discard the cap.
- Using a transfer pipette supplied with the kit, immediately add 100 μ L (first mark from the tip of the pipette) of Sample Diluent/Negative Control to the Conjugate Tube. Dispense directly into the center of the tube. Vortex or swirl the contents of the Conjugate Tube for 10 seconds.*
- Mix patient sample thoroughly. Use one of the transfer pipettes supplied with the kit to mix the sample gently but thoroughly by squeezing the pipette bulb three times. Alternatively, mix for at least 10 seconds using a vortex mixer.
- Using the same pipette, draw 100 μ L of specimen (first mark from the end of the pipette) and add it to the Conjugate Tube.
- Using the same pipette, mix the sample and conjugate thoroughly but gently by squeezing the pipette bulb three times. Alternatively, mix for at least 10 seconds using a vortex mixer. Discard the pipette.

***NOTE: Dilution errors may affect test performance. Failure to add sufficient sample to the Sample Diluent/Negative Control may result in falsely negative tests. Failure to add the full amount of Sample Diluent/Negative Control may result in falsely positive tests. Addition of too much sample may result in invalid test results due to the inhibition of proper sample flow.**


TEST PROCEDURE

- Remove the Test Strip from its foil pouch and discard the pouch. Visually inspect Test Strip. Do not use Test Strip if a green line is present.
- Insert the narrow end of the Test Strip into the Conjugate Tube and firmly press down on the cap to close the tube.
- Incubate at 19-27 C in an upright position in a rack for 20 minutes.
- Read the results on the Test Strip within 1 minute. Do not read results beyond this period. (NOTE: Remove the Test Strip from the Conjugate Tube if Test or Control Lines are difficult to read. Recap the Conjugate Tube with the Test Strip holder and discard when testing is completed.)

EXTERNAL CONTROL TESTS

- Bring all test components, reagents and samples to room temperature (19-27 C) before testing.
- Use 1 Conjugate Tube and 1 Test Strip for Positive Control testing and 1 Conjugate Tube and 1 Test Strip for Sample Diluent/Negative Control testing.
- Remove the Conjugate Tubes from their foil pouches and label accordingly. Discard the pouches.
- Remove the caps from the Conjugate Tubes and discard the caps.
- Add exactly 5 drops of the Positive Control reagent to the Conjugate Tube marked for the Positive Control. The drops should be dispensed directly into the center of the tube.
- Using 1 of the transfer pipettes supplied with the kit, add 200 μ L (second mark from the end of the pipette tip) of Sample Diluent/Negative Control to the Conjugate Tube marked for the Negative Control. Add directly to the center of the tube.
- Vortex or swirl the contents of the tubes for 10 seconds.
- Remove 2 Test Strips from their foil pouches and discard the pouches.
- Insert the narrow end of a Test Strip to each Conjugate Tube and firmly press down on the caps to close each tube.
- Incubate both tubes at 19-27 C for 20 minutes.
- Read the results on the test strip within 1 minute. Do not read results beyond this period. (NOTE: Remove the Test Strip from the Conjugate Tube if Test or Control Lines are difficult to read. Recap the Conjugate Tube with the Test Strip holder and discard when testing is completed.)

INTERPRETATION OF RESULTS

Negative test: A PINK-RED band at the Control Line position. No other bands are present.

Positive test: PINK-RED band of any color intensity at the Control and *Streptococcus pneumoniae* Test Line positions. The color of the Test Line can be lighter than that of the Control Line. Test Lines may appear strongly visible or may appear less strongly visible.

Invalid test results:

- No band at the designated position for the Control Line. The test is invalid since the absence of a Control band indicates the test procedure was performed improperly or that deterioration of reagents has occurred.
- A PINK-RED band appearing at the Test Line position of the device after 21 minutes of incubation, or a band of any color other than PINK-RED. Falsely positive results may occur if tests are incubated too long. Bands with colors other than PINK-RED may indicate reagent deterioration.

If any result is difficult to interpret, the test should be repeated with the same sample to eliminate the potential for error. Obtain a new sample and retest when the original sample repeatedly produces unreadable results.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

At the time of each use, kit components should be visually examined for obvious signs of microbial contamination, freezing or leakage. Do not use contaminated or suspect reagents.

Internal procedural controls: Internal procedural controls are contained within the Test Strip and therefore are evaluated with each test.

- A PINK-RED band appearing at the Control Line serves as a procedural control and indicates the test has been performed correctly, that proper flow occurred and that the test reagents were active at the time of use.
- A clean background around the Control or Test Lines also serves as a procedural control. Control or Test Lines that are obscured by heavy background color may invalidate the test and may be an indication of reagent deterioration, use of an inappropriate sample or improper test performance.

External Control reagents should be tested according to the requirements of the laboratory or applicable local, state or accrediting agencies.

- See section EXTERNAL CONTROL TESTS for instructions on performing these control tests.
- The reactivity of each new lot and each new shipment of Tru Strep Pneumo should be verified on receipt using External Positive and Negative Control reagents. The number of additional tests performed with External Controls will be determined by the requirements of local, state or federal regulations or accrediting agencies.
- The External Controls are used to monitor reagent reactivity. Failure of the controls to produce the expected results can mean that one of the reagents or components is no longer reactive at the time of use, the test was not performed correctly, or that reagents or samples were not added. If the Positive or Negative External Controls fail, do not report test results to the clinician.
- The results expected with the Controls are described in the section on INTERPRETATION OF RESULTS.
- The kit should not be used if control tests do not produce the correct results. **If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If the control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.**
- Positive and Negative Control reagents manufactured for this assay are prepared in the matrix of the Sample Diluent/Negative Control, which may not mimic test specimens. If control materials that are identical in composition to test specimens are preferred, the user can prepare those by diluting known positive and negative specimens in Sample Diluent/Negative Control according to the SPECIMEN PREPARATION section of this insert.

EXPECTED VALUES

Streptococcus pneumoniae is the leading cause of community-acquired pneumonia and bacterial meningitis in adults worldwide. There were approximately 20,843 confirmed cases of invasive pneumococcal disease within 26 European countries in 2011; approximately 5.59 per 100,000 people.²

LIMITATIONS OF THE PROCEDURE

- The test is qualitative and no quantitative interpretation should be made with respect to the intensity of the positive Test Line when reporting the result.
- Tru Strep Pneumo cannot be used as the sole means of determining *Streptococcus pneumoniae* infection. Test results must be used in conjunction with information available from the patient clinical evaluation and other diagnostic procedures.
- Overincubation of tests may lead to an increase in false-positive test results. Conversely, incubation for periods less than those defined in this insert can result in an increase in false-negative tests. Follow incubation times defined in this insert.
- Tru Strep Pneumo detects both viable and non-viable *Streptococcus pneumoniae*. The appearance of Tru Strep Pneumo positive tests depends on the antigen load in the specimen; therefore a Tru Strep Pneumo positive test may not always correlate with the results of bacterial culture.
- A negative test result does not exclude infection with *Streptococcus pneumoniae* nor does it rule out other microbial-caused respiratory infections or invasive disease caused by *S. pneumoniae*. A positive test result does not rule out coinfection with other microbes.
- The Tru Strep Pneumo assay is intended for use with urine or CSF specimens. Use of this assay is not recommended with other specimens.
- A positive result implies the presence of antigen to *Streptococcus pneumoniae*, however, all test results should be reviewed in light of other clinical data by the physician. Young children are known asymptomatic carriers of *S. pneumoniae* that may produce urinary antigen.⁵
- The Tru Strep Pneumo assay is not intended for monitoring therapy.
- The performance of the Tru Strep Pneumo assay has only been established for manual reading and/or visual result determination.
- The Tru Strep Pneumo assay is reactive with *Streptococcus mitis*. *S. mitis* is a normal commensal species of the human mouth that is known to cause rare opportunistic infections, such as infective endocarditis, bacteraemia and septicemia, primarily in immune-compromised patients.
- Tru Strep Pneumo has not been validated for use with either boiled or concentrated urine specimens.

SPECIFIC PERFORMANCE CHARACTERISTICS

Performance characteristics for the Tru Strep Pneumo assay were established in 2014 by the manufacturer using retrospective urine clinical samples from two geographically distinct regions in the United States and contrived positive cerebrospinal fluid (CSF) samples.

A total of one hundred and ten (110) leftover, de-identified, retrospective frozen urine samples were evaluated. Ninety (90) samples were collected from symptomatic patients and frozen and twenty (20) samples were collected prospectively from otherwise normal donors and frozen. There were no exclusions for patient age or gender and demographic information was not collected. There is no expected difference in assay performance due to patient age or gender. Table 1 compares the performance of the Tru Strep Pneumo assay for all urine specimens to the comparator assay.

Table 1. Percent agreement of Tru Strep Pneumo to the comparator assay for urine specimens

Comparator Assay	Tru Strep Pneumo		
	Positive	Negative	Total
Positive	39	1	40
Negative	1	69	70
Total	40	70	110
			95% CI
Positive Agreement	39/40	97.5%	87.1 - 99.6%
Negative Agreement	69/70	98.6%	92.3 - 99.7%
Overall Agreement	108/110	98.2%	93.6 - 99.5%
Tru Invalid Rate	0/110	0.0%	0.0 - 3.4%
Comparator Assay Invalid Rate*	3/110	2.7%	0.9 - 7.7%

*NOTE: Three invalid results were produced by the comparator assay for an initial invalid rate of 2.7%. The invalid samples were repeated and produced acceptable results to disposition samples.

Performance of the Tru Strep Pneumo assay with CSF samples was established using retrospective, frozen CSF samples from patients with symptoms of meningitis or from patients on whom a lumbar puncture was otherwise indicated. There were no exclusions for patient age or gender and demographic information was not collected. There is no expected difference in assay performance due to patient age or gender. Fifty (50) samples were collected and confirmed negative for *Streptococcus pneumoniae* by the comparator method. Therefore, the fifty (50) confirmed negative samples were manufactured as contrived positive specimens using *S. pneumoniae* strain ATCC 6301 spiked at 2 times the limit of detect (LoD), 10X LoD, and 100X LoD. A total of 100 CSF specimens were evaluated with the Tru Strep Pneumo assay in a blind panel. Table 2 compares the performance of the Tru Strep Pneumo assay to the comparator assay for CSF specimens.

Table 2. Percent agreement of Tru Strep Pneumo to comparator assay for CSF specimens

Comparator Assay	Tru Strep Pneumo		
	Positive	Negative	Total
Positive	50	0	50
Negative	0	50	50
Total	50	50	100
			95% CI
Positive Agreement	50/50	100.0%	92.9 – 100.0%
Negative Agreement	50/50	100.0%	92.9 – 100.0%
Overall Agreement	100/100	100.0%	96.3 – 100.0%
Tru Invalid Rate	0/100	0.0%	0.0 - 3.7%
Comparator Invalid Rate	0/100	0.0%	0.0 - 3.7%

REPRODUCIBILITY

Assay precision, intra-assay variability, and inter-assay variability were assessed with a reference panel prepared from negative urine and CSF samples spiked with *S. pneumoniae*. The reproducibility panel consisted of moderately positive (n=3), low positive (n=3), high negative (n=3), and negative (n=1) urine and CSF specimens. Each sample was evaluated by different operators on the same day for five days by three different internal laboratories. Three lots of Tru Strep Pneumo were used in this study. Positive and Negative Controls were tested each day of testing. The results are provided in the table below:

Sample Type	Site 1		Site 2		Site 3		Total	
Urine	Percent Agreement		Percent Agreement		Percent Agreement		Percent Agreement	
Moderate Positive	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
Low Positive	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
High Negative	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
Negative	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
CSF								
Moderate Positive	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
Low Positive	30/30	100.0%	30/30	100.0%	30/30	100.0%	90/90	100.0%
High Negative	30/30	100.0%	30/30	100.0%	28/30	93.3%	88/90	97.8%
Negative	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
Negative Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%
Positive Control	10/10	100.0%	10/10	100.0%	10/10	100.0%	30/30	100.0%

CROSSREACTIVITY STUDIES

Crossreactivity studies were performed with positive and negative urine and CSF specimens inoculated with bacteria or fungi at 10⁵ CFU/mL or viruses at or greater than 10⁵ TCID₅₀/mL.

None of the following organisms in urine reacted with Tru Strep Pneumo:

Acinetobacter baumannii, *Acinetobacter calcoaceticus*, *Acinetobacter lwoffii*, *Bacillus subtilis*, *Corynebacterium diphtheriae*, *Cryptococcus neoformans*, *Escherichia coli* (2 strains), *Enterobacter aerogenes*, *Enterobacter cloacae*, *Enterococcus avium*, *Enterococcus durans*, *Enterococcus faecalis*, *Flavobacterium meningosepticum*, *Gardnerella vaginalis*, *Haemophilus influenzae*, *Haemophilus parainfluenzae*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Lactobacillus acidophilus*, *Lactobacillus brevis*, *Lactobacillus casei*, *Lactobacillus delbruekii*, *Lactobacillus jensenii*, *Legionella pneumophila*, *Micrococcus luteus*, *Moraxella osloensis*, *Morganella morganii*, *Neisseria meningitidis*, *Neisseria polysaccharea*, *Neisseria subflava*, *Nocardia asteroides*, *Proteus mirabilis*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Salmonella* (Groups B, C, D and E), *Serratia marcescens*, *Staphylococcus aureus* (2 strains), *Staphylococcus epidermidis*, *Staphylococcus saprophyticus*, *Streptococcus agalactiae* (Group B), *Streptococcus anginosus* (Groups C, F, G), *Streptococcus bovis* (Group D), *Streptococcus equinus* (Group D), *Streptococcus dysgalactiae ssp. equisimilis* (Group C), *Streptococcus intermedius* (Groups C, F, G), *Streptococcus mutans*, *Streptococcus oralis*, *Streptococcus pyogenes* (Group A), *Streptococcus sanguinis* (Viridans group).

None of the following organisms in CSF reacted with Tru Strep Pneumo:

Acinetobacter baumannii, *Acinetobacter calcoaceticus*, *Acinetobacter lwoffii*, *Bacillus subtilis*, *Corynebacterium diphtheriae*, *Cryptococcus neoformans*, *Escherichia coli* (2 strains), *Enterobacter aerogenes*, *Enterobacter cloacae*, *Enterococcus avium*, *Enterococcus durans*, *Enterococcus faecalis*, *Flavobacterium meningosepticum*, *Haemophilus influenzae*, *Klebsiella pneumoniae*, *Legionella pneumophila*, *Listeria monocytogenes*, *Micrococcus luteus*, *Morganella morganii*, *Neisseria meningitidis* (serotypes A, B, C, W135, Y), *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Staphylococcus saprophyticus*, *Streptococcus agalactiae* (Group B), *Streptococcus anginosus* (Groups C, F, G), *Streptococcus bovis* (Group D), *Streptococcus equinus* (Group D), *Streptococcus dysgalactiae ssp. equisimilis* (Group C), *Streptococcus intermedius* (Groups C, F, G), *Streptococcus mutans*, *Streptococcus oralis*, *Streptococcus pyogenes* (Group A), *Streptococcus sanguinis* (Viridans group), *Coxsackievirus* (type B4), *Herpes simplex virus* (type 1), *Herpes simplex virus* (type 2).

Streptococcus mitis was tested at 1.2 x 10⁸ CFU/mL in both urine and CSF and was found to be reactive with the Tru Strep Pneumo assay.

ASSAY SPECIFICITY

There are 93 known serotypes of *Streptococcus pneumoniae*; as many as 30 different serotypes are commonly associated with disease worldwide. Twenty serotypes associated with invasive pneumococcal disease were tested and produced positive results at the limit of detect for urine (3.4 x 10⁴ CFU/mL) and CSF (6.7 x 10⁴ CFU/mL) with the Tru Strep Pneumo assay. The following serotypes were evaluated: 3, 4, 5, 6A, 6B, 7F, 8, 9V, 10A, 11A, 14, 15B, 16F, 18C, 19A, 19F, 22F, 23A, 23F, 33F.

ANALYTICAL SENSITIVITY

The analytical sensitivity for the Tru Strep Pneumo assay was assessed using *Streptococcus pneumoniae* strain ATCC 6301 (Serotype 1) in two sample matrices, urine and CSF. The analytical sensitivities are listed in the table below.

Strain ID	Sample Type	Limit of Detection (LoD)
<i>S. pneumoniae</i> ATCC 6301	Urine	3.4 x 10 ⁴ CFU/mL
	CSF	6.7 x 10 ⁴ CFU/mL

TESTS FOR INTERFERING SUBSTANCES

The following substances, at the specified saturated solvent/diluent concentrations in urine and CSF, do not interfere with the test results in the final concentrations listed: Acetylsalicylic acid (4.0 mg/mL), Amphotericin B (0.22 mg/mL), Erythromycin (0.067 mg/mL), Prednisone (0.22 mg/mL), Rifampicin (0.09 mg/mL).

The following biological substances, at the specified concentrations in urine and CSF, do not interfere with the test results in the final concentrations listed: Bilirubin (unconjugated, 0.1 mg/mL), Glucose (20 mg/mL), Hemoglobin (64 mg/mL), Protein (30 mg/mL), White blood cells (5.84 x 10⁴ cells/mL).

The following biological substance, at the specified concentration in urine, does not interfere with the test results in the final concentration listed: Urea (20 mg/mL).

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REFERENCES

- World Health Organization. *Weekly Epidemiological Record*. No. 14, 2012, 87, 129-144.
- European Centre for Disease Prevention and Control. 2011. Surveillance of invasive bacterial diseases in Europe: Invasive pneumococcal disease, invasive Haemophilus influenza disease, and invasive meningococcal disease, 2011.
- European Centre for Disease Prevention and Control. Pneumococcal Diseases: Factsheet for health professionals. www.ecdc.europa.eu
- Marchisio, P et al. Nasopharyngeal Carriage of *Streptococcus pneumoniae* in Healthy Children: Implications for the Use of Heptavalent Pneumococcal Conjugate Vaccine. *Emerging Infectious Diseases* 2002, Vol. 8 (5): 479-484.
- Sinclair, A et al. Systematic Review and Meta-Analysis of a Urine-Based Pneumococcal Antigen Test for Diagnosis of Community-Acquired Pneumonia Caused by *Streptococcus pneumoniae*. *JCM* 2013, 51: 2303-2310.
- Navarro D., et al. Performance of Binax NOW Streptococcus pneumoniae Urinary Antigen Assay for Diagnosis of Pneumonia in Children with Underlying Pulmonary Diseases in the Absence of Acute Pneumococcal Infection. *JCM* 2004, 42(10): 4853-4855.

SN11132

REV. 08/23

 Manufactured By	Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com <u>Contacts:</u> Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com
	Meridian Bioscience Europe, SRL Via Dell'Industria 7, 20035 Villa Cortese (Milano) ITALY www.meridianbioscience.com <u>Contacts:</u> Main Telephone (+39) 0331.433636 E-mail: info@meridianbioscience.eu Technical Support: MBE-TechService@meridianbioscience.eu Customer Service/Orders: • For Italian Customers: ordini@meridianbioscience.com • For Distributors / International Customers: Export.CustomerService@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE CHART

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department 800-543-1980.