



ImmunoCard STAT! hCG



A simple, rapid test for detection of hCG



Better patient care starts with STAT!

The ImmunoCard STAT! hCG Combo Test is a rapid immunoassay for the qualitative detection of human chorionic gonadotropin (hCG) in urine or serum as an aid in the early determination of pregnancy.

ALSO AVAILABLE


CONTROLS	URINE	SERUM
Catalog #	755625	755525
# of tests per kit	60	60
Shelf Life	15 months	12 months
Storage Temp	Refrigerated	Refrigerated

- Sensitivity: 10 mIU/ML serum, 20 mIU/ML urine
- Simple, easy to read results
- Rapid results: 3 minutes for urine; 5 minutes for serum
- Room temperature storage




Meridian
Bioscience, Inc.
Inspired Science. Trusted Solutions.™


Test Procedure



1. Squeeze the bulb of the pipette and insert the barrel into the patient sample. Release the bulb and draw up enough sample to fill the barrel to the line indicated on the pipette. Do not overfill.



2. Expel the entire contents of the barrel (135µL) into the sample well of the Test Device. No drop counting required. Read results: 3 minutes for urine or 5 minutes for serum.



POSITIVE

Two separate black or gray bands - one at "T: Test" and the other at "C: Control" - are visible in the results window, indicating that the specimen contains detectable levels of hCG. While the intensity of the test band may vary with different specimens, the appearance of 2 distinct bands should be interpreted as a positive result.

NEGATIVE

If no band appears at "T" and a black or gray line is visible at the "C: Control" position the test can be considered negative, indicating that a detectable level of hCG is not present.

INVALID

If no band appears at the "C: Control" position, the test is invalid. The test is also invalid if incomplete or beaded bands appear at either the stated read time regardless of whether the sample has migrated all the way to the end of the sample window.

**This chart does not contain complete instructions for use. For further information, please refer to the package insert.*

Ordering Information

	Test Size	Meridian#	Cardinal#	Fisher#
<i>ImmunoCard STAT!® hCG Combo</i>	25 Test Cards	755425	B1041-77	23030150

Other Meridian Products

	Test Size	Meridian#	Cardinal#	Fisher#
<i>hCG Combo Control - Urine</i>	60 Tests	755625	B1041-79	
<i>hCG Combo Control - Serum</i>	60 Tests	755525	B1041-78	23030152



USA/Corporate Office:
3471 River Hills Drive
Cincinnati, Ohio 45244
Telephone: 513.271.3700
Orders/Customer Service: 1.800.543.1980
Technical Support: 1.800.343.3858
Information Fax: 513.272.5432
Ordering Fax: 513.271.0124
www.meridianbioscience.com

Belgium / France / Italy / The Netherlands

