



curian®

HpSA™

TEST PROCEDURE

How to perform the test


1



- Prior to sampling, mix stool thoroughly with wooden applicator.
- Remove the top cap from the Aioprep device and set aside.



- Using the Curian HpSA sample collection brush, insert the brush into the stool specimen to coat only the bristles of the brush.
- Insert the Curian HpSA sample collection brush into the Aioprep device.
- Snap off the sample collection brush at the 'break' point and discard handle.
- Re-cap the Aioprep, leaving the sample brush in the device.
- Mix the sample by vortexing for approximately 5 seconds. Do not invert or shake to mix.

2



- After mixing the sample, remove 1 test card from its foil pouch.
- Remove white tip cap from bottom of the Aioprep and discard.
- Holding vertically, with both hands, squeeze near the bottom to add 7 drops into the 'SAMPLE' port of the test card.
- Discard the Aioprep device immediately.
- Start the test read by navigating to the 'TEST' Menu. Follow the on-screen instructions.

3



- **Analyze Now:** Incubate the test card at 19-27 C on the bench top for 20 minutes. Insert test card into the drawer of the Curian Analyzer and read test within 2 minutes of the end of incubation.

- **Incubate and Analyze:** insert the test card immediately into the analyzer and the analyzer will time the incubation for 20 minutes.

- The analyzer will read the test card and automatically report the test result.



This illustration is representative of the current Package Insert at the time of publication. Please refer to the most current version of the Package Insert for complete instructions.
www.meridianbioscience.com/pi

USA/CORPORATE OFFICE

Meridian Bioscience
3471 River Hills Drive
Cincinnati, Ohio 45244 USA
Telephone: 513-271-3700

Orders/Customer Service: 1-800-543-1980
Technical Support: 1-800-343-3858
Information Fax: 513-272-5432
Ordering Fax: 513-271-0124
meridianbioscience.com

meridian BIOSCIENCE®