For the detection of *Helicobacter pylori* antigen in human stool

REF 760130

IVD

Rx Only  
For Professional Use Only

#### INTENDED USE

Curian HpSA, for use with the Curian Analyzer is a rapid, qualitative, fluorescent immunoassay for the detection of *Helicobacter pylori* antigen in human stool. Test results are intended to aid in the diagnosis of *H. pylori* infection and to demonstrate loss of *H. pylori* antigen following treatment. Accepted medical practice recommends that testing by any current method, to confirm eradication, be done at least four weeks following completion of therapy. Test results should be taken into consideration by the physician in conjunction with the patient history and symptoms.

#### SUMMARY AND EXPLANATION OF TEST

Recent studies have shown that *Helicobacter pylori* infections are a major public health concern globally.<sup>1</sup> It is estimated that over half of the global population is infected with *Helicobacter pylori*.<sup>1, 3</sup> Individuals infected with the *H. pylori* organism fall under two groups. The first group, "colonized", make up the majority of the global population and exhibit no gastrointestinal signs and symptoms. The second group, "infected", exhibit gastrointestinal signs and symptoms which includes gastritis, peptic or gastric ulcers, or gastric or duodenal cancer.<sup>5, 9</sup> *H. pylori* infection is considered the leading cause for chronic gastritis and is also considered an etiological agent for gastric cancer (adenocarcinoma) and peptic ulcer disease.<sup>1-3</sup> *H. pylori* infection has also been associated with mucosa-associated lymphoid tissue (MALT) and has been attributed to approximately 89% of gastric cancer diagnoses.<sup>1, 4</sup>

There are multiple methods in which to diagnose an *H. pylori* infection. Biopsy based tests are invasive methods which include culture, polymerase chain reaction (PCR) and the rapid urease test (RUT). The urea breath test (UBT), serology and stool antigen test are non-invasive methods of detection.<sup>1, 6</sup> It is recommended for clinicians to use non-invasive testing methods to aid in the diagnosis of patients who are: (i) without alarming symptoms (e.g., unexplained weight loss, progressive dysphagia, odynophagia), (ii) under the age of 55, and (iii) have a low risk of gastric cancer.<sup>7, 8</sup> Curian HpSA can detect *H. pylori* antigen in human stool and is a rapid, qualitative, non-invasive test designed to aid in the diagnosis of *H. pylori* infection.

#### BIOLOGICAL PRINCIPLES

The Curian HpSA assay consists of a test strip enclosed in a plastic frame (test card), positive control reagent, and Aioprep™ sample preparation device. Curian HpSA is a lateral flow-based immunoassay for the direct detection of *H. pylori* antigens in human stool. Curian HpSA utilizes monoclonal anti-*H. pylori* antibodies as the capture and detector antibodies. A sample of the patient's stool specimen is collected using the sample collection brush that is included in the Curian HpSA assay kit. The brush is pushed through a metering insert, which is integral to the Aioprep device, to remove excess stool. The brush is then pushed directly into the pre-filled sample diluent. The diluted sample is mixed and dispensed drop-wise into the sample port of the Curian HpSA test card. The *H. pylori* antigen (if present) binds to the anti-*H. pylori* antibody conjugated to fluorescent particles, forming a complex. As the sample moves through the test strip, the anti-*H. pylori* capture antibody, bound to the assay membrane at the test position of the strip, binds the complex and yields a test line. When antigen is not present, a complex is not formed and a test line will not occur. As the conjugate continues to move further up the test strip, a detector antibody forms another complex in the control area of the test strip. A line at the control position of the test strip should be present each time a sample or external control is tested. If the control line is not generated, adequate sample flow has not occurred, and the Curian Analyzer will consider the test invalid.

#### REAGENTS/MATERIALS PROVIDED

*The maximum number of tests obtained from this kit is listed on the outer box.*

1. Curian HpSA Test Card: A test strip enclosed in a plastic frame which is in a foil pouch with a desiccant. The test card is supplied ready to use.
2. Curian HpSA Aioprep Sample Preparation Device/ Negative Control: A buffered protein solution containing blue dye and 0.1% ProClin. The Aioprep device is fitted with a metering insert and dropper tip. Supplied ready to use.
3. Curian HpSA Positive Control: Inactivated *H. pylori* purified antigen in a phosphate buffered solution containing 0.094% sodium azide. Supplied ready to use.
4. Curian HpSA Sample Collection Brushes
5. Wooden Applicators

#### MATERIALS NOT PROVIDED

1. Disposable gloves, powder free.

#### EQUIPMENT NOT PROVIDED

1. Vortex mixer
2. Interval timer (Optional)
3. Curian Analyzer, Meridian Bioscience, Inc. Catalog 610190

#### PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Store the kit at the temperature indicated on labeling when not in use.
3. Do not interchange sample collection brushes between assays. Specimens must be sampled and prepared with the Curian HpSA sample collection brush and used with the Curian HpSA Aioprep device. These components are not interchangeable with other assays.
4. Handle and dispose of all human specimens as if they are biologically hazardous.
5. Inspect foil pouch before removing the test card. Do not use test cards that have holes in the foil pouch or where the pouch has not been completely sealed.
6. Do not use test cards where the desiccant indicator has changed from blue to pink.
7. Inspect Curian HpSA Aioprep / Negative Control prior to use. Tap the Aioprep to ensure the liquid is in the main chamber of the Aioprep.
8. Do not interchange positive control or test cards between kit lots.
9. The Curian HpSA test card must be incubated outside of a laminar flow hood, this can be done either inside the Analyzer or on the lab bench.
10. Do not mark over or near the barcode on the test card.

## HAZARD and PRECAUTIONARY STATEMENTS

    Sample Prep	<p><b>Signal word</b> Danger <b>Hazard statements</b> H302 - Harmful if swallowed H315 - Causes skin irritation H317 - May cause an allergic skin reaction H319 - Causes serious eye irritation H411 - Toxic to aquatic life with long lasting effects H311 - Toxic in contact with skin H331 - Toxic if inhaled <b>Precautionary Statements - EU (§28, 1272/2008)</b> P321 - See SDS Section 4 or Section 11 for specific medical treatment information P280 - Wear eye protection/ face protection P264 - Wash face, hands and any exposed skin thoroughly after handling P270 - Do not eat, drink or smoke when using this product P280 - Wear protective gloves/ protective clothing P261 - Avoid breathing dust/ fume/ gas/ mist/ vapors/ spray P271 - Use only outdoors or in a well-ventilated area P272 - Contaminated work clothing should not be allowed out of the workplace P302 + P352 - IF ON SKIN: Wash with plenty of soap and water P312 - Call a POISON CENTER or doctor/ physician if you feel unwell P362 - Take off contaminated clothing and wash before reuse P332 + P313 - If skin irritation occurs: Get medical advice/ attention P363 - Wash contaminated clothing before reuse P304 + P340 - IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing P310 - Immediately call a POISON CENTER or doctor/ physician P301 + P312 - IF SWALLOWED: Call a POISON CENTER or doctor/ physician if you feel unwell P330 - Rinse mouth P405 - Store locked up P403 + P233 - Store in a well-ventilated place. Keep container tightly closed P501 - Dispose of contents/ container to an approved waste disposal plant. P320 - Specific treatment is urgent (see Section 8 on this label)</p>
--	--

## SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-8 C as indicated on the labeling. Return all kit components to the indicated storage temperature after use.

## OPERATION OF THE CURIAN ANALYZER

The Curian Analyzer is an easy to use, menu driven analyzer system. Instructions to complete testing are provided on the analyzer touchscreen and in the Curian Analyzer operator's manual.

1. Power-on the Curian Analyzer by pressing the power button located on the left side of the analyzer. The Curian Analyzer will initialize and perform self-checks.
2. Instrument Checks, Quality Control testing, and patient specimen testing are completed by navigating the Curian Analyzer menus. User ID and Sample ID entries are required for Quality Control and specimen testing.
3. Complete use and operation of the Curian Analyzer by following the instructions in the Curian operator's manual and on-screen prompts. **Refer to operator's manual for detailed instructions for the Curian Analyzer.**

## QUALITY CONTROL

There are four controls for the Curian HpSA test system. Two internal controls, one for the assay and one for the analyzer, and two external controls, one for the assay and one for the analyzer.

1. Analyzer Internal Control: Self Test check (automatic)
2. Analyzer External Control: Instrument Check (IC Check)
3. Assay External Controls: Positive and Negative Controls (QC test)
4. Assay Internal Control: Control line

## INSTRUMENT CHECK TEST PROCEDURE (Analyzer External Control Test)

Note: The Analyzer's External Control (IC Check) test has a default schedule of 30 days. The IC Check test will be required to be performed prior to running a QC test or a patient test if this has not been performed within this timeframe. The user will be locked-out of running a patient test if this test is overdue.

1. Power-on the Curian by pressing the power button located on the side of the analyzer. The Curian will initialize and perform Self Test checks.
2. From the HOME screen, select "TEST", and then select "Instrument Check".
3. Insert the Fluorescent IC Card into the drawer. A pass or fail result will be displayed once the analysis is complete.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. **Refer to operator's manual for detailed instructions for the Analyzer.**

## QUALITY CONTROL TEST PROCEDURE (Assay External Quality Control Test)

Note: Good Laboratory Practice<sup>10</sup> guidelines recommend the use of control material. Users should follow the appropriate federal, state, and local guidelines concerning the running of external quality controls.

1. Bring all test components and reagents to 19-27 C before testing.
2. Use 1 Curian test card for each QC sample control. When ready to perform testing, remove the test card from its foil pouch.
3. Discard the pouch and desiccant. Do not use test cards where the desiccant indicator has changed from blue to pink.
4. Label the test card with control identification information. Do not mark over or near the barcode on the test card

## Negative Control

1. Start the test read by navigating to the "TEST" Menu, selecting the "QC Test". Select "Analyze Now" or "Incubate and Analyze". Follow the on-screen instructions.
2. Remove the white tip cap from the bottom of the Aioprep and discard.
3. While holding the Aioprep vertically, with both hands, squeeze near the bottom of the barrel to add 7 drops of Negative Control reagent (Sample Diluent) into the sample port of the test card.
4. Discard the Aioprep device immediately.
5. **Analyze Now:** Incubate the test card at 19-27 C on the benchtop for 20 minutes. Insert test card into the drawer of the Curian Analyzer and select "Curian HpSA QC Negative". Initiate read within 2 minutes after incubation is complete. Analysis of the reaction must be initiated by the Curian Analyzer. The analyzer will read the test card and automatically report the QC test result.
6. **Incubate and Analyze:** Insert test card immediately into the Curian Analyzer and select "Curian HpSA QC Negative". The Curian will time the incubation for 20 minutes then read the test card and automatically report the QC test result.

### **Positive Control**

1. Start the test read by navigating to the "TEST" Menu and select "QC Test". Select "Analyze Now" or "Incubate and Analyze". Follow the on-screen instructions
2. Invert the positive control bottle to mix.
3. Remove the positive control tip cap and hold vertically to add 4 drops of positive control into the sample port of the test card.
4. **Analyze Now:** Incubate the test card at 19-27 C on the benchtop for 20 minutes. Insert test card into the drawer of the Curian Analyzer and select "Curian HpSA QC Positive". Initiate read within 2 minutes after incubation is complete. Analysis of the reaction must be initiated by the Curian Analyzer. The Curian will read the test card and automatically report the QC test result.
5. **Incubate and Analyze:** If "Incubate and Analyze" is selected, insert the test card immediately into the Curian Analyzer and select "Curian HpSA QC Positive". The Curian Analyzer will time the incubation for 20 minutes then read the test card and automatically report the test result.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

### **SPECIMEN PREPARATION AND TEST PROCEDURE**

This procedure is designed to be used with unpreserved stool. DO NOT USE stool in transport media, on swabs, or mixed with preservatives. The patient specimen should be received in an airtight transport container and stored at 2-8 C until tested. The specimen should be tested as soon as possible but may be held up to 72 hours at 2-8 C prior to testing. If testing cannot be performed within this time frame, specimens should be frozen immediately upon receipt and stored frozen for up to 14 days at -20 C until tested. Specimens may be frozen and thawed twice.

**Note:** Handle all waste and specimens as biohazardous.

**Note:** Mix stool specimen thoroughly with the included wooden applicator (regardless of consistency) before testing.

1. Bring all test components, reagents and specimens to 19-27 C before testing.
2. Prior to sampling, mix stool specimen thoroughly with wooden applicator.
3. Remove the top cap from the Aioprep device and set aside.
4. Using the Curian HpSA sample collection brush, insert the brush into the stool specimen to coat the bristles of the brush. It is important that only the brush is covered with stool. **Do not get stool on the shaft or handle of the brush.** Refer to the examples below.

5. Insert the Curian HpSA sample collection brush into the Aioprep device, pushing the brush through the metering insert.
6. Snap off the sample collection brush at the marked 'Break' point and discard the handle.
7. Recap the Aioprep device, leaving the Curian HpSA sample collection brush in the device.
8. Mix the sample by vortexing for approximately 5 seconds. **Do not invert or shake to mix.**
9. Use 1 Curian HpSA test card for each sample. After mixing the sample, remove 1 test card from its foil pouch.
10. Discard the pouch and desiccant. **Do not use test cards where the desiccant indicator has changed from blue to pink**
11. Label the test card with sample identification information. **Do not mark over or near the barcode on the test card**
12. Remove the white tip cap from the bottom of the Aioprep and discard.
13. While holding the Aioprep vertically, with both hands, squeeze near the bottom of the barrel to add 7 drops into the sample port of the test card.
14. Discard the Aioprep device immediately.
15. Start the test read by navigating to the "TEST" Menu, selecting the appropriate test and following the on-screen instructions.
16. **Analyze Now:** Incubate the test card at 19-27 C on the benchtop for 20 minutes. Insert test card into the drawer of the Curian Analyzer and read test within 2 minutes of the end of incubation. Analysis of the reaction must be initiated by the Curian Analyzer. The analyzer will read the test card and automatically report the test result.
17. **Incubate and Analyze:** If "Incubate and Analyze" is selected, insert the test card immediately into the Curian Analyzer and the analyzer will time the incubation for 20 minutes then read the test card and automatically report the test result.

### **INTERPRETATION OF RESULTS**

Result interpretation is completed automatically by the Curian Analyzer system. The result will be shown on screen. Results can be retrieved from instrument storage, printed and/or exported.

#### **Instrument Check:**

One of the following result interpretations will be generated by the Curian Analyzer for IC Check.

1. **Pass:** This indicates that the essential and critical components of the Curian Analyzer are working correctly.
2. **Fail:** This indicates that there may be an issue with the Curian Analyzer. Refer to the Curian operator's manual if an error message has occurred. Repeat testing to assist in troubleshooting the error. If repeat testing results in a failed output, please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

#### **External Quality Controls:**

One of the following result interpretations will be generated by the Curian Analyzer for External Controls.

1. **Pass:** This indicates the test card and reagents are performing as intended.
2. **Fail:** This result indicates the test card and/ or reagents are not performing as intended or a user error has occurred. Test should be repeated to assist in trouble-shooting the error. If repeat testing results in a failed output, please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.
3. **Control Invalid/ Fail:** This result indicates inadequate sample flow or incorrect sample preparation. Repeat the test. If repeat testing yields the same results, please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor for additional assistance.

#### **Patient Test Results:**

One of the following result interpretations will be generated by the Curian Analyzer for patient specimens.

1. **Positive:** *H. pylori* stool antigen is detected in sample.
2. **Negative:** *H. pylori* stool antigen is not detected in sample.
3. **Invalid:** Result indicates inadequate flow of sample or incorrect sample preparation. Repeat the test from the original stool specimen. If repeat testing yields the same results, please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor for additional assistance.

## EXPECTED VALUES

Recent studies have shown that *Helicobacter pylori* infections are a major public health concern globally.<sup>1</sup> It is estimated that over half of the global population is infected with *Helicobacter pylori*.<sup>1, 3</sup> In patients diagnosed with duodenal ulcers, however, it has been shown to be 90%-95%.<sup>12</sup> Currently, recommended eradication treatments have an overall eradication rate of 84.3% (82.1% to 86.4%).<sup>11</sup>

The Curian HpSA assay detects the presence of *H. pylori* antigens in human stool. Expected values for a given population should be determined for each laboratory. The rate of positivity may vary depending on geographic location, method of specimen collection, handling and transportation, test employed and general health environment of patient population under study. The observed prevalence of *H. pylori* in this study was 13.5% (73/542). The prevalence of *H. pylori* by gender is provided below.

Gender	Samples Positive by Comparator Method and Curian HpSA	Prevalence in Tested Population
Female	42	7.7%
Male	31	5.7%
Total	73	13.5%

## LIMITATIONS OF THE PROCEDURE:

1. The Curian HpSA assay must only be used with the Curian Analyzer.
2. Curian HpSA is a qualitative, in vitro diagnostic test. The Curian Analyzer will only provide qualitative results. This test is not intended to provide quantitative values.
3. The Curian HpSA assay is designed to be used with unpreserved stool samples. DO NOT USE stool in transport media, on swabs, or mixed with preservatives.
4. Antimicrobials, proton pump inhibitors and bismuth preparations are known to suppress *H. pylori* and ingestion of these prior to *H. pylori* testing (culture, histology, rapid urease, UBT, antigen) may give a false negative result. If a negative result is obtained for a patient ingesting these compounds within two weeks prior to performing the Curian HpSA assay, a false-negative result may be obtained. The test should be repeated on a new specimen two weeks after discontinuing treatment. A positive result for a patient ingesting these compounds within two weeks prior to performing the Curian HpSA test, should be considered accurate.
5. Performance characteristics have not been established in asymptomatic populations.
6. A negative test result does not preclude the possibility of the presence of *H. pylori* antigen in the specimen which may occur if the level of antigen is below the detection limit of the test.
7. No data exists on the effects of colonic washes, barium enemas, laxatives, or bowel preparations on the performance of the Curian HpSA assay. These procedures can result in extensive dilution or the presence of additives that may affect test performance.

## SPECIFIC PERFORMANCE CHARACTERISTICS

### CLINICAL PERFORMANCE

The Curian HpSA assay was evaluated at three study sites representing geographically distinct regions throughout the United States. Five hundred forty-two specimens from patients suspected of having an *H. pylori* infection were prospectively collected and enrolled into the study.

All specimens were tested at the study sites with the Curian HpSA assay and an FDA-cleared *H. pylori* stool antigen EIA that was previously evaluated relative to the endoscopy biopsy composite reference method (i.e., culture, histology, and RUT) for initial *H. pylori* diagnosis with a demonstrated sensitivity and specificity greater than or equal to 95% and a lower bound of the two-sided 95% confidence interval (CI) greater than 89%. Specimens with discordant results between the Curian HpSA assay and the FDA-cleared *H. pylori* stool antigen EIA assay were further tested in-house using a validated PCR assay.

The overall clinical performance (positive percent agreement (PPA) and negative percent agreement (NPA)) of the Curian HpSA assay compared to the FDA-cleared *H. pylori* stool antigen EIA is presented in the table below. PCR results of discordant testing are footnoted below.

Positive and Negative Curian HpSA Results vs. FDA-cleared <i>H. pylori</i> stool antigen EIA				
		FDA-cleared <i>H. pylori</i> stool antigen EIA		
		Positive	Negative	Total
Curian HpSA Assay	Positive	73	14 <sup>b</sup>	87
	Negative	3 <sup>a</sup>	452	455
	Total	76	466	542
Agreement			95% CI (%)	
PPA	96.1%	(73/76)	89.0%, 98.6%	
NPA	97.0%	(452/466)	95.0%, 98.2%	

<sup>a</sup> 2/3 Curian HpSA false negatives were dispositioned as negative by PCR

<sup>b</sup> 8/14 Curian HpSA false positives were dispositioned as positive by PCR

### Invalid Results by Curian HpSA assay

Three samples produced invalid results during initial testing with the Curian HpSA assay. The rate of initial invalid results was 0.5% (3/542). After re-testing, the 3 samples had valid negative results. The rate of final invalid results was 0.0%.

### ANALYTICAL SENSITIVITY

The Limit of Detection (LoD) of the Curian HpSA assay was determined for *Helicobacter pylori* stool antigens in a negative diluted natural stool matrix (70% natural stool/ 30% physiological saline). Diluted natural stool was used because of difficulties preparing dilutions with neat stool for analytical testing. Diluted stool was used for all analytical studies. The LoD is defined as the lowest concentration of measurand which produced positive results ≥ 95% of the time. The Curian HpSA assay LoD concentration is 2.0 ng/mL.

### ASSAY REACTIVITY/ INCLUSIVITY

A total of 5 strains of *H. pylori* were evaluated for reactivity with the Curian HpSA assay. The final reactive concentrations observed for each strain are shown in the table below.

<i>H. pylori</i> strain tested	Geographic origin & other information	Reactive Concentration
ATCC 43504	Australia	1.0 x 10 <sup>5</sup> CFU/mL
CCUG 38771	Unknown	3.0 x 10 <sup>5</sup> CFU/mL
CCUG 19087	South Africa	3.0 x 10 <sup>5</sup> CFU/mL
ATCC 700392	UK; clade hpEurope	6.0 x 10 <sup>5</sup> CFU/mL
ATCC 700824	US; clade hpAfrica1	3.0 x 10 <sup>5</sup> CFU/mL

**REPRODUCIBILITY**

Ten blinded-coded panels consisting of 16 samples per panel were supplied to three study sites for this study. Each panel consisted of 5 moderately positive samples, 5 low positive samples, 5 high negative samples, and 1 true negative sample. Panel members were generated by spiking *H. pylori* purified flagellar antigen into diluted, natural stool matrix.

Each site used 2 unique operators and tested one panel of samples each day for five days. Positive and negative external controls were tested with each panel. Three Curian HpSA kit lots and 3 Curian Analyzers were used in this study. Each site tested 2 Curian HpSA kit lots using one analyzer. The results of the study are provided in the table below.

	Overall Agreement							
	Site 1 Percent Agreement		Site 2 Percent Agreement		Site 3 Percent Agreement		Total Percent Agreement	
Sample Type	Rate	%	Rate	%	Rate	%	Rate	%
<b>High Negative</b>	44/50	88.0	43/50	86.0	46/50	92.0	133/150	88.7
<b>Low Positive</b>	50/50	100.0	47/50	94.0	50/50	100.0	147/150	98.0
<b>Moderate Positive</b>	49/50	98.0	50/50	100.0	50/50	100.0	149/150	99.3
<b>Negative Control</b>	10/10	100.0	10/10	100.0	10/10	100.0	30/30	100.0
<b>Positive Control</b>	10/10	100.0	10/10	100.0	10/10	100.0	30/30	100.0
<b>True Negative</b>	10/10	100.0	9/10	90.0	10/10	100.0	29/30	96.7

**CROSSREACTIVITY/ MICROBIAL INTERFERENCE**

The specificity of Curian HpSA was tested utilizing the following bacterial, fungal and viral strains. Each potentially cross-reactive microorganism was added at minimum concentrations of  $1.0 \times 10^7$  CFU/mL (bacteria/fungi) or  $1.0 \times 10^5$  TCID<sub>50</sub>/mL (for viruses) to a diluted, natural negative stool matrix and a contrived positive matrix sample. No crossreactivity or microbial interference with the Curian HpSA assay was observed.

Organism Name	Strain ID	Organism Name	Strain ID
Adenovirus 40	Dugan	<i>Klebsiella pneumoniae</i>	ATCC 13883
<i>Aeromonas hydrophila</i>	ATCC 35654	<i>Proteus vulgaris</i>	CCUG 6380
<i>Bacillus cereus</i>	CCLUG 52704	<i>Pseudomonas aeruginosa</i>	ATCC 39324
<i>Borrelia burgdorferi</i>	B31.5A19	<i>Rotavirus</i>	WA
<i>Campylobacter coli</i>	ATCC 10956	<i>Salmonella spp. Dublin</i>	ATCC 15480
<i>Campylobacter jejuni</i>	ATCC 29411	<i>Salmonella spp. Hilversum</i>	ATCC 15784
<i>Candida albicans</i>	ATCC 18804	<i>Salmonella spp. Minnesota</i>	ATCC 9700
<i>Citrobacter freundii</i>	ATCC 8090	<i>Salmonella typhimurium Group B</i>	ATCC 14028
<i>Clostridium difficile</i>	ATCC 43255	<i>Shigella bovis</i>	ATCC 9207
<i>Clostridium perfringens</i>	ATCC 12915	<i>Shigella dysenteriae</i>	ATCC 9361
<i>Enterobacter cloacae</i>	ATCC 15337	<i>Shigella flexneri</i>	ATCC 12022
<i>Enterococcus faecalis</i>	ATCC 49532	<i>Shigella sonnei</i>	ATCC 25931
<i>E. coli</i> O157:H7	ATCC 43895	<i>Staphylococcus aureus</i>	ATCC 6538
<i>E. coli</i>	ATCC 9637	<i>Staphylococcus aureus</i> Cowan I	ATCC 12598
<i>Escherichia fergusonii</i>	ATCC 35469	<i>Staphylococcus epidermidis</i>	ATCC 51625
<i>Haemophilus influenzae</i>	ATCC 9006	<i>Yersinia enterocolitica</i>	ATCC 23715

**TESTS FOR INTERFERING SUBSTANCES**

Interference testing was performed in the presence of chemical and biological substances introduced directly into contrived HpSA low positive and negative samples. No interference was observed with the Curian HpSA assay for any of the substances tested. The substances tested and concentrations evaluated are shown in the table below.

Substance (active ingredient(s))	Test Concentration
Barium Sulfate	5% w/v (50 mg/mL)
Benzalkonium chloride	1% v/v
Ciprofloxacin	0.25% w/v (2.5 mg/mL)
Ethanol	1% v/v
Hog gastric mucin	3.5% w/v (35 mg/mL)
Human blood (whole)	40% v/v
Human hemoglobin	12.5% w/v (125 mg/ml)
Human urine	5% v/v
Hydrocortisone	1% w/v (10 mg/mL)
Imodium® (Loperamide HCl, 1 mg/7.5 mL)	5% v/v
Kaopectate® (Bismuth subsalicylate 262 mg/15 mL)	5% v/v
Leukocytes	0.05% v/v
Mesalazine (5-Aminosalicylic acid)	10% w/v (100 mg/mL)
Metronidazole	0.25% w/v (2.5 mg/mL)
MiralAX® (Polivethylene Glycol 3350, 17 g/dose)	7% w/v (70 mg/mL)
Mineral Oil	10% v/v
Mylanta® (per 10 mL: Aluminum hydroxide 800 mg, Magnesium hydroxide 800 mg)	4.2 mg/mL (2.5% v/v)
Naproxen Sodium	5% w/v (50 mg/mL)
Nonoxynol-9	1% v/v
Nystatin	1% w/v (10 mg/mL)
Palmitic acid (fecal fat)	20% w/v (200 mg/mL)
Pepto-Bismol® (Bismuth subsalicylate 525 mg/30 mL)	5% v/v
Phenylephrine	1% w/v (10 mg/mL)
Prilosec OTC® (Omeprazole 20 mg/tablet)	5 mg/mL
Sennosides	1% w/v (10 mg/mL)
Simethicone	10% v/v
Stearic acid (fecal fat)	20% w/v (200 mg/mL)
Tagamet HB 200® (Cimetidine 200 mg/tablet)	5 mg/mL
TUMS®	5 mg/mL
Vancomycin	0.25% w/v (2.5 mg/mL)

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## REFERENCES

1. Hooi James KY, Lai Wan Ying, Ng Wee Khoon, Suen Michael MY. Global Prevalence of *Helicobacter pylori* Infection: Systematic Review and Meta-Analysis, Gastroenterology. 2017; Vol 153, No. 2.
2. Malfertheiner P, Megraud F, O'Morain C, et al. Current concepts in the management of *Helicobacter pylori* infection: the Maastricht III Consensus Report. Gut 2007;56:772–781.
3. Wang C, Yuan Y, Hunt RH. The association between *Helicobacter pylori* infection and early gastric cancer: a meta-analysis. Am J Gastroenterol. 2007;102:1789–1798.
4. International Agency for Research on Cancer *Helicobacter pylori* Working Group. *Helicobacter pylori* Eradication as a Strategy for Preventing Gastric Cancer. (IARC Working Group Reports, No. 8). Lyon, France: International Agency for Research on Cancer, 2014.
5. Guadalupe Ayala, Escobedo-Hinojosa, Wendy Itzel, de la Cruz-Herrera Carlos Felipe, Romero Irma. Exploring alternative treatments for *Helicobacter pylori* infection, World J Gastroenterology. 2014;February 14; 20(6): 1450-1469
6. Garza-González Elvira, Perez-Perez Guillermo Ignacio, Maldonado-Garza Héctor Jesús, Bosques-Padilla Francisco Javier. A review of *Helicobacter pylori* diagnosis, treatment, and methods to detect eradication, World J Gastroenterology. 2014;February 14; 20(6): 1438-1449
7. Moayyedi Paul M, MB, ChB, PhD, MPH, FACG 1 , Lacy Brian E, MD, PhD, FACG 2, Andrews Christopher N, MD 3, Enns Robert A, MD 4, Howden Colin W, MD, FACG 5 and Vakil Nimish, MD, FACG, ACG and CAG Clinical Guideline: Management of Dyspepsia, AM J Gastroenterol advance online publication., 2017;20 June doi 10.1038/ajg.2017.154.
8. Fashner Julia, MD, Gitu Alfred C. MD, Diagnosis and Treatment of Peptic Ulcer Disease and *H. pylori* Infection, Am Fam Physician. 2015;91(4):236-242.
9. Logan Robert PH, Walker Marjorie M. Epidemiology and diagnosis of *Helicobacter pylori* infection, BMJ. 2001 Oct 20; 323(7318): 920–922.
10. US Department of Health and Human Services PHS/CDC/NIH. Biosafety in microbiology and biomedical laboratories, Washington DC: US Government Printing Office. 2007.
11. Gatta, Luigi, Vakil, Nimish, Vaira, Dino, et al. Global eradication rates for *Helicobacter pylori* infection: systematic review and meta-analysis of sequential therapy BMJ 2013;347:f4587 doi: 10.1136/bmj.f4587, 7 August 2013.
12. Testerman, Traci L., Morris, James, Beyond the stomach: An updated view of *Helicobacter pylori* pathogenesis, diagnosis, and treatment, World J Gastroenterol. 2014 Sep 28; 20(36): 12781–12808.

REV. 04/24

 <b>Manufactured By</b>	<p><b>Meridian Bioscience, Inc.</b>  3471 River Hills Drive  Cincinnati, OHIO - 45244 USA  <a href="http://www.meridianbioscience.com">www.meridianbioscience.com</a></p> <p><b>Contacts:</b>  Main Telephone (+1) 513.271.3700  Customer Service/Orders 800.543.1980  Technical Support Center 800.343.3858  Information Fax: 513.272.5432  Ordering Fax: 513.271.0124  E-mail: <a href="mailto:info@meridianbioscience.com">info@meridianbioscience.com</a></p>
	<p><b>Meridian Bioscience Europe, SRL</b>  Via Dell'Industria 7, 20035 Villa Cortese  (Milano) ITALY  <a href="http://www.meridianbioscience.com">www.meridianbioscience.com</a></p> <p><b>Contacts:</b>  Main Telephone (+39) 0331.433636  E-mail: <a href="mailto:info@meridianbioscience.eu">info@meridianbioscience.eu</a>  Technical Support: <a href="mailto:MBE-TechService@meridianbioscience.eu">MBE-TechService@meridianbioscience.eu</a>  Customer Service/Orders:  • For Italian Customers:  <a href="mailto:ordini@meridianbioscience.com">ordini@meridianbioscience.com</a>  • For Distributors / International Customers:  <a href="mailto:Export.CustomerService@meridianbioscience.eu">Export.CustomerService@meridianbioscience.eu</a></p>
<b>AUSTRALIAN SPONSOR</b>	<p><b>Emergo Australia</b>  Level 20, Tower II  Darling Park  201 Sussex Street  Sydney, NSW 2000  Australia</p>
	<p><b>MedEnvoy Switzerland</b>  Gotthardstrasse 28  6302 Zug  Switzerland</p>

## INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

**Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de símbolos, Zeichenerklärung, Sleutelgids voor symbolen, Guia-chave para símbolos)**

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis / Uiterste gebruiksdatum / Prazo de validade		Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle / Positieve controle / Controlo positivo
	Batch Code / Codice del lotto / Code du lot / Código de lote / chargebenzelnung / Batchcode / Código do lote		Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle / Negatieve controle / Controlo negativo
	In vitro diagnostic medical device / Dispositivo medico-diagnóstico In vitro / Dispositif médical de diagnostic in vitro . Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum / Medisch hulpmiddel voor in-vitrodiagnosek / Dispositivo médico de diagnóstico in vitro		Authorized representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandatario dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft / Gevolmachtigd vertegenwoordiger in de Europese Gemeenschap / Representante autorizado na Comunidade Europeia
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sul dispositivo medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostika 98/79/EG / CE-markering / Dit product voldoet aan de eisen van Richtlijn 98/79/EG met betrekking tot medische hulpmiddelen voor in-vitrodiagnosek / Este producto cumple os requisitos da Diretiva 98/79/CE relativa a dispositivos médicos de diagnóstico in vitro		Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungsbuffer befindet / Monsterpreparatieapparaat met monsterooplosmiddel / Aparelho de preparação de amostras com diluente de amostras
	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer / catalogusnummer / Referência		CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risques de danger / PRECAUCIÓN: Peligroso / WARNUNG: Risikogefahr / LET OP: Risico van gevaar / ATENÇÃO: potencial perigo
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisungsbeachten / Raadpleeg de gebruiksaanwijzing / Consultar as Instruções de Utilização		Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer / Reactiebuffer / Tampão de reação
	Manufacturer / Fabbricante / Fabricant / Fabricante / Hersteller / Fabrikant / Fabricante		For IVD Performance Evaluation Only / Solamente para evaluación de las prestaciones / Réactifs IVD reservados à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung / Uitsluitend voor prestatieverlating IVD / Apenas para avaliação do desempenho IVD
	Contains sufficient for <n> tests / Contiene suficiente per "n" saggi / Contenu suffisant pour "n" test / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen / Bevat voldoende voor <n> testen / Conteúdo suficiente para <n> testes		Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung / Stopvloeistof / Solução de paragem
	Temperature limitation / Limiti di temperatura / Limites de température / Límite de temperatura / Temperaturbegrenzung / Temperatuurgrens / Limites de temperatura		Enzyme Conjugate / Conjugato enzimatico / Conjugué enzymatique / Conjulado enzimático / enzymkonjugat / Enzymconjugaat / Conjulado Enzimático
	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer / Serienummer / Número de série		Assay Control / Control del test / Test de contrôle / Control de Ensayo / Kontrolltest / Assaycontrol / Controlo de Ensaio
	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / testgarät / Testinstrument / Dispositivo de teste		Reagent / Reagente / Réactifs / Reactivos / Reagenzien / Reagens / Reagente
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum / Datum van vervaardiging / Data de fábrica		Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer / Concentratie washbuffer: 20X / Tampão de lavagem
	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer / Buffer / Tampão		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise / Waarschuwing / Aviso
	Conjugate / Conjugato / Conjugué / Conjulado / Konjugat / Conjugat / Conjulado		Specimen Diluent (or Sample Diluent) / Diluente del Campione / diluant échantillons / Diluyente de muestra / Probenverdünnungsbuffer / Specimenverdünnungsmittel (of monsterverdünnungsmittel) / Diluente de amostra
	Substrate / Substrato / Substrat / Substrato / Substrat / Substrato		Wash Buffer Concentration 20X / Soluzione dil lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat / Bufferconcentrate 20x wassen / Concentração do tampão de lavagem: 20X
	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso por Receta / verschreiberpflichtig / Gebrauk uitsluitend op voorschrijf / Apenas para utilização mediante prescrição		Detection Reagent / Reagente Diretto / Réactif de Detection / Reactivo de Deteción / Nachweis Reagenz / Detectiereagens / Reagente de deteção
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / ne pas utiliser si le paquet est endommagé / No use si el paquete esta dañado / Nicht verwenden, wenn die Verpackung beschädigt ist / Niet gebruiken als de verpakking beschadigd is / Não utilizar se a embalagem estiver danificada		Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß / Leeg huisje / Tubo vazio
	Single Use Only / Prodotto Monouso / A usage unique / Para Un Solo Uso / nur für die einmalige Anwendung / Uitsluitend voor eenmalig gebruik / Apenas para uma única utilização		Swiss Authorized Representative / Mandatario svizzero / Mandataire Suisse / Representante Autorizado Suizo / Schweizer Bevollmächtigter / Gevolmachtigd vertegenwoordiger in Zwitserland / Representante autorizado na Suíça

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.