MONOSPOT[®] Latex

REF

776150

IVD**Rx Only****INTENDED USE**

MONOSPOT Latex is a one step rapid latex particle agglutination test for the qualitative and semiquantitative determination of infectious mononucleosis heterophile antibodies in serum or plasma. MONOSPOT Latex aids in the diagnosis of infectious mononucleosis.

SUMMARY AND EXPLANATION OF THE TEST

Infectious mononucleosis is an acute infectious disease of viral etiology. The most frequent symptoms are fever, sore throat, tender lymphadenopathy, anorexia, malaise, headache and myalgia. Splenomegaly occurs in most patients. A macular, maculopapular or petechial rash occurs in up to 50% of the cases, but such rashes occur most commonly in patients who have been treated with ampicillin.

The complications of infectious mononucleosis include secondary bacterial pharyngitis, rupture of the spleen, autoimmune hemolytic anemia, autoimmune thrombocytopenia, myocarditis, hepatitis and central nervous system involvement with meningoencephalitis or transverse myelitis. Fatal fulminant infectious mononucleosis or acquired hypogammaglobulinemia is rarely seen.

The diagnosis made on clinical history and symptomatology alone is difficult. Numerous cases in which infectious mononucleosis has been misidentified with other non-related viral and bacterial diseases have been cited.¹ For this reason, hematologic and serologic tests are very helpful in diagnosis. In 1932, Paul and Bunnell² noted that sera from patients with infectious mononucleosis have heterophile antibodies to sheep erythrocytes. Also described were agglutinins to red blood cells from other mammals.^{3,4}

The proteins responsible for this agglutination are glycoproteins from red cell membranes called Paul-Bunnell antigen by several authors. Studies made on these glycoproteins show that those purified from bovine red blood cells are the most sensitive to infectious mononucleosis heterophile antibodies. Heterophile antibodies to sheep erythrocytes (which are different from those present during infectious mononucleosis), may also be detected in sera from normal people, from individuals who have received injections of serum, and others.^{3,5}

Traditionally the infectious mononucleosis heterophile antibodies have been distinguished from other heterophile antibodies by a "differential" absorption test^{6,7} with bovine red blood cells and guinea pig kidney tissue. Now, the use of the purified Paul-Bunnell antigen attached to latex particles provides a simple method with improved sensitivity for the specific detection of heterophile antibodies associated with infectious mononucleosis.

In 1968 the etiologic agent of infectious mononucleosis was described.⁸ It was called the Epstein-Barr virus (EBV), a member of the herpes virus group. Subsequently, several serologic techniques involving EBV-related antigens have been developed.

The mode of transmission of infectious mononucleosis appears to be intimate salivary contact, salivary contamination of eating and drinking vessels and airborne dissemination of EBV.⁹

BIOLOGICAL PRINCIPLES

The MONOSPOT Latex reagent is a suspension of polystyrene latex particles of uniform size coated with highly purified Paul-Bunnell antigen from bovine red cell membranes. The degree of purity of the antigen is such that MONOSPOT Latex only reacts with infectious mononucleosis heterophile antibodies. For this reason, "differential" absorptions are not necessary. Latex particles allow visual observation of the antigen-antibody reaction. If infectious mononucleosis heterophile antibodies are present in either serum or plasma, the latex suspension changes its uniform appearance and a clear agglutination becomes evident.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. **Latex Reagent** – Suspension of polystyrene latex particles coated with Paul-Bunnell antigen in a buffer. Contains < 0.1% sodium azide.
2. **Positive Control** – Rabbit IgG anti-Paul-Bunnell antigen diluted in a buffer. Contains < 0.1% sodium azide.
3. **Negative Control** – Nonreactive diluted human serum. Contains < 0.1% sodium azide.
4. **Pipettestirs**
5. **Test Slides**

MATERIALS REQUIRED BUT NOT PROVIDED

1. Normal saline (0.9% NaCl, only for semiquantitative technique)
2. Automatic pipettes
3. Timer
4. Rotator
5. Stirrers

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. To assure proper delivery the reagent dropper must be held vertically and a single drop allowed to fall.
3. The reagents in this kit contain sodium azide as a preservative. Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion. Flush drains with water thoroughly after disposing of fluids containing sodium azide.

WARNING: POTENTIALLY BIOHAZARDOUS MATERIAL.

Each human source material used in the preparation of this material was tested by an FDA approved method for the presence of the antibody to HIV I and II and HCV antibodies as well as for HBsAg and found to be negative.

Because no test method can offer complete assurance of the absence of infectious agents, the reagents should be handled carefully.¹⁰ Dispose all used materials in a suitable biohazardous waste container.

HAZARD and PRECAUTIONARY STATEMENTS

 MONOSPOT Detection Latex MONOSPOT Negative Control MONOSPOT Positive Control	Signal Word Danger Precautionary Statements - EU (§28, 1272/2008) P301 + P310 - IF SWALLOWED: Immediately call a POISON CENTER or doctor/ physician
---	---

SHELF LIFE AND STORAGE

1. The Latex Reagent and the controls will remain stable through the expiration date shown on the label if stored between 2 and 8 C. Do not freeze. The reagents can be damaged by improper handling, especially temperature extremes.
2. Checking the Latex Reagent with the positive and negative controls provided will permit detection of reagent deterioration. The reagents should not be used after the expiration date shown on the label.
3. The Latex Reagent, once shaken, must be uniform without visible clumping.
4. When stored, a slight sedimentation may occur and should be considered normal.
5. Do not use the Latex Reagent or controls if they become contaminated.
6. The reagent dropper dispenses drops of 28 µL ± 10%. The dropper must be held perpendicular to the slide surface and a single drop allowed to fall. Do not use another dropper without previously checking the volume of the drop.

SPECIMEN COLLECTION AND PREPARATION**SERUM:**

Use fresh serum collected by centrifuging clotted blood. If the test cannot be carried out on the same day, serum may be stored between 2 and 8 C for no longer than 48 hours after collection. For longer periods the sample must be frozen (-20 C or below).

PLASMA:

Collect blood into a tube containing anticoagulant (EDTA). Other anticoagulants should be evaluated before use. Centrifuge to separate plasma from cellular elements. Test the specimen within 24 hours of blood collection. Do not use hemolyzed or contaminated samples.

TEST PROCEDURE**QUALITATIVE TECHNIQUE**

1. Allow the Latex Reagent and controls to reach room temperature (20 to 30 C).
2. Gently shake the Latex Reagent vial to disperse and suspend the latex particles in the buffer solution. Vigorous shaking should be avoided.
3. Place 50 µL of the sample (or one drop of control) on one section of the disposable slide.
4. Shake the reagent vial and add one drop of Latex Reagent next to the drop of sample.
5. Mix both drops with a stirrer covering the whole surface of the slide section.
6. Gently rotate the slide for 3 minutes manually or on a rotary shaker set at 80-100 rpm.
7. Look for the presence or absence of agglutination after the aforementioned period of time.

SEMIQUANTITATIVE TECHNIQUE

Allow the Latex Reagent to reach room temperature (20 to 30 °C). Prepare the sample dilutions on 2 slides (see descriptive diagram):

1. Place 50 µL of normal saline on slide sections 2 through 6.
2. Using an automatic pipette, place 50 µL of the sample on slide sections 1 and 2.
3. Using the same pipette, take in and release the sample and the normal saline in section 2 several times until they are well mixed.
4. Take 50 µL of the mixture made on section 2 and transfer it to section 3.
5. Repeat the aforementioned operations to obtain a thorough mixing of reagents, transferring 50 µL from section 4 and so on, in succession, through section 6, thereafter discarding 50 µL.

SECTION	1	2	3	4	5	6
SALINE µL	—	50	50	50	50	50
SAMPLE µL	50	50	50	50	50	50
MIX and TRANSFER		↓ 50	↑ 50	↓ 50	↑ 50	↓ 50
DILUTION	1:1	1:2	1:4	1:8	1:16	1:32

6. Gently shake the Latex Reagent vial and place a drop on section 1 through 6 of the slides containing the different sample dilution.
7. Mix both drops using a stirrer covering the whole surface of the slide section.
8. Gently rotate the slide for 3 minutes manually or on a rotary shaker set at 80-100 rpm.
9. Look for the presence of agglutination after the aforementioned period of time.

INTERPRETATION OF RESULTS

QUALITATIVE TECHNIQUE

The presence of agglutination indicates a clinically significant concentration of infectious mononucleosis heterophile antibodies in the sample.

SEMIQUANTITATIVE TECHNIQUE

The approximate titer will correspond to the highest sample dilution that still presents a clearly visible agglutination (see diagram).

POSITIVE REACTIONS

- 3+ Large clumping with clear background.
- 2+ Moderate clumping with fluid slightly opaque in background.
- 1+ Small clumping with opaque fluid in background.

NEGATIVE REACTIONS

No visible clumping, uniform suspension.

QUALITY CONTROL

This test should be performed by applicable local, state, or federal regulations or accrediting agencies.

Control of the latex reagent:

1. Before performing a set of determinations it is advisable to test the latex reagent with each of the controls, positive and negative, included in the kit.
2. Both controls should be used following the steps outlined in the QUALITATIVE TECHNIQUE.
3. The reaction between the Positive Control and the Latex Reagent should show a clear agglutination, different from the uniform appearance of the negative control.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

EXPECTED VALUES

Studies^{11, 12} on the presence of infectious mononucleosis heterophile antibodies in normal donors show the incidence of the disease in from 0.9 to 1.7% of the population. As the presence of antibodies indicates a relatively recent infection these results suggest that the true incidence of the disease is higher than the diagnosed cases.

LIMITATIONS OF THE PROCEDURE

1. As with all diagnostic assays, the results of the MONOSPOT Latex assay should be interpreted in light of the clinical symptoms shown by the patient.
2. Occasionally detectable levels of heterophile antibodies are late in developing in patients symptomatic for infectious mononucleosis. If symptoms persist it is recommended to repeat the assay in several days. Some patients may remain persistently negative, especially children and adolescents. It has been reported that only 80 to 90% of adults and less than 50% of young children develop heterophile antibodies.
3. Detectable levels of heterophile antibodies may persist for months, and more rarely for years, in some individuals.
4. Although titers of heterophile antibodies have little relation with the severity of infection, the semiquantitative procedure can be used to follow the evolution of the disease.

SPECIFIC PERFORMANCE CHARACTERISTICS

Evaluations comparing MONOSPOT Latex to a commercially available differential hemagglutination test were performed to determine the sensitivity and specificity of the reagent.¹³ A differential hemagglutination slide test was used to identify 106 positive and 114 negative sera for the study. Discrepancies between the results given by MONOSPOT Latex and the hemagglutination test were resolved using Epstein-Barr Virus (EBV) specific serological assays.

In these assays, the titers of specific antibodies to the EBV capsid antigen (both IgG and IgM), EBV early antigen (both diffused -D- and restricted -R-) and EBV nuclear antigen were determined. The results of these assays specified whether EBV infections were recent or acute (in which case the sera was considered positive), or if antibodies were absent or relics of an old infection (in which case the serum was considered negative). Twelve of the 100 nondiscrepant positive sera and five of the 106 nondiscrepant negative sera were also analyzed using these same EBV specific serological assays. In all cases the EBV specific assay confirmed the positivity or negativity of the samples.

Compared with hemagglutination, MONOSPOT Latex was found to have a sensitivity of 94% and a specificity of 93%. Assuming that the concordant results of the EBV serology performed on nondiscrepant sera applies to all the samples tested, it can be inferred that the sensitivity of MONOSPOT Latex relative to EBV specific tests is 99% and its specificity relative to the same is 93%. There was only one result negative with MONOSPOT Latex, and positive with EBV serology. It was determined to be a recent infection and not an acute infection due to the absence of anti-VCA IgM and nuclear antigen. For this reason this negative result cannot be interpreted as prozone effect.

In a separate study MONOSPOT Latex was compared to a qualitative horse red cell slide test involving a total of 224 EDTA plasma samples. There was complete agreement in test results which included 51 positive and 173 negative samples. Overall, the results of this study indicate clearly that MONOSPOT Latex is highly sensitive and specific for the diagnosis of infectious mononucleosis.

A panel of 10 positive serum samples was tested on three consecutive days using the semiquantitative technique. The results of the study indicate that MONOSPOT Latex has 100% precision. The error of repeated titrations was expected to be only one doubling dilution.

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REFERENCES

- 1. Davidsohn I, Stern K and Kashiwagi C. The differential test for infectious mononucleosis. Amer J Clin Pathol 1951;21:1101-1113.
- 2. Paul JR and Bunnell WW. The presence of heterophile antibodies in infectious mononucleosis. Amer J Med Sci 1932;183:90-104.
- 3. Beer P. The heterophile antibodies in infectious mononucleosis and after the injection of serum. J Clin Invest 1936;15:591-599.
- 4. Bailey GH and Raffel S. Hemolytic antibodies for sheep and ox erythrocytes in infectious mononucleosis. J Clin Invest 1935;14:228-244.
- 5. Forssman J. Die Herstellung hochwertiger spezifischer schaffamolsine ohne Verwendung von schaffblut. Biochem Ztschr. 1911;37:78-115.
- 6. Davidsohn I. Serologic diagnosis of infectious mononucleosis J.A.M.A. 1937;108:289-295.
- 7. Stuart CA, Burgess AM, Lawson HA and Wellman HE, Some cytologic and serologic aspects of infectious mononucleosis Arch Intern Med 1934;54:199-214.
- 8. Henle G, Henle W and Diehl V. Relation of burkitt tumor associated herpes-type virus to infectious Mononucleosis Proc Nat Acad Sci USA 1968;59:94-101.
- 9. Hoagland RJ. The transmission of infectious mononucleosis Amer J Med Sci 1955;229:262-272.
- 10. Biosafety in microbiological and biomedical laboratories CDC/NIH Manual 1984.
- 11. Barret AM. The serological diagnosis of glandular fever (infectious mononucleosis): a new technique J Hyg 1941;41:330.
- 12. Virtanen S. Incidence of infectious mononucleosis antibodies in blood donors Acta Pathol Microbiol Immunol Scand 1962;56:53-56.
- 13. Tilton RC, Dias F and Ryan RW. Comparative evaluation of three commercial tests for detection of heterophile antibody in patients with infectious mononucleosis J Clin Microbiol 1988;26:275-278.

Manufactured By

Meridian Bioscience, Inc.
Corporate Office
3471 River Hills Drive
Cincinnati, Ohio 45244 USA
Telephone: 513.271.3700
Orders/Customer Service:
800.543.1980
Technical Support Center:
800.343.3858
Information Fax: 513.272.5432
Ordering Fax: 513.271.0124

Authorized Representative

Meridian Bioscience Europe S. r. L
Via dell' Industria, 7
20020 Villa Cortese, Milano
ITALY
Tel: +39 0331 43 36 36
Fax: +39 0331 43 36 16
Email: info@meridianbioscience.eu
WEB: www.meridianbioscience.com/eu

Meridian Bioscience Europe s.a./n.v.
2 Avenue du Japon - 1420 Braine l'Alleud
BELGIUM
Tel: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
Email: info.bnl@meridianbioscience.eu

Meridian Bioscience Europe France
34 rue de Ponthieu - 75008 Paris
FRANCE
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 9 70 06 62 10
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
Postbus 301 - 5460 AH Veghel
NETHERLANDS
Tel: +31 (0) 411 62 11 66
Fax: +31 (0) 411 62 48 41
Email: info.bnl@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfils the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet
			CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risques de danger / Precaución: Peligroso / WARNUNG: Risikogefahr
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Numero de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Eingrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabbricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" test / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	SOLN STOP	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limit / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Coniugato enzimatico / Conjugué enzymatique / Conjugado enzimático / Enzymkonjugat
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / testgerät	REAG	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
CONJ	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
SUBS	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentration 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
Rx Only	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig	DET REAG	Detection Reagent / Reagente Diretto / Réactif de Detection / Reactivo de Detección / Nachweis Reagenz
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / ne pas utiliser si le paquet est endommagé / No use si el paquete esta dañado / Nicht verwenden, wenn die Verpackung beschädigt ist	TUBE	Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.