REF 901612, 981612

IVD

Rx Only

INTENDED USE

Para-Pak SVT (Single Vial Transport) is a convenient one-vial system for the routine collection, transportation, preservation and examination of stool specimens for intestinal parasites. Kit systems are designed for easy use by individuals not trained in microbiological procedures and afford an excellent means of minimizing the adverse effects of delay in specimen transportation. Stool specimens preserved in Para-Pak SVT can be used for direct examination, concentration, permanent stain and in fecal immunoassays. The fixative affords an excellent means of minimizing the adverse effects of delay in specimen transport.

SUMMARY AND EXPLANATION OF THE TEST

Diagnosis of intestinal parasitic disease is confirmed by recovery and identification of helminth eggs and larvae, or protozoan trophozoites and cysts in the clinical parasitology laboratory. Timely collection and transportation of "fresh" stool specimens to the laboratory cannot always be assured. Workload conditions and priorities in clinical laboratories frequently do not permit immediate examination of "fresh" specimens. Procedures such as incubation, refrigeration, or freezing will not guarantee the recovery of all diagnostic stages of all parasites.

Proper use of Para-Pak SVT thus assures the parasitologist that diagnostic stages of intestinal parasites will be preserved if present in fecal material.

BIOLOGICAL PRINCIPLES

Para-Pak SVT provides a multi-purpose PVA-free, mercury-free and formalin-free fixative/preservative suitable for a variety of parasitological procedures. Para-Pak SVT preserved specimens may be examined directly, concentrated, permanently stained, or used in fecal immunoassays.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

Each vial contains 15 mL Para-Pak SVT with a built-in collection spoon. Simple directions for patients and nursing personnel are also provided.

Catalog# 901612

SVT (item# 9016)

Catalog# 981612

ULTRA SVT (item# 9016)

Surfactant (item# 9601)

MATERIALS PROVIDED SEPARATELY

1. Para-Pak® CON-Trate® Concentration System, Meridian Bioscience Catalog 960200, 960500
2. Para-Pak® SpinCon® Stool Concentration System, Meridian Bioscience Catalog 972200, 972500
3. Para-Pak® Trichrome Stain, Meridian Bioscience Catalog 400101, 400112
4. Para-Pak® EcoStain®, Meridian Bioscience Catalog 801400
5. Para-Pak® Mayer's Albumin, Meridian Bioscience Catalog 800700

MATERIALS NOT PROVIDED

1. Physiological Saline
2. Lugol's Iodine
3. 70% Ethanol
4. Cotton Tipped Applicator Sticks
5. Transfer Pipettes
6. Microscope Slides and Cover Slips
7. Fecal Immunoassay Products

EQUIPMENT NOT PROVIDED

1. Centrifuge
2. Microscope

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. For best results, stool specimens should be added to preservative as soon as possible.²
3. Avoid contact of Para-Pak SVT solution with the skin and eyes. Avoid breathing fumes. Should contact occur, flush with running water. If irritation develops, contact a physician.
4. Para-Pak SVT solution is poisonous. If ingested, call the local poison center or a physician immediately.
5. Follow Biosafety Level 2 and Good Laboratory practices. Treat all specimens as capable of transmitting infectious agents. Do not eat, drink or smoke in areas where specimens or reagents are being handled.
6. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
7. Para-Pak SVT may be considered a hazardous waste in some areas. Consult local authorities for the regulations on disposal in your area.
8. Any serious incident that has occurred in relation to the device should be reported to Meridian Bioscience, Inc., 3471 River Hills Drive, Cincinnati, Ohio 45244 USA or Technical Support Center 800.343-3858 and competent authority of the EU Member State in which the clinician and/or patient is established.
9. IMPORTANT: See SDS for additional safety and hazard information.

HAZARD AND PRECAUTIONARY STATEMENTS

 Para-Pak Single Vial Transport	<p>Signal word Danger</p> <p>Hazard Statements</p> <p>H302 - Harmful if swallowed H332 - Harmful if inhaled H412 - Harmful to aquatic life with long lasting effects H311 - Toxic in contact with skin H370 - Causes damage to organs H226 - Flammable liquid and vapor Contains Isopropyl alcohol, Methyl alcohol</p> <p>Precautionary Statements - EU (§28, 1272/2008)</p> <p>P210 - Keep away from heat/sparks/open flames/hot surfaces. - No smoking P370 + P378 - In case of fire: Use dry sand, dry chemical or alcohol-resistant foam for extinction P312 - Call a POISON CENTER or doctor/ physician if you feel unwell P303 + P361 + P353 - IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/shower P304 + P340 - IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing P301 + P312 - IF SWALLOWED: Call a POISON CENTER or doctor/ physician if you feel unwell P330 - Rinse mouth P378 - Use foam for extinction P405 - Store locked up P403 + P233 - Store in a well-ventilated place. Keep container tightly closed P501 - Dispose of contents/ container to an approved waste disposal plant.</p>
---	--

 Reagent A (MucoPenX), CON-Trate Reagent A and Macro-CON Surfactant	<p>Signal word Warning</p> <p>Hazard Statements</p> <p>H302 - Harmful if swallowed Contains Polyethylene glycol 1,1,3,3-tetramethylbutylphenyl ether</p>
---	---

SHELF LIFE AND STORAGE

Expiration dating is indicated on the outer packaging, where appropriate. Store at room temperature (15-30 C). Excessive heat and cold should be avoided.

SPECIMEN COLLECTION AND PREPARATION

1. The patient should be cautioned against the use of antacids, barium, bismuth, antidiarrheal medication, or oily laxatives prior to collection of the specimen.
2. To assure recovery of parasitic elements which are passed intermittently and in fluctuating numbers, three specimens spaced a few days apart must be examined. In the case of hospitalized patients it is suggested that all fecal passages be collected for a designated length of time to avoid prolonging the hospital stay.^{3,4}
3. The specimen is ideally passed into a bedpan but must not be contaminated with urine. Alternatively, a large plastic bag or "saran wrap" may be placed over the toilet seat opening and the specimen passed into the bag. A thoroughly cleaned and dried milk carton, cut so as to remove the upper two thirds of the carton may also be used. It will be easier to collect the specimen if the water supply to the toilet is shut off and water drained from the bowl by flushing twice.
4. An appropriate (i.e. bloody, mucoid, watery) area of stool should be selected and sampled with the collection spoon provided in the cap of the vial. Sufficient stool is added to each container to bring the liquid level up to the "Fill to Here" line. This will result in approximately 5 mL of sample. To insure adequate sampling of a formed stool, material should be removed from the sides, ends, and middle of the bolus.
5. Agitate each specimen with the spoon along the sides of the vial, tighten the cap and shake firmly to insure that the specimen is adequately mixed. When mixing is completed the specimen should appear homogenous. Allow the specimen to stand for at least 30 minutes at room temperature prior to processing to allow adequate fixation.
6. Return the vial(s) to their container, seal the container, and label appropriately.

TEST PROCEDURE

- Para-Pak SVT lends itself to a wide variety of procedures in common use. The following discussion is not exhaustive and alternatives may be found. While variations exist from lab to lab, a thorough examination should include the following steps:
1. Gross examination: Record the presence of blood, worms, mucus, or proglottids.
 2. Direct microscopic examination from the Para-Pak SVT preserved specimen:
 - a. Place a clean glass slide on a sheet of newsprint.
 - b. Add a drop of saline or Lugol's Iodine to the slide.
 - c. Add a representative sample of Para-Pak SVT preserved specimen to the drop of saline or Lugol's Iodine and mix thoroughly with the collecting spoon. The newsprint must be just legible through the slide.
 - d. Place a double width coverslip on the suspension and examine immediately.
 3. Concentration Procedures (for permanent slide preparation):
 - a. **Para-Pak SpinCon Stool Concentration System**
 1. Thoroughly mix the Para Pak SVT preserved sample in the vial until it appears homogenous. Follow the Specimen Collection and Specimen Processing sections of the Para-Pak® SpinCon® Stool Concentration System Package Insert for product catalog 972200 and 972500.
 2. **Proceed with Permanent Slide Preparation.**
 - b. **Para-Pak® CON-Trate® System**
Refer to the Para-Pak® CON-Trate® System Package Insert for product catalog 960200, 960500 for additional information.
 1. Thoroughly mix the Para-Pak SVT preserved sample in the vial until it appears homogenous.
 2. Add 4 drops of CON-Trate® Reagent A to the specimen in the fixative vial. (Up to 8 drops of Reagent A may be added if the specimen is highly mucoid)
 3. Cap and mix the contents thoroughly shaking.
 4. Insert one of the CON-Trate® Filtering devices into the top of one of the disposable centrifuge tubes provided.
 5. Pour fecal suspension through the filtering device into the centrifuge tube. Usually 3 mL in the centrifuge tube is sufficient unless the fecal suspension is thin.
 6. Discard filtering device, add 10 mL saline and centrifuge at 500 xg (1800-2200 rpm) for 10 minutes. Decant the supernatant fluids, retaining the sediment. **Proceed with Permanent Slide Preparation.**
 7. Re-suspend the sediment in 9.0 mL physiological saline.
 8. Add approximately 3 mL of CON-Trate® Reagent B, cap the tube and shake for 30 seconds. Invert the tube while shaking. **CAUTION:** Pressure may build up within the tube during shaking. Carefully release the pressure by opening the cap on the centrifuge tube away from your person.
 9. Centrifuge the tube at 500 xg (1800-2200 rpm) for 10 minutes. Examination of the tube after centrifugation should reveal four distinct layers from the top down:
 - a. Top Layer: Reagent B
 - b. Second Layer: a "plug" of fecal debris
 - c. Third Layer: a discolored aqueous layer
 - d. Bottom Layer: a sediment layer, containing the parasites. The final sediment remaining should be approximately 0.25 mL.
 10. Hold the tube in a vertical position. Free the plug of debris by ringing with a wooden applicator stick. Decant the upper layers, leaving the sediment. DO NOT TURN THE TUBE UPRIGHT UNTIL THE SIDES OF THE TUBE HAVE BEEN CLEANED WITH COTTON TIPPED SWABS. **NOTE: The prepared sediment should be re-suspended in approximately an equal volume of un-used Para-Pak SVT if the sample is not to be examined immediately.**
 11. Transfer a portion of the sediment to a clean glass microscope slide and prepare the mount of choice. Examine by microscope. Consult an appropriate reference for proper preparation and examination.^{2, 4, 5}
NOTE: If the prepared slides are too dense (newsprint should be legible through them), additional saline or Para-Pak SVT may be added.
 - c. **Alternate Specimen Concentration:**
If the Para-Pak® CON-Trate® or SpinCon® Systems are not available, the following procedure for specimen concentration prior to permanent slide preparation may be utilized:
 1. Thoroughly mix the Para-Pak SVT preserved specimen in the vial until it appears homogenous.
 2. Strain approximately 2-3 mL (amount will vary in direct proportion to the density of the suspension) of the suspended material through 1 layer of narrow mesh or 2 layers of wide mesh gauze into a suitable centrifuge tube.
 3. Centrifuge 10 minutes at 500 xg (1800-2200 rpm).
 4. Decant the supernatant fluid. Approximately 0.5 to 1.0 mL of sediment should remain. If necessary, remove sediment or add more strained suspension and repeat steps 2-4.
 5. Mix the sediment and prepare the permanent slide as described in **Permanent Slide Preparation.**

Refer to Instructions for Use, Para-Pak Trichrome Stain Set product catalog 400101 and 400112 for complete instructions.

	Staining Procedure	Time
1	Trichrome Stain	6-8 minutes
2	90% Alcohol Acidified	2 quick dips
3	95% Alcohol	Dip twice
4	95% Alcohol	5 minutes
5	100% Ethanol or Carbol-xylene	10 minutes
6	Xylene	10 minutes
7	Mount with coverslip using the mounting medium of choice	

5. Fecal Immunoassays

When using Para-Pak SVT-preserved specimens in fecal immunoassay procedures, shake the Para-Pak SVT vial and allow to stand undisturbed for a minimum of 5 minutes. Remove specimen for immunoassay testing from the top of the vial, taking care to avoid stool particulate material.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

1. The Para-Pak SVT vial should contain approximately 15 mL of fluid to insure 1:3 stool to preservative ratio.
2. When a Para-Pak SVT-fixed film of stock trophozoite, stock cyst or human buffy coat is stained, the organisms or cells should appear well-fixed and defined.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor.

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REFERENCES

- Garcia LS. Diagnostic medical parasitology. 5th Edition. 2007. ASM Press.
- Garcia LS, Voge M. Diagnostic clinical parasitology: I. proper specimen collection and processing. Am J Med Technol 1980;46:459-467.
- Brooke MM. "Intestinal and urogenital protozoa", Manual of Clinical Microbiology, ASM, Washington, D.C., Second Edition, 1974;582-601.
- Melvin DM, Brooke MM. "Laboratory procedures for the diagnosis of intestinal parasites", USDHEW, 1980;79:8282, CDC, Atlanta, GA., 23-65.
- Burrows RB. "Microscopic diagnosis of the parasites of man", Yale University Press, New Haven, 1965;319-381.

SN10750

REV. 03/23

 Manufactured By	<p>Meridian Bioscience, Inc. 3471 River Hills Drive Cincinnati, OHIO - 45244 USA www.meridianbioscience.com</p> <p>Contacts: Main Telephone (+1) 513.271.3700 Customer Service/Orders 800.543.1980 Technical Support Center 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124 E-mail: info@meridianbioscience.com</p>
 Authorized Representative	<p>Meridian Bioscience Europe, SRL Via Dell'Industria 7, 20035 Villa Cortese (Milano) ITALY www.meridianbioscience.com</p> <p>Contacts: Main Telephone (+39) 0331.433636 E-mail: info@meridianbioscience.eu Technical Support: MBE-TechService@meridianbioscience.eu Customer Service/Orders: <ul style="list-style-type: none"> • For Italian Customers: ordini@meridianbioscience.com • For Distributors / International Customers: Export.CustomerService@meridianbioscience.eu </p>
UK Authorised Representative	<p>Launch Diagnostics Ash House Ash Road Longfield DA3 8JD UK</p>
 Authorized Representative	<p>MedEnvoy Switzerland Gotthardstrasse 28 6302 Zug Switzerland</p>

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de símbolos, Zeichenerklärung)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.