# merifluor®

## EBV-EA IgG IFA

### Immunofluorescence assay for detection of IgG antibodies to EBV-EA in human serum

REF EA101

IVD

R<sub>x</sub> Only

#### INTENDED USE

The MERIFLUOR EBV-EA IgG IFA Test is a sensitive and rapid immunofluorescence method for the qualitative and quantitative detection of antibodies to early antigen (EA) of Epstein-Barr virus (EBV) in human serum. When performed according to instructions, the MERIFLUOR EBV-EA IgG IFA Test detects IgG antibodies to the diffuse (D) and restricted (R) components of the EBV-EA complex. It is of value in providing supportive information for the diagnosis of EBV-associated disease.

#### SUMMARY AND EXPLANATION OF THE TEST

Epstein-Barr virus is a widely-disseminated human pathogen which commonly causes infection in early childhood, resulting in asymptomatic or mild disease. Individuals who contract infection with EBV in early adult life frequently develop classic infectious mononucleosis (IM). Approximately half of the patients with IM develop splenomegaly with a smaller percentage having hepatomegaly or both.<sup>1</sup> Abnormal liver function is more marked with EBV mononucleosis than in CMV-associated mononucleosis and both disease conditions must be considered in the differential diagnosis of hepatitis.<sup>2</sup>

Infections with EBV can result in complications involving the neurologic, cardiac, ocular, respiratory, hematologic, digestive and renal systems.<sup>3</sup> Associated neurologic syndromes include Bell's palsy, Guillain-Barre syndrome, meningoencephalitis, Reye's syndrome, myelitis, cranial nerve neuritis, and psychotic disorders.<sup>4,5</sup> Bulbar involvement with ensuing respiratory paralysis can be fatal.<sup>5</sup> Burkitt's lymphoma, nasopharyngeal carcinoma and neoplasias of the thymus, parotid gland and supraglottic larynx have also been associated with EBV infections.<sup>6,7,8,9</sup>

Like other herpes viruses, EBV establishes a latent infection.<sup>10</sup> The genome resides in B cells of the lymphoreticular system and epithelial cells of the oropharynx.<sup>11</sup> Blood transfusion from an immune donor to a non-immune recipient may produce a primary infection in the recipient known as IM post-perfusion syndrome.<sup>12</sup> Reactivation of latent infection has been implicated in persistent illness referred to as chronic mononucleosis or chronic EBV infection. Persistent infection with EBV may be a possible cofactor in the development of AIDS in asymptomatic individuals infected with human immunodeficiency virus (HIV).<sup>13</sup>

Infection with EBV results in the expression of viral capsid antigen (VCA), early antigen (EA), and nuclear antigen (NA) with corresponding antibody responses. The synthesis of EA occurs early in viral replication, prior to DNA synthesis.<sup>14</sup> In acute disease, IgM and IgG antibodies to VCA are present and IgG antibodies to EA may be present.<sup>5</sup> The presence of antibodies to EA in conjunction with antibodies to NA helps substantiate reactivation of EBV infection.<sup>13,15</sup> Elevated levels of IgG to EA and VCA frequently occur in patients with chronic mononucleosis and in immunocompromised or elderly patients upon reactivation of latent EBV infection.<sup>13,16</sup> A persistent elevated antibody response to EA in the presence of antibodies to NA is a serologic feature indicative of a reactivation or persistently active infection with EBV.<sup>13</sup>

The indirect fluorescent antibody (IFA) technique is a valuable serological procedure for the determination of antibody status to the EA of EBV. The MERIFLUOR EBV-EA IgG IFA Test detects antibodies to the D and R components of EA and is not intended to differentiate between antibodies to these antigens. Screening for the presence of antibodies to EA and other EBV-associated antigens can provide important diagnostic information.

#### BIOLOGICAL PRINCIPLES

The indirect fluorescent antibody (IFA) method is used in the MERIFLUOR EBV-EA IgG IFA Test. Patient serum is reacted with the virus-infected cell substrate and if antibodies to the EA of EBV are present they will bind to the antigen substrate and not rinse off. Subsequently, when fluoresceinated antihuman IgG is added to the reaction site, it will bind to the IgG antibodies, causing the EA-antibody complex to fluoresce when viewed through the fluorescence microscope.

#### REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

- EBV-EA Substrate Slides:** With Raji cells from a Burkitt's lymphoma fixed onto each well. Approximately 15 to 25% of the cells express EBV-EA to permit easy reading and optimal contrast. Positive cells are chemically induced in the presence of inhibitors of DNA synthesis to ensure the expression of EBV-EA and prevent the production of infectious EBV virions. Negative cells are Raji cells which have not been chemically induced. The individually packaged slides are ready for use when removed from the foil packets. Slides stored at 2-8 C are stable until the date stated on the slide package. Do not use a slide if the desiccant indicator (line in center of desiccant) has changed from blue to pink.
- Phosphate Buffered Saline (PBS powder):** Rehydrate to 1 liter with distilled water. When rehydrated, the PBS is a 0.01 M phosphate buffer with a pH of 7.5 and contains sodium azide. Rehydrated PBS stored at 2-8 C is stable indefinitely.
- Serum Diluent:** Rehydrate with 20 mL PBS solution. Allow at least 10 minutes for contents to dissolve. Mix thoroughly by inversion before each use. The reconstituted diluent contains caprine proteins and stabilizers and is stable for 6 months at 2-8 C.
- Positive and Negative Control Sera (See WARNINGS step 4):** Lyophilized. Each human control serum is reconstituted with 1 mL PBS to give a 1:10 working dilution. The titer is stated on the label of the Positive Control. The reconstituted sera are stable for 6 weeks at 2-8 C or 8 months at  $\leq -20$  C.
- Conjugate:** Lyophilized. Fluorescein-labeled antihuman IgG (goat) is reconstituted with 3 mL of PBS. The conjugate contains less than 0.1% Evans blue counterstain and is prefiltered for use with each kit. The reconstituted conjugate can be stored up to 2 weeks at 2-8 C or aliquoted and stored up to 8 months at  $\leq -20$  C. Thawed aliquots should not be refrozen.
- Mounting Fluid:** The mounting fluid, which is a 0.2 M phosphate-buffered glycerol formulation at pH 8.0, is designed to minimize elution of counterstain. It is ready for use.

#### MATERIALS REQUIRED BUT NOT PROVIDED

- Serology laboratory supplies: 12 x 75 mm test tubes, test tube rack, serological and Pasteur pipettes
- Volumetric flask, 1 liter
- Wash bottle
- Distilled water
- Staining dish
- Cotton-tipped swabs, absorbent towels
- Moist chamber
- Cover slips, no. 1, 22 x 50 mm
- 37 C Incubator
- Fluorescence microscope. A FITC blue light excitation filter and a 515 nm barrier filter, or any comparable filter system, is suggested for transmitted-light microscopy using a darkfield condenser and for incident-light microscopy using a 500 nm dichroic mirror.

#### PRECAUTIONS

All reagents are for in vitro diagnostic use only.

#### WARNINGS

- The conjugate contains Evans Blue dye as a counterstain. This dye may be carcinogenic. Avoid contact with the skin.
- The PBS Powder contains sodium azide, which may react with lead or copper plumbing to form highly explosive metal azide compounds. When disposing of reagents through plumbing fixtures, flush with a large volume of water to prevent azide build-up in drains.
- Sodium azide is very toxic in case of ingestion. In contact with acid it emits a very toxic gas. In case of contact with the skin wash immediately and abundantly with water.
- Handle specimens, Controls and the materials that contact them as potential biohazards. The source materials used to prepare the control reagents were tested at the donor level and found to be nonreactive for the antigens or nuclear materials associated with, or antibodies to, the following viruses, as defined by required FDA-licensed tests: Human immunodeficiency viruses Types 1 and 2, Hepatitis B virus and Hepatitis C virus. However, since no test method can offer complete assurance that infectious agents are absent, these materials should be handled at the Biosafety Level 2 as recommended in the CDC/NIH "Biosafety and Microbiological and Biomedical Laboratories".

#### HAZARD and PRECAUTIONARY STATEMENTS

  MERIFLUOR PBS Powder	<b>Signal Word</b> Danger <b>Hazard Statements</b> H301 – Toxic if swallowed H310 – Fatal in contact with skin H411 – Toxic to aquatic life with long lasting effects Contains Sodium azide <b>Precautionary Statements - EU (§28, 1272/2008)</b> P301 + P310 – IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician P280 – Wear eye protection/ face protection P361 – Remove/Take off immediately all contaminated clothing
--	---

#### SHELF LIFE AND STORAGE

- The expiration date is indicated on the kit label. Store the kit at 2-8 C and return the kit promptly to the refrigerator after each use.
- The kit is ready for use after reconstitution of reagents.
- If reconstituted conjugate and control sera are to be stored at  $\leq -20$  C, self-defrosting freezers are not recommended.
- Reagents should not be refrozen once thawed.
- Individual slide packets should remain sealed until just before use.
- Reagents and antigen slides should not be used beyond stated expiration dates.
- If lyophilized reagents show evidence of rehydration, they should not be used.
- Incubation times and temperatures other than those specified may give erroneous results.
- The components of this kit have been tested and standardized as a unit. Use of components from other lots or other manufacturers may yield unsatisfactory results.

## SPECIMEN COLLECTION AND PREPARATION

Blood obtained aseptically by venipuncture should be allowed to clot at room temperature and then be centrifuged at 1500 x g for 10 minutes at room temperature. The serum should be separated as soon as possible and refrigerated (2-8 C), or stored frozen ( $\leq -20$  C) if not tested within one week. Serum samples should not be repeatedly frozen and thawed. Self-defrosting freezers are not recommended as storage units. The use of sera exhibiting hemolysis, lipemia or microbial growth is not recommended.

Sera should be screened at a 1:10 dilution. The screening dilution must be prepared in Serum Diluent. Sera which are positive at the screening dilution can be titrated using serial two-fold dilutions. The titration dilutions should be made from the 1:10 screening dilution and may be prepared using either Serum Diluent or PBS.

## TEST PROCEDURE

1. Prepare a 1:10 dilution of the patient serum with Serum Diluent.
2. Remove slide from foil packet.
3. Using a Pasteur pipette, add just enough appropriately diluted serum (approximately 15  $\mu$ L) to cover each reaction site.
4. Place the slide in a moist chamber and incubate at 37 C for 30 minutes.
5. Rinse slide briefly with a gentle stream of PBS (avoid directing PBS at wells) and then immerse in PBS for 5 minutes.
6. Air dry the slide by standing on end on absorbent toweling.
7. Wipe between wells with a cotton-tipped swab moistened with distilled water.
8. Add just enough conjugate (approximately 15  $\mu$ L) to cover each reaction site.
9. Incubate the slide at 37 C in a moist chamber for 30 minutes.
10. Rinse and dry the slide as described in steps 5 and 6.
11. Just before reading, add a sufficient amount of mounting fluid to each well to entirely cover the substrate when the cover slip is applied. Cover with a 22 x 50 mm, no. 1 cover slip.  
**NOTE:** Addition of too much mounting fluid may cause excessive movement of the cover slip. Addition of insufficient mounting fluid to a well may make part of the substrate unreadable.
12. Read slide at 150X to 200X magnification within 24 hours. Store in the dark until read and then discard.

## INTERPRETATION OF RESULTS

The reaction is POSITIVE when 15 to 25% of the cells in each field exhibit a yellow to greenish-yellow cytoplasmic or whole cell fluorescence. The remaining cells provide a contrasting red background. The EBV-EA antibody titer is the reciprocal of the highest dilution of serum which produces a 1+ fluorescence. A positive reaction at a dilution of 1:10 indicates the presence of IgG antibody to EBV-EA.

The reaction is NEGATIVE when cells do not fluoresce yellow to greenish-yellow but appear red due to the counterstain or exhibit less than a 1+ reaction at the screening dilution. A negative reaction at the 1:10 dilution indicates no detectable IgG antibody to EBV-EA.

Test results for antibodies to EBV-EA should be evaluated in relation to patient symptoms, clinical history, and antibody responses to VCA and NA of EBV to establish a diagnosis. The positive/negative cut-off has been set at the 1:10 dilution as an indicator of the presence of IgG antibody to EBV-EA. It is recognized that the presence of antibody at a titer of 10 is not necessarily related to current symptoms due to infection with EBV. Some normal healthy individuals may have titers of 10.

Some sera produce reactions that cannot be determined to be either positive or negative due to the high levels of non-specific fluorescence (NSF) that they produce. These sera generally cause fluorescence in >50% of the cells at the 1:10 dilution, obscuring the true test reaction. Further dilutions of these sera may alleviate the NSF problem and allow the reading of a positive reaction; however, negative reactions at these dilutions cannot be interpreted as a true negative specimen for IgG antibody to EBV-EA since the 1:10 dilution was unreadable. A serum which produces a negative reaction at higher dilutions but is unreadable at the 1:10 dilution due to NSF must be retested using a freshly drawn serum sample or be reported as "indeterminate for IgG to EBV-EA due to NSF".

## QUALITY CONTROL

**This test should be performed per applicable local, state, or federal regulations or accrediting agencies.**

Each kit contains positive and negative control sera which should be incorporated into each test run.

The Negative Control must exhibit less than a 1+ reaction at the 1:10 working dilution.

The Positive Control serum is titrated to provide a standard for checking test sensitivity. It should exhibit a strong positive reaction at the screening dilution and diminish to a 1+ reaction at its stated titer. A day-to-day variance of one two-fold dilution on either side of the stated titer is considered acceptable performance.

**If the test run is invalid because the Positive or Negative Control values do not fall within the limits specified and if the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If the control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor. Test results should not be reported until the condition(s) is corrected.**

Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organizations. Refer to Clinical Laboratory Medicine – In Vitro Diagnostic Medical Devices – Validation of User Quality Control Procedures by the Manufacturer (ISO/FDIS 15198) for further guidance regarding Quality Control Practices.

## EXPECTED VALUES

In acute primary infections with EBV the IgG response to EA closely follows that to VCA but is more transient, declining to low but detectable levels months or, in some cases, years after an IM episode. Antibody titers to EA tend to be lower than those to VCA, and usually do not exceed 640 with this test. The presence of antibodies to EA concurrent with IgM antibody to VCA is evidence of primary infection; however, antibodies to EA do not occur in approximately 10 to 20% of adults and children with acute IM.<sup>13</sup>

Antibodies to EA are characteristically absent or present in low titer in individuals with latent infections. Immunosuppressed and immunodeficient patients and some persons of advanced age often have high titers of IgG antibodies to both EA and VCA, which most likely can be attributed to reactivation of latent virus.<sup>16</sup>

Elevated IgG titers to EA and VCA concurrent with stable levels of antibody to NA is serologic evidence for reactivation of EBV. Fourfold or greater increases in IgG titers to EA and VCA may occur as a nonspecific reactivation of latent virus in some renal transplant patients and AIDS, immunosuppressed or immunodeficient patients.<sup>13</sup>

## LIMITATIONS OF THE PROCEDURE

1. Test results should not be used if nonspecific reactants in the test sample bind to the substrate cells. Nonspecific binding is easily discerned since the staining produced by this reaction involves a majority of the cells whereas specific staining for EBV-EA involves only 15-25% of the cells (those expressing early antigens to EBV). Care must be taken in the interpretation of results when more than 30% of the cells exhibit fluorescence. Under these conditions the positive result must demonstrate the characteristic whole-cell fluorescence present in the positive control as opposed to nonspecific fluorescence that is restricted to only the nucleus or cell membrane.
2. Test results for sera with high antibody titers to nuclear antigen of EBV may be difficult to interpret. These sera typically demonstrate a dull yellow nuclear fluorescence in more than 50% of the cells which is not typical of the greenish-yellow cytoplasmic or whole cell fluorescence characteristic of EA. As indicated above, close observation of the characteristic fluorescence staining pattern of the EBV-EA positive control will help in MERIFLUOR EBV-EA IgG IFA Test interpretation.
3. Fourfold or greater increases in IgG titers to EA and VCA may occur as a nonspecific reactivation of latent virus in some AIDS patients and other immunosuppressed and immunodeficient patients.<sup>13</sup> Elevated levels of IgG antibody to EBV-EA are frequently noted in patients with Burkitt's lymphoma, nasopharyngeal carcinoma, and other malignancies.<sup>16</sup>
4. Endpoint titers in the MERIFLUOR EBV-EA IgG IFA Test can vary depending on the type of microscope, light source, filter system, and age of bulb used to read the reactions.

## SPECIFIC PERFORMANCE CHARACTERISTICS

The MERIFLUOR EBV-EA IgG IFA Test is a sensitive and rapid method for the qualitative and quantitative determination of IgG antibodies to D and R components of EBV-EA in human serum. The sensitivity of each kit lot is controlled by testing with reference sera to ensure reproducible results from lot to lot. The reproducibility of the MERIFLUOR EBV-EA IgG IFA Test was evaluated and found to be within an acceptable range of plus or minus one two-fold dilution. Test results from clinical studies and the testing of reference sera indicate the MERIFLUOR EBV-EA IgG IFA Test is specific for antibodies to EA of EBV and does not cross-react with IgG antibodies to the other herpes viruses. Reference sera with EBV-EA specific IgG titers as high as 640 have been tested without evidence of a prozone phenomenon.

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- REFERENCES**
1. Evans, AS and Niederman, JC Epstein-Barr Virus. In *Viral Infections of Humans: Epidemiology and Control*. 2<sup>nd</sup> ed. Evans, ed. Plenum Medical. New York, 1982.
  2. Horwitz, C.A., et al. Heterophil-Negative Infectious Mononucleosis and Mononucleosis-Like Illnesses. *Am. J. Med.* 1977; 63:947-57.
  3. Sumaya, C.V. and Ench, Y. Epstein-Barr Virus Infectious Mononucleosis in Children. I. Clinical and General Laboratory Findings. *Pediatrics* 1985;75:1003-10.
  4. Fleisher, G., et al. Primary Epstein-Barr Virus Infection in Association with Reye Syndrome. *J. Pediatr.* 1980;97:935-37.
  5. Gottlieb-Stematsky, T. and Glaser, R. Association of Epstein-Barr Virus with Neurologic Diseases. In *Human Herpesvirus Infections: Clinical Aspects*. Glaser and Gottlieb-Stematsky, eds. Marcel Dekker. New York, 1982.
  6. Henle W. and Henle G. Epstein-Barr Virus and Infectious Mononucleosis. In *Human Herpesvirus Infections: Clinical Aspects*. Glaser and Gottlieb-Stematsky, eds. Marcel Dekker. New York, 1982.
  7. Leyvraz, S., et al. Association of Epstein-Barr Virus with Thymic Carcinoma. *N. Engl. J. Med.* 1985;312:1296-99.
  8. Saemundsen, A.K. et al. Epstein-Barr Virus in Nasopharyngeal and Salivary Gland Carcinomas of Greenland Eskimoes. *Br. J. Cancer.* 1982;46:721-28.
  9. Brichacek, B., et al. Association of Some Supraglottic Laryngeal Carcinomas with EB Virus. *Int. J. Cancer.* 1983;32:193-97.
  10. Henle, W. and Henle, G. Epstein-Barr Virus and Infectious Mononucleosis. *N. Engl. J. Med.* 1973;288:263-64.
  11. Sixbey, J.W., et al. Epstein-Barr Virus Replication in Oropharyngeal Epithelial Cells. *N. Engl. J. Med.* 1984;310:1225-30.
  12. Henle, W., et al. Antibodies to Early Antigens Induced by Epstein-Barr Virus in Infectious Mononucleosis. *J. Infect. Dis.* 1971;124:58-67.
  13. Sumaya, C.V. Infectious Mononucleosis and Other EBV Infections: Diagnostic Factors. *Laboratory Management.* 1986;24:37-46.
  14. Ooka, T. The Molecular Biology of Epstein-Barr Virus. *Biomedicine & Pharmacotherapy.* 1985;39:59-66.
  15. Sumaya, C.V. Epstein-Barr Virus Serologic Testing: Diagnostic Indications and Interpretations. *Pediatr. Infect. Dis.* 1986;5:337-42.
  16. Sumaya, C.V. Seriological Testing for Epstein-Barr Virus: Developments in Interpretation. *J. Infect. Dis.* 1985;151:984-87.

SN11746

REV. 09/20

 <p>Manufactured By</p>	<p><b>Meridian Bioscience, Inc.</b>  <b>Corporate Office</b>  <b>3471 River Hills Drive</b>  <b>Cincinnati, Ohio 45244 USA</b>  <b>Telephone: 513.271.3700</b>  <b>Orders/ Customer Service:</b>  <b>800.543.1980</b>  <b>Technical Support Center:</b>  <b>800.343.3858</b>  <b>Information Fax: 513.272.5432</b>  <b>Ordering Fax: 513.271.0124</b></p>
 <p>Authorized Representative</p>	<p><b>Meridian Bioscience Europe S. r. L</b>  <b>Via dell' Industria, 7</b>  <b>20035 Villa Cortese, Milano</b>  <b>ITALY</b>  <b>Tel: +39 0331 43 36 36</b>  <b>Fax: +39 0331 43 36 16</b>  <b>Email: info@meridianbioscience.eu</b>  <b>WEB: www.meridianbioscience.com/eu</b></p>

**Meridian Bioscience Europe s.a./n.v.**  
**2 Avenue du Japon - 1420 Braine l'Alleud**  
**BELGIUM**  
**Tel: +32 (0) 67 89 59 59**  
**Fax: +32 (0) 67 89 59 58**  
**Email: info.bnl@meridianbioscience.eu**

**Meridian Bioscience Europe France**  
**34 rue de Ponthieu - 75008 Paris**  
**FRANCE**  
**Tel: +33 (0) 1 42 56 04 40**  
**Fax: +33 (0) 9 70 06 62 10**  
**Email: info.fr@meridianbioscience.eu**

**Meridian Bioscience Europe b.v.**  
**Postbus 301 - 5460 AH Veghel**  
**NETHERLANDS**  
**Tel: +31 (0) 411 62 11 66**  
**Fax: +31 (0) 411 62 48 41**  
**Email: info.bnl@meridianbioscience.eu**

**INTERNATIONAL SYMBOL USAGE**

You may see one or more of these symbols on the labeling/packaging of this product:

**Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)**

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	<b>CONTROL +</b>	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
<b>LOT</b>	Batch Code / Codice del lotto / Code du lot / Código de lote / chargenbezeichnung	<b>CONTROL -</b>	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
<b>IVD</b>	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	<b>EC REP</b>	Authorized representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfils the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	<b>SMP PREP DIL SPE</b>	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet
			<b>CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risques de danger / Precaución: Peligroso / WARNUNG: Risikogefahr</b>
<b>REF</b>	Catalogue number / Numero di catalogo / Référence du catalogue / Numero de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Eingrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	<b>BUF RXN</b>	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" test / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	<b>SOLN STOP</b>	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	<b>CONJ ENZ</b>	Enzyme Conjugate / Coniugato enzimatico / Conjugué enzymatique / Conjugado enzimático / Enzymkonjugat
<b>SN</b>	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	<b>CONTROL</b>	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
<b>TEST</b>	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / testgerät	<b>REAG</b>	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	<b>BUF WASH</b>	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
<b>BUF</b>	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
<b>CONJ</b>	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	<b>DIL SPE</b>	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
<b>SUBS</b>	Substrate / Substrato / Substrat / Substrato / Substrat	<b>BUF WASH 20X</b>	Wash Buffer Concentration 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
<b>R. Only</b>	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig	<b>DET REAG</b>	Detection Reagent / Reagente Diretto / Réactif de Detection / Reactivo de Detección / Nachweis Reagenz
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / ne pas utiliser si le paquet est endommagé / No use si el paquete esta dañado / Nicht verwenden, wenn die Verpackung beschädigt ist	<b>TUBE</b>	Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß

For technical assistance, call Technical Support Services at (800) 343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at (800) 543-1980.