MERIFLUOR[®]

EBV VCA IgG IFA

An immunofluorescence test for the detection of EBV VCA IgG Antibodies to Epstein-Barr Virus

REF

EB100

IVD

In vitro diagnostic medical device

R. Only

INTENDED USE

An immunofluorescence test for the detection of EBV VCA IgG Antibodies to Epstein-Barr Virus.

SUMMARY AND EXPLANATION OF THE TEST

Epstein-Barr virus (EBV) is recognized as the cause of infectious mononucleosis (IM). Infection with EBV is common in early childhood. Most of these childhood infections result in asymptomatic or mild disease with the development of EBV antibodies and permanent immunity. Unprotected persons who become infected in early adult life, especially those in high socioeconomic groups, generally develop full-blown, symptomatic IM.¹ IM is one of the most important disabling diseases of young adults.²

Infections produced by EBV appear to be restricted to lymphoid cells.³ Like other herpes viruses, EBV establishes a persistent carrier state.⁴ The EBV genome is carried in the peripheral leukocytes indefinitely. Blood transfusion from an immune donor to a nonimmune recipient produces a primary infection referred to as IM post-perfusion syndrome.⁵ Antibodies to viral capsid antigens (VCA) are the predominant and earliest appearing antibodies in EBV infection.³ They develop in all patients with EBV IM, Burkitt's lymphoma, and nasopharyngeal carcinoma.⁶ Elevated antibody titers to EBV-VCA are also found in some patients with Hodgkin's disease, chronic lymphocytic leukemia, systemic lupus erythematosus, sarcoidosis, and Iizumi fever.^{6,7}

While EBV is the etiologic agent of IM, cytomegalovirus and *Toxoplasma gondii* also produce IM-like disease which is often clinically indistinguishable from EBV IM. The differential diagnosis of EBV infection is especially important in those patients who are heterophil negative or have atypical clinical and hematological features or prolonged or recurrent illness.⁸

The indirect fluorescent antibody (IFA) technique is a valuable serological procedure for the determination of immune status as well as active disease. EBV-VCA seroconversion, the presence of EBV capsid IgM, or a significant rise in EBV capsid antibodies are considered evidence of active EBV infection.⁹

BIOLOGICAL PRINCIPLES

The indirect fluorescent antibody (IFA) method is used in the MERIFLUOR EBV VCA IgG IFA Test. Lymphocytes from a Burkitt's patient are used as the VCA substrate. Patient serum is reacted with the substrate and if antibodies to EBV are present they will bind to the antigen substrate and not rise off. Subsequently, when fluoresceinated antihuman IgG is added to the reaction site, it will bind to the IgG antibodies, causing the EBV-infected cells to fluoresce when viewed through the fluorescence microscope.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

- EBV Antigen Substrate Slides:** With HR1 Burkitt's lymphocytic cells fixed onto each well. Approximately 10-20% of the cells express the viral antigen to permit easy reading and optimal contrast. The cells are rendered noninfectious by the PsoralSafe™ process. The individually packaged slides are ready for use when removed from the foil packets. Slides stored at 2-8 C are stable until the date stated on the slide package. Do not use a slide if the desiccant indicator (line in center of desiccant) has changed from blue to pink.
- Phosphate Buffered Saline (PBS) Powder:** Rehydrate to 1 liter with distilled water. The rehydrated solution is a 0.01 M phosphate buffer with a pH of 7.5 and contains sodium azide. Rehydrated PBS stored at 2-8 C is stable indefinitely.
- Positive and Negative Control Sera (See WARNINGS step 4):** Lyophilized. Each human control serum is reconstituted with 1 mL PBS to give a 1:10 working dilution. The titer is stated on the label of the Positive Control. The reconstituted sera are stable for 6 weeks at 2-8 C or 8 months at ≤ -20 C.
- Conjugate:** Lyophilized. Fluorescein-labeled antihuman IgG (goat) is reconstituted with 3 mL of PBS. The conjugate contains counterstain and is preterified for use with each kit. The reconstituted conjugate can be stored up to 2 weeks at 2-8 C or aliquoted and stored up to 8 months at ≤ -20 C. Thawed aliquots should not be refrozen.
- Mounting Fluid:** The mounting fluid, which is buffered to pH 8.0, is a glycerol-water combination formulated to minimize elution of counterstain. It is ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- Serology laboratory supplies: 12 x 75 mm test tubes, test tube rack, serological and Pasteur pipets
- Liter volumetric flask
- Wash bottle
- Distilled water
- Staining dish
- Cotton-tipped swabs, absorbent towels
- Moist chamber
- Cover slips, No. 1, 22 x 50 mm
- Fluorescence microscope. A FITC blue light excitation filter and a 515 nm barrier filter, or any comparable filter system is suggested for transmitted-light microscopy using a darkfield condenser and for incident-light microscopy using a 500 nm dichroic mirror.

PRECAUTIONS

All reagents are for in vitro diagnostic use only.

WARNINGS

- The conjugate contains Evans Blue dye as a counterstain. This dye may be carcinogenic. Avoid contact with the skin.
- The PBS Powder contains sodium azide, which may react with lead or copper plumbing to form highly explosive metal azide compounds. When disposing of reagents through plumbing fixtures, flush with a large volume of water to prevent azide build-up in drains.
- Sodium azide is very toxic in case of ingestion. In contact with acid it emits a very toxic gas. In case of contact with the skin wash immediately and abundantly with water.
- Handle specimens, Controls and the materials that contact them as potential biohazards. The source materials used to prepare the control reagents were tested at the donor level and found to be nonreactive for the antigens or nuclear materials associated with, or antibodies to, the following viruses, as defined by required FDA-licensed tests: Human immunodeficiency viruses Types 1 and 2, Hepatitis B virus and Hepatitis C virus. However, since no test method can offer complete assurance that infectious agents are absent, these materials should be handled at the Biosafety Level 2 as recommended in the CDC/NIH "Biosafety and Microbiological and Biomedical Laboratories".

HAZARD and PRECAUTIONARY STATEMENTS

Refer to the SDS, available at www.meridianbioscience.com for Hazards and Precautionary Statements.

SHELF LIFE AND STORAGE

- The entire kit should be stored in the refrigerator (2-8 C).
- The kit is ready for use after reconstitution of reagents.
- If reconstituted conjugate and sera are to be stored at ≤ -20 C, self-defrosting freezers are not recommended.
- Individual slide packets should remain sealed until just before use.
- Reagents and antigen slides should not be used beyond stated expiration dates.
- The components of this kit have been tested and standardized as a unit. Use of components from other lots or other manufacturers may yield unsatisfactory results.

SPECIMEN COLLECTION AND PREPARATION

Blood obtained by venipuncture should be allowed to clot at room temperature and then be centrifuged. The serum should be separated as soon as possible and refrigerated (2-8 C), or be stored frozen (≤ -20 C) if not tested within one week. Self-defrosting freezers are not recommended as storage units. The use of sera exhibiting hemolysis, lipemia or microbial growth is not recommended.

If the MERIFLUOR EBV VCA IgG IFA Test is used to diagnose active disease, paired sera should be collected and then tested together. The first serum should be obtained as close as possible to the onset of illness and a convalescent serum obtained about 10 days later.

Sera should be screened at a 1:10 dilution. All serum dilutions should be prepared in PBS. Sera which are positive at the screening dilution can be titrated using serial two-fold dilutions.

TEST PROCEDURE

- Prepare a 1:10 dilution of the patient serum with PBS.**
- Remove slide from foil packet.
- Using a Pasteur pipet, add just enough appropriately diluted serum (approximately 15 μ L) to cover each reaction site.
- Place the slide in a moist chamber and incubate at room temperature for 30 minutes.
- Rinse slide briefly with a gentle stream of PBS (avoid directing PBS at wells) and then immerse in PBS for 5 minutes.
- Air dry the slide by standing on end on absorbent toweling.
- Wipe between wells with a cotton-tipped swab moistened with distilled water.
- Add just enough conjugate (approximately 15 μ L) to cover each reaction site.
- Incubate the slide at room temperature in a moist chamber for 30 minutes.
- Rinse and dry the slide as described in Steps 5 and 6.
- Just before reading, add a sufficient amount of mounting fluid to each well to entirely cover the substrate when the cover slip is applied. Cover with a 22 x 50 mm, no. 1 cover slip.**
NOTE: Addition of too much mounting fluid may cause excessive movement of the cover slip. Addition of insufficient mounting fluid to a well may make part of the substrate unreadable.
- Read slide as soon as possible at 150X to 200X magnification.

INTERPRETATION OF RESULTS

The reaction is POSITIVE when 10-20% of the cells in each field exhibit a greenish-yellow fluorescence. The remaining cells provide a contrasting red background. The EBV IgG titer is the reciprocal of the highest dilution of serum which produces a 1+ fluorescence. A positive reaction at a dilution of 1:10 of greater indicates prior exposure to EBV.

The reaction is NEGATIVE when cells do not fluoresce greenish-yellow but appear red due to the counterstain or exhibit less than a 1+ reaction at the screening dilution. A negative reaction at the screening dilution indicates no detectable IgG antibody to EBV.

QUALITY CONTROL

This test should be performed by applicable local, state, or federal regulations or accrediting agencies. Each kit contains positive and negative control sera which should be incorporated into each test run.

The negative control must exhibit less than a 1+ reaction at the 1:10 working dilution.

The positive control serum is titered to provide a standard for checking test sensitivity. It should exhibit a strong positive reaction at the screening dilution and diminish to a 1+ reaction at its stated titer. A day-to-day variance of one two-fold dilution on either side of the stated titer is considered acceptable performance.

If the test run is invalid because the Positive or Negative Control values do not fall within the limits specified and if the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If the control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor. Test results should not be reported until the condition(s) is corrected.

Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organizations. Refer to Clinical Laboratory Medicine – In Vitro Diagnostic Medical Devices – Validation of User Quality Control Procedures by the Manufacturer (ISO/FDIS 15198) for further guidance regarding Quality Control Practices.

EXPECTED VALUES

In classical infectious mononucleosis, IgG antibodies to EBV develop early, reach peak titers within a few weeks, and then slowly decline to lower levels which persist indefinitely.¹¹ Measuring a four-fold or greater rise in titer is diagnostic of active disease. However, a significant titer rise is observed in less than 20% of patients, since peak titers are often reached before the acute phase specimen is drawn.³

Females and tonsillectomized children develop higher titers to EBV than their respective male and nontonsillectomized counterparts. Higher titer levels are also developed by persons with histories of pneumonia or urinary tract infections.¹² A later rise in titer, which can exceed 2560, is the result of secondary disease, such as Burkitt's lymphoma or nasopharyngeal carcinoma.^{13,14} The titer of antibodies to EBV, as with heterophil antibodies, does not reflect the severity of clinical symptoms in IM.¹¹

LIMITATIONS OF THE PROCEDURE

- Antinuclear antibodies present in serum may interfere with the MERIFLUOR EBV VCA IgG IFA Test. However, the diffuse, reticular or peripheral staining produced by the antinuclear antibodies is easily distinguished since it is dull and involves most of the cells.⁹ The homogeneous staining of EBV IgG antibody involves only 10-20% of the cells.
- Endpoint titers can vary depending on the type of microscope, light source, age of bulb, and filter system used to read reactions.
- A single serum titer is of limited diagnostic value in determining active disease. Persisting antibody titers can be misinterpreted in patients later infected with other agents, especially when the original EBV infection was not diagnosed.¹⁰ Also, the upper range of persisting convalescent titers overlaps with the lower range of maximal titers seen in the acute phase.⁹ Meridian's MERIFLUOR EBV VCA IgG IFA Test is recommended when the diagnosis of active disease is to be based on a single serum specimen.

SPECIFIC PERFORMANCE CHARACTERISTICS

The MERIFLUOR EBV VCA IgG IFA Test is a sensitive and rapid method for the determination of immune status. It is also of value in diagnosing active disease by testing acute and convalescent sera to determine sero-conversion or a significant rise in antibody titer. The sensitivity of each kit lot is controlled by testing with reference sera to ensure results which are reproducible within one tube dilution. The test is specific for antibody to EBV and titers are not affected by other febrile viruses, such as coxsackievirus, adenovirus, myxovirus, and other herpes viruses.¹¹

Results of IFA tests, such as the MERIFLUOR EBV VCA IgG IFA Test, are comparable to those of the complement fixation test¹⁵ but do not correlate with results of the heterophil antibody test. In cases of EBV IM, heterophil antibody fails to develop in about 10-15% of adults and in an even greater percentage of children. Heterophil antibody is seldom detected in infants with primary EBV IM.³ By contrast, antibody to EBV-VCA develops in all EBV infections and is not affected by age.¹³ Heterophil antibody is short-lived; whereas EBV IgG antibody persists indefinitely and can be detected by the MERIFLUOR EBV VCA IgG IFA Test as a measure of immune status.

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REFERENCES

- Epstein MA and Achong BG. The Eb Virus. Annual Rev. Microbiol. 1973;27: 413-36.
- Henle G and Henle W. EB Virus in the Etiology of Infectious Mononucleosis. Hospital Practice, July 1970; pp 33-41.
- Henle W, Henle G and Horowitz CA. Epstein-Barr Virus Specific Diagnostic Tests in Infectious Mononucleosis. Human Pathol. 1974;5: 551-65.
- Henle W and Henle G. Epstein-Barr Virus and Infectious Mononucleosis. N. Engl. J. Med. 1973;288: 263-64.
- Henle W, et al. Antibodies to Early Antigens Induced by Epstein-Barr Virus in Infectious Mononucleosis. J. Infect. Dis. 1971;124: 58-67.
- Henle W and Henle G. Epstein-Barr Virus-Related Serology in Hodgkin's Disease. Natl. Cancer Inst. Monogr. 1973;36: 79-84.
- Nishioka K, et al. An Outbreak of a Febrile Disease (Izumi Fever) Associated with High Anti-EB Virus Antibody Titers. Japan J. Exp. Med. 1973;43: 201-07.
- Chang RS. Infectious Mononucleosis. G. K. Hall. Boston. 1980: p93.
- Evans AS. E. B. Virus Antibody in Systemic Lupus Erythematosus. Lancet 1971;1: 1023-24.
- Nikoskelainen J, Klemola E, and Leikola J. Epstein-Barr Virus-IgM Antibody Test in Infectious Mononucleosis. J. Infect. Dis. 1976;134: 313-14.
- Niederman JC, et al. Infectious Mononucleosis: Clinical Manifestations in Relation to EB Virus Antibodies. JAMA 1968;203: 205-09.
- Sumaya CV, et al. Seroepidemiologic Study of Epstein-Barr Virus Infections in a Rural Community. J. Infect. Dis. 1975;131: 403-08.
- Evans AS. The Spectrum of Infections with Epstein-Barr Virus: A Hypothesis. J. Infect. Dis. 1971;124: 330-35.
- Goldman JM, Goodman ML and Miller D. Antibody to Epstein-Barr Virus in American Patients with Carcinoma of the Nasopharynx. JAMA 1971;216: 1618-22.
- Gerber P, et al. Infectious Mononucleosis: Complement-Fixing Antibodies to Herpes-like Virus Associated with Burkitt Lymphoma. Science 1968;161: 173-75.

 Manufactured By	Meridian Bioscience, Inc. USA/Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 Telephone: 513.271.3700 Orders/Customer Service: 800.543.1980 Technical Support Center: 800.343.3858 Information Fax: 513.272.5432 Ordering Fax: 513.271.0124
	Meridian Bioscience Europe S. r. L. Via dell'Industria, 7 20020 Villa Cortese, Milano ITALY Tel: +39 0331 43 36 36 Fax: +39 0331 43 36 16 Email: info@meridianbioscience.eu WEB: www.meridianbioscience.eu

Meridian Bioscience Europe s.a./n.v.
2 Avenue du Japon - 1420 Braine l'Alleud
BELGIUM
Tel: +33 (0) 67 89 59 59
Fax: +33 (0) 67 89 59 58
Email: info.bn@meridianbioscience.eu










Meridian Bioscience Europe France
34 rue de Ponthieu - 75008 Paris
FRANCE
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 1 97 06 62 10
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
Postbus 301 - 5460 AH Veghel
NETHERLANDS
Tel: +31 (0) 411 62 11 66
Fax: +31 (0) 411 62 48 41
Email: info.bn@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis		Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung		Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo medico para diagnóstico in vitro / In-Vitro-Diagnostikum		Authorized Representative in the European Community / Rappresentante autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostika 98/79/EG.		Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenzubereitung, in dem sich Probenverdünnungspuffer befindet
	Catalogue number / Numero di catalogo / Référence ou catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten		Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur für IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen		Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung		Enzyme Conjugate / Coniugato enzimatico / Coniugado enzimático / Coniugado enzimático / Enzymkonjugat
	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer		Assay Control / Contrôle del test / Test de contrôle / Control de Ensayo / Kontrolltest
	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät		Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstelldatum		Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenza / Mise En Garde / Advertencia / Warnhinweise
	Conjugate / Coniugato / Coniugado / Coniugado / Konjugat		Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
	Substrate / Substrato / Substrat / Substrato / Substrat		Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
	Caution: US Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner. / La legge Federale negli Stati Uniti prevede che questo test sia venduto solo a un operatore abilitato ad operare nella salute. / Attention: La vente ou la prescription de ce dispositif par un médecin agréé est soumise à des restrictions par la loi fédérale des US. / Cuidado: Leyes federales de los EU limitan la venta de dispositivos por o por el orden de licenciados en la práctica del cuidado de la salud. / Vorsicht: Der Verkauf oder die Verschreibung dieses Geräts durch einen zugelassenen Arzt unterliegt den Beschränkungen des US Bundesgesetzes.		Detection Reagent / Reagente Diretto / Réactif de Détection / Reactivo de Detección / Nachweis Reagenz

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.

SN11740

REV. 02/16