EBV VCA IgM IFA

An Immunofluorescence Test for the Detection of
IgM Antibodies to Epstein-Barr Virus

REF EB150

IVD

Rx Only

INTENDED USE

The MERIFLUOR EBV VCA IgM IFA Test is a sensitive and rapid immunofluorescence method for the qualitative detection of antibodies to the Viral Capsid Antigen (VCA) of Epstein-Barr virus (EBV) in human serum. It is of value in providing supportive information for the diagnosis of active infection with EBV.

SUMMARY AND EXPLANATION OF THE TEST

Epstein-Barr virus is perhaps the most common and widely disseminated human virus. Virtually everyone acquires antibodies to EBV with advancing age.¹ Responses to primary EBV infection range from subclinical infection in children less than two years to severe debilitating disease in young adults. The predominant symptomatic disease of EBV infection is infectious mononucleosis (IM). IM has an insidious onset with headache, sore throat, fever, enlarged cervical lymph nodes, splenomegaly, malaise, and some form of pharyngitis. Abnormal liver function is more marked with EBV mononucleosis than CMV mononucleosis and must be considered in the differential diagnosis of hepatitis.² Other EBV-induced disease can involve the neurologic, cardiac, ocular, respiratory, hematologic, digestive, and renal systems. Neurologic syndromes associated with EBV infection include meningitis, encephalitis, Guillain-Barré syndrome, Bell's palsy, myelitis, cranial nerve neuritis, and psychotic disorders. Bulbar involvement with ensuing respiratory paralysis can be fatal.³ EBV is also associated with Burkitt's lymphoma and nasopharyngeal carcinoma.¹

Many viral antigens are released in EBV infection. Viral capsid antigens (VCA) not only appear first but also develop in all cases of EBV infection. Since IgM to EBV capsid antigen is the acute phase antibody, it is the antibody of choice for detection of active disease. Serology is the only means to identify EBV infection in the absence of typical diagnostic symptoms of IM. Likewise, EBV serology is important in diagnosing infection in young children since they usually do not develop typical IM symptoms or heterophil antibodies.⁴

BIOLOGICAL PRINCIPLES

The indirect fluorescent antibody (IFA) method is used in the MERIFLUOR EBV VCA IgM IFA Test. Patient serum is reacted with the Burkitt's lymphocytic-cell substrate and if IgM antibodies to the viral capsid antigen of EBV are present they will bind to the antigen substrate and not rinse off. Subsequently, when fluoresceinated antihuman IgM is added to the reaction site, it will bind to the IgM antibodies, causing the EBV infected cells to fluoresce when viewed through the fluorescence microscope.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

- EBV Antigen Substrate Slides:** With HR1 Burkitt's lymphocytic cells fixed onto each well. Approximately 10 to 20% of the cells express the viral antigen to permit easy reading and optimal contrast. The cells are rendered noninfectious by the PsoralSafe™ process. The individually packaged slides are ready for use when removed from the foil packets. Slides stored at 2-8 C are stable until the date stated on the slide package. Do not use a slide if the desiccant indicator (line in center of desiccant) has changed from blue to pink.
- Phosphate Buffered Saline (PBS) Powder:** Rehydrate to 1 liter with distilled water. When rehydrated, the PBS is a 0.01 M phosphate buffer with pH of 7.5 and contains sodium azide. Rehydrated PBS stored at 2-8 C is stable indefinitely.
- Control Sera (See WARNINGS, step 4):** Lyophilized. Each control serum is reconstituted with 1 mL PBS to give a 1:10 working dilution. The reconstituted sera are stable for 6 weeks at 2-8 C or 8 months at ≤ -20 C. The Positive Control is human serum containing IgM antibodies to EBV. The titer is stated on the label.
- Conjugate:** Lyophilized. Fluorescein-labeled anti-human IgM, (goat), is reconstituted with 2 mL of PBS. The conjugate contains counterstain and is pretitered for use with each kit. The reconstituted conjugate can be stored up to 2 weeks at 2-8 C or aliquoted and stored up to 8 months at ≤ -20 C. Thawed aliquots should not be refrozen.
- Mounting Fluid:** The mounting fluid, which is buffered to pH 8.0, is a glycerol-water combination formulated to minimize elution of counterstain. It is ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- Serology laboratory supplies: 12 x 75 mm test tubes, test tube rack, serological and Pasteur pipettes
- Volumetric flask, 1 liter
- Wash bottle
- Distilled water
- Staining dish
- Cotton-tipped swabs, absorbent towels
- Moist chamber
- Cover slips, No. 1, 22 x 50 mm
- 37 C Incubator
- Fluorescence microscope. A FITC blue light excitation filter and a 515 nm barrier filter, or any comparable filter system, is suggested for transmitted-light microscopy using a darkfield condenser and for incident-light microscopy using a 500 nm dichroic mirror.

PRECAUTIONS

All reagents are for in vitro diagnostic use only.

WARNINGS

- The conjugate contains Evans Blue dye as a counterstain. This dye may be carcinogenic. Avoid contact with the skin.
- The PBS Powder contains sodium azide, which may react with lead or copper plumbing to form highly explosive metal azide compounds. When disposing of reagents through plumbing fixtures, flush with a large volume of water to prevent azide build-up in drains.
- Sodium azide is very toxic in case of ingestion. In contact with acid it emits a very toxic gas. In case of contact with the skin wash immediately and abundantly with water.
- Handle specimens, Controls and the materials that contact them as potential biohazards. The source materials used to prepare the control reagents were tested at the donor level and found to be nonreactive for the antigens or nuclear materials associated with, or antibodies to, the following viruses, as defined by required FDA-licensed tests: Human immunodeficiency viruses Types 1 and 2, Hepatitis B virus and Hepatitis C virus. However, since no test method can offer complete assurance that infectious agents are absent, these materials should be handled at the Biosafety Level 2 as recommended in the CDC/NIH "Biosafety and Microbiological and Biomedical Laboratories".

HAZARD and PRECAUTIONARY STATEMENTS

	Signal Word
  PBS Powder	Danger Hazard Statements H301 – Toxic if swallowed H310 – Fatal in contact with skin H411 – Toxic to aquatic life with long lasting effects Contains Sodium azide Precautionary Statements - EU (§28, 1272/2008) P301 + P310 – IF SWALLOWED: Immediately call a POISON CENTER or doctor/ physician P321 – Specific treatment (see .? on this label) P280 – Wear eye protection/ face protection P361 – Remove/Take off immediately all contaminated clothing P322 – Specific measures (see .? on this label)

SHELF LIFE AND STORAGE

- The expiration date is indicated on the kit label. Store the kit at 2-8 C and return the kit promptly to the refrigerator after each use.
- The kit is ready for use after reconstitution of reagents.
- If reconstituted conjugate and sera are to be stored at ≤ -20 C, self-defrosting freezers are not recommended.
- Individual slide packets should remain sealed until just before use.
- Reagents and antigen slides should not be used beyond stated expiration dates.
- Incubation times and temperatures other than those specified, may give erroneous results.
- The components of this kit have been tested and standardized as a unit. Use of components from other lots or other manufacturers may yield unsatisfactory results.

SPECIMEN COLLECTION AND PREPARATION

Blood obtained by venipuncture should be allowed to clot at room temperature and then be centrifuged. The serum should be separated as soon as possible and refrigerated (2-8 C), or stored frozen (≤ -20 C) if not tested within one week. Self-defrosting freezers are not recommended as storage units. The use of sera exhibiting hemolysis, lipemia or microbial growth is not recommended.

Screening dilutions of 1:10 and 1:40 are recommended for testing unless a serum pretreatment procedure is employed to remove interfering IgG (e.g., GullSORB™ IgG Inactivation Reagent, Product No. XX715). If the serum is pretreated with GullSORB, only the resulting 1:10 dilution need be tested. Other methods of serum pretreatment may require additional screening dilutions.

TEST PROCEDURE

1. Prepare screening dilution(s) of the patient serum with PBS or use pretreated serum.
2. Remove slide from foil packet.
3. Using a Pasteur pipette, add just enough appropriately diluted serum (approximately 15 µL) to cover each reaction site.
4. Place the slide in a moist chamber and incubate at 37 C for 90 minutes.
5. Rinse slide briefly with a gentle stream of PBS (avoid directing PBS at wells) and then immerse in PBS for 5 minutes.
6. Air dry the slide by standing on end on absorbent toweling.
7. Wipe between wells with a cotton-tipped swab moistened with distilled water.
8. Add just enough conjugate (approximately 15 µL) to cover each reaction site.
9. Incubate the slide at 37 C in a moist chamber for 30 minutes.
10. Rinse and dry the slide as described in Steps 5 and 6.
11. Just before reading, add a sufficient amount of mounting fluid to each well to entirely cover the substrate when the cover slip is applied. Cover with a 22 x 50 mm, no. 1 cover slip.
NOTE: Addition of too much mounting fluid may cause excessive movement of the cover slip. Addition of insufficient mounting fluid to a well may make part of the substrate unreadable.
12. Read the slide as soon as possible at 150X to 200X magnification.

INTERPRETATION OF RESULTS

The reaction is POSITIVE when 10 to 20% of the cells in each field exhibit a greenish-yellow fluorescence. The remaining cells provide a contrasting red background. A positive reaction at a 1:10 or greater dilution is indicative of a current EBV infection. It is not necessary to titrate positive sera in the MERIFLUOR EBV VCA IgM IFA Test. (See LIMITATIONS OF THE PROCEDURE to rule out false positives due to rheumatoid factor).

The reaction is NEGATIVE when cells do not fluoresce greenish-yellow but appear red due to the counterstain or exhibit less than a 1+ reaction at the screening dilution. A negative reaction at the screening dilution(s) indicates no detectable IgM antibody to EBV.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

Each kit contains positive and negative control sera which should be incorporated into each test run.

The negative control must exhibit less than a 1+ reaction at the 1:10 working dilution.

The positive control serum should be titered to provide a standard for checking test sensitivity. It should exhibit a 1+ reaction at its stated titer and increase to a strong positive reaction at the 1:10 dilution. A day-to-day variance of one two-fold dilution on either side of the stated titer is considered acceptable performance. A reaction of less than 2+ at the 1:10 dilution may indicate the test is not functioning within specifications.

If the test run is invalid because the Positive or Negative Control values do not fall within the limits specified, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated, please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor. Test results should not be reported until the condition(s) is corrected.

Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organizations. Refer to Clinical Laboratory Medicine – In Vitro Diagnostic Medical Devices – Validation of User Quality Control Procedures by the Manufacturer (ISO/FDIS 15198) for further guidance regarding Quality Control Practices.

EXPECTED VALUES

A positive result at a screening dilution of 1:10 (after RF interference has been ruled out) indicates an active (primary) EBV infection. IgM to EBV capsid is produced in all primary EBV infections, whether symptomatic or asymptomatic. Peak IgM titers in patients with diagnosed EBV infection range from a low of 40 to as high as 1280 with the MERIFLUOR EBV VCA IgM IFA Test. IgM antibodies decline rapidly to undetectable levels in 8-10 weeks.⁴

LIMITATIONS OF THE PROCEDURE

1. Immune IgG in the test sample can cause false negatives by competing with EBV-specific IgM for substrate binding sites. Immune IgG can also cause false positives by forming immune complexes with the antigen substrate which can then bind IgM-class rheumatoid factor (IgM-RF). Routine RF tests may not detect small amounts of IgM-RF which are recognized by the more sensitive IA technique. Therefore, serum samples intended for EBV IgM IFA Testing should be pretreated (e.g., GullSORB™ reagent, Product No. XX715) to eliminate immune IgG and RF interference.
2. IgM anticell antibodies, if present in the serum, may interfere with the MERIFLUOR EBV VCA IgM IFA Test. However, this interference is easily discerned since the staining produced by EBV-specific antibodies involves only the EBV infected cells (approximately 15%).
3. Endpoint titers can vary depending on the type of microscope, light source, age of bulb, and filter system used to read the reactions.
4. Due to the innate ability of HR1 lymphocytes to synthesize IgM in vitro, some lots of slides will exhibit a faint background fluorescence. This staining consists of unevenly distributed fluorescing dots seen within the majority of cells. This stippled staining, when it appears, is easily distinguished from the relatively uniform, whole-cell fluorescence produced by immune IgM to EBV, because the number of fluorescing cells and the fluorescent intensity, do not look like the positive control at the stated titer.
5. Serum samples with high IgM antibody titers to Cytomegalovirus and Herpes Simplex Virus may demonstrate a non-specific reaction at the 1:10 dilution.

SPECIFIC PERFORMANCE CHARACTERISTICS

When performed according to instructions, the MERIFLUOR EBV VCA IgM IFA Test is a sensitive and rapid method for diagnosing EBV infection. The sensitivity of each kit lot is controlled by testing with reference sera to ensure reproducible results. Clinical studies of the MERIFLUOR EBV VCA IgM IFA Test showed a 96% sensitivity and 100% specificity for IgM antibody to EBV capsid antigen.⁵

Tests for heterophil antibody and IgG antibody to EBV have inherent limitations in the rapid diagnosis of active disease, since they measure antibodies which develop more slowly and persist for longer periods of time. This antibody persistence can create confusion in patients later infected with other agents, especially when the original EBV infection was not diagnosed.⁶ IgM testing resolves the heterophil false positive responses and the false negative problem in children.⁴ Because IgG VCA titers are frequently at their peak at the time the first serum specimen is taken, diagnosis of IM by seroconversion or by demonstrating a rise in titer is unlikely. Likewise, diagnosing active disease by the IgG titer of a single serum specimen is difficult since the upper range of persisting convalescent titers overlaps with the lower range of maximal titers obtained in the acute phase of IM.⁷ Since IgM is present only in current disease, the MERIFLUOR EBV VCA IgM IFA Test requires only a single serum specimen, positive sera do not need to be titered, and results can be simply reported in terms of presence or absence of immune IgM.

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REFERENCES

1. Henle, W. and Henle, G. Epstein-Barr Virus and Infectious Mononucleosis, *Human Herpesvirus Infections*. Glaser and Gotlieb-Stematsky, ed. Marcel Dekker, Inc. New York. Basel. 1982.
2. Horwitz, C.A., et al. Heterophil-Negative Infectious Mononucleosis and Mononucleosis-Like Illnesses. Am. J. Med. 1977;63: 947-957.
3. Gotlieb-Stematsky, T. and Glaser, R. Association of Epstein-Barr Virus with Neurologic Diseases. *Human Herpesvirus Infections*. Glaser and Gotlieb-Stematsky, ed. Marcel Dekker, Inc. New York. Basel. 1982.
4. Nikoskelainen, J. and Hanninen, P. Antibody Response to Epstein-Barr Virus in Infectious Mononucleosis. Infect. Immun. 1975;11: 42-51.
5. Swenson, P.D., Clogston, A.G. and Fung, J.C. Detection of Immunoglobulin M Antibodies to Epstein-Barr Virus Capsid Antigen. Abstract of the Annual Meeting of the American Society for Microbiology. (New Orleans, Louisiana, 1983), p. 330.
6. Nikoskelainen, J., Klemola, E., and Leikola, J. Epstein-Barr Virus-IgM Antibody Test in Infectious Mononucleosis. J. Infect. Dis. 1976;134: 313-314.
7. Henle, W., Henle, G. and Horwitz, C.A. Epstein-Barr Virus Specific Diagnostic Tests in Infectious Mononucleosis. Human Path. 1974;5: 551-565..

SN11742

REV. 09/20

Manufactured By

Meridian Bioscience, Inc.
Corporate Office
3471 River Hills Drive
Cincinnati, Ohio 45244 USA
Telephone: 513.271.3700
Orders/Customer Service:
800.543.1980
Technical Support Center:
800.343.3858
Information Fax: 513.272.5432
Ordering Fax: 513.271.0124

EC REP

Authorized Representative

Meridian Bioscience Europe S. r. L
Via dell' Industria, 7
20035 Villa Cortese, Milano
ITALY
Tel: +39 0331 43 36 36
Fax: +39 0331 43 36 16
Email: info@meridianbioscience.eu
WEB: www.meridianbioscience.com/eu

Meridian Bioscience Europe s.a./n.v.
2 Avenue du Japon - 1420 Braine l'Alleud
BELGIUM
Tel: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
Email: info.bn1@meridianbioscience.eu

Meridian Bioscience Europe France
34 rue de Ponthieu - 75008 Paris
FRANCE
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 9 70 06 62 10
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
Postbus 301 - 5460 AH Veghel
NETHERLANDS
Tel: +31 (0) 411 62 11 66
Fax: +31 (0) 411 62 48 41
Email: info.bn1@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simblos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositivo medical de diagnóstico in vitro / Dispositivo medico para diagnostico in vitro / In-Vitro-Diagnosiskum	EC REP	Authorized representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandatario dans la Communauté Européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79/CE sur les dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon / Borsa, diluente inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenverarbeitung, in dem sich Probenverdünnschlüsse befindet
	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht eingefrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabbricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contiene suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	SOLN STOP	Stopping Solution / Soluzione di stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limite di temperatura / Limites de température / Límite de temperatura / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Conjugato enzimatico / Conjugé enzymatique / Conjugado enzimático / enzymkonzjugat
	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / testgarät	REAG	Reagent / Reagente / Réactifs / Reactivos / Reagenzen
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
	Buffer / Soluzione tamponi / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
	Conjugate / Conjugato / Conjugué / Conjugado / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentrate / Soluzione diluente di lavaggio / Solution de lavage concentrée 20X / Solución de tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
	Prescription Use Only / Per l'uso su prescrizione medica / Uniquement sur prescription / Solo Para Uso Por Receta / verschreibungspflichtig	DET REAG	Detection Reagent / Reagente Directo / Réactif de Détection / Reactivo de Deteción / Nachweis Reagens
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / ne pas utiliser si le paquet est endommagé / No use si el paquete está dañado / Nicht verwenden, wenn die Verpackung beschädigt ist	TUBE	Empty Tube / Provetta vuota / Tube vide / Tubo vacío / Leeres Gefäß

For technical assistance, call Technical Support Services at (800) 343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at (800) 543-1980.