MERIFLUOR[®]

EBV-NA

Immunofluorescence assay for antibodies to EBV-NA

REF EN100

IVD

In vitro diagnostic medical device

R. Only

INTENDED USE

The MERIFLUOR EBV-NA Test is a rapid, sensitive anti-complement immunofluorescence (ACIF) test for the qualitative or quantitative detection of antibodies to the nuclear antigen of Epstein-Barr virus in human serum. The test can provide valuable diagnostic information to facilitate patient management when used in conjunction with other EBV antibody tests.

SUMMARY AND EXPLANATION OF THE TEST

Epstein-Barr virus (EBV) is a widely-disseminated human pathogen which commonly causes infection in early childhood, resulting in asymptomatic or mild disease. However, individuals who contract infection with EBV in early adult life frequently develop classic infectious mononucleosis with symptoms of fever, pharyngitis, and cervical lymphadenopathy. Splenomegaly occurs in 50% of cases.¹ A variety of factors, in particular the immune competence of the individual, determine the clinical and pathologic expression of the disease. EBV infections can result in complications involving the neurologic, cardiac, ocular, respiratory, hematologic, digestive and renal systems.² Reactivation of latent infection has also been implicated in a persistent illness referred to in recent literature as the EBV-Associated Fatigue Syndrome.³

Infection with EBV results in the expression of multiple antigens with corresponding antibody responses. Of these, viral capsid antigen (VCA), early antigen (EA), and nuclear antigen (NA) are most useful diagnostically. At the onset of primary infection with EBV, VCA, IgG, VCA IgM, and EA antibody levels are elevated. Anti-NA levels begin to rise several weeks or months thereafter and remain at measurable titers, normally 1:10-1:160, indefinitely due to the persistent viral carrier state established following primary EBV infection. The lack of anti-NA (or very low levels of anti-NA) in the presence of IgG, IgM, and anti-EA antibodies is useful in establishing the presence of an acute primary infection.^{4,5}

Patients suffering from Burkitt's lymphoma and nasopharyngeal carcinoma have been shown to exhibit high levels of anti-NA antibodies.^{6,7} Conversely, anti-NA antibodies may be absent in immunosuppressed patients or those with severe immunologic defects.⁴

Indirect immunofluorescence (IFA) is the method most commonly used to test for EBV-specific antibodies. Testing for antibodies to EBV-NA requires greater sensitivity, and the ACIF technique is an easily-visualized means of determining anti-NA titers in human serum.⁸

BIOLOGICAL PRINCIPLES

The anti-complement immunofluorescence (ACIF) procedure is used in the MERIFLUOR EBV-NA Test. Patient serum is reacted with EBV-NA substrate and if antibodies to NA of EBV are present, they will bind to the substrate and not rinse off. Guinea pig complement (C') is then added and is bound by the antigen-antibody complex. Subsequently, when fluorescein-conjugated antibody to guinea pig C3 is added to the reaction site, it will bind to the C3 component of complement, causing the EBV-infected cells to have a granular pattern of nuclear fluorescence when viewed through a fluorescence microscope.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

- EBV-NA Substrate Slides:** Each with EBV-NA positive and negative cells fixed onto each well. Approximately 20-40% of the cells in each well are a Burkitt's lymphoma (Raji) cell line which is EBV-NA positive but cannot produce infectious EBV virions. The remaining cells do not express NA and are present to permit easy reading and optimal contrast with the positive cells. The individually packaged slides are ready for use when removed from the foil packets. Slides stored at 2-8 C are stable until the date stated on the slide package. Do not use a slide if the desiccant indicator (line in center of desiccant) has changed from blue to pink.
- Guinea Pig Complement (C'):** Lyophilized. Each vial sufficient for approximately 100 tests. Upon initial use, rehydrate one vial of C' with 1 mL of cold (2-8 C) sterile distilled water and divide into 100 µL aliquots. The aliquots should be tightly capped and immediately frozen at ≤ -20 C. Reconstituted C' is stable for 6 weeks at ≤ -20 C. Thawed aliquots should not be refrozen and reused.
- Complement Diluent:** Contains PBS with at least 0.01% each Ca++ and Mg++. Store at 2-8 C.
- Anti-Guinea Pig C3 (fluorescein-labeled):** Lyophilized, (goat). Reconstitute with 3 mL of PBS. This conjugate contains less than 0.1% Evans blue counterstain and is pretiered for use with each kit. The reconstituted conjugate can be stored for up to 2 weeks at 2-8 C or can be divided into aliquots and stored for up to 8 months at ≤ -20 C. Thawed aliquots should not be refrozen.
- Phosphate Buffered Saline (PBS) Powder:** Rehydrate to 1 liter with distilled water. The rehydrated solution is a 0.01 M phosphate buffer at pH 7.5 and contains sodium azide. Rehydrated PBS stored at 2-8 C is stable indefinitely.
- Positive and Negative Control Sera (See Warnings Step 4):** Heat inactivated and lyophilized. Each human control serum is reconstituted with 1 mL PBS to give a 1:10 working dilution. The titer is stated on the label of the Positive Control. The reconstituted sera are stable for 6 weeks at 2-8 C or for 8 months at ≤ -20 C. Thawed aliquots should not be refrozen. The control sera have been heat inactivated at 56 C for 30 minutes and should not be reheated before use.
- Mounting Fluid:** The mounting fluid, which is buffered to pH 8.0, is a glycerol-water combination formulated to minimize elution of counterstain. It is ready for use.

MATERIALS REQUIRED BUT NOT PROVIDED

- Serology laboratory supplies: 12 x 75 mm test tubes, test tube rack, serological and Pasteur pipettes
- Volumetric flask, 1 liter
- Ice bath
- Distilled water
- Sterile distilled water
- Staining dish
- Disposable tubes with tight-fitting caps
- Moist chamber
- Cover slips, No. 1, 22 x 50 mm
- 37 C Incubator
- 56 C heat bath
- Fluorescence microscope. A FITC blue light excitation filter and a 515 nm barrier filter, or any comparable filter system is suggested for transmitted-light microscopy using a darkfield condenser and for incident-light microscopy using a 500 nm dichroic mirror.

PRECAUTIONS

All reagents are for in vitro diagnostic use only.

WARNINGS

- The conjugate contains Evans Blue dye as a counterstain. This dye may be carcinogenic. Avoid contact with the skin.
- The PBS Powder contains sodium azide, which may react with lead or copper plumbing to form highly explosive metal azide compounds. When disposing of reagents through plumbing fixtures, flush with a large volume of water to prevent azide build-up in drains.
- Sodium azide is very toxic in case of ingestion. In contact with acid it emits a very toxic gas. In case of contact with the skin wash immediately and abundantly with water.

- Handle specimens, Controls and the materials that contact them as potential biohazards. The source materials used to prepare the control reagents were tested at the donor level and found to be nonreactive for the antigens or nuclear materials associated with, or antibodies to, the following viruses, as defined by required FDA-licensed tests: Human immunodeficiency viruses Types 1 and 2, Hepatitis B virus and Hepatitis C virus. However, since no test method can offer complete assurance that infectious agents are absent, these materials should be handled at the Biosafety Level 2 as recommended in the CDC/NIH "Biosafety and Microbiological and Biomedical Laboratories".

HAZARD and PRECAUTIONARY STATEMENTS

Refer to the SDS, available at www.meridianbioscience.com for Hazards and Precautionary Statements.

SHELF LIFE AND STORAGE

- Upon receipt of the kit, the two vials of complement should be detached and stored at ≤ -20 C. The remaining kit components should be stored in the refrigerator (2-8 C).
- The kit is ready for use after reconstitution of reagents.
- Individual slide packets should remain sealed until just before use.
- Reagents and antigen slides should not be used beyond stated expiration dates.
- If lyophilized reagents show evidence of rehydration, they should not be used.
- Self-defrosting freezers are not recommended for storage of complement, conjugate or control sera.
- The components of this kit have been tested and standardized as a unit. Use of components from other lots or other manufacturers may yield unsatisfactory results.

SPECIMEN COLLECTION AND PREPARATION

Blood obtained by venipuncture should be allowed to clot at room temperature and then be centrifuged. The serum should be separated as soon as possible and refrigerated (2-8 C), or be stored frozen (≤ -20 C) if not tested within one week. Self-defrosting freezers are not recommended as storage units. The use of sera exhibiting hemolysis, lipemia or microbial growth is not recommended.

Before testing, an aliquot of the serum sample must be heated for 30 minutes at 56 C to inactivate endogenous complement.

For qualitative results, sera should be screened at a 1:10 dilution. For quantitative results, sera can be titrated using serial two-fold dilutions. All serum dilutions should be prepared with PBS.

TEST PROCEDURE

- Prepare dilution(s) of heat-inactivated patient serum with PBS.
- Remove slide from foil packet.
- Using a Pasteur pipette cover each reaction site with approximately 15 µL of diluted serum.
- Place the slide in a moist chamber and incubate at 37 C for 30 minutes. **Important: From this point through Step 16, slides must be kept moist.**
- Five to ten minutes before the end of the incubation period, thaw an aliquot of guinea pig complement and immediately add the appropriate volume of cold Complement Diluent to obtain a working complement solution (See chart below). **Keep this reagent on ice during test procedure.**

# of slides in run	C' (µL) (freshly thawed)	C' diluent (mL)
1	10	0.4
2	20	0.8
3	30	1.2
4	40	1.6
5	50	2.0
6	60	2.4
7	70	2.8
8	80	3.2
9	90	3.6
10	100	4.0

- At the end of the incubation period, rinse slide briefly 3 times with a gentle stream of PBS, gently flicking slide between rinses (avoid directing stream at wells). Rinse again, and allow PBS to bead up on wells. Let slide sit in a moist chamber for 5-10 minutes.
- Holding slide upright, tap gently on a folded piece of paper towel to remove the PBS. **Do not allow the wells to dry.**
- Quickly add just enough cold working complement to cover each reaction site (approximately 15 µL).
- Repeat steps 7 and 8 until all slides are processed.
- Incubate the slides at 37 C in a moist chamber for 30 minutes.
- Wash as in step 6.
- Working with one slide at a time, hold slide and remove excess PBS with a gentle flick of the wrist. **Do not allow the wells to dry.**
- Add approximately 15 µL of Anti Guinea-Pig C3 to each reaction site.
- Repeat steps 12 and 13 until all slides are processed.
- Incubate the slide at 37 C in a moist chamber for 30 minutes.
- Rinse slide briefly 3 times with a gentle stream of PBS (avoid directing stream at wells).
- Air dry the slide by standing on end on absorbent toweling.
- Just before reading, add a sufficient amount of mounting fluid to each well to entirely cover the substrate when the cover slip is applied. Cover with a 22 x 50 mm, no. 1 cover slip. **NOTE:** Addition of too much mounting fluid may cause excessive movement of the cover slip. Addition of insufficient mounting fluid to a well may make part of the substrate unreadable.
- Read the slide at 200X magnification. If the slide cannot be read immediately, it may be stored in the dark at room temperature for up to 24 hours.

INTERPRETATION OF RESULTS

The reaction is POSITIVE when 20-40% of the cells in each field exhibit a greenish-yellow granular pattern of nuclear fluorescence against a contrasting red background. The EBV-NA antibody titer is the highest dilution of serum which produces a 1+ fluorescence. A positive reaction at a dilution of 1:10 or greater indicates the presence of antibody to EBV-NA. If desired, positive sera may be titered using two-fold serial dilutions.

The reaction is NEGATIVE when cells do not fluoresce greenish-yellow but appear red due to the counterstain or exhibit less than a 1+ reaction at the screening dilution. A negative reaction at the 1:10 dilution indicates no detectable antibody to EBV-NA.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies. Each kit contains positive and negative control sera which should be incorporated into each test run.

The Negative Control must exhibit less than a 1+ reaction at the 1:10 working dilution.

The Positive Control provides a standard for checking test sensitivity. When titrated, it should exhibit a strong positive reaction at the 1:10 working dilution and diminish to a 1+ reaction at its stated titer. A day-to-day variance of one two-fold dilution on either side of the stated titer is considered acceptable performance.

If the test run is invalid because the Positive or Negative Control values do not fall within the limits specified, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800 343-3858 (US) or your local distributor. Test results should not be reported until the condition(s) is corrected. Repeat the control tests as the first step in determining the root cause of the failure.

Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organizations. Refer to Clinical Laboratory Medicine – In Vitro Diagnostic Medical Devices – Validation of User Quality Control Procedures by the Manufacturer. (ISO/FDIS 15198) for further guidance regarding Quality Control Practices.

EXPECTED VALUES

EBV-NA antibodies are almost always present in sera containing IgG antibodies to VCA of EBV. However, NA is generally absent during the acute phase of IM. NA may be absent for several months post-acute infection. Patients with severe immunologic defects or immunosuppressive diseases may not have EBV-NA antibodies, even if antibodies to VCA are present.⁴

During convalescence, antibody titers to NA gradually increase and persist at moderate levels for life. Most healthy individuals previously exposed to EBV have antibody titers to EBV-NA that range from 1:10 to 1:160.⁴

In EBV-associated malignancies, levels of antibodies to NA of EBV are usually high in patients with nasopharyngeal carcinoma and can range from barely detectable to very high in patients with Burkitt's lymphoma.^{6,7}

MERIFLUOR EBV-NA Test results can be especially useful when correlated with results from other serological tests for EBV. For example, a patient with an IM-like illness caused by reactivation of a persistent EBV infection due to an immunosuppressive malignant or nonmalignant disease may demonstrate high titers of both IgM and IgG to EBV-VCA. However, if antibody to EBV-NA is also present, a diagnosis of primary EBV infection can be excluded.⁹

LIMITATIONS OF THE PROCEDURE

1. **Allowing wells to dry during the test procedure will invalidate the test results, usually giving a false negative reaction, and require the test to be repeated.**
2. Storage of complement under conditions other than those specified may result in decreased titers or negative reactions.
3. Anti-cell antibodies, if present in the serum, may interfere with the MERIFLUOR EBV-NA Test. However, this interference is easily discerned since the staining produced by these antibodies involves all of the cells, whereas antibody staining specific for EBV-NA involves only the 20-40% of the cells expressing EBV-NA.
4. Test results for antibodies to EBV-NA should be evaluated by clinicians in relation to patient symptoms, clinical history, and antibody response patterns to EBV-VCA and EA to establish a diagnosis.
5. Endpoint titers can vary depending on the type of microscope, light source, age of bulb, and filter system used.

SPECIFIC PERFORMANCE CHARACTERISTICS

The MERIFLUOR EBV-NA Test is a sensitive and rapid method for the qualitative and quantitative determination of antibodies to EBV-NA in human serum. It is of value in providing supportive information for the diagnosis of EBV-associated disease.

The sensitivity of each kit lot is controlled by testing with reference sera to ensure reproducible results from lot to lot. The reproducibility of the MERIFLUOR EBV-NA Test was evaluated and found to be within an acceptable range of one two-fold dilution. Test results from clinical studies and the testing of reference sera indicate the MERIFLUOR EBV-NA Test is specific for antibodies to NA of EBV and does not cross react with antibodies to the other herpesviruses. Clinical specimens from 106 patients were tested using the MERIFLUOR EBV-NA Test and another ACIF method. The MERIFLUOR EBV-NA Test demonstrated a relative sensitivity (70/71) of 98.6% and a specificity (35/35) of 100% compared to the reference method. The positive predictive value of the MERIFLUOR EBV-NA Test was 100% and the negative predictive value was 97.2%, for an overall predictive value of 99.1%.

Other ACIF	MERIFLUOR EBV-NA TEST	
	Neg	Pos
	Neg	70
Pos	1	35

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REFERENCES

- Evans AS and Niederman JC. Epstein-Barr Virus. In *Viral Infections of Humans: Epidemiology and Control*. 2nd ed. Evans, ed. Plenum Medical. New York, 1982.
- Gotlieb-Stematsky T. and Glaser R. Association of Epstein-Barr Virus with Neurologic Diseases. In *Human Herpesvirus Infections: Clinical Aspects*. Glaser and Gotlieb-Stematsky, eds. Marcel Dekker. New York, 1982.
- Sumaya CV. Infectious Mononucleosis and Other EBV Infections: Diagnostic Factors. *Laboratory Management* 1986;24:37-46.
- Henle W, Henle G, and Horwitz CA. Infectious Mononucleosis and Epstein-Barr Virus Associated Malignancies. In *Diagnostic Procedures for Viral Rickettsial and Chlamydial Infection*. 5th ed. Lenette and Schmidt, eds. American Public Health Association, Inc. Washington, D.C., 1979.
- Henle G, et al. Rheumatoid Factor as a Cause of Positive Reactions in Tests for Epstein-Barr Virus-Specific IgM Antibodies. *Clin. Exp. Immunol.* 1979;36:415-22.
- Lindahl T, et al. Relationship Between Epstein-Barr Virus (EBV) DNA and the EBV-Determined Nuclear Antigen (EBNA) in Burkitt Lymphoma Biopsies and Other Lymphoproliferative Malignancies. *Int. J. Cancer* 1974;13:764-72.
- Huang DP, et al. Demonstration of Epstein-Barr Virus-Associated Nuclear Antigen in Nasopharyngeal Carcinoma Cells from Fresh Biopsies. *Int. J. Cancer* 1974;14:580-88.
- Ernberg I and Klein G. EB Virus-Induced Antigens. In *The Epstein-Barr Virus*. Epstein and Achong, eds. Springer-Verlag, Berlin. 1979.
- Henle W and Henle G. The Virus as the Etiologic Agent of Infectious Mononucleosis. In *The Epstein-Barr Virus*. Epstein and Achong, eds. Springer Verlag, Berlin. 1979.

SN11748

REV. 02/16

Manufactured By

Meridian Bioscience, Inc.
USA/Corporate Office
3471 River Hills Drive
Cincinnati, Ohio 45244
Telephone: 513.271.3700
Orders/Customer Service:
800.543.1980
Technical Support Center:
800.343.3858
Information Fax: 513.272.5432
Ordering Fax: 513.271.0124

Authorized Representative

Meridian Bioscience Europe S. r. l.
Via dell'Industria, 7
20020 Villa Cortese, Milano
ITALY
Tel: +39 0331 43 36 36
Fax: +39 0331 43 36 16
Email: info@meridianbioscience.eu
WEB: www.meridianbioscience.eu

Meridian Bioscience Europe s.a./n.v.
2 Avenue du Japon - 1420 Braine l'Alleud
BELGIUM
Tel: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
Email: info.bn@meridianbioscience.eu


Meridian Bioscience Europe France
34 rue de Ponthieu - 75008 Paris
FRANCE
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 67 89 59 58
Email: info.fr@meridianbioscience.eu

Meridian Bioscience Europe b.v.
Postbus 301 - 5460 AH Veghel
NETHERLANDS
Tel: +31 (0) 411 62 11 66
Fax: +31 (0) 411 62 48 41
Email: info.bn@meridianbioscience.eu

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symboles, Guía de símbolos, Zeichenerklärung)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostika 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller		For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur IVD Leistungsbewertung
	Contains sufficient for <n> tests / Contenuto sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	SOLN STOP	Stopping Solution / Soluzione di Stop / Solution d'arrêt / Solución de parada / Stopplösung
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung	CONJ ENZ	Enzyme Conjugate / Coniugato enzimatico / Conjugué enzymatique / Conjugado enzimático / Enzymkonjugat
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	REAG	Reagent / Reagente / Réactifs / Reactivos / Reagenzien
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum	BUF WASH	Wash Buffer / Soluzione di lavaggio / Solution de lavage / Tampón de lavado / Waschpuffer
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
CONJ	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat	DIL SPE	Specimen Diluent (or Sample Diluent) / Diluente del Campione / Diluant échantillons / Diluyente de muestra / Probenverdünnungspuffer
SUBS	Substrate / Substrato / Substrat / Substrato / Substrat	BUF WASH 20X	Wash Buffer Concentration: 20X / Soluzione di lavaggio 20X / Solution de lavage concentrée 20X / Solución tampón de lavado 20X / 20fach konzentriertes Waschkonzentrat
Rx Only	Caution: US Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner. / La legge Federale negli Stati Uniti prevede che questo test sia venduto solo a un operatore abilitato ad operare nella salute. / Attention: La vente ou la prescription de ce dispositif par un médecin agréé est soumise à des restrictions par la loi fédérale des US. / Cuidado: leyes federales de los EU limitan la venta de dispositivos por o por la orden de licenciados en la práctica del cuidado de la salud. / Vorsicht: Der Verkauf oder die Verschreibung dieses Geräts durch einen zugelassenen Arzt unterliegt den Beschränkungen des US-Bundesgesetzes.	DET REAG	Detection Reagent / Reagente Diretto / Réactif de Detection / Reactivo de Detección / Nachweise Reagenz

For technical assistance, call Technical Support Services at (800) 343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at (800) 543-1980.